

# Research Collection

## Steel Laboratory Furniture



KEWAUNEE®

encouraging new discovery...*Worldwide*

# RESEARCH COLLE

## *Steel Furniture*

- Fully-welded, heavy-gauge steel construction
- Constructed to meet SEFA standards
- Powder coat, lab grade finish
- Stainless steel construction option
- Customize with multiple styles, pulls, and hardware



# CTION



## Contents

Steel Furniture Styles & Construction .....	2
Catalog Number Explanation .....	4
Steel Cabinet Options .....	6
Sitting Height Cabinets.....	11
ADA Height Cabinets.....	14
Standing Height Cabinets.....	22
Mobile Cabinets.....	34
Apron Assemblies .....	38
Wall Mounted Storage Cabinets.....	46
Full Height Storage Cabinets.....	53
Fillers, Scribes & Panels .....	58
Sloping Tops & Enclosures.....	66
Fume Hoods & Ventilated Enclosures ...	68
Work Tops .....	70
Center Bench Shelving.....	72
Wall Shelving.....	74
Reagent Racks & Shelving.....	76
Pegboards .....	78
Miscellaneous Accessories.....	79
General Plumbing Information.....	82
Service Fittings .....	83
Drain Accessories .....	91
Sinks & Troughs .....	92
Electrical Fittings .....	96
Seating .....	98
Technical Information.....	100
Alphabetical Index.....	106
Numerical Index .....	108
Typical Installations .....	112

Beautifully finished, smooth, clean, easy to maintain surfaces. Rugged, strong, heavy gauge steel, full-frame, four sided, flush front construction. That's Kewaunee's Steel Laboratory Furniture.

# Kewaunee Steel Cabinet Construction

Kewaunee's fully-welded, heavy-gauge, steel cabinets are designed to SEFA 8 performance standards and built to rigorous benchmarks. The Kewaunee classic "unit" design allows doors

and drawers to be reconfigured in place. The configurable design provides a rich offering of hardware and aesthetic options.

## DRAWER OPTIONS



Standard Drawers feature one-piece drawer body construction with full-extension, side mounted, soft-close, slides. (model F-8400-00)




Optional Concealed Slide Drawers feature integral concealed, full-extension, soft-close slides



### Additional Cabinet Features

- Four sided, flush front, frame construction
- Solid finished back with optional removable back
- Full-depth adjustable shelf
- Leveling bolts accessible through cabinet top

## PULL STYLE OPTIONS




## DOOR CATCH



Positive Catch

Catch furnished with 5-knuckle hinges only

## DOOR HINGE OPTIONS



# Kewaunee Steel Cabinet Styles



*shown with five knuckle hinge and rounded aluminum pull style 1*

## Inset - Square Edge

Doors and drawer heads are flush with face of cabinet. Available with Positive Catch and a wide selection of pulls.



*shown with maple veneer door and drawer head, five knuckle hinge and square aluminum pull style 8*

## Overlay - Wood on Steel

Oak or maple veneer doors and drawer heads overlay face of cabinet. Available with 5-knuckle, semi-concealed or concealed hinges and a wide selection of pulls. *(Flush Pull 9 not available on wood)*



*shown with five knuckle hinge and chrome plated wire pull style 6*

## Overlay - Square Edge

Doors and drawer heads overlay face of cabinet. Available with 5-knuckle, semi-concealed or concealed hinges and a wide selection of pulls.



*shown with concealed hinge and beveled label holder pull style 7*

## Overlay - Bevel Edge

Doors and drawer heads, with beveled top edge, overlay face of cabinet. Available with 5-knuckle, semi-concealed or concealed hinges and a wide selection of pulls.

# Steel Cabinet Information

## Catalog Number Explanation

**E40C352236-1611-B11,C...**

### Cabinet Type

A = Apron Assemblies and Legs  
 B = Open Cupboard Base Cabinet  
 C = All Cupboard Base Cabinet  
 D = All Drawer Base Cabinet  
 E = Drawer & Cupboard Base Cabinet  
 G = Special Purpose Base Cabinet  
 S = Full Height Storage Cabinet  
 W = Wall Storage Cabinet

### Door & Drawer Configuration

### Cabinet Construction & Design

### Cabinet Height

#### Base Cabinets

29 = 29" High    32 = 32" High    35 = 35" High

#### Aprons

4 = 4" High

#### Full Height Cabinets

84 = 84" High

#### Wall Cabinets

24 = 24" High    30 = 30" High    36 = 36" High    48 = 48" High

### Cabinet Options

(See Cabinet Option list page 5)

### Style – Hardware – Material

(See Selection list below)

### Special Configurations

L = Reverse Construction/Left Hand Swing

### Cabinet Length

Standard Lengths: 12" – 15" – 18" – 21"

24" – 30" – 36" – 42"

48" – 52" – 58" – 60"

### Cabinet Depth

13 = 12<sup>3</sup>/<sub>4</sub>" Deep

20 = 20" Deep

16 = 16" Deep

22 = 22" Deep

28 = 28" Deep

Note: Full Overlay Styles add 7/8" to cabinet depth

## Door & Drawer Head Style – Hardware – Material Selections

**E40C352236-1611-B11,C...**

### Cabinet Style - Hinge Type - Catch

0 = Inset            - 5-knuckle Hinge - Positive Catch  
 1 = Full Overlay   - 5-knuckle Hinge - Positive Catch  
 4 = Full Overlay   - Semi-concealed Hinge  
 6 = Full Overlay   - Concealed Hinge (170° Swing)

### Door & Drawer Head Style

1 = Square Edge Steel  
 6 = Wood – available on Full Overlay Cabinet Styles only  
 7 = Beveled Edge Steel – available on Full Overlay Cabinet Styles only

### Cabinet Material / D&D Material

0 = Cold Rolled Steel / Cold Rolled Steel    \*3 = Stainless Steel / Stainless Steel  
 1 = Cold Rolled Steel / Red Oak                5 = Cold Rolled Steel / FSC® Red Oak  
 2 = Cold Rolled Steel / White Maple           6 = Cold Rolled Steel / FSC® White Maple  
 7 = Cold Rolled Steel / Steamed Beech        8 = Cold Rolled Steel / FSC® Steamed Beech


\* Stainless Steel Cabinets available in D&D Style 1 only

### Pull Type

0 = No Option  
 1 = Rounded Aluminum  
 2 = Black Wire  
 3 = White Wire  
 4 = Stainless Steel Wire  
 6 = Chrome Plated Wire  
 ‡7 = Beveled Label Holder  
 8 = Square Aluminum  
 †9 = Flush Aluminum  
 ‡ mounts horizontally  
 † not available on wood doors & drawer heads

# Steel Cabinet Information

## Drawer and Cupboard Size Coding



Drawer Code

Cupboard Code

<b>E40C352215L</b> _____	15"L
<b>E40C352218L</b> _____	18"L
<b>E40C352221L</b> _____	21"L
<b>E40C352224L</b> _____	24"L

Base cabinets are designed on a "unit" system to allow maximum flexibility in the selection of drawer and/or cupboard door sizes.

Standing height cabinets (35" high) have an opening 60 "units" high. Sitting height cabinets (29" high) have an opening 48 "units" high. ADA height cabinets (32" high) have an opening 54 "units" high. This means that any combination of drawer and cupboard door heights equaling 60, 54 or 48 units can be used in a base cabinet.

Mobile standing height cabinets (33" high) have an opening 56 "units" high. Mobile sitting height cabinets (26 1/4" high) have an opening 42 "units" high.

The adjacent tables indicate the "unit" height and inside clearance of the Inset style base cabinets shown in this catalog. Simply match the unit number shown on the cabinet drawers and cupboard to the appropriate number in the tables.

### Drawer Size Codes

"Units" High	Drawer Head Height	Inside Dimensions	
		Clear/Height	Length (Front-to-Back)
8*	3 3/4"	2 3/8"	18"
9*	4 1/4"	2 3/8"	18"
10*	4 3/4"	3 7/16"	18"
11	5 3/16"	3 15/16"	18"
12	5 11/16"	4 7/16"	18"
14	6 5/8"	5 3/8"	18"
15	7 1/8"	5 7/8"	18"
16	7 5/8"	6 3/8"	18"
20	9 9/16"	8 5/16"	18"
24	11 1/2"	10 1/4"	18"
3*	1 3/16"	† Pull Out Shelf	

\* Not available with Option Z Concealed Slide Drawer  
 † Pull Out Shelf mounted on full extension slides only

### Drawer Width I.D. (side-to-side)

Base Cabinet Length	Number of Drawers Across		
	1	2	3
12"	8 3/8"		
15"	11 3/8"		
18"	14 3/8"		
21"	17 3/8"		
24"	20 3/8"	8 3/8"	
30"	26 3/8"	11 3/8"	
36"	32 3/8"	14 3/8"	8 3/8"
42"	38 3/8"	17 3/8"	
48"	44 3/8"	20 3/8"	12 3/8"

### Cupboard Size Codes

"Units" High	Door Height	Inside Dimensions	
		Height (Top-to-Bottom)	Depth (Front-to-Back)
30	14 9/16"	13 9/16"	20 1/2"
36	17 7/16"	16 7/16"	20 1/2"
45	21 13/16"	20 13/16"	20 1/2"
48	23 3/16"	22 3/16"	20 1/2"
54	26 1/16"	25 1/16"	20 1/2"
60	29"	28"	20 1/2"

### Door Opening Clearance

#### Width I.D. (side-to-side)

Base Cabinet Length	Number of Drawers Across		
	1	2	3
12"	9 1/2"		
15"	12 1/2"		
18"	15 1/2"		
21"	18 1/2"		
24"	21 1/2"	9 1/2"	
30"	27 1/2"	12 1/2"	
36"	33 1/2"	15 1/2"	9 1/2"
42"	39 1/2"	18 1/2"	
48"	45 1/2"	21 1/2"	13 1/2"

Cupboard and Drawer sizes listed in these charts are for Inset style cabinets only. Sizes for Overlay styles are generally slightly larger. Contact your Kewaunee Dealer for more information.

## Cabinet Options – (Add to End of Catalog Number)

### B – Base Cabinet Caster Option

- B11** – 4" Swivel Caster w/ brakes & with counter balance
- B21** – 4" Swivel Caster w/o brakes & with counter balance
- B31** – 4" Swivel Caster w/ brakes at front & counter balance
- B41** – 3" Swivel Caster w/ brakes & with counter balance
- B51** – 3" Swivel Caster w/o brakes & with counter balance
- B61** – 3" Swivel Caster w/ brakes at front & counter balance
- B71** – Fixed rear – Swivel front – with counter balance
- B81** – Fixed rear – Glides front – with counter balance

**C** – Toe space depth increased by 7/8" (for B00 & G07)

**F** – Removable Back (not available on Mobile Cabinets)

**H** – Label Holders \*

\* **J** – Locks w/ Lock Panels    **J**=Disk Locks †    **J2**=Pin Locks †

\* **K** – Locks w/o Lock Panels    **K**=Disk Locks †    **K2**=Pin Locks †

**O** – Shelf Lips

### P – Pullout Trays & Pullout Pans

**P1** – One Pullout Tray - at bottom of cupboard

**P2** – Two Pullout Trays - center & bottom

**P3** – One Pullout Pan - at bottom of cupboard

**P4** – Two Pullout Pans - center & bottom

**I** – 1/4" Laminated Safety Glass

**R** – 1/8" Tempered Safety Glass

**T** – Heavy Duty Full Extension Slide

**U** – Toe Space on Base Cabinet Ends

**U1** – Toe Space on Left End

**U2** – Toe Space on Right End

**U3** – Toe Space on Both Ends

**Z** – Concealed Slide Drawer

**Z1** – Cold Rolled Steel Drawer Sides

**Z2** – Stainless Steel Drawer Sides

\* Furnished on all cupboards and drawers on cabinet    † Lock keying must be specified.

# Steel Cabinet Options

## Base Cabinet Caster Option

Base Cabinet Caster Option adds casters to the base cabinet. Eight caster styles are available and all are furnished with counter balance weight for cabinets that may be top heavy when loaded and/or prone to tipping when opened. (**Note:** Casters ship loose)

### Swivel Casters

adds four (4) swivel casters and counter balance weight for stability.

- B11** – \*4" Swivel Casters with brakes
- B21** – \*4" Swivel Casters without brakes
- B31** – \*4" Swivel Casters with brakes at front only
- B41** – 3" Swivel Casters with brakes
- B51** – 3" Swivel Casters without brakes
- B61** – 3" Swivel Casters with brakes at front only

\* Adds 1" to height of cabinet

### Fixed and Swivel Casters

adds two (2) fixed 4" casters at the rear, two locking swivel casters at the front, and counter balance weight for stability.

- B71** – Fixed rear casters, locking front swivel casters

### Fixed Casters and Glides

adds two (2) fixed 4" casters at the rear, two (2) adjustable glides at the front, and counter balance weight for stability

- B81** – Fixed rear casters, front adjustable glides.

To specify Casters, add **B**\_\_ to the end of the cabinet number.

Example: **E40C352236-0101-B71**  
(Available on base cabinets only.)

## Toe Space Depth Increase by 7/8"

Option **C** adds 7/8" to the depth of a base cabinet toe space. This option is only used on base cabinets without doors or drawers (B00 & G07) to

allow the face of the cabinet to be pulled forward to align with the face of the doors and drawers on adjacent cabinets.

To specify increased toe space depth add **C** to the end of the

cabinet number.  
Example: **B00C352236-0000-C**

## Removable Back

All cabinets are furnished with a solid back for improved cleanliness. A removable back is available when required for access to the pipe space behind the cabinet.

**Note:** Removable Backs are not available on Mobile Cabinets.

To specify Removable Back, add **F** to the end of the cabinet number.

Example: **D30C352236-0101-F**

## Label Holders

Option **H** adds the F-0260-00 aluminum Label Holder to each drawer and door on the cabinet. Overall size is 3" x 1<sup>5</sup>/<sub>8</sub>" and holds a 2<sup>1</sup>/<sub>2</sub>" x 1<sup>1</sup>/<sub>8</sub>" card. Double doors


receive only one Label Holder per pair.

**Shipped loose for field installation.**


To specify Label Holders, add **H** to the end of the cabinet number.

Example: **E40C352236-0101-H**  
(Added to all cupboards and drawers.)


## Option B



**B11 – B21 – B31 Caster Option**




**B71 Caster Option**




**B81 Caster Option**

## Option C

## Option F



## Option H



One-half Actual Size



## Steel Cabinet Options

### Shelf Lips

### Option O

Option **O** adds a 3/4" high steel lip to the front edge of all steel shelves in the cabinet.

To specify Shelf Lips, add **O** to the end of the cabinet number.

Example: **W00C361336-0101-O**  
(Added to all shelves in cabinet.)

### Pullout Trays & Pans

### Option P

**P1** – One Pullout Tray - at bottom of cupboard

**P2** – Two Pullout Trays - center & bottom

**P3** – One Pullout Pan - at bottom of cupboard

**P4** – Two Pullout Pans - center & bottom

Pullout Trays and Pans are fabricated from 18 gauge cold rolled steel and are mounted on full-extension slides. Both have 2<sup>3</sup>/<sub>8</sub>" high sides and back. Trays are open at the front while pans are enclosed on all four sides. Each is available in two configurations: one tray or pan mounted at bottom of cupboard; or

two trays or pans mounted at bottom and in center of cupboard. When pullout trays or pans are added, the adjustable shelf is omitted.


Pullout Trays and Pans are available for base cabinets only.

To specify a Pullout Tray, add **P1** or **P2** to the end of the cabinet number.


To specify a Pullout Pan, add **P3** or **P4** to the end of the cabinet number.

Example: **E40C352236-0101-P2**  
(Available for base cabinets only.)

Example: **E40C352236-0101-P3**  
(Available for base cabinets only.)



Pullout Tray



Pullout Pan

### Laminated Safety Glass

### Option I

### Tempered Safety Glass

### Option R

Option **I** replaces the standard float glass in plate glass or glazed door wall, counter mounted and full height cabinets with 1/4" Laminated Safety Glass.

Option **R** replaces the standard 1/8" float glass in glazed door wall, counter mounted and full height cabinets with 1/8" Tempered Safety Glass.

To specify Laminated Safety Glass, add **I** to the end of the cabinet number.  
Example: **W20C351336-0101-I**

To specify Tempered Glass, add **R** to the end of the cabinet number.  
Example: **W20C351336-0101-R**

### Full Extension Heavy Duty Slides


### Option T

Full extension ball bearing drawer slide with cushioned in and outstops and positive drawer hold-in detent. Loading capacity of 200 lbs. Replaces standard 100 lbs full-

extension slides on all drawers in cabinet when **T** option is specified.

To specify Heavy Duty Slides, add **T** to the end of the cabinet number.

Example: **E40C352236-1011-T**



### Toe Space on Cabinet Ends

### Option U

When cabinets are used at the end of center benches and in some other special applications, they may require a toe space on one or

both ends of the cabinets. These toe spaces match the toe space at the front of the cabinets and may be added to the left, right, or both ends.

- U1** – Toe Space on Left End
- U2** – Toe Space on Right End
- U3** – Toe Space on Both Ends

To specify End Toe Spaces, add **U\_** to the end of the cabinet number.

Example: **E40C352236-0101-U2**  
(Available for base cabinets only.)

# Steel Cabinet Options

## Locks with Lock Panels Locks without Lock Panels

**Option J**  
**Option K**


All cabinets may be ordered with locks by specifying the K option. In addition, base cabinets may be specified with the J option to add lock panels. Lock panels are

installed between drawers and/or drawers and cupboards to provide additional security when the locks on a cabinet are keyed different.

### Locks Types

Cabinet locks are available in two types. Disk Tumbler which offers medium security, and Pin Tumbler which offer a higher level of security. Pin Tumbler also allows more keying arrangements. Both are fabricated with die cast sleeves and cylinders

with a satin nickel face and zinc plated steel cams. Both locks use the same sleeve and cam assembly allowing the lock cylinder to be removed and switched after the lock has been installed. All locks are furnished with two keys unless otherwise specified.



**Lock Panel**

### Disk Tumbler Locks

**Security:** Medium

**Keying:** Different, Alike, Alike in Sets, Masterkeyed

**Key**

**Changes:** Up to 2000

**Master Key**

**Changes:** Up to 5

### Pin Tumbler Locks

**Security:** Medium-High

**Keying:** Different, Alike, Alike in Sets, Masterkeyed, Grand-masterkeyed, Sub-masterkeyed, Mason Keyed

**Key**


**Changes:** Up to 2000

**Master Key**


**Changes:** Up to 52

**J** – Disk Locks with Lock Panels  
**K** – Disk Locks w/o Lock Panels

**J2** – Pin Locks with Lock Panels  
**K2** – Pin Locks w/o Lock Panels



**Disk Tumbler Lock**  
(Medium Security)



**Pin Tumbler Lock**  
(Higher Security)

**Lock schedule with keying arrangement must be furnished when locks are specified.** See page 9 & 10 for keying arrangement diagrams, and page 104 & 105 for sample lock schedules.

To specify Locks without Lock Panels, add **K** or **K2** to the end of the cabinet number.

Example: **E40C352236-0101K2**  
(Added to all cupboards and drawers.)

To specify Locks with Lock Panels, add **J** or **J2** to the end of the cabinet number.

Example: **E40C352236-0101J**  
(Added to all cupboards and drawers.)

## Concealed Slide Drawers

**Option Z**

Overlay cabinets may be upgraded to a Concealed Slide Drawer with integral, concealed, full-extension, soft-close, slides. Concealed Slide Drawers are available with either powder-coated

steel drawer sides or with stainless steel drawer sides.

Concealed Slide Drawers are available in Overlay cabinets only and are not available for all drawer sizes.

To specify Concealed Slide Drawers with Powder-coated Steel sides, add **Z1** to the end of the cabinet number.

Example: **E40C352236-1101Z1**  
(Added to all drawers.)

To specify Concealed Slide Drawers with Stainless Steel Sides, add **Z2** to the end of the cabinet number.


Example: **E40C352236-1101Z2**  
(Added to all drawers.)

# Lock Keying Options

## Keyed Different

**Disk Tumbler Locks:** Available  
**Pin Tumbler Locks:** Available


All locks are keyed different from each other with individual keys for each lock. No key will open more than one lock.



## Keyed Alike

**Disk Tumbler Locks:** Available  
**Pin Tumbler Locks:** Available


All locks are keyed alike. One key will open all locks.



## Keyed Alike in Sets

**Disk Tumbler Locks:** Available  
**Pin Tumbler Locks:** Available


Locks are grouped into sets, with as many locks as needed in each set. All locks in a set can be opened by a single key, but that key can open no other locks. No key will open all locks.



## Master Keyed

**Disk Tumbler Locks:** Available  
**Pin Tumbler Locks:** Available

All locks are keyed different from each other with individual keys for each lock. Master key provided that can open all locks in the set. Multiple sets are available, but no key can open all locks in all sets.




# Lock Keying Options

## Grand Master Keyed


- Disk Tumbler Locks:** N/A
- Pin Tumbler Locks:** Available

All locks are keyed different from each other with individual keys for each lock. Master key provided that can open all locks in the set. Multiple sets are available, with Grand master key that can open all locks in all sets.



## Sub-master Keyed & Mason Keyed

- Disk Tumbler Locks:** N/A
- Pin Tumbler Locks:** Available




### Sub-master Keyed

All locks are keyed different from each other with individual keys for each lock. Sub-master key provided that can open all locks in sub-set. Multiple sub-sets are available, with master key that can open all locks in all sub-sets. Multiple sets are available, with grand master key that can open all locks in all sets and sub-sets.

### Mason Keyed

Locks are keyed different from each other with individual keys for each lock. Each individual key can also open all common locks in set. Master key provided that can open all locks in set. Multiple sets are available with grand master key that can open all locks in all sets.

# Steel Sitting Height Base Cabinets





## Specifications:

Sitting Height Cabinets are 29" high and 22" deep. Lengths as shown.

Open and Cupboard Cabinets are furnished with one adjustable shelf.

**Note:** Last four digits of part number designate door and drawer head style, material, and pull selection. See page 4 for full explanation. Actual product may vary slightly from catalog illustration due to style and hardware selection.

## Open Cabinets






<b>B00C292212-00_0</b>	12"L	<b>B00C292230-00_0</b>	30"L
<b>B00C292215-00_0</b>	15"L	<b>B00C292236-00_0</b>	36"L
<b>B00C292218-00_0</b>	18"L	<b>B00C292242-00_0</b>	42"L
<b>B00C292221-00_0</b>	21"L	<b>B00C292248-00_0</b>	48"L
<b>B00C292224-00_0</b>	24"L		

Open Cabinets may be ordered with the topspace 7/8" deeper so that cabinet fronts line up with the doors and drawer fronts on adjacent Full Overlay cabinets. To specify the deeper toe space add a **C** to the end of the part number.




Example: **B00C292236-0000-C**

## Cupboard Cabinets







<b>C03C292215L</b> _____	15"L	<b>C03C292215-</b> _____	15"L	<b>C03C292230-</b> _____	30"L
<b>C03C292218L</b> _____	18"L	<b>C03C292218-</b> _____	18"L	<b>C03C292236-</b> _____	36"L
<b>C03C292221L</b> _____	21"L	<b>C03C292221-</b> _____	21"L	<b>C03C292242-</b> _____	42"L
<b>C03C292224L</b> _____	24"L	<b>C03C292224-</b> _____	24"L	<b>C03C292248-</b> _____	48"L

## Drawer Cabinets




<b>D00C292215-</b> _____	15"L	<b>D00C292230-</b> _____	30"L
<b>D00C292218-</b> _____	18"L	<b>D00C292236-</b> _____	36"L
<b>D00C292221-</b> _____	21"L	<b>D00C292242-</b> _____	42"L
<b>D00C292224-</b> _____	24"L	<b>D00C292248-</b> _____	48"L

**D26C292218-**\_\_\_\_\_ 18"L  
1 Pull Out Shelf



<b>D30C292215-</b> _____	15"L	<b>D30C292230-</b> _____	30"L	<b>D31C292230-</b> _____	30"L	<b>D32C292230-</b> _____	30"L
<b>D30C292218-</b> _____	18"L	<b>D30C292236-</b> _____	36"L	<b>D31C292236-</b> _____	36"L	<b>D32C292236-</b> _____	36"L
<b>D30C292221-</b> _____	21"L	<b>D30C292242-</b> _____	42"L	<b>D31C292242-</b> _____	42"L	<b>D32C292242-</b> _____	42"L
<b>D30C292224-</b> _____	24"L	<b>D30C292248-</b> _____	48"L	<b>D31C292248-</b> _____	48"L	<b>D32C292248-</b> _____	48"L

# Steel Sitting Height Base Cabinets




## Specifications:

Sitting Height Cabinets are 29" high and 22" deep. Lengths as shown.


**Note:** Last four digits of part number designate door and drawer head style, material, and pull selection. See page 4 for full explanation. Actual product may vary slightly from catalog illustration due to style and hardware selection.

## Drawer Cabinets




**D81C292218-** 18"L  
**D81C292221-** 21"L  
**D81C292224-** 24"L

2 File Drawer with hangers



**D84C292218-** 18"L  
**D84C292221-** 21"L  
**D84C292224-** 24"L


Pull Out Shelf,  
1 File Drawer with hangers




**D85C292218-** 18"L  
**D85C292221-** 21"L  
**D85C292224-** 24"L

1 File Drawer with hangers


## Cupboard and Drawer Cabinets




**E05C292230L** 30"L  
**E05C292236L** 36"L  
**E05C292242L** 42"L  
**E05C292248L** 48"L




**E05C292230-** 30"L  
**E05C292236-** 36"L  
**E05C292242-** 42"L  
**E05C292248-** 48"L



**E40C292215L** 15"L  
**E40C292218L** 18"L  
**E40C292221L** 21"L  
**E40C292224L** 24"L




**E40C292215-** 15"L  
**E40C292218-** 18"L  
**E40C292221-** 21"L  
**E40C292224-** 24"L



**E41C292230-** 30"L  
**E41C292236-** 36"L  
**E41C292242-** 42"L  
**E41C292248-** 48"L

# Steel Sitting Height Base Cabinets





### Specifications:

Sitting Height Cabinets are 29" high and 22" deep. Lengths as shown.

Sink cabinets are provided with lowered back for plumbing access. Sink support kits are available when needed to support under-mount sinks.




**Note:** Last four digits of part number designate door and drawer head style, material, and pull selection. See page 4 for full explanation. Actual product may vary slightly from catalog illustration due to style and hardware selection.

### Cupboard and Drawer Cabinets




<b>E45C292230L</b> -----	30"L	<b>E45C292230-</b> -----	30"L
<b>E45C292236L</b> -----	36"L	<b>E45C292236-</b> -----	36"L
<b>E45C292242L</b> -----	42"L	<b>E45C292242-</b> -----	42"L
<b>E45C292248L</b> -----	48"L	<b>E45C292248-</b> -----	48"L

### Sink Cabinets




<b>G00C292215L</b> -----	15"L	<b>G00C292215-</b> -----	15"L	<b>G00C292230-</b> -----	30"L
<b>G00C292218L</b> -----	18"L	<b>G00C292218-</b> -----	18"L	<b>G00C292236-</b> -----	36"L
<b>G00C292221L</b> -----	21"L	<b>G00C292221-</b> -----	21"L	<b>G00C292242-</b> -----	42"L
<b>G00C292224L</b> -----	24"L	<b>G00C292224-</b> -----	24"L	<b>G00C292248-</b> -----	48"L


### Sink Support Kits



### Diagonal Corner Cabinets



**G14C292214-**-----




Part Number	Fits
<b>K-MSS-15</b>	15" Cabinet
<b>K-MSS-18</b>	18" Cabinet
<b>K-MSS-24</b>	24" Cabinet
<b>K-MSS-30</b>	30" Cabinet
<b>K-MSS-36</b>	36" Cabinet
<b>K-MSS-42</b>	42" Cabinet
<b>K-MSS-48</b>	48" Cabinet

Sink support kit includes: four (4) threaded hanging rods and two (2) support rails and are designed to support under-mount tub sinks up to 14" deep. For sinks from 14" to 20" deep add a **-D** to the end of the part number.

Example: **K-MSS-36-D**

# Steel ADA Height Base Cabinets




## Specifications:

ADA Height Cabinets are 32" high and 22" deep. Lengths as shown.


Open and Cupboard Cabinets are furnished with one adjustable shelf.

**Note:** Last four digits of part number designate door and drawer head style, material, and pull selection. See page 4 for full explanation. Actual product may vary slightly from catalog illustration due to style and hardware selection.

## Open Cabinets



**B00C322212-00\_0** 12"L  
**B00C322215-00\_0** 15"L  
**B00C322218-00\_0** 18"L  
**B00C322221-00\_0** 21"L  
**B00C322224-00\_0** 24"L




**B00C322230-00\_0** 30"L  
**B00C322236-00\_0** 36"L  
**B00C322242-00\_0** 42"L  
**B00C322248-00\_0** 48"L


Open Cabinets may be ordered with the topspace 7/8" deeper so that cabinet fronts line up with the doors and drawer fronts on adjacent Full Overlay cabinets. To specify the deeper toe space add a **C** to the end of the part number.

Example: **G00C322236-0000-C**


## Cupboard Cabinets



**C03C322215L** 15"L  
**C03C322218L** 18"L  
**C03C322221L** 21"L  
**C03C322224L** 24"L




**C03C322215-** 15"L  
**C03C322218-** 18"L  
**C03C322221-** 21"L  
**C03C322224-** 24"L



**C03C322230-** 30"L  
**C03C322236-** 36"L  
**C03C322242-** 42"L  
**C03C322248-** 48"L



# Steel ADA Height Base Cabinets




## Specifications:


ADA Height Cabinets are 32" high and 22" deep. Lengths as shown.

**Note:** Last four digits of part number designate door and drawer head style, material, and pull selection. See page 4 for full explanation. Actual product may vary slightly from catalog illustration due to style and hardware selection.


## Drawer Cabinets




**D00C322215-**----- 15"L  
**D00C322218-**----- 18"L  
**D00C322221-**----- 21"L  
**D00C322224-**----- 24"L




**D00C322230-**----- 30"L  
**D00C322236-**----- 36"L  
**D00C322242-**----- 42"L  
**D00C322248-**----- 48"L




**D02C322230-**----- 30"L  
**D02C322236-**----- 36"L  
**D02C322242-**----- 42"L  
**D02C322248-**----- 48"L




**D30C322215-**----- 15"L  
**D30C322218-**----- 18"L  
**D30C322221-**----- 21"L  
**D30C322224-**----- 24"L




**D30C322230-**----- 30"L  
**D30C322236-**----- 36"L  
**D30C322242-**----- 42"L  
**D30C322248-**----- 48"L




**D31C322230-**----- 30"L  
**D31C322236-**----- 36"L  
**D31C322242-**----- 42"L  
**D31C322248-**----- 48"L




**D32C322230-**----- 30"L  
**D32C322236-**----- 36"L  
**D32C322242-**----- 42"L  
**D32C322248-**----- 48"L




**D36C322230-**----- 30"L  
**D36C322236-**----- 36"L  
**D36C322242-**----- 42"L  
**D36C322248-**----- 48"L




**D50C322215-**----- 15"L  
**D50C322218-**----- 18"L  
**D50C322221-**----- 21"L  
**D50C322224-**----- 24"L




**D51C322230-**----- 30"L  
**D51C322236-**----- 36"L  
**D51C322242-**----- 42"L  
**D51C322248-**----- 48"L



**D52C322230-**----- 30"L  
**D52C322236-**----- 36"L  
**D52C322242-**----- 42"L  
**D52C322248-**----- 48"L




**D83C322215-**----- 15"L  
**D83C322218-**----- 18"L  
**D83C322221-**----- 21"L  
**D83C322224-**----- 24"L



**D84C322215-**----- 15"L  
**D84C322218-**----- 18"L  
**D84C322221-**----- 21"L  
**D84C322224-**----- 24"L

# Steel ADA Height Base Cabinets




## Specifications:


ADA Height Cabinets are 32" high and 22" deep. Lengths as shown.

**Note:** Last four digits of part number designate door and drawer head style, material, and pull selection. See page 4 for full explanation. Actual product may vary slightly from catalog illustration due to style and hardware selection.


## Cupboard and Drawer Cabinets




**E05C322230L** \_\_\_\_\_ 30"L  
**E05C322236L** \_\_\_\_\_ 36"L  
**E05C322242L** \_\_\_\_\_ 42"L  
**E05C322248L** \_\_\_\_\_ 48"L




**E05C322230-** \_\_\_\_\_ 30"L  
**E05C322236-** \_\_\_\_\_ 36"L  
**E05C322242-** \_\_\_\_\_ 42"L  
**E05C322248-** \_\_\_\_\_ 48"L




**E40C322215L** \_\_\_\_\_ 15"L  
**E40C322218L** \_\_\_\_\_ 18"L  
**E40C322221L** \_\_\_\_\_ 21"L  
**E40C322224L** \_\_\_\_\_ 24"L




**E40C322215-** \_\_\_\_\_ 15"L  
**E40C322218-** \_\_\_\_\_ 18"L  
**E40C322221-** \_\_\_\_\_ 21"L  
**E40C322224-** \_\_\_\_\_ 24"L




**E40C322230-** \_\_\_\_\_ 30"L  
**E40C322236-** \_\_\_\_\_ 36"L  
**E40C322242-** \_\_\_\_\_ 42"L  
**E40C322248-** \_\_\_\_\_ 48"L




**E41C322230-** \_\_\_\_\_ 30"L  
**E41C322236-** \_\_\_\_\_ 36"L  
**E41C322242-** \_\_\_\_\_ 42"L  
**E41C322248-** \_\_\_\_\_ 48"L




**E42C322230L** \_\_\_\_\_ 30"L  
**E42C322236L** \_\_\_\_\_ 36"L  
**E42C322242L** \_\_\_\_\_ 42"L  
**E42C322248L** \_\_\_\_\_ 48"L




**E42C322230-** \_\_\_\_\_ 30"L  
**E42C322236-** \_\_\_\_\_ 36"L  
**E42C322242-** \_\_\_\_\_ 42"L  
**E42C322248-** \_\_\_\_\_ 48"L




**E43C322230L** \_\_\_\_\_ 30"L  
**E43C322236L** \_\_\_\_\_ 36"L  
**E43C322242L** \_\_\_\_\_ 42"L  
**E43C322248L** \_\_\_\_\_ 48"L




**E43C322230-** \_\_\_\_\_ 30"L  
**E43C322236-** \_\_\_\_\_ 36"L  
**E43C322242-** \_\_\_\_\_ 42"L  
**E43C322248-** \_\_\_\_\_ 48"L



**E45C322230L** \_\_\_\_\_ 30"L  
**E45C322236L** \_\_\_\_\_ 36"L  
**E45C322242L** \_\_\_\_\_ 42"L  
**E45C322248L** \_\_\_\_\_ 48"L




**E45C322230-** \_\_\_\_\_ 30"L  
**E45C322236-** \_\_\_\_\_ 36"L  
**E45C322242-** \_\_\_\_\_ 42"L  
**E45C322248-** \_\_\_\_\_ 48"L



**E56C322230-** \_\_\_\_\_ 30"L  
**E56C322236-** \_\_\_\_\_ 36"L  
**E56C322242-** \_\_\_\_\_ 42"L  
**E56C322248-** \_\_\_\_\_ 48"L

# Steel ADA Height Base Cabinets




**Specifications:**


ADA Height Cabinets are 32" high and 22" deep. Lengths as shown.

**Note:** Last four digits of part number designate door and drawer head style, material, and pull selection. See page 4 for full explanation. Actual product may vary slightly from catalog illustration due to style and hardware selection.


**Cupboard and Drawer Cabinets**




- E62C322230L \_\_\_\_\_ 30"L
- E62C322236L \_\_\_\_\_ 36"L
- E62C322242L \_\_\_\_\_ 42"L
- E62C322248L \_\_\_\_\_ 48"L




- E62C322230- \_\_\_\_\_ 30"L
- E62C322236- \_\_\_\_\_ 36"L
- E62C322242- \_\_\_\_\_ 42"L
- E62C322248- \_\_\_\_\_ 48"L




- E63C322230L \_\_\_\_\_ 30"L
- E63C322236L \_\_\_\_\_ 36"L
- E63C322242L \_\_\_\_\_ 42"L
- E63C322248L \_\_\_\_\_ 48"L



- E63C322230- \_\_\_\_\_ 30"L
- E63C322236- \_\_\_\_\_ 36"L
- E63C322242- \_\_\_\_\_ 42"L
- E63C322248- \_\_\_\_\_ 48"L




- E65C322230L \_\_\_\_\_ 30"L
- E65C322236L \_\_\_\_\_ 36"L
- E65C322242L \_\_\_\_\_ 42"L
- E65C322248L \_\_\_\_\_ 48"L



- E65C322230- \_\_\_\_\_ 30"L
- E65C322236- \_\_\_\_\_ 36"L
- E65C322242- \_\_\_\_\_ 42"L
- E65C322248- \_\_\_\_\_ 48"L

# Steel ADA Height Base Cabinets




**Specifications:**


ADA Height Cabinets are 32" high and 22" deep. Lengths as shown.

**Note:** Last four digits of part number designate door and drawer head style, material, and pull selection. See page 4 for full explanation. Actual product may vary slightly from catalog illustration due to style and hardware selection.


**Blind Corner Cabinets**



G12C322248-\_\_\_\_\_ 48"L




PLAN VIEW OF  
G12C322248- & G13C322248-




G13C322248-\_\_\_\_\_ 48"L

**Acid Storage Cabinets**




20" Deep

G21C322018L \_\_\_\_\_ 18"L  
G21C322024L \_\_\_\_\_ 24"L



20" Deep

G21C320218- \_\_\_\_\_ 18"L  
G21C320224- \_\_\_\_\_ 24"L



20" Deep

G21C322030- \_\_\_\_\_ 30"L  
G21C322036- \_\_\_\_\_ 36"L  
G21C322042- \_\_\_\_\_ 42"L  
G21C322048- \_\_\_\_\_ 48"L

22" Deep

G21C322218L \_\_\_\_\_ 18"L  
G21C322224L \_\_\_\_\_ 24"L

22" Deep


G21C322218- \_\_\_\_\_ 18"L  
G21C322224- \_\_\_\_\_ 24"L

22" Deep


G21C322230- \_\_\_\_\_ 30"L  
G21C322236- \_\_\_\_\_ 36"L  
G21C322242- \_\_\_\_\_ 42"L  
G21C322248- \_\_\_\_\_ 48"L

**Note:** G21 Acid Storage Cabinets are available with five-knuckle hinges only


**Diagonal Corner Cabinet**




G14C322214-\_\_\_\_\_



G15C322214-\_\_\_\_\_



# Steel ADA Height Base Cabinets




### Specifications:

ADA Height Cabinets are 32" high and 22" deep. Lengths as shown.


Sink cabinets are provided with lowered back for plumbing access. Sink support kits are available when needed to support under-mount sinks.

**Note:** Last four digits of part number designate door and drawer head style, material, and pull selection. See page 4 for full explanation. Actual product may vary slightly from catalog illustration due to style and hardware selection.


### Sink Cabinets



- G00C322212L \_\_\_\_\_ 12"L
- G00C322215L \_\_\_\_\_ 15"L
- G00C322218L \_\_\_\_\_ 18"L
- G00C322221L \_\_\_\_\_ 21"L
- G00C322224L \_\_\_\_\_ 24"L




- G00C322212- \_\_\_\_\_ 12"L
- G00C322215- \_\_\_\_\_ 15"L
- G00C322218- \_\_\_\_\_ 18"L
- G00C322221- \_\_\_\_\_ 21"L
- G00C322224- \_\_\_\_\_ 24"L



- G00C322230- \_\_\_\_\_ 30"L
- G00C322236- \_\_\_\_\_ 36"L
- G00C322242- \_\_\_\_\_ 42"L
- G00C322248- \_\_\_\_\_ 48"L

### Sink Support Kits




Part Number	Fits
K-MSS-12	12" Cabinet
K-MSS-15	15" Cabinet
K-MSS-18	18" Cabinet
K-MSS-24	24" Cabinet
K-MSS-30	30" Cabinet
K-MSS-36	36" Cabinet
K-MSS-42	42" Cabinet
K-MSS-48	48" Cabinet

Sink support kit includes: four (4) threaded hanging rods and two (2) support rails and are designed to support under-mount tub sinks up to 14" deep. For sinks from 14" to 20" deep add a **-D** to the end of the part number.


Example: **K-MSS-36-D**

### Solvent Storage Cabinets




20" Deep

- G63C322018L\_100 18"L
- G63C322024L\_100 24"L



20" Deep

- G63C322018-100 18"L
- G63C322024-100 24"L



20" Deep

- G63C322030-100 30"L
- G63C322036-100 36"L
- G63C322048-100 48"L

22" Deep

- G63C322218L\_100 18"L
- G63C322224L\_100 24"L

22" Deep

- G63C322218-100 18"L
- G63C322224-100 24"L

22" Deep

- G63C322230-100 30"L
- G63C322236-100 36"L
- G63C322248-100 48"L

G63 Solvent Storage Cabinets are available as mobile.

To add casters to the cabinets above, add **-B7** to the end of the part number.

Example: **G63C322230-0100-B7** — Mobile, Inset style cabinet, 22" deep by 30" long.


### Bin Storage Cabinet



- G23C322218- \_\_\_\_\_ 18"L

Furnished with 11" high drawer body with full height drawer front mounted on full extension slides.

# Steel ADA Height Base Cabinets




### Specifications:

ADA Height Fume Hood Cabinets are 32" high and either 20" or 22" deep. Lengths as shown.

**Note:** Last four digits of part number designate door and drawer head style, material, and pull selection. See page 4 for full explanation. Actual product may vary slightly from catalog illustration due to style and hardware selection.

## Fume Hood Base Cabinets




20" Deep

G08C322012L	12"L
G08C322015L	15"L
G08C322018L	18"L
G08C322021L	21"L
G08C322024L	24"L

22" Deep

G08C322212L	12"L
G08C322215L	15"L
G08C322218L	18"L
G08C322221L	21"L
G08C322224L	24"L




20" Deep

G08C322012-	12"L
G08C322015-	15"L
G08C322018-	18"L
G08C322021-	21"L
G08C322024-	24"L

22" Deep

G08C322212-	12"L
G08C322215-	15"L
G08C322218-	18"L
G08C322221-	21"L
G08C322224-	24"L




20" Deep

G08C322030-	30"L
G08C322036-	36"L
G08C322042-	42"L
G08C322048-	48"L

22" Deep

G08C322230-	30"L
G08C322236-	36"L
G08C322242-	42"L
G08C322248-	48"L

## Fume Hood Solvent Storage Base Cabinets




20" Deep

G68C322024L_100	24"L
-----------------	------

22" Deep

G68C322224L_100	24"L
-----------------	------




20" Deep

G68C322024-_100	24"L
-----------------	------

22" Deep

G68C322224-_100	24"L
-----------------	------



20" Deep


G68C322030-_100	30"L
G68C322036-_100	36"L
G68C322048-_100	48"L

22" Deep

G68C322230-_100	30"L
G68C322236-_100	36"L
G68C322248-_100	48"L

G68 Solvent Storage Cabinets are available as mobile. To add casters to the cabinets above, add **-B7** to the end of the part number. Example: **G68C322018-0100-B7** — Mobile, Inset style cabinet, 22" deep by 18" long.

## Fume Hood Vacuum Storage Base Cabinets




20" Deep

G35C322024L	24"L
-------------	------

22" Deep

G35C322224L	24"L
-------------	------




20" Deep

G35C322024-	24"L
-------------	------

22" Deep

G35C322224-	24"L
-------------	------



20" Deep


G35C322030-	30"L
G35C322036-	36"L
G35C322048-	48"L

22" Deep

G35C322230-	30"L
G35C322236-	36"L
G35C322248-	48"L


# Steel ADA Height Base Cabinets

## Fume Hood Acid Storage Base Cabinets




20" Deep

<b>G80C322018L</b> _____	24"L
<b>G80C322024L</b> _____	24"L



20" Deep

<b>G80C322018-</b> _____	18"L
<b>G80C322024-</b> _____	24"L



20" Deep

<b>G80C322030-</b> _____	30"L
<b>G80C322036-</b> _____	36"L
<b>G80C322048-</b> _____	48"L

22" Deep

<b>G80C322218L</b> _____	18"L
<b>G80C322224L</b> _____	24"L

22" Deep


<b>G80C322218-</b> _____	18"L
<b>G80C322224-</b> _____	24"L

22" Deep

<b>G80C322230-</b> _____	30"L
<b>G80C322236-</b> _____	36"L
<b>G80C322248-</b> _____	48"L


**Note:** G80 Acid Storage Cabinets are available with five-knuckle hinges only

## Wheelchair Accessible Sink Cabinet with Doors




The Wheelchair Accessible Base Cabinet with Doors is designed to be used as a sink unit, but may be used anywhere wheelchair accessibility is needed. The cabinet is constructed without a bottom, and with the toe-

space rails attached to the doors so it will match adjacent cabinets when the doors are closed. The cabinet incorporates a 2<sup>3</sup>/<sub>4</sub>" top rail to give a minimum of 29<sup>1</sup>/<sub>4</sub>" of clearance from the floor to the bottom of the rail.



<b>G03C322236-</b> _____	36"L
<b>G03C322242-</b> _____	42"L
<b>G03C322248-</b> _____	48"L

## Wheelchair Accessible Base Cabinet without Doors




The Wheelchair Accessible Base Cabinet without Doors is designed to be used as a sink unit, but may be used anywhere wheelchair accessibility is needed. The cabinet is constructed without a bottom, and without doors for easy accessibility.

The cabinet incorporates a 4<sup>1</sup>/<sub>2</sub>" top rail to give a minimum of 27<sup>1</sup>/<sub>2</sub>" of clearance from the floor to the bottom of the rail and a removable sloping knee space panel to enclose piping and drain lines.


Wheelchair Accessible Base Cabinets without Doors may be ordered with the topspace 7/8" deeper so that cabinet fronts line up with the doors and drawer fronts on adjacent Full Overlay cabinets. To specify the deeper toe space, add a "C" to the end of the part number.

Example: **G07C322236-0000-C**



<b>G07C322236-00_0</b>	36"L
<b>G07C322242-00_0</b>	42"L
<b>G07C322248-00_0</b>	48"L

# Steel Standing Height Base Cabinets




## Specifications:

Standing Height Cabinets are 35" high and 22" deep. Lengths as shown.


Open and Cupboard Cabinets are furnished with one adjustable shelf.

**Note:** Last four digits of part number designate door and drawer head style, material, and pull selection. See page 4 for full explanation. Actual product may vary slightly from catalog illustration due to style and hardware selection.

## Open Cabinets



**B00C352212-00\_0** 12"L  
**B00C352215-00\_0** 15"L  
**B00C352218-00\_0** 18"L  
**B00C352221-00\_0** 21"L  
**B00C352224-00\_0** 24"L




**B00C352230-00\_0** 30"L  
**B00C352236-00\_0** 36"L  
**B00C352242-00\_0** 42"L  
**B00C352248-00\_0** 48"L


Open Cabinets may be ordered with the topspace 7/8" deeper so that cabinet fronts line up with the doors and drawer fronts on adjacent Full Overlay cabinets. To specify the deeper toe space add a **C** to the end of the part number.

Example: **G00C352236-0000-C**


## Cupboard Cabinets



**C03C352212L** 12"L  
**C03C352215L** 15"L  
**C03C352218L** 18"L  
**C03C352221L** 21"L  
**C03C352224L** 24"L




**C03C352212-** 12"L  
**C03C352215-** 15"L  
**C03C352218-** 18"L  
**C03C352221-** 21"L  
**C03C352224-** 24"L




**C03C352230-** 30"L  
**C03C352236-** 36"L  
**C03C352242-** 42"L  
**C03C352248-** 48"L


## Drawer Cabinets




**D00C352215-** 15"L  
**D00C352218-** 18"L  
**D00C352221-** 21"L  
**D00C352224-** 24"L




**D00C352230-** 30"L  
**D00C352236-** 36"L  
**D00C352242-** 42"L  
**D00C352248-** 48"L



**D02C352230-** 30"L  
**D02C352236-** 36"L  
**D02C352242-** 42"L  
**D02C352248-** 48"L




**D08C352215-** 15"L  
**D08C352218-** 18"L  
**D08C352221-** 21"L  
**D08C352224-** 24"L



**D08C352230-** 30"L  
**D08C352236-** 36"L  
**D08C352242-** 42"L  
**D08C352248-** 48"L



# Steel Standing Height Base Cabinets




**Specifications:**


Standing Height Cabinets are 35" high and 22" deep. Lengths as shown.

**Note:** Last four digits of part number designate door and drawer head style, material, and pull selection. See page 4 for full explanation. Actual product may vary slightly from catalog illustration due to style and hardware selection.


**Drawer Cabinets**




- D30C352212- 12"L
- D30C352215- 15"L
- D30C352218- 18"L
- D30C352221- 21"L
- D30C352224- 24"L




- D30C352230- 30"L
- D30C352236- 36"L
- D30C352242- 42"L
- D30C352248- 48"L




- D31C352230- 30"L
- D31C352236- 36"L
- D31C352242- 42"L
- D31C352248- 48"L




- D32C352230- 30"L
- D32C352236- 36"L
- D32C352242- 42"L
- D32C352248- 48"L




- D36C352230- 30"L
- D36C352236- 36"L
- D36C352242- 42"L
- D36C352248- 48"L




- D50C352215- 15"L
- D50C352218- 18"L
- D50C352221- 21"L
- D50C352224- 24"L




- D51C352230- 30"L
- D51C352236- 36"L
- D51C352242- 42"L
- D51C352248- 48"L



- D52C352230- 30"L
- D52C352236- 36"L
- D52C352242- 42"L
- D52C352248- 48"L



- D81C352218- 18"L
- D81C352221- 21"L
- D81C352224- 24"L




- D84C352218- 18"L
- D84C352221- 21"L
- D84C352224- 24"L

2 File Drawers with hangers

Pull Out Shelf  
1 File Drawer with hangers

# Steel Standing Height Base Cabinets




**Specifications:**

Standing Height Cabinets are 35" high and 22" deep. Lengths as shown.


Cupboards have one adjustable shelf.

**Note:** Last four digits of part number designate door and drawer head style, material, and pull selection. See page 4 for full explanation. Actual product may vary slightly from catalog illustration due to style and hardware selection.


**Cupboard and Drawer Cabinets**




- E05C352230L \_\_\_\_\_ 30"L
- E05C352236L \_\_\_\_\_ 36"L
- E05C352242L \_\_\_\_\_ 42"L
- E05C352248L \_\_\_\_\_ 48"L




- E05C352230- \_\_\_\_\_ 30"L
- E05C352236- \_\_\_\_\_ 36"L
- E05C352242- \_\_\_\_\_ 42"L
- E05C352248- \_\_\_\_\_ 48"L




- E40C352212L \_\_\_\_\_ 12"L
- E40C352215L \_\_\_\_\_ 15"L
- E40C352218L \_\_\_\_\_ 18"L
- E40C352221L \_\_\_\_\_ 21"L
- E40C352224L \_\_\_\_\_ 24"L




- E40C352212- \_\_\_\_\_ 12"L
- E40C352215- \_\_\_\_\_ 15"L
- E40C352218- \_\_\_\_\_ 18"L
- E40C352221- \_\_\_\_\_ 21"L
- E40C352224- \_\_\_\_\_ 24"L




- E40C352230- \_\_\_\_\_ 30"L
- E40C352236- \_\_\_\_\_ 36"L
- E40C352242- \_\_\_\_\_ 42"L
- E40C352248- \_\_\_\_\_ 48"L




- E41C352230- \_\_\_\_\_ 30"L
- E41C352236- \_\_\_\_\_ 36"L
- E41C352242- \_\_\_\_\_ 42"L
- E41C352248- \_\_\_\_\_ 48"L




- E42C352230L \_\_\_\_\_ 30"L
- E42C352236L \_\_\_\_\_ 36"L
- E42C352242L \_\_\_\_\_ 42"L
- E42C352248L \_\_\_\_\_ 48"L




- E42C352230- \_\_\_\_\_ 30"L
- E42C352236- \_\_\_\_\_ 36"L
- E42C352242- \_\_\_\_\_ 42"L
- E42C352248- \_\_\_\_\_ 48"L




- E43C352230L \_\_\_\_\_ 30"L
- E43C352236L \_\_\_\_\_ 36"L
- E43C352242L \_\_\_\_\_ 42"L
- E43C352248L \_\_\_\_\_ 48"L




- E43C352230- \_\_\_\_\_ 30"L
- E43C352236- \_\_\_\_\_ 36"L
- E43C352242- \_\_\_\_\_ 42"L
- E43C352248- \_\_\_\_\_ 48"L



- E45C352230L \_\_\_\_\_ 30"L
- E45C352236L \_\_\_\_\_ 36"L
- E45C352242L \_\_\_\_\_ 42"L
- E45C352248L \_\_\_\_\_ 48"L




- E45C352230- \_\_\_\_\_ 30"L
- E45C352236- \_\_\_\_\_ 36"L
- E45C352242- \_\_\_\_\_ 42"L
- E45C352248- \_\_\_\_\_ 48"L



- E56C352230- \_\_\_\_\_ 30"L
- E56C352236- \_\_\_\_\_ 36"L
- E56C352242- \_\_\_\_\_ 42"L
- E56C352248- \_\_\_\_\_ 48"L

# Steel Standing Height Base Cabinets




### Specifications:

Standing Height Cabinets are 35" high and 22" deep. Lengths as shown.


Cupboards have one adjustable shelf.

**Note:** Last four digits of part number designate door and drawer head style, material, and pull selection. See page 4 for full explanation. Actual product may vary slightly from catalog illustration due to style and hardware selection.


### Cupboard and Drawer Cabinets




- E60C352212L \_\_\_\_\_ 12"L
- E60C352215L \_\_\_\_\_ 15"L
- E60C352218L \_\_\_\_\_ 18"L
- E60C352221L \_\_\_\_\_ 21"L
- E60C352224L \_\_\_\_\_ 24"L




- E60C352212- \_\_\_\_\_ 12"L
- E60C352215- \_\_\_\_\_ 15"L
- E60C352218- \_\_\_\_\_ 18"L
- E60C352221- \_\_\_\_\_ 21"L
- E60C352224- \_\_\_\_\_ 24"L




- E61C352230- \_\_\_\_\_ 30"L
- E61C352236- \_\_\_\_\_ 36"L
- E61C352242- \_\_\_\_\_ 42"L
- E61C352248- \_\_\_\_\_ 48"L




- E62C352230L \_\_\_\_\_ 30"L
- E62C352236L \_\_\_\_\_ 36"L
- E62C352242L \_\_\_\_\_ 42"L
- E62C352248L \_\_\_\_\_ 48"L




- E62C352230- \_\_\_\_\_ 30"L
- E62C352236- \_\_\_\_\_ 36"L
- E62C352242- \_\_\_\_\_ 42"L
- E62C352248- \_\_\_\_\_ 48"L




- E63C352230L \_\_\_\_\_ 30"L
- E63C352236L \_\_\_\_\_ 36"L
- E63C352242L \_\_\_\_\_ 42"L
- E63C352248L \_\_\_\_\_ 48"L



- E63C352230- \_\_\_\_\_ 30"L
- E63C352236- \_\_\_\_\_ 36"L
- E63C352242- \_\_\_\_\_ 42"L
- E63C352248- \_\_\_\_\_ 48"L




- E65C352230L \_\_\_\_\_ 30"L
- E65C352236L \_\_\_\_\_ 36"L
- E65C352242L \_\_\_\_\_ 42"L
- E65C352248L \_\_\_\_\_ 48"L



- E65C352230- \_\_\_\_\_ 30"L
- E65C352236- \_\_\_\_\_ 36"L
- E65C352242- \_\_\_\_\_ 42"L
- E65C352248- \_\_\_\_\_ 48"L

# Steel Standing Height Base Cabinets




## Specifications:

Standing Height Cabinets are 35" high and 22" deep. Lengths as shown.


Sink cabinets are provided with lowered back for plumbing access except where noted. Sink support kits are available when needed to support under-mount sinks.

**Note:** Last four digits of part number designate door and drawer head style, material, and pull selection. See page 4 for full explanation. Actual product may vary slightly from catalog illustration due to style and hardware selection.


## Sink Cabinets




<b>G00C352212L</b> -----	12"L
<b>G00C352215L</b> -----	15"L
<b>G00C352218L</b> -----	18"L
<b>G00C352221L</b> -----	21"L
<b>G00C352224L</b> -----	24"L




<b>G00C352212-</b> -----	12"L
<b>G00C352215-</b> -----	15"L
<b>G00C352218-</b> -----	18"L
<b>G00C352221-</b> -----	21"L
<b>G00C352224-</b> -----	24"L



<b>G00C352230-</b> -----	30"L
<b>G00C352236-</b> -----	36"L
<b>G00C352242-</b> -----	42"L
<b>G00C352248-</b> -----	48"L




<b>G00C352260-</b> -----	60"L
--------------------------	------



<b>G03C352230-</b> -----	30"L
<b>G03C352236-</b> -----	36"L
<b>G03C352242-</b> -----	42"L
<b>G03C352248-</b> -----	48"L

G03 Cabinets are constructed without a bottom and with the toe-space rail attached to the doors so it will match adjacent cabinets when the doors are closed.


Furnished with full height back



<b>G04C352252-</b> -----	52"L
<b>G04C352258-</b> -----	58"L
<b>G04C352260-</b> -----	60"L

Center Bench End Sink Cabinet  
Furnished with full height back

## Sink Support Kits



Part Number	Fits
<b>K-MSS-12</b>	12" Cabinet
<b>K-MSS-15</b>	15" Cabinet
<b>K-MSS-18</b>	18" Cabinet
<b>K-MSS-21</b>	21" Cabinet
<b>K-MSS-24</b>	24" Cabinet
<b>K-MSS-30</b>	30" Cabinet
<b>K-MSS-36</b>	36" Cabinet
<b>K-MSS-42</b>	42" Cabinet
<b>K-MSS-48</b>	48" Cabinet
<b>K-MSS-52</b>	52" Cabinet
<b>K-MSS-58</b>	58" Cabinet
<b>K-MSS-60</b>	60" Cabinet

Sink support kit includes: four (4) threaded hanging rods and two (2) support rails and are designed to support under-mount tub sinks up to 14" deep. For sinks from 14" to 20" deep add a **-D** to the end of the part number.

Example: **K-MSS-36-D**

# Steel Standing Height Base Cabinets


## Specifications

Acid Storage Base Cabinets are specifically designed for the storage of corrosive chemicals. These cabinets are lined with 1/8" thick corrosion resistant polyethylene sheets.

The cabinet doors are provided with louvers and are latched using a nylon roller catch.

An optional 1" deep polyethylene pan is available to contain spills.

Each cabinet is furnished with a PVC vent pipe for venting. (Requires a 2" hole in the counter top above.)




## Optional Polyethylene Pan

Seamless 1" deep polyethylene pan designed to contain spills.

Part Number	Fits 20" Deep
F-8440-15	15" Cabinet
F-8440-18	18" Cabinet
F-8440-21	21" Cabinet
F-8440-24	24" Cabinet
F-8440-30	30" Cabinet
F-8440-36	36" Cabinet
F-8440-42	42" Cabinet
F-8440-48	48" Cabinet


Part Number	Fits 22" Deep
F-8441-15	15" Cabinet
F-8441-18	18" Cabinet
F-8441-21	21" Cabinet
F-8441-24	24" Cabinet
F-8441-30	30" Cabinet
F-8441-36	36" Cabinet
F-8441-42	42" Cabinet
F-8441-48	48" Cabinet

## Acid Storage Cabinet




20" Deep

G21C352015L	15"L
G21C352018L	18"L
G21C352021L	21"L
G21C352024L	24"L



20" Deep

G21C352015-	15"L
G21C352018-	18"L
G21C352021-	21"L
G21C352024-	24"L



20" Deep

G21C352030-	30"L
G21C352036-	36"L
G21C352042-	42"L
G21C352048-	48"L

22" Deep	
G21C352215L	15"L
G21C352218L	18"L
G21C352221L	21"L
G21C352224L	24"L

22" Deep	
G21C352215-	15"L
G21C352218-	18"L
G21C352221-	21"L
G21C352224-	24"L

22" Deep	
G21C352230-	30"L
G21C352236-	36"L
G21C352242-	42"L
G21C352248-	48"L

Note: G21 Acid Storage Cabinets are available with five-knuckle hinges only

# Steel Standing Height Base Cabinets

## Specifications

Solvent Storage Cabinets are UL listed for the storage of flammable and combustible liquids, meeting UFC, OSHA and NFPA No. 30-1993 standards.


A 2" deep, steel, liquid-tight pan covers the entire bottom to contain liquid leaks and spills. The cabinet is provided with a full depth adjustable shelf and two

diametrically opposed vents with spark screens in the back of the cabinet for cases when ventilation is required.


Steel Solvent Storage cabinets are all 18 gauge steel, double panel construction with self-closing doors, synchronized so that both doors will always fully close. The right hand door is equipped

with a three-point latching system that automatically engages the cabinet frame. Each door is equipped with a fusible-link hold-open feature that ensures the doors close should the temperature outside the cabinet exceed 165°F. Cabinets are provided with a grounding screw at the rear.


All Solvent Storage Cabinets are labeled: CAUTION FLAMMABLE – KEEP FIRE AWAY in English, Spanish, and French.




## Vent Locations




## Solvent Storage Cabinets



**20" Deep** Capacity  
**G63C352024L\_100** 30"L 12 gal



**22" Deep** Capacity  
**G63C352224L\_100** 30"L 12 gal



**20" Deep** Capacity  
**G63C352030-\_100** 30"L 12 gal  
**G63C352036-\_100** 36"L 15 gal  
**G63C352048-\_100** 48"L 20 gal


**22" Deep** Capacity  
**G63C352230-\_100** 30"L 12 gal  
**G63C352236-\_100** 36"L 15 gal  
**G63C352248-\_100** 48"L 20 gal

G63 Solvent Storage Cabinets are available as mobile.

To add casters to the cabinets above, add **-B7** to the end of the part number.

Example: **G63C352036-0100-B7** — Mobile, Inset style cabinet, 20" deep by 36" long.

# Steel Standing Height Base Cabinets




**Specifications:**

Standing Height Cabinets are 35" high and 22" deep. Lengths as shown.


Blind Corner Cabinets have one adjustable shelf.

**Note:** Last four digits of part number designate door and drawer head style, material, and pull selection. See page 4 for full explanation. Actual product may vary slightly from catalog illustration due to style and hardware selection.


**Blind Corner Cabinets**



G12C352248-\_\_\_\_\_ 48"L




PLAN VIEW OF  
G12C352248- & G13C352248-




G13C352248-\_\_\_\_\_ 48"L

**Diagonal Corner Cabinets**




G14C352214-\_\_\_\_\_




G15C352214-\_\_\_\_\_

Provided with Lazy Susan



**Waste Storage Cabinet      Bin Storage Cabinet**



- G20C352215-\_\_\_\_\_ 15"L
- G20C352218-\_\_\_\_\_ 18"L
- G20C352221-\_\_\_\_\_ 21"L
- G20C352224-\_\_\_\_\_ 21"L


Cabinet furnished with 13" deep, cold rolled steel bin, mounted on full extension drawer slides for ease of access. Includes 6" x 5" (opening) stainless steel waste chute



- G23C352218-\_\_\_\_\_ 18"L

Furnished with 11" high drawer body with full height drawer front mounted on full extension slides.

# Fume Hood Base Cabinets




## Specifications

Standing Height Steel Fume Hood Base Cabinets are 35" high and 22" deep; ADA Height Steel Fume Hood Base Cabinets are 32" high and 22" deep; lengths as shown.

## Cabinet Style Option:


Specify cabinet door and hardware style when ordering by replacing blanks in last four digits of catalog number with style numbers.

## Standing Height Fume Hood Base Cabinets




20" Deep

<b>G08C352012L</b> -----	12"L
<b>G08C352015L</b> -----	15"L
<b>G08C352018L</b> -----	18"L
<b>G08C352021L</b> -----	21"L
<b>G08C352024L</b> -----	24"L




20" Deep

<b>G08C352012-</b> -----	12"L
<b>G08C352015-</b> -----	15"L
<b>G08C352018-</b> -----	18"L
<b>G08C352021-</b> -----	21"L
<b>G08C352024-</b> -----	24"L



20" Deep

<b>G08C352030-</b> -----	30"L
<b>G08C352036-</b> -----	36"L
<b>G08C352042-</b> -----	42"L
<b>G08C352048-</b> -----	48"L



20" Deep

<b>G08C352060-</b> -----	60"L
--------------------------	------

22" Deep

<b>G08C352212L</b> -----	12"L
<b>G08C352215L</b> -----	15"L
<b>G08C352218L</b> -----	18"L
<b>G08C352221L</b> -----	21"L
<b>G08C352224L</b> -----	24"L

22" Deep

<b>G08C352212-</b> -----	12"L
<b>G08C352215-</b> -----	15"L
<b>G08C352218-</b> -----	18"L
<b>G08C352221-</b> -----	21"L
<b>G08C352224-</b> -----	24"L

22" Deep

<b>G08C352230-</b> -----	30"L
<b>G08C352236-</b> -----	36"L
<b>G08C352242-</b> -----	42"L
<b>G08C352248-</b> -----	48"L

22" Deep

<b>G08C352260-</b> -----	60"L
--------------------------	------




# Acid Storage Fume Hood Base Cabinets

## Specifications

Acid Storage Fume Hood Base Cabinets are specifically designed for the storage of corrosive chemicals. These cabinets are lined with a molded one piece linear low density polyethylene tub with coved corners

and a 1" lip at the bottom front. The cabinet doors are lined with 1/8" sheet polyethylene and the doors are latched using a nylon roller catch.


Each cabinet is furnished with a PVC vent pipe for venting to the fume hood above. (Requires a 2" hole in the hood worksurface.)




## Optional Phenolic Resin Removable Half-Depth Shelf

Part No.	
019623	Shelf for 18" cabinet
019624	Shelf for 24" cabinet
019625	Shelf for 30" cabinet
019626	Shelf for 36" cabinet
019627	Shelf for 48" cabinet


## Standing Height Acid Storage Fume Hood Cabinets




**Vent Location**  
On Center  
11 5/8" from Top of Cabinet



**Vent Location**  
On Center  
11 5/8" from Top of Cabinet



**Vent Location**  
30" - 7 5/8" Right of Center  
36" & 48" - 5" Right of Center  
11 5/8" from Top of Cabinet



**Vent Location**  
On Center  
11 5/8" from Top of Cabinet

20" Deep	
G80C352018L	18"L
G80C352024L	24"L

20" Deep	
G80C352018-	18"L
G80C352024-	24"L

20" Deep	
G80C352030-	30"L
G80C352036-	36"L
G80C352048-	48"L

20" Deep	
G80C352060-	60"L

22" Deep	
G80C352218L	18"L
G80C352224L	24"L

22" Deep	
G80C352218-	18"L
G80C352224-	24"L

22" Deep	
G80C352230-	30"L
G80C352236-	36"L
G80C352248-	48"L

22" Deep	
G80C352260-	60"L

**Note:** G80 Acid Storage Cabinets are available with five-knuckle hinges only

# Solvent Storage Fume Hood Base Cabinets

## Specifications

Solvent Storage Cabinets are UL listed for the storage of flammable and combustible liquids, meeting UFC, OSHA and NFPA No. 30-1993 standards.


A 2" deep, steel, liquid-tight pan covers the entire bottom to contain liquid leaks

and spills. The cabinet is provided with a full depth adjustable shelf and two diametrically opposed vents with spark screens in the back of the cabinet for cases when ventilation is required.


Steel Solvent Storage cabinets are all 18

gauge steel, double panel construction with self-closing doors, synchronized so that both doors will always fully close. The right hand door is equipped with a three-point latching system that automatically engages the cabinet frame. Each door is equipped with a fusible-link hold-open feature that ensures the doors close should the temperature outside the cabinet exceed 165°F. Cabinets are provided with a grounding screw at the rear.


All Solvent Storage Cabinets are labeled: CAUTION FLAMMABLE – KEEP FIRE AWAY in English, Spanish, and French.




## Vent Locations




## Fume Hood Solvent Storage Base Cabinets



**20" Deep** Capacity  
**G68C352024L\_100** 30"L 12 gal



**20" Deep** Capacity  
**G68C352024-\_100** 30"L 12 gal



**20" Deep** Capacity  
**G68C352030-\_100** 30"L 12 gal  
**G68C352036-\_100** 36"L 15 gal  
**G68C352048-\_100** 48"L 20 gal

**22" Deep** Capacity  
**G68C352224L\_100** 30"L 12 gal

**22" Deep** Capacity  
**G68C352224-\_100** 30"L 12 gal

**22" Deep** Capacity  
**G68C352230-\_100** 30"L 12 gal  
**G68C352236-\_100** 36"L 15 gal  
**G68C352248-\_100** 48"L 20 gal

G68 Fume Hood Solvent Storage Cabinets are available as mobile.  
To add casters to the cabinets above, add **-B7** to the end of the part number.  
Example: **G68C352224-1100-B7** — Mobile, Overlay style cabinet, 22" deep by 24" long.


# Vacuum Pump Storage Base Cabinets

## Specifications


Vacuum Pump Storage Fume Hood Cabinets are designed without a bottom to allow vacuum pumps and other equipment to be rolled in or out of the cabinets. The interior is lined with an easy to clean sound absorbant material. Each cabinet is furnished with

a 120 VAC, 20 amp, duplex receptacle mounted on the inside cabinet back and a pilot lighted toggle switch mounted in the top front rail. (Wiring is not included.) Cabinets are also furnished with a 1 1/2" diameter PVC vent pipe for venting or access to the

hood above. (Requires a 2" hole in the hood worksurface.) The toespace rail is attached to the door to allow total access to the cabinet. Cabinet inside clearance at the floor is 14 1/2" or 16 1/2" front-to-back, 25" high, and 3" less than the overall cabinet length.




## Fume Hood Vacuum Pump Storage Cabinets




20" Deep

G35C352024L \_\_\_\_\_ 24"L



20" Deep

G35C352024- \_\_\_\_\_ 24"L



20" Deep

G35C352030- \_\_\_\_\_ 30"L  
G35C352036- \_\_\_\_\_ 36"L  
G35C352048- \_\_\_\_\_ 48"L

22" Deep

G35C352224L \_\_\_\_\_ 24"L

22" Deep


G35C352224- \_\_\_\_\_ 24"L

22" Deep

G35C352230- \_\_\_\_\_ 30"L  
G35C352236- \_\_\_\_\_ 36"L  
G35C352248- \_\_\_\_\_ 48"L

## Vacuum Pump Carts

Vacuum Pump Carts are designed for use with Kewaunee's Vacuum Pump cabinets shown on this page. The cart consists of a 1" deep pan, fabricated from powder coated or stainless steel and mounted on 4" swivel casters.




### Carts for 20" Deep Cabinets

Part Number	Length
VPCC061414-00_0	14"
VPCC061420-00_0	20"
VPCC061426-00_0	26"
VPCC061432-00_0	32"
VPCC061438-00_0	38"
VPCC061444-00_0	44"

### Carts for 22" Deep Cabinets

Part Number	Length
VPCC061614-00_0	14"
VPCC061620-00_0	20"
VPCC061626-00_0	26"
VPCC061632-00_0	32"
VPCC061638-00_0	38"
VPCC061644-00_0	44"

# Steel Sitting Height Mobile Cabinets





## Specifications

Steel Sitting Height Mobile Cabinets are 26 1/4" high and 22" deep. Lengths as shown. Designed to fit under 30" high Workstation frames.

Drawer cabinets are furnished with counter weights, multiple-drawer cabinets include an interlock system. All mobile cabinets have finished backs.

**Note:** Last four digits of part number designates door and drawer head style, hinge, hardwood species and pull selection. See page 4 for full explanation. Actual product may vary slightly from catalog illustration due to style and hardware selection.





## Cupboard Cabinets



<b>M10C262215L0000</b>	15"L	<b>M10C262215-0000</b>	15"L
<b>M10C262218L0000</b>	18"L	<b>M10C262218-0000</b>	18"L
<b>M10C262221L0000</b>	21"L	<b>M10C262221-0000</b>	21"L
<b>M10C262224L0000</b>	24"L	<b>M10C262224-0000</b>	24"L

## Drawer/Cupboard Cabinets




## Drawer/Cupboard Cabinets with Pullboards





<b>M20C262215L</b>	15"L	<b>M20C262215-</b>	15"L	<b>M25C262215L</b>	15"L	<b>M25C262215-</b>	15"L
<b>M20C262218L</b>	18"L	<b>M20C262218-</b>	18"L	<b>M25C262218L</b>	18"L	<b>M25C262218-</b>	18"L
<b>M20C262221L</b>	21"L	<b>M20C262221-</b>	21"L	<b>M25C262221L</b>	21"L	<b>M25C262221-</b>	21"L
<b>M20C262224L</b>	24"L	<b>M20C262224-</b>	24"L	<b>M25C262224L</b>	24"L	<b>M25C262224-</b>	24"L

## Drawer Cabinets

## Drawer Cabinets with Pullboards




<b>M30C262215-</b>	15"L	<b>M35C262215-</b>	15"L	<b>M45C262215-</b>	15"L
<b>M30C262218-</b>	18"L	<b>M35C262218-</b>	18"L	<b>M45C262218-</b>	18"L
<b>M30C262221-</b>	21"L	<b>M35C262221-</b>	21"L	<b>M45C262221-</b>	21"L
<b>M30C262224-</b>	24"L	<b>M35C262224-</b>	24"L	<b>M45C262224-</b>	24"L



<b>M62C262215-</b>	15"L	<b>M66C262215-</b>	15"L
<b>M62C262218-</b>	18"L	<b>M66C262218-</b>	18"L
<b>M62C262221-</b>	21"L	<b>M66C262221-</b>	21"L
<b>M62C262224-</b>	24"L	<b>M66C262224-</b>	24"L

# Steel Standing Height Mobile Cabinets





## Specifications

Steel Standing Height Mobile Cabinets are 33" high and 22" deep. Lengths as shown. Designed to fit under 36" high Workstation frames.

Drawer cabinets are furnished with counter weights, multiple-drawer cabinets include an interlock system. All mobile cabinets have finished backs.

**Note:** Last four digits of part number designates door and drawer head style, hinge, hardwood species and pull selection. See page 4 for full explanation. Actual product may vary slightly from catalog illustration due to style and hardware selection.




### Cupboard Cabinets



<b>M10C332215L0000</b>	15"L	<b>M10C332215-0000</b>	15"L
<b>M10C332218L0000</b>	18"L	<b>M10C332218-0000</b>	18"L
<b>M10C332221L0000</b>	21"L	<b>M10C332221-0000</b>	21"L
<b>M10C332224L0000</b>	24"L	<b>M10C332224-0000</b>	24"L

### Drawer/Cupboard Cabinets




### Drawer/Cupboard Cabinets with Pullboards





<b>M20C332215L</b>	15"L	<b>M20C332215-</b>	15"L	<b>M25C332215L</b>	15"L	<b>M25C332215-</b>	15"L
<b>M20C332218L</b>	18"L	<b>M20C332218-</b>	18"L	<b>M25C332218L</b>	18"L	<b>M25C332218-</b>	18"L
<b>M20C332221L</b>	21"L	<b>M20C332221-</b>	21"L	<b>M25C332221L</b>	21"L	<b>M25C332221-</b>	21"L
<b>M20C332224L</b>	24"L	<b>M20C332224-</b>	24"L	<b>M25C332224L</b>	24"L	<b>M25C332224-</b>	24"L

### Drawer Cabinets

### Drawer Cabinets with Pullboards



<b>M40C332215-</b>	15"L	<b>M46C332215-</b>	15"L	<b>M56C332215-</b>	15"L
<b>M40C332218-</b>	18"L	<b>M46C332218-</b>	18"L	<b>M56C332218-</b>	18"L
<b>M40C332221-</b>	21"L	<b>M46C332221-</b>	21"L	<b>M56C332221-</b>	21"L
<b>M40C332224-</b>	24"L	<b>M46C332224-</b>	24"L	<b>M56C332224-</b>	24"L



<b>M71C332215-</b>	15"L	<b>M68C332215-</b>	15"L
<b>M71C332218-</b>	18"L	<b>M68C332218-</b>	18"L
<b>M71C332221-</b>	21"L	<b>M68C332221-</b>	21"L
<b>M71C332224-</b>	24"L	<b>M68C332224-</b>	24"L


# Mobile Cabinet Tops & Accessories

## Specifications:

**1" Thick Mobile Cabinet Tops** are available in three depths, 22", 23", or 24" and three materials: Kemresin, Phenolic Resin or Plastic Laminate. They may be used on any Mobile


Cabinet and add 1" to the overall height of the cabinets. All edges of top are finished.

### 1" Thick — Mobile Cabinet Top (adds 1" to height of Mobile Cabinet)




#### 22" Deep flush to front of case

Length	Kemresin®	Kemresin Colors	Black Phenolic Resin	Black Plastic Laminate
15"	<b>F94R2215-00-</b>	<b>BK = Black</b>	<b>F94T2215-00</b>	<b>F94P2215-00</b>
18"	<b>F94R2218-00-</b>	<b>GR = Grey</b>	<b>F94T2218-00</b>	<b>F94P2218-00</b>
21"	<b>F94R2221-00-</b>	<b>PT = Putty</b>	<b>F94T2221-00</b>	<b>F94P2221-00</b>
24"	<b>F94R2224-00-</b>	<b>SL = Slate</b>	<b>F94T2224-00</b>	<b>F94P2224-00</b>
30"	<b>F94R2230-00-</b>		<b>F94T2230-00</b>	<b>F94P2230-00</b>
36"	<b>F94R2236-00-</b>		<b>F94T2236-00</b>	<b>F94P2236-00</b>
42"	<b>F94R2242-00-</b>		<b>F94T2242-00</b>	<b>F94P2242-00</b>
48"	<b>F94R2248-00-</b>		<b>F94T2248-00</b>	<b>F94P2248-00</b>




Inset Cabinet  
End View




#### 23" Deep flush to overlay door/drawer fronts - 1" overhang of case

Length	Kemresin®	Kemresin Colors	Black Phenolic Resin	Black Plastic Laminate
15"	<b>F94R2315-00-</b>	<b>BK = Black</b>	<b>F94T2315-00</b>	<b>F94P2315-00</b>
18"	<b>F94R2318-00-</b>	<b>GR = Grey</b>	<b>F94T2318-00</b>	<b>F94P2318-00</b>
21"	<b>F94R2321-00-</b>	<b>PT = Putty</b>	<b>F94T2321-00</b>	<b>F94P2321-00</b>
24"	<b>F94R2324-00-</b>	<b>SL = Slate</b>	<b>F94T2324-00</b>	<b>F94P2324-00</b>
30"	<b>F94R2330-00-</b>		<b>F94T2330-00</b>	<b>F94P2330-00</b>
36"	<b>F94R2336-00-</b>		<b>F94T2336-00</b>	<b>F94P2336-00</b>
42"	<b>F94R2342-00-</b>		<b>F94T2342-00</b>	<b>F94P2342-00</b>
48"	<b>F94R2348-00-</b>		<b>F94T2348-00</b>	<b>F94P2348-00</b>




Inset Cabinet  
End View



#### 24" Deep 1" overhang of overlay door/drawer fronts

Length	Kemresin®	Kemresin Colors	Black Phenolic Resin	Black Plastic Laminate
15"	<b>F94R2415-00-</b>	<b>BK = Black</b>	<b>F94T2415-00</b>	<b>F94P2415-00</b>
18"	<b>F94R2418-00-</b>	<b>GR = Grey</b>	<b>F94T2418-00</b>	<b>F94P2418-00</b>
21"	<b>F94R2421-00-</b>	<b>PT = Putty</b>	<b>F94T2421-00</b>	<b>F94P2421-00</b>
24"	<b>F94R2424-00-</b>	<b>SL = Slate</b>	<b>F94T2424-00</b>	<b>F94P2424-00</b>
30"	<b>F94R2430-00-</b>		<b>F94T2430-00</b>	<b>F94P2430-00</b>
36"	<b>F94R2436-00-</b>		<b>F94T2436-00</b>	<b>F94P2436-00</b>
42"	<b>F94R2442-00-</b>		<b>F94T2442-00</b>	<b>F94P2442-00</b>
48"	<b>F94R2448-00-</b>		<b>F94T2448-00</b>	<b>F94P2448-00</b>



Overlay Cabinet  
End View

# Mobile Cabinet Tops & Accessories



## Specifications:

**Steel Mobile Cabinet Tops** are fabricated from 18 gauge steel and finished to match the cabinet. They are designed for Steel Mobile Cabinets only.

**Mobile Cabinet Grabbar** allows easier movement and relocation of mobile cabinets. Fabricated from powder coated steel or stainless steel, they may be used on any


Mobile Cabinet, with or without Tops. Bars must be attached in field with supplied screws and/or bolts.

### 22" Deep — Steel — Mobile Cabinet Top (adds less than 1/4" to height of Mobile Cabinet)




*Inset Cabinet*  
End View



22" Deep Top for Steel Mobile Cabinets only		
Lg	Painted Steel	Stainless Steel
15"	<b>F94M2215-00</b>	<b>F94S2215-00</b>
18"	<b>F94M2218-00</b>	<b>F94S2218-00</b>
21"	<b>F94M2221-00</b>	<b>F94S2221-00</b>
24"	<b>F94M2224-00</b>	<b>F94S2224-00</b>
30"	<b>F94M2230-00</b>	<b>F94S2230-00</b>
36"	<b>F94M2236-00</b>	<b>F94S2236-00</b>
42"	<b>F94M2242-00</b>	<b>F94S2242-00</b>
48"	<b>F94M2248-00</b>	<b>F94S2248-00</b>



*Overlay Cabinet*  
End View




### 23" Deep — Steel — Mobile Cabinet Top (adds 1/2" to height of Mobile Cabinet)




*Inset Cabinet*  
End View


23" Deep Top for Steel Mobile Cabinets only		
Lg	Painted Steel	Stainless Steel
15"	<b>F94M2315-00</b>	<b>F94S2315-00</b>
18"	<b>F94M2318-00</b>	<b>F94S2318-00</b>
21"	<b>F94M2321-00</b>	<b>F94S2321-00</b>
24"	<b>F94M2324-00</b>	<b>F94S2324-00</b>
30"	<b>F94M2330-00</b>	<b>F94S2330-00</b>
36"	<b>F94M2336-00</b>	<b>F94S2336-00</b>
42"	<b>F94M2342-00</b>	<b>F94S2342-00</b>
48"	<b>F94M2348-00</b>	<b>F94S2348-00</b>




*Overlay Cabinet*  
End View



### Mobile Cabinet Grabbar




Length	Painted Steel	Stainless Steel
15"	<b>F73M0015-00</b>	<b>F73S0015-00</b>
18"	<b>F73M0018-00</b>	<b>F73S0018-00</b>
21"	<b>F73M0021-00</b>	<b>F73S0021-00</b>
24"	<b>F73M0024-00</b>	<b>F73S0024-00</b>



# Steel Apron Assemblies

## Components of an Apron Assembly




## Specifications:

Box Aprons are 4" high in lengths and depths as shown. Dimensions shown are sizes of the apron without legs. Each leg or pedestal adds 2" to the overall length of the apron assembly. A free-standing table with legs on both ends will therefore be 4" longer than the apron length.


All aprons can be used between cabinets, with leg assemblies or with pedestal legs.

For aprons 42" long and longer, spreaders are available to connect the two leg assemblies. When a spreader is used, order the A08 or A28 leg assembly which are designed to accept the spreader.

## 22" Deep Plain Rail Box Apron



**Inset Style Apron**



**Full Overlay Style Apron**


All Plain Rail Box Aprons available in either Inset or Full Overlay Style

22" Deep Apron		22" Deep Apron		22" Deep Apron	
Part Number	Apron Length	Part Number	Apron Length	Part Number	Apron Length
A00C042215-0	15"	A00C042236-0	36"	A00C042257-0	57"
A00C042216-0	16"	A00C042237-0	37"	A00C042258-0	58"
A00C042217-0	17"	A00C042238-0	38"	A00C042259-0	59"
A00C042218-0	18"	A00C042239-0	39"	A00C042260-0	60"
A00C042219-0	19"	A00C042240-0	40"	A00C042261-0	61"
A00C042220-0	20"	A00C042241-0	41"	A00C042262-0	62"
A00C042221-0	21"	A00C042242-0	42"	A00C042263-0	63"
A00C042222-0	22"	A00C042243-0	43"	A00C042264-0	64"
A00C042223-0	23"	A00C042244-0	44"	A00C042265-0	65"
A00C042224-0	24"	A00C042245-0	45"	A00C042266-0	66"
A00C042225-0	25"	A00C042246-0	46"	A00C042267-0	67"
A00C042226-0	26"	A00C042247-0	47"	A00C042268-0	68"
A00C042227-0	27"	A00C042248-0	48"	A00C042219-0	69"
A00C042228-0	28"	A00C042249-0	49"	A00C042270-0	70"
A00C042229-0	29"	A00C042250-0	50"	A00C042271-0	71"
A00C042230-0	30"	A00C042251-0	51"	A00C042272-0	72"
A00C042231-0	31"	A00C042252-0	52"		
A00C042232-0	32"	A00C042253-0	53"		
A00C042233-0	33"	A00C042254-0	54"		
A00C042234-0	34"	A00C042255-0	55"		
A00C042235-0	35"	A00C042256-0	56"		



# Steel Apron Assemblies

## Components of an Apron Assembly




## Specifications:

Box Aprons are 4" high in lengths and depths as shown. Dimensions shown are sizes of the apron without legs. Each leg or pedestal adds 2" to the overall length of the apron assembly. A free-standing table with legs on both ends will therefore be 4" longer than the apron length.


All aprons can be used between cabinets, with leg assemblies or with pedestal legs.

For aprons 42" long and longer, spreaders are available to connect the two leg assemblies. When a spreader is used, order the A08 or A28 leg assembly which are designed to accept the spreader.

## 28" Deep Plain Rail Box Apron



Inset Style Apron




Full Overlay Style Apron

All Plain Rail Box Aprons available in either Inset or Full Overlay Style

28" Deep Apron		28" Deep Apron		28" Deep Apron	
	Apron Length		Apron Length		Apron Length
A00C042815-___0	15"	A00C042836-___0	36"	A00C042857-___0	57"
A00C042816-___0	16"	A00C042837-___0	37"	A00C042858-___0	58"
A00C042817-___0	17"	A00C042838-___0	38"	A00C042859-___0	59"
A00C042818-___0	18"	A00C042839-___0	39"	A00C042860-___0	60"
A00C042819-___0	19"	A00C042840-___0	40"	A00C042861-___0	61"
A00C042820-___0	20"	A00C042841-___0	41"	A00C042862-___0	62"
A00C042821-___0	21"	A00C042842-___0	42"	A00C042863-___0	63"
A00C042822-___0	22"	A00C042843-___0	43"	A00C042864-___0	64"
A00C042823-___0	23"	A00C042844-___0	44"	A00C042865-___0	65"
A00C042824-___0	24"	A00C042845-___0	45"	A00C042866-___0	66"
A00C042825-___0	25"	A00C042846-___0	46"	A00C042867-___0	67"
A00C042826-___0	26"	A00C042847-___0	47"	A00C042868-___0	68"
A00C042827-___0	27"	A00C042848-___0	48"	A00C042819-___0	69"
A00C042828-___0	28"	A00C042849-___0	49"	A00C042870-___0	70"
A00C042829-___0	29"	A00C042850-___0	50"	A00C042871-___0	71"
A00C042830-___0	30"	A00C042851-___0	51"	A00C042872-___0	72"
A00C042831-___0	31"	A00C042852-___0	52"		
A00C042832-___0	32"	A00C042853-___0	53"		
A00C042833-___0	33"	A00C042854-___0	54"		
A00C042834-___0	34"	A00C042855-___0	55"		
A00C042835-___0	35"	A00C042856-___0	56"		

# Steel Apron Assemblies

## Components of an Apron Assembly




## Specifications:

Box Aprons are 4" high in lengths and depths as shown. Dimensions shown are sizes of the apron without legs. Each leg or pedestal adds 2" to the overall length of the apron assembly. A free-standing table with legs on both ends will therefore be 4" longer than the apron length.


All aprons can be used between cabinets, with leg assemblies or with pedestal legs.

For aprons 42" long and longer, spreaders are available to connect the two leg assemblies. When a spreader is used, order the A08 or A28 leg assembly which are designed to accept the spreader.

## 22" Deep Box Apron with Drawers



1 Drawer




22" Deep Apron

2 Drawers



22" Deep Apron

3 Drawers



22" Deep Apron


	Apron Length	Dwr Size
A51C042220-1	20"	18"
A51C042224-1	24"	21"
A51C042226-1	26"	24"
A51C042230-1	30"	27"
A51C042232-1	32"	30"
A51C042236-1	36"	33"

	Apron Length	Dwr Size
A52C042226-1	26"	12"
A52C042230-1	30"	12"
A52C042232-1	32"	15"
A52C042236-1	36"	15"
A52C042238-1	38"	18"
A52C042242-1	42"	18"
A52C042244-1	44"	21"
A52C042248-1	48"	21"
A52C042250-1	50"	24"
A52C042254-1	54"	24"
A52C042256-1	56"	27"
A52C042260-1	60"	27"
A52C042262-1	62"	30"
A52C042266-1	66"	30"
A52C042268-1	68"	33"
A52C042272-1	72"	33"

	Apron Length	Dwr Size
A53C042250-1	50"	15"
A53C042254-1	54"	15"
A53C042256-1	56"	17"
A53C042260-1	60"	17"
A53C042262-1	62"	19"
A53C042266-1	66"	19"
A53C042268-1	68"	21"
A53C042272-1	72"	21"

# Steel Apron Assemblies

## Components of an Apron Assembly




## Specifications:

Box Aprons are 4" high in lengths and depths as shown. Dimensions shown are sizes of the apron without legs. Each leg or pedestal adds 2" to the overall length of the apron assembly. A free-standing table with legs on both ends will therefore be 4" longer than the apron length.


All aprons can be used between cabinets, with leg assemblies or with pedestal legs.

For aprons 42" long and longer, spreaders are available to connect the two leg assemblies. When a spreader is used, order the A08 or A28 leg assembly which are designed to accept the spreader.

## 28" Deep Box Apron with Drawers



1 Drawer




28" Deep Apron

2 Drawers



28" Deep Apron

3 Drawers



28" Deep Apron

	Apron Length	Dwr Size
A51C042820-1	20"	18"
A51C042824-1	24"	21"
A51C042826-1	26"	24"
A51C042830-1	30"	27"
A51C042832-1	32"	30"
A51C042836-1	36"	33"

	Apron Length	Dwr Size
A52C042826-1	26"	12"
A52C042830-1	30"	12"
A52C042832-1	32"	15"
A52C042836-1	36"	15"
A52C042838-1	38"	18"
A52C042842-1	42"	18"
A52C042844-1	44"	21"
A52C042848-1	48"	21"
A52C042850-1	50"	24"
A52C042854-1	54"	24"
A52C042856-1	56"	27"
A52C042860-1	60"	27"
A52C042862-1	62"	30"
A52C042866-1	66"	30"
A52C042868-1	68"	33"
A52C042872-1	72"	33"

	Apron Length	Dwr Size
A53C042850-1	50"	15"
A53C042854-1	54"	15"
A53C042856-1	56"	17"
A53C042860-1	60"	17"
A53C042862-1	62"	19"
A53C042866-1	66"	19"
A53C042868-1	68"	21"
A53C042872-1	72"	21"

# Steel Apron Assemblies


## Specifications:

Corner Box Aprons provide support for corner work tops as well as the extra storage space of an apron drawer. The blind end may either be attached to the wall or supported by a leg or pedestal assembly.


Corner Box Aprons are furnished with steel panels when specified with wood drawer heads.

2" High Box Aprons are designed for use when knee space height is critical. They are available 22" deep and in lengths from 15" to 72" on 1" increments. 2" High Box Aprons are not designed for use with legs.

## Corner Box Apron



*Furnished with steel panels when specified with wood drawer heads.*



**1 Drawer on Right**



**22" Deep**

**1 Drawer on Left**



**22" Deep**

	Apron Length	Dwr Size
<b>A05C042248-1</b> _____	48"	21"
<b>A05C042254-1</b> _____	54"	21"
<b>A05C042260-1</b> _____	60"	21"
<b>A05C042266-1</b> _____	66"	21"
<b>A05C042272-1</b> _____	72"	21"

**28" Deep**

	Apron Length	Dwr Size
<b>A05C042848-1</b> _____	48"	21"
<b>A05C042854-1</b> _____	54"	21"
<b>A05C042860-1</b> _____	60"	21"
<b>A05C042866-1</b> _____	66"	21"
<b>A05C042872-1</b> _____	72"	21"

	Apron Length	Dwr Size
<b>A04C042248-1</b> _____	48"	21"
<b>A04C042254-1</b> _____	54"	21"
<b>A04C042260-1</b> _____	60"	21"
<b>A04C042266-1</b> _____	66"	21"
<b>A04C042272-1</b> _____	72"	21"

**28" Deep**

	Apron Length	Dwr Size
<b>A04C042848-1</b> _____	48"	21"
<b>A04C042854-1</b> _____	54"	21"
<b>A04C042860-1</b> _____	60"	21"
<b>A04C042866-1</b> _____	66"	21"
<b>A04C042872-1</b> _____	72"	21"

# Steel Apron Assemblies

## Specifications:


Apron Rails may be used in place of Box Aprons to provide worksurface support between two units. They are available 2" or 4" high and from 15" to 72" long in 1" increments. They may be installed vertically to provide a 2"

or 4" high rail, or turned 90 degrees to provide a 1 1/2" high rail. Apron rails should be ordered in pairs to provide both front and rear top support.

Work Top Support Rails and Assemblies may be used anywhere

maximum knee space is wanted but top support is required. Fabricated of heavy gauge steel channels, it is only 7/8" high. The Keyboard Tray may be added to any of the frames or may be used independently.

### Inset 4" High Apron Rail




**A06C0402**□□-00\_0 4" High  
Inset Apron Rails are available from 15" to 72" long in 1" increments. To specify apron rail, insert the length required into the part number.

**Example:** A06C040248-0000 is a 48" long cold rolled steel apron rail.

Apron rails should be ordered in pairs to provide both front and rear worksurface support.

### Overlay 4" High Apron Rail




**A06C0403**□□-1\_\_0 4" High  
Overlay Apron Rails are available from 15" to 72" long in 1" increments. To specify apron rail, insert the length required into the part number.

**Example:** A06C040248-1610 is a 48" long apron rail with Oak overlay panel.

Apron rails should be ordered in pairs to provide both front and rear worksurface support.

### 4" High Adjustable Rail




**A56C040236-00\_0** 36"  
**A56C040248-00\_0** 48"  
**A56C040260-00\_0** 60"  
**A56C040272-00\_0** 72"

Same as A06 Apron Rails except furnished with a loose angle at one end to allow the rail to be cut in the field to fit available space.

Apron rails should be ordered in pairs to provide both front and rear worksurface support.

### 2" High Plain Rail Box Apron




**A00C0222**□□-00\_0  
2" high Box Aprons are available from 15" to 72" long in 1" increments. To specify apron, insert the length required into the part number.

**Example:** A00C020248-0000 is a 48" long 2" high cold rolled steel box apron.

**Not for use with legs**

### 2" High Apron Rail




**A06C0202**□□-00\_0 2" High  
Apron rails are available from 15" to 72" long in 1" increments. To specify apron rail, insert the length required into the part number.

**Example:** A06C020248-0000 is a 48" long cold rolled steel apron rail.

Apron rails should be ordered in pairs to provide both front and rear worksurface support.

### Support Rails



**A81C0102**□□-0000  
Support rails may be used individually or assembled to form a support framework. Rails are available from 12" to 72" long in 1" increments and are furnished with angles attached to each end. To specify rail, insert the length required in the part number.

**Example:** A81010248-0000 is a 48" long support rail.

Support rails should be ordered in pairs to provide both front and rear worksurface support.

# Steel Leg Assemblies & Pedestals

## Specifications:

Legs are 2" square and all except the Top Mounted Leg are joined in pairs with a 2" high stretcher and a 4" high apron rail. Dimensions are as shown.

For aprons 42" long and longer, spreaders are available to connect the two leg assemblies. When a

spreader is used, order the A08 or A28 leg assembly which are designed for use with a spreader.





All legs may be ordered with 4" swivel locking casters by adding **-B10** to the end of the part numbers.

Example: A07C292202-0000-**B10**

Adjustable Height Legs are adjustable in 1" increments from 29" to 35" high.

### End Leg Assemblies

### Intermediate Leg Assemblies



For use with Spreader

For use with Spreader

**22<sup>3</sup>/<sub>16</sub>" Deep**  
**A07C292202-00\_0** 29"H  
**A07C322202-00\_0** 32"H  
**A07M352202-00\_0** 35"H

**22<sup>3</sup>/<sub>16</sub>" Deep**  
**A08C292202-00\_0** 29"H  
**A08C322202-00\_0** 32"H  
**A08C352202-00\_0** 35"H

**22<sup>3</sup>/<sub>16</sub>" Deep**  
**A27C292202-00\_0** 29"H  
**A27C322202-00\_0** 32"H  
**A27C352202-00\_0** 35"H

**22<sup>3</sup>/<sub>16</sub>" Deep**  
**A28C292202-00\_0** 29"H  
**A28C322202-00\_0** 32"H  
**A28C352202-00\_0** 35"H

**28<sup>3</sup>/<sub>16</sub>" Deep**  
**A07C292802-00\_0** 29"H  
**A07C322802-00\_0** 32"H  
**A07C352802-00\_0** 35"H

**28<sup>3</sup>/<sub>16</sub>" Deep**  
**A08C292802-00\_0** 29"H  
**A08C322802-00\_0** 32"H  
**A08C352802-00\_0** 35"H





**28<sup>3</sup>/<sub>16</sub>" Deep**  
**A27C292802-00\_0** 29"H  
**A27C322802-00\_0** 32"H  
**A27C352802-00\_0** 35"H

**28<sup>3</sup>/<sub>16</sub>" Deep**  
**A28C292802-00\_0** 29"H  
**A28C322802-00\_0** 32"H  
**A28C352802-00\_0** 35"H

Each Leg Assembly adds 2" to the overall length of the apron assembly to which it is attached..

### Adjustable Height End Leg Assemblies

### Adjustable Height Intermediate Leg Ass'y



For use with Spreader

For use with Spreader

**22<sup>3</sup>/<sub>16</sub>" Deep**  
**A57C352202-00\_0** Adj.  
 Adjustable Height 29"-35"

**22<sup>3</sup>/<sub>16</sub>" Deep**  
**A58C352202-00\_0** Adj.  
 Adjustable Height 29"-35"

**22<sup>3</sup>/<sub>16</sub>" Deep**  
**A77C352202-00\_0** Adj.  
 Adjustable Height 29"-35"

**22<sup>3</sup>/<sub>16</sub>" Deep**  
**A78C352202-00\_0** Adj.  
 Adjustable Height 29"-35"

**28<sup>3</sup>/<sub>16</sub>" Deep**  
**A57C352802-00\_0** Adj.  
 Adjustable Height 29"-35"

**28<sup>3</sup>/<sub>16</sub>" Deep**  
**A58C352802-00\_0** Adj.  
 Adjustable Height 29"-35"

**28<sup>3</sup>/<sub>16</sub>" Deep**  
**A77C352802-00\_0** Adj.  
 Adjustable Height 29"-35"

**28<sup>3</sup>/<sub>16</sub>" Deep**  
**A78C352802-00\_0** Adj.  
 Adjustable Height 29"-35"

# Steel Apron Rails & Panels

## Specifications:


Knee Space panels are designed to close off the space between the floor and the bottom of the apron or countertop. They are sized for 29" high, 32" high, and 35" high cabinets. Each Knee Space Panel is made with a separate base rail, to allow the removal

of the panel without removing the base molding.

Knee Space Panels are furnished with two angles to attach to panel to the cabinets.

Pedestal Legs are 2" wide, depths and heights as shown. They are designed to be used with or without Box Aprons to form knee spaces at the end of assemblies.

## Knee Space Panels



Knee space panels are available from 15" to 72" long in 1" increments. To specify knee space panel, insert the length required into the part number.

**Example:** A10C310048-0000 is a 48" long knee space panel.

### Sitting Height Knee Space Panels

Apron Part Number	Height	Length	for use with:
A10C2500 □□-00_0	25"	15"-72"L	4" high apron
A10C2800 □□-00_0	28"	15"-72"L	A80 support frame
A10C2900 □□-00_0	29"	15"-72"L	no apron


### ADA Height Knee Space Panels

Apron Part Number	Height	Length	for use with:
A10C2800 □□-00_0	28"	15"-72"L	4" high apron
A10C3100 □□-00_0	31"	15"-72"L	A80 support frame
A10C3200 □□-00_0	32"	15"-72"L	no apron or frame

### Standing Height Knee Space Panels


Apron Part Number	Height	Length	for use with:
A10C3100 □□-00_0	31"	15"-72"L	4" high apron
A10C3400 □□-00_0	34"	15"-72"L	A80 support frame
A10C3500 □□-00_0	35"	15"-72"L	no apron or frame

## Spreader



A09C020230-00_0	30"L
A09C020236-00_0	36"L
A09C020242-00_0	42"L
A09C020248-00_0	48"L
A09C020254-00_0	54"L
A09C020260-00_0	60"L
A09C020266-00_0	66"L
A09C020272-00_0	72"L


## Top Mounted Leg



A15C290202-0000	29"H
A15C320202-00_0	32"H
A15C350202-00_0	35"H
A65C350202-00_0	Adj. Adjustable Height 29"-35"

2" square leg welded to 4" square steel plate

## Wall Pedestal Leg



### Left Side 22" Deep

A11C292202-00_0	29"H
A11C322202-00_0	32"H
A11C352202-00_0	35"H

### Left Side 28" Deep

A11C292802-00_0	29"H
A11C322802-00_0	32"H
A11C352802-00_0	35"H


### Right Side 22" Deep

A12C292202-00_0	29"H
A12C322202-00_0	32"H
A12C352202-00_0	35"H

### Right Side 28" Deep

A12C292802-00_0	29"H
A12C322802-00_0	32"H
A12C352802-00_0	35"H

## Center Pedestal Leg



### 52" Deep


A13C295202-00_0	29"H
A13C325202-00_0	32"H
A13C355202-00_0	35"H

Each Pedestal Leg adds 2" to the overall length of the assembly to which it is attached.

Pedestal Legs may be ordered with the topspace 7/8" deeper so that the front edge lines up with the doors and drawer fronts on adjacent Full Overlay cabinets. To specify the deeper toe space, add a **C** to the end of the part number.

Example: **A11C352202-0000-C**

# Steel Wall Mounted Storage Cabinets



## Specifications:

Open Wall Mounted Storage Cabinets are available in four heights, 24", 30", 36" and 48" and two depths; 12<sup>3</sup>/<sub>4</sub>" and 16". Lengths are as shown. All are furnished with one, two, or three adjustable shelves as shown.


Wall cabinets are designed to be mounted on hangers rails. Shown on page 47 and 81.

## Open Storage

Height	Depth	Product Code	Length
24" High	12 <sup>3</sup> / <sub>4</sub> " Deep	W00C241315-00_0	15"L
		W00C241318-00_0	18"L
		W00C241321-00_0	21"L
		W00C241324-00_0	24"L
	16" Deep	W00C241615-00_0	15"L
		W00C241618-00_0	18"L
		W00C241621-00_0	21"L
		W00C241624-00_0	24"L
30" High	12 <sup>3</sup> / <sub>4</sub> " Deep	W00C301315-00_0	15"L
		W00C301318-00_0	18"L
		W00C301321-00_0	21"L
		W00C301324-00_0	24"L
	16" Deep	W00C301615-00_0	15"L
		W00C301618-00_0	18"L
		W00C301621-00_0	21"L
		W00C301624-00_0	24"L
36" High	12 <sup>3</sup> / <sub>4</sub> " Deep	W00C361315-00_0	15"L
		W00C361318-00_0	18"L
		W00C361321-00_0	21"L
		W00C361324-00_0	24"L
	16" Deep	W00C361615-00_0	15"L
		W00C361618-00_0	18"L
		W00C361621-00_0	21"L
		W00C361624-00_0	24"L
48" High	12 <sup>3</sup> / <sub>4</sub> " Deep	W00C481315-00_0	15"L
		W00C481318-00_0	18"L
		W00C481321-00_0	21"L
		W00C481324-00_0	24"L
	16" Deep	W00C481615-00_0	15"L
		W00C481618-00_0	18"L
		W00C481621-00_0	21"L
		W00C481624-00_0	24"L
24" High	12 <sup>3</sup> / <sub>4</sub> " Deep	W00C241330-00_0	30"L
		W00C241336-00_0	36"L
		W00C241342-00_0	42"L
		W00C241348-00_0	48"L
	16" Deep	W00C241630-00_0	30"L
		W00C241636-00_0	36"L
		W00C241642-00_0	42"L
		W00C241648-00_0	48"L
30" High	12 <sup>3</sup> / <sub>4</sub> " Deep	W00C301330-00_0	30"L
		W00C301336-00_0	36"L
		W00C301342-00_0	42"L
		W00C301348-00_0	48"L
	16" Deep	W00C301630-00_0	30"L
		W00C301636-00_0	36"L
		W00C301642-00_0	42"L
		W00C301648-00_0	48"L
36" High	12 <sup>3</sup> / <sub>4</sub> " Deep	W00C361330-00_0	30"L
		W00C361336-00_0	36"L
		W00C361342-00_0	42"L
		W00C361348-00_0	48"L
	16" Deep	W00C361630-00_0	30"L
		W00C361636-00_0	36"L
		W00C361642-00_0	42"L
		W00C361648-00_0	48"L
48" High	12 <sup>3</sup> / <sub>4</sub> " Deep	W00C481330-00_0	30"L
		W00C481336-00_0	36"L
		W00C481342-00_0	42"L
		W00C481348-00_0	48"L
	16" Deep	W00C481630-00_0	30"L
		W00C481636-00_0	36"L
		W00C481642-00_0	42"L
		W00C481648-00_0	48"L



# Steel Wall Mounted Storage Cabinets



### Specifications:

Plate Glass Door Wall Mounted Storage Cabinets are available in three heights, 24", 30", 36" and 48" and two depths; 12<sup>3</sup>/<sub>4</sub>" and 16". Lengths are as shown. Cabinets are furnished with one, two, or three adjustable shelves as shown.


Doors are 1/4" float glass and are mounted on friction-free roller glides.

Wall cabinets are designed to be mounted on hangers rails. Shown on page 47 and 81.

### Sliding Plate Glass Doors — Clear Glass

### Wall Mounted Cabinet Hanging Rails

#### 24" High




#### 12<sup>3</sup>/<sub>4</sub>" Deep

W10C241330-00_0	30"L
W10C241336-00_0	36"L
W10C241342-00_0	42"L
W10C241348-00_0	48"L

#### 16" Deep

W10C241630-00_0	30"L
W10C241636-00_0	36"L
W10C241642-00_0	42"L
W10C241648-00_0	48"L

#### 30" High




#### 12<sup>3</sup>/<sub>4</sub>" Deep

W10C301330-00_0	30"L
W10C301336-00_0	36"L
W10C301342-00_0	42"L
W10C301348-00_0	48"L

#### 16" Deep

W10C301630-00_0	30"L
W10C301636-00_0	36"L
W10C301642-00_0	42"L
W10C301648-00_0	48"L

#### 36" High




#### 12<sup>3</sup>/<sub>4</sub>" Deep

W10C361330-00_0	30"L
W10C361336-00_0	36"L
W10C361342-00_0	42"L
W10C361348-00_0	48"L

#### 16" Deep

W10C361630-00_0	30"L
W10C361636-00_0	36"L
W10C361642-00_0	42"L
W10C361648-00_0	48"L

#### 48" High




#### 12<sup>3</sup>/<sub>4</sub>" Deep

W10C481330-00_0	30"L
W10C481336-00_0	36"L
W10C481342-00_0	42"L
W10C481348-00_0	48"L

#### 16" Deep


W10C481630-00_0	30"L
W10C481636-00_0	36"L
W10C481642-00_0	42"L
W10C481648-00_0	48"L



- L01-020109-00 for 15" wall cabinet
- L01-020112-00 for 18" wall cabinet
- L01-020115-00 for 21" wall cabinet
- L01-020118-00 for 24" wall cabinet
- L01-020124-00 for 30" wall cabinet
- L01-020130-00 for 36" wall cabinet
- L01-020136-00 for 42" wall cabinet
- L01-020142-00 for 48" wall cabinet

Formed of heavy gauge cold-rolled steel. Designed to capture lip on top rear of wall cabinets for easy installation.

# Steel Wall Mounted Storage Cabinets



## Specifications:

Sliding Door Wall Mounted Storage Cabinets are available in four heights, 24", 30", 36" and 48" and two depths; 12<sup>3</sup>/<sub>4</sub>" and 16". Lengths are as shown. Cabinets are glazed with 1/8" float glass unless otherwise specified and are furnished with one, two or three adjustable shelves as shown. Doors are

mounted on friction-free roller glides.


Sliding Door Wall Mounted Storage Cabinets are available in Door Style 1 and 6 only.

Wall cabinets are designed to be mounted on hangers rails. Shown on page 47 and 81.

### Sliding Glazed Doors

### Sliding Panel Doors

#### 24" High




#### 12<sup>3</sup>/<sub>4</sub>" Deep

W41C241330-0__0	30"
W41C241336-0__0	36"
W41C241342-0__0	42"
W41C241348-0__0	48"

#### 16" Deep

W41C241630-0__0	30"
W41C241636-0__0	36"
W41C241642-0__0	42"
W41C241648-0__0	48"

#### 24" High




#### 12<sup>3</sup>/<sub>4</sub>" Deep

W12C241330-0__0	30"
W12C241336-0__0	36"
W12C241342-0__0	42"
W12C241348-0__0	48"

#### 16" Deep

W12C241630-0__0	30"
W12C241636-0__0	36"
W12C241642-0__0	42"
W12C241648-0__0	48"

#### 30" High




#### 12<sup>3</sup>/<sub>4</sub>" Deep

W41C301330-0__0	30"
W41C301336-0__0	36"
W41C301342-0__0	42"
W41C301348-0__0	48"

#### 16" Deep

W41C301630-0__0	30"
W41C301636-0__0	36"
W41C301642-0__0	42"
W41C301648-0__0	48"

#### 30" High




#### 12<sup>3</sup>/<sub>4</sub>" Deep

W12C301330-0__0	30"
W12C301336-0__0	36"
W12C301342-0__0	42"
W12C301348-0__0	48"

#### 16" Deep

W12C301630-0__0	30"
W12C301636-0__0	36"
W12C301642-0__0	42"
W12C301648-0__0	48"

#### 36" High




#### 12<sup>3</sup>/<sub>4</sub>" Deep

W41C361330-0__0	30"
W41C361336-0__0	36"
W41C361342-0__0	42"
W41C361348-0__0	48"

#### 16" Deep

W41C361630-0__0	30"
W41C361636-0__0	36"
W41C361642-0__0	42"
W41C361648-0__0	48"

#### 36" High




#### 12<sup>3</sup>/<sub>4</sub>" Deep

W12C361330-0__0	30"
W12C361336-0__0	36"
W12C361342-0__0	42"
W12C361348-0__0	48"

#### 16" Deep

W12C361630-0__0	30"
W12C361636-0__0	36"
W12C361642-0__0	42"
W12C361648-0__0	48"

#### 48" High




#### 12<sup>3</sup>/<sub>4</sub>" Deep

W41C481330-0__0	30"
W41C481336-0__0	36"
W41C481342-0__0	42"
W41C481348-0__0	48"

#### 16" Deep

W41C481630-0__0	30"
W41C481636-0__0	36"
W41C481642-0__0	42"
W41C481648-0__0	48"

#### 48" High




#### 12<sup>3</sup>/<sub>4</sub>" Deep

W12C481330-0__0	30"
W12C481336-0__0	36"
W12C481342-0__0	42"
W12C481348-0__0	48"

#### 16" Deep

W12C481630-0__0	30"
W12C481636-0__0	36"
W12C481642-0__0	42"
W12C481648-0__0	48"

# Steel Wall Mounted Storage Cabinets



## Specifications:


Double Sided Plate Glass Pass-Thru Cabinets are available 24", 30", and 36" high, 12<sup>3</sup>/<sub>4</sub>" and 16" deep and in lengths shown. Double Sided Pass-Thru Cabinets have doors on both sides and are designed to be mounted on center bench reagent racks.

Doors are 1/4" float glass and are mounted on friction-free roller glides.

Wall cabinets are designed to be mounted on hangers rails. Shown on page 47 and 81.

## Double Sided Pass-Thru Cabinet - Clear Glass

### 24" High




#### 12<sup>3</sup>/<sub>4</sub>" Deep

W13C241330-00_0	30"L
W13C241336-00_0	36"L
W13C241342-00_0	42"L
W13C241348-00_0	48"L

#### 16" Deep

W13C241630-00_0	30"L
W13C241636-00_0	36"L
W13C241642-00_0	42"L
W13C241648-00_0	48"L

### 30" High




#### 12<sup>3</sup>/<sub>4</sub>" Deep

W13C301330-00_0	30"L
W13C301336-00_0	36"L
W13C301342-00_0	42"L
W13C301348-00_0	48"L

#### 16" Deep

W13C301630-00_0	30"L
W13C301636-00_0	36"L
W13C301642-00_0	42"L
W13C301648-00_0	48"L

### 36" High




#### 12<sup>3</sup>/<sub>4</sub>" Deep

W13C361330-00_0	30"L
W13C361336-00_0	36"L
W13C361342-00_0	42"L
W13C361348-00_0	48"L

#### 16" Deep

W13C361630-00_0	30"L
W13C361636-00_0	36"L
W13C361642-00_0	42"L
W13C361648-00_0	48"L

# Steel Wall Mounted Storage Cabinets





## Specifications:


Swinging Glazed Door Wall Mounted Storage Cabinets are available in four heights, 24", 30", 36" and 48" and two depths; 12<sup>3</sup>/<sub>4</sub>" and 16". Lengths are as shown. Cabinets are glazed with 1/8" float glass unless otherwise specified and are furnished with one, two or three adjustable shelves as shown.


**Note:** Last four digits of part number designate door style, material, and pull selection. See page 4 for full explanation. Beveled Edge, Door Style 7 is not available for wall cabinets. Wall cabinets are designed to be mounted on hanging rails. Shown on page 47 and 81.


## Swinging Glazed Doors


24" High	Left Hand		Right Hand		12 <sup>3</sup> / <sub>4</sub> " D
		W50C241315L	W50C241318L	W50C241315-	W50C241318-
	W50C241321L	W50C241324L	W50C241321-	W50C241324-	21"L
			W50C241321-	W50C241324-	24"L
					16" D
	W50C241615L	W50C241618L	W50C241615-	W50C241618-	15"L
	W50C241621L	W50C241624L	W50C241621-	W50C241624-	21"L
			W50C241621-	W50C241624-	24"L


24" High	12 <sup>3</sup> / <sub>4</sub> " Deep		16" Deep
		W50C241330-	
	W50C241342-	W50C241348-	42"L
			48"L
	W50C241630-	W50C241636-	30"L
	W50C241642-	W50C241648-	42"L
			48"L


30" High	Left Hand		Right Hand		12 <sup>3</sup> / <sub>4</sub> " D
		W50C301315L	W50C301318L	W50C301315-	W50C301318-
	W50C301321L	W50C301324L	W50C301321-	W50C301324-	21"L
			W50C301321-	W50C301324-	24"L
					16" D
	W50C301615L	W50C301618L	W50C301615-	W50C301618-	15"L
	W50C301621L	W50C301624L	W50C301621-	W50C301624-	21"L
			W50C301621-	W50C301624-	24"L

30" High	12 <sup>3</sup> / <sub>4</sub> " Deep		16" Deep
		W50C301330-	
	W50C301342-	W50C301348-	42"L
			48"L
	W50C301630-	W50C301636-	30"L
	W50C301642-	W50C301648-	42"L
			48"L


36" High	Left Hand		Right Hand		12 <sup>3</sup> / <sub>4</sub> " D
		W50C361315L	W50C361318L	W50C361315-	W50C361318-
	W50C361321L	W50C361324L	W50C361321-	W50C361324-	21"L
			W50C361321-	W50C361324-	24"L
					16" D
	W50C361615L	W50C361618L	W50C361615-	W50C361618-	15"L
	W50C361621L	W50C361624L	W50C361621-	W50C361624-	21"L
			W50C361621-	W50C361624-	24"L

36" High	12 <sup>3</sup> / <sub>4</sub> " Deep		16" Deep
		W50C361330-	
	W50C361342-	W50C361348-	42"L
			48"L
	W50C361630-	W50C361636-	30"L
	W50C361642-	W50C361648-	42"L
			48"L

48" High	Left Hand		Right Hand		12 <sup>3</sup> / <sub>4</sub> " D
		W50C481315L	W50C481318L	W50C481315-	W50C481318-
	W50C481321L	W50C481324L	W50C481321-	W50C481324-	21"L
			W50C481321-	W50C481324-	24"L
					16" D
	W50C481615L	W50C481618L	W50C481615-	W50C481618-	15"L
	W50C481621L	W50C481624L	W50C481621-	W50C481624-	21"L
			W50C481621-	W50C481624-	24"L

48" High	12 <sup>3</sup> / <sub>4</sub> " Deep		16" Deep
		W50C481330-	
	W50C481342-	W50C481348-	42"L
			48"L
	W50C481630-	W50C481636-	30"L
	W50C481642-	W50C481648-	42"L
			48"L

# Steel Wall Mounted Storage Cabinets





### Specifications:


Swinging Panel Door Wall Mounted Storage Cabinets are available in four heights, 24", 30", 36" and 48" and two depths; 12<sup>3</sup>/<sub>4</sub>" and 16". Lengths are as shown. Cabinets are furnished with one, two or three adjustable shelves as shown.


**Note:** Last four digits of part number designate door style, material, and pull selection. See page 4 for full explanation. Beveled Edge, Door Style 7 is not available for wall cabinets. Wall cabinets are designed to be mounted on hanging rails. Shown on page 47 and 81.


### Swinging Panel Doors


24" High	Left Hand		Right Hand		12 <sup>3</sup> / <sub>4</sub> " D
		W25C241315L	W25C241318L	W25C241321L	W25C241324L
					16" D
	W25C241615L	W25C241618L	W25C241621L	W25C241624L	15" 18" 21" 24"


24" High	12 <sup>3</sup> / <sub>4</sub> " Deep		16" Deep
		W25C241330-	
	W25C241630-	W25C241636-	30" 36" 42" 48"


30" High	Left Hand		Right Hand		12 <sup>3</sup> / <sub>4</sub> " D
		W25C301315L	W25C301318L	W25C301321L	W25C301324L
					16" D
	W25C301615L	W25C301618L	W25C301621L	W25C301624L	15" 18" 21" 24"

30" High	12 <sup>3</sup> / <sub>4</sub> " Deep		16" Deep
		W25C301330-	
	W25C301630-	W25C301636-	30" 36" 42" 48"


36" High	Left Hand		Right Hand		12 <sup>3</sup> / <sub>4</sub> " D
		W25C361315L	W25C361318L	W25C361321L	W25C361324L
					16" D
	W25C361615L	W25C361618L	W25C361621L	W25C361624L	15" 18" 21" 24"

36" High	12 <sup>3</sup> / <sub>4</sub> " Deep		16" Deep
		W25C361330-	
	W25C361630-	W25C361636-	30" 36" 42" 48"

48" High	Left Hand		Right Hand		12 <sup>3</sup> / <sub>4</sub> " D
		W25C481315L	W25C481318L	W25C481321L	W25C481324L
					16" D
	W25C481615L	W25C481618L	W25C481621L	W25C481624L	15" 18" 21" 24"

48" High	12 <sup>3</sup> / <sub>4</sub> " Deep		16" Deep
		W25C481330-	
	W25C481630-	W25C481636-	30" 36" 42" 48"

# Steel Wall Mounted Storage Cabinets




**Specifications:**

Corner Wall Mounted Storage Cabinets are available in four heights, 24", 30", 36" and 48" and two depths; 12<sup>3</sup>/<sub>4</sub>" and 16" and with either glazed or panel swinging doors. Cabinets are furnished with one, two or three adjustable shelves as shown.


**Note:** Last four digits of part number designate door style, material, and pull selection. See page 4 for full explanation. Beveled Edge, Door Style 7 is not available for wall cabinets. Wall cabinets are designed to be mounted on hanging rails. Shown on page 47 and 81.

**Swinging Door Corner Cabinets**

Glazed Door	Panel Door	Glazed Door	Panel Door
<b>24" High</b>	<b>12<sup>3</sup>/<sub>4</sub>" Deep</b>	<b>24" High</b>	<b>16" Deep</b>
W51C241320-_____	W26C241320-_____	W51C241615-_____	W26C241615-_____
<b>30" High</b>	<b>12<sup>3</sup>/<sub>4</sub>" Deep</b>	<b>30" High</b>	<b>16" Deep</b>
W51C301320-_____	W26C301320-_____	W51C301615-_____	W26C301615-_____
<b>36" High</b>	<b>12<sup>3</sup>/<sub>4</sub>" Deep</b>	<b>36" High</b>	<b>16" Deep</b>
W51C361320-_____	W26C361320-_____	W51C361615-_____	W26C361615-_____
<b>48" High</b>	<b>12<sup>3</sup>/<sub>4</sub>" Deep</b>	<b>48" High</b>	<b>16" Deep</b>
W51C481320-_____	W26C481320-_____	W51C481615-_____	W26C481615-_____




PLAN VIEW OF  
W51C241320  
W51C301320  
W51C361320  
W51C481320  
W26C241320  
W26C301320  
W26C361320  
W26C481320



PLAN VIEW OF  
W51C241615  
W51C301615  
W51C361615  
W51C481615  
W26C241615  
W26C301615  
W26C361615  
W26C481615

# Steel Full Height Storage Cabinets




**Specifications:**

Full Height Storage Cabinets are 84" high and 12<sup>3</sup>/<sub>4</sub>", 16", or 22" deep. Lengths as shown. They are furnished with five adjustable shelves.

**CAUTION:** Full height cabinets may present a potential tipping over problem. Proper installation requires these cabinets to be anchored to a wall or other cabinets to improve stability and ensure smooth door movement.

**NOTE:** Full height storage cabinets are not available in Overlay – Beveled Edge Door Style 7.

**Open Storage**



**12<sup>3</sup>/<sub>4</sub>" Deep**


S00C841318-00_0	18"L
S00C841324-00_0	24"L

**16" Deep**

S00C841618-00_0	18"L
S00C841624-00_0	24"L

**22" Deep**

S00C842218-00_0	18"L
S00C842224-00_0	24"L



**12<sup>3</sup>/<sub>4</sub>" Deep**

S00C841330-00_0	30"L
S00C841336-00_0	36"L
S00C841342-00_0	42"L
S00C841348-00_0	48"L


**16" Deep**

S00C841630-00_0	30"L
S00C841636-00_0	36"L
S00C841642-00_0	42"L
S00C841648-00_0	48"L

**22" Deep**

S00C842230-00_0	30"L
S00C842236-00_0	36"L
S00C842242-00_0	42"L
S00C842248-00_0	48"L

**Sliding Glazed Doors**



**12<sup>3</sup>/<sub>4</sub>" Deep**

S41C841330-0__0	30"L
S41C841336-0__0	36"L
S41C841342-0__0	42"L
S41C841348-0__0	48"L


**16" Deep**

S41C841630-0__0	30"L
S41C841636-0__0	36"L
S41C841642-0__0	42"L
S41C841648-0__0	48"L

**22" Deep**

S41C842230-0__0	30"L
S41C842236-0__0	36"L
S41C842242-0__0	42"L
S41C842248-0__0	48"L

**Sliding Panel Doors**



**12<sup>3</sup>/<sub>4</sub>" Deep**

S12C841330-0__0	30"L
S12C841336-0__0	36"L
S12C841342-0__0	42"L
S12C841348-0__0	48"L


**16" Deep**

S12C841630-0__0	30"L
S12C841636-0__0	36"L
S12C841642-0__0	42"L
S12C841648-0__0	48"L

**22" Deep**

S12C842230-0__0	30"L
S12C842236-0__0	36"L
S12C842242-0__0	42"L
S12C842248-0__0	48"L

# Steel Full Height Storage Cabinets




### Specifications:

Full Height Storage Cabinets are 84" high and 12<sup>3</sup>/<sub>4</sub>", 16", or 22" deep. Lengths as shown. They are furnished with five adjustable shelves.


**CAUTION:** Full height cabinets may present a potential tipping over problem. Proper installation requires these cabinets to be anchored to a wall or other cabinets to improve stability and ensure smooth door movement.

**NOTE:** Full height storage cabinets are not available in Overlay – Beveled Edge Door Style 7.

### Swinging Glazed Doors




Left Hand	Right Hand	12 <sup>3</sup> / <sub>4</sub> " D
S50C841318L	S50C841318-	18"L
S50C841324L	S50C841324-	24"L
<b>16" D</b>		
S50C841618L	S50C841618-	18"L
S50C841624L	S50C841624-	24"L
<b>22" D</b>		
S50C842218L	S50C8422-18	18"L
S50C842224L	S50C8422-24	24"L




12 <sup>3</sup> / <sub>4</sub> " Deep	
S50C841330-	30"L
S50C841336-	36"L
S50C841342-	42"L
S50C841348-	48"L
<b>16" Deep</b>	
S50C841630-	30"L
S50C841636-	36"L
S50C841642-	42"L
S50C841648-	48"L
<b>22" Deep</b>	
S50C842230-	30"L
S50C842236-	36"L
S50C842242-	42"L
S50C842248-	48"L

### Swinging Panel Doors




Left Hand	Right Hand	12 <sup>3</sup> / <sub>4</sub> " D
S25C841318L	S25C841318-	18"L
S25C841324L	S25C841324-	24"L
<b>16" D</b>		
S25C841618L	S25C841618-	18"L
S25C841624L	S25C841624-	24"L
<b>22" D</b>		
S25C842218L	S25C842218-	18"L
S25C842224L	S25C842224-	24"L



12 <sup>3</sup> / <sub>4</sub> " Deep	
S25C841330-	30"L
S25C841336-	36"L
S25C841342-	42"L
S25C841348-	48"L
<b>16" Deep</b>	
S25C841630-	30"L
S25C841636-	36"L
S25C841642-	42"L
S25C841648-	48"L
<b>22" Deep</b>	
S25C842230-	30"L
S25C842236-	36"L
S25C842242-	42"L
S25C842248-	48"L



# Steel Full Height Storage Cabinets



**Specifications:**


Microscope Storage Cabinets are 84" high and 16" or 22" deep. Lengths as shown. They are furnished with five adjustable shelves.

Full height swinging doors are hung on three hinges per door.

**CAUTION:** Full height cabinets may present a potential tipping over problem. Proper installation requires these cabinets to be anchored to a wall or other cabinets to improve stability and ensure smooth door movement.

**NOTE:** Full height storage cabinets are not available in Overlay – Beveled Edge Door Style 7.

**Microscope Storage Cabinets**




**16" Deep**

**S51C841624L** \_\_\_\_\_ 24"L

**22" Deep**

**S51C842224L** \_\_\_\_\_ 24"L

Provides 15 storage compartments




**16" Deep**

**S51C841624-** \_\_\_\_\_ 24"L

**22" Deep**

**S51C842224-** \_\_\_\_\_ 24"L

Provides 15 storage compartments




**16" Deep**

**S51C841636-** \_\_\_\_\_ 36"L

**22" Deep**

**S51C842236-** \_\_\_\_\_ 36"L

Provides 20 storage compartments



**16" Deep**


**S51C841648-** \_\_\_\_\_ 48"L

**22" Deep**

**S51C842248-** \_\_\_\_\_ 48"L

Provides 30 storage compartments

# Steel Full Height Storage Cabinets



**Specifications:**


Full Height Storage Cabinets are 84" high and 12<sup>3</sup>/<sub>4</sub>", 16", or 22" deep. Lengths as shown. Split door cabinets are furnished with a fixed center shelf and four adjustable shelves.

Wardrobe Storage Cabinets are 84" high and 22" deep. Lengths as shown. They are provided with a hat shelf and a horizontal bar for hanging garments.

**CAUTION:** Full height cabinets may present a potential tipping over problem. Proper installation requires these cabinets to be anchored to a wall or other cabinets to improve stability and ensure smooth door movement.

**NOTE:** Full height storage cabinets are not available in Overlay – Beveled Edge Door Style 7.

**Split Swinging Glazed Doors**



**12<sup>3</sup>/<sub>4</sub>" Deep**

S60C841330-	30"L
S60C841336-	36"L
S60C841342-	42"L
S60C841348-	48"L


**16" Deep**

S60C841630-	30"L
S60C841636-	36"L
S60C841642-	42"L
S60C841648-	48"L

**22" Deep**

S60C842230-	30"L
S60C842236-	36"L
S60C842242-	42"L
S60C842248-	48"L

**Split Swinging Panel Doors**



**12<sup>3</sup>/<sub>4</sub>" Deep**

S35C841330-	30"L
S35C841336-	36"L
S35C841342-	42"L
S35C841348-	48"L


**16" Deep**

S35C841630-	30"L
S35C841636-	36"L
S35C841642-	42"L
S35C841648-	48"L

**22" Deep**

S35C842230-	30"L
S35C842236-	36"L
S35C842242-	42"L
S35C842248-	48"L

**Wardrobe Storage Cabinets**




**Left Hand — 22" Deep**

S40C842218L	18"L
S40C842224L	24"L

**Right Hand — 22" Deep**


S40C842218-	18"L
S40C842224-	24"L



**22" Deep**

S40C842236-	36"L
S40C842248-	48"L

# Steel Full Height Storage Cabinets



**Specifications:**

Full Height Storage Cabinets are 84" high. Lengths as shown.


Acid Storage Cabinets are furnished with a polyethylene corrosion resistant lining and five full-width corrosion resistant shelves. Ventilation is provided by louvers at the bottom of the door and a 3 1/2" O.D. vent outlet at the top of the cabinet.

Full height swinging doors are hung on three hinges per door.


**CAUTION:** Full height cabinets may present a potential tipping over problem. Proper installation requires these cabinets to be anchored to a wall or other cabinets to improve stability and ensure smooth door movement.

**NOTE:** Full height storage cabinets are not available in Overlay – Beveled Edge Door Style 7.

**Acid Storage Cabinets**




<b>Left Hand — 16" Deep</b>	
S70C841624L	24"L
<b>Left Hand — 22" Deep</b>	
S70C842224L	24"L
<b>Right Hand — 16" Deep</b>	
S70C841624-	24"L
<b>Right Hand — 22" Deep</b>	
S70C842224-	24"L



<b>16" Deep</b>	
S70C841636-	36"L
S70C841648-	48"L
<b>22" Deep</b>	
S70C842236-	36"L
S70C842248-	48"L

# Base Cabinet Scribes & Fillers

## End Scribes



To close opening between wall and front of adjacent base cabinets. Includes detached toe space scribe and wall angle.

Available in widths from 1" - 24" in one inch increments

**Sitting Height — 29" High**

**BESC2900--00\_0**


**ADA Height — 32" High**

**BESC3200--00\_0**

**Standing Height — 35" High**

**BESC3500--00\_0**

## End Scribe Applique Panels



End Applique Panels may be used with End Scribes to match the face of doors and drawer heads on full overlay steel cabinets.

**End Scribe must be ordered separately**

**Steel Square Edge**

Widths from 3" - 24" in one inch increments

**Sitting Height — 29" High**

**BEAC2900--11\_0**

**Steel Beveled Edge**

Widths from 3" - 24" in one inch increments

**BEAC2900--1700**

**Wood Square Edge**

Widths from 3" - 24" in one inch increments

**BEAC2900--16\_0**

**ADA Height — 32" High — Steel**

**BEAC3200--11\_0**

**BEAC3200--1700**

**BEAC3200--16\_0**

**Standing Height — 35" High — Steel**


**BEAC3500--11\_0**

**BEAC3500--1700**

**BEAC3500--16\_0**

Insert width required before dash in part number to specify width of End Scribe or Applique Panel.  
Example: BESC3500**06**-0000 for 6" wide End Scribe or BEAC3500**06**-1100 for 6" wide Steel Square Edge Applique Panel

## Intermediate Fillers



To close opening between two adjacent base cabinets. Includes detached toe space filler.

Available in widths from 6" - 48" in one inch increments

**Sitting Height — 29" High**

**BIFC2900--00\_0**


**ADA Height — 32" High**

**BIFC3200--00\_0**

**Standing Height — 35" High**

**BIFC3500--00\_0**

## Intermediate Applique Panels



Intermediate Applique Panels may be used with Intermediate Fillers to match the face of doors and drawer heads on full overlay steel cabinets.

**Intermediate Filler must be ordered separately**

**Steel Square Edge**

Widths from 6" - 48" in one inch increments

**Sitting Height — 29" High**

**BIAC2900--11\_0**

**Steel Beveled Edge**

Widths from 6" - 48" in one inch increments

**BIAC2900--1700**

**Wood Square Edge**

Widths from 6" - 48" in one inch increments

**BIAC2900--16\_0**

**ADA Height — 32" High**

**BIAC3200--11\_0**

**BIAC3200--1700**

**BIAC3200--16\_0**

**Standing Height — 35" High**

**BIAC3500--11\_0**


**BIAC3500--1700**

**BIAC3500--16\_0**

Insert width required before dash in part number to specify width of Intermediate Filler or Applique Panel.  
Example: BIFC3500**12**-0000 for 12" wide Intermediate Filler or BIAC3500**12**-1100 for 12" wide Steel Square Edge Applique Panel

# Base Cabinet Scribes & Fillers

## Corner Fillers



To close opening between two adjacent base cabinets joined at right angles. Includes detached toe space filler.

Available in widths from 1" - 12" in one inch increments

**Sitting Height — 29" High**

**BCFC2900--00\_0**


**ADA Height — 32" High**

**BCFC3200--00\_0**

**Standing Height — 35" High**

**BCFC3500--00\_0**

## Corner Applique Panels



Corner Applique Panels may be used with Corner Fillers to match the face of doors and drawer heads on full overlay steel cabinets.

**Corner Filler must be ordered separately**

**Steel Square Edge**

Widths from 3" - 12" in one inch increments

**Sitting Height — 29" High**

**BCAC2900--11\_0**

**ADA Height — 32" High**

**BCAC3200--11\_0**

**Standing Height — 35" High**

**BCAC3500--11\_0**

**Steel Beveled Edge**

Widths from 3" - 12" in one inch increments

**BCAC2900--1700**

**BCAC3200--1700**

**BCAC3500--1700**

**Wood Square Edge**

Widths from 3" - 12" in one inch increments


**BCAC2900--16\_0**

**BCAC3200--16\_0**

**BCAC3500--16\_0**

Insert width required before dash in part number to specify size of Corner Filler or Applique Panel. Example: BCFC3500**04**-0000 for 4" Corner Filler or BCAC3500**04**-1100 for 4" Steel Square Edge Applique Panel

## Rear Scribes



To close opening between wall and rear of adjacent base cabinet. Includes wall mounting angle.

Widths from 1" - 24" in one inch increments

**Sitting Height — 29" High**

**BRSC2900--00\_0**


**ADA Height — 32" High**

**BRSC3200--00\_0**

**Standing Height — 35" High**

**BRSC3500--00\_0**

## Pipe Chase Fillers



To close opening between two adjacent base cabinets installed back-to-back or a cabinet and a finished back.

Widths from 1" - 24" in one inch increments

**Sitting Height — 29" High**

**BPCC2900--00\_0**


**ADA Height — 32" High**

**BPCC3200--00\_0**

**Standing Height — 35" High**

**BPCC3500--00\_0**

## Cabinet End Applique Panel



Wood Applique Panels to match exposed ends of base cabinets to cabinet faces.

**Wood Square Edge**

**Sitting Height — 29" High**

**BPAC290023-16\_0** 23" W

**BPAC290029-16\_0** 29" W

**ADA Height — 32" High**

**BPAC320023-16\_0** 23" W

**BPAC320029-16\_0** 29" W

**Standing Height — 35" High**


**BPAC350023-16\_0** 23" W

**BPAC350029-16\_0** 29" W

Insert width required before dash in part number to specify width of Rear Scribe or Pipe Chase Filler. Example: BRSC3500**08**-0000 for 8" wide Rear Scribe or BPCC3500**08**-0000 for 8" wide Pipe Chase Filler.

# Wall Cabinet Scribes & Fillers

## End Scribes




To close opening between wall and front of adjacent wall cabinet. Includes wall mounting angles.

**12<sup>3</sup>/<sub>4</sub>" Deep**      **16" Deep**  
Available in widths from 1" - 24" in one inch increments

24" High	30" High	36" High	48" High
WESC2413□□-00_0	WESC2416□□-00_0	WESC3013□□-00_0	WESC3016□□-00_0
WESC3613□□-00_0	WESC3616□□-00_0	WESC4813□□-00_0	WESC4816□□-00_0

## End Scribe Applique Panels




End Applique Panels may be used with End Scribes to match the face of doors and drawer heads on full overlay steel cabinets.

**End Scribe must be ordered separately**

Steel Square Edge	Steel Beveled Edge	Wood Square Edge
Widths from 3" - 24" in one inch increments		Widths from 3" - 24" in one inch increments
<b>24" High</b>		
WEAC2400□□-11_0	not available	WEAC2400□□-16_0
<b>30" High</b>		
WEAC3000□□-11_0	not available	WEAC3000□□-16_0
<b>36" High</b>		
WEAC3600□□-11_0	not available	WEAC3600□□-16_0
<b>48" High</b>		
WEAC4800□□-11_0	not available	WEAC4800□□-16_0

Insert width required before dash in part number to specify width of End Scribe or Applique Panel.  
Example: WESC361306-0000 for 6" wide End Scribe or WEAC360006-1100 for 6" wide Steel Square Edge Applique Panel

## Intermediate Fillers




To close opening between two adjacent wall cabinets.

**12<sup>3</sup>/<sub>4</sub>" Deep**      **16" Deep**  
Available in widths from 6" - 48" in one inch increments

24" High	30" High	36" High	48" High
WIFC2413□□-00_0	WIFC2416□□-00_0	WIFC3013□□-00_0	WIFC3016□□-00_0
WIFC3613□□-00_0	WIFC3616□□-00_0	WIFC4813□□-00_0	WIFC4816□□-00_0

## Intermediate Applique Panels



Intermediate Applique Panels may be used with Intermediate Fillers to match the face of doors and drawer heads on full overlay steel cabinets.


**Intermediate Filler must be ordered separately**

Steel Square Edge	Steel Beveled Edge	Wood Square Edge
Widths from 6" - 48" in one inch increments		Widths from 6" - 48" in one inch increments
<b>24" High</b>		
WIAC2400□□-11_0	not available	WIAC2400□□-16_0
<b>30" High</b>		
WIAC3000□□-11_0	not available	WIAC3000□□-16_0
<b>36" High</b>		
WIAC3600□□-11_0	not available	WIAC3600□□-16_0
<b>48" High</b>		
WIAC4800□□-11_0	not available	WIAC4800□□-16_0

Insert width required before dash in part number to specify width of Intermediate Filler or Applique Panel.  
Example: WIFC361312-0000 for 12" wide Intermediate Filler or WIAC360012-1100 for 12" wide Steel Square Edge Applique Panel

# Wall Cabinet Scribes & Fillers


## Corner Fillers



To close opening between two adjacent wall cabinets joined at right angles.

<b>12<sup>3</sup>/<sub>4</sub>" Deep</b>	<b>16" Deep</b>
Available in depth and widths from 1" - 12" in one inch increments.	
<b>24" High</b>	
WCFC2413□□-00_0	WCFC2416□□-00_0
<b>30" High</b>	
WCFC3013□□-00_0	WCFC3016□□-00_0
<b>36" High</b>	
WCFC3613□□-00_0	WCFC3616□□-00_0
<b>48" High</b>	
WCFC4813□□-00_0	WCFC4816□□-00_0

## Corner Applique Panels




Corner Applique Panels may be used with Corner Fillers to match the face of doors and drawer heads on full overlay steel cabinets.

**Corner Filler must be ordered separately**

<b>Steel Square Edge</b>	<b>Steel Beveled Edge</b>	<b>Wood Square Edge</b>
Widths from 3" - 12" in one inch increments		Widths from 3" - 12" in one inch increments
<b>24" High</b>		
WCAC2400□□-11_0	not available	WCAC2400□□-16_0
<b>30" High</b>		
WCAC3000□□-11_0	not available	WCAC3000□□-16_0
<b>36" High</b>		
WCAC3600□□-11_0	not available	WCAC3600□□-16_0
<b>48" High</b>		
WCAC4800□□-11_0	not available	WCAC4800□□-16_0

Insert width required before dash in part number to specify width of Corner Filler or Applique Panel.  
 Example: WCFC3613**06**-0000 for 6" wide Corner Filler or WCAC3600**06**-1100 for 6" wide Steel Square Edge Applique Panel

## Cabinet End Applique Panel




Wood Applique Panels to match exposed ends of wall cabinets to cabinet faces.

<b>Wood Square Edge</b>	
<b>24" High</b>	
WPAC240013-16_0	12 <sup>3</sup> / <sub>4</sub> " W
WPAC240016-16_0	16" W
<b>30" High</b>	
WPAC300013-16_0	12 <sup>3</sup> / <sub>4</sub> " W
WPAC300016-16_0	16" W
<b>36" High</b>	
WPAC360013-16_0	12 <sup>3</sup> / <sub>4</sub> " W
WPAC360016-16_0	16" W
<b>48" High</b>	
WPAC480013-16_0	12 <sup>3</sup> / <sub>4</sub> " W
WPAC480016-16_0	16" W

# Full Height Cabinet Scribes & Fillers

## Full Height End Scribes



To close opening between wall and front of adjacent full height cabinets. Includes detached toe space scribe and wall mounting angles.

Available in widths from 1" - 24" in one inch increments

12<sup>3</sup>/<sub>4</sub>" Deep

**FESC8413**□□-00\_0


16" Deep

**FESC8416**□□-00\_0

22" Deep

**FESC8422**□□-00\_0

## Full Height End Applique Panel



End Applique Panels may be used with End Scribes to match the face of doors and drawer heads on full overlay steel cabinets.

**End Scribe must be ordered separately**

**Steel Square Edge**

Widths from 3" - 24" in one inch increments

**FEAC8400**□□-11\_0

**Steel Beveled Edge**

not available


**Wood Square Edge**

Widths from 3" - 24" in one inch increments

**FEAC8400**□□-16\_0

Insert width required before dash in part number to specify width of End Scribe or Applique Panel.  
Example: FESC8413**06**-0000 for 6" wide End Scribe or FEAC8400**06**-1100 for 6" wide Steel Square Edge Applique Panel

## Full Height Intermediate Fillers



To close opening between two adjacent full height cabinets. Includes detached toe space filler.

Available in widths from 6" - 48" in one inch increments

12<sup>3</sup>/<sub>4</sub>" Deep

**FIFC8413**□□-00\_0


16" Deep

**FIFC8416**□□-00\_0

22" Deep

**FIFC8422**□□-00\_0

## Full Height Intermediate Applique Panel



Intermediate Applique Panels may be used with Intermediate Fillers to match the face of doors and drawer heads on full overlay steel cabinets.

**Intermediate Fillers must be ordered separately**

**Steel Square Edge**

Widths from 6" - 48" in one inch increments

**FIAC8400**□□-11\_0

**Steel Beveled Edge**

not available

**Wood Square Edge**

Widths from 6" - 48" in one inch increments


**FIAC8400**□□-16\_0

Insert width required before dash in part number to specify width of Intermediate Filler or Applique Panel.  
Example: FIFC8413**06**-0000 for 6" wide Intermediate Filler or FIAC8400**06**-1100 for 6" wide Steel Square Edge Applique Panel



# Full Height Cabinet Scribes & Fillers

## Corner Filler



To close opening between two adjacent full height cabinets joined at right angles. Includes detached toe space filler.

Available in widths from 1" - 12" in one inch increments.

12<sup>3</sup>/<sub>4</sub>" Deep

**FCFC8413**□□-00\_0


16" Deep

**FCFC8416**□□-00\_0

22" Deep

**FCFC8422**□□-00\_0

## Corner Fillers Applique Panel



Corner Applique Panels may be used with Corner Fillers to match the face of doors and drawer heads on full overlay steel cabinets.

**Corner Filler must be ordered separately**

### Steel Square Edge

Widths from 3" - 12" in one inch increments

**FCAC8400**□□-11\_0

### Steel Beveled Edge

not available

### Wood Square Edge


Widths from 3" - 12" in one inch increments

**FCAC8400**□□-16\_0

Insert width required before dash in part number to specify width of Corner Filler or Applique Panel.

Example: FCFC8413**06**-0000 for 6" wide Corner Filler or FCAC8400**06**-1100 for 6" wide Steel Square Edge Applique Panel

## Full Height Rear Scribe



To close opening between wall and rear of adjacent full height cabinets.

Available in widths from 1" - 24" in one inch increments

56" High

84" High


**FRSC5600**□□-00\_0

**FRSC8400**□□-00\_0

Insert width required before dash in part number to specify width of Rear Scribe.

Example: FRSC8400**06**-0000 for 6" wide Rear Scribe.

## Cabinet End Applique Panel



Wood Applique Panels to match exposed ends of full height cabinets to cabinet faces.

### Wood Square Edge

84" High

**FPAC840013-16\_0**

12<sup>3</sup>/<sub>4</sub>" W

**FPAC840016-16\_0**


16" W

**FPAC840022-16\_0**

22" W

# Steel Cabinet Finished Backs

## Base Cabinet Finished Back



To cover backs of base cabinets exposed to view. Back adds 1" to the depth of the base cabinet.

### Sitting Height — 29" High

BFBC290112-00_0	12"W
BFBC290115-00_0	18"W
BFBC290118-00_0	18"W
BFBC290121-00_0	21"W
BFBC290124-00_0	24"W
BFBC290130-00_0	30"W
BFBC290136-00_0	36"W
BFBC290142-00_0	42"W
BFBC290148-00_0	48"W
BFBC290160-00_0	60"W

### ADA Height — 32" High


BFBC320112-00_0	12"W
BFBC320115-00_0	15"W
BFBC320118-00_0	18"W
BFBC320121-00_0	21"W
BFBC320124-00_0	24"W
BFBC320130-00_0	30"W
BFBC320136-00_0	36"W
BFBC320142-00_0	42"W
BFBC320148-00_0	48"W
BFBC320160-00_0	60"W

### Standing Height — 35" High

BFBC350112-00_0	12"W
BFBC350115-00_0	15"W
BFBC350118-00_0	18"W
BFBC350121-00_0	21"W
BFBC350124-00_0	24"W
BFBC350130-00_0	30"W
BFBC350136-00_0	36"W
BFBC350142-00_0	42"W
BFBC350148-00_0	48"W
BFBC350152-00_0U3*	52"W
BFBC350158-00_0U3*	58"W
BFBC350160-00_0U3*	60"W
BFBC350160-00_0	60"W

\* Includes toespace notches for End Sink Cabinets

## Wall Cabinet Finished Back



To cover backs of wall cabinets exposed to view. Adds 1" to depth of cabinet.

### 24" High

WFBC240115-00_0	15"W
WFBC240118-00_0	18"W
WFBC240121-00_0	21"W
WFBC240124-00_0	24"W
WFBC240130-00_0	30"W
WFBC240136-00_0	36"W
WFBC240142-00_0	42"W
WFBC240148-00_0	48"W

### 30" High

WFBC300115-00_0	15"W
WFBC300118-00_0	18"W
WFBC300121-00_0	21"W
WFBC300124-00_0	24"W
WFBC300130-00_0	30"W
WFBC300136-00_0	36"W
WFBC300142-00_0	42"W
WFBC300148-00_0	48"W


### 36" High

WFBC360115-00_0	15"W
WFBC360118-00_0	18"W
WFBC360121-00_0	21"W
WFBC360124-00_0	24"W
WFBC360130-00_0	30"W
WFBC360136-00_0	36"W
WFBC360142-00_0	42"W
WFBC360148-00_0	48"W

### 48" High

WFBC480115-00_0	15"W
WFBC480118-00_0	18"W
WFBC480121-00_0	21"W
WFBC480124-00_0	24"W
WFBC480130-00_0	30"W
WFBC480136-00_0	36"W
WFBC480142-00_0	42"W
WFBC480148-00_0	48"W

## Full Height Finished Back




To cover backs of full height cabinets exposed to view. Adds 1" to depth of cabinet.

### 84" High

FFBC840118-00_0	18"W
FFBC840124-00_0	24"W
FFBC840130-00_0	30"W
FFBC840136-00_0	36"W
FFBC840142-00_0	42"W
FFBC840148-00_0	48"W

# Steel Adjustable Slip Scribes




## Specifications:

Steel Adjustable Slip Scribes are designed the same as standard scribes with the addition of a Slip Scribe Angle at the wall edge to dress-up the edge of the scribe.

The Steel Adjustable Slip Scribe is designed to be cut to width in the field to more easily fit the space between the wall and the end of the furniture assembly.

### Adjustable Base Cabinet Scribe



To close opening between wall and front of adjacent base cabinets. Includes detached toe space scribe, end scribe and slip scribe angle.

#### Sitting Height — 29" High

**ABSC290012-00\_0** 12"W


#### ADA Height — 32" High

**ABSC320012-00\_0** 12"W

#### Standing Height — 35" High

**ABSC350012-00\_0** 12"W

### Adjustable Wall Cabinet Scribe



To close opening between wall and front of adjacent wall cabinet. Includes end scribe and slip scribe angle.

#### 24" High

**AWSC241312-00\_0** 12<sup>3</sup>/<sub>4</sub>" D 12"W

**AWSC241612-00\_0** 16" D 12"W

#### 30" High

**AWSC301312-00\_0** 12<sup>3</sup>/<sub>4</sub>" D 12"W

**AWSC301612-00\_0** 16" D 12"W

#### 36" High

**AWSC361312-00\_0** 12<sup>3</sup>/<sub>4</sub>" D 12"W


**AWSC361612-00\_0** 16" D 12"W

#### 48" High

**AWSC481312-00\_0** 12<sup>3</sup>/<sub>4</sub>" D 12"W

**AWSC481612-00\_0** 116" D 12"W

### Adjustable Full Height Scribe



To close opening between wall and front of adjacent full height cabinets. Includes detached toe space scribe, end scribe, and slip scribe angle.

#### 84" High

**AFSC841312-00\_0** 12<sup>3</sup>/<sub>4</sub>" D 12"W

**AFSC841612-00\_0** 16" D 12"W

**AFSC842212-00\_0** 22" D 12"W

### Slip Scribe Angle




May be used with any standard scribe or filler to make it adjustable. Available in lengths from 3" to 96" in one inch increments.

**SCAC0101□-00\_0**

Insert length required before dash in part number to specify length of angle.

### Spring Clips




18 gauge spring steel clip, 3/4" x 1 1/2". May be used with any scribe or filler to make them easily removable.

**F-4995-PK** package of 10

# Steel Sloping Tops


## Sloping Tops



May be used on wall cabinets and full height cabinets to keep the tops dust free.

8 <sup>3</sup> / <sub>8</sub> " High 12 <sup>3</sup> / <sub>4</sub> " Deep	10 <sup>3</sup> / <sub>4</sub> " High 16" Deep	13 <sup>3</sup> / <sub>4</sub> " High 22" Deep	
SLTC081315-00_0	SLTC111615-00_0	SLTC142215-00_0	15"L
SLTC081318-00_0	SLTC111618-00_0	SLTC142218-00_0	18"L
SLTC081321-00_0	SLTC111621-00_0	SLTC142221-00_0	21"L
SLTC081324-00_0	SLTC111624-00_0	SLTC142224-00_0	24"L
SLTC081330-00_0	SLTC111630-00_0	SLTC142230-00_0	30"L
SLTC081336-00_0	SLTC111636-00_0	SLTC142236-00_0	36"L
SLTC081342-00_0	SLTC111642-00_0	SLTC142242-00_0	42"L
SLTC081348-00_0	SLTC111648-00_0	SLTC142248-00_0	48"L
SLTC081360-00_0	SLTC111660-00_0	SLTC142260-00_0	60"L
SLTC081372-00_0	SLTC111672-00_0	SLTC142272-00_0	72"L
SLTC081384-00_0	SLTC111684-00_0	SLTC142284-00_0	84"L
SLTC081396-00_0	SLTC111696-00_0	SLTC142296-00_0	96"L


## Sloping Top Fillers



To close opening between wall and front of adjacent sloping top. Includes wall mounting angle.

8 <sup>3</sup> / <sub>8</sub> " High 12 <sup>3</sup> / <sub>4</sub> " Deep	10 <sup>3</sup> / <sub>4</sub> " High 16" Deep	13 <sup>3</sup> / <sub>4</sub> " High 22" Deep	
STFC081301-00_0	STFC111601-00_0	STFC142201-00_0	1"W
STFC081306-00_0	STFC111606-00_0	STFC142206-00_0	6"W
STFC081312-00_0	STFC111612-00_0	STFC142212-00_0	12"W


## Corner Cabinet Sloping Tops



To be used on corner wall cabinets.

8 <sup>3</sup> / <sub>8</sub> " High 12 <sup>3</sup> / <sub>4</sub> " Deep	10 <sup>3</sup> / <sub>4</sub> " High 16" Deep	
CSTC081327-00_0	CSTC111627-00_0	26 <sup>7</sup> / <sub>8</sub> "L

## Sloping Top Corner Fillers




May be used on wall cabinets and full height cabinets when cabinets are joined at right angles.

8 <sup>3</sup> / <sub>8</sub> " High 12 <sup>3</sup> / <sub>4</sub> " Deep	10 <sup>3</sup> / <sub>4</sub> " High 16" Deep	13 <sup>3</sup> / <sub>4</sub> " High 22" Deep	
SCFC081301-00_0	SCFC111601-00_0	SCFC142201-00_0	1"W
SCFC081306-00_0	SCFC111606-00_0	SCFC142206-00_0	6"W

# Steel Ceiling & Pipe Enclosures

## Ceiling Enclosures


## Corner Ceiling Enclosures



To close opening between ceiling and top of wall or full height cabinets.

### Front Ceiling Enclosure

12" High	24" High	
FCEC120018-00_0	FCEC240018-00_0	18"L
FCEC120024-00_0	FCEC240024-00_0	24"L
FCEC120030-00_0	FCEC240030-00_0	30"L
FCEC120036-00_0	FCEC240036-00_0	36"L
FCEC120042-00_0	FCEC240042-00_0	42"L
FCEC120048-00_0	FCEC240048-00_0	48"L
FCEC120060-00_0	FCEC240060-00_0	60"L
FCEC120072-00_0	FCEC240072-00_0	72"L
FCEC120084-00_0	FCEC240084-00_0	84"L
FCEC120096-00_0	FCEC240096-00_0	96"L



To close opening between ceiling and top of corner wall cabinets.


12" High	24" High	Depth	Length
CCEC121320-00_0	CCEC241320-00_0	12 <sup>3</sup> / <sub>4</sub> "	20"
CCEC121615-00_0	CCEC241615-00_0	16"	15 <sup>3</sup> / <sub>8</sub> "

### End Ceiling Enclosure

12" High	24" High	
ECEC120013-00_0	ECEC240013-00_0	12 <sup>3</sup> / <sub>4</sub> "D
ECEC120016-00_0	ECEC240016-00_0	16"D
ECEC120022-00_0	ECEC240022-00_0	22"D


## Pipe Enclosures

### Wall Bench Enclosure



Pipe Enclosures conceal vertical service lines coming down from the ceiling to the worksurface. One side is removable for easy access with no exposed fasteners. Designed to sit on top of curb on worksurface and penetrate drop ceiling. Available in powder coated steel or stainless steel.

### Center Bench Enclosure



**U01C** [ ] [ ] [ ] [ ] -00\_0 Wall Bench Enclosure  
**U02C** [ ] [ ] [ ] [ ] -00\_0 Center Bench Enclosure

Available Heights	Available Depths	Available Lengths	Available Materials
60=60"	05=5"	12=12"	0 = Powder Coated Steel
66=66"	06=6"	18=18"	3 = Stainless Steel
72=72"	07=7"	24=24"	
78=78"	08=8"	30=30"	
84=84"	09=9"	36=36"	
90=90"	10=10"	42=42"	
96=96"		48=48"	

# Fume Hoods & Ventilated Enclosures

## Ducted Ventilated Products

Kewaunee offers a full line of Laboratory Fume Hoods to meet the needs of every lab application.

In addition to our newest line of **Supreme Air Venturi Fume Hoods**, we also offer uniquely designed

hoods for education and specialty hoods for custom applications.

### Supreme Air – General Purpose Fume Hoods



The design is based on rigid frame construction that assures a solid installation and low vibration and sound levels. Access panels are easily removable to access service lines and fittings. Radius corner posts and airfoils, plus, easy operator control of interior baffle settings assure a high level of comfort, safety and efficiency.

### Supreme Air Venturi – Next Generation Fume Hoods

The Supreme Air Venturi utilizes proven containment enhancing techniques coupled with innovative laboratory convenience features, including: adjustable LED lights, electromechanical sash stop, ultra-smooth sash system, and push-button sash management. The Supreme Air Venturi provides the safest user friendly fume hood manufactured today.

### TruView – Teaching Fume Hoods

Designed for today's classroom allowing clear line-of-sight through multiple hoods to allow easier supervision of students by the instructor. May be used as a stand-alone hood or in assemblies, both single-sided for use against the wall, or double-sided for use in the center of the laboratory.



For more information on Kewaunee's **Supreme Air Fume Hoods** contact your local Kewaunee sales representative or visit the Kewaunee website: <http://www.kewaunee.com/LaboratoryHoods/FumeHoods/SupremeAirHoods/>

For more information on Kewaunee's **Supreme Air Venturi Fume Hoods** contact your local Kewaunee sales representative or visit the Kewaunee website: <http://www.kewaunee.com/LaboratoryHoods/FumeHoods/SupremeAirVenturi/>

For more information on Kewaunee's **TruView Teaching Fume Hoods** contact your local Kewaunee sales representative or visit the Kewaunee website: <http://www.kewaunee.com/LaboratoryHoods/FumeHoods/TruViewHood/>


## Fume Hoods & Ventilated Enclosures

### Ductless Ventilated Products

Kewaunee offers a full line of filtered products to meet every protection requirement. **DetectAir** ductless fume hoods provide a safe alternative to conventional chemical fume hoods

#### **DetectAir – Ductless/Filtered Fume Hoods**

With the introduction of DetectAir, the filtered fume hood has come of age. DetectAir combines Kewaunee's traditional fume hood expertise with the filtering technology of AirClean Systems.




For more information on Kewaunee's **DetectAir - Ductless/Filtered Fume Hoods** contact your local Kewaunee sales representative or visit the Kewaunee website: <http://www.kewaunee.com/LaboratoryHoods/DuctlessHoods/>

when total exhaust is not an option. The **Interceptor Biological Safety Cabinets** are offered in both class II A2 and B2 models, (NSF Certified) for both recirculating and direct exhaust

#### **the Interceptor – Biological Safety Cabinet**

The Interceptor Class II, A2 and Class II, B2 Biological Safety Cabinets are Kewaunee's latest design. Developed with comprehensive features created to enhance both safety and ergonomics... the Interceptor leads its class.




For more information on **the Interceptor Biological Safety Cabinet** from Kewaunee, contact your local Kewaunee sales representative for visit the Kewaunee website: <http://www.kewaunee.com/LaboratoryHoods/BioSafetyCabinets/>

requirements. The **Interceptor VLF** provides ISO 14644-1 Standard Class 5 conditions, providing a sterile environment for product protection.

#### **the Interceptor – Vertical Laminar Flow Workstation**

Interceptor Vertical Laminar Flow Workstations feature a front loading HEPA filter, hinged front panel with self-support brackets, 10° angled safety glass, Minihelic gauge and T5 florescent lamp.



For more information on **the Interceptor - Vertical Laminar Flow Workstation** from Kewaunee, contact your local Kewaunee sales representative for visit the Kewaunee website: <http://www.kewaunee.com/LaboratoryHoods/LaminarFlowHoods/>

# Work Tops


Work Top material selections are among the most critical of the laboratory furniture decisions. The intended applications for your work tops dictate material choice. Surfaces could

be exposed to reagents and caustic substances. Heavy equipment may need to be supported or only light duty tasks may be performed.

Kewaunee can recommend and furnish materials fabricated to your individual needs, complete with appropriate sinks, drains and accessories.

## Kemresin® Work Tops

*manufactured exclusively by Kewaunee*



**Appearance:** Ultra-smooth non-porous surface.

**Available Colors:** Black, Grey, Putty, and Slate

**Composition:** 1" thick molded epoxy resin.

**Physical Properties:** Excellent resistance to practically all acids, alkalis and solvents. Highly resistant to heat and normal physical abuse.

**Applications:** Use wherever optimum physical and chemical resisting properties are required.

## Other Work Top Materials

MATERIAL	APPEARANCE	COMPOSITION	CHEMICAL RESISTANCE	APPLICATIONS
STAINLESS STEEL	No.4 Finish	Type 304 furnished as standard, but other types are available when specified. Type 316 supplied for Perchloric Acid fume hoods.	Good chemical resistance to a wide range of reagents.	Recommended for sterile and radiological laboratories.
HARDWOOD, OAK OR MAPLE	Natural varnish, or penetrating oil finish	Solid edge grain Maple or Oak strips, glued together electronically, and finished with baked-on clear synthetic varnish or penetrating oil.	Fair chemical resistance, but not recommended for use in chemical laboratories.	Excellent top material for Physics tables, Electronic benches, and other tables where mechanical or physical work levels are prevalent.
COMPOSITE RESIN	Hard, smooth non-glare surface	Cellulose fiber reinforced phenolic resin core with a polyurethane copolymer surface. Available in multiple colors and patterns.	Good chemical resistance to wide range of acids, alkalis and solvents.	Use wherever a high degree of chemical resistance is required. (Not recommended for fume hood worksurfaces.)
HIGH-PRESSURE LAMINATES	Hard, smooth surface	Plastic laminated surface attached with water resistant adhesive to solid particle board. Available in choice of colors and patterns.	Fair chemical resistance to wide range of acids, alkalis and solvents.	Recommended for light duty purposes.



## Work Tops


### Specifications:

The configurations shown are for Kewaunee base cabinets. Standard top thickness is 1", standard curb height is 4".

To order, select configuration and top material desired. Specify number of linear or square feet per product number. Dimensions are nominal. Not all con-


figurations are available in all materials. Indented tops available on special order.

### Wall Bench Tops




**Flat Top no curb** **Style B**

Standard	18"
Sizes:	23" 24" 30" 31"




**Top with Applied Curb** **Style C**

Standard	18"
Sizes:	23" 24" 30" 31"




**Corner Top no curbs** **Style K**

	<b>A</b>	<b>B</b>	<b>C</b>
Standard	33"	24"	9"
Sizes:	39" 40"	30" 31"	9" 9"



**Flat Top & Ledge no curb** **Style E**

	<b>A</b>	<b>B</b>	<b>C</b>
Standard	30"	23"	7"
Sizes:	31"	23"	8"




**Top & Ledge with Applied curb** **Style F**

	<b>A</b>	<b>B</b>	<b>C</b>
Standard	30"	23"	7"
Sizes:	31"	23"	8"

**Corner Top with Applied Curbs** **Style L**


	<b>A</b>	<b>B</b>	<b>C</b>
Standard	33"	24"	9"
Sizes:	39" 40"	30" 31"	9" 9"

### Center Bench Tops




**Top without Center Joint** **Style H**

Standard	54"
Sizes:	60" 62"



**Top with Center Joint** **Style I**


	<b>A</b>	<b>B</b>	<b>C</b>
Standard	54"	27"	27"
Sizes:	54" 60" 62"	23" 30" 31"	31" 30" 31"



**Top with Center Ledge** **Style J**


	<b>A</b>	<b>B</b>	<b>C</b>	<b>D</b>
Standard	54"	23"	8"	23"
Sizes:	60" 62"	23" 23"	14" 16"	23" 23"

### Ledges and Shelves



**Ledge Unfinished Edges** **Style A**

Standard	8"
Sizes:	10" 12" 14" 16"



**Shelf Finished Edges** **Style N**

Standard	8"
Sizes:	10" 12" 14" 16"

# Center Bench Shelving Assemblies


## Specifications:

Shelving Support Assemblies are designed to add shelving to center benches. Each assembly consists of four uprights with slots on 1" centers and a top frame. They are 12" deep in lengths as shown.

Each assembly is furnished with mounting brackets for attaching to the counter top and cabinets below but no fasteners. Fasteners must be furnished in field depending on top material and cabinet configuration.


Shelving Assembly Support Struts are especially designed to provide added support to the Shelving Support Assemblies while allowing space for connection of plumbing or electrical connections.

### 12" Deep Shelving Support Assembly




Length	Part Number
24"	<b>F17-4824-00</b>
30"	<b>F17-4830-00</b>
36"	<b>F17-4836-00</b>
42"	<b>F17-4842-00</b>
48"	<b>F17-4848-00</b>
54"	<b>F17-4854-00</b>
60"	<b>F17-4860-00</b>
66"	<b>F17-4866-00</b>
72"	<b>F17-4872-00</b>

### Shelving Assembly Support Strut



Width	Part Number
12"	<b>K17-3612-00</b>
14"	<b>K17-3614-00</b>

### Shelving Uprights



#### Powder Coated Steel Uprights

Length	Single Slot Upright	Double Slot Upright	
24"	<b>F13M000224-00</b>	<b>F12M000224-00</b>	each
30"	<b>F13M000230-00</b>	<b>F12M000230-00</b>	each
36"	<b>F13M000236-00</b>	<b>F12M000236-00</b>	each
42"	<b>F13M000242-00</b>	<b>F12M000242-00</b>	each
48"	<b>F13M000248-00</b>	<b>F12M000248-00</b>	each
72"	<b>F13M000272-00</b>	<b>F12M000272-00</b>	each
79"	<b>F13M000279-00</b>	<b>F12M000279-00</b>	each
84"	<b>F13M000284-00</b>	<b>F12M000284-00</b>	each

#### Stainless Steel Uprights

Length	Single Slot Upright	Double Slot Upright	
24"	<b>F13S000224-00</b>	<b>F12S000224-00</b>	each
30"	<b>F13S000230-00</b>	<b>F12S000230-00</b>	each
36"	<b>F13S000236-00</b>	<b>F12S000236-00</b>	each
42"	<b>F13S000242-00</b>	<b>F12S000242-00</b>	each
48"	<b>F13S000248-00</b>	<b>F12S000248-00</b>	each
72"	<b>F13S000272-00</b>	<b>F12S000272-00</b>	each
79"	<b>F13S000279-00</b>	<b>F12S000279-00</b>	each
84"	<b>F13S000284-00</b>	<b>F12S000284-00</b>	each

Designed to be used with Adjustable Module Shelving Assemblies to create wall mounted shelving.

# Center Bench Shelving Assemblies

## Specifications:

Pass-Thru Shelf Assemblies for 12" Island Modules are six inches deeper than the Shelf Bracket to allow the Shelf to pass through the module core and form a continuous surface. These shelves may also be moved forward and mounted to form


a 2" cord gap between the shelves. Shelf Brackets are 11 gauge steel, adjustable on 1" centers.

**Painted Steel Shelves** – 1" thick, formed from 16 gauge cold rolled steel.


**Veneer Plywood Shelves** – 1" thick hardwood plywood with matching veneer on all exposed surfaces. Available in Oak, Maple or Beech.

**Phenolic Resin Shelves** – 1" thick phenolic resin.

## Pass-Thru Shelf Assemblies



Pass-Thru Shelves



Pass-Thru Shelves with Cord Gap

	Length	Painted Steel	Veneer Plywood	Phenolic Resin
<b>15" DEEP SHELF</b> 9" Deep Bracket	24"	<b>F21M0924-60</b>	<b>F21W0924-6_</b>	<b>F21T0924-60</b>
	30"	<b>F21M0930-60</b>	<b>F21W0930-6_</b>	<b>F21T0930-60</b>
	36"	<b>F21M0936-60</b>	<b>F21W0936-6_</b>	<b>F21T0936-60</b>
	42"	<b>F21M0942-60</b>	<b>F21W0942-6_</b>	<b>F21T0942-60</b>
	48"	<b>F21M0948-60</b>	<b>F21W0948-6_</b>	<b>F21T0948-60</b>
	54"	<b>F21M0954-60</b>		
	60"	<b>F21M0960-60</b>		
	66"	<b>F21M0966-60</b>		
<b>18" DEEP SHELF</b> 12" Deep Bracket	24"	<b>F21M1224-60</b>	<b>F21W1224-6_</b>	<b>F21T1224-60</b>
	30"	<b>F21M1230-60</b>	<b>F21W1230-6_</b>	<b>F21T1230-60</b>
	36"	<b>F21M1236-60</b>	<b>F21W1236-6_</b>	<b>F21T1236-60</b>
	42"	<b>F21M1242-60</b>	<b>F21W1242-6_</b>	<b>F21T1242-60</b>
	48"	<b>F21M1248-60</b>	<b>F21W1248-6_</b>	<b>F21T1248-60</b>
	54"	<b>F21M1254-60</b>		
	60"	<b>F21M1260-60</b>		
	66"	<b>F21M1266-60</b>		
<b>24" DEEP SHELF</b> 18" Deep Bracket	24"	<b>F21M1824-60</b>	<b>F21W1824-6_</b>	<b>F21T1824-60</b>
	30"	<b>F21M1830-60</b>	<b>F21W1830-6_</b>	<b>F21T1830-60</b>
	36"	<b>F21M1836-60</b>	<b>F21W1836-6_</b>	<b>F21T1836-60</b>
	42"	<b>F21M1842-60</b>	<b>F21W1842-6_</b>	<b>F21T1842-60</b>
	48"	<b>F21M1848-60</b>	<b>F21W1848-6_</b>	<b>F21T1848-60</b>
	60"	<b>F21M1860-60</b>		

To specify Veneer Plywood Species replace blank at end of part number with number that indicates material required.

### Veneer Plywood Material Options

- 1 = Red Oak                      -5 = FSC® Red Oak
- 2 = White Maple                -6 = FSC® White Maple
- 7 = Steamed Beech              -8 = FSC® St'md Beech

## Shelf Retaining Lips



Length	Shelf Lip	
24"	<b>F63M0224-00</b>	each
30"	<b>F63M0230-00</b>	each
36"	<b>F63M0236-00</b>	each
42"	<b>F63M0242-00</b>	each
48"	<b>F63M0248-00</b>	each
54"	<b>F63M0254-00</b>	each
60"	<b>F63M0260-00</b>	each
66"	<b>F63M0266-00</b>	each
72"	<b>F63M0272-00</b>	each

Designed to be added to Pass-thru and Wall Shelf Assemblies. Steel Shelf Lips may be mounted to steel,

veneer plywood, or phenolic resin shelves with supplied screws.

# Wall Shelving Assemblies

## Specifications:

**Vertical Uprights** are 1 1/4" deep and either 79" or 42" high. They are manufactured with a double-slotted face with slots on 1" centers. They are designed to be mounted between two Support Rails to hold accurate spacing.


**Support Rails** are available in 8 foot lengths and may be cut to the required length. Support Rails are notched on 6" centers to accept the Vertical Uprights.

**Support Rail Covers** are required to finish the assembly and must be

ordered in pairs to match the Vertical Upright spacing.

**Optional Center Rails** may be ordered to provide a horizontal member at counter top height. They should be sized to match the Vertical Upright spacing.

## Full Height Wall Suspension System



### Vertical Upright

**F10M0279-00** 79"H pair

### Support Rail

**F11M0196-00** 96"L each

### Support Rail Covers

**F11M0324-00** 24"L pair

**F11M0330-00** 30"L pair

**F11M0336-00** 36"L pair

**F11M0342-00** 42"L pair

**F11M0348-00** 48"L pair

**F11M0354-00** 54"L pair

**F11M0360-00** 60"L pair

**F11M0366-00** 66"L pair

**F11M0372-00** 72"L pair

### Optional Center Rail Assembly

**F11M0224-00** 24"L each

**F11M0230-00** 30"L each

**F11M0236-00** 36"L each

**F11M0242-00** 42"L each

**F11M0248-00** 48"L each


**F11M0254-00** 54"L each

**F11M0260-00** 60"L each

**F11M0266-00** 66"L each

**F11M0272-00** 72"L each

## Above Countertop Wall Suspension System



### Vertical Upright

**F10M0242-00** 42"H pair

### Support Rail

**F11M0196-00** 96"L each

### Support Rail Covers

**F11M0324-00** 24"L pair

**F11M0330-00** 30"L pair

**F11M0336-00** 36"L pair

**F11M0342-00** 42"L pair

**F11M0348-00** 48"L pair

**F11M0354-00** 54"L pair

**F11M0360-00** 60"L pair

**F11M0366-00** 66"L pair

**F11M0372-00** 72"L pair

# Wall Shelving Assemblies

## Specifications:

Adjustable Wall Shelf Assemblies consist of a pair of 11 gauge steel Shelf Brackets and a Shelf. The Shelf Bracket engages the inner slot of the double slotted Wall Vertical Upright and is adjustable on 1" centers. Shelves are available in Painted



Steel, Wood veneer plywood, or Phenolic Resin.

**Painted Steel Shelves** – 1" thick, formed from 16 gauge cold rolled steel.

**Veneer Plywood Shelves** – 1" thick hardwood plywood with matching veneer on all exposed surfaces. Available in Oak, Maple or Beech.

**Phenolic Resin Shelves** – 1" thick phenolic resin.

## Adjustable Wall Shelf Assemblies



### Optional Shelf Locking Tab

Sliding tab mounted on Shelf Bracket locks Shelf into upright to keep it from being accidentally removed.

To specify, add an **K** to the end of the Shelf Assembly part number.

Example: **F21M1236-00,K**

	Length	Painted Steel	Veneer Plywood	Phenolic Resin
<b>6" DEEP SHELF</b> 6" Deep Bracket	24"	<b>F21M0624-00</b>	<b>F21W0624-0_</b>	<b>F21T0624-00</b>
	30"	<b>F21M0630-00</b>	<b>F21W0630-0_</b>	<b>F21T0630-00</b>
	36"	<b>F21M0636-00</b>	<b>F21W0636-0_</b>	<b>F21T0636-00</b>
	42"	<b>F21M0642-00</b>	<b>F21W0642-0_</b>	<b>F21T0642-00</b>
	48"	<b>F21M0648-00</b>	<b>F21W0648-0_</b>	<b>F21T0648-00</b>
	54"	<b>F21M0654-00</b>		
	60"	<b>F21M0660-00</b>		
	66"	<b>F21M0666-00</b>		
<b>9" DEEP SHELF</b> 9" Deep Bracket	24"	<b>F21M0924-00</b>	<b>F21W0924-0_</b>	<b>F21T0924-00</b>
	30"	<b>F21M0930-00</b>	<b>F21W0930-0_</b>	<b>F21T0930-00</b>
	36"	<b>F21M0936-00</b>	<b>F21W0936-0_</b>	<b>F21T0936-00</b>
	42"	<b>F21M0942-00</b>	<b>F21W0942-0_</b>	<b>F21T0942-00</b>
	48"	<b>F21M0948-00</b>	<b>F21W0948-0_</b>	<b>F21T0948-00</b>
	54"	<b>F21M0954-00</b>		
	60"	<b>F21M0960-00</b>		
	66"	<b>F21M0966-00</b>		
<b>12" DEEP SHELF</b> 12" Deep Bracket	24"	<b>F21M1224-00</b>	<b>F21W1224-0_</b>	<b>F21T1224-00</b>
	30"	<b>F21M1230-00</b>	<b>F21W1230-0_</b>	<b>F21T1230-00</b>
	36"	<b>F21M1236-00</b>	<b>F21W1236-0_</b>	<b>F21T1236-00</b>
	42"	<b>F21M1242-00</b>	<b>F21W1242-0_</b>	<b>F21T1242-00</b>
	48"	<b>F21M1248-00</b>	<b>F21W1248-0_</b>	<b>F21T1248-00</b>
	54"	<b>F21M1254-00</b>		
	60"	<b>F21M1260-00</b>		
	66"	<b>F21M1266-00</b>		
<b>18" DEEP SHELF</b> 18" Deep Bracket	24"	<b>F21M1824-00</b>	<b>F21W1824-0_</b>	<b>F21T1824-00</b>
	30"	<b>F21M1830-00</b>	<b>F21W1830-0_</b>	<b>F21T1830-00</b>
	36"	<b>F21M1836-00</b>	<b>F21W1836-0_</b>	<b>F21T1836-00</b>
	42"	<b>F21M1842-00</b>	<b>F21W1842-0_</b>	<b>F21T1842-00</b>
	48"	<b>F21M1848-00</b>	<b>F21W1848-0_</b>	<b>F21T1848-00</b>
	54"	<b>F21M1854-00</b>		
	60"	<b>F21M1860-00</b>		
	<b>24" DEEP SHELF</b> 24" Deep Bracket	24"	<b>F21M2424-00</b>	<b>F21W2424-0_</b>
30"		<b>F21M2430-00</b>	<b>F21W2430-0_</b>	<b>F21T2430-00</b>
36"		<b>F21M2436-00</b>	<b>F21W2436-0_</b>	<b>F21T2436-00</b>
42"		<b>F21M2442-00</b>	<b>F21W2442-0_</b>	<b>F21T2442-00</b>
48"		<b>F21M2448-00</b>	<b>F21W2448-0_</b>	<b>F21T2448-00</b>
60"		<b>F21M2460-00</b>		

To specify Veneer Plywood Species replace blank at end of part number with number that indicates material required.

### Veneer Plywood Material Options

- 1 = Red Oak
- 2 = White Maple
- 7 = Steamed Beech
- 5 = FSC® Red Oak
- 6 = FSC® White Maple
- 8 = FSC® St'md Beech

# Reagent Rack & Shelving Systems

**Reagent Rack Uprights** are available in two heights and two depths. All are 1½" thick and available in either powder coated steel to match or complement the casework color or in type 304 stainless steel. 12" deep uprights are designed for


use on Center Benches with a 14" deep shelf, 10" deep uprights are designed for use on Wall Benches with a 12" deep shelf, but both can be used on either bench.

Installation hardware is included with each upright.


Spacing between upright should not exceed 48".

**Plastic Laminate Shelves** are 1" thick particle board faced on top and all four edges with black plastic laminate. A backer sheet covers the bottom.


## Typical Reagent Rack Assemblies



**Single Tier Center Bench Reagent Rack** using 18" high by 12" deep uprights and 14" deep shelving.




**Double Tier Center Bench Reagent Rack** using 18" high by 12" deep and 12" by 12" deep uprights and 14" deep shelving.



**Single Tier Wall Bench Reagent Rack** using 18" high by 10" deep uprights and 12" deep shelving.

## Reagent Rack Uprights


18" High x 12" Deep



18" H x 12" D x 1½" T

**RRU-1812-0000** P.C. Steel  
**RRU-1812-0030** St. Steel


12" High x 12" Deep



12" H x 12" D x 1½" T

**RRU-1212-0000** P.C. Steel  
**RRU-1212-0030** St. Steel


18" High x 10" Deep



18" H x 10" D x 1½" T

**RRU-1810-0000** P.C. Steel  
**RRU-1810-0030** St. Steel

12" High x 10" Deep




12" H x 10" D x 1½" T

**RRU-1210-0000** P.C. Steel  
**RRU-1210-0030** St. Steel


## Plastic Laminate Shelves

14" Deep Shelf



**PL-4814-00** 14" Deep x 48" Long  
**PL-7214-00** 14" Deep x 72" Long

12" Deep Shelf



**PL-4812-00** 12" Deep x 48" Long  
**PL-7212-00** 12" Deep x 72" Long

# Reagent Rack & Shelving Systems

**Fixed Wall Mounted Shelf Brackets** are designed to mount directly to the wall and furnished without mounting hardware.

**Adjustable Wall Mounted Shelf Brackets** are designed to mount to KF-0170-04 Shelving Uprights and


include mounting bolts and nuts.

**KF-0170-04 Shelving Uprights** are designed to hold a maximum of three shelves.


Brackets and Uprights are finished in Kewaunee's unique corrosion

resistant powder coat finish to match or complement the casework color.


## Typical Reagent Rack Assemblies



Fixed Wall Shelving




Adjustable Wall Shelving



Adjustable Wall Shelving


## Flat Top Shelf Bracket    Lipped Shelf Bracket    Book End Bracket    Shelving Upright



**Fixed Shelf Brackets**

	bracket size	shelf size
<b>M-0031-0A</b>	11" x 8"	12"
<b>M-0030-0A</b>	15" x 10"	16"
<b>M-0032-0A</b>	17" x 12"	18"


part number is for single bracket



**Fixed Shelf Brackets**

	bracket size	shelf size
<b>M-0038-0A</b>	12 <sup>1</sup> / <sub>8</sub> " x 6"	12"
<b>M-0044-0A</b>	16 <sup>1</sup> / <sub>8</sub> " x 6"	16"
<b>M-0039-0A</b>	18 <sup>1</sup> / <sub>8</sub> " x 6"	18"

part number is for single bracket



**Fixed Shelf Brackets**


	bracket size	shelf size
<b>M-0060-0A</b>	12" x 6"	12"
<b>M-0061-0A</b>	16" x 6"	16"
<b>M-0062-0A</b>	18" x 6"	18"

part number is for pair of brackets



**KF-0170-04**                      48" L


High-strength, cold-rolled steel shelf support upright allows adjustable shelf bracket to be positioned along 48" length for maximum use of space. A pair of uprights holds a maximum of three shelf brackets.



**Adjustable Shelf Brackets**

	bracket size	shelf size
<b>M-0031-AD</b>	11" x 8"	12"
<b>M-0030-AD</b>	15" x 10"	16"
<b>M-0032-AD</b>	17" x 12"	18"


part number is for single bracket



**Adjustable Shelf Brackets**

	bracket size	shelf size
<b>M-0038-AD</b>	12 <sup>1</sup> / <sub>8</sub> " x 6"	12"
<b>M-0044-AD</b>	16 <sup>1</sup> / <sub>8</sub> " x 6"	16"
<b>M-0039-AD</b>	18 <sup>1</sup> / <sub>8</sub> " x 6"	18"

part number is for single bracket




**Adjustable Shelf Brackets**

	bracket size	shelf size
<b>M-0060-AD</b>	12" x 6"	12"
<b>M-0061-AD</b>	16" x 6"	16"
<b>M-0062-AD</b>	18" x 6"	18"

part number is for pair of brackets

# Laboratory Pegboards


## Kemresin Laboratory Pegboards



**Pegboard without Drip Trough**  
29" High x 1" Thick x 20" Long

	Pegboard Color	Peg Color
<b>X-020014-BK</b>	Black	Black
<b>X-020014-GR</b>	Grey	White
<b>X-020014-PT</b>	Putty	White
<b>X-020014-SL</b>	Slate	White


Peg Size: 1/2" dia. x 6 1/8" long  
Peg Material: Polypropylene  
Number of Pegs: 32



**Pegboard without Drip Trough**  
29" High x 1" Thick x 32" Long

	Pegboard Color	Peg Color
<b>X-020015-BK</b>	Black	Black
<b>X-020015-GR</b>	Grey	White
<b>X-020015-PT</b>	Putty	White
<b>X-020015-SL</b>	Slate	White


Peg Size: 1/2" dia. x 6 1/8" long  
Peg Material: Polypropylene  
Number of Pegs: 53



**Pegboard with Drip Trough**  
30" High x 1" Thick x 20" Long

	Pegboard Color	Peg Color
<b>X-020016-BK</b>	Black	Black
<b>X-020016-GR</b>	Grey	White
<b>X-020016-PT</b>	Putty	White
<b>X-020016-SL</b>	Slate	White

Peg Size: 1/2" dia. x 6 1/8" long  
Peg Material: Polypropylene  
Number of Pegs: 32  
Trough Size: 1" high x 2 1/4" deep



**Pegboard with Drip Trough**  
30" High x 1" Thick x 32" Long

	Pegboard Color	Peg Color
<b>X-020017-BK</b>	Black	Black
<b>X-020017-GR</b>	Grey	White
<b>X-020017-PT</b>	Putty	White
<b>X-020017-SL</b>	Slate	White

Peg Size: 1/2" dia. x 6 1/8" long  
Peg Material: Polypropylene  
Number of Pegs: 53  
Trough Size: 1" high x 2 1/4" deep



**Stainless Steel Drip Trough**  
for use on  
**F-2939-00** X-020014 Pegboard


Size: 1 1/2"H x 2 1/2"W x 20"L  
Material: Type 304 Stainless Steel  
Includes: 1/2" OD drain outlet



**Stainless Steel Drip Trough**  
for use on  
**F-2940-00** X-020015 Pegboard

Size: 1 1/2"H x 2 1/2"W x 32"L  
Material: Type 304 Stainless Steel  
Includes: 1/2" OD drain outlet


## Type 304 Stainless Steel Laboratory Pegboards



**Pegboard without Drip Trough**  
30" High x 3/4" Thick x 24" Long

<b>X-020008</b>	with Open Back
<b>X-020010</b>	with Finished Back


Peg Size: 1/2" diameter x 6" long  
Peg Material: Polypropylene  
Peg Color: White  
Number of Pegs: 32



**Pegboard without Drip Trough**  
30" High x 3/4" Thick x 30" Long

<b>X-020009</b>	with Open Back
<b>X-020011</b>	with Finished Back


Peg Size: 1/2" diameter x 6" long  
Peg Material: Polypropylene  
Peg Color: White  
Number of Pegs: 50



**Pegboard with Drip Trough**  
30" High x 3/4" Thick x 24" Long

<b>X-020004</b>	with Open Back
<b>X-020006</b>	with Finished Back

Peg Size: 1/2" diameter x 6" long  
Peg Material: Polypropylene  
Peg Color: White  
Number of Pegs: 32



**Pegboard with Drip Trough**  
30" High x 3/4" Thick x 30" Long


<b>X-020005</b>	with Open Back
<b>X-020007</b>	with Finished Back

Peg Size: 1/2" diameter x 6" long  
Peg Material: Polypropylene  
Peg Color: White  
Number of Pegs: 50



# Miscellaneous Accessories


## Writing Boards & Tackboards



**Black Chalkboard**  
**F-5336-00** 18<sup>1</sup>/<sub>4</sub>"H x 24<sup>1</sup>/<sub>4</sub>"L  
**F-5337-00** 24<sup>1</sup>/<sub>4</sub>"H x 36<sup>1</sup>/<sub>4</sub>"L  
**F-5338-00** 36<sup>1</sup>/<sub>4</sub>"H x 48<sup>1</sup>/<sub>4</sub>"L

---


**Surface:** washable Duracite  
**Frame Material:** anodized aluminum  
**Trough Material:** anodized aluminum



**White Magnetic Markerboard**  
**F-5342-00** 18<sup>1</sup>/<sub>4</sub>"H x 24<sup>1</sup>/<sub>4</sub>"L  
**F-5343-00** 24<sup>1</sup>/<sub>4</sub>"H x 36<sup>1</sup>/<sub>4</sub>"L  
**F-5344-00** 36<sup>1</sup>/<sub>4</sub>"H x 48<sup>1</sup>/<sub>4</sub>"L

---


**Surface:** 1/4" white magnetic  
**Frame Material:** grey mylar vinyl



**Green Chalkboard**  
**F-5333-00** 18<sup>1</sup>/<sub>4</sub>"H x 24<sup>1</sup>/<sub>4</sub>"L  
**F-5340-00** 24<sup>1</sup>/<sub>4</sub>"H x 36<sup>1</sup>/<sub>4</sub>"L  
**F-5341-00** 36<sup>1</sup>/<sub>4</sub>"H x 48<sup>1</sup>/<sub>4</sub>"L

---

**Surface:** washable Duracite  
**Frame Material:** anodized aluminum  
**Trough Material:** anodized aluminum




**Cork Tackboard**  
**F-5418-24** 18"H x 24"L  
**F-5424-36** 24"H x 36"L  
**F-5436-48** 36"H x 48"L

---

**Surface:** Resilient natural cork  
**Frame Material:** anodized aluminum

## Kemresin Balance Table




30" High x 24" Deep x 35" Long  
**K7-302435-BK** Black  
**K7-302435-GR** Grey  
**K7-302435-PT** Putty  
**K7-302435-SL** Slate

---

**Work Top:** 2" thick Kemresin  
**Work Top Size:** 2"T x 24"D x 35"L  
**Pedestal Legs:** 2" thick Kemresin  
**Total Weight:** 300 pounds

## Marvel Refrigerators & Freezers



**L-110** Refrigerator/Freezer  
**Size:** 34<sup>1</sup>/<sub>2</sub>" H x 24" D x 23<sup>7</sup>/<sub>8</sub>" L  
 6.1 cubic feet  
**Door Swing:** Right Hand  
 (User changeable with swing reversing kit)  
**Shelves:** (3) Adjustable  
**Service:** Single-phase,  
 60 hertz, 120 volt AC

---

**L-105** Refrigerator  
 Same as above except  
 without freezer compartment.

---


**L-115** Freezer  
 Same as above except  
 entire unit is a freezer.

## Gillis Stainless Steel Cart



**F-5349-00**  
**Overall Size:** 33<sup>3</sup>/<sub>8</sub>"H x 19"D x 31"L  
**Shelf Size:** (3) 18"D x 27"L  
 (1" high lip on 3 sides of  
 top & middle shelf)  
**Capacity:** 400 pounds  
**Casters:** (4) Black rubber  
 4" Diameter – Swivel

## Lund Key Cases




**F-5332-00**  
**Height:** 11"  
**Depth:** 2<sup>1</sup>/<sub>2</sub>"  
**Length:** 17<sup>5</sup>/<sub>16</sub>"  
**Key Capacity:** 90 Keys

---

**F-5332-01**  
**Height:** 24<sup>3</sup>/<sub>4</sub>"  
**Depth:** 2<sup>1</sup>/<sub>2</sub>"  
**Length:** 17<sup>5</sup>/<sub>16</sub>"  
**Key Capacity:** 200 Keys

# Miscellaneous Accessories

## Pipe Space Kemstruts




Floor mounted support to provide stability to rear overhang of fume hood tops or worksurfaces, and for mounting piping and drainlines. Adjustable to match variances in height.

**K12-**    **-0A**

Available Heights	Available Depths	
<i>adjustable</i>	<i>depth</i>	<i>depth</i>
<b>29</b> = 27" - 31"	<b>05</b> = 5"	<b>13</b> = 13"
<b>32</b> = 31" - 35"	<b>06</b> = 6"	<b>14</b> = 14"
<b>35</b> = 34" - 38"	<b>07</b> = 7"	<b>15</b> = 15"
	<b>08</b> = 8"	<b>16</b> = 16"
	<b>09</b> = 9"	<b>17</b> = 17"
	<b>10</b> = 10"	<b>18</b> = 18"
	<b>11</b> = 11"	<b>19</b> = 19"
	<b>12</b> = 12"	<b>20</b> = 20"

## Corner Kemstruts




**K-0020-1A** 28" - 30" H  
**K-0016-1A** 33" - 36" H

Floor and wall mounted support system for adding stability to corner tops in L-shaped furniture assemblies. Adjustable to match variances in height. Includes two support angles to attach to adjoining base cabinets.


**NOTE:** Kemstrut Corner Top Support Assemblies are not required when a corner base cabinet is used.

## Kemstrut Pipe Hangers



**K-0005-0A** Drainline Support  
**M-0051-0A** Pipe Clamp for  $\frac{3}{8}$ " -  $\frac{1}{2}$ " IPS  
**M-0052-0A** Pipe Clamp for  $\frac{3}{4}$ " - 1" IPS


## Support Angles



**L-2034-06** 6" L  
**L-2034-12** 12" L  
**L-2034-18** 18" L  
**L-2034-24** 24" L  
**L-2034-48** 48" L  
**L-2034-96** 96" L

Formed from 16 gauge cold-rolled steel, with  $1\frac{1}{2}$ " sides.

## LED Task Lighting




Part Number	Length
<b>X-030001-00</b>	300mm - 11 $\frac{7}{8}$ "
<b>X-030002-00</b>	500mm - 19 $\frac{5}{8}$ "
<b>X-030003-00</b>	1000mm - 39 $\frac{3}{8}$ "

**LED Task Lighting** consists of a thin-profile, white, LED strip, 9.5mm high by 33mm deep with a touch dimmer switch. It is UL listed and furnished with 120 VAC to 24 VDC power supply with ten foot power cord. The light includes a magnetic mounting strip for mounting to steel shelving and integral tabs for mounting to wood and plastic laminate.

# Miscellaneous Accessories


## Rod Assembly



**4740-00** Rod Assembly

Consists of two 0553-00 upright rods, one 0557-00 cross rod and two 0550-00 rod clamps.

## Support Rods



- 0549-18** 1/2" Dia. x 18" H
- 0549-24** 1/2" Dia. x 24" H
- 0549-36** 1/2" Dia. x 36" H
- 0553-00** 3/4" Dia. x 36" H

Duraluminum support rods are available in 1/2" and 3/4" diameter and have a tapered end to fit 0548-00 rod sockets.


## Cross Rods



- 0554-00** 1/2" Dia. x 39" L
- 0556-00** 1/2" Dia. x 18" L
- 0557-00** 3/4" Dia. x 39" L

Duraluminum cross rods are available in 1/2" and 3/4" diameter and have rounded ends


## Rod Clamp



**0550-00** Rod Clamp

Made of cast aluminum. Designed to clamp horizontal cross rods on vertical upright rods. Can be used for either 1/2" or 3/4" diameter rods


## Rod Sockets



**0548-00** Rod Socket

Made of chrome-plated brass. Use with 1/2" and 3/4" diameter upright rods. 1" diameter x 1 7/8" long threaded shank. Complete with locknut and washer.

## Label Holder



**F-0260-00** Label Holder

Aluminum with adhesive back for easy installation. Holds 1 1/8" H x 2 1/2" L identification card.

## Base Molding




**F-0862-04** 48" L

4" high black base molding. Made of pliable material to match variances in floor surface. Comes with high-strength adhesive applied to back.

**Note: Specify number of 4' sections required.**


## Corner Guard



**F-1368-00** Corner Guard

4" high black corner guard to finish off corner joints of base molding.

## Wall Cabinet Hanger Rail



- L01-020109-00** for 15" wall cabinet
- L01-020112-00** for 18" wall cabinet
- L01-020115-00** for 21" wall cabinet
- L01-020118-00** for 24" wall cabinet
- L01-020124-00** for 30" wall cabinet
- L01-020130-00** for 36" wall cabinet
- L01-020136-00** for 42" wall cabinet
- L01-020142-00** for 48" wall cabinet

Formed of heavy gauge cold-rolled steel. Designed to capture lip on top rear of wall cabinets for easy installation.

# General Plumbing Information

### Recommended Supply Line Pressures

**Gas** ..... 0.216 psi  
**Air** (Ordinary) ..... 15 psi  
**Air** (Air Powered Stirrer Operation) ..... 40 psi  
**Vacuum** (Ordinary) ..... 18"-22" of Mercury  
**Vacuum** (High) ..... 28.5" of Mercury  
**Hot and Cold Water** ..... Local Pressure  
**Steam** (Low Pressure) ..... 0-15 psi as required  
**Steam** (High Pressure) ..... 15-30 psi as required

### Service & Rough-in Abbreviations

**A** – Air  
**E** – Electric - 120 volt AC  
**G** – Gas  
**N** – Nitrogen  
**O** – Oxygen  
**S** – Steam  
**V** – Vacuum  
**W** – Waste  
**CW** – Cold Water  
**DW** – Distilled Water  
**HV** – Hood Vent  
**HW** – Hot Water  
**WV** – Waste Vent

### Service Fitting Index Color Coding

When ordering service fittings, specify service required, so that the proper index button(s) will be supplied.

**Air** – Orange  
**Cold Water** – Dark Green  
**Distilled Water** – White  
**Gas** – Dark Blue  
**Hot Water** – Red  
**Hydrogen** – Pink  
**Nitrogen** – Brown  
**Oxygen** – Light Green  
**Special Gases** – Light Blue  
**Steam** – Black  
**Vacuum** – Yellow

### Recommended Pipe & Drain Line Sizes

Service	Number of Outlets	Pipe Size	Pipe Material
Steam	1-8	3/8"	In accordance with governing codes
	9-16	1/2"	
	17-24	3/4"	
*Water	1-10	1/2"	In accordance with governing codes
	11-18	3/4"	
	19-24	1"	
Gas	1-12	3/8"	In accordance with governing codes
	13-18	1/2"	
	19-24	3/4"	
	25-30	1"	
Air or Vacuum	1-8	3/8"	In accordance with governing codes
	9-16	1/2"	
	17-24	3/4"	
	<b>Fittings or Filter Pumps</b>	<b>Pipe Size</b>	<b>Pipe Material</b>
Drain Lines	1-6 Fittings or 4 Pumps	1 1/2"	In accordance with governing codes
	7-10 Fittings or 8 Pumps	2"	


\* If table is equipped with a series of filter pumps (aspiration), next larger pipe size should be used.

**NOTE:** For information not supplied, or if additional details are required, contact your local Kewaunee Representative.

### Open Reagent Rack Horizontal Pipe Spacing


Water pipe lines should always be the lowest and a minimum of 9" above the working surface. Electric conduit should always be the uppermost service line. (Insulated pipe will require a greater spacing than indicated below.)

Pipe Diameters	Minimum Vertical Spacing Between Centerlines of 2 Service Pipe Lines
3/8" and 3/8"	2"
3/8" and 1/2"	2"
1/2" and 1/2"	2"
1/2" and 3/4"	2 1/2"
3/4" and 3/4"	2 1/2"
3/4" and 1"	3"
1" and 1"	3"




### Minimum Spacing Between Fittings

These drawings show the allowable minimum spacing between service and electrical fittings on laboratory bench tops and ledges.




# Laboratory Service Fittings


## Laboratory Ball Valves (for use with air, vacuum or gases – not to be used for liquids, steam or oxygen)




**W-0242-00 Gas, Air or Vacuum Valve** is chrome plated brass. Valve provides on/off control of gas, air or vacuum at pressures up to 75 PSI.  $\frac{3}{8}$ " I.P.S. male inlet.  
**Specify service desired.**




**W-0286-02 Double Gas, Air or Vacuum Fitting (curb or panel mounted)** is chrome plated brass. Consists of two-way 70° outlet panel flange and two laboratory ball valves. Complete with  $\frac{3}{8}$ " I.P.S. tank nipple, locknut and washer.  
**Specify service desired.**




**W-0287-00 Gas, Air or Vacuum Fitting (curb or panel mounted)** is chrome plated brass. Consists of a single outlet panel flange and a laboratory ball valve. Complete with a  $\frac{3}{8}$ " I.P.S. tank nipple, locknut and washer.  
**Specify service desired.**




**W-0260-00 Gas, Air or Vacuum Deck Mounted Fitting** is chrome plated brass. Consists of a single outlet turret with a laboratory ball valve. Complete with a  $\frac{3}{8}$ " I.P.S. tank nipple, locknut and washer.  
**Specify service desired.**



**W-0263-00 Double Gas, Air or Vacuum Deck Mounted Fitting** Same as No. W-0260-00 except consists of a 2-way 90° outlet turret and two laboratory ball valves. Complete with a  $\frac{3}{8}$ " I.P.S. tank nipple, locknut and washer.  
**Specify service desired.**




**W-0266-00 Double Gas, Air or Vacuum Deck Mounted Fitting** Same as No. W-0260-00, except consists of a 2-way 180° outlet turret and two laboratory ball valves. Complete with a  $\frac{3}{8}$ " I.P.S. tank nipple, locknut and washer.  
**Specify service desired.**




**W-0272-00 Quadruple Gas, Air or Vacuum Deck Mounted Fitting** Same as W-0260-00, except consists of a 4-way outlet turret and four valves. Complete with a  $\frac{3}{8}$ " I.P.S. tank nipple, locknut and washer.  
**Specify service desired.**


## Needle Point Cock Valves (for use with air, vacuum or gases – not to be used for liquids or steam)



**W-0243-00 Gas, Air or Vacuum Valve** is chrome plated brass. Single needle point valve with integral straight serrated hose ends. Brass floating cone and integral brass seat. Valve provides metering control of laboratory gases at pressures up to 60 PSI.  $\frac{3}{8}$ " I.P.S. male inlet.  
**Specify service desired.**



**W-0400-00 Double Gas, Air or Vacuum Fitting (curb or panel mounted)** with two needle point valves with integral serrated hose ends. Brass floating cone and integral brass seat. Complete with a  $\frac{3}{8}$ " I.P.S. tank nipple, locknut and washer.  
**Specify service desired.**




**W-0401-00 Gas, Air or Vacuum Fitting (curb or panel mounted)** Same as No. W-0400-00 except single needle point valve. Complete with a  $\frac{3}{8}$ " I.P.S. tank nipple, locknut and washer.  
**Specify service desired.**

**NOTE:** If valve is to be used for oxygen service, order W-0243-01.


# Laboratory Service Fittings

## Needle Point Cock Valves


(for use with air, vacuum or gases – not to be used for liquids or steam)




**W-0259-00 Single Needle Cock** Deck Mounted Fitting  
Turret Base is chrome plated brass. Single needle point valve with integral straight serrated hose end. Complete with a  $\frac{3}{8}$ " I.P.S. tank nipple, locknut and washer.  
**Specify service desired.**



**W-0402-00 Double Needle Cock** Deck Mounted Fitting  
Turret Base is chrome plated brass. Two needle point valves with integral straight serrated hose ends at 90°. Complete with a  $\frac{3}{8}$ " I.P.S. tank nipple, locknut and washer.  
**Specify service desired.**




**W-0403-00 Double Needle Cock** Deck Mounted Fitting  
Turret Base is chrome plated brass. Same as No. W-0402-00 except integral straight serrated hose ends are at 180°. Complete with a  $\frac{3}{8}$ " I.P.S. tank nipple, locknut and washer.  
**Specify service desired.**




**W-0404-00 Quadruple Needle Cock** Deck Mounted Fitting  
Turret Base is chrome plated brass. Four needle point valves with integral straight serrated hose ends. Complete with a  $\frac{3}{8}$ " I.P.S. tank nipple, locknut and washer.  
**Specify service desired.**

## Water Valve




**W-0306-00 Water Fitting (deck or panel mounted)**  
Made of chrome plated brass. Consists of a single outlet turret and a straight water valve with removable straight serrated hose connector. Complete with a  $\frac{3}{8}$ " I.P.S. tank nipple, locknut and washer.

## Steam Valves




**W-0300-00 Steam Fitting (deck or panel mounted)**  
Made of chrome plated brass. Consists of a single outlet turret and a straight steam valve. Designed for maximum steam pressure of 20 psi. Complete with a  $\frac{3}{8}$ " I.P.S. tank nipple, locknut and washer.




**W-0301-00 Steam Fitting (curb or panel mounted)**  
Made of chrome plated brass. Consists of a single outlet panel flange and a straight steam valve. Designed for maximum steam pressure of 20 psi. Complete with a  $\frac{3}{8}$ " I.P.S. tank nipple, locknut and washer.

## Multiple Service Fitting



**W-0273-02 Deck Mounted Combination Water/Gas Fitting** Made of chrome plated brass. Water faucet has a 6" spread gooseneck with removable aerator outlet and disc type valve. Includes a two-way turret with two laboratory ball valves set at 180°. Complete with two  $\frac{3}{8}$ " I.P.S. tank nipples, locknuts and washers.  
**(Specify service desired for ball valves.)**

**W-0273-VB Deck Mounted Combination Water/Gas Fitting**  
Same as W-0273-02 except furnished with vacuum breaker type gooseneck and removable straight serrated hose connector.  
**(Specify service desired for ball valves.)**




**W-0274-VB Deck Mounted Combination Water/Gas Fitting** Made of chrome plated brass. Water faucet has a 8" spread vacuum breaker type gooseneck with removable straight serrated hose connector and disc type valve. Includes a two-way turret with two laboratory ball valves set at 180°. Complete with two  $\frac{3}{8}$ " I.P.S. tank nipples, locknuts and washers.  
**(Specify service desired for valves.)**


**W-0274-02 Deck Mounted Combination Water/Gas Fitting**  
Same as W-0274-VB except furnished with standard gooseneck and removable aerator outlet.  
**(Specify service desired for valves.)**

# Laboratory Service Fittings


## Single Water Faucet




**W-0337-00 Deck Mounted Gooseneck Water Faucet**  
 Chrome plated brass with an 8" spread rigid gooseneck. Removable straight serrated hose connector. Valve is disc type with a renewable internal assembly. Faucet body is furnished with a  $\frac{3}{8}$ " I.P.S. threaded mounting shank with locknut and washer.




**W-0337-0V Deck Mounted Vacuum Breaker Gooseneck Water Faucet**  
 Same as W-0337-00 except furnished with 8" spread rigid vacuum breaker type gooseneck.




**W-0337-06 Deck Mounted Gooseneck Water Faucet**  
 Same as W-0337-00 except furnished with 6" spread rigid gooseneck




**W-0337-6V Deck Mounted Vacuum Breaker Gooseneck Water Faucet**  
 Same as W-0337-0V except furnished with 6" spread rigid vacuum breaker type gooseneck.




**W-0350-00 Right Hand Curb Mounted Gooseneck Water Faucet**  
 Made of chrome plated brass. Has an 8" spread rigid gooseneck with removable straight serrated hose connector. Valve is disc type with renewable internal assembly. Flanged inlet is complete with  $\frac{1}{2}$ " I.P.S. tank nipple, locknut and washer.




**W-0350-0V Right Hand Curb Mounted Vacuum Breaker Gooseneck Water Faucet**  
 Same as W-0350-00 except supplied with vacuum breaker type gooseneck.



**W-0351-00 Left Hand Curb Mounted Gooseneck Water Faucet**  
 Same as W-0350-00 except supplied with valve handle on left.




**W-0351-0V Left Hand Curb Mounted Vacuum Breaker Gooseneck Water Faucet**  
 Same as W-0350-0V except supplied with valve handle on left.




**W-0405-00 Portable Pump Gooseneck Faucet** for laboratory demonstration table. Made of chrome plated brass with  $7\frac{5}{8}$ " spread, rigid,  $\frac{1}{4}$ " I.P.S. gooseneck. Faucet/pump body is complete with an integral threaded mounting shank and flanged locknut.  $\frac{5}{8}$ " O.D. hose connector inlet.

# Laboratory Service Fittings


## Mixing Faucet




**W-0340-00 Deck Mounted Hot and Cold Water Mixing Gooseneck Faucet** Made of chrome plated brass with an 8" spread rigid gooseneck. Removable straight serrated hose connector. Valves are disc type with renewable internal assemblies. Faucet body is complete with an integral threaded mounting shank and a flanged locknut. Furnished with 1/2" O.D. slip joint coupling nuts for 3/8" tubing.




**W-0340-0V Deck Mounted Hot and Cold Water Mixing Vacuum Breaker Gooseneck Faucet** Same as W-0340-00 except supplied with vacuum breaker type gooseneck.




**W-0340-06 Deck Mounted Hot and Cold Water Mixing Gooseneck Faucet** Same as W-0340-00 except supplied with 6" spread rigid gooseneck.




**W-0340-6V Deck Mounted Hot and Cold Water Mixing Vacuum Breaker Gooseneck Faucet** Same as W-0340-06 except supplied with vacuum breaker type gooseneck.




**W-0333-00 Deck Mounted Hot and Cold Water Mixing Swing Spout Faucet** Made of chrome plated brass. Has a 9" spread swinging spout with anti-splash outlet. Valves are disc type with renewable internal assemblies. Faucet body is complete with an integral threaded shank and a flanged locknut. Furnished with 1/2" O.D. slip joint coupling nuts for 3/8" tubing.



**W-0354-00 Curb Mounted Hot and Cold Water Mixing Gooseneck Faucet** Made of chrome plated brass. Has an 8" spread rigid gooseneck with removable straight serrated hose connector. Valves are disc type with renewable internal assemblies. Flanged inlets are complete with 3/8" I.P.S. tank nipples, locknuts and washers on 8" centers.



**W-0354-0V Curb Mounted Hot and Cold Water Mixing Vacuum Breaker Gooseneck Faucet** Same as No. W-0354-00 except supplied with vacuum breaker type gooseneck.




**W-0355-00 Curb Mounted Hot and Cold Water Mixing Swing Spout Faucet** Made of chrome plated brass. Has a 9" spread swinging spout with anti-splash outlet. Valves are disc type with renewable internal assemblies. Flanged inlets are complete with 3/8" I.P.S. tank nipples, locknuts and washers on 8" centers.



# Laboratory Service Fittings

## Mixing Faucet




**W-0361-00 Panel Mounted Hot and Cold Water Mixing Gooseneck Faucet** Chrome plated brass. Designed for mounting on fascia type reagent rack plumbing box uprights. It has an 8" spread rigid gooseneck equipped with a removable straight serrated hose connector. Valves are disc type with renewable internal assemblies. Flanged inlets on 3" centers are complete with  $\frac{3}{8}$ " I.P.S. tank nipples, locknuts and washers.

**W-0361-0V Panel Mounted Hot and Cold Water Mixing Vacuum Breaker Gooseneck Faucet** Same as W-0361-00 except supplied with vacuum breaker type gooseneck.

**W-0825-00 Floor Mounted Foot Control Double Pedal Mixing Valve** is chrome plated brass. Has twin pedals with self-closing valves for mixing hot and cold water.  $\frac{1}{2}$ " I.P.S. female thread inlets and outlet at the rear. Inlets are on 2" centers,  $\frac{7}{8}$ " above floor. Outlet is on center  $1\frac{1}{2}$ " above floor.

**W-0833-00 Deck Mounted Water Gooseneck Fitting** is chrome plated brass. Consisting of a 8" gooseneck with mounting flange and W-0532-00 aerator.  $\frac{3}{8}$ " I.P.S. female thread inlet complete with tank nipple, locknut and washer. **Intended for use as an outlet for W-0825-00.**

## Faucet Accessories





**W-0532-00 Aerator** Anti-splash aerator can be used in lieu of the removable straight serrated hose connectors supplied with the gooseneck water faucet.  $\frac{3}{8}$ " I.P.S. male thread inlet.

**W-0510-00 Filter Pump (Aspirator)**  $5\frac{1}{2}$ " high overall. Has serrated hose connector side and bottom outlets.  $\frac{3}{8}$ " I.P.S. male thread inlet.

**W-0535-00 In-Line Vacuum Breaker** Chrome plated. Mounts between the outlet end of a fitting (gooseneck, flange, etc.) and a removable serrated hose connector.  $\frac{3}{8}$ " I.P.S. male thread inlet and female thread outlet. (Check plumbing codes for required location of vacuum breakers before ordering.)

# Laboratory Service Fittings




## Distilled Water Faucet



**W-0367-00 Right Hand Distilled Water (Curb or Panel Mounted) Faucet** Chrome plated with a hot tin lined interior. It has a self-closing handle and a disc type valve with a renewable internal assembly. The 8" spread rigid gooseneck has a removable serrated hose end. 1/2" I.P.S. flanged female thread inlet complete with tank nipple, locknut and washer.

**W-0368-00 Deck Mounted Distilled Water Gooseneck Faucet** Made of chrome plated brass. Has a hot tin lined interior. It has a self-closing handle and a disc type valve with renewable internal assembly. The 8" spread rigid gooseneck has a removable serrated hose end. Faucet body is complete with an integral threaded mounting shank and a flanged locknut. 1/2" I.P.S. male thread inlet with coupling nut and 1/4" I.P.S. tailpiece.

## Ultra-Pure Deionized Water Faucet




**W-0407-00 Distilled, Deionized and Ultra-Pure Water Faucet (deck mounted)** made of chrome plated brass. All polypropylene internal construction with 6" spread gooseneck. Furnished with 1/2" I.P.S. shank with locknut and washer. Inlet is 3/8" O.D. compression.

**W-0408-00 Distilled, Deionized and Ultra-Pure Water Faucet (panel mounted)** made of chrome plated brass. All polypropylene internal construction with 6" spread gooseneck. Furnished with 1/2" I.P.S. shank with locknut and washer. Inlet is 3/8" O.D. compression.


**W-0409-00 Distilled Deionized and Ultra-Pure Water Faucet (panel mounted)** made of chrome plated brass. All polypropylene internal construction with 8" spread gooseneck. Furnished with 1/2" I.P.S. shank with locknut and washer. Inlet is 3/8" O.D. compression.

# Laboratory Service Fittings


## Vandal-Resistant Fittings




**W-0410-VR Vandal Resistant Deck Mounted Gas, Air or Vacuum Fitting** with single vandal-resistant laboratory ball valve. Complete with a  $\frac{3}{8}$ " I.P.S. tank nipple, locknut and washer. **Specify service desired.**




**W-0411-VR Vandal Resistant Deck Mounted Double Gas, Air or Vacuum Fitting** with two vandal resistant laboratory ball valves with integral straight serrated hose ends at 90°. Complete with a  $\frac{3}{8}$ " I.P.S. tank nipple, locknut and washer. **Specify service desired.**




**W-0412-VR Vandal Resistant Deck Mounted Double Gas, Air or Vacuum Fitting** Same as W-0411-VR except at 180°. Complete with a  $\frac{3}{8}$ " I.P.S. tank nipple, locknut and washer. **Specify service desired.**




**W-0413-VR Vandal Resistant Deck Mounted Hot and Cold Mixing Faucet** with 8" rigid vacuum breaker gooseneck. Valves are disc type with renewable internal assemblies. Faucet body is complete with an integral threaded mounting shank and a flanged locknut.  $\frac{1}{2}$ " I.P.S. male thread inlets with coupling nuts for slip joint connections.




**W-0414-VR Vandal Resistant Panel Mounted Hot and Cold Mixing Faucet** with 8" rigid vacuum breaker gooseneck and mounting shanks.




**W-0415-VR Vandal Resistant Deck Mounted Water Faucet** with 8" rigid vacuum breaker gooseneck. Faucet body is furnished with a  $\frac{3}{8}$ " I.P.S. threaded mounting shank with locknut and washer.



**W-0417-VR Vandal Resistant Deck Mounted Hand-Held Drench Hose** with single spray type outlet heads, squeeze-lever operated valve and flexible stainless steel hose.  $\frac{3}{8}$ " I.P.S. male inlet.




**W-0416-VB Vandal Resistant Deck Mounted Combination Vacuum Breaker Water/Gas Fitting** Made of chrome plated base. Rigid gooseneck with disc type valve and in-line vacuum breaker. Includes two-way turret with two laboratory ball valves. Complete with two  $\frac{3}{8}$ " I.P.S. tank nipples, locknuts and washers. **Specify service desired.**




**W-0416-VR Vandal Resistant Deck Mounted Combination Water/Gas Fitting** Same as W-0416-VB except furnished without in-line vacuum breakers. **Specify service desired.**

# Laboratory Service Fittings


## Emergency Eye Wash/Shower Equipment




**W-0927-00 Panel Mounted Eye Wash Fixture** Consists of stainless steel bowl, two Gentle Spray outlet heads, a "PUSH TO OPERATE" ball valve and mounting bracket. 1/2" I.P.S. female thread inlet and 1 1/4" I.P.S. female thread outlet.




**W-0933-00 Eye/Face Wash (wall mounted)** the same as No. W-0927-00, except with four Gentle Spray outlet heads for complete coverage of eyes and face.




**W-0930-00 BARRIER-FREE (ADA) Eye/Face Wash (wall mounted)** with four Gentle Spray outlet heads, stainless steel skirt and bowl and stay-open ball valve. 1/2" I.P.S. female thread inlet and 1 1/4" I.P.S. female thread outlet.




**W-0931-00 Eye Wash (deck mounted)** with two Gentle Spray outlet heads, stay-open ball valve and 90° swivel feature (right hand mounting).




**W-0932-00 Eye Wash (deck mounted)** the same as No. W-0931-00 but (left hand mounting).




**W-0923-00 Deck Mounted Hose Spray Fitting** Consists of 8 ft. reinforced PVC hose with a squeeze-lever operated valve and spray-type outlet head with protective rubber collar. Flange, shank and locknut for attaching to table top are included. 3/8" I.P.S. male thread swivel inlet.




**W-0934-00 Deck Mounted Dual Purpose Eye Wash/Drench Hose** with two Gentle Spray outlet heads angled at 45° and squeeze valve with locking clip. Flange, shank and locknut for attaching to table top are included. 3/8" I.P.S. male thread swivel inlet.




**W-0935-00 Wall Mounted Dual Purpose Eye Wash/Drench Hose** the same as No. W-0934-00 except wall mounted.



**W-0374-01 Emergency Shower Head Assembly (horizontally mounted).** Includes 10" diameter orange ABS plastic shower head. Operated by a rigid pull-rod and easy to grip handle. Complete with an elbow, connecting pipe and ball valve. 1" I.P.S. female thread inlet.



**W-0928-00 Emergency Shower Head Assembly (vertically mounted).** The same as No. W-0374-01 except vertically mounted.




**W-0926-00 Drench Shower and Eye Wash Combination**  
A necessity for chemical stockrooms and other areas handling a quantity of corrosive chemicals. Free-standing 10" dia. drench shower made of orange ABS plastic. Has a 1" ball valve operated by a rigid pull-rod and easy to grip handle. Eye wash fixture with stainless steel bowl, two Gentle Spray heads and a "PUSH TO OPERATE" ball valve. Eye wash and shower operate independently. Overall height is 80". 1 1/4" I.P.S. female thread inlet and outlet.

**W-0936-00 BARRIER-FREE (ADA) Drench Shower and Eye Wash Combination** Similar to No. W-0926-00 except this unit is barrier-free and eye wash furnished with four Gentle Spray heads.

# Cupsinks & Drain Accessories

## Cupsinks



**0490-BE** Black Epoxy  
**0490-GE** Grey Epoxy  
**0490-PE** Putty Epoxy  
**0490-SE** Slate Epoxy  
 Complete with removable strainer. 6"x 3" inside dimension. 1 1/2" I.P.S. male straight thread outlet. (Overall height is 7 1/2")



**0491-BP** Black Poly  
**0491-GP** Grey Poly  
 Complete with removable strainer. 6"x 3" inside dimension. 1 1/2" I.P.S. male straight thread outlet. (Overall height is 8 1/2")



**0492-BP** Black Poly  
**0492-GP** Grey Poly  
 Complete with removable strainer. 9"x 3" inside dimension. 1 1/2" I.P.S. male straight thread outlet. (Overall height is 7")



**0493-BE** Black Epoxy  
**0493-GE** Grey Epoxy  
**0493-PE** Putty Epoxy  
**0493-SE** Slate Epoxy  
 Complete with removable type 316 stainless steel wire screen. 13 3/8" x 4 1/8" x 5 1/4" I.D. 1 1/2" I.P.S. male straight thread outlet. (Overall height is 8")

## Stainless Steel Cupsinks




**0950-0F** Stainless Steel  
 Made of type 304 stainless steel and has integral cross bars. 5" I.D. 1 1/2" I.P.S. male straight thread outlet.




**0975-0F** Stainless Steel  
 Made of type 316 stainless steel and has integral cross bars. 6"x 3" I.D. 1 1/2" I.P.S. male straight thread outlet.

## Drain Accessories



**0464-00** Polyethylene Trap and Sink Outlet  
 Adjustable riser (tailpiece) can be cut to the desired length. Sink outlet is complete with plug strainer, gasket and locknut. Has a 1 1/2" I.P.S. male straight thread outlet. Running type trap can be used for either "P" or "S" type drain connections. 1 1/2" I.P.S. male straight thread outlet.



**0469-00** Polyethylene Trap and Sink Outlet  
 Adjustable riser (tailpiece) with 4 1/2" adjustment that does not require cutting. The sink outlet is complete with plug, strainer, gasket and locknut. Has a 1 1/2" I.P.S. male straight thread outlet. Trap is the bottle "P" type, with removable bottom for cleanout. 1 1/2" I.P.S. male straight thread outlet.

## Sink Outlet




**0482-BP** Black Poly  
**0482-GP** Grey Poly  
 Complete with locknut. 1 1/2" I.P.S. male straight thread outlet. For use with epoxy resin tub sinks and troughs.

## Sink Overflow



**0494BP-04** 4" High  
**0494BP-06** 6" High  
**0494BP-08** 8" High  
**0494BP-10** 10" High  
 Sink Overflows are designed to be used with 0482-BP sink outlets. 1 1/2" diameter.

## St. Steel Sink Drain



**0995-00** Stainless Steel  
 Crumb Cup Strainer  
 Complete with gaskets, locknut and basket strainer. 1 1/2" male straight thread outlet.

# Kemresin Tub Sinks & Drain Troughs

Kemresin® Tub Sinks and Drain Troughs are made from a combination of modified epoxy resins, upgraded for maximum properties with carefully selected curing resins and inert fillers to produce a thermosetting material. They are formed in metal molds to provide one-piece units with coved corners and bottoms pitched to the drain outlet. After being removed










from the molds, they are cured at high temperatures for several hours to complete the chemical reaction, and ensure maximum chemical, impact and thermal shock resistance. Drain outlets are sized for 1½" sink outlets unless otherwise specified.

Kemresin Drain Troughs are molded with no longitudinal joints. A slope of 1/8" per foot assures complete drainage, and a 7/8" cove radius

in the bottom corners is provided for ease of cleaning. Five trough sections, or any part of them, are available to be cemented together to provide Open End Troughs up to 21' long, Closed End Troughs with an end drain outlet up to 21'-8" long and Closed End Troughs with a center drain outlet up to 42'-8" long. Section TR-4 is provided with a center drain recess for a standard 1½" sink outlet.

## Epoxy Resin Drain Troughs

## Epoxy Resin Tub Sinks



# Kemresin Tub Sinks & Drain Troughs

## Kemresin Colors

Kemresin® Tub Sinks and Drain Troughs are available in four standard colors: Black, Grey, Putty, and Slate.

To specify color, add color code to the end of the part number.  
**Example: 1005-00-BK**

## Color Codes

**BK** = Black      **PT** = Putty  
**GR** = Grey      **SL** = Slate

## Kemresin Tub Sinks

Catalog Number	Inside Dimensions (at Top of Sink)			Drain Location		Overall Sink Height
	Length	Width	Height	A	B	
<b>1000-00-__</b>	14"	10"	6"	4"	5"	6 <sup>7</sup> / <sub>8</sub> "
<b>1001-00-__</b>	16"	8"	7"	4"	4"	8"
<b>1003-00-__</b>	18"	15"	11 <sup>13</sup> / <sub>16</sub> "	4"	4"	12 <sup>7</sup> / <sub>8</sub> "
(1) <b>1003-DT-__</b>	37"	15"	11"	4 <sup>1</sup> / <sub>2</sub> "	4 <sup>1</sup> / <sub>2</sub> "	12"
<b>1005-00-__</b>	25"	15"	10"	4"	4"	11 <sup>1</sup> / <sub>4</sub> "
<b>1006-00-__</b>	18"	15"	8"	4"	4"	9"
(2) <b>1006-DT-__</b>	36 <sup>3</sup> / <sub>4</sub> "	15"	8 <sup>1</sup> / <sub>4</sub> "	4 <sup>1</sup> / <sub>2</sub> "	4 <sup>1</sup> / <sub>2</sub> "	9"
<b>1007-00-__</b>	24"	16"	8"	4"	8"	9 <sup>1</sup> / <sub>8</sub> "
<b>1008-00-__</b>	16"	12"	8"	4"	6"	9 <sup>1</sup> / <sub>8</sub> "
<b>1011-00-__</b>	18"	6 <sup>1</sup> / <sub>2</sub> "	6"	9"	3 <sup>1</sup> / <sub>4</sub> "	7"
(3) <b>1014-00-__</b>	25"	15"	10"	4"	4"	11 <sup>1</sup> / <sub>4</sub> "
(4) <b>1015-00-__</b>	23 <sup>3</sup> / <sub>8</sub> "	13 <sup>3</sup> / <sub>4</sub> "	13 <sup>1</sup> / <sub>4</sub> "	11 <sup>11</sup> / <sub>16</sub> "	9 <sup>1</sup> / <sub>2</sub> "	17 <sup>3</sup> / <sub>4</sub> "
(5) <b>1016-00-__</b>	23 <sup>3</sup> / <sub>8</sub> "	13 <sup>3</sup> / <sub>4</sub> "	13 <sup>1</sup> / <sub>4</sub> "	11 <sup>11</sup> / <sub>16</sub> "	9 <sup>1</sup> / <sub>2</sub> "	17 <sup>3</sup> / <sub>4</sub> "
<b>1021-00-__</b>	16"	16"	7 <sup>1</sup> / <sub>2</sub> "	4"	4"	8 <sup>3</sup> / <sub>4</sub> "
<b>1022-00-__</b>	21"	18"	10"	4 <sup>1</sup> / <sub>2</sub> "	9"	11"
<b>1023-00-__</b>	21 <sup>1</sup> / <sub>2</sub> "	15 <sup>1</sup> / <sub>2</sub> "	7"	4"	4"	8 <sup>3</sup> / <sub>8</sub> "
<b>1024-00-__</b>	21 <sup>1</sup> / <sub>2</sub> "	15 <sup>1</sup> / <sub>2</sub> "	11"	4"	4"	12 <sup>1</sup> / <sub>8</sub> "
(6) <b>1024-DT-__</b>	44"	15 <sup>1</sup> / <sub>2</sub> "	11"	4"	4"	12 <sup>1</sup> / <sub>2</sub> "
<b>1025-00-__</b>	28"	20"	9 <sup>1</sup> / <sub>8</sub> "	14"	6"	10 <sup>1</sup> / <sub>8</sub> "
<b>1026-00-__</b>	35 <sup>1</sup> / <sub>2</sub> "	19 <sup>1</sup> / <sub>2</sub> "	10"	17 <sup>3</sup> / <sub>4</sub> "	9 <sup>3</sup> / <sub>4</sub> "	11 <sup>1</sup> / <sub>8</sub> "
<b>1027-00-__</b>	42"	20"	10"	21"	10"	11 <sup>1</sup> / <sub>4</sub> "
<b>1030-00-__</b>	16"	12"	8"	8"	6"	9"
<b>1031-00-__</b>	24"	16"	8"	12"	8"	9"
<b>1032-00-__</b>	29 <sup>3</sup> / <sub>4</sub> "	14 <sup>3</sup> / <sub>4</sub> "	6"	14 <sup>7</sup> / <sub>8</sub> "	7 <sup>3</sup> / <sub>8</sub> "	7"
<b>1034-00-__</b>	21 <sup>1</sup> / <sub>2</sub> "	19 <sup>1</sup> / <sub>4</sub> "	8"	10 <sup>3</sup> / <sub>4</sub> "	9 <sup>5</sup> / <sub>8</sub> "	9 <sup>1</sup> / <sub>8</sub> "
<b>1035-00-__</b>	45 <sup>3</sup> / <sub>4</sub> "	17 <sup>3</sup> / <sub>8</sub> "	12 <sup>1</sup> / <sub>2</sub> "	22 <sup>7</sup> / <sub>8</sub> "	8 <sup>3</sup> / <sub>4</sub> "	13 <sup>3</sup> / <sub>8</sub> "
<b>1036-00-__</b>	48 <sup>5</sup> / <sub>8</sub> "	24 <sup>1</sup> / <sub>2</sub> "	5 <sup>3</sup> / <sub>8</sub> "	3 <sup>3</sup> / <sub>4</sub> "	12 <sup>1</sup> / <sub>4</sub> "	6"
(7) <b>1037-00-__</b>	23 <sup>1</sup> / <sub>4</sub> "	14 <sup>1</sup> / <sub>4</sub> "	12 <sup>3</sup> / <sub>4</sub> "	11 <sup>5</sup> / <sub>8</sub> "	9 <sup>1</sup> / <sub>4</sub> "	22 <sup>3</sup> / <sub>8</sub> "
(8) <b>1038-00-__</b>	12" diameter		7 <sup>5</sup> / <sub>8</sub> "	center		8 <sup>3</sup> / <sub>8</sub> "
(9) <b>1039-00-__</b>	18" diameter		8 <sup>1</sup> / <sub>8</sub> "	center		8 <sup>7</sup> / <sub>8</sub> "
<b>SM-01-00-__</b>	12"	6"	6"	6"	3"	7"
<b>SM-02-00-__</b>	19 <sup>1</sup> / <sub>2</sub> "	16 <sup>1</sup> / <sub>2</sub> "	8"	9 <sup>3</sup> / <sub>4</sub> "	8 <sup>1</sup> / <sub>4</sub> "	9"
<b>SM-03-00-__</b>	24"	21"	8"	6"	4 <sup>1</sup> / <sub>2</sub> "	9 <sup>1</sup> / <sub>8</sub> "
<b>SM-04-00-__</b>	25 <sup>1</sup> / <sub>2</sub> "	16 <sup>1</sup> / <sub>2</sub> "	19"	4 <sup>1</sup> / <sub>2</sub> "	4 <sup>1</sup> / <sub>2</sub> "	20 <sup>1</sup> / <sub>8</sub> "
<b>SM-05-00-__</b>	28"	18"	10"	4"	4"	11 <sup>1</sup> / <sub>8</sub> "
<b>SM-06-00-__</b>	28"	18"	10"	24"	4"	11 <sup>1</sup> / <sub>8</sub> "
<b>SM-07-00-__</b>	30"	16"	10"	4"	4"	11 <sup>1</sup> / <sub>8</sub> "
<b>SM-08-00-__</b>	30"	16"	10"	26"	4"	11 <sup>1</sup> / <sub>8</sub> "
<b>SM-09-00-__</b>	30"	18"	18"	4 <sup>1</sup> / <sub>2</sub> "	4 <sup>3</sup> / <sub>4</sub> "	18 <sup>3</sup> / <sub>4</sub> "
<b>SM-10-00-__</b>	36"	21"	4 <sup>3</sup> / <sub>4</sub> "	4"	4"	5 <sup>3</sup> / <sub>4</sub> "
<b>SM-11-00-__</b>	48"	13"	12 <sup>1</sup> / <sub>4</sub> "	24"	6 <sup>1</sup> / <sub>2</sub> "	13 <sup>1</sup> / <sub>4</sub> "
<b>SM-12-00-__</b>	8"	8"	8"	4"	4"	8 <sup>7</sup> / <sub>8</sub> "

- (1) Double tub sink. Distance between centerlines of drain outlets is 19<sup>3</sup>/<sub>8</sub>".
- (2) Double tub sink. Distance between centerlines of drain outlets is 19".
- (3) Sink has 7<sup>1</sup>/<sub>8</sub>" long by 5<sup>5</sup>/<sub>8</sub>" high trough cutout.
- (4) Center table end sink has 7<sup>1</sup>/<sub>8</sub>" long by 5<sup>5</sup>/<sub>8</sub>" high trough cutout and 24<sup>1</sup>/<sub>2</sub>" long by 3<sup>3</sup>/<sub>4</sub>" high back splash. Table side of sink is not tapered.
- (5) Center table end sink has 24<sup>1</sup>/<sub>2</sub>" long by 3<sup>3</sup>/<sub>4</sub>" high back splash without trough cutout. Table side of sink is not tapered.
- (6) Double tub sink. Distance between centerlines of drain outlets is 22<sup>1</sup>/<sub>2</sub>".
- (7) Assembly of tub sink, back splash, ledge and end fillers. Used for T69W4525-00 Scrub Sink Assembly.
- (8) Cylindrical sink –12" inside diameter at top of sink with center drain outlet.
- (9) Cylindrical sink –18" inside diameter at top of sink with center drain outlet.
- (10) Sides have 1<sup>1</sup>/<sub>2</sub>" taper per foot in height.
- (11) Exterior of one 21" long side is not tapered.

Catalog Number	Inside Dimensions (at Top of Sink)			Drain Location		Overall Sink Height
	Length	Width	Height	A	B	
<b>SM-13-00-__</b>	14"	10"	8"	4"	5"	9"
<b>SM-14-00-__</b>	23 <sup>1</sup> / <sub>4</sub> "	14 <sup>1</sup> / <sub>4</sub> "	12 <sup>3</sup> / <sub>4</sub> "	11 <sup>5</sup> / <sub>8</sub> "	9 <sup>1</sup> / <sub>4</sub> "	13 <sup>5</sup> / <sub>8</sub> "
<b>SM-15-00-__</b>	19 <sup>3</sup> / <sub>4</sub> "	18"	24"	9 <sup>7</sup> / <sub>8</sub> "	9"	25"
<b>SM-16-00-__</b>	40"	18"	18 <sup>1</sup> / <sub>4</sub> "	20"	9"	19 <sup>1</sup> / <sub>4</sub> "
<b>SM-17-00-__</b>	29 <sup>3</sup> / <sub>4</sub> "	23 <sup>1</sup> / <sub>2</sub> "	6"	14 <sup>7</sup> / <sub>8</sub> "	11 <sup>3</sup> / <sub>4</sub> "	7"
<b>SM-18-00-__</b>	18"	15"	15 <sup>1</sup> / <sub>4</sub> "	4 <sup>3</sup> / <sub>4</sub> "	4 <sup>3</sup> / <sub>4</sub> "	16 <sup>1</sup> / <sub>2</sub> "
<b>SM-19-00-__</b>	35 <sup>1</sup> / <sub>2</sub> "	19 <sup>1</sup> / <sub>2</sub> "	15 <sup>3</sup> / <sub>8</sub> "	3 <sup>1</sup> / <sub>2</sub> "	4 <sup>1</sup> / <sub>2</sub> "	16 <sup>1</sup> / <sub>2</sub> "
<b>SM-20-00-__</b>	21"	18"	15 <sup>1</sup> / <sub>4</sub> "	4 <sup>1</sup> / <sub>2</sub> "	4 <sup>1</sup> / <sub>2</sub> "	16 <sup>1</sup> / <sub>2</sub> "
<b>SM-21-00-__</b>	60 <sup>1</sup> / <sub>2</sub> "	22"	7 <sup>1</sup> / <sub>4</sub> "	6"	6"	8 <sup>1</sup> / <sub>8</sub> "
<b>SM-22-00-__</b>	24"	16"	15 <sup>3</sup> / <sub>8</sub> "	12"	8"	16 <sup>1</sup> / <sub>8</sub> "
<b>SM-23-00-__</b>	28"	18"	15 <sup>1</sup> / <sub>2</sub> "	5 <sup>1</sup> / <sub>2</sub> "	4"	16 <sup>1</sup> / <sub>2</sub> "
<b>SM-24-00-__</b>	48"	8"	8"	24"	4"	9"
<b>SM-25-00-__</b>	21"	18"	15 <sup>1</sup> / <sub>8</sub> "	10 <sup>1</sup> / <sub>2</sub> "	9"	16 <sup>1</sup> / <sub>8</sub> "
<b>SM-26-00-__</b>	21"	11 <sup>3</sup> / <sub>8</sub> "	15 <sup>1</sup> / <sub>8</sub> "	10 <sup>1</sup> / <sub>2</sub> "	5 <sup>1</sup> / <sub>16</sub> "	16 <sup>1</sup> / <sub>8</sub> "
<b>SM-27-00-__</b>	36"	20"	16 <sup>1</sup> / <sub>4</sub> "	18"	10"	17 <sup>1</sup> / <sub>4</sub> "
<b>SM-28-00-__</b>	30"	18"	4 <sup>1</sup> / <sub>8</sub> "	15"	9"	5 <sup>1</sup> / <sub>4</sub> "
<b>SM-29-00-__</b>	36"	24"	12 <sup>1</sup> / <sub>4</sub> "	18"	12"	13 <sup>1</sup> / <sub>4</sub> "
<b>SM-30-00-__</b>	40"	30"	8 <sup>1</sup> / <sub>8</sub> "	20"	15"	9"
<b>SM-31-00-__</b>	24"	15 <sup>3</sup> / <sub>4</sub> "	12 <sup>1</sup> / <sub>8</sub> "	4 <sup>3</sup> / <sub>4</sub> "	8"	13"
<b>SM-32-00-__</b>	59"	21"	12 <sup>5</sup> / <sub>8</sub> "	29 <sup>1</sup> / <sub>2</sub> "	10 <sup>1</sup> / <sub>2</sub> "	13 <sup>3</sup> / <sub>8</sub> "
<b>SM-33-00-__</b>	24"	18"	11"	4"	9"	11 <sup>7</sup> / <sub>8</sub> "
(10) <b>SM-34-00-__</b>	15"	14"	10 <sup>1</sup> / <sub>2</sub> "	7 <sup>1</sup> / <sub>2</sub> "	7"	11 <sup>1</sup> / <sub>4</sub> "
<b>SM-35-00-__</b>	51"	17"	12 <sup>1</sup> / <sub>2</sub> "	25 <sup>1</sup> / <sub>2</sub> "	8 <sup>1</sup> / <sub>2</sub> "	13 <sup>1</sup> / <sub>4</sub> "
<b>SM-36-00-__</b>	14"	10"	6"	7"	5"	6 <sup>7</sup> / <sub>8</sub> "
(11) <b>SM-37-00-__</b>	24"	20"	10 <sup>5</sup> / <sub>8</sub> "	6 <sup>3</sup> / <sub>4</sub> "	10"	11 <sup>1</sup> / <sub>4</sub> "
<b>SM-39-00-__</b>	35"	6 <sup>5</sup> / <sub>8</sub> "	4 <sup>3</sup> / <sub>4</sub> "	17 <sup>1</sup> / <sub>2</sub> "	3 <sup>5</sup> / <sub>16</sub> "	5 <sup>1</sup> / <sub>2</sub> "
<b>SM-40-00-__</b>	33 <sup>3</sup> / <sub>8</sub> "	29 <sup>1</sup> / <sub>8</sub> "	15 <sup>3</sup> / <sub>4</sub> "	6"	14 <sup>9</sup> / <sub>16</sub> "	16 <sup>1</sup> / <sub>2</sub> "
<b>SM-41-00-__</b>	18"	15"	8 <sup>1</sup> / <sub>4</sub> "	9"	7 <sup>1</sup> / <sub>2</sub> "	9"
<b>SM-42-00-__</b>	25"	15"	10 <sup>1</sup> / <sub>2</sub> "	12 <sup>1</sup> / <sub>2</sub> "	7 <sup>1</sup> / <sub>2</sub> "	11 <sup>1</sup> / <sub>4</sub> "
<b>SM-43-00-__</b>	21"	18"	10 <sup>3</sup> / <sub>8</sub> "	10 <sup>1</sup> / <sub>2</sub> "	9"	11 <sup>1</sup> / <sub>8</sub> "
<b>SM-44-00-__</b>	36"	16"	11"	12"	4"	12"
<b>SM-45-00-__</b>	27 <sup>3</sup> / <sub>4</sub> "	14 <sup>3</sup> / <sub>4</sub> "	8 <sup>1</sup> / <sub>2</sub> "	10 <sup>1</sup> / <sub>2</sub> "	7 <sup>1</sup> / <sub>2</sub> "	9 <sup>3</sup> / <sub>8</sub> "
<b>SM-47-00-__</b>	43"	15 <sup>1</sup> / <sub>2</sub> "	11 <sup>1</sup> / <sub>8</sub> "	21 <sup>1</sup> / <sub>2</sub> "	7 <sup>3</sup> / <sub>4</sub> "	12"
<b>SM-48-00-__</b>	56"	24"	8 <sup>1</sup> / <sub>4</sub> "	28"	12"	9"
<b>SM-49-00-__</b>	36"	16"	8 <sup>3</sup> / <sub>8</sub> "	6"	8"	9"
<b>SM-50-00-__</b>	24"	24"	12 <sup>1</sup> / <sub>2</sub> "	12"	12"	13 <sup>1</sup> / <sub>4</sub> "
<b>SM-51-00-__</b>	30"	15 <sup>7</sup> / <sub>8</sub> "	16 <sup>3</sup> / <sub>8</sub> "	6"	7 <sup>15</sup> / <sub>16</sub> "	17 <sup>1</sup> / <sub>8</sub> "
<b>SM-52-00-__</b>	45 <sup>1</sup> / <sub>8</sub> "	11 <sup>1</sup> / <sub>2</sub> "	8 <sup>1</sup> / <sub>8</sub> "	22 <sup>9</sup> / <sub>16</sub> "	5 <sup>3</sup> / <sub>4</sub> "	9"
<b>SM-53-00-__</b>	31 <sup>1</sup> / <sub>2</sub> "	16"	16 <sup>1</sup> / <sub>2</sub> "	7 <sup>1</sup> / <sub>8</sub> "	8"	17 <sup>1</sup> / <sub>4</sub> "
<b>SM-55-00-__</b>	16"	16"	12 <sup>3</sup> / <sub>8</sub> "	4 <sup>1</sup> / <sub>4</sub> "	4 <sup>1</sup> / <sub>4</sub> "	13"
<b>SM-56-00-__</b>	42 <sup>1</sup> / <sub>4</sub> "	18 <sup>1</sup> / <sub>8</sub> "	12 <sup>5</sup> / <sub>8</sub> "	21 <sup>1</sup> / <sub>8</sub> "	9 <sup>1</sup> / <sub>16</sub> "	13 <sup>3</sup> / <sub>8</sub> "

### General Notes:

- Sinks with center drain outlets, except 1011-00, can be furnished with a drain recess sized for a crumb cup strainer or garbage disposal when specified.
- Any sink may be furnished with trough cutout. Specify location and size. (Trough location must be at least 1<sup>1</sup>/<sub>2</sub>" from inside corner and 2" from bottom.)
- Most sinks can be joined together to form multiple tub sinks provided adjoining sides of sinks are the same size. Sinks longer than 30", wider than 24" or deeper than 12" cannot be factory joined, but they can be factory prepared for field assembly. Assemblies of three or more tub sinks should be no longer than 60" overall. For longer assemblies, one factory joint can be made with a second joint prepared for field assembly.





# Kemresin Drop-in Sinks

Kemresin Drop-in Sinks are formed in metal molds to provide a lipped one-piece tub with covered corners and bottoms pitched to the drain outlet. They are made from a combination of modified epoxy resins, upgraded for maximum properties with carefully selected curing resins and inert fillers to produce a thermosetting


material. After being removed from the molds, they are cured at high temperatures for several hours to complete the chemical reaction, and ensure maximum chemical, impact and thermal shock resistance. Drain outlets are sized for 1/2" sink outlets unless otherwise specified.

They are designed to be mounted in Kemresin tops in a specially rabbetted cutout so that no support from below is required.

## Drop-in Kemresin Sinks



PLAN VIEW  
**1048-DI-**



SECTION VIEW  
**1049-DI-**

Part Number	I.D.		Depth	O.D.		Drain Location	
	A	B	C	D	E	G	H
<b>1000-DI-</b> ___	13 <sup>7</sup> / <sub>8</sub> "	9 <sup>7</sup> / <sub>8</sub> "	6"	16 <sup>1</sup> / <sub>8</sub> "	12 <sup>1</sup> / <sub>8</sub> "	7"	5"
<b>1001-DI-</b> ___	8"	16"	7"	10 <sup>1</sup> / <sub>4</sub> "	18 <sup>1</sup> / <sub>4</sub> "	4"	4"
<b>1003-DI-</b> ___	17 <sup>15</sup> / <sub>16</sub> "	14 <sup>15</sup> / <sub>16</sub> "	11"	20 <sup>3</sup> / <sub>16</sub> "	17 <sup>3</sup> / <sub>16</sub> "	3 <sup>1</sup> / <sub>2</sub> "	3 <sup>1</sup> / <sub>2</sub> "
<b>1005-DI-</b> ___	25"	15"	10"	27 <sup>1</sup> / <sub>4</sub> "	17 <sup>1</sup> / <sub>4</sub> "	4"	4"
<b>1006-DI-</b> ___	17 <sup>7</sup> / <sub>8</sub> "	14 <sup>7</sup> / <sub>8</sub> "	8"	20 <sup>1</sup> / <sub>8</sub> "	17 <sup>1</sup> / <sub>8</sub> "	3 <sup>1</sup> / <sub>2</sub> "	3 <sup>1</sup> / <sub>2</sub> "
<b>1008-DI-</b> ___	16"	12"	8"	18 <sup>1</sup> / <sub>4</sub> "	14 <sup>1</sup> / <sub>4</sub> "	4"	6"
<b>1011-DI-</b> ___	18"	6 <sup>1</sup> / <sub>2</sub> "	6"	20 <sup>1</sup> / <sub>4</sub> "	8 <sup>3</sup> / <sub>4</sub> "	9"	3 <sup>1</sup> / <sub>4</sub> "
<b>1021-DI-</b> ___	16"	16"	7 <sup>1</sup> / <sub>2</sub> "	18 <sup>1</sup> / <sub>4</sub> "	18 <sup>1</sup> / <sub>4</sub> "	4"	4"
<b>1022-DI-</b> ___	21"	18"	10"	23 <sup>1</sup> / <sub>4</sub> "	20 <sup>1</sup> / <sub>4</sub> "	4"	9"
(1) <b>1038-DI-</b> ___	12" diameter		7 <sup>5</sup> / <sub>8</sub> "	14 <sup>1</sup> / <sub>4</sub> " diameter		center center	
(2) <b>1039-DI-</b> ___	18" diameter		8 <sup>1</sup> / <sub>8</sub> "	20 <sup>1</sup> / <sub>4</sub> " diameter		center center	
<b>1040-DI-</b> ___	15 <sup>3</sup> / <sub>8</sub> "	11 <sup>3</sup> / <sub>8</sub> "	8"	17.6"	13.6"	4"	4"
<b>1041-DI-</b> ___	23 <sup>3</sup> / <sub>8</sub> "	15 <sup>3</sup> / <sub>8</sub> "	8.3"	25.6"	17.6"	4"	4"
<b>1042-DI-</b> ___	28"	15"	11.8"	30 <sup>1</sup> / <sub>4</sub> "	17 <sup>1</sup> / <sub>4</sub> "	5"	5"
<b>1043-DI-</b> ___	30"	16"	17.8"	32 <sup>1</sup> / <sub>4</sub> "	18 <sup>1</sup> / <sub>4</sub> "	5"	8"
<b>1044-DI-</b> ___	12"	8"	6"	14 <sup>1</sup> / <sub>4</sub> "	10 <sup>1</sup> / <sub>4</sub> "	6"	4"
<b>1045-DI-</b> ___	21"	17"	10"	23 <sup>1</sup> / <sub>4</sub> "	19 <sup>1</sup> / <sub>4</sub> "	4 <sup>1</sup> / <sub>2</sub> "	8 <sup>1</sup> / <sub>2</sub> "
<b>1046-DI-</b> ___	21 <sup>1</sup> / <sub>2</sub> "	15 <sup>1</sup> / <sub>2</sub> "	11"	23 <sup>3</sup> / <sub>4</sub> "	17 <sup>3</sup> / <sub>4</sub> "	4"	4"
<b>1047-DI-</b> ___	24"	18"	11"	26 <sup>1</sup> / <sub>4</sub> "	20 <sup>1</sup> / <sub>4</sub> "	4"	4"
(3) <b>1048-DI-</b> ___	17 <sup>1</sup> / <sub>2</sub> "	12 <sup>3</sup> / <sub>8</sub> "	6 <sup>1</sup> / <sub>2</sub> "	19 <sup>1</sup> / <sub>2</sub> "	14 <sup>5</sup> / <sub>8</sub> "	8 <sup>3</sup> / <sub>4</sub> "	6 <sup>3</sup> / <sub>16</sub> "
(4) <b>1049-DI-</b> ___	18"	15"	11"	20 <sup>1</sup> / <sub>4</sub> "	17 <sup>1</sup> / <sub>4</sub> "	4 <sup>1</sup> / <sub>2</sub> "	5"
<b>1070-DI-</b> ___	24"	16"	15 <sup>1</sup> / <sub>2</sub> "	26 <sup>1</sup> / <sub>4</sub> "	18 <sup>1</sup> / <sub>4</sub> "	12"	8"
ADA Compliant							
<b>1000-ADA-DI-</b> ___	13 <sup>7</sup> / <sub>8</sub> "	9 <sup>7</sup> / <sub>8</sub> "	5"	16 <sup>1</sup> / <sub>8</sub> "	12 <sup>1</sup> / <sub>8</sub> "	4"	4"
<b>1005-ADA-DI-</b> ___	25"	15"	4.8"	27 <sup>1</sup> / <sub>4</sub> "	17 <sup>1</sup> / <sub>4</sub> "	4"	4"
<b>1006-ADA-DI-</b> ___	17 <sup>7</sup> / <sub>8</sub> "	14 <sup>7</sup> / <sub>8</sub> "	5"	20 <sup>1</sup> / <sub>8</sub> "	17 <sup>1</sup> / <sub>8</sub> "	3 <sup>1</sup> / <sub>2</sub> "	3 <sup>1</sup> / <sub>2</sub> "

(1) Cylindrical sink, 12" inside diameter  
(2) Cylindrical sink, 18" inside diameter

(3) Triangular sink used in X92 and X94 assemblies  
(4) Stepped sink, 4<sup>3</sup>/<sub>4</sub>" deep at front

## Kemresin Colors

Kemresin® Drop-in Sinks are available in four standard colors: Black, Grey, Putty, and Slate.

To specify color, add color code to the end of the part number.  
**Example: 1005-DI-BK**

## Color Codes

**BK** = Black      **PT** = Putty  
**GR** = Grey      **SL** = Slate



# Research Collection

## Stainless Steel Drop-in Sinks


Seamless die-drawn from nickel bearing type 304, 18-8, 18 gauge stainless steel. Allows easy cleaning and stain resistance. Interior and top surfaces polished to a gleaming finish with a highlighted bowl rim. Corners and bottom turns are radiused  $1\frac{3}{4}$ ". Sinks are thoroughly sound deadened with environmentally

safe compounds, to ensure proper sound insulation and to reduce condensation.

All sinks are self-rimming with stainless steel mounting channels for easy, top mount, snap-in clip installation. Sinks are compatible with any type work-surface material from  $\frac{3}{4}$ " to  $2\frac{1}{2}$ " thick.


Each sink is center punched to receive a standard cup strainer. ADA compliant sinks are punched in the center,  $4\frac{1}{2}$ " from the rear of the bowl for extra clearance.

### Stainless Steel Sinks



Part Number	O.D.		I.D.		Sink Depth
	Front to Back	Side to Side	Front to Back	Side to Side	
<b>F-5353-08</b>	15"	15"	12"	12"	7 $\frac{1}{2}$ "
<b>F-5354-08</b>	17"	19"	14"	16"	7 $\frac{1}{2}$ "
<b>F-5355-08</b>	19"	22"	16"	19"	7 $\frac{1}{2}$ "
<b>F-5356-08</b>	19"	25"	16"	22"	7 $\frac{1}{2}$ "
<b>F-5357-08</b>	19"	31"	16"	28"	7 $\frac{1}{2}$ "
<b>F-5353-09</b>	15"	15"	12"	12"	9"
<b>F-5354-11</b>	17"	19"	14"	16"	10 $\frac{1}{2}$ "
<b>F-5355-11</b>	19"	22"	16"	19"	10 $\frac{1}{2}$ "
<b>F-5356-11</b>	19"	25"	16"	22"	10 $\frac{1}{2}$ "
<b>F-5357-11</b>	19"	31"	16"	28"	10 $\frac{1}{2}$ "
<b>F-5354-12</b>	17"	19"	14"	16"	12"
<b>F-5355-12</b>	19"	22"	16"	19"	12"
<b>F-5356-12</b>	19"	25"	16"	22"	12"
<b>F-5357-12</b>	19"	31"	16"	28"	12"


### Stainless Steel Sinks with Fixture Ledge



Part Number	O.D.		I.D.		Sink Depth	Fixture Centers
	Front to Back	Side to Side	Front to Back	Side to Side		
<b>F-5358-08</b>	18"	15"	12"	12"	7 $\frac{1}{2}$ "	4"
<b>F-5359-08</b>	20"	19"	14"	16"	7 $\frac{1}{2}$ "	6"
<b>F-5360-08</b>	22"	22"	16"	19"	7 $\frac{1}{2}$ "	6"
<b>F-5361-08</b>	22"	25"	16"	22"	8"	8"
<b>F-5362-08</b>	22"	31"	16"	28"	7 $\frac{1}{2}$ "	8"
<b>F-5358-09</b>	18"	15"	12"	12"	9"	4"
<b>F-5359-11</b>	20"	19"	14"	16"	10 $\frac{1}{2}$ "	6"
<b>F-5360-11</b>	22"	22"	16"	19"	10 $\frac{1}{2}$ "	6"
<b>F-5361-11</b>	22"	25"	16"	22"	10 $\frac{1}{2}$ "	8"
<b>F-5362-11</b>	22"	31"	16"	28"	10 $\frac{1}{2}$ "	8"
<b>F-5360-12</b>	22"	22"	16"	19"	12"	6"
<b>F-5361-12</b>	22"	25"	16"	22"	12"	8"
<b>F-5362-12</b>	22"	31"	16"	28"	12"	8"


### ADA Compliant

Part Number	O.D.		I.D.		Sink Depth
	Front to Back	Side to Side	Front to Back	Side to Side	
<b>F-5353-06</b>	15"	15"	12"	12"	5 $\frac{1}{2}$ "
<b>F-5354-06</b>	17"	19"	14"	16"	5 $\frac{1}{2}$ "
<b>F-5355-06</b>	19"	22"	16"	19"	5 $\frac{1}{2}$ "
<b>F-5356-06</b>	19"	25"	16"	22"	5 $\frac{1}{2}$ "
<b>F-5357-06</b>	19"	31"	16"	28"	5 $\frac{1}{2}$ "



### ADA Compliant with Fixture Ledge





Part Number	O.D.		I.D.		Sink Depth	Fixture Centers
	Front to Back	Side to Side	Front to Back	Side to Side		
<b>F-5358-06</b>	18"	15"	12"	12"	5 $\frac{1}{2}$ "	4"
<b>F-5359-06</b>	20"	19"	14"	16"	5 $\frac{1}{2}$ "	6"
<b>F-5360-06</b>	22"	22"	16"	19"	5 $\frac{1}{2}$ "	6"
<b>F-5361-06</b>	22"	25"	16"	22"	5 $\frac{1}{2}$ "	8"
<b>F-5362-06</b>	22"	31"	16"	28"	5 $\frac{1}{2}$ "	8"










# Laboratory Electrical Fixtures

## Electrical Receptacles

All flush mount fixtures furnished complete with device, stainless steel face plate and electrical box. Pedestals are furnished with device, wrinkled black epoxy powder coated pedestal box, stainless steel face plate, tank nipple and locknut. All devices and boxes are UL listed. (Ship loose)









Type 2 and Type 3 pedestal boxes furnished with 3/4" IPS tank nipple, Type 7 and Type 8 pedestal boxes furnished with 1" IPS tank nipple.		Flush mount single gang box	Single Wide Pedestal Type 2 – Single Faced Type 3 – Double Faced			Flush mount double gang box	Double Wide Pedestal Type 7 – Single Faced Type 8 – Double Faced	
Electrical Device	Color	Type 1	Type 2	Type 3	Type 6	Type 7	Type 8	
<b>Specification Grade</b> 120 Volt AC, 20 amp, 3-wire polarized duplex receptacle NEMA 5-20 	Black	<b>S581-1K</b>	<b>S581-2K</b>	<b>S581-3K</b>	<b>S581-6K</b>	<b>S581-7K</b>	<b>S581-8K</b>	
	Ivory	<b>S581-1V</b>	<b>S581-2V</b>	<b>S581-3V</b>	<b>S581-6V</b>	<b>S581-7V</b>	<b>S581-8V</b>	
	White	<b>S581-1W</b>	<b>S581-2W</b>	<b>S581-3W</b>	<b>S581-6W</b>	<b>S581-7W</b>	<b>S581-8W</b>	
	Grey	<b>S581-1G</b>	<b>S581-2G</b>	<b>S581-3G</b>	<b>S581-6G</b>	<b>S581-7G</b>	<b>S581-8G</b>	
	Red	<b>S581-1R</b>	<b>S581-2R</b>	<b>S581-3R</b>	<b>S581-6R</b>	<b>S581-7R</b>	<b>S581-8R</b>	
<b>Specification Grade</b> 120 Volt AC, 20 amp, 3-wire polarized duplex receptacle NEMA 5-20 	Black	<b>S681-1K</b>	<b>S681-2K</b>	<b>S681-3K</b>	<b>S681-6K</b>	<b>S681-7K</b>	<b>S681-8K</b>	
	Ivory	<b>S681-1V</b>	<b>S681-2V</b>	<b>S681-3V</b>	<b>S681-6V</b>	<b>S681-7V</b>	<b>S681-8V</b>	
	White	<b>S681-1W</b>	<b>S681-2W</b>	<b>S681-3W</b>	<b>S681-6W</b>	<b>S681-7W</b>	<b>S681-8W</b>	
	Grey	<b>S681-1G</b>	<b>S681-2G</b>	<b>S681-3G</b>	<b>S681-6G</b>	<b>S681-7G</b>	<b>S681-8G</b>	
	Red	<b>S681-1R</b>	<b>S681-2R</b>	<b>S681-3R</b>	<b>S681-6R</b>	<b>S681-7R</b>	<b>S681-8R</b>	
<b>Hospital Grade</b> 120 Volt AC, 20 amp, 3-wire polarized duplex receptacle NEMA 5-20 	Black	<b>H581-1K</b>	<b>H581-2K</b>	<b>H581-3K</b>	<b>H581-6K</b>	<b>H581-7K</b>	<b>H581-8K</b>	
	Ivory	<b>H581-1V</b>	<b>H581-2V</b>	<b>H581-3V</b>	<b>H581-6V</b>	<b>H581-7V</b>	<b>H581-8V</b>	
	White	<b>H581-1W</b>	<b>H581-2W</b>	<b>H581-3W</b>	<b>H581-6W</b>	<b>H581-7W</b>	<b>H581-8W</b>	
	Grey	<b>H581-1G</b>	<b>H581-2G</b>	<b>H581-3G</b>	<b>H581-6G</b>	<b>H581-7G</b>	<b>H581-8G</b>	
	Red	<b>H581-1R</b>	<b>H581-2R</b>	<b>H581-3R</b>	<b>H581-6R</b>	<b>H581-7R</b>	<b>H581-8R</b>	
<b>Specification Grade</b> 240 Volt AC, 20 amp, 3-wire polarized duplex receptacle NEMA 6-20 	Black	<b>S582-1K</b>	<b>S582-2K</b>	<b>S582-3K</b>	<b>S582-6K</b>	<b>S582-7K</b>	<b>S582-8K</b>	
	Ivory	<b>S582-1V</b>	<b>S582-2V</b>	<b>S582-3V</b>	<b>S582-6V</b>	<b>S582-7V</b>	<b>S582-8V</b>	
	White	<b>S582-1W</b>	<b>S582-2W</b>	<b>S582-3W</b>	<b>S582-6W</b>	<b>S582-7W</b>	<b>S582-8W</b>	
	Grey	<b>S582-1G</b>	<b>S582-2G</b>	<b>S582-3G</b>	<b>S582-6G</b>	<b>S582-7G</b>	<b>S582-8G</b>	
	Red	<b>S582-1R</b>	<b>S582-2R</b>	<b>S582-3R</b>	<b>S582-6R</b>	<b>S582-7R</b>	<b>S582-8R</b>	
<b>Ground Fault Protected</b> 120 Volt AC, 20 amp, specification grade, duplex receptacle NEMA 5-20 	Black	<b>S656-1K</b>	<b>S656-2K</b>	<b>S656-3K</b>				
	Ivory	<b>S656-1V</b>	<b>S656-2V</b>	<b>S656-3V</b>				
	White	<b>S656-1W</b>	<b>S656-2W</b>	<b>S656-3W</b>				
	Grey	<b>S656-1G</b>	<b>S656-2G</b>	<b>S656-3G</b>				
	Red	<b>S656-1R</b>	<b>S656-2R</b>	<b>S656-3R</b>				
<b>Ground Fault Protected</b> 120 Volt AC, 20 amp, duplex receptacle and 3-wire polarized duplex receptacle NEMA 5-20 	Black				<b>S657-6K</b>	<b>S657-7K</b>	<b>S657-8K</b>	
	Ivory				<b>S657-6V</b>	<b>S657-7V</b>	<b>S657-8V</b>	
	White				<b>S657-6W</b>	<b>S657-7W</b>	<b>S657-8W</b>	
	Grey				<b>S657-6G</b>	<b>S657-7G</b>	<b>S657-8G</b>	
	Red				<b>S657-6R</b>	<b>S657-7R</b>	<b>S657-8R</b>	
<b>Ground Fault Protected</b> 120 Volt AC, 20 amp, duplex receptacle and 3-wire polarized duplex receptacle NEMA 5-20 	Black				<b>S658-6K</b>	<b>S658-7K</b>	<b>S658-8K</b>	
	Ivory				<b>S658-6V</b>	<b>S658-7V</b>	<b>S658-8V</b>	
	White				<b>S658-6W</b>	<b>S658-7W</b>	<b>S658-8W</b>	
	Grey				<b>S658-6G</b>	<b>S658-7G</b>	<b>S658-8G</b>	
	Red				<b>S658-6R</b>	<b>S658-7R</b>	<b>S658-8R</b>	

# Laboratory Electrical Fixtures

## Electrical Switches and Data Ports

All flush mount fixtures furnished complete with device, stainless steel face plate and electrical box. All devices and boxes are UL listed. (Ship loose)



Electrical Device		Color	Flush mount single gang box		Flush mount double gang box	
			Cat-6	Cat-6a	Cat-6	Cat-6a
<b>Specification Grade</b> Single pole, 120/240 Volt AC, 20 amp toggle switch. 	Black Ivory White Grey Red	<b>S655-1K</b> <b>S655-1V</b> <b>S655-1W</b> <b>S655-1G</b> <b>S655-1R</b>				
<b>RJ-45 Data Port Duplex</b> Cat-6 Data Port 	Black Ivory White Grey Red		<b>0720-1K-6</b> <b>0720-1V-6</b> <b>0720-1W-6</b> <b>0720-1G-6</b> <b>0720-1R-6</b>	<b>0720-1K-6a</b> <b>0720-1V-6a</b> <b>0720-1W-6a</b> <b>0720-1G-6a</b> <b>0720-1R-6a</b>		
<b>RJ-45 Data Port Duplex</b> Cat-6 Data Port with 120 Volt AC, 20 amp, specification grade, duplex receptacle. NEMA 5-20 	Black Ivory White Grey Red				<b>S721-6K-6</b> <b>S721-6V-6</b> <b>S721-6W-6</b> <b>S721-6G-6</b> <b>S721-6R-6</b>	<b>S721-6K-6a</b> <b>S721-6V-6a</b> <b>S721-6W-6a</b> <b>S721-6G-6a</b> <b>S721-6R-6a</b>
<b>RJ-45 Data Port Duplex</b> Cat-6 Data Port with 120 Volt AC, 20 amp, specification grade, GFI duplex receptacle. NEMA 5-20 	Black Ivory White Grey Red				<b>S751-6K-6</b> <b>S751-6V-6</b> <b>S751-6W-6</b> <b>S751-6G-6</b> <b>S751-6R-6</b>	<b>S751-6K-6a</b> <b>S751-6V-6a</b> <b>S751-6W-6a</b> <b>S751-6G-6a</b> <b>S751-6R-6a</b>

# Laboratory Seating

**Specifications:**

The seating shown on this page features chrome-plated metal parts. Upholstery is No. 684 Black Vinyl. Soft

Touch pneumatic mechanism allows easy adjustment to the seating height. Although not shown in the photographs

below, models F-4940-00 thru F-4949-00 come with 2" soft wheel ball bearing casters for hard floors.



- Upholstered seat/backrest
- 5-legged base with 18" fixed footring
- Tubular base

Model No.	Seat Height Adjustment	Footring Height
<b>F-4940-00</b>	16 <sup>3</sup> / <sub>4</sub> "-22"	8 <sup>1</sup> / <sub>2</sub> "
<b>F-4941-00</b>	22 <sup>1</sup> / <sub>4</sub> "-27 <sup>1</sup> / <sub>2</sub> "	8 <sup>1</sup> / <sub>2</sub> "
<b>F-4942-00</b>	26 <sup>1</sup> / <sub>4</sub> "-31 <sup>1</sup> / <sub>2</sub> "	12 <sup>1</sup> / <sub>2</sub> "



- Upholstered seat/backrest
- 5-legged base with 18" fixed foot ring
- Tubular base
- Backrest with lumbar support

Model No.	Seat Height Adjustment	Footring Height
<b>F-4943-00</b>	16 <sup>3</sup> / <sub>4</sub> "-22"	8 <sup>1</sup> / <sub>2</sub> "
<b>F-4944-00</b>	22 <sup>1</sup> / <sub>4</sub> "-27 <sup>1</sup> / <sub>2</sub> "	8 <sup>1</sup> / <sub>2</sub> "
<b>F-4945-00</b>	26 <sup>1</sup> / <sub>4</sub> "-31 <sup>1</sup> / <sub>2</sub> "	12 <sup>1</sup> / <sub>2</sub> "

# Laboratory Seating

**Specifications:**

The seating shown on this page features chrome-plated metal parts. Upholstery is No. 684 Black Vinyl. Soft

Touch pneumatic mechanism allows easy adjustment to the seating height. Although not shown in the photographs

below, models F-4940-00 thru F-4949-00 come with 2" soft wheel ball bearing casters for hard floors.



- 13" round chrome-plated seat
- 5-legged base with 18" fixed foot ring
- Tubular base

Model No.	Seat Height Adjustment	Footring Height
<b>F-4946-00</b>	18 <sup>1</sup> / <sub>2</sub> "-22 <sup>3</sup> / <sub>4</sub> "	8 <sup>1</sup> / <sub>2</sub> "
<b>F-4947-00</b>	22 <sup>1</sup> / <sub>2</sub> "-27 <sup>3</sup> / <sub>4</sub> "	12 <sup>1</sup> / <sub>2</sub> "



- 15" round upholstered seat
- 5-legged base with 18" fixed foot ring
- Tubular base

Model No.	Seat Height Adjustment	Footring Height
<b>F-4948-00</b>	21 <sup>1</sup> / <sub>2</sub> "-24 <sup>3</sup> / <sub>4</sub> "	8 <sup>1</sup> / <sub>2</sub> "
<b>F-4949-00</b>	25 <sup>1</sup> / <sub>2</sub> "-20 <sup>3</sup> / <sub>4</sub> "	12 <sup>1</sup> / <sub>2</sub> "




- 13" round wood seat
- Tubular 4-legged base with 16" fixed foot ring
- Black tubular steel base

Model No.	Seat Height
<b>F-5350-00</b>	18"
<b>F-5351-00</b>	24"
<b>F-5352-00</b>	30"

# Steel Cabinet Technical Information

## Inset Style Wall Bench End Views




**24" Deep Wall Bench**

**30" Deep Wall Bench**

**31" Deep Wall Bench**

## Inset Style Center Bench End Views




**Center Bench End/Section View**


**Center Bench End Sink Cabinet**

# Steel Cabinet Technical Information


## Inset Style Full Height Storage Cabinet End Views



12<sup>3</sup>/<sub>4</sub>" Deep Full Height Cabinet




16" Deep Full Height Cabinet



22" Deep Full Height Cabinet


## Inset Style Center Bench Plan View




Plan View  
Center Bench End Sink Cabinet

# Steel Cabinet Technical Information


## Full Overlay Style Wall Bench End Views



**24" Deep Wall Bench**




**30" Deep Wall Bench**




**31" Deep Wall Bench**

## Full Overlay Style Center Bench End Views



**Center Bench End/Section View**




**Center Bench End Sink Cabinet**



# Steel Cabinet Technical Information

## Full Overlay Style Full Height Storage Cabinet End Views




12 3/4" Deep Full Height Cabinet

16" Deep Full Height Cabinet

22" Deep Full Height Cabinet

## Full Overlay Style Center Bench Plan View



Plan View  
Center Bench End Sink Cabinet

# Lock Schedule Sample and Examples

## Lock Schedule Examples

### Keyed Different & Master-keyed

An example of a simple lock schedule with disk tumbler locks. All locks on the project are keyed different with a single master key that can open all locks. Either disk tumbler or pin tumbler locks could be used for this arrangement.

### LOCK SCHEDULE

Sheet \_\_\_\_ of \_\_\_\_

Project Name: SAMPLE PROJECT				Project/Quote Number		Date:		for Order Dept. use	
Lock Type	<input checked="" type="checkbox"/> Disk Tumbler <input type="checkbox"/> Pin Tumbler	Number of Keys per Lock	2	Number of Master Keys per change	4	Number of Grand-Master Keys per change			
Room Number	Item Number	Cabinet (all cabinets with locks must be listed)	Number of Locks	Keyed Different	Keyed Alike	Key Group	Masterkey Group	Grand-Masterkey Group	
100	1	D30C3518-0101-J1	4	<input checked="" type="checkbox"/>	<input type="checkbox"/>	A			
100	1	D36C3536-0101-J1	6	<input checked="" type="checkbox"/>	<input type="checkbox"/>	A			
100	2	E40C3524-0101-J1	2	<input checked="" type="checkbox"/>	<input type="checkbox"/>	A			
105	6	D30C3524-0101-J1	4	<input checked="" type="checkbox"/>	<input type="checkbox"/>	A			
105	6	E41C3536-0101-J1	3	<input checked="" type="checkbox"/>	<input type="checkbox"/>	A			
105	6	W20C2936-0101-J1	1	<input checked="" type="checkbox"/>	<input type="checkbox"/>	A			

### Keyed Alike in Sets

An example of another simple lock schedule with disk tumbler locks. On this project all the locks in each room are keyed alike, but each room is different from the other. A master key could be used in this type of keying arrangement, but in this example it is not. Either disk tumbler or pin tumbler locks could be used for this arrangement.

### LOCK SCHEDULE

Sheet \_\_\_\_ of \_\_\_\_

Project Name: SAMPLE PROJECT				Project/Quote Number:		Date:		for Order Dept. use	
Lock Type	<input checked="" type="checkbox"/> Disk Tumbler <input type="checkbox"/> Pin Tumbler	Number of Keys per Lock	2	Number of Master Keys per change		Number of Grand-Master Keys per change			
Room Number	Item Number	Cabinet (all cabinets with locks must be listed)	Number of Locks	Keyed Different	Keyed Alike	Key Group	Masterkey Group	Grand-Masterkey Group	
115	5	D30C3518-0101-K1	4	<input type="checkbox"/>	<input checked="" type="checkbox"/>	A			
115	5	E56C3536-0101-K1	4	<input type="checkbox"/>	<input checked="" type="checkbox"/>	A			
115	5	E41C3536-0101-K1	3	<input type="checkbox"/>	<input checked="" type="checkbox"/>	A			
126	16	E40C3524-0101-K1	2	<input type="checkbox"/>	<input checked="" type="checkbox"/>	B			
126	16	E42C3536-0101-K1	5	<input type="checkbox"/>	<input checked="" type="checkbox"/>	B			
126	17	C03C3536-0101-K1	1	<input type="checkbox"/>	<input checked="" type="checkbox"/>	B			
126	17	D30C3524-0101-K1	4	<input type="checkbox"/>	<input checked="" type="checkbox"/>	B			

### Keyed Different – Master-keyed & Grand-Masterkeyed

An example of a more complex lock schedule with pin tumbler locks. All locks are keyed different, with a separate master key for each room, and a grand-master key for the entire project. This keying arrangement could not be used with disk tumbler locks.

### LOCK SCHEDULE

Sheet \_\_\_\_ of \_\_\_\_

Project Name: SAMPLE PROJECT				Project/Quote Number:		Date:		for Order Dept. use	
Lock Type	<input type="checkbox"/> Disk Tumbler <input checked="" type="checkbox"/> Pin Tumbler	Number of Keys per Lock	2	Number of Master Keys per change	5	Number of Grand-Master Keys per change	5		
Room Number	Item Number	Cabinet (all cabinets with locks must be listed)	Number of Locks	Keyed Different	Keyed Alike	Key Group	Masterkey Group	Grand-Masterkey Group	
102	1	E43C3536-0109-J2	6	<input checked="" type="checkbox"/>	<input type="checkbox"/>	A	A		
102	1	E40C3518-0109-J2	2	<input checked="" type="checkbox"/>	<input type="checkbox"/>	A	A		
102	1	D30C3524-0109-J2	4	<input checked="" type="checkbox"/>	<input type="checkbox"/>	A	A		
105	3	D40C3524-0109-J2	2	<input checked="" type="checkbox"/>	<input type="checkbox"/>	B	A		
105	3	D30C3524-0109-J2	4	<input checked="" type="checkbox"/>	<input type="checkbox"/>	B	A		
1113	9	D32C3536-0109-J2	8	<input checked="" type="checkbox"/>	<input type="checkbox"/>	C	A		
1113	9	C03C3536-0109-J2	1	<input checked="" type="checkbox"/>	<input type="checkbox"/>	C	A		
1113	10	D30C3524-0109-J2	4	<input checked="" type="checkbox"/>	<input type="checkbox"/>	C	A		



# Alphabetical Index

<b>A</b>	Acid Storage Cabinet–ADA Height . . . . .	18 & 21		
	Acid Storage Cabinet–Standing Height. . . . .	27 & 32		
	Acid Storage Fume Hood Cabinet. . . . .	21 & 32		
	Acid Storage Cabinet - Full Height . . . . .	57		
	<b>ADA Height Cabinets</b> . . . . .	14-21		
	ADA Compliant Sink . . . . .	94-95		
	Adjustable Slip Scribes . . . . .	65		
	Aerator . . . . .	87		
	Air Valves . . . . .	83-84		
	Applique Panels. . . . .	58-63		
	Apron Rails . . . . .	43		
	Aprons . . . . .	38-43		
	Aspirator . . . . .	87		
<b>B</b>				
	Ball Valves . . . . .	83		
	Blind Corner Apron . . . . .	42		
	Blind Corner Cabinet . . . . .	18 & 29		
	Balance Table . . . . .	79		
	Base Molding. . . . .	81		
<b>C</b>				
	Cabinet Dimensions . . . . .	5 & 100		
	Cabinet Styles . . . . .	3		
	Cart - Stainless Steel . . . . .	79		
	Caster Options. . . . .	6		
	Catalog Number Explanation. . . . .	4		
	Center Bench Shelving . . . . .	72-73		
	Ceiling Enclosure. . . . .	67		
	Chalkboard . . . . .	79		
	Cold Water Faucet. . . . .	85		
	Combination Faucet. . . . .	84		
	Construction Details. . . . .	2		
	Corner Box Apron . . . . .	42		
	Corner Cabinets–ADA Height . . . . .	18		
	Corner Cabinets–Sitting Height. . . . .	13		
	Corner Cabinets–Standing Height. . . . .	29		
	Corner Cabinets–Wall Mounted. . . . .	52		
	Corner Fillers . . . . .	59, 61 & 63		
	Corner Guard. . . . .	81		
	Crumb Cup Strainer. . . . .	91		
	Cupboard/Drawer Cabinets–ADA Height . . . . .	16-17		
	Cupboard/Drawer Cabinets–Sitting Height. . . . .	12-13		
	Cupboard/Drawer Cabinets–Standing Height. . . . .	24-25		
	Cupboard Cabinets–ADA Height . . . . .	14		
	Cupboard Cabinets–Sitting Height . . . . .	11		
	Cupboard Cabinets–Standing Height . . . . .	22		
	Cupsink . . . . .	91		
<b>D</b>				
	Deionized Water Faucet . . . . .	88		
	Diagonal Corner Cabinets–ADA Height . . . . .	18		
	Diagonal Corner Cabinets–Sitting Height . . . . .	13		
	Diagonal Corner Cabinets–Standing Height . . . . .	29		
	Distilled Water Faucet . . . . .	88		
	Double-Sided Wall Cabinet . . . . .	49		
	Drain Outlet . . . . .	91		
	Drain Trough . . . . .	92		
	Drawer Cabinets–ADA Height. . . . .	15		
	Drawer Cabinets–Sitting Height. . . . .	11-12		
	Drawer Cabinets–Standing Height . . . . .	22-23		
	Drench Shower . . . . .	90		
<b>E</b>				
	Electrical Fixtures. . . . .	96-97		
	Emergency Eye Wash/Shower. . . . .	90		
	End Leg Assembly. . . . .	44		
	End View Dimensions . . . . .	100-103		
	Epoxy Resin Cupsink. . . . .	91		
	Epoxy Resin Drain Trough. . . . .	92		
	Epoxy Resin Top . . . . .	70-71		
	Epoxy Resin Tub Sink . . . . .	92-93		
	Eye Wash. . . . .	90		
<b>F</b>				
	Face Wash . . . . .	90		
	Faucet Accessories . . . . .	87		
	Filler Panel . . . . .	58-65		
	Filter Pump. . . . .	87		
	Finished Back . . . . .	64		
	Fixtures–Electrical . . . . .	96-97		
	Foot Valve . . . . .	87		
	Freezer. . . . .	79		
	<b>Full Height Storage Cabinets</b> . . . . .	53-57		
	Full Overlay Applique Panels. . . . .	59, 61 & 63		
	Fume Hood Acid Storage Cabinet–ADA Ht. . . . .	20		
	Fume Hood Acid Storage Cabinet–Standing Ht. . . . .	31		
	Fume Hood Base Cabinet–ADA Height. . . . .	20		
	Fume Hood Base Cabinet–Standing Height . . . . .	30		
	Fume Hoods . . . . .	68-69		
<b>G</b>				
	Gas Valves. . . . .	83-89		
	Gooseneck. . . . .	84-87		
	Gooseneck Coupling. . . . .	87		
<b>H</b>				
	Hoods . . . . .	68-69		
	Hose Spray . . . . .	90		
	Hot & Cold Water Mixing Faucet . . . . .	86-87		
<b>I</b>				
	Intermediate Leg Assembly. . . . .	44		
<b>K</b>				
	Kemresin Work Tops . . . . .	70-71		
	Kemresin Trough . . . . .	92		
	Kemresin Tub Sink. . . . .	92-94		
	Kemstruts. . . . .	80		
	Key Case . . . . .	79		
	Keying Options . . . . .	9-10		
	Knee Space Panel. . . . .	45		
<b>L</b>				
	Label Holder . . . . .	81		
	Laboratory Cart . . . . .	79		
	Laboratory Service Fittings . . . . .	83-90		
	LED Lighting . . . . .	80		
	Leg Assemblies . . . . .	44		
	Lock Keying. . . . .	9-10		
	Lock Panels . . . . .	8		
	Lock Schedule Examples . . . . .	104		
	Lock Schedule Sample . . . . .	105		
	Locks . . . . .	8		

## Alphabetical Index (continued)

<b>M</b>	Markerboard . . . . .	79	Stainless Steel Sink . . . . .	95
	Mobile Cabinets. . . . .	34-37	Stainless Steel Sink Outlet . . . . .	91
	Mixing Faucet . . . . .	86-87	<b>Standing Height Base Cabinets</b> . . . . .	22-33
	Multiple Service Fittings. . . . .	84	Steam Valves . . . . .	84
<b>N</b>			Stools . . . . .	98-99
	Needle Point Cock Valves . . . . .	83-84	Styles . . . . .	3
<b>O</b>			Support Angle . . . . .	80
	Open Cabinets–Sitting Height. . . . .	11	Support Strut . . . . .	72 & 80
	Open Cabinets–ADA Height . . . . .	14	Swinging Glazed Door Full Height Cabinet. . . . .	54
	Open Cabinets–Standing Height. . . . .	22	Swinging Glazed Door Wall Cabinet . . . . .	50
	Open Cabinets–Full Height . . . . .	53	Swinging Panel Door Full Height Cabinet . . . . .	54
	Open Cabinets–Wall Mounted. . . . .	46	Swinging Panel Door Wall Cabinet . . . . .	51
<b>P</b>			Switch . . . . .	97
	Part Number Explanation. . . . .	4	<b>T</b>	
	Pass-Thru Wall Cabinet . . . . .	49	Tackboard . . . . .	79
	Pedestal Legs . . . . .	45	Task Lighting . . . . .	80
	Pegboard . . . . .	78	Tops . . . . .	68-69
	Pipe Enclosure. . . . .	67	Trap . . . . .	91
	Plastic Laminate Work Tops . . . . .	70-71	Trough . . . . .	92
	Plastic Laminate Ledges & Shelves . . . . .	76	Tub Sink. . . . .	92-93
	Plumbing Information. . . . .	82	<b>U</b>	
	Portable Pump Faucet. . . . .	85	Ultra-Pure Water Faucet . . . . .	88
<b>R</b>			Under Counter Refrigerator . . . . .	79
	Reagent Rack Assemblies. . . . .	76	Upright Rod Assembly. . . . .	81
	Reagent Rack Ledges. . . . .	76	<b>V</b>	
	Rear Scribes & Filler . . . . .	59 & 63	Vacuum Breaker . . . . .	87
	Receptacle - Electric . . . . .	96-97	Vacuum Pump Storage Fume Hood Cabinet 20&33	
	Refrigerator . . . . .	79	Vacuum Pump Cart. . . . .	33
	Rods & Accessories . . . . .	81	Vacuum Valves . . . . .	83-84
<b>S</b>			Vandal-Resistant Fittings . . . . .	89
	Safety Shower . . . . .	90	Ventilated Enclosures . . . . .	68-69
	Scribe . . . . .	58-67	<b>W</b>	
	Seating. . . . .	98-99	<b>Wall Cabinets</b> . . . . .	46-52
	Service Fittings . . . . .	82-90	Wall Shelving . . . . .	75-77
	Shelf Brackets . . . . .	77	Wardrobe Storage Cabinet . . . . .	56
	Shelf Lips . . . . .	6 & 73	Waste Storage Cabinet . . . . .	29
	Shelving Assemblies . . . . .	72-75	Water Faucet . . . . .	84-89
	Sink . . . . .	91-95	Water Mixing Faucet . . . . .	86-89
	Sink Cabinets–ADA Height . . . . .	19	Water Valves . . . . .	84
	Sink Cabinets–Sitting Height . . . . .	13	Water/Gas Fitting . . . . .	84
	Sink Cabinets–Standing Height. . . . .	26	Work Tops . . . . .	70-71
	Sink Outlet . . . . .	91	Work Top Support Frame. . . . .	43
	Sink Overflow. . . . .	91	Wrist Blade Handle . . . . .	87
	Sink Trap . . . . .	91		
	<b>Sitting Height Base Cabinets</b> . . . . .	11-13		
	Sloping Tops . . . . .	66		
	Sliding Glazed Door Full Height Cabinet. . . . .	53		
	Sliding Glazed Door Wall Cabinet . . . . .	48		
	Sliding Panel Door Full Height Cabinet . . . . .	53		
	Sliding Panel Door Wall Cabinet . . . . .	48		
	Sliding Plate Glass Door Wall Cabinet. . . . .	47 & 49		
	Solvent Storage Cabinet . . . . .	19-20, 28 & 32		
	Spreader . . . . .	45		
	Spring Clips . . . . .	65		
	Stainless Steel Cart . . . . .	79		
	Stainless Steel Cupsink . . . . .	91		

# Numerical Index

<b>Part Number</b>	<b>Page</b>
019623 - 019627 .....	31
0464-00 - .....	91
0469-00 - .....	91
0482-BP - 0482-GP .....	91
0490-BE - 0490-SE .....	91
0491-BP - 0491-GP .....	91
0492-BP - 0492-GP .....	91
0493-BE - 0493-SE .....	91
0494BP-04 - 0494BP-10 .....	91
0548-00 - 0557-00 .....	81
0720-1K-6 - 0720-1R-6a .....	97
0950-0F - 0975-0F .....	91
0995-00 - .....	91
1000-00-__ - 1039-00-__ .....	93
1000-DI-__ - 1070-DI-__ .....	94
1000-ADA-DI-__ - 1006-ADA-DI-__ .....	94
4740-00 .....	81
A00C022215-00_0 - A00C022272-00_0 .....	43
A00C042215-__0 - A00C042272-__0 .....	38
A00C042815-__0 - A00C042872-__0 .....	39
A04C042248-1__ - A05C042872-1__ .....	42
A06C020215-00_0 - A06C040272-00_0 .....	43
A06C040315-1__0 - A06C040372-1__0 .....	43
A07C292202-00_0 - A08C352802-00_0 .....	44
A09C020230-00_0 - A09C020272-00_0 .....	45
A15C290202-00_0 - A15C350202-00_0 .....	45
A10C250015-00_0 - A10C350072-00_0 .....	45
A11C292202-00_0 - A13C355202-00_0 .....	45
A27C292202-00_0 - A28C352802-00_0 .....	44
A51C042220-1__ - A53C042272-1__ .....	40
A51C042820-1__ - A53C042872-1__ .....	41
A56C040236-00_0 - A56C040272-00_0 .....	43
A57C352202-00_0 - A58C352802-00_0 .....	44
A65C350202-00_0 .....	45
A77C352202-00_0 - A78C352802-00_0 .....	44
ABSC290012-00_0 - ABSC350012-00_0 .....	65
AWSC241312-00_0 - AWSC481612-00_0 .....	65
AFSC841312-00_0 - AFSC842212-00_0 .....	65
B00C292215-00_0 - B00C292248-00_0 .....	11
B00C322215-00_0 - B00C322248-00_0 .....	14
B00C352215-00_0 - B00C352248-00_0 .....	22
BCAC290003-11_0 - BCAC350012-17_0 .....	59
BCFC290001-00_0 - BCFC350012-00_0 .....	59
BEAC290003-11_0 - BEAC350024-17_0 .....	58
BESC290001-00_0 - BESC350024-00_0 .....	58
BFBC290112-00_0 - BFBC350160-00_0 .....	64
BIAC290006-11_0 - BIAC350048-17_0 .....	58
BIFC290006-11_0 - BIFC350048-00_0 .....	58
BPC290001-00_0 - BPC350024-00_0 .....	59
BRSC290001-00_0 - BRSC350024-00_0 .....	59
BPAC290023-16_0 - BPAC350029-16_0 .....	59

<b>Part Number</b>	<b>Page</b>
C03C292215L____ - C03C292248-____ .....	11
C03C322215L____ - C03C322248-____ .....	14
C03C352212L____ - C03C352248-____ .....	22
CCEC121320-00_0 - CCEC241615-00_0 .....	67
CSTC081327-00_0 - CSTC111627-00_0 .....	66
D00C292215-____ - D00C292248-____ .....	11
D00C322215-____ - D00C322248-____ .....	15
D00C352215-____ - D00C352248-____ .....	22
D02C322230-____ - D02C322248-____ .....	15
D02C352230-____ - D02C352248-____ .....	22
D08C352215-____ - D08C352248-____ .....	22
D26C292218-____ - .....	11
D30C292215-____ - D30C292248-____ .....	11
D30C322215-____ - D30C322248-____ .....	15
D30C352212-____ - D30C352248-____ .....	23
D31C292230-____ - D31C292248-____ .....	11
D31C322230-____ - D31C322248-____ .....	15
D31C352230-____ - D31C352248-____ .....	23
D32C292230-____ - D32C292248-____ .....	11
D32C322230-____ - D32C322248-____ .....	15
D32C352230-____ - D32C352248-____ .....	23
D36C322230-____ - D36C322248-____ .....	15
D36C352230-____ - D36C352248-____ .....	23
D50C322215-____ - D50C322248-____ .....	15
D50C352215-____ - D50C352248-____ .....	23
D51C322230-____ - D51C322248-____ .....	15
D51C352230-____ - D51C352248-____ .....	23
D52C322230-____ - D52C322248-____ .....	15
D52C352230-____ - D52C352248-____ .....	23
D81C292218-____ - D81C292248-____ .....	12
D81C352218-____ - D81C352248-____ .....	23
D83C322215-____ - D83C322248-____ .....	15
D84C292218-____ - D84C292248-____ .....	12
D84C322215-____ - D84C322248-____ .....	15
D84C352218-____ - D84C352248-____ .....	23
D85C292218-____ - D85C292248-____ .....	12
E05C292230L____ - E05C292248-____ .....	12
E05C322230L____ - E05C322248-____ .....	16
E05C352230L____ - E05C352248-____ .....	24
E40C292215L____ - E40C292248-____ .....	12
E40C322215L____ - E40C322248-____ .....	16
E40C352212L____ - E40C352248-____ .....	24
E41C292230-____ - E41C292248-____ .....	12
E41C322230-____ - E41C322248-____ .....	16
E41C352230-____ - E41C352248-____ .....	24
E42C322230L____ - E42C322248-____ .....	16
E42C352230L____ - E42C352248-____ .....	24
E43C322230L____ - E43C322248-____ .....	16
E43C352230L____ - E43C352248-____ .....	24
E45C292230L____ - E45C292248-____ .....	13
E45C322230L____ - E45C322248-____ .....	16
E45C352230L____ - E45C352248-____ .....	24
E56C322236-____ - E56C352248-____ .....	16
E56C352236-____ - E56C352248-____ .....	24
E60C352212L____ - E60C352248-____ .....	25
E61C352230-____ - E61C352248-____ .....	25
E62C322230L____ - E62C322248-____ .....	17
E62C352230L____ - E62C352248-____ .....	25

# Research Collection

## Numerical Index (continued)

Part Number	Page	Part Number	Page
E63C322230L_____ - E63C322248-_____	17	G35C322024L_____ - G35C322248-_____	20
E63C352230L_____ - E63C352248-_____	25	G35C352024L_____ - G35C352248-_____	33
E65C322230L_____ - E65C322248-_____	17	G63C322018L_100 - G63C322248-_____	19
E65C352230L_____ - E65C352248-_____	25	G63C352024L_100 - G63C352248-_____	28
ECEC120013-00_0 - ECEC240022-00_0	67	G68C322024L_100 - G68C322248-_____	20
F-0260-00 - _____	81	G68C352024L_100 - G68C352248-_____	32
F-0862-04 - _____	81	G80C322018L_____ - G80C322248-_____	21
F-1368-00 - _____	81	G80C352018L_____ - G80C352260-_____	31
F-2939-00 - F-2940-00	78	H581-1K - H581-8R	96
F-4940-00 - F-4945-00	98	K-0005-0A - K-0020-1A	80
F-4946-00 - F-5352-00	99	K7-302435-BK - K7-302435-SL	79
F-4995-PK - _____	65	K12-2905-0A - K12-3520-0A	80
F-5332-00 - F-5332-01	79	K17-3612-00 - K17-3614-00	73
F-5333-00 - F-5348-00	79	KF-0170-04 - _____	77
F-5349-00 - _____	79	K-MSS-12 - K-MSS-60	13, 19 & 26
F-5353-06 - F-5362-12	95	L-105 - L-115	79
F-8440-15 - F-8441-48	27	L-2034-06 - L-2034-96	80
F10M0242-00 - F10M0279-00	74	L01-020109-00 - L01-020142-00	47 & 81
F11M0196-00 - F11M0372-00	74	M-0051-0A - M0052-0A	80
F12M000224-00 - F13S000284-00	72	M-0030-0A - M-0044-AD	77
F17-4824-00 - F17-4872-00	72	M-0060-0A - M-0062-AD	77
F21M0624-00 - F21W2448-0_	75	M10C262215L_____ - M10C262224-_____	34
F21M0924-60 - F21W1848-6_	73	M10C332215L_____ - M10C332224-_____	35
F63M0224-00 - F63M0272-00	73	M20C262215L_____ - M20C262224-_____	34
F73M0015-00 - F73S0024-00	37	M20C332215L_____ - M20C332224-_____	35
F94M2215-00 - F94M2348-00	37	M25C262215L_____ - M25C262224-_____	34
F94P2215-00 - F94T2448-00	36	M25C332215L_____ - M25C332224-_____	35
FCAC840003-11_0 - FCAC840012-16_0	63	M30C262215-_____ - M30C262224-_____	34
FCEC120018-00_0 - FCEC240096-00_0	67	M35C262215-_____ - M35C263224-_____	34
FCCFC841301-00_0 - FRSC840024-00_0	63	M40C332215-_____ - M40C332224-_____	35
FEAC840003-11_0 - FEAC840024-16_0	62	M45C262215-_____ - M45C263224-_____	34
FESC841301-11_0 - FESC842224-00_0	62	M46C332215-_____ - M46C332224-_____	35
FFBC840118-00_0 - FFBC840148-00_0	64	M56C332215-_____ - M56C332224-_____	35
FIAC840006-11_0 - FIAC840024-16_0	62	M62C262215-_____ - M62C262224-_____	34
FIFC841301-00_0 - FIFC842224-00_0	62	M66C262215-_____ - M66C262224-_____	34
FPAC840013-16_0 - FPAC840022-16_0	63	M68C332215-_____ - M68C332224-_____	35
FRSC560001-00_0 - FRSC840024-00_0	63	M71C332215-_____ - M71C332224-_____	35
G00C292215L_____ - G00C292248-_____	13	PL-4812-00 - PL-7214-00	76
G00C322212L_____ - G00C322248-_____	19	RRU-1210-0000 - RRU-1812-0030	76
G00C352212L_____ - G00C352260-_____	26	S00C841318-00-0 - S00C842248-00_0	53
G03C322236-_____ - G03C322248-_____	21	S12C841330-_____0 - S12C842248-_____0	53
G03C352230-_____ - G03C352248-_____	26	S25C841318L_____ - S25C842248-_____	54
G04C352252-_____ - G04C352260-_____	26	S35C841330-_____ - S35C842248-_____	56
G07C322236-_____ - G07C322248-_____	21	S40C842218L_____ - S40C842248-_____	56
G08C322012L_____ - G08C322248-_____	20	S41C841330-_____0 - S41C842248-_____0	53
G08C352012-_____ - G08C352260-_____	30	S50C841318L_____ - S50C842248-_____	54
G12C322248-_____ - G13C322248-_____	18	S51C841624L_____ - S51C842248-_____	55
G12C352248-_____ - G13C352248-_____	29	S60C841330-_____ - S60C842248-_____	56
G14C292214-_____ - _____	13	S70C841624L_____ - S70C842248-_____	57
G14C322214-_____ - G15C352214-_____	18		
G14C352214-_____ - G15C352224-_____	29		
G20C352215-_____ - G20C352224-_____	29		
G23C352218-_____ - _____	29		
G21C322018L_____ - G21C322248-_____	18		
G21C352015L_____ - G21C352248-_____	27		
G23C322218-_____ - _____	19		

# Numerical Index (continued)

<b>Part Number</b>	<b>Page</b>	<b>Part Number</b>	<b>Page</b>
S581-1K	- S581-8R ..... 96	W00C241315-00_0	- W00C481648-00_0 ..... 46
S582-1K	- S582-8R ..... 96	W10C241330-00_0	- W10C481648-00_0 ..... 47
S655-1K	- S655-1R ..... 97	W12C241330-__0	- W12C481648-__0 ..... 48
S656-1K	- S656-3R ..... 96	W13C241330-00_0	- W13C361648-00-0 ..... 49
S657-6K	- S657-8R ..... 96	W25C241315-__	- W25C481648-__ ..... 51
S658-6K	- S658-8R ..... 96	W26C241320-__	- W26C481615-__ ..... 52
S681-1K	- S681-8R ..... 96	W41C241330-__0	- W41C481648-__0 ..... 48
S721-6K-6	- S721-6R-6a ..... 97	W50C241315-__	- W50C481648-__ ..... 50
S751-6K-6	- S751-6R-6a ..... 97	W51C241320-__	- W51C481615-__ ..... 52
SCAC010103-00_0	- SCAC010196-00_0 ..... 65	WCAC240003-11_0	- WCAC480012-16_0 ..... 61
SCFC081301-00_0	- SCFC142206-00_0 ..... 66	WCFC241301-00_0	- WCFC481612-00_0 ..... 61
SLTC081315-00_0	- SLTC142296-00_0 ..... 66	WEAC240003-11_0	- WEAC470024-16_0 ..... 60
SM-01-00-__	- SM-56-00-__ ..... 93	WESC241301-00_0	- WESC481624-00_0 ..... 60
STFC081301-00_0	- STFC142212-00_0 ..... 66	WFBC240115-00_0	- WFBC480148-00_0 ..... 64
TR-01-00-__	- TR-05-00-__ ..... 92	WIAC240006-11_0	- WIAC480048-16_0 ..... 60
U01C600512-00_0	- U02C961048-00_0 ..... 67	WIFC241301-00_0	- WIFC481648-00_0 ..... 60
VPCC061414-00_0	- VPCC061644-00_0 ..... 33	WPAC240013-16_0	- WPAC480016-16_0 ..... 61
W-0242-00	- W-0243-00 ..... 83	X-020004	- X-020017 ..... 78
W-0259-00	..... 84	X-030001-00	- X-030003-00 ..... 80
W-0260-00	- W-0272-00 ..... 83		
W-0273-02	- W-0274-VB ..... 84		
W-0286-02	- W-0287-00 ..... 83		
W-0300-00	- W-0306-00 ..... 84		
W-0333-00	..... 86		
W-0337-00	- W-0337-6V ..... 85		
W-0340-00	- W-0340-0V ..... 86		
W-0350-00	- W-0350-0V ..... 85		
W-0351-00	- W-0351-0V ..... 85		
W-0354-00	- W-0355-00 ..... 86		
W-0358-00	..... 87		
W-0361-00	- W-361-0V ..... 87		
W-0367-00	- W-0368-00 ..... 88		
W-0374-01	..... 90		
W-0400-00	- W-0401-00 ..... 83		
W-0402-00	- W-0404-00 ..... 84		
W-0405-00	..... 84		
W-0407-00	- W-0409-00 ..... 88		
W-0410-VR	- W-0417-VR ..... 89		
W-0510-00	..... 87		
W-0520-00	..... 87		
W-0532-00	..... 87		
W-0535-00	..... 87		
W-0825-00	..... 87		
W-0833-00	..... 87		
W-0923-00	..... 90		
W-0926-00	..... 90		
W-0927-00	..... 90		
W-0928-00	..... 90		
W-0930-00	..... 90		
W-0931-00	..... 90		
W-0932-00	..... 90		
W-0933-00	..... 90		
W-0934-00	..... 90		
W-0935-00	..... 90		
W-0936-00	..... 90		



# Research Collection

# Typical Installations



## Trademarks

Kewaunee, Alpha, Basic Bench, EarthLine, Enterprise, Evolution, Interceptor, Kemresin, Kewaunee Matrix, Research Collection, Signature, Sturdillite, and Supreme Air are registered trademarks of Kewaunee Scientific Corporation.

## Patents:

Products in this catalog are protected by one or more of the following patents:

4,856,420	5,241,788	5,302,030	5,407,389
5,582,472	5,728,626	D388,986	DES. 346,437

Additional patents are pending.

Part numbers, illustrations, descriptions, specifications and other information contained in this catalog are based on the most current design and production information available

at the time of publication. Kewaunee Scientific Corporation reserves the right to make changes at any time without notice.

# Typical Installations



We are pleased to have been selected to supply laboratory furniture for the facilities of some of the best companies, organizations and institutions in the world. We believe the repeat opportunities, over time, to participate in their construction projects is good evidence of our product quality and experience in successfully completing their laboratory projects.

Whether your laboratory projects are large or small, new construction or remodeling, Kewaunee has the qualified people and quality products to transform your visions into successful reality.

From pre-planning to final completion, experienced Kewaunee representatives are available to help you through the intricacies of laboratory design and layout, product selection, ordering, scheduling and installation.

Kewaunee Scientific Corporation is dedicated to manufacturing high quality products for the laboratory marketplace. We offer **Total Laboratory Solutions, Innovative Designs, and Technological Expertise** with a worldwide distribution network to provide laboratories that are truly World Class.

**Laboratory Products from Kewaunee®:**

ALPHA® – Adaptable Casework System

ALPHA Columns – Expandable Workstations

ALPHA Overhead Carrier Systems

BasikBench® – Welded C-Leg Bench

DetectAir – Filtered/Ductless Fume Hoods

EarthLine® – Environmentally Friendly Casework

Enterprise® – Moveable Workstations

Evolution® – Column Based Workstations

Interceptor® – Biological Safety Cabinet

Interceptor® – Laminar Flow Cabinet

Kemresin® – Epoxy Resin Counters & Sinks

Kewaunee Matrix® – Educational Stations

Research Collection® – Steel Casework

Signature® Series – Wood Casework

Spektrum – Laminate Architectural Millwork

Sturdilite® – Flexible Pedestal System

Supreme Air® – Fume Hoods

For our complete product offering, specifications, catalogs, brochures, case studies, and dealer list, go to [www.kewaunee.com](http://www.kewaunee.com), or contact your local sales representative directly.



Kewaunee, Alpha, Basic Bench, EarthLine, Enterprise, Evolution, Interceptor, Kemresin, Kewaunee Matrix, Research Collection, Signature, Sturdilite, and Supreme Air are registered trademarks of Kewaunee Scientific Corporation



P.O. Box 1842 • Statesville, NC 28687-1842  
Phone: (704) 873-7202 • Fax: (704) 873-5160  
E-mail: [kscmarketing@kewaunee.com](mailto:kscmarketing@kewaunee.com)  
Web Site: [www.kewaunee.com](http://www.kewaunee.com)

**Kewaunee Labway India Pvt. Ltd.**

No. CA-9A, 2nd Floor, Jigani Link Road  
Jigani, Bangalore 562 106, India  
Phone: 9180-27826725 • Fax: 9180-27826724  
Web Site: [www.kewaunee.in](http://www.kewaunee.in)

**Kewaunee Labway Asia Pte. Ltd.**

194 Pandan Loop, #06-22  
Pantech Business Hub, Singapore 128 383  
Phone: 65-6773-0288 • Fax: 65-6773-2322

**Kewaunee Scientific (Suzhou) Co Ltd**

333 Xingpu Road, Modern Industrial Square  
Building 4, Unit #2-A • Suzhou Industrial Park • China  
Phone: (0512) 6262 5090 • Fax: (0512) 6262 5089  
Web Site: [china.kewaunee.com](http://china.kewaunee.com)

# Research Collection

Kewaunee Steel Cabinet Construction	2
Kewaunee Steel Cabinet Styles	3
Steel Cabinet Information	4
Steel Cabinet Options	6
Lock Keying Options	9
Steel Sitting Height Base Cabinets	11
Steel ADA Height Base Cabinets	14
Steel Standing Height Base Cabinets	22
Fume Hood Base Cabinets	30
Steel Sitting Height Mobile Cabinets	34
Steel Standing Height Mobile Cabinets	35
Mobile Cabinet Tops & Accessories	36
Steel Apron Assemblies	38
22" Deep Apron	38
22" Deep Apron	38
Steel Box Apron	38
22" Deep Apron	38
Applique Panel	38
Steel Box Apron	38
28" Deep Apron	39
28" Deep Apron	39
Steel Box Apron	39
28" Deep Apron	39
Applique Panel	39
Steel Box Apron	39
Steel Leg Assemblies & Pedestals	44
Steel Apron Rails & Panels	45
Steel Wall Mounted Storage Cabinets	46
Steel Full Height Storage Cabinets	53
Base Cabinet Scribes & Fillers	58
Wall Cabinet Scribes & Fillers	60
Full Height Cabinet Scribes & Fillers	62
Steel Cabinet Finished Backs	64
Steel Adjustable Slip Scribes	65
Steel Sloping Tops	66
Steel Ceiling & Pipe Enclosures	67
Fume Hoods & Ventilated Enclosures	68
Work Tops	70
Center Bench Shelving Assemblies	72
Wall Shelving Assemblies	74
Reagent Rack & Shelving Systems	76
Pegboards	78
Miscellaneous Accessories	79
Miscellaneous Accessories	80
K12-	80
-0A	80
General Plumbing Information	82
Laboratory Service Fittings	83
Cupsinks & Drain Accessories	91
Kemresin Tub Sinks & Drain Troughs	92
Kemresin Drop-in Sinks	94
Stainless Steel Drop-in Sinks	95
Laboratory Electrical Fixtures	96
Laboratory Seating	98
Steel Cabinet Technical Information	100
Lock Schedule Sample and Examples	104
Alphabetical Index	106
Numerical Index	108
Typical Installations	112

# Research Collection

## Steel Laboratory Furniture



KEWAUNEE<sup>®</sup>

encouraging new discovery...*Worldwide*

# RESEARCH COLLE

## *Steel Furniture*

- Fully-welded, heavy-gauge steel construction
- Constructed to meet SEFA standards
- Powder coat, lab grade finish
- Stainless steel construction option
- Customize with multiple styles, pulls, and hardware



# CTION



## Contents

Steel Furniture Styles & Construction .....	2
Catalog Number Explanation .....	4
Steel Cabinet Options .....	6
Sitting Height Cabinets.....	11
ADA Height Cabinets.....	14
Standing Height Cabinets.....	22
Mobile Cabinets.....	34
Apron Assemblies .....	38
Wall Mounted Storage Cabinets.....	46
Full Height Storage Cabinets.....	53
Fillers, Scribes & Panels .....	58
Sloping Tops & Enclosures.....	66
Fume Hoods & Ventilated Enclosures ...	68
Work Tops .....	70
Center Bench Shelving.....	72
Wall Shelving.....	74
Reagent Racks & Shelving.....	76
Pegboards .....	78
Miscellaneous Accessories.....	79
General Plumbing Information.....	82
Service Fittings .....	83
Drain Accessories .....	91
Sinks & Troughs .....	92
Electrical Fittings .....	96
Seating .....	98
Technical Information.....	100
Alphabetical Index.....	106
Numerical Index .....	108
Typical Installations .....	112

Beautifully finished, smooth, clean, easy to maintain surfaces. Rugged, strong, heavy gauge steel, full-frame, four sided, flush front construction. That's Kewaunee's Steel Laboratory Furniture.



# Kewaunee Steel Cabinet Construction

Kewaunee's fully-welded, heavy-gauge, steel cabinets are designed to SEFA 8 performance standards and built to rigorous benchmarks. The Kewaunee classic "unit" design allows doors

and drawers to be reconfigured in place. The configurable design provides a rich offering of hardware and aesthetic options.

## DRAWER OPTIONS



Standard Drawers feature one-piece drawer body construction with full-extension, side mounted, soft-close, slides. (model F-8400-00)




Optional Concealed Slide Drawers feature integral concealed, full-extension, soft-close slides



### Additional Cabinet Features

- Four sided, flush front, frame construction
- Solid finished back with optional removable back
- Full-depth adjustable shelf
- Leveling bolts accessible through cabinet top

## PULL STYLE OPTIONS




## DOOR CATCH



Positive Catch

Catch furnished with 5-knuckle hinges only

## DOOR HINGE OPTIONS



# Kewaunee Steel Cabinet Styles



*shown with five knuckle hinge and rounded aluminum pull style 1*

## Inset - Square Edge

Doors and drawer heads are flush with face of cabinet. Available with Positive Catch and a wide selection of pulls.



*shown with maple veneer door and drawer head, five knuckle hinge and square aluminum pull style 8*

## Overlay - Wood on Steel

Oak or maple veneer doors and drawer heads overlay face of cabinet. Available with 5-knuckle, semi-concealed or concealed hinges and a wide selection of pulls. *(Flush Pull 9 not available on wood)*



*shown with five knuckle hinge and chrome plated wire pull style 6*

## Overlay - Square Edge

Doors and drawer heads overlay face of cabinet. Available with 5-knuckle, semi-concealed or concealed hinges and a wide selection of pulls.



*shown with concealed hinge and beveled label holder pull style 7*

## Overlay - Bevel Edge

Doors and drawer heads, with beveled top edge, overlay face of cabinet. Available with 5-knuckle, semi-concealed or concealed hinges and a wide selection of pulls.

# Steel Cabinet Information

## Catalog Number Explanation

**E40C352236-1611-B11,C...**

### Cabinet Type

A = Apron Assemblies and Legs  
 B = Open Cupboard Base Cabinet  
 C = All Cupboard Base Cabinet  
 D = All Drawer Base Cabinet  
 E = Drawer & Cupboard Base Cabinet  
 G = Special Purpose Base Cabinet  
 S = Full Height Storage Cabinet  
 W = Wall Storage Cabinet

### Door & Drawer Configuration

### Cabinet Construction & Design

### Cabinet Height

#### Base Cabinets

29 = 29" High    32 = 32" High    35 = 35" High

#### Aprons

4 = 4" High

#### Full Height Cabinets

84 = 84" High

#### Wall Cabinets

24 = 24" High    30 = 30" High    36 = 36" High    48 = 48" High

### Cabinet Options

(See Cabinet Option list page 5)

### Style – Hardware – Material

(See Selection list below)

### Special Configurations

L = Reverse Construction/Left Hand Swing

### Cabinet Length

Standard Lengths: 12" – 15" – 18" – 21"

24" – 30" – 36" – 42"

48" – 52" – 58" – 60"

### Cabinet Depth

13 = 12<sup>3</sup>/<sub>4</sub>" Deep

20 = 20" Deep

16 = 16" Deep

22 = 22" Deep

28 = 28" Deep

Note: Full Overlay Styles add 7/8" to cabinet depth

## Door & Drawer Head Style – Hardware – Material Selections

**E40C352236-1611-B11,C...**

### Cabinet Style - Hinge Type - Catch

0 = Inset            - 5-knuckle Hinge - Positive Catch  
 1 = Full Overlay   - 5-knuckle Hinge - Positive Catch  
 4 = Full Overlay   - Semi-concealed Hinge  
 6 = Full Overlay   - Concealed Hinge (170° Swing)

### Door & Drawer Head Style

1 = Square Edge Steel  
 6 = Wood – available on Full Overlay Cabinet Styles only  
 7 = Beveled Edge Steel – available on Full Overlay Cabinet Styles only

### Cabinet Material / D&D Material

0 = Cold Rolled Steel / Cold Rolled Steel    \*3 = Stainless Steel / Stainless Steel  
 1 = Cold Rolled Steel / Red Oak                5 = Cold Rolled Steel / FSC® Red Oak  
 2 = Cold Rolled Steel / White Maple           6 = Cold Rolled Steel / FSC® White Maple  
 7 = Cold Rolled Steel / Steamed Beech        8 = Cold Rolled Steel / FSC® Steamed Beech


\* Stainless Steel Cabinets available in D&D Style 1 only

### Pull Type

0 = No Option  
 1 = Rounded Aluminum  
 2 = Black Wire  
 3 = White Wire  
 4 = Stainless Steel Wire  
 6 = Chrome Plated Wire  
 ‡7 = Beveled Label Holder  
 8 = Square Aluminum  
 †9 = Flush Aluminum  
 ‡ mounts horizontally  
 † not available on wood doors & drawer heads

# Steel Cabinet Information

## Drawer and Cupboard Size Coding



**Drawer Code**

**Cupboard Code**

<b>E40C352215L</b> _____	15"L
<b>E40C352218L</b> _____	18"L
<b>E40C352221L</b> _____	21"L
<b>E40C352224L</b> _____	24"L

Base cabinets are designed on a "unit" system to allow maximum flexibility in the selection of drawer and/or cupboard door sizes.

Standing height cabinets (35" high) have an opening 60 "units" high. Sitting height cabinets (29" high) have an opening 48 "units" high. ADA height cabinets (32" high) have an opening 54 "units" high. This means that any combination of drawer and cupboard door heights equaling 60, 54 or 48 units can be used in a base cabinet.

Mobile standing height cabinets (33" high) have an opening 56 "units" high. Mobile sitting height cabinets (26 1/4" high) have an opening 42 "units" high.

The adjacent tables indicate the "unit" height and inside clearance of the Inset style base cabinets shown in this catalog. Simply match the unit number shown on the cabinet drawers and cupboard to the appropriate number in the tables.

### Drawer Size Codes

"Units" High	Drawer Head Height	Inside Dimensions	
		Clear/Height	Length (Front-to-Back)
8*	3 3/4"	2 3/8"	18"
9*	4 1/4"	2 3/8"	18"
10*	4 3/4"	3 7/16"	18"
11	5 3/16"	3 15/16"	18"
12	5 11/16"	4 7/16"	18"
14	6 5/8"	5 3/8"	18"
15	7 1/8"	5 7/8"	18"
16	7 5/8"	6 3/8"	18"
20	9 9/16"	8 5/16"	18"
24	11 1/2"	10 1/4"	18"
3*	1 3/16"	† Pull Out Shelf	

\* Not available with Option Z Concealed Slide Drawer  
 † Pull Out Shelf mounted on full extension slides only

### Drawer Width I.D. (side-to-side)

Base Cabinet Length	Number of Drawers Across		
	1	2	3
12"	8 3/8"		
15"	11 3/8"		
18"	14 3/8"		
21"	17 3/8"		
24"	20 3/8"	8 3/8"	
30"	26 3/8"	11 3/8"	
36"	32 3/8"	14 3/8"	8 3/8"
42"	38 3/8"	17 3/8"	
48"	44 3/8"	20 3/8"	12 3/8"

### Cupboard Size Codes

"Units" High	Door Height	Inside Dimensions	
		Height (Top-to-Bottom)	Depth (Front-to-Back)
30	14 9/16"	13 9/16"	20 1/2"
36	17 7/16"	16 7/16"	20 1/2"
45	21 13/16"	20 13/16"	20 1/2"
48	23 3/16"	22 3/16"	20 1/2"
54	26 1/16"	25 1/16"	20 1/2"
60	29"	28"	20 1/2"

### Door Opening Clearance

#### Width I.D. (side-to-side)

Base Cabinet Length	Number of Drawers Across		
	1	2	3
12"	9 1/2"		
15"	12 1/2"		
18"	15 1/2"		
21"	18 1/2"		
24"	21 1/2"	9 1/2"	
30"	27 1/2"	12 1/2"	
36"	33 1/2"	15 1/2"	9 1/2"
42"	39 1/2"	18 1/2"	
48"	45 1/2"	21 1/2"	13 1/2"

Cupboard and Drawer sizes listed in these charts are for Inset style cabinets only. Sizes for Overlay styles are generally slightly larger. Contact your Kewaunee Dealer for more information.

## Cabinet Options – (Add to End of Catalog Number)

### B – Base Cabinet Caster Option

- B11** – 4" Swivel Caster w/ brakes & with counter balance
- B21** – 4" Swivel Caster w/o brakes & with counter balance
- B31** – 4" Swivel Caster w/ brakes at front & counter balance
- B41** – 3" Swivel Caster w/ brakes & with counter balance
- B51** – 3" Swivel Caster w/o brakes & with counter balance
- B61** – 3" Swivel Caster w/ brakes at front & counter balance
- B71** – Fixed rear – Swivel front – with counter balance
- B81** – Fixed rear – Glides front – with counter balance

**C** – Toe space depth increased by 7/8" (for B00 & G07)

**F** – Removable Back (not available on Mobile Cabinets)

**H** – Label Holders \*

\* **J** – Locks w/ Lock Panels    **J**=Disk Locks†    **J2**=Pin Locks†

\* **K** – Locks w/o Lock Panels    **K**=Disk Locks†    **K2**=Pin Locks†

**O** – Shelf Lips

### P – Pullout Trays & Pullout Pans

**P1** – One Pullout Tray - at bottom of cupboard

**P2** – Two Pullout Trays - center & bottom

**P3** – One Pullout Pan - at bottom of cupboard

**P4** – Two Pullout Pans - center & bottom

**I** – 1/4" Laminated Safety Glass

**R** – 1/8" Tempered Safety Glass

**T** – Heavy Duty Full Extension Slide

**U** – Toe Space on Base Cabinet Ends

**U1** – Toe Space on Left End

**U2** – Toe Space on Right End

**U3** – Toe Space on Both Ends

**Z** – Concealed Slide Drawer

**Z1** – Cold Rolled Steel Drawer Sides

**Z2** – Stainless Steel Drawer Sides

\* Furnished on all cupboards and drawers on cabinet    † Lock keying must be specified.

# Steel Cabinet Options

## Base Cabinet Caster Option

Base Cabinet Caster Option adds casters to the base cabinet. Eight caster styles are available and all are furnished with counter balance weight for cabinets that may be top heavy when loaded and/or prone to tipping when opened. (**Note:** Casters ship loose)

### Swivel Casters

adds four (4) swivel casters and counter balance weight for stability.

- B11** – \*4" Swivel Casters with brakes
- B21** – \*4" Swivel Casters without brakes
- B31** – \*4" Swivel Casters with brakes at front only
- B41** – 3" Swivel Casters with brakes
- B51** – 3" Swivel Casters without brakes
- B61** – 3" Swivel Casters with brakes at front only

\* Adds 1" to height of cabinet

### Fixed and Swivel Casters

adds two (2) fixed 4" casters at the rear, two locking swivel casters at the front, and counter balance weight for stability.

- B71** – Fixed rear casters, locking front swivel casters

### Fixed Casters and Glides

adds two (2) fixed 4" casters at the rear, two (2) adjustable glides at the front, and counter balance weight for stability

- B81** – Fixed rear casters, front adjustable glides.

To specify Casters, add **B**\_\_ to the end of the cabinet number.

Example: **E40C352236-0101-B71**  
(Available on base cabinets only.)

## Toe Space Depth Increase by 7/8"

Option **C** adds 7/8" to the depth of a base cabinet toe space. This option is only used on base cabinets without doors or drawers (B00 & G07) to

allow the face of the cabinet to be pulled forward to align with the face of the doors and drawers on adjacent cabinets.

To specify increased toe space depth add **C** to the end of the

cabinet number.  
Example: **B00C352236-0000-C**

## Removable Back

All cabinets are furnished with a solid back for improved cleanliness. A removable back is available when required for access to the pipe space behind the cabinet.

**Note:** Removable Backs are not available on Mobile Cabinets.

To specify Removable Back, add **F** to the end of the cabinet number.

Example: **D30C352236-0101-F**

## Label Holders

Option **H** adds the F-0260-00 aluminum Label Holder to each drawer and door on the cabinet. Overall size is 3" x 1<sup>5</sup>/<sub>8</sub>" and holds a 2<sup>1</sup>/<sub>2</sub>" x 1<sup>1</sup>/<sub>8</sub>" card. Double doors


receive only one Label Holder per pair.

**Shipped loose for field installation.**


To specify Label Holders, add **H** to the end of the cabinet number.

Example: **E40C352236-0101-H**  
(Added to all cupboards and drawers.)


## Option B



**B11 – B21 – B31 Caster Option**




**B71 Caster Option**




**B81 Caster Option**

## Option C

## Option F



## Option H



One-half Actual Size

## Steel Cabinet Options

### Shelf Lips

### Option O

Option **O** adds a 3/4" high steel lip to the front edge of all steel shelves in the cabinet.

To specify Shelf Lips, add **O** to the end of the cabinet number.

Example: **W00C361336-0101-O**  
(Added to all shelves in cabinet.)

### Pullout Trays & Pans

### Option P

**P1** – One Pullout Tray - at bottom of cupboard

**P2** – Two Pullout Trays - center & bottom

**P3** – One Pullout Pan - at bottom of cupboard

**P4** – Two Pullout Pans - center & bottom

Pullout Trays and Pans are fabricated from 18 gauge cold rolled steel and are mounted on full-extension slides. Both have 2<sup>3</sup>/<sub>8</sub>" high sides and back. Trays are open at the front while pans are enclosed on all four sides. Each is available in two configurations: one tray or pan mounted at bottom of cupboard; or

two trays or pans mounted at bottom and in center of cupboard. When pullout trays or pans are added, the adjustable shelf is omitted.


Pullout Trays and Pans are available for base cabinets only.

To specify a Pullout Tray, add **P1** or **P2** to the end of the cabinet number.


To specify a Pullout Pan, add **P3** or **P4** to the end of the cabinet number.

Example: **E40C352236-0101-P2**  
(Available for base cabinets only.)

Example: **E40C352236-0101-P3**  
(Available for base cabinets only.)



Pullout Tray



Pullout Pan

### Laminated Safety Glass

### Option I

### Tempered Safety Glass

### Option R

Option **I** replaces the standard float glass in plate glass or glazed door wall, counter mounted and full height cabinets with 1/4" Laminated Safety Glass.

Option **R** replaces the standard 1/8" float glass in glazed door wall, counter mounted and full height cabinets with 1/8" Tempered Safety Glass.

To specify Laminated Safety Glass, add **I** to the end of the cabinet number.  
Example: **W20C351336-0101-I**

To specify Tempered Glass, add **R** to the end of the cabinet number.  
Example: **W20C351336-0101-R**

### Full Extension Heavy Duty Slides


### Option T

Full extension ball bearing drawer slide with cushioned in and outstops and positive drawer hold-in detent. Loading capacity of 200 lbs. Replaces standard 100 lbs full-

extension slides on all drawers in cabinet when **T** option is specified.

To specify Heavy Duty Slides, add **T** to the end of the cabinet number.

Example: **E40C352236-1011-T**



### Toe Space on Cabinet Ends

### Option U

When cabinets are used at the end of center benches and in some other special applications, they may require a toe space on one or

both ends of the cabinets. These toe spaces match the toe space at the front of the cabinets and may be added to the left, right, or both ends.

- U1** – Toe Space on Left End
- U2** – Toe Space on Right End
- U3** – Toe Space on Both Ends

To specify End Toe Spaces, add **U\_** to the end of the cabinet number.

Example: **E40C352236-0101-U2**  
(Available for base cabinets only.)


# Steel Cabinet Options

## Locks with Lock Panels Locks without Lock Panels

**Option J**  
**Option K**

All cabinets may be ordered with locks by specifying the K option. In addition, base cabinets may be specified with the J option to add lock panels. Lock panels are

installed between drawers and/or drawers and cupboards to provide additional security when the locks on a cabinet are keyed different.



**Lock Panel**

### Locks Types

Cabinet locks are available in two types. Disk Tumbler which offers medium security, and Pin Tumbler which offer a higher level of security. Pin Tumbler also allows more keying arrangements. Both are fabricated with die cast sleeves and cylinders

with a satin nickel face and zinc plated steel cams. Both locks use the same sleeve and cam assembly allowing the lock cylinder to be removed and switched after the lock has been installed. All locks are furnished with two keys unless otherwise specified.

### Disk Tumbler Locks

**Security:** Medium

**Keying:** Different, Alike, Alike in Sets, Masterkeyed

**Key Changes:** Up to 2000

**Master Key Changes:** Up to 5

**J** – Disk Locks with Lock Panels  
**K** – Disk Locks w/o Lock Panels

### Pin Tumbler Locks


**Security:** Medium-High

**Keying:** Different, Alike, Alike in Sets, Masterkeyed, Grand-masterkeyed, Sub-masterkeyed, Mason Keyed


**Key Changes:** Up to 2000

**Master Key Changes:** Up to 52

**J2** – Pin Locks with Lock Panels  
**K2** – Pin Locks w/o Lock Panels



**Disk Tumbler Lock**  
(Medium Security)



**Pin Tumbler Lock**  
(Higher Security)

**Lock schedule with keying arrangement must be furnished when locks are specified.** See page 9 & 10 for keying arrangement diagrams, and page 104 & 105 for sample lock schedules.

To specify Locks without Lock Panels, add **K** or **K2** to the end of the cabinet number.

Example: **E40C352236-0101K2**  
(Added to all cupboards and drawers.)

To specify Locks with Lock Panels, add **J** or **J2** to the end of the cabinet number.

Example: **E40C352236-0101J**  
(Added to all cupboards and drawers.)

## Concealed Slide Drawers

**Option Z**

Overlay cabinets may be upgraded to a Concealed Slide Drawer with integral, concealed, full-extension, soft-close, slides. Concealed Slide Drawers are available with either powder-coated

steel drawer sides or with stainless steel drawer sides.

Concealed Slide Drawers are available in Overlay cabinets only and are not available for all drawer sizes.

To specify Concealed Slide Drawers with Powder-coated Steel sides, add **Z1** to the end of the cabinet number.

Example: **E40C352236-1101Z1**  
(Added to all drawers.)

To specify Concealed Slide Drawers with Stainless Steel Sides, add **Z2** to the end of the cabinet number.


Example: **E40C352236-1101Z2**  
(Added to all drawers.)

# Lock Keying Options

## Keyed Different

- Disk Tumbler Locks:** Available
- Pin Tumbler Locks:** Available


All locks are keyed different from each other with individual keys for each lock. No key will open more than one lock.



## Keyed Alike

- Disk Tumbler Locks:** Available
- Pin Tumbler Locks:** Available


All locks are keyed alike. One key will open all locks.



## Keyed Alike in Sets

- Disk Tumbler Locks:** Available
- Pin Tumbler Locks:** Available


Locks are grouped into sets, with as many locks as needed in each set. All locks in a set can be opened by a single key, but that key can open no other locks. No key will open all locks.



## Master Keyed

- Disk Tumbler Locks:** Available
- Pin Tumbler Locks:** Available

All locks are keyed different from each other with individual keys for each lock. Master key provided that can open all locks in the set. Multiple sets are available, but no key can open all locks in all sets.






# Lock Keying Options

## Grand Master Keyed


- Disk Tumbler Locks:** N/A
- Pin Tumbler Locks:** Available

All locks are keyed different from each other with individual keys for each lock. Master key provided that can open all locks in the set. Multiple sets are available, with Grand master key that can open all locks in all sets.



## Sub-master Keyed & Mason Keyed

- Disk Tumbler Locks:** N/A
- Pin Tumbler Locks:** Available




### Sub-master Keyed

All locks are keyed different from each other with individual keys for each lock. Sub-master key provided that can open all locks in sub-set. Multiple sub-sets are available, with master key that can open all locks in all sub-sets. Multiple sets are available, with grand master key that can open all locks in all sets and sub-sets.

### Mason Keyed

Locks are keyed different from each other with individual keys for each lock. Each individual key can also open all common locks in set. Master key provided that can open all locks in set. Multiple sets are available with grand master key that can open all locks in all sets.

# Steel Sitting Height Base Cabinets





## Specifications:

Sitting Height Cabinets are 29" high and 22" deep. Lengths as shown.

Open and Cupboard Cabinets are furnished with one adjustable shelf.

**Note:** Last four digits of part number designate door and drawer head style, material, and pull selection. See page 4 for full explanation. Actual product may vary slightly from catalog illustration due to style and hardware selection.

## Open Cabinets






<b>B00C292212-00_0</b>	12"L	<b>B00C292230-00_0</b>	30"L
<b>B00C292215-00_0</b>	15"L	<b>B00C292236-00_0</b>	36"L
<b>B00C292218-00_0</b>	18"L	<b>B00C292242-00_0</b>	42"L
<b>B00C292221-00_0</b>	21"L	<b>B00C292248-00_0</b>	48"L
<b>B00C292224-00_0</b>	24"L		

Open Cabinets may be ordered with the topspace 7/8" deeper so that cabinet fronts line up with the doors and drawer fronts on adjacent Full Overlay cabinets. To specify the deeper toe space add a **C** to the end of the part number.




Example: **B00C292236-0000-C**

## Cupboard Cabinets







<b>C03C292215L</b> _____	15"L	<b>C03C292215-</b> _____	15"L	<b>C03C292230-</b> _____	30"L
<b>C03C292218L</b> _____	18"L	<b>C03C292218-</b> _____	18"L	<b>C03C292236-</b> _____	36"L
<b>C03C292221L</b> _____	21"L	<b>C03C292221-</b> _____	21"L	<b>C03C292242-</b> _____	42"L
<b>C03C292224L</b> _____	24"L	<b>C03C292224-</b> _____	24"L	<b>C03C292248-</b> _____	48"L

## Drawer Cabinets




<b>D00C292215-</b> _____	15"L	<b>D00C292230-</b> _____	30"L
<b>D00C292218-</b> _____	18"L	<b>D00C292236-</b> _____	36"L
<b>D00C292221-</b> _____	21"L	<b>D00C292242-</b> _____	42"L
<b>D00C292224-</b> _____	24"L	<b>D00C292248-</b> _____	48"L

**D26C292218-**\_\_\_\_\_ 18"L  
1 Pull Out Shelf



<b>D30C292215-</b> _____	15"L	<b>D30C292230-</b> _____	30"L	<b>D31C292230-</b> _____	30"L	<b>D32C292230-</b> _____	30"L
<b>D30C292218-</b> _____	18"L	<b>D30C292236-</b> _____	36"L	<b>D31C292236-</b> _____	36"L	<b>D32C292236-</b> _____	36"L
<b>D30C292221-</b> _____	21"L	<b>D30C292242-</b> _____	42"L	<b>D31C292242-</b> _____	42"L	<b>D32C292242-</b> _____	42"L
<b>D30C292224-</b> _____	24"L	<b>D30C292248-</b> _____	48"L	<b>D31C292248-</b> _____	48"L	<b>D32C292248-</b> _____	48"L

# Steel Sitting Height Base Cabinets




## Specifications:

Sitting Height Cabinets are 29" high and 22" deep. Lengths as shown.


**Note:** Last four digits of part number designate door and drawer head style, material, and pull selection. See page 4 for full explanation. Actual product may vary slightly from catalog illustration due to style and hardware selection.

## Drawer Cabinets




**D81C292218-** 18"L  
**D81C292221-** 21"L  
**D81C292224-** 24"L

2 File Drawer with hangers



**D84C292218-** 18"L  
**D84C292221-** 21"L  
**D84C292224-** 24"L


Pull Out Shelf,  
1 File Drawer with hangers




**D85C292218-** 18"L  
**D85C292221-** 21"L  
**D85C292224-** 24"L

1 File Drawer with hangers


## Cupboard and Drawer Cabinets




**E05C292230L** 30"L  
**E05C292236L** 36"L  
**E05C292242L** 42"L  
**E05C292248L** 48"L




**E05C292230-** 30"L  
**E05C292236-** 36"L  
**E05C292242-** 42"L  
**E05C292248-** 48"L



**E40C292215L** 15"L  
**E40C292218L** 18"L  
**E40C292221L** 21"L  
**E40C292224L** 24"L




**E40C292215-** 15"L  
**E40C292218-** 18"L  
**E40C292221-** 21"L  
**E40C292224-** 24"L



**E41C292230-** 30"L  
**E41C292236-** 36"L  
**E41C292242-** 42"L  
**E41C292248-** 48"L

# Steel Sitting Height Base Cabinets





### Specifications:

Sitting Height Cabinets are 29" high and 22" deep. Lengths as shown.

Sink cabinets are provided with lowered back for plumbing access. Sink support kits are available when needed to support under-mount sinks.




**Note:** Last four digits of part number designate door and drawer head style, material, and pull selection. See page 4 for full explanation. Actual product may vary slightly from catalog illustration due to style and hardware selection.

### Cupboard and Drawer Cabinets




<b>E45C292230L</b> -----	30"L	<b>E45C292230-</b> -----	30"L
<b>E45C292236L</b> -----	36"L	<b>E45C292236-</b> -----	36"L
<b>E45C292242L</b> -----	42"L	<b>E45C292242-</b> -----	42"L
<b>E45C292248L</b> -----	48"L	<b>E45C292248-</b> -----	48"L

### Sink Cabinets





<b>G00C292215L</b> -----	15"L	<b>G00C292215-</b> -----	15"L	<b>G00C292230-</b> -----	30"L
<b>G00C292218L</b> -----	18"L	<b>G00C292218-</b> -----	18"L	<b>G00C292236-</b> -----	36"L
<b>G00C292221L</b> -----	21"L	<b>G00C292221-</b> -----	21"L	<b>G00C292242-</b> -----	42"L
<b>G00C292224L</b> -----	24"L	<b>G00C292224-</b> -----	24"L	<b>G00C292248-</b> -----	48"L

### Sink Support Kits



### Diagonal Corner Cabinets




**G14C292214-**-----

Part Number	Fits
<b>K-MSS-15</b>	15" Cabinet
<b>K-MSS-18</b>	18" Cabinet
<b>K-MSS-24</b>	24" Cabinet
<b>K-MSS-30</b>	30" Cabinet
<b>K-MSS-36</b>	36" Cabinet
<b>K-MSS-42</b>	42" Cabinet
<b>K-MSS-48</b>	48" Cabinet

Sink support kit includes: four (4) threaded hanging rods and two (2) support rails and are designed to support under-mount tub sinks up to 14" deep. For sinks from 14" to 20" deep add a **-D** to the end of the part number.

Example: **K-MSS-36-D**

# Steel ADA Height Base Cabinets




## Specifications:

ADA Height Cabinets are 32" high and 22" deep. Lengths as shown.


Open and Cupboard Cabinets are furnished with one adjustable shelf.

**Note:** Last four digits of part number designate door and drawer head style, material, and pull selection. See page 4 for full explanation. Actual product may vary slightly from catalog illustration due to style and hardware selection.

## Open Cabinets



**B00C322212-00\_0** 12"L  
**B00C322215-00\_0** 15"L  
**B00C322218-00\_0** 18"L  
**B00C322221-00\_0** 21"L  
**B00C322224-00\_0** 24"L




**B00C322230-00\_0** 30"L  
**B00C322236-00\_0** 36"L  
**B00C322242-00\_0** 42"L  
**B00C322248-00\_0** 48"L


Open Cabinets may be ordered with the topspace 7/8" deeper so that cabinet fronts line up with the doors and drawer fronts on adjacent Full Overlay cabinets. To specify the deeper toe space add a **C** to the end of the part number.

Example: **G00C322236-0000-C**


## Cupboard Cabinets



**C03C322215L** 15"L  
**C03C322218L** 18"L  
**C03C322221L** 21"L  
**C03C322224L** 24"L




**C03C322215-** 15"L  
**C03C322218-** 18"L  
**C03C322221-** 21"L  
**C03C322224-** 24"L



**C03C322230-** 30"L  
**C03C322236-** 36"L  
**C03C322242-** 42"L  
**C03C322248-** 48"L

# Steel ADA Height Base Cabinets




### Specifications:


ADA Height Cabinets are 32" high and 22" deep. Lengths as shown.

**Note:** Last four digits of part number designate door and drawer head style, material, and pull selection. See page 4 for full explanation. Actual product may vary slightly from catalog illustration due to style and hardware selection.


### Drawer Cabinets




**D00C322215-** 15"  
**D00C322218-** 18"  
**D00C322221-** 21"  
**D00C322224-** 24"




**D00C322230-** 30"  
**D00C322236-** 36"  
**D00C322242-** 42"  
**D00C322248-** 48"




**D02C322230-** 30"  
**D02C322236-** 36"  
**D02C322242-** 42"  
**D02C322248-** 48"




**D30C322215-** 15"  
**D30C322218-** 18"  
**D30C322221-** 21"  
**D30C322224-** 24"




**D30C322230-** 30"  
**D30C322236-** 36"  
**D30C322242-** 42"  
**D30C322248-** 48"




**D31C322230-** 30"  
**D31C322236-** 36"  
**D31C322242-** 42"  
**D31C322248-** 48"




**D32C322230-** 30"  
**D32C322236-** 36"  
**D32C322242-** 42"  
**D32C322248-** 48"




**D36C322230-** 30"  
**D36C322236-** 36"  
**D36C322242-** 42"  
**D36C322248-** 48"




**D50C322215-** 15"  
**D50C322218-** 18"  
**D50C322221-** 21"  
**D50C322224-** 24"




**D51C322230-** 30"  
**D51C322236-** 36"  
**D51C322242-** 42"  
**D51C322248-** 48"



**D52C322230-** 30"  
**D52C322236-** 36"  
**D52C322242-** 42"  
**D52C322248-** 48"




**D83C322215-** 15"  
**D83C322218-** 18"  
**D83C322221-** 21"  
**D83C322224-** 24"



**D84C322215-** 15"  
**D84C322218-** 18"  
**D84C322221-** 21"  
**D84C322224-** 24"

# Steel ADA Height Base Cabinets




## Specifications:


ADA Height Cabinets are 32" high and 22" deep. Lengths as shown.

**Note:** Last four digits of part number designate door and drawer head style, material, and pull selection. See page 4 for full explanation. Actual product may vary slightly from catalog illustration due to style and hardware selection.


## Cupboard and Drawer Cabinets




**E05C322230L** \_\_\_\_\_ 30"L  
**E05C322236L** \_\_\_\_\_ 36"L  
**E05C322242L** \_\_\_\_\_ 42"L  
**E05C322248L** \_\_\_\_\_ 48"L




**E05C322230-** \_\_\_\_\_ 30"L  
**E05C322236-** \_\_\_\_\_ 36"L  
**E05C322242-** \_\_\_\_\_ 42"L  
**E05C322248-** \_\_\_\_\_ 48"L




**E40C322215L** \_\_\_\_\_ 15"L  
**E40C322218L** \_\_\_\_\_ 18"L  
**E40C322221L** \_\_\_\_\_ 21"L  
**E40C322224L** \_\_\_\_\_ 24"L




**E40C322215-** \_\_\_\_\_ 15"L  
**E40C322218-** \_\_\_\_\_ 18"L  
**E40C322221-** \_\_\_\_\_ 21"L  
**E40C322224-** \_\_\_\_\_ 24"L




**E40C322230-** \_\_\_\_\_ 30"L  
**E40C322236-** \_\_\_\_\_ 36"L  
**E40C322242-** \_\_\_\_\_ 42"L  
**E40C322248-** \_\_\_\_\_ 48"L




**E41C322230-** \_\_\_\_\_ 30"L  
**E41C322236-** \_\_\_\_\_ 36"L  
**E41C322242-** \_\_\_\_\_ 42"L  
**E41C322248-** \_\_\_\_\_ 48"L




**E42C322230L** \_\_\_\_\_ 30"L  
**E42C322236L** \_\_\_\_\_ 36"L  
**E42C322242L** \_\_\_\_\_ 42"L  
**E42C322248L** \_\_\_\_\_ 48"L




**E42C322230-** \_\_\_\_\_ 30"L  
**E42C322236-** \_\_\_\_\_ 36"L  
**E42C322242-** \_\_\_\_\_ 42"L  
**E42C322248-** \_\_\_\_\_ 48"L



**E43C322230L** \_\_\_\_\_ 30"L  
**E43C322236L** \_\_\_\_\_ 36"L  
**E43C322242L** \_\_\_\_\_ 42"L  
**E43C322248L** \_\_\_\_\_ 48"L




**E43C322230-** \_\_\_\_\_ 30"L  
**E43C322236-** \_\_\_\_\_ 36"L  
**E43C322242-** \_\_\_\_\_ 42"L  
**E43C322248-** \_\_\_\_\_ 48"L



**E45C322230L** \_\_\_\_\_ 30"L  
**E45C322236L** \_\_\_\_\_ 36"L  
**E45C322242L** \_\_\_\_\_ 42"L  
**E45C322248L** \_\_\_\_\_ 48"L




**E45C322230-** \_\_\_\_\_ 30"L  
**E45C322236-** \_\_\_\_\_ 36"L  
**E45C322242-** \_\_\_\_\_ 42"L  
**E45C322248-** \_\_\_\_\_ 48"L



**E56C322230-** \_\_\_\_\_ 30"L  
**E56C322236-** \_\_\_\_\_ 36"L  
**E56C322242-** \_\_\_\_\_ 42"L  
**E56C322248-** \_\_\_\_\_ 48"L

# Steel ADA Height Base Cabinets




**Specifications:**


ADA Height Cabinets are 32" high and 22" deep. Lengths as shown.

**Note:** Last four digits of part number designate door and drawer head style, material, and pull selection. See page 4 for full explanation. Actual product may vary slightly from catalog illustration due to style and hardware selection.


**Cupboard and Drawer Cabinets**




- E62C322230L \_\_\_\_\_ 30"L
- E62C322236L \_\_\_\_\_ 36"L
- E62C322242L \_\_\_\_\_ 42"L
- E62C322248L \_\_\_\_\_ 48"L




- E62C322230- \_\_\_\_\_ 30"L
- E62C322236- \_\_\_\_\_ 36"L
- E62C322242- \_\_\_\_\_ 42"L
- E62C322248- \_\_\_\_\_ 48"L




- E63C322230L \_\_\_\_\_ 30"L
- E63C322236L \_\_\_\_\_ 36"L
- E63C322242L \_\_\_\_\_ 42"L
- E63C322248L \_\_\_\_\_ 48"L



- E63C322230- \_\_\_\_\_ 30"L
- E63C322236- \_\_\_\_\_ 36"L
- E63C322242- \_\_\_\_\_ 42"L
- E63C322248- \_\_\_\_\_ 48"L




- E65C322230L \_\_\_\_\_ 30"L
- E65C322236L \_\_\_\_\_ 36"L
- E65C322242L \_\_\_\_\_ 42"L
- E65C322248L \_\_\_\_\_ 48"L



- E65C322230- \_\_\_\_\_ 30"L
- E65C322236- \_\_\_\_\_ 36"L
- E65C322242- \_\_\_\_\_ 42"L
- E65C322248- \_\_\_\_\_ 48"L



# Steel ADA Height Base Cabinets




**Specifications:**


ADA Height Cabinets are 32" high and 22" deep. Lengths as shown.

**Note:** Last four digits of part number designate door and drawer head style, material, and pull selection. See page 4 for full explanation. Actual product may vary slightly from catalog illustration due to style and hardware selection.


**Blind Corner Cabinets**



G12C322248-\_\_\_\_\_ 48"L




PLAN VIEW OF  
G12C322248- & G13C322248-




G13C322248-\_\_\_\_\_ 48"L

**Acid Storage Cabinets**




20" Deep

G21C322018L \_\_\_\_\_ 18"L  
G21C322024L \_\_\_\_\_ 24"L



20" Deep

G21C320218- \_\_\_\_\_ 18"L  
G21C320224- \_\_\_\_\_ 24"L



20" Deep

G21C322030- \_\_\_\_\_ 30"L  
G21C322036- \_\_\_\_\_ 36"L  
G21C322042- \_\_\_\_\_ 42"L  
G21C322048- \_\_\_\_\_ 48"L

22" Deep

G21C322218L \_\_\_\_\_ 18"L  
G21C322224L \_\_\_\_\_ 24"L

22" Deep


G21C322218- \_\_\_\_\_ 18"L  
G21C322224- \_\_\_\_\_ 24"L

22" Deep


G21C322230- \_\_\_\_\_ 30"L  
G21C322236- \_\_\_\_\_ 36"L  
G21C322242- \_\_\_\_\_ 42"L  
G21C322248- \_\_\_\_\_ 48"L

**Note:** G21 Acid Storage Cabinets are available with five-knuckle hinges only


**Diagonal Corner Cabinet**




G14C322214-\_\_\_\_\_



G15C322214-\_\_\_\_\_



# Steel ADA Height Base Cabinets




### Specifications:

ADA Height Cabinets are 32" high and 22" deep. Lengths as shown.


Sink cabinets are provided with lowered back for plumbing access. Sink support kits are available when needed to support under-mount sinks.

**Note:** Last four digits of part number designate door and drawer head style, material, and pull selection. See page 4 for full explanation. Actual product may vary slightly from catalog illustration due to style and hardware selection.


### Sink Cabinets



- G00C322212L \_\_\_\_\_ 12"L
- G00C322215L \_\_\_\_\_ 15"L
- G00C322218L \_\_\_\_\_ 18"L
- G00C322221L \_\_\_\_\_ 21"L
- G00C322224L \_\_\_\_\_ 24"L




- G00C322212- \_\_\_\_\_ 12"L
- G00C322215- \_\_\_\_\_ 15"L
- G00C322218- \_\_\_\_\_ 18"L
- G00C322221- \_\_\_\_\_ 21"L
- G00C322224- \_\_\_\_\_ 24"L



- G00C322230- \_\_\_\_\_ 30"L
- G00C322236- \_\_\_\_\_ 36"L
- G00C322242- \_\_\_\_\_ 42"L
- G00C322248- \_\_\_\_\_ 48"L

### Sink Support Kits




Part Number	Fits
K-MSS-12	12" Cabinet
K-MSS-15	15" Cabinet
K-MSS-18	18" Cabinet
K-MSS-24	24" Cabinet
K-MSS-30	30" Cabinet
K-MSS-36	36" Cabinet
K-MSS-42	42" Cabinet
K-MSS-48	48" Cabinet


Sink support kit includes: four (4) threaded hanging rods and two (2) support rails and are designed to support under-mount tub sinks up to 14" deep. For sinks from 14" to 20" deep add a **-D** to the end of the part number.

Example: **K-MSS-36-D**


### Solvent Storage Cabinets



- 20" Deep**
- G63C322018L\_100 \_\_\_\_\_ 18"L
  - G63C322024L\_100 \_\_\_\_\_ 24"L



- 20" Deep**
- G63C322018-100 \_\_\_\_\_ 18"L
  - G63C322024-100 \_\_\_\_\_ 24"L



- 20" Deep**
- G63C322030-100 \_\_\_\_\_ 30"L
  - G63C322036-100 \_\_\_\_\_ 36"L
  - G63C322048-100 \_\_\_\_\_ 48"L

- 22" Deep**
- G63C322218L\_100 \_\_\_\_\_ 18"L
  - G63C322224L\_100 \_\_\_\_\_ 24"L


- 22" Deep**
- G63C322218-100 \_\_\_\_\_ 18"L
  - G63C322224-100 \_\_\_\_\_ 24"L

- 22" Deep**
- G63C322230-100 \_\_\_\_\_ 30"L
  - G63C322236-100 \_\_\_\_\_ 36"L
  - G63C322248-100 \_\_\_\_\_ 48"L

G63 Solvent Storage Cabinets are available as mobile.

To add casters to the cabinets above, add **-B7** to the end of the part number.  
Example: **G63C322230-0100-B7** — Mobile, Inset style cabinet, 22" deep by 30" long.


### Bin Storage Cabinet



- G23C322218- \_\_\_\_\_ 18"L

Furnished with 11" high drawer body with full height drawer front mounted on full extension slides.

# Steel ADA Height Base Cabinets




## Specifications:

ADA Height Fume Hood Cabinets are 32" high and either 20" or 22" deep. Lengths as shown.

**Note:** Last four digits of part number designate door and drawer head style, material, and pull selection. See page 4 for full explanation. Actual product may vary slightly from catalog illustration due to style and hardware selection.

## Fume Hood Base Cabinets




20" Deep

G08C322012L	12"L
G08C322015L	15"L
G08C322018L	18"L
G08C322021L	21"L
G08C322024L	24"L

22" Deep

G08C322212L	12"L
G08C322215L	15"L
G08C322218L	18"L
G08C322221L	21"L
G08C322224L	24"L




20" Deep

G08C322012-	12"L
G08C322015-	15"L
G08C322018-	18"L
G08C322021-	21"L
G08C322024-	24"L

22" Deep

G08C322212-	12"L
G08C322215-	15"L
G08C322218-	18"L
G08C322221-	21"L
G08C322224-	24"L




20" Deep

G08C322030-	30"L
G08C322036-	36"L
G08C322042-	42"L
G08C322048-	48"L

22" Deep

G08C322230-	30"L
G08C322236-	36"L
G08C322242-	42"L
G08C322248-	48"L

## Fume Hood Solvent Storage Base Cabinets




20" Deep

G68C322024L_100	24"L
-----------------	------

22" Deep

G68C322224L_100	24"L
-----------------	------




20" Deep

G68C322024-_100	24"L
-----------------	------

22" Deep

G68C322224-_100	24"L
-----------------	------



20" Deep


G68C322030-_100	30"L
G68C322036-_100	36"L
G68C322048-_100	48"L

22" Deep

G68C322230-_100	30"L
G68C322236-_100	36"L
G68C322248-_100	48"L

G68 Solvent Storage Cabinets are available as mobile. To add casters to the cabinets above, add **-B7** to the end of the part number. Example: **G68C322018-0100-B7** — Mobile, Inset style cabinet, 22" deep by 18" long.

## Fume Hood Vacuum Storage Base Cabinets




20" Deep

G35C322024L	24"L
-------------	------

22" Deep

G35C322224L	24"L
-------------	------




20" Deep

G35C322024-	24"L
-------------	------

22" Deep

G35C322224-	24"L
-------------	------



20" Deep


G35C322030-	30"L
G35C322036-	36"L
G35C322048-	48"L

22" Deep

G35C322230-	30"L
G35C322236-	36"L
G35C322248-	48"L


# Steel ADA Height Base Cabinets

## Fume Hood Acid Storage Base Cabinets




20" Deep

<b>G80C322018L</b> _____	24"L
<b>G80C322024L</b> _____	24"L



20" Deep

<b>G80C322018-</b> _____	18"L
<b>G80C322024-</b> _____	24"L



20" Deep

<b>G80C322030-</b> _____	30"L
<b>G80C322036-</b> _____	36"L
<b>G80C322048-</b> _____	48"L

22" Deep

<b>G80C322218L</b> _____	18"L
<b>G80C322224L</b> _____	24"L

22" Deep


<b>G80C322218-</b> _____	18"L
<b>G80C322224-</b> _____	24"L

22" Deep

<b>G80C322230-</b> _____	30"L
<b>G80C322236-</b> _____	36"L
<b>G80C322248-</b> _____	48"L


**Note:** G80 Acid Storage Cabinets are available with five-knuckle hinges only

## Wheelchair Accessible Sink Cabinet with Doors




The Wheelchair Accessible Base Cabinet with Doors is designed to be used as a sink unit, but may be used anywhere wheelchair accessibility is needed. The cabinet is constructed without a bottom, and with the toe-

space rails attached to the doors so it will match adjacent cabinets when the doors are closed. The cabinet incorporates a 2<sup>3</sup>/<sub>4</sub>" top rail to give a minimum of 29<sup>1</sup>/<sub>4</sub>" of clearance from the floor to the bottom of the rail.



<b>G03C322236-</b> _____	36"L
<b>G03C322242-</b> _____	42"L
<b>G03C322248-</b> _____	48"L

## Wheelchair Accessible Base Cabinet without Doors




The Wheelchair Accessible Base Cabinet without Doors is designed to be used as a sink unit, but may be used anywhere wheelchair accessibility is needed. The cabinet is constructed without a bottom, and without doors for easy accessibility.

The cabinet incorporates a 4<sup>1</sup>/<sub>2</sub>" top rail to give a minimum of 27<sup>1</sup>/<sub>2</sub>" of clearance from the floor to the bottom of the rail and a removable sloping knee space panel to enclose piping and drain lines.


Wheelchair Accessible Base Cabinets without Doors may be ordered with the topspace 7/8" deeper so that cabinet fronts line up with the doors and drawer fronts on adjacent Full Overlay cabinets. To specify the deeper toe space, add a "C" to the end of the part number.

Example: **G07C322236-0000-C**



<b>G07C322236-00_0</b>	36"L
<b>G07C322242-00_0</b>	42"L
<b>G07C322248-00_0</b>	48"L

# Steel Standing Height Base Cabinets




**Specifications:**

Standing Height Cabinets are 35" high and 22" deep. Lengths as shown.


Open and Cupboard Cabinets are furnished with one adjustable shelf.

**Note:** Last four digits of part number designate door and drawer head style, material, and pull selection. See page 4 for full explanation. Actual product may vary slightly from catalog illustration due to style and hardware selection.

**Open Cabinets**



**B00C352212-00\_0** 12"L  
**B00C352215-00\_0** 15"L  
**B00C352218-00\_0** 18"L  
**B00C352221-00\_0** 21"L  
**B00C352224-00\_0** 24"L




**B00C352230-00\_0** 30"L  
**B00C352236-00\_0** 36"L  
**B00C352242-00\_0** 42"L  
**B00C352248-00\_0** 48"L


Open Cabinets may be ordered with the topspace 7/8" deeper so that cabinet fronts line up with the doors and drawer fronts on adjacent Full Overlay cabinets. To specify the deeper toe space add a **C** to the end of the part number.

Example: **G00C352236-0000-C**


**Cupboard Cabinets**



**C03C352212L** 12"L  
**C03C352215L** 15"L  
**C03C352218L** 18"L  
**C03C352221L** 21"L  
**C03C352224L** 24"L




**C03C352212-** 12"L  
**C03C352215-** 15"L  
**C03C352218-** 18"L  
**C03C352221-** 21"L  
**C03C352224-** 24"L




**C03C352230-** 30"L  
**C03C352236-** 36"L  
**C03C352242-** 42"L  
**C03C352248-** 48"L


**Drawer Cabinets**




**D00C352215-** 15"L  
**D00C352218-** 18"L  
**D00C352221-** 21"L  
**D00C352224-** 24"L




**D00C352230-** 30"L  
**D00C352236-** 36"L  
**D00C352242-** 42"L  
**D00C352248-** 48"L



**D02C352230-** 30"L  
**D02C352236-** 36"L  
**D02C352242-** 42"L  
**D02C352248-** 48"L




**D08C352215-** 15"L  
**D08C352218-** 18"L  
**D08C352221-** 21"L  
**D08C352224-** 24"L



**D08C352230-** 30"L  
**D08C352236-** 36"L  
**D08C352242-** 42"L  
**D08C352248-** 48"L

# Steel Standing Height Base Cabinets




**Specifications:**


Standing Height Cabinets are 35" high and 22" deep. Lengths as shown.

**Note:** Last four digits of part number designate door and drawer head style, material, and pull selection. See page 4 for full explanation. Actual product may vary slightly from catalog illustration due to style and hardware selection.


**Drawer Cabinets**




- D30C352212- 12"L
- D30C352215- 15"L
- D30C352218- 18"L
- D30C352221- 21"L
- D30C352224- 24"L




- D30C352230- 30"L
- D30C352236- 36"L
- D30C352242- 42"L
- D30C352248- 48"L




- D31C352230- 30"L
- D31C352236- 36"L
- D31C352242- 42"L
- D31C352248- 48"L




- D32C352230- 30"L
- D32C352236- 36"L
- D32C352242- 42"L
- D32C352248- 48"L




- D36C352230- 30"L
- D36C352236- 36"L
- D36C352242- 42"L
- D36C352248- 48"L




- D50C352215- 15"L
- D50C352218- 18"L
- D50C352221- 21"L
- D50C352224- 24"L




- D51C352230- 30"L
- D51C352236- 36"L
- D51C352242- 42"L
- D51C352248- 48"L



- D52C352230- 30"L
- D52C352236- 36"L
- D52C352242- 42"L
- D52C352248- 48"L



- D81C352218- 18"L
- D81C352221- 21"L
- D81C352224- 24"L




- D84C352218- 18"L
- D84C352221- 21"L
- D84C352224- 24"L

2 File Drawers with hangers

Pull Out Shelf  
1 File Drawer with hangers

# Steel Standing Height Base Cabinets




**Specifications:**

Standing Height Cabinets are 35" high and 22" deep. Lengths as shown.


Cupboards have one adjustable shelf.

**Note:** Last four digits of part number designate door and drawer head style, material, and pull selection. See page 4 for full explanation. Actual product may vary slightly from catalog illustration due to style and hardware selection.


**Cupboard and Drawer Cabinets**




- E05C352230L \_\_\_\_\_ 30"L
- E05C352236L \_\_\_\_\_ 36"L
- E05C352242L \_\_\_\_\_ 42"L
- E05C352248L \_\_\_\_\_ 48"L




- E05C352230- \_\_\_\_\_ 30"L
- E05C352236- \_\_\_\_\_ 36"L
- E05C352242- \_\_\_\_\_ 42"L
- E05C352248- \_\_\_\_\_ 48"L




- E40C352212L \_\_\_\_\_ 12"L
- E40C352215L \_\_\_\_\_ 15"L
- E40C352218L \_\_\_\_\_ 18"L
- E40C352221L \_\_\_\_\_ 21"L
- E40C352224L \_\_\_\_\_ 24"L




- E40C352212- \_\_\_\_\_ 12"L
- E40C352215- \_\_\_\_\_ 15"L
- E40C352218- \_\_\_\_\_ 18"L
- E40C352221- \_\_\_\_\_ 21"L
- E40C352224- \_\_\_\_\_ 24"L




- E40C352230- \_\_\_\_\_ 30"L
- E40C352236- \_\_\_\_\_ 36"L
- E40C352242- \_\_\_\_\_ 42"L
- E40C352248- \_\_\_\_\_ 48"L




- E41C352230- \_\_\_\_\_ 30"L
- E41C352236- \_\_\_\_\_ 36"L
- E41C352242- \_\_\_\_\_ 42"L
- E41C352248- \_\_\_\_\_ 48"L




- E42C352230L \_\_\_\_\_ 30"L
- E42C352236L \_\_\_\_\_ 36"L
- E42C352242L \_\_\_\_\_ 42"L
- E42C352248L \_\_\_\_\_ 48"L




- E42C352230- \_\_\_\_\_ 30"L
- E42C352236- \_\_\_\_\_ 36"L
- E42C352242- \_\_\_\_\_ 42"L
- E42C352248- \_\_\_\_\_ 48"L




- E43C352230L \_\_\_\_\_ 30"L
- E43C352236L \_\_\_\_\_ 36"L
- E43C352242L \_\_\_\_\_ 42"L
- E43C352248L \_\_\_\_\_ 48"L




- E43C352230- \_\_\_\_\_ 30"L
- E43C352236- \_\_\_\_\_ 36"L
- E43C352242- \_\_\_\_\_ 42"L
- E43C352248- \_\_\_\_\_ 48"L



- E45C352230L \_\_\_\_\_ 30"L
- E45C352236L \_\_\_\_\_ 36"L
- E45C352242L \_\_\_\_\_ 42"L
- E45C352248L \_\_\_\_\_ 48"L




- E45C352230- \_\_\_\_\_ 30"L
- E45C352236- \_\_\_\_\_ 36"L
- E45C352242- \_\_\_\_\_ 42"L
- E45C352248- \_\_\_\_\_ 48"L



- E56C352230- \_\_\_\_\_ 30"L
- E56C352236- \_\_\_\_\_ 36"L
- E56C352242- \_\_\_\_\_ 42"L
- E56C352248- \_\_\_\_\_ 48"L

# Steel Standing Height Base Cabinets




### Specifications:

Standing Height Cabinets are 35" high and 22" deep. Lengths as shown.


Cupboards have one adjustable shelf.

**Note:** Last four digits of part number designate door and drawer head style, material, and pull selection. See page 4 for full explanation. Actual product may vary slightly from catalog illustration due to style and hardware selection.


### Cupboard and Drawer Cabinets




- E60C352212L \_\_\_\_\_ 12"L
- E60C352215L \_\_\_\_\_ 15"L
- E60C352218L \_\_\_\_\_ 18"L
- E60C352221L \_\_\_\_\_ 21"L
- E60C352224L \_\_\_\_\_ 24"L




- E60C352212- \_\_\_\_\_ 12"L
- E60C352215- \_\_\_\_\_ 15"L
- E60C352218- \_\_\_\_\_ 18"L
- E60C352221- \_\_\_\_\_ 21"L
- E60C352224- \_\_\_\_\_ 24"L




- E61C352230- \_\_\_\_\_ 30"L
- E61C352236- \_\_\_\_\_ 36"L
- E61C352242- \_\_\_\_\_ 42"L
- E61C352248- \_\_\_\_\_ 48"L




- E62C352230L \_\_\_\_\_ 30"L
- E62C352236L \_\_\_\_\_ 36"L
- E62C352242L \_\_\_\_\_ 42"L
- E62C352248L \_\_\_\_\_ 48"L




- E62C352230- \_\_\_\_\_ 30"L
- E62C352236- \_\_\_\_\_ 36"L
- E62C352242- \_\_\_\_\_ 42"L
- E62C352248- \_\_\_\_\_ 48"L




- E63C352230L \_\_\_\_\_ 30"L
- E63C352236L \_\_\_\_\_ 36"L
- E63C352242L \_\_\_\_\_ 42"L
- E63C352248L \_\_\_\_\_ 48"L



- E63C352230- \_\_\_\_\_ 30"L
- E63C352236- \_\_\_\_\_ 36"L
- E63C352242- \_\_\_\_\_ 42"L
- E63C352248- \_\_\_\_\_ 48"L




- E65C352230L \_\_\_\_\_ 30"L
- E65C352236L \_\_\_\_\_ 36"L
- E65C352242L \_\_\_\_\_ 42"L
- E65C352248L \_\_\_\_\_ 48"L



- E65C352230- \_\_\_\_\_ 30"L
- E65C352236- \_\_\_\_\_ 36"L
- E65C352242- \_\_\_\_\_ 42"L
- E65C352248- \_\_\_\_\_ 48"L



# Steel Standing Height Base Cabinets




## Specifications:

Standing Height Cabinets are 35" high and 22" deep. Lengths as shown.


Sink cabinets are provided with lowered back for plumbing access except where noted. Sink support kits are available when needed to support under-mount sinks.

**Note:** Last four digits of part number designate door and drawer head style, material, and pull selection. See page 4 for full explanation. Actual product may vary slightly from catalog illustration due to style and hardware selection.


## Sink Cabinets




<b>G00C352212L</b> -----	12"L
<b>G00C352215L</b> -----	15"L
<b>G00C352218L</b> -----	18"L
<b>G00C352221L</b> -----	21"L
<b>G00C352224L</b> -----	24"L




<b>G00C352212-</b> -----	12"L
<b>G00C352215-</b> -----	15"L
<b>G00C352218-</b> -----	18"L
<b>G00C352221-</b> -----	21"L
<b>G00C352224-</b> -----	24"L



<b>G00C352230-</b> -----	30"L
<b>G00C352236-</b> -----	36"L
<b>G00C352242-</b> -----	42"L
<b>G00C352248-</b> -----	48"L




<b>G00C352260-</b> -----	60"L
--------------------------	------



<b>G03C352230-</b> -----	30"L
<b>G03C352236-</b> -----	36"L
<b>G03C352242-</b> -----	42"L
<b>G03C352248-</b> -----	48"L

G03 Cabinets are constructed without a bottom and with the toe-space rail attached to the doors so it will match adjacent cabinets when the doors are closed.


Furnished with full height back



<b>G04C352252-</b> -----	52"L
<b>G04C352258-</b> -----	58"L
<b>G04C352260-</b> -----	60"L

Center Bench End Sink Cabinet  
Furnished with full height back

## Sink Support Kits



Part Number	Fits
<b>K-MSS-12</b>	12" Cabinet
<b>K-MSS-15</b>	15" Cabinet
<b>K-MSS-18</b>	18" Cabinet
<b>K-MSS-21</b>	21" Cabinet
<b>K-MSS-24</b>	24" Cabinet
<b>K-MSS-30</b>	30" Cabinet
<b>K-MSS-36</b>	36" Cabinet
<b>K-MSS-42</b>	42" Cabinet
<b>K-MSS-48</b>	48" Cabinet
<b>K-MSS-52</b>	52" Cabinet
<b>K-MSS-58</b>	58" Cabinet
<b>K-MSS-60</b>	60" Cabinet

Sink support kit includes: four (4) threaded hanging rods and two (2) support rails and are designed to support under-mount tub sinks up to 14" deep. For sinks from 14" to 20" deep add a **-D** to the end of the part number.

Example: **K-MSS-36-D**

# Steel Standing Height Base Cabinets


## Specifications

Acid Storage Base Cabinets are specifically designed for the storage of corrosive chemicals. These cabinets are lined with 1/8" thick corrosion resistant polyethylene sheets.

The cabinet doors are provided with louvers and are latched using a nylon roller catch.

An optional 1" deep polyethylene pan is available to contain spills.

Each cabinet is furnished with a PVC vent pipe for venting. (Requires a 2" hole in the counter top above.)




## Optional Polyethylene Pan

Seamless 1" deep polyethylene pan designed to contain spills.

Part Number	Fits 20" Deep
<b>F-8440-15</b>	15" Cabinet
<b>F-8440-18</b>	18" Cabinet
<b>F-8440-21</b>	21" Cabinet
<b>F-8440-24</b>	24" Cabinet
<b>F-8440-30</b>	30" Cabinet
<b>F-8440-36</b>	36" Cabinet
<b>F-8440-42</b>	42" Cabinet
<b>F-8440-48</b>	48" Cabinet


Part Number	Fits 22" Deep
<b>F-8441-15</b>	15" Cabinet
<b>F-8441-18</b>	18" Cabinet
<b>F-8441-21</b>	21" Cabinet
<b>F-8441-24</b>	24" Cabinet
<b>F-8441-30</b>	30" Cabinet
<b>F-8441-36</b>	36" Cabinet
<b>F-8441-42</b>	42" Cabinet
<b>F-8441-48</b>	48" Cabinet

## Acid Storage Cabinet




20" Deep

<b>G21C352015L</b>	15"L
<b>G21C352018L</b>	18"L
<b>G21C352021L</b>	21"L
<b>G21C352024L</b>	24"L



20" Deep

<b>G21C352015-</b>	15"L
<b>G21C352018-</b>	18"L
<b>G21C352021-</b>	21"L
<b>G21C352024-</b>	24"L



20" Deep

<b>G21C352030-</b>	30"L
<b>G21C352036-</b>	36"L
<b>G21C352042-</b>	42"L
<b>G21C352048-</b>	48"L

22" Deep

<b>G21C352215L</b>	15"L
<b>G21C352218L</b>	18"L
<b>G21C352221L</b>	21"L
<b>G21C352224L</b>	24"L

22" Deep

<b>G21C352215-</b>	15"L
<b>G21C352218-</b>	18"L
<b>G21C352221-</b>	21"L
<b>G21C352224-</b>	24"L

22" Deep

<b>G21C352230-</b>	30"L
<b>G21C352236-</b>	36"L
<b>G21C352242-</b>	42"L
<b>G21C352248-</b>	48"L

**Note:** G21 Acid Storage Cabinets are available with five-knuckle hinges only

# Steel Standing Height Base Cabinets

## Specifications

Solvent Storage Cabinets are UL listed for the storage of flammable and combustible liquids, meeting UFC, OSHA and NFPA No. 30-1993 standards.


A 2" deep, steel, liquid-tight pan covers the entire bottom to contain liquid leaks and spills. The cabinet is provided with a full depth adjustable shelf and two

diametrically opposed vents with spark screens in the back of the cabinet for cases when ventilation is required.


Steel Solvent Storage cabinets are all 18 gauge steel, double panel construction with self-closing doors, synchronized so that both doors will always fully close. The right hand door is equipped

with a three-point latching system that automatically engages the cabinet frame. Each door is equipped with a fusible-link hold-open feature that ensures the doors close should the temperature outside the cabinet exceed 165°F. Cabinets are provided with a grounding screw at the rear.


All Solvent Storage Cabinets are labeled: CAUTION FLAMMABLE – KEEP FIRE AWAY in English, Spanish, and French.




## Vent Locations




## Solvent Storage Cabinets



**20" Deep** Capacity  
**G63C352024L\_100** 30"L 12 gal



**20" Deep** Capacity  
**G63C352024-\_100** 30"L 12 gal



**20" Deep** Capacity  
**G63C352030-\_100** 30"L 12 gal  
**G63C352036-\_100** 36"L 15 gal  
**G63C352048-\_100** 48"L 20 gal


**22" Deep** Capacity  
**G63C352224L\_100** 30"L 12 gal

**22" Deep** Capacity  
**G63C352224-\_100** 30"L 12 gal

**22" Deep** Capacity  
**G63C352230-\_100** 30"L 12 gal  
**G63C352236-\_100** 36"L 15 gal  
**G63C352248-\_100** 48"L 20 gal

G63 Solvent Storage Cabinets are available as mobile.  
To add casters to the cabinets above, add **-B7** to the end of the part number.  
Example: **G63C352036-0100-B7** — Mobile, Inset style cabinet, 20" deep by 36" long.

# Steel Standing Height Base Cabinets




**Specifications:**

Standing Height Cabinets are 35" high and 22" deep. Lengths as shown.


Blind Corner Cabinets have one adjustable shelf.

**Note:** Last four digits of part number designate door and drawer head style, material, and pull selection. See page 4 for full explanation. Actual product may vary slightly from catalog illustration due to style and hardware selection.


**Blind Corner Cabinets**



G12C352248-\_\_\_\_\_ 48"L




PLAN VIEW OF  
G12C352248- & G13C352248-




G13C352248-\_\_\_\_\_ 48"L

**Diagonal Corner Cabinets**




G14C352214-\_\_\_\_\_




G15C352214-\_\_\_\_\_

Provided with Lazy Susan



**Waste Storage Cabinet      Bin Storage Cabinet**



- G20C352215-\_\_\_\_\_ 15"L
- G20C352218-\_\_\_\_\_ 18"L
- G20C352221-\_\_\_\_\_ 21"L
- G20C352224-\_\_\_\_\_ 21"L


Cabinet furnished with 13" deep, cold rolled steel bin, mounted on full extension drawer slides for ease of access. Includes 6" x 5" (opening) stainless steel waste chute



- G23C352218-\_\_\_\_\_ 18"L

Furnished with 11" high drawer body with full height drawer front mounted on full extension slides.

# Fume Hood Base Cabinets




## Specifications

Standing Height Steel Fume Hood Base Cabinets are 35" high and 22" deep; ADA Height Steel Fume Hood Base Cabinets are 32" high and 22" deep; lengths as shown.

## Cabinet Style Option:


Specify cabinet door and hardware style when ordering by replacing blanks in last four digits of catalog number with style numbers.

## Standing Height Fume Hood Base Cabinets




20" Deep

G08C352012L	12"L
G08C352015L	15"L
G08C352018L	18"L
G08C352021L	21"L
G08C352024L	24"L




20" Deep

G08C352012-	12"L
G08C352015-	15"L
G08C352018-	18"L
G08C352021-	21"L
G08C352024-	24"L



20" Deep

G08C352030-	30"L
G08C352036-	36"L
G08C352042-	42"L
G08C352048-	48"L



20" Deep

G08C352060-	60"L
-------------	------

22" Deep

G08C352212L	12"L
G08C352215L	15"L
G08C352218L	18"L
G08C352221L	21"L
G08C352224L	24"L

22" Deep

G08C352212-	12"L
G08C352215-	15"L
G08C352218-	18"L
G08C352221-	21"L
G08C352224-	24"L

22" Deep

G08C352230-	30"L
G08C352236-	36"L
G08C352242-	42"L
G08C352248-	48"L

22" Deep

G08C352260-	60"L
-------------	------


# Acid Storage Fume Hood Base Cabinets

## Specifications

Acid Storage Fume Hood Base Cabinets are specifically designed for the storage of corrosive chemicals. These cabinets are lined with a molded one piece linear low density polyethylene tub with coved corners

and a 1" lip at the bottom front. The cabinet doors are lined with 1/8" sheet polyethylene and the doors are latched using a nylon roller catch.


Each cabinet is furnished with a PVC vent pipe for venting to the fume hood above. (Requires a 2" hole in the hood worksurface.)




## Optional Phenolic Resin Removable Half-Depth Shelf

Part No.	
019623	Shelf for 18" cabinet
019624	Shelf for 24" cabinet
019625	Shelf for 30" cabinet
019626	Shelf for 36" cabinet
019627	Shelf for 48" cabinet


## Standing Height Acid Storage Fume Hood Cabinets




**Vent Location**  
On Center  
11<sup>5</sup>/<sub>8</sub>" from Top of Cabinet



**Vent Location**  
On Center  
11<sup>5</sup>/<sub>8</sub>" from Top of Cabinet



**Vent Location**  
30" - 7<sup>5</sup>/<sub>8</sub>" Right of Center  
36" & 48" - 5" Right of Center  
11<sup>5</sup>/<sub>8</sub>" from Top of Cabinet



**Vent Location**  
On Center  
11<sup>5</sup>/<sub>8</sub>" from Top of Cabinet

20" Deep	
G80C352018L	18"L
G80C352024L	24"L

20" Deep	
G80C352018-	18"L
G80C352024-	24"L

20" Deep	
G80C352030-	30"L
G80C352036-	36"L
G80C352048-	48"L

20" Deep	
G80C352060-	60"L

22" Deep	
G80C352218L	18"L
G80C352224L	24"L

22" Deep	
G80C352218-	18"L
G80C352224-	24"L

22" Deep	
G80C352230-	30"L
G80C352236-	36"L
G80C352248-	48"L

22" Deep	
G80C352260-	60"L

**Note:** G80 Acid Storage Cabinets are available with five-knuckle hinges only

# Solvent Storage Fume Hood Base Cabinets

## Specifications

Solvent Storage Cabinets are UL listed for the storage of flammable and combustible liquids, meeting UFC, OSHA and NFPA No. 30-1993 standards.


A 2" deep, steel, liquid-tight pan covers the entire bottom to contain liquid leaks

and spills. The cabinet is provided with a full depth adjustable shelf and two diametrically opposed vents with spark screens in the back of the cabinet for cases when ventilation is required.


Steel Solvent Storage cabinets are all 18

gauge steel, double panel construction with self-closing doors, synchronized so that both doors will always fully close. The right hand door is equipped with a three-point latching system that automatically engages the cabinet frame. Each door is equipped with a fusible-link hold-open feature that ensures the doors close should the temperature outside the cabinet exceed 165°F. Cabinets are provided with a grounding screw at the rear.


All Solvent Storage Cabinets are labeled: CAUTION FLAMMABLE – KEEP FIRE AWAY in English, Spanish, and French.




## Vent Locations




## Fume Hood Solvent Storage Base Cabinets



**20" Deep** Capacity  
**G68C352024L\_100** 30"L 12 gal



**20" Deep** Capacity  
**G68C352024-\_100** 30"L 12 gal



**20" Deep** Capacity  
**G68C352030-\_100** 30"L 12 gal  
**G68C352036-\_100** 36"L 15 gal  
**G68C352048-\_100** 48"L 20 gal

**22" Deep** Capacity  
**G68C352224L\_100** 30"L 12 gal

**22" Deep** Capacity  
**G68C352224-\_100** 30"L 12 gal

**22" Deep** Capacity  
**G68C352230-\_100** 30"L 12 gal  
**G68C352236-\_100** 36"L 15 gal  
**G68C352248-\_100** 48"L 20 gal

G68 Fume Hood Solvent Storage Cabinets are available as mobile.  
To add casters to the cabinets above, add **-B7** to the end of the part number.  
Example: **G68C352224-1100-B7** — Mobile, Overlay style cabinet, 22" deep by 24" long.


# Vacuum Pump Storage Base Cabinets

## Specifications


Vacuum Pump Storage Fume Hood Cabinets are designed without a bottom to allow vacuum pumps and other equipment to be rolled in or out of the cabinets. The interior is lined with an easy to clean sound absorbant material. Each cabinet is furnished with

a 120 VAC, 20 amp, duplex receptacle mounted on the inside cabinet back and a pilot lighted toggle switch mounted in the top front rail. (Wiring is not included.) Cabinets are also furnished with a 1 1/2" diameter PVC vent pipe for venting or access to the

hood above. (Requires a 2" hole in the hood worksurface.) The toespace rail is attached to the door to allow total access to the cabinet. Cabinet inside clearance at the floor is 14 1/2" or 16 1/2" front-to-back, 25" high, and 3" less than the overall cabinet length.




## Fume Hood Vacuum Pump Storage Cabinets




20" Deep

G35C352024L \_\_\_\_\_ 24"L



20" Deep

G35C352024- \_\_\_\_\_ 24"L



20" Deep

G35C352030- \_\_\_\_\_ 30"L  
G35C352036- \_\_\_\_\_ 36"L  
G35C352048- \_\_\_\_\_ 48"L

22" Deep

G35C352224L \_\_\_\_\_ 24"L

22" Deep


G35C352224- \_\_\_\_\_ 24"L

22" Deep

G35C352230- \_\_\_\_\_ 30"L  
G35C352236- \_\_\_\_\_ 36"L  
G35C352248- \_\_\_\_\_ 48"L

## Vacuum Pump Carts

Vacuum Pump Carts are designed for use with Kewaunee's Vacuum Pump cabinets shown on this page. The cart consists of a 1" deep pan, fabricated from powder coated or stainless steel and mounted on 4" swivel casters.



### Carts for 20" Deep Cabinets


Part Number	Length
VPCC061414-00_0	14"
VPCC061420-00_0	20"
VPCC061426-00_0	26"
VPCC061432-00_0	32"
VPCC061438-00_0	38"
VPCC061444-00_0	44"

### Carts for 22" Deep Cabinets

Part Number	Length
VPCC061614-00_0	14"
VPCC061620-00_0	20"
VPCC061626-00_0	26"
VPCC061632-00_0	32"
VPCC061638-00_0	38"
VPCC061644-00_0	44"



# Steel Sitting Height Mobile Cabinets





## Specifications

Steel Sitting Height Mobile Cabinets are 26 1/4" high and 22" deep. Lengths as shown. Designed to fit under 30" high Workstation frames.

Drawer cabinets are furnished with counter weights, multiple-drawer cabinets include an interlock system. All mobile cabinets have finished backs.

**Note:** Last four digits of part number designates door and drawer head style, hinge, hardwood species and pull selection. See page 4 for full explanation. Actual product may vary slightly from catalog illustration due to style and hardware selection.





## Cupboard Cabinets



<b>M10C262215L0000</b>	15"L	<b>M10C262215-0000</b>	15"L
<b>M10C262218L0000</b>	18"L	<b>M10C262218-0000</b>	18"L
<b>M10C262221L0000</b>	21"L	<b>M10C262221-0000</b>	21"L
<b>M10C262224L0000</b>	24"L	<b>M10C262224-0000</b>	24"L

## Drawer/Cupboard Cabinets




## Drawer/Cupboard Cabinets with Pullboards





<b>M20C262215L</b> _____	15"L	<b>M20C262215-</b> _____	15"L	<b>M25C262215L</b> _____	15"L	<b>M25C262215-</b> _____	15"L
<b>M20C262218L</b> _____	18"L	<b>M20C262218-</b> _____	18"L	<b>M25C262218L</b> _____	18"L	<b>M25C262218-</b> _____	18"L
<b>M20C262221L</b> _____	21"L	<b>M20C262221-</b> _____	21"L	<b>M25C262221L</b> _____	21"L	<b>M25C262221-</b> _____	21"L
<b>M20C262224L</b> _____	24"L	<b>M20C262224-</b> _____	24"L	<b>M25C262224L</b> _____	24"L	<b>M25C262224-</b> _____	24"L

## Drawer Cabinets

## Drawer Cabinets with Pullboards




<b>M30C262215-</b> _____	15"L	<b>M35C262215-</b> _____	15"L	<b>M45C262215-</b> _____	15"L
<b>M30C262218-</b> _____	18"L	<b>M35C262218-</b> _____	18"L	<b>M45C262218-</b> _____	18"L
<b>M30C262221-</b> _____	21"L	<b>M35C262221-</b> _____	21"L	<b>M45C262221-</b> _____	21"L
<b>M30C262224-</b> _____	24"L	<b>M35C262224-</b> _____	24"L	<b>M45C262224-</b> _____	24"L



<b>M62C262215-</b> _____	15"L	<b>M66C262215-</b> _____	15"L
<b>M62C262218-</b> _____	18"L	<b>M66C262218-</b> _____	18"L
<b>M62C262221-</b> _____	21"L	<b>M66C262221-</b> _____	21"L
<b>M62C262224-</b> _____	24"L	<b>M66C262224-</b> _____	24"L

# Steel Standing Height Mobile Cabinets





## Specifications

Steel Standing Height Mobile Cabinets are 33" high and 22" deep. Lengths as shown. Designed to fit under 36" high Workstation frames.

Drawer cabinets are furnished with counter weights, multiple-drawer cabinets include an interlock system. All mobile cabinets have finished backs.

**Note:** Last four digits of part number designates door and drawer head style, hinge, hardwood species and pull selection. See page 4 for full explanation. Actual product may vary slightly from catalog illustration due to style and hardware selection.




### Cupboard Cabinets



<b>M10C332215L0000</b>	15"L	<b>M10C332215-0000</b>	15"L
<b>M10C332218L0000</b>	18"L	<b>M10C332218-0000</b>	18"L
<b>M10C332221L0000</b>	21"L	<b>M10C332221-0000</b>	21"L
<b>M10C332224L0000</b>	24"L	<b>M10C332224-0000</b>	24"L

### Drawer/Cupboard Cabinets




### Drawer/Cupboard Cabinets with Pullboards





<b>M20C332215L</b>	15"L	<b>M20C332215-</b>	15"L	<b>M25C332215L</b>	15"L	<b>M25C332215-</b>	15"L
<b>M20C332218L</b>	18"L	<b>M20C332218-</b>	18"L	<b>M25C332218L</b>	18"L	<b>M25C332218-</b>	18"L
<b>M20C332221L</b>	21"L	<b>M20C332221-</b>	21"L	<b>M25C332221L</b>	21"L	<b>M25C332221-</b>	21"L
<b>M20C332224L</b>	24"L	<b>M20C332224-</b>	24"L	<b>M25C332224L</b>	24"L	<b>M25C332224-</b>	24"L

### Drawer Cabinets

### Drawer Cabinets with Pullboards



<b>M40C332215-</b>	15"L	<b>M46C332215-</b>	15"L	<b>M56C332215-</b>	15"L
<b>M40C332218-</b>	18"L	<b>M46C332218-</b>	18"L	<b>M56C332218-</b>	18"L
<b>M40C332221-</b>	21"L	<b>M46C332221-</b>	21"L	<b>M56C332221-</b>	21"L
<b>M40C332224-</b>	24"L	<b>M46C332224-</b>	24"L	<b>M56C332224-</b>	24"L



<b>M71C332215-</b>	15"L	<b>M68C332215-</b>	15"L
<b>M71C332218-</b>	18"L	<b>M68C332218-</b>	18"L
<b>M71C332221-</b>	21"L	<b>M68C332221-</b>	21"L
<b>M71C332224-</b>	24"L	<b>M68C332224-</b>	24"L


# Mobile Cabinet Tops & Accessories

## Specifications:

**1" Thick Mobile Cabinet Tops** are available in three depths, 22", 23", or 24" and three materials: Kemresin, Phenolic Resin or Plastic Laminate. They may be used on any Mobile


Cabinet and add 1" to the overall height of the cabinets. All edges of top are finished.

### 1" Thick — Mobile Cabinet Top (adds 1" to height of Mobile Cabinet)




#### 22" Deep flush to front of case

Length	Kemresin®	Kemresin Colors	Black Phenolic Resin	Black Plastic Laminate
15"	<b>F94R2215-00-__</b>	<b>BK = Black</b>	<b>F94T2215-00</b>	<b>F94P2215-00</b>
18"	<b>F94R2218-00-__</b>	<b>GR = Grey</b>	<b>F94T2218-00</b>	<b>F94P2218-00</b>
21"	<b>F94R2221-00-__</b>	<b>PT = Putty</b>	<b>F94T2221-00</b>	<b>F94P2221-00</b>
24"	<b>F94R2224-00-__</b>	<b>SL = Slate</b>	<b>F94T2224-00</b>	<b>F94P2224-00</b>
30"	<b>F94R2230-00-__</b>		<b>F94T2230-00</b>	<b>F94P2230-00</b>
36"	<b>F94R2236-00-__</b>		<b>F94T2236-00</b>	<b>F94P2236-00</b>
42"	<b>F94R2242-00-__</b>		<b>F94T2242-00</b>	<b>F94P2242-00</b>
48"	<b>F94R2248-00-__</b>		<b>F94T2248-00</b>	<b>F94P2248-00</b>




*Inset Cabinet End View*




#### 23" Deep flush to overlay door/drawer fronts - 1" overhang of case

Length	Kemresin®	Kemresin Colors	Black Phenolic Resin	Black Plastic Laminate
15"	<b>F94R2315-00-__</b>	<b>BK = Black</b>	<b>F94T2315-00</b>	<b>F94P2315-00</b>
18"	<b>F94R2318-00-__</b>	<b>GR = Grey</b>	<b>F94T2318-00</b>	<b>F94P2318-00</b>
21"	<b>F94R2321-00-__</b>	<b>PT = Putty</b>	<b>F94T2321-00</b>	<b>F94P2321-00</b>
24"	<b>F94R2324-00-__</b>	<b>SL = Slate</b>	<b>F94T2324-00</b>	<b>F94P2324-00</b>
30"	<b>F94R2330-00-__</b>		<b>F94T2330-00</b>	<b>F94P2330-00</b>
36"	<b>F94R2336-00-__</b>		<b>F94T2336-00</b>	<b>F94P2336-00</b>
42"	<b>F94R2342-00-__</b>		<b>F94T2342-00</b>	<b>F94P2342-00</b>
48"	<b>F94R2348-00-__</b>		<b>F94T2348-00</b>	<b>F94P2348-00</b>




*Inset Cabinet End View*



#### 24" Deep 1" overhang of overlay door/drawer fronts

Length	Kemresin®	Kemresin Colors	Black Phenolic Resin	Black Plastic Laminate
15"	<b>F94R2415-00-__</b>	<b>BK = Black</b>	<b>F94T2415-00</b>	<b>F94P2415-00</b>
18"	<b>F94R2418-00-__</b>	<b>GR = Grey</b>	<b>F94T2418-00</b>	<b>F94P2418-00</b>
21"	<b>F94R2421-00-__</b>	<b>PT = Putty</b>	<b>F94T2421-00</b>	<b>F94P2421-00</b>
24"	<b>F94R2424-00-__</b>	<b>SL = Slate</b>	<b>F94T2424-00</b>	<b>F94P2424-00</b>
30"	<b>F94R2430-00-__</b>		<b>F94T2430-00</b>	<b>F94P2430-00</b>
36"	<b>F94R2436-00-__</b>		<b>F94T2436-00</b>	<b>F94P2436-00</b>
42"	<b>F94R2442-00-__</b>		<b>F94T2442-00</b>	<b>F94P2442-00</b>
48"	<b>F94R2448-00-__</b>		<b>F94T2448-00</b>	<b>F94P2448-00</b>



*Overlay Cabinet End View*

# Mobile Cabinet Tops & Accessories



**Specifications:**

**Steel Mobile Cabinet Tops** are fabricated from 18 gauge steel and finished to match the cabinet. They are designed for Steel Mobile Cabinets only.

**Mobile Cabinet Grabbar** allows easier movement and relocation of mobile cabinets. Fabricated from powder coated steel or stainless steel, they may be used on any

Mobile Cabinet, with or without Tops. Bars must be attached in field with supplied screws and/or bolts.


**22" Deep — Steel — Mobile Cabinet Top (adds less than 1/4" to height of Mobile Cabinet)**




*Inset Cabinet*  
**End View**

**22" Deep Top**  
*for Steel Mobile Cabinets only*



Lg	Painted Steel	Stainless Steel
15"	<b>F94M2215-00</b>	<b>F94S2215-00</b>
18"	<b>F94M2218-00</b>	<b>F94S2218-00</b>
21"	<b>F94M2221-00</b>	<b>F94S2221-00</b>
24"	<b>F94M2224-00</b>	<b>F94S2224-00</b>
30"	<b>F94M2230-00</b>	<b>F94S2230-00</b>
36"	<b>F94M2236-00</b>	<b>F94S2236-00</b>
42"	<b>F94M2242-00</b>	<b>F94S2242-00</b>
48"	<b>F94M2248-00</b>	<b>F94S2248-00</b>



*Overlay Cabinet*  
**End View**




**23" Deep — Steel — Mobile Cabinet Top (adds 1/2" to height of Mobile Cabinet)**




*Inset Cabinet*  
**End View**

**23" Deep Top**  
*for Steel Mobile Cabinets only*


Lg	Painted Steel	Stainless Steel
15"	<b>F94M2315-00</b>	<b>F94S2315-00</b>
18"	<b>F94M2318-00</b>	<b>F94S2318-00</b>
21"	<b>F94M2321-00</b>	<b>F94S2321-00</b>
24"	<b>F94M2324-00</b>	<b>F94S2324-00</b>
30"	<b>F94M2330-00</b>	<b>F94S2330-00</b>
36"	<b>F94M2336-00</b>	<b>F94S2336-00</b>
42"	<b>F94M2342-00</b>	<b>F94S2342-00</b>
48"	<b>F94M2348-00</b>	<b>F94S2348-00</b>




*Overlay Cabinet*  
**End View**



**Mobile Cabinet Grabbar**




Length	Painted Steel	Stainless Steel
15"	<b>F73M0015-00</b>	<b>F73S0015-00</b>
18"	<b>F73M0018-00</b>	<b>F73S0018-00</b>
21"	<b>F73M0021-00</b>	<b>F73S0021-00</b>
24"	<b>F73M0024-00</b>	<b>F73S0024-00</b>



# Steel Apron Assemblies

## Components of an Apron Assembly





## Specifications:

Box Aprons are 4" high in lengths and depths as shown. Dimensions shown are sizes of the apron without legs. Each leg or pedestal adds 2" to the overall length of the apron assembly. A free-standing table with legs on both ends will therefore be 4" longer than the apron length.

All aprons can be used between cabinets, with leg assemblies or with pedestal legs.

For aprons 42" long and longer, spreaders are available to connect the two leg assemblies. When a spreader is used, order the A08 or A28 leg assembly which are designed to accept the spreader.

## 22" Deep Plain Rail Box Apron




All Plain Rail Box Aprons available in either Inset or Full Overlay Style

22" Deep Apron		22" Deep Apron		22" Deep Apron	
Part Number	Apron Length	Part Number	Apron Length	Part Number	Apron Length
A00C042215-0	15"	A00C042236-0	36"	A00C042257-0	57"
A00C042216-0	16"	A00C042237-0	37"	A00C042258-0	58"
A00C042217-0	17"	A00C042238-0	38"	A00C042259-0	59"
A00C042218-0	18"	A00C042239-0	39"	A00C042260-0	60"
A00C042219-0	19"	A00C042240-0	40"	A00C042261-0	61"
A00C042220-0	20"	A00C042241-0	41"	A00C042262-0	62"
A00C042221-0	21"	A00C042242-0	42"	A00C042263-0	63"
A00C042222-0	22"	A00C042243-0	43"	A00C042264-0	64"
A00C042223-0	23"	A00C042244-0	44"	A00C042265-0	65"
A00C042224-0	24"	A00C042245-0	45"	A00C042266-0	66"
A00C042225-0	25"	A00C042246-0	46"	A00C042267-0	67"
A00C042226-0	26"	A00C042247-0	47"	A00C042268-0	68"
A00C042227-0	27"	A00C042248-0	48"	A00C042219-0	69"
A00C042228-0	28"	A00C042249-0	49"	A00C042270-0	70"
A00C042229-0	29"	A00C042250-0	50"	A00C042271-0	71"
A00C042230-0	30"	A00C042251-0	51"	A00C042272-0	72"
A00C042231-0	31"	A00C042252-0	52"		
A00C042232-0	32"	A00C042253-0	53"		
A00C042233-0	33"	A00C042254-0	54"		
A00C042234-0	34"	A00C042255-0	55"		
A00C042235-0	35"	A00C042256-0	56"		

# Steel Apron Assemblies

## Components of an Apron Assembly




## Specifications:

Box Aprons are 4" high in lengths and depths as shown. Dimensions shown are sizes of the apron without legs. Each leg or pedestal adds 2" to the overall length of the apron assembly. A free-standing table with legs on both ends will therefore be 4" longer than the apron length.


All aprons can be used between cabinets, with leg assemblies or with pedestal legs.

For aprons 42" long and longer, spreaders are available to connect the two leg assemblies. When a spreader is used, order the A08 or A28 leg assembly which are designed to accept the spreader.

## 28" Deep Plain Rail Box Apron



Inset Style Apron




Full Overlay Style Apron

All Plain Rail Box Aprons available in either Inset or Full Overlay Style

28" Deep Apron		28" Deep Apron		28" Deep Apron	
	Apron Length		Apron Length		Apron Length
A00C042815-___0	15"	A00C042836-___0	36"	A00C042857-___0	57"
A00C042816-___0	16"	A00C042837-___0	37"	A00C042858-___0	58"
A00C042817-___0	17"	A00C042838-___0	38"	A00C042859-___0	59"
A00C042818-___0	18"	A00C042839-___0	39"	A00C042860-___0	60"
A00C042819-___0	19"	A00C042840-___0	40"	A00C042861-___0	61"
A00C042820-___0	20"	A00C042841-___0	41"	A00C042862-___0	62"
A00C042821-___0	21"	A00C042842-___0	42"	A00C042863-___0	63"
A00C042822-___0	22"	A00C042843-___0	43"	A00C042864-___0	64"
A00C042823-___0	23"	A00C042844-___0	44"	A00C042865-___0	65"
A00C042824-___0	24"	A00C042845-___0	45"	A00C042866-___0	66"
A00C042825-___0	25"	A00C042846-___0	46"	A00C042867-___0	67"
A00C042826-___0	26"	A00C042847-___0	47"	A00C042868-___0	68"
A00C042827-___0	27"	A00C042848-___0	48"	A00C042819-___0	69"
A00C042828-___0	28"	A00C042849-___0	49"	A00C042870-___0	70"
A00C042829-___0	29"	A00C042850-___0	50"	A00C042871-___0	71"
A00C042830-___0	30"	A00C042851-___0	51"	A00C042872-___0	72"
A00C042831-___0	31"	A00C042852-___0	52"		
A00C042832-___0	32"	A00C042853-___0	53"		
A00C042833-___0	33"	A00C042854-___0	54"		
A00C042834-___0	34"	A00C042855-___0	55"		
A00C042835-___0	35"	A00C042856-___0	56"		

# Steel Apron Assemblies

## Components of an Apron Assembly




## Specifications:

Box Aprons are 4" high in lengths and depths as shown. Dimensions shown are sizes of the apron without legs. Each leg or pedestal adds 2" to the overall length of the apron assembly. A free-standing table with legs on both ends will therefore be 4" longer than the apron length.


All aprons can be used between cabinets, with leg assemblies or with pedestal legs.

For aprons 42" long and longer, spreaders are available to connect the two leg assemblies. When a spreader is used, order the A08 or A28 leg assembly which are designed to accept the spreader.

## 22" Deep Box Apron with Drawers



1 Drawer




22" Deep Apron

2 Drawers



22" Deep Apron

3 Drawers



22" Deep Apron


	Apron Length	Dwr Size
A51C042220-1	20"	18"
A51C042224-1	24"	21"
A51C042226-1	26"	24"
A51C042230-1	30"	27"
A51C042232-1	32"	30"
A51C042236-1	36"	33"

	Apron Length	Dwr Size
A52C042226-1	26"	12"
A52C042230-1	30"	12"
A52C042232-1	32"	15"
A52C042236-1	36"	15"
A52C042238-1	38"	18"
A52C042242-1	42"	18"
A52C042244-1	44"	21"
A52C042248-1	48"	21"
A52C042250-1	50"	24"
A52C042254-1	54"	24"
A52C042256-1	56"	27"
A52C042260-1	60"	27"
A52C042262-1	62"	30"
A52C042266-1	66"	30"
A52C042268-1	68"	33"
A52C042272-1	72"	33"

	Apron Length	Dwr Size
A53C042250-1	50"	15"
A53C042254-1	54"	15"
A53C042256-1	56"	17"
A53C042260-1	60"	17"
A53C042262-1	62"	19"
A53C042266-1	66"	19"
A53C042268-1	68"	21"
A53C042272-1	72"	21"

# Steel Apron Assemblies

## Components of an Apron Assembly




## Specifications:

Box Aprons are 4" high in lengths and depths as shown. Dimensions shown are sizes of the apron without legs. Each leg or pedestal adds 2" to the overall length of the apron assembly. A free-standing table with legs on both ends will therefore be 4" longer than the apron length.


All aprons can be used between cabinets, with leg assemblies or with pedestal legs.

For aprons 42" long and longer, spreaders are available to connect the two leg assemblies. When a spreader is used, order the A08 or A28 leg assembly which are designed to accept the spreader.

## 28" Deep Box Apron with Drawers



1 Drawer




28" Deep Apron

2 Drawers



28" Deep Apron

3 Drawers



28" Deep Apron

	Apron Length	Dwr Size
A51C042820-1	20"	18"
A51C042824-1	24"	21"
A51C042826-1	26"	24"
A51C042830-1	30"	27"
A51C042832-1	32"	30"
A51C042836-1	36"	33"

	Apron Length	Dwr Size
A52C042826-1	26"	12"
A52C042830-1	30"	12"
A52C042832-1	32"	15"
A52C042836-1	36"	15"
A52C042838-1	38"	18"
A52C042842-1	42"	18"
A52C042844-1	44"	21"
A52C042848-1	48"	21"
A52C042850-1	50"	24"
A52C042854-1	54"	24"
A52C042856-1	56"	27"
A52C042860-1	60"	27"
A52C042862-1	62"	30"
A52C042866-1	66"	30"
A52C042868-1	68"	33"
A52C042872-1	72"	33"

	Apron Length	Dwr Size
A53C042850-1	50"	15"
A53C042854-1	54"	15"
A53C042856-1	56"	17"
A53C042860-1	60"	17"
A53C042862-1	62"	19"
A53C042866-1	66"	19"
A53C042868-1	68"	21"
A53C042872-1	72"	21"



# Steel Apron Assemblies


## Specifications:

Corner Box Aprons provide support for corner work tops as well as the extra storage space of an apron drawer. The blind end may either be attached to the wall or supported by a leg or pedestal assembly.


Corner Box Aprons are furnished with steel panels when specified with wood drawer heads.

2" High Box Aprons are designed for use when knee space height is critical. They are available 22" deep and in lengths from 15" to 72" on 1" increments. 2" High Box Aprons are not designed for use with legs.


## Corner Box Apron



*Furnished with steel panels when specified with wood drawer heads.*




**1 Drawer on Right**



**22" Deep**

**1 Drawer on Left**



**22" Deep**

	Apron Length	Dwr Size
<b>A05C042248-1</b> _____	48"	21"
<b>A05C042254-1</b> _____	54"	21"
<b>A05C042260-1</b> _____	60"	21"
<b>A05C042266-1</b> _____	66"	21"
<b>A05C042272-1</b> _____	72"	21"

	Apron Length	Dwr Size
<b>A04C042248-1</b> _____	48"	21"
<b>A04C042254-1</b> _____	54"	21"
<b>A04C042260-1</b> _____	60"	21"
<b>A04C042266-1</b> _____	66"	21"
<b>A04C042272-1</b> _____	72"	21"

**28" Deep**

	Apron Length	Dwr Size
<b>A05C042848-1</b> _____	48"	21"
<b>A05C042854-1</b> _____	54"	21"
<b>A05C042860-1</b> _____	60"	21"
<b>A05C042866-1</b> _____	66"	21"
<b>A05C042872-1</b> _____	72"	21"

**28" Deep**

	Apron Length	Dwr Size
<b>A04C042848-1</b> _____	48"	21"
<b>A04C042854-1</b> _____	54"	21"
<b>A04C042860-1</b> _____	60"	21"
<b>A04C042866-1</b> _____	66"	21"
<b>A04C042872-1</b> _____	72"	21"

# Steel Apron Assemblies

## Specifications:


Apron Rails may be used in place of Box Aprons to provide worksurface support between two units. They are available 2" or 4" high and from 15" to 72" long in 1" increments. They may be installed vertically to provide a 2"

or 4" high rail, or turned 90 degrees to provide a 1 1/2" high rail. Apron rails should be ordered in pairs to provide both front and rear top support.

Work Top Support Rails and Assemblies may be used anywhere

maximum knee space is wanted but top support is required. Fabricated of heavy gauge steel channels, it is only 7/8" high. The Keyboard Tray may be added to any of the frames or may be used independently.

### Inset 4" High Apron Rail




**A06C0402**□□-00\_0 4" High  
Inset Apron Rails are available from 15" to 72" long in 1" increments. To specify apron rail, insert the length required into the part number.

**Example:** A06C040248-0000 is a 48" long cold rolled steel apron rail.

Apron rails should be ordered in pairs to provide both front and rear worksurface support.

### Overlay 4" High Apron Rail




**A06C0403**□□-1\_\_0 4" High  
Overlay Apron Rails are available from 15" to 72" long in 1" increments. To specify apron rail, insert the length required into the part number.

**Example:** A06C040248-1610 is a 48" long apron rail with Oak overlay panel.

Apron rails should be ordered in pairs to provide both front and rear worksurface support.

### 4" High Adjustable Rail




**A56C040236-00\_0** 36"  
**A56C040248-00\_0** 48"  
**A56C040260-00\_0** 60"  
**A56C040272-00\_0** 72"

Same as A06 Apron Rails except furnished with a loose angle at one end to allow the rail to be cut in the field to fit available space.

Apron rails should be ordered in pairs to provide both front and rear worksurface support.

### 2" High Plain Rail Box Apron




**A00C0222**□□-00\_0  
2" high Box Aprons are available from 15" to 72" long in 1" increments. To specify apron, insert the length required into the part number.

**Example:** A00C020248-0000 is a 48" long 2" high cold rolled steel box apron.

**Not for use with legs**

### 2" High Apron Rail




**A06C0202**□□-00\_0 2" High  
Apron rails are available from 15" to 72" long in 1" increments. To specify apron rail, insert the length required into the part number.

**Example:** A06C020248-0000 is a 48" long cold rolled steel apron rail.

Apron rails should be ordered in pairs to provide both front and rear worksurface support.

### Support Rails



**A81C0102**□□-0000  
Support rails may be used individually or assembled to form a support framework. Rails are available from 12" to 72" long in 1" increments and are furnished with angles attached to each end. To specify rail, insert the length required in the part number.

**Example:** A81010248-0000 is a 48" long support rail.

Support rails should be ordered in pairs to provide both front and rear worksurface support.

# Steel Leg Assemblies & Pedestals

## Specifications:

Legs are 2" square and all except the Top Mounted Leg are joined in pairs with a 2" high stretcher and a 4" high apron rail. Dimensions are as shown.

For aprons 42" long and longer, spreaders are available to connect the two leg assemblies. When a

spreader is used, order the A08 or A28 leg assembly which are designed for use with a spreader.





All legs may be ordered with 4" swivel locking casters by adding **-B10** to the end of the part numbers.

Example: A07C292202-0000-**B10**

Adjustable Height Legs are adjustable in 1" increments from 29" to 35" high.

### End Leg Assemblies

### Intermediate Leg Assemblies



For use with Spreader

For use with Spreader

**22<sup>3</sup>/<sub>16</sub>" Deep**  
**A07C292202-00\_0** 29"H  
**A07C322202-00\_0** 32"H  
**A07M352202-00\_0** 35"H

**22<sup>3</sup>/<sub>16</sub>" Deep**  
**A08C292202-00\_0** 29"H  
**A08C322202-00\_0** 32"H  
**A08C352202-00\_0** 35"H

**22<sup>3</sup>/<sub>16</sub>" Deep**  
**A27C292202-00\_0** 29"H  
**A27C322202-00\_0** 32"H  
**A27C352202-00\_0** 35"H

**22<sup>3</sup>/<sub>16</sub>" Deep**  
**A28C292202-00\_0** 29"H  
**A28C322202-00\_0** 32"H  
**A28C352202-00\_0** 35"H

**28<sup>3</sup>/<sub>16</sub>" Deep**  
**A07C292802-00\_0** 29"H  
**A07C322802-00\_0** 32"H  
**A07C352802-00\_0** 35"H

**28<sup>3</sup>/<sub>16</sub>" Deep**  
**A08C292802-00\_0** 29"H  
**A08C322802-00\_0** 32"H  
**A08C352802-00\_0** 35"H





**28<sup>3</sup>/<sub>16</sub>" Deep**  
**A27C292802-00\_0** 29"H  
**A27C322802-00\_0** 32"H  
**A27C352802-00\_0** 35"H

**28<sup>3</sup>/<sub>16</sub>" Deep**  
**A28C292802-00\_0** 29"H  
**A28C322802-00\_0** 32"H  
**A28C352802-00\_0** 35"H

Each Leg Assembly adds 2" to the overall length of the apron assembly to which it is attached..

### Adjustable Height End Leg Assemblies

### Adjustable Height Intermediate Leg Ass'y



For use with Spreader

For use with Spreader

**22<sup>3</sup>/<sub>16</sub>" Deep**  
**A57C352202-00\_0** Adj.  
 Adjustable Height 29"-35"

**22<sup>3</sup>/<sub>16</sub>" Deep**  
**A58C352202-00\_0** Adj.  
 Adjustable Height 29"-35"

**22<sup>3</sup>/<sub>16</sub>" Deep**  
**A77C352202-00\_0** Adj.  
 Adjustable Height 29"-35"

**22<sup>3</sup>/<sub>16</sub>" Deep**  
**A78C352202-00\_0** Adj.  
 Adjustable Height 29"-35"

**28<sup>3</sup>/<sub>16</sub>" Deep**  
**A57C352802-00\_0** Adj.  
 Adjustable Height 29"-35"

**28<sup>3</sup>/<sub>16</sub>" Deep**  
**A58C352802-00\_0** Adj.  
 Adjustable Height 29"-35"

**28<sup>3</sup>/<sub>16</sub>" Deep**  
**A77C352802-00\_0** Adj.  
 Adjustable Height 29"-35"

**28<sup>3</sup>/<sub>16</sub>" Deep**  
**A78C352802-00\_0** Adj.  
 Adjustable Height 29"-35"

# Steel Apron Rails & Panels

## Specifications:


Knee Space panels are designed to close off the space between the floor and the bottom of the apron or countertop. They are sized for 29" high, 32" high, and 35" high cabinets. Each Knee Space Panel is made with a separate base rail, to allow the removal

of the panel without removing the base molding.

Knee Space Panels are furnished with two angles to attach to panel to the cabinets.

Pedestal Legs are 2" wide, depths and heights as shown. They are designed to be used with or without Box Aprons to form knee spaces at the end of assemblies.

## Knee Space Panels



Knee space panels are available from 15" to 72" long in 1" increments. To specify knee space panel, insert the length required into the part number.

**Example:** A10C310048-0000 is a 48" long knee space panel.

### Sitting Height Knee Space Panels

Apron Part Number	Height	Length	for use with:
A10C2500 □□-00_0	25"	15"-72"L	4" high apron
A10C2800 □□-00_0	28"	15"-72"L	A80 support frame
A10C2900 □□-00_0	29"	15"-72"L	no apron


### ADA Height Knee Space Panels

Apron Part Number	Height	Length	for use with:
A10C2800 □□-00_0	28"	15"-72"L	4" high apron
A10C3100 □□-00_0	31"	15"-72"L	A80 support frame
A10C3200 □□-00_0	32"	15"-72"L	no apron or frame

### Standing Height Knee Space Panels


Apron Part Number	Height	Length	for use with:
A10C3100 □□-00_0	31"	15"-72"L	4" high apron
A10C3400 □□-00_0	34"	15"-72"L	A80 support frame
A10C3500 □□-00_0	35"	15"-72"L	no apron or frame

## Spreader



A09C020230-00_0	30"L
A09C020236-00_0	36"L
A09C020242-00_0	42"L
A09C020248-00_0	48"L
A09C020254-00_0	54"L
A09C020260-00_0	60"L
A09C020266-00_0	66"L
A09C020272-00_0	72"L


## Top Mounted Leg



A15C290202-0000	29"H
A15C320202-00_0	32"H
A15C350202-00_0	35"H
A65C350202-00_0	Adj. Adjustable Height 29"-35"

2" square leg welded to 4" square steel plate

## Wall Pedestal Leg



### Left Side 22" Deep

A11C292202-00_0	29"H
A11C322202-00_0	32"H
A11C352202-00_0	35"H

### Left Side 28" Deep

A11C292802-00_0	29"H
A11C322802-00_0	32"H
A11C352802-00_0	35"H


### Right Side 22" Deep

A12C292202-00_0	29"H
A12C322202-00_0	32"H
A12C352202-00_0	35"H

### Right Side 28" Deep

A12C292802-00_0	29"H
A12C322802-00_0	32"H
A12C352802-00_0	35"H

## Center Pedestal Leg



### 52" Deep


A13C295202-00_0	29"H
A13C325202-00_0	32"H
A13C355202-00_0	35"H

Each Pedestal Leg adds 2" to the overall length of the assembly to which it is attached.

Pedestal Legs may be ordered with the topspace 7/8" deeper so that the front edge lines up with the doors and drawer fronts on adjacent Full Overlay cabinets. To specify the deeper toe space, add a **C** to the end of the part number.

Example: A11C352202-0000-C

# Steel Wall Mounted Storage Cabinets











## Specifications:


Open Wall Mounted Storage Cabinets are available in four heights, 24", 30", 36" and 48" and two depths; 12<sup>3</sup>/<sub>4</sub>" and 16". Lengths are as shown. All are furnished with one, two, or three adjustable shelves as shown.

Wall cabinets are designed to be mounted on hangers rails. Shown on page 47 and 81.

## Open Storage

<b>24" High</b>  	<b>12<sup>3</sup>/<sub>4</sub>" Deep</b> <hr/> W00C241315-00_0 15"L W00C241318-00_0 18"L W00C241321-00_0 21"L W00C241324-00_0 24"L	<b>24" High</b>  	<b>12<sup>3</sup>/<sub>4</sub>" Deep</b> <hr/> W00C241330-00_0 30"L W00C241336-00_0 36"L W00C241342-00_0 42"L W00C241348-00_0 48"L
	<b>16" Deep</b> <hr/> W00C241615-00_0 15"L W00C241618-00_0 18"L W00C241621-00_0 21"L W00C241624-00_0 24"L		<b>16" Deep</b> <hr/> W00C241630-00_0 30"L W00C241636-00_0 36"L W00C241642-00_0 42"L W00C241648-00_0 48"L
<b>30" High</b>  	<b>12<sup>3</sup>/<sub>4</sub>" Deep</b> <hr/> W00C301315-00_0 15"L W00C301318-00_0 18"L W00C301321-00_0 21"L W00C301324-00_0 24"L	<b>30" High</b>  	<b>12<sup>3</sup>/<sub>4</sub>" Deep</b> <hr/> W00C301330-00_0 30"L W00C301336-00_0 36"L W00C301342-00_0 42"L W00C301348-00_0 48"L
	<b>16" Deep</b> <hr/> W00C301615-00_0 15"L W00C301618-00_0 18"L W00C301621-00_0 21"L W00C301624-00_0 24"L		<b>16" Deep</b> <hr/> W00C301630-00_0 30"L W00C301636-00_0 36"L W00C301642-00_0 42"L W00C301648-00_0 48"L
<b>36" High</b>  	<b>12<sup>3</sup>/<sub>4</sub>" Deep</b> <hr/> W00C361315-00_0 15"L W00C361318-00_0 18"L W00C361321-00_0 21"L W00C361324-00_0 24"L	<b>36" High</b>  	<b>12<sup>3</sup>/<sub>4</sub>" Deep</b> <hr/> W00C361330-00_0 30"L W00C361336-00_0 36"L W00C361342-00_0 42"L W00C361348-00_0 48"L
	<b>16" Deep</b> <hr/> W00C361615-00_0 15"L W00C361618-00_0 18"L W00C361621-00_0 21"L W00C361624-00_0 24"L		<b>16" Deep</b> <hr/> W00C361630-00_0 30"L W00C361636-00_0 36"L W00C361642-00_0 42"L W00C361648-00_0 48"L
<b>48" High</b>  	<b>12<sup>3</sup>/<sub>4</sub>" Deep</b> <hr/> W00C481315-00_0 15"L W00C481318-00_0 18"L W00C481321-00_0 21"L W00C481324-00_0 24"L	<b>48" High</b>  	<b>12<sup>3</sup>/<sub>4</sub>" Deep</b> <hr/> W00C481330-00_0 30"L W00C481336-00_0 36"L W00C481342-00_0 42"L W00C481348-00_0 48"L
	<b>16" Deep</b> <hr/> W00C481615-00_0 15"L W00C481618-00_0 18"L W00C481621-00_0 21"L W00C481624-00_0 24"L		<b>16" Deep</b> <hr/> W00C481630-00_0 30"L W00C481636-00_0 36"L W00C481642-00_0 42"L W00C481648-00_0 48"L

# Steel Wall Mounted Storage Cabinets



### Specifications:

Plate Glass Door Wall Mounted Storage Cabinets are available in three heights, 24", 30", 36" and 48" and two depths; 12<sup>3</sup>/<sub>4</sub>" and 16". Lengths are as shown. Cabinets are furnished with one, two, or three adjustable shelves as shown.


Doors are 1/4" float glass and are mounted on friction-free roller glides.

Wall cabinets are designed to be mounted on hangers rails. Shown on page 47 and 81.

## Sliding Plate Glass Doors — Clear Glass

## Wall Mounted Cabinet Hanging Rails

### 24" High




#### 12<sup>3</sup>/<sub>4</sub>" Deep

W10C241330-00_0	30"L
W10C241336-00_0	36"L
W10C241342-00_0	42"L
W10C241348-00_0	48"L

#### 16" Deep

W10C241630-00_0	30"L
W10C241636-00_0	36"L
W10C241642-00_0	42"L
W10C241648-00_0	48"L

### 30" High




#### 12<sup>3</sup>/<sub>4</sub>" Deep

W10C301330-00_0	30"L
W10C301336-00_0	36"L
W10C301342-00_0	42"L
W10C301348-00_0	48"L

#### 16" Deep

W10C301630-00_0	30"L
W10C301636-00_0	36"L
W10C301642-00_0	42"L
W10C301648-00_0	48"L

### 36" High




#### 12<sup>3</sup>/<sub>4</sub>" Deep

W10C361330-00_0	30"L
W10C361336-00_0	36"L
W10C361342-00_0	42"L
W10C361348-00_0	48"L

#### 16" Deep

W10C361630-00_0	30"L
W10C361636-00_0	36"L
W10C361642-00_0	42"L
W10C361648-00_0	48"L

### 48" High




#### 12<sup>3</sup>/<sub>4</sub>" Deep

W10C481330-00_0	30"L
W10C481336-00_0	36"L
W10C481342-00_0	42"L
W10C481348-00_0	48"L

#### 16" Deep


W10C481630-00_0	30"L
W10C481636-00_0	36"L
W10C481642-00_0	42"L
W10C481648-00_0	48"L



- L01-020109-00** for 15" wall cabinet
- L01-020112-00** for 18" wall cabinet
- L01-020115-00** for 21" wall cabinet
- L01-020118-00** for 24" wall cabinet
- L01-020124-00** for 30" wall cabinet
- L01-020130-00** for 36" wall cabinet
- L01-020136-00** for 42" wall cabinet
- L01-020142-00** for 48" wall cabinet

Formed of heavy gauge cold-rolled steel. Designed to capture lip on top rear of wall cabinets for easy installation.

# Steel Wall Mounted Storage Cabinets



## Specifications:

Sliding Door Wall Mounted Storage Cabinets are available in four heights, 24", 30", 36" and 48" and two depths; 12<sup>3</sup>/<sub>4</sub>" and 16". Lengths are as shown. Cabinets are glazed with 1/8" float glass unless otherwise specified and are furnished with one, two or three adjustable shelves as shown. Doors are

mounted on friction-free roller glides.


Sliding Door Wall Mounted Storage Cabinets are available in Door Style 1 and 6 only.

Wall cabinets are designed to be mounted on hangers rails. Shown on page 47 and 81.

### Sliding Glazed Doors

### Sliding Panel Doors

#### 24" High




#### 12<sup>3</sup>/<sub>4</sub>" Deep

W41C241330-0__0	30"
W41C241336-0__0	36"
W41C241342-0__0	42"
W41C241348-0__0	48"

#### 16" Deep

W41C241630-0__0	30"
W41C241636-0__0	36"
W41C241642-0__0	42"
W41C241648-0__0	48"

#### 24" High




#### 12<sup>3</sup>/<sub>4</sub>" Deep

W12C241330-0__0	30"
W12C241336-0__0	36"
W12C241342-0__0	42"
W12C241348-0__0	48"

#### 16" Deep

W12C241630-0__0	30"
W12C241636-0__0	36"
W12C241642-0__0	42"
W12C241648-0__0	48"

#### 30" High




#### 12<sup>3</sup>/<sub>4</sub>" Deep

W41C301330-0__0	30"
W41C301336-0__0	36"
W41C301342-0__0	42"
W41C301348-0__0	48"

#### 16" Deep

W41C301630-0__0	30"
W41C301636-0__0	36"
W41C301642-0__0	42"
W41C301648-0__0	48"

#### 30" High




#### 12<sup>3</sup>/<sub>4</sub>" Deep

W12C301330-0__0	30"
W12C301336-0__0	36"
W12C301342-0__0	42"
W12C301348-0__0	48"

#### 16" Deep

W12C301630-0__0	30"
W12C301636-0__0	36"
W12C301642-0__0	42"
W12C301648-0__0	48"

#### 36" High




#### 12<sup>3</sup>/<sub>4</sub>" Deep

W41C361330-0__0	30"
W41C361336-0__0	36"
W41C361342-0__0	42"
W41C361348-0__0	48"

#### 16" Deep

W41C361630-0__0	30"
W41C361636-0__0	36"
W41C361642-0__0	42"
W41C361648-0__0	48"

#### 36" High




#### 12<sup>3</sup>/<sub>4</sub>" Deep

W12C361330-0__0	30"
W12C361336-0__0	36"
W12C361342-0__0	42"
W12C361348-0__0	48"

#### 16" Deep

W12C361630-0__0	30"
W12C361636-0__0	36"
W12C361642-0__0	42"
W12C361648-0__0	48"

#### 48" High




#### 12<sup>3</sup>/<sub>4</sub>" Deep

W41C481330-0__0	30"
W41C481336-0__0	36"
W41C481342-0__0	42"
W41C481348-0__0	48"

#### 16" Deep

W41C481630-0__0	30"
W41C481636-0__0	36"
W41C481642-0__0	42"
W41C481648-0__0	48"

#### 48" High




#### 12<sup>3</sup>/<sub>4</sub>" Deep

W12C481330-0__0	30"
W12C481336-0__0	36"
W12C481342-0__0	42"
W12C481348-0__0	48"

#### 16" Deep

W12C481630-0__0	30"
W12C481636-0__0	36"
W12C481642-0__0	42"
W12C481648-0__0	48"

# Steel Wall Mounted Storage Cabinets



### Specifications:


Double Sided Plate Glass Pass-Thru Cabinets are available 24", 30", and 36" high, 12<sup>3</sup>/<sub>4</sub>" and 16" deep and in lengths shown. Double Sided Pass-Thru Cabinets have doors on both sides and are designed to be mounted on center bench reagent racks.

Doors are 1/4" float glass and are mounted on friction-free roller glides.

Wall cabinets are designed to be mounted on hangers rails. Shown on page 47 and 81.

## Double Sided Pass-Thru Cabinet - Clear Glass

### 24" High




#### 12<sup>3</sup>/<sub>4</sub>" Deep

W13C241330-00_0	30"L
W13C241336-00_0	36"L
W13C241342-00_0	42"L
W13C241348-00_0	48"L

#### 16" Deep

W13C241630-00_0	30"L
W13C241636-00_0	36"L
W13C241642-00_0	42"L
W13C241648-00_0	48"L

### 30" High




#### 12<sup>3</sup>/<sub>4</sub>" Deep

W13C301330-00_0	30"L
W13C301336-00_0	36"L
W13C301342-00_0	42"L
W13C301348-00_0	48"L

#### 16" Deep

W13C301630-00_0	30"L
W13C301636-00_0	36"L
W13C301642-00_0	42"L
W13C301648-00_0	48"L

### 36" High



#### 12<sup>3</sup>/<sub>4</sub>" Deep


W13C361330-00_0	30"L
W13C361336-00_0	36"L
W13C361342-00_0	42"L
W13C361348-00_0	48"L

#### 16" Deep

W13C361630-00_0	30"L
W13C361636-00_0	36"L
W13C361642-00_0	42"L
W13C361648-00_0	48"L



# Steel Wall Mounted Storage Cabinets




## Specifications:

Swinging Glazed Door Wall Mounted Storage Cabinets are available in four heights, 24", 30", 36" and 48" and two depths; 12<sup>3</sup>/<sub>4</sub>" and 16". Lengths are as shown. Cabinets are glazed with 1/8" float glass unless otherwise specified and are furnished with one, two or three adjustable shelves as shown.

**Note:** Last four digits of part number designate door style, material, and pull selection. See page 4 for full explanation. Beveled Edge, Door Style 7 is not available for wall cabinets. Wall cabinets are designed to be mounted on hanging rails. Shown on page 47 and 81.


## Swinging Glazed Doors

### 24" High




Left Hand	Right Hand	12 <sup>3</sup> / <sub>4</sub> " D
W50C241315L	W50C241315-	15"L
W50C241318L	W50C241318-	18"L
W50C241321L	W50C241321-	21"L
W50C241324L	W50C241324-	24"L
		<b>16" D</b>
W50C241615L	W50C241615-	15"L
W50C241618L	W50C241618-	18"L
W50C241621L	W50C241621-	21"L
W50C241624L	W50C241624-	24"L

### 24" High




12 <sup>3</sup> / <sub>4</sub> " Deep	16" Deep
W50C241330-	W50C241630-
W50C241336-	W50C241636-
W50C241342-	W50C241642-
W50C241348-	W50C241648-

### 30" High




Left Hand	Right Hand	12 <sup>3</sup> / <sub>4</sub> " D
W50C301315L	W50C301315-	15"L
W50C301318L	W50C301318-	18"L
W50C301321L	W50C301321-	21"L
W50C301324L	W50C301324-	24"L
		<b>16" D</b>
W50C301615L	W50C301615-	15"L
W50C301618L	W50C301618-	18"L
W50C301621L	W50C301621-	21"L
W50C301624L	W50C301624-	24"L

### 30" High




12 <sup>3</sup> / <sub>4</sub> " Deep	16" Deep
W50C301330-	W50C301630-
W50C301336-	W50C301636-
W50C301342-	W50C301642-
W50C301348-	W50C301648-

### 36" High




Left Hand	Right Hand	12 <sup>3</sup> / <sub>4</sub> " D
W50C361315L	W50C361315-	15"L
W50C361318L	W50C361318-	18"L
W50C361321L	W50C361321-	21"L
W50C361324L	W50C361324-	24"L
		<b>16" D</b>
W50C361615L	W50C361615-	15"L
W50C361618L	W50C361618-	18"L
W50C361621L	W50C361621-	21"L
W50C361624L	W50C361624-	24"L

### 36" High




12 <sup>3</sup> / <sub>4</sub> " Deep	16" Deep
W50C361330-	W50C361630-
W50C361336-	W50C361636-
W50C361342-	W50C361642-
W50C361348-	W50C361648-

### 48" High




Left Hand	Right Hand	12 <sup>3</sup> / <sub>4</sub> " D
W50C481315L	W50C481315-	15"L
W50C481318L	W50C481318-	18"L
W50C481321L	W50C481321-	21"L
W50C481324L	W50C481324-	24"L
		<b>16" D</b>
W50C481615L	W50C481615-	15"L
W50C481618L	W50C481618-	18"L
W50C481621L	W50C481621-	21"L
W50C481624L	W50C481624-	24"L

### 48" High



12 <sup>3</sup> / <sub>4</sub> " Deep	16" Deep
W50C481330-	W50C481630-
W50C481336-	W50C481636-
W50C481342-	W50C481642-
W50C481348-	W50C481648-

# Steel Wall Mounted Storage Cabinets




### Specifications:

Swinging Panel Door Wall Mounted Storage Cabinets are available in four heights, 24", 30", 36" and 48" and two depths; 12<sup>3</sup>/<sub>4</sub>" and 16". Lengths are as shown. Cabinets are furnished with one, two or three adjustable shelves as shown.

**Note:** Last four digits of part number designate door style, material, and pull selection. See page 4 for full explanation. Beveled Edge, Door Style 7 is not available for wall cabinets. Wall cabinets are designed to be mounted on hanging rails. Shown on page 47 and 81.


### Swinging Panel Doors

#### 24" High




Left Hand	Right Hand	12 <sup>3</sup> / <sub>4</sub> " D
W25C241315L	W25C241315-	15"L
W25C241318L	W25C241318-	18"L
W25C241321L	W25C241321-	21"L
W25C241324L	W25C241324-	24"L
		<b>16" D</b>
W25C241615L	W25C241615-	15"L
W25C241618L	W25C241618-	18"L
W25C241621L	W25C241621-	21"L
W25C241624L	W25C241624-	24"L

#### 24" High




12 <sup>3</sup> / <sub>4</sub> " Deep	16" Deep
W25C241330-	W25C241630-
W25C241336-	W25C241636-
W25C241342-	W25C241642-
W25C241348-	W25C241648-

#### 30" High




Left Hand	Right Hand	12 <sup>3</sup> / <sub>4</sub> " D
W25C301315L	W25C301315-	15"L
W25C301318L	W25C301318-	18"L
W25C301321L	W25C301321-	21"L
W25C301324L	W25C301324-	24"L
		<b>16" D</b>
W25C301615L	W25C301615-	15"L
W25C301618L	W25C301618-	18"L
W25C301621L	W25C301621-	21"L
W25C301624L	W25C301624-	24"L

#### 30" High




12 <sup>3</sup> / <sub>4</sub> " Deep	16" Deep
W25C301330-	W25C301630-
W25C301336-	W25C301636-
W25C301342-	W25C301642-
W25C301348-	W25C301648-

#### 36" High




Left Hand	Right Hand	12 <sup>3</sup> / <sub>4</sub> " D
W25C361315L	W25C361315-	15"L
W25C361318L	W25C361318-	18"L
W25C361321L	W25C361321-	21"L
W25C361324L	W25C361324-	24"L
		<b>16" D</b>
W25C361615L	W25C361615-	15"L
W25C361618L	W25C361618-	18"L
W25C361621L	W25C361621-	21"L
W25C361624L	W25C361624-	24"L

#### 36" High




12 <sup>3</sup> / <sub>4</sub> " Deep	16" Deep
W25C361330-	W25C361630-
W25C361336-	W25C361636-
W25C361342-	W25C361642-
W25C361348-	W25C361648-

#### 48" High




Left Hand	Right Hand	12 <sup>3</sup> / <sub>4</sub> " D
W25C481315L	W25C481315-	15"L
W25C481318L	W25C481318-	18"L
W25C481321L	W25C481321-	21"L
W25C481324L	W25C481324-	24"L
		<b>16" D</b>
W25C481615L	W25C481615-	15"L
W25C481618L	W25C481618-	18"L
W25C481621L	W25C481621-	21"L
W25C481624L	W25C481624-	24"L

#### 48" High



12 <sup>3</sup> / <sub>4</sub> " Deep	16" Deep
W25C481330-	W25C481630-
W25C481336-	W25C481636-
W25C481342-	W25C481642-
W25C481348-	W25C481648-

# Steel Wall Mounted Storage Cabinets




**Specifications:**

Corner Wall Mounted Storage Cabinets are available in four heights, 24", 30", 36" and 48" and two depths; 12<sup>3</sup>/<sub>4</sub>" and 16" and with either glazed or panel swinging doors. Cabinets are furnished with one, two or three adjustable shelves as shown.


**Note:** Last four digits of part number designate door style, material, and pull selection. See page 4 for full explanation. Beveled Edge, Door Style 7 is not available for wall cabinets. Wall cabinets are designed to be mounted on hanging rails. Shown on page 47 and 81.

**Swinging Door Corner Cabinets**

Glazed Door	Panel Door	Glazed Door	Panel Door
<b>24" High</b>	<b>12<sup>3</sup>/<sub>4</sub>" Deep</b>	<b>24" High</b>	<b>16" Deep</b>
W51C241320-_____	W26C241320-_____	W51C241615-_____	W26C241615-_____
<b>30" High</b>	<b>12<sup>3</sup>/<sub>4</sub>" Deep</b>	<b>30" High</b>	<b>16" Deep</b>
W51C301320-_____	W26C301320-_____	W51C301615-_____	W26C301615-_____
<b>36" High</b>	<b>12<sup>3</sup>/<sub>4</sub>" Deep</b>	<b>36" High</b>	<b>16" Deep</b>
W51C361320-_____	W26C361320-_____	W51C361615-_____	W26C361615-_____
<b>48" High</b>	<b>12<sup>3</sup>/<sub>4</sub>" Deep</b>	<b>48" High</b>	<b>16" Deep</b>
W51C481320-_____	W26C481320-_____	W51C481615-_____	W26C481615-_____




PLAN VIEW OF  
W51C241320  
W51C301320  
W51C361320  
W51C481320  
W26C241320  
W26C301320  
W26C361320  
W26C481320



PLAN VIEW OF  
W51C241615  
W51C301615  
W51C361615  
W51C481615  
W26C241615  
W26C301615  
W26C361615  
W26C481615

# Steel Full Height Storage Cabinets




### Specifications:

Full Height Storage Cabinets are 84" high and 12<sup>3</sup>/<sub>4</sub>", 16", or 22" deep. Lengths as shown. They are furnished with five adjustable shelves.

**CAUTION:** Full height cabinets may present a potential tipping over problem. Proper installation requires these cabinets to be anchored to a wall or other cabinets to improve stability and ensure smooth door movement.

**NOTE:** Full height storage cabinets are not available in Overlay – Beveled Edge Door Style 7.

### Open Storage



#### 12<sup>3</sup>/<sub>4</sub>" Deep


S00C841318-00_0	18"L
S00C841324-00_0	24"L

#### 16" Deep

S00C841618-00_0	18"L
S00C841624-00_0	24"L

#### 22" Deep

S00C842218-00_0	18"L
S00C842224-00_0	24"L



#### 12<sup>3</sup>/<sub>4</sub>" Deep

S00C841330-00_0	30"L
S00C841336-00_0	36"L
S00C841342-00_0	42"L
S00C841348-00_0	48"L


#### 16" Deep

S00C841630-00_0	30"L
S00C841636-00_0	36"L
S00C841642-00_0	42"L
S00C841648-00_0	48"L

#### 22" Deep

S00C842230-00_0	30"L
S00C842236-00_0	36"L
S00C842242-00_0	42"L
S00C842248-00_0	48"L

### Sliding Glazed Doors



#### 12<sup>3</sup>/<sub>4</sub>" Deep

S41C841330-0__0	30"L
S41C841336-0__0	36"L
S41C841342-0__0	42"L
S41C841348-0__0	48"L


#### 16" Deep

S41C841630-0__0	30"L
S41C841636-0__0	36"L
S41C841642-0__0	42"L
S41C841648-0__0	48"L

#### 22" Deep

S41C842230-0__0	30"L
S41C842236-0__0	36"L
S41C842242-0__0	42"L
S41C842248-0__0	48"L

### Sliding Panel Doors



#### 12<sup>3</sup>/<sub>4</sub>" Deep

S12C841330-0__0	30"L
S12C841336-0__0	36"L
S12C841342-0__0	42"L
S12C841348-0__0	48"L


#### 16" Deep

S12C841630-0__0	30"L
S12C841636-0__0	36"L
S12C841642-0__0	42"L
S12C841648-0__0	48"L

#### 22" Deep

S12C842230-0__0	30"L
S12C842236-0__0	36"L
S12C842242-0__0	42"L
S12C842248-0__0	48"L

# Steel Full Height Storage Cabinets




### Specifications:

Full Height Storage Cabinets are 84" high and 12<sup>3</sup>/<sub>4</sub>", 16", or 22" deep. Lengths as shown. They are furnished with five adjustable shelves.


**CAUTION:** Full height cabinets may present a potential tipping over problem. Proper installation requires these cabinets to be anchored to a wall or other cabinets to improve stability and ensure smooth door movement.

**NOTE:** Full height storage cabinets are not available in Overlay – Beveled Edge Door Style 7.

### Swinging Glazed Doors




Left Hand	Right Hand	12 <sup>3</sup> / <sub>4</sub> " D
S50C841318L	S50C841318-	18"L
S50C841324L	S50C841324-	24"L
<b>16" D</b>		
S50C841618L	S50C841618-	18"L
S50C841624L	S50C841624-	24"L
<b>22" D</b>		
S50C842218L	S50C8422-18	18"L
S50C842224L	S50C8422-24	24"L




12 <sup>3</sup> / <sub>4</sub> " Deep	
S50C841330-	30"L
S50C841336-	36"L
S50C841342-	42"L
S50C841348-	48"L
<b>16" Deep</b>	
S50C841630-	30"L
S50C841636-	36"L
S50C841642-	42"L
S50C841648-	48"L
<b>22" Deep</b>	
S50C842230-	30"L
S50C842236-	36"L
S50C842242-	42"L
S50C842248-	48"L

### Swinging Panel Doors




Left Hand	Right Hand	12 <sup>3</sup> / <sub>4</sub> " D
S25C841318L	S25C841318-	18"L
S25C841324L	S25C841324-	24"L
<b>16" D</b>		
S25C841618L	S25C841618-	18"L
S25C841624L	S25C841624-	24"L
<b>22" D</b>		
S25C842218L	S25C842218-	18"L
S25C842224L	S25C842224-	24"L



12 <sup>3</sup> / <sub>4</sub> " Deep	
S25C841330-	30"L
S25C841336-	36"L
S25C841342-	42"L
S25C841348-	48"L
<b>16" Deep</b>	
S25C841630-	30"L
S25C841636-	36"L
S25C841642-	42"L
S25C841648-	48"L
<b>22" Deep</b>	
S25C842230-	30"L
S25C842236-	36"L
S25C842242-	42"L
S25C842248-	48"L

# Steel Full Height Storage Cabinets



**Specifications:**


Microscope Storage Cabinets are 84" high and 16" or 22" deep. Lengths as shown. They are furnished with five adjustable shelves.

Full height swinging doors are hung on three hinges per door.

**CAUTION:** Full height cabinets may present a potential tipping over problem. Proper installation requires these cabinets to be anchored to a wall or other cabinets to improve stability and ensure smooth door movement.

**NOTE:** Full height storage cabinets are not available in Overlay – Beveled Edge Door Style 7.

**Microscope Storage Cabinets**




**16" Deep**

**S51C841624L** \_\_\_\_\_ 24"L

**22" Deep**

**S51C842224L** \_\_\_\_\_ 24"L

Provides 15 storage compartments




**16" Deep**

**S51C841624-** \_\_\_\_\_ 24"L

**22" Deep**

**S51C842224-** \_\_\_\_\_ 24"L

Provides 15 storage compartments




**16" Deep**

**S51C841636-** \_\_\_\_\_ 36"L

**22" Deep**

**S51C842236-** \_\_\_\_\_ 36"L

Provides 20 storage compartments



**16" Deep**


**S51C841648-** \_\_\_\_\_ 48"L

**22" Deep**

**S51C842248-** \_\_\_\_\_ 48"L

Provides 30 storage compartments

# Steel Full Height Storage Cabinets



**Specifications:**


Full Height Storage Cabinets are 84" high and 12<sup>3</sup>/<sub>4</sub>", 16", or 22" deep. Lengths as shown. Split door cabinets are furnished with a fixed center shelf and four adjustable shelves.

Wardrobe Storage Cabinets are 84" high and 22" deep. Lengths as shown. They are provided with a hat shelf and a horizontal bar for hanging garments.

**CAUTION:** Full height cabinets may present a potential tipping over problem. Proper installation requires these cabinets to be anchored to a wall or other cabinets to improve stability and ensure smooth door movement.

**NOTE:** Full height storage cabinets are not available in Overlay – Beveled Edge Door Style 7.

**Split Swinging Glazed Doors**



**12<sup>3</sup>/<sub>4</sub>" Deep**

S60C841330-	30"L
S60C841336-	36"L
S60C841342-	42"L
S60C841348-	48"L


**16" Deep**

S60C841630-	30"L
S60C841636-	36"L
S60C841642-	42"L
S60C841648-	48"L

**22" Deep**

S60C842230-	30"L
S60C842236-	36"L
S60C842242-	42"L
S60C842248-	48"L

**Split Swinging Panel Doors**



**12<sup>3</sup>/<sub>4</sub>" Deep**

S35C841330-	30"L
S35C841336-	36"L
S35C841342-	42"L
S35C841348-	48"L


**16" Deep**

S35C841630-	30"L
S35C841636-	36"L
S35C841642-	42"L
S35C841648-	48"L

**22" Deep**

S35C842230-	30"L
S35C842236-	36"L
S35C842242-	42"L
S35C842248-	48"L

**Wardrobe Storage Cabinets**




**Left Hand — 22" Deep**

S40C842218L	18"L
S40C842224L	24"L

**Right Hand — 22" Deep**


S40C842218-	18"L
S40C842224-	24"L



**22" Deep**

S40C842236-	36"L
S40C842248-	48"L

# Steel Full Height Storage Cabinets



**Specifications:**

Full Height Storage Cabinets are 84" high. Lengths as shown.


Acid Storage Cabinets are furnished with a polyethylene corrosion resistant lining and five full-width corrosion resistant shelves. Ventilation is provided by louvers at the bottom of the door and a 3 1/2" O.D. vent outlet at the top of the cabinet.

Full height swinging doors are hung on three hinges per door.


**CAUTION:** Full height cabinets may present a potential tipping over problem. Proper installation requires these cabinets to be anchored to a wall or other cabinets to improve stability and ensure smooth door movement.

**NOTE:** Full height storage cabinets are not available in Overlay – Beveled Edge Door Style 7.

**Acid Storage Cabinets**



<b>Left Hand — 16" Deep</b>	
S70C841624L	24"L
<b>Left Hand — 22" Deep</b>	
S70C842224L	24"L
<b>Right Hand — 16" Deep</b>	
S70C841624-	24"L
<b>Right Hand — 22" Deep</b>	
S70C842224-	24"L




<b>16" Deep</b>	
S70C841636-	36"L
S70C841648-	48"L
<b>22" Deep</b>	
S70C842236-	36"L
S70C842248-	48"L



# Base Cabinet Scribes & Fillers

## End Scribes



To close opening between wall and front of adjacent base cabinets. Includes detached toe space scribe and wall angle.

Available in widths from 1" - 24" in one inch increments

**Sitting Height — 29" High**

**BESC2900--00\_0**


**ADA Height — 32" High**

**BESC3200--00\_0**

**Standing Height — 35" High**

**BESC3500--00\_0**

## End Scribe Applique Panels



End Applique Panels may be used with End Scribes to match the face of doors and drawer heads on full overlay steel cabinets.

**End Scribe must be ordered separately**

**Steel Square Edge**

Widths from 3" - 24" in one inch increments

**Sitting Height — 29" High**

**BEAC2900--11\_0**

**Steel Beveled Edge**

Widths from 3" - 24" in one inch increments

**BEAC2900--1700**

**Wood Square Edge**

Widths from 3" - 24" in one inch increments

**BEAC2900--16\_0**

**ADA Height — 32" High — Steel**

**BEAC3200--11\_0**

**BEAC3200--1700**

**BEAC3200--16\_0**

**Standing Height — 35" High — Steel**


**BEAC3500--11\_0**

**BEAC3500--1700**

**BEAC3500--16\_0**

Insert width required before dash in part number to specify width of End Scribe or Applique Panel.  
Example: BESC3500**06**-0000 for 6" wide End Scribe or BEAC3500**06**-1100 for 6" wide Steel Square Edge Applique Panel

## Intermediate Fillers



To close opening between two adjacent base cabinets. Includes detached toe space filler.

Available in widths from 6" - 48" in one inch increments

**Sitting Height — 29" High**

**BIFC2900--00\_0**


**ADA Height — 32" High**

**BIFC3200--00\_0**

**Standing Height — 35" High**

**BIFC3500--00\_0**

## Intermediate Applique Panels



Intermediate Applique Panels may be used with Intermediate Fillers to match the face of doors and drawer heads on full overlay steel cabinets.

**Intermediate Filler must be ordered separately**

**Steel Square Edge**

Widths from 6" - 48" in one inch increments

**Sitting Height — 29" High**

**BIAC2900--11\_0**

**Steel Beveled Edge**

Widths from 6" - 48" in one inch increments

**BIAC2900--1700**

**Wood Square Edge**

Widths from 6" - 48" in one inch increments

**BIAC2900--16\_0**

**ADA Height — 32" High**

**BIAC3200--11\_0**

**BIAC3200--1700**

**BIAC3200--16\_0**

**Standing Height — 35" High**

**BIAC3500--11\_0**


**BIAC3500--1700**

**BIAC3500--16\_0**

Insert width required before dash in part number to specify width of Intermediate Filler or Applique Panel.  
Example: BIFC3500**12**-0000 for 12" wide Intermediate Filler or BIAC3500**12**-1100 for 12" wide Steel Square Edge Applique Panel

# Base Cabinet Scribes & Fillers

## Corner Fillers



To close opening between two adjacent base cabinets joined at right angles. Includes detached toe space filler.

Available in widths from 1" - 12" in one inch increments

**Sitting Height — 29" High**

**BCFC2900--00\_0**


**ADA Height — 32" High**

**BCFC3200--00\_0**

**Standing Height — 35" High**

**BCFC3500--00\_0**

## Corner Applique Panels



Corner Applique Panels may be used with Corner Fillers to match the face of doors and drawer heads on full overlay steel cabinets.

**Corner Filler must be ordered separately**

**Steel Square Edge**

Widths from 3" - 12" in one inch increments

**Sitting Height — 29" High**

**BCAC2900--11\_0**

**Steel Beveled Edge**

Widths from 3" - 12" in one inch increments

**BCAC2900--1700**

**Wood Square Edge**

Widths from 3" - 12" in one inch increments

**BCAC2900--16\_0**

**ADA Height — 32" High**

**BCAC3200--11\_0**

**BCAC3200--1700**

**BCAC3200--16\_0**

**Standing Height — 35" High**


**BCAC3500--11\_0**

**BCAC3500--1700**

**BCAC3500--16\_0**

Insert width required before dash in part number to specify size of Corner Filler or Applique Panel. Example: BCFC3500**04**-0000 for 4" Corner Filler or BCAC3500**04**-1100 for 4" Steel Square Edge Applique Panel

## Rear Scribes



To close opening between wall and rear of adjacent base cabinet. Includes wall mounting angle.

Widths from 1" - 24" in one inch increments

**Sitting Height — 29" High**

**BRSC2900--00\_0**


**ADA Height — 32" High**

**BRSC3200--00\_0**

**Standing Height — 35" High**

**BRSC3500--00\_0**

## Pipe Chase Fillers



To close opening between two adjacent base cabinets installed back-to-back or a cabinet and a finished back.

Widths from 1" - 24" in one inch increments

**Sitting Height — 29" High**

**BPCC2900--00\_0**


**ADA Height — 32" High**

**BPCC3200--00\_0**

**Standing Height — 35" High**

**BPCC3500--00\_0**

## Cabinet End Applique Panel



Wood Applique Panels to match exposed ends of base cabinets to cabinet faces.

**Wood Square Edge**

**Sitting Height — 29" High**

**BPAC290023-16\_0** 23" W

**BPAC290029-16\_0** 29" W

**ADA Height — 32" High**

**BPAC320023-16\_0** 23" W

**BPAC320029-16\_0** 29" W

**Standing Height — 35" High**


**BPAC350023-16\_0** 23" W

**BPAC350029-16\_0** 29" W

Insert width required before dash in part number to specify width of Rear Scribe or Pipe Chase Filler. Example: BRSC3500**08**-0000 for 8" wide Rear Scribe or BPCC3500**08**-0000 for 8" wide Pipe Chase Filler.

# Wall Cabinet Scribes & Fillers

## End Scribes




To close opening between wall and front of adjacent wall cabinet. Includes wall mounting angles.

**12<sup>3</sup>/<sub>4</sub>" Deep**      **16" Deep**  
Available in widths from 1" - 24" in one inch increments

24" High	30" High	36" High	48" High
WESC2413□□-00_0	WESC2416□□-00_0	WESC3013□□-00_0	WESC3016□□-00_0
WESC3613□□-00_0	WESC3616□□-00_0	WESC4813□□-00_0	WESC4816□□-00_0

## End Scribe Applique Panels




End Applique Panels may be used with End Scribes to match the face of doors and drawer heads on full overlay steel cabinets.

**End Scribe must be ordered separately**

Steel Square Edge	Steel Beveled Edge	Wood Square Edge
Widths from 3" - 24" in one inch increments		Widths from 3" - 24" in one inch increments
<b>24" High</b>		
WEAC2400□□-11_0	not available	WEAC2400□□-16_0
<b>30" High</b>		
WEAC3000□□-11_0	not available	WEAC3000□□-16_0
<b>36" High</b>		
WEAC3600□□-11_0	not available	WEAC3600□□-16_0
<b>48" High</b>		
WEAC4800□□-11_0	not available	WEAC4800□□-16_0

Insert width required before dash in part number to specify width of End Scribe or Applique Panel.  
Example: WESC361306-0000 for 6" wide End Scribe or WEAC360006-1100 for 6" wide Steel Square Edge Applique Panel

## Intermediate Fillers




To close opening between two adjacent wall cabinets.

**12<sup>3</sup>/<sub>4</sub>" Deep**      **16" Deep**  
Available in widths from 6" - 48" in one inch increments

24" High	30" High	36" High	48" High
WIFC2413□□-00_0	WIFC2416□□-00_0	WIFC3013□□-00_0	WIFC3016□□-00_0
WIFC3613□□-00_0	WIFC3616□□-00_0	WIFC4813□□-00_0	WIFC4816□□-00_0

## Intermediate Applique Panels



Intermediate Applique Panels may be used with Intermediate Fillers to match the face of doors and drawer heads on full overlay steel cabinets.


**Intermediate Filler must be ordered separately**

Steel Square Edge	Steel Beveled Edge	Wood Square Edge
Widths from 6" - 48" in one inch increments		Widths from 6" - 48" in one inch increments
<b>24" High</b>		
WIAC2400□□-11_0	not available	WIAC2400□□-16_0
<b>30" High</b>		
WIAC3000□□-11_0	not available	WIAC3000□□-16_0
<b>36" High</b>		
WIAC3600□□-11_0	not available	WIAC3600□□-16_0
<b>48" High</b>		
WIAC4800□□-11_0	not available	WIAC4800□□-16_0

Insert width required before dash in part number to specify width of Intermediate Filler or Applique Panel.  
Example: WIFC361312-0000 for 12" wide Intermediate Filler or WIAC360012-1100 for 12" wide Steel Square Edge Applique Panel

# Wall Cabinet Scribes & Fillers


## Corner Fillers



To close opening between two adjacent wall cabinets joined at right angles.

<b>12<sup>3</sup>/<sub>4</sub>" Deep</b>	<b>16" Deep</b>
Available in depth and widths from 1"- 12" in one inch increments.	
<b>24" High</b>	
WCFC2413□□-00_0	WCFC2416□□-00_0
<b>30" High</b>	
WCFC3013□□-00_0	WCFC3016□□-00_0
<b>36" High</b>	
WCFC3613□□-00_0	WCFC3616□□-00_0
<b>48" High</b>	
WCFC4813□□-00_0	WCFC4816□□-00_0

## Corner Applique Panels




Corner Applique Panels may be used with Corner Fillers to match the face of doors and drawer heads on full overlay steel cabinets.

**Corner Filler must be ordered separately**

<b>Steel Square Edge</b>	<b>Steel Beveled Edge</b>	<b>Wood Square Edge</b>
Widths from 3"- 12" in one inch increments		Widths from 3"- 12" in one inch increments
<b>24" High</b>		
WCAC2400□□-11_0	not available	WCAC2400□□-16_0
<b>30" High</b>		
WCAC3000□□-11_0	not available	WCAC3000□□-16_0
<b>36" High</b>		
WCAC3600□□-11_0	not available	WCAC3600□□-16_0
<b>48" High</b>		
WCAC4800□□-11_0	not available	WCAC4800□□-16_0

Insert width required before dash in part number to specify width of Corner Filler or Applique Panel.  
 Example: WCFC3613**06**-0000 for 6" wide Corner Filler or WCAC3600**06**-1100 for 6" wide Steel Square Edge Applique Panel

## Cabinet End Applique Panel




Wood Applique Panels to match exposed ends of wall cabinets to cabinet faces.

<b>Wood Square Edge</b>	
<b>24" High</b>	
WPAC240013-16_0	12 <sup>3</sup> / <sub>4</sub> " W
WPAC240016-16_0	16" W
<b>30" High</b>	
WPAC300013-16_0	12 <sup>3</sup> / <sub>4</sub> " W
WPAC300016-16_0	16" W
<b>36" High</b>	
WPAC360013-16_0	12 <sup>3</sup> / <sub>4</sub> " W
WPAC360016-16_0	16" W
<b>48" High</b>	
WPAC480013-16_0	12 <sup>3</sup> / <sub>4</sub> " W
WPAC480016-16_0	16" W

# Full Height Cabinet Scribes & Fillers

## Full Height End Scribes



To close opening between wall and front of adjacent full height cabinets. Includes detached toe space scribe and wall mounting angles.

Available in widths from 1" - 24" in one inch increments

12<sup>3</sup>/<sub>4</sub>" Deep

**FESC8413**□□-00\_0


16" Deep

**FESC8416**□□-00\_0

22" Deep

**FESC8422**□□-00\_0

## Full Height End Applique Panel



End Applique Panels may be used with End Scribes to match the face of doors and drawer heads on full overlay steel cabinets.

**End Scribe must be ordered separately**

**Steel Square Edge**

Widths from 3" - 24" in one inch increments

**FEAC8400**□□-11\_0

**Steel Beveled Edge**

not available


**Wood Square Edge**

Widths from 3" - 24" in one inch increments

**FEAC8400**□□-16\_0

Insert width required before dash in part number to specify width of End Scribe or Applique Panel.  
Example: FESC8413**06**-0000 for 6" wide End Scribe or FEAC8400**06**-1100 for 6" wide Steel Square Edge Applique Panel

## Full Height Intermediate Fillers



To close opening between two adjacent full height cabinets. Includes detached toe space filler.

Available in widths from 6" - 48" in one inch increments

12<sup>3</sup>/<sub>4</sub>" Deep

**FIFC8413**□□-00\_0


16" Deep

**FIFC8416**□□-00\_0

22" Deep

**FIFC8422**□□-00\_0

## Full Height Intermediate Applique Panel



Intermediate Applique Panels may be used with Intermediate Fillers to match the face of doors and drawer heads on full overlay steel cabinets.

**Intermediate Fillers must be ordered separately**

**Steel Square Edge**

Widths from 6" - 48" in one inch increments

**FIAC8400**□□-11\_0

**Steel Beveled Edge**

not available

**Wood Square Edge**


Widths from 6" - 48" in one inch increments

**FIAC8400**□□-16\_0

Insert width required before dash in part number to specify width of Intermediate Filler or Applique Panel.  
Example: FIFC8413**06**-0000 for 6" wide Intermediate Filler or FIAC8400**06**-1100 for 6" wide Steel Square Edge Applique Panel

# Full Height Cabinet Scribes & Fillers

## Corner Filler



To close opening between two adjacent full height cabinets joined at right angles. Includes detached toe space filler.

Available in widths from 1" - 12" in one inch increments.

12<sup>3</sup>/<sub>4</sub>" Deep

**FCFC8413**□□-00\_0


16" Deep

**FCFC8416**□□-00\_0

22" Deep

**FCFC8422**□□-00\_0

## Corner Fillers Applique Panel



Corner Applique Panels may be used with Corner Fillers to match the face of doors and drawer heads on full overlay steel cabinets.

**Corner Filler must be ordered separately**

### Steel Square Edge

Widths from 3" - 12" in one inch increments

**FCAC8400**□□-11\_0

### Steel Beveled Edge

not available


### Wood Square Edge

Widths from 3" - 12" in one inch increments

**FCAC8400**□□-16\_0

Insert width required before dash in part number to specify width of Corner Filler or Applique Panel.  
 Example: FCFC8413**06**-0000 for 6" wide Corner Filler or FCAC8400**06**-1100 for 6" wide Steel Square Edge Applique Panel

## Full Height Rear Scribe



To close opening between wall and rear of adjacent full height cabinets.

Available in widths from 1" - 24" in one inch increments

56" High


84" High

**FRSC5600**□□-00\_0

**FRSC8400**□□-00\_0

Insert width required before dash in part number to specify width of Rear Scribe.  
 Example: FRSC8400**06**-0000 for 6" wide Rear Scribe.

## Cabinet End Applique Panel



Wood Applique Panels to match exposed ends of full height cabinets to cabinet faces.

### Wood Square Edge

84" High

**FPAC840013-16\_0**

12<sup>3</sup>/<sub>4</sub>" W

**FPAC840016-16\_0**


16" W

**FPAC840022-16\_0**

22" W

# Steel Cabinet Finished Backs

## Base Cabinet Finished Back



To cover backs of base cabinets exposed to view. Back adds 1" to the depth of the base cabinet.

### Sitting Height — 29" High

BFBC290112-00_0	12"W
BFBC290115-00_0	18"W
BFBC290118-00_0	18"W
BFBC290121-00_0	21"W
BFBC290124-00_0	24"W
BFBC290130-00_0	30"W
BFBC290136-00_0	36"W
BFBC290142-00_0	42"W
BFBC290148-00_0	48"W
BFBC290160-00_0	60"W

### ADA Height — 32" High


BFBC320112-00_0	12"W
BFBC320115-00_0	15"W
BFBC320118-00_0	18"W
BFBC320121-00_0	21"W
BFBC320124-00_0	24"W
BFBC320130-00_0	30"W
BFBC320136-00_0	36"W
BFBC320142-00_0	42"W
BFBC320148-00_0	48"W
BFBC320160-00_0	60"W

### Standing Height — 35" High

BFBC350112-00_0	12"W
BFBC350115-00_0	15"W
BFBC350118-00_0	18"W
BFBC350121-00_0	21"W
BFBC350124-00_0	24"W
BFBC350130-00_0	30"W
BFBC350136-00_0	36"W
BFBC350142-00_0	42"W
BFBC350148-00_0	48"W
BFBC350152-00_0U3*	52"W
BFBC350158-00_0U3*	58"W
BFBC350160-00_0U3*	60"W
BFBC350160-00_0	60"W

\* Includes toespace notches for End Sink Cabinets

## Wall Cabinet Finished Back



To cover backs of wall cabinets exposed to view. Adds 1" to depth of cabinet.

### 24" High

WFBC240115-00_0	15"W
WFBC240118-00_0	18"W
WFBC240121-00_0	21"W
WFBC240124-00_0	24"W
WFBC240130-00_0	30"W
WFBC240136-00_0	36"W
WFBC240142-00_0	42"W
WFBC240148-00_0	48"W

### 30" High

WFBC300115-00_0	15"W
WFBC300118-00_0	18"W
WFBC300121-00_0	21"W
WFBC300124-00_0	24"W
WFBC300130-00_0	30"W
WFBC300136-00_0	36"W
WFBC300142-00_0	42"W
WFBC300148-00_0	48"W


### 36" High

WFBC360115-00_0	15"W
WFBC360118-00_0	18"W
WFBC360121-00_0	21"W
WFBC360124-00_0	24"W
WFBC360130-00_0	30"W
WFBC360136-00_0	36"W
WFBC360142-00_0	42"W
WFBC360148-00_0	48"W

### 48" High

WFBC480115-00_0	15"W
WFBC480118-00_0	18"W
WFBC480121-00_0	21"W
WFBC480124-00_0	24"W
WFBC480130-00_0	30"W
WFBC480136-00_0	36"W
WFBC480142-00_0	42"W
WFBC480148-00_0	48"W

## Full Height Finished Back




To cover backs of full height cabinets exposed to view. Adds 1" to depth of cabinet.

### 84" High

FFBC840118-00_0	18"W
FFBC840124-00_0	24"W
FFBC840130-00_0	30"W
FFBC840136-00_0	36"W
FFBC840142-00_0	42"W
FFBC840148-00_0	48"W

# Steel Adjustable Slip Scribes




## Specifications:

Steel Adjustable Slip Scribes are designed the same as standard scribes with the addition of a Slip Scribe Angle at the wall edge to dress-up the edge of the scribe.

The Steel Adjustable Slip Scribe is designed to be cut to width in the field to more easily fit the space between the wall and the end of the furniture assembly.

### Adjustable Base Cabinet Scribe



To close opening between wall and front of adjacent base cabinets. Includes detached toe space scribe, end scribe and slip scribe angle.

#### Sitting Height — 29" High

**ABSC290012-00\_0** 12"W


#### ADA Height — 32" High

**ABSC320012-00\_0** 12"W

#### Standing Height — 35" High

**ABSC350012-00\_0** 12"W

### Adjustable Wall Cabinet Scribe



To close opening between wall and front of adjacent wall cabinet. Includes end scribe and slip scribe angle.

#### 24" High

**AWSC241312-00\_0** 12<sup>3</sup>/<sub>4</sub>" D 12"W

**AWSC241612-00\_0** 16" D 12"W

#### 30" High

**AWSC301312-00\_0** 12<sup>3</sup>/<sub>4</sub>" D 12"W

**AWSC301612-00\_0** 16" D 12"W

#### 36" High

**AWSC361312-00\_0** 12<sup>3</sup>/<sub>4</sub>" D 12"W


**AWSC361612-00\_0** 16" D 12"W

#### 48" High

**AWSC481312-00\_0** 12<sup>3</sup>/<sub>4</sub>" D 12"W

**AWSC481612-00\_0** 116" D 12"W

### Adjustable Full Height Scribe



To close opening between wall and front of adjacent full height cabinets. Includes detached toe space scribe, end scribe, and slip scribe angle.

#### 84" High

**AFSC841312-00\_0** 12<sup>3</sup>/<sub>4</sub>" D 12"W

**AFSC841612-00\_0** 16" D 12"W

**AFSC842212-00\_0** 22" D 12"W

### Slip Scribe Angle




May be used with any standard scribe or filler to make it adjustable. Available in lengths from 3" to 96" in one inch increments.

**SCAC0101□-00\_0**

Insert length required before dash in part number to specify length of angle.

### Spring Clips




18 gauge spring steel clip, 3/4" x 1 1/2". May be used with any scribe or filler to make them easily removable.

**F-4995-PK** package of 10



# Steel Sloping Tops


## Sloping Tops



May be used on wall cabinets and full height cabinets to keep the tops dust free.

8 <sup>3</sup> / <sub>8</sub> " High 12 <sup>3</sup> / <sub>4</sub> " Deep	10 <sup>3</sup> / <sub>4</sub> " High 16" Deep	13 <sup>3</sup> / <sub>4</sub> " High 22" Deep	
SLTC081315-00_0	SLTC111615-00_0	SLTC142215-00_0	15"L
SLTC081318-00_0	SLTC111618-00_0	SLTC142218-00_0	18"L
SLTC081321-00_0	SLTC111621-00_0	SLTC142221-00_0	21"L
SLTC081324-00_0	SLTC111624-00_0	SLTC142224-00_0	24"L
SLTC081330-00_0	SLTC111630-00_0	SLTC142230-00_0	30"L
SLTC081336-00_0	SLTC111636-00_0	SLTC142236-00_0	36"L
SLTC081342-00_0	SLTC111642-00_0	SLTC142242-00_0	42"L
SLTC081348-00_0	SLTC111648-00_0	SLTC142248-00_0	48"L
SLTC081360-00_0	SLTC111660-00_0	SLTC142260-00_0	60"L
SLTC081372-00_0	SLTC111672-00_0	SLTC142272-00_0	72"L
SLTC081384-00_0	SLTC111684-00_0	SLTC142284-00_0	84"L
SLTC081396-00_0	SLTC111696-00_0	SLTC142296-00_0	96"L


## Sloping Top Fillers



To close opening between wall and front of adjacent sloping top. Includes wall mounting angle.

8 <sup>3</sup> / <sub>8</sub> " High 12 <sup>3</sup> / <sub>4</sub> " Deep	10 <sup>3</sup> / <sub>4</sub> " High 16" Deep	13 <sup>3</sup> / <sub>4</sub> " High 22" Deep	
STFC081301-00_0	STFC111601-00_0	STFC142201-00_0	1"W
STFC081306-00_0	STFC111606-00_0	STFC142206-00_0	6"W
STFC081312-00_0	STFC111612-00_0	STFC142212-00_0	12"W


## Corner Cabinet Sloping Tops



To be used on corner wall cabinets.

8 <sup>3</sup> / <sub>8</sub> " High 12 <sup>3</sup> / <sub>4</sub> " Deep	10 <sup>3</sup> / <sub>4</sub> " High 16" Deep	
CSTC081327-00_0	CSTC111627-00_0	26 <sup>7</sup> / <sub>8</sub> "L

## Sloping Top Corner Fillers




May be used on wall cabinets and full height cabinets when cabinets are joined at right angles.

8 <sup>3</sup> / <sub>8</sub> " High 12 <sup>3</sup> / <sub>4</sub> " Deep	10 <sup>3</sup> / <sub>4</sub> " High 16" Deep	13 <sup>3</sup> / <sub>4</sub> " High 22" Deep	
SCFC081301-00_0	SCFC111601-00_0	SCFC142201-00_0	1"W
SCFC081306-00_0	SCFC111606-00_0	SCFC142206-00_0	6"W

# Steel Ceiling & Pipe Enclosures

## Ceiling Enclosures


## Corner Ceiling Enclosures



To close opening between ceiling and top of wall or full height cabinets.

### Front Ceiling Enclosure

12" High	24" High	
FCEC120018-00_0	FCEC240018-00_0	18"L
FCEC120024-00_0	FCEC240024-00_0	24"L
FCEC120030-00_0	FCEC240030-00_0	30"L
FCEC120036-00_0	FCEC240036-00_0	36"L
FCEC120042-00_0	FCEC240042-00_0	42"L
FCEC120048-00_0	FCEC240048-00_0	48"L
FCEC120060-00_0	FCEC240060-00_0	60"L
FCEC120072-00_0	FCEC240072-00_0	72"L
FCEC120084-00_0	FCEC240084-00_0	84"L
FCEC120096-00_0	FCEC240096-00_0	96"L



To close opening between ceiling and top of corner wall cabinets.


12" High	24" High	Depth	Length
CCEC121320-00_0	CCEC241320-00_0	12 <sup>3</sup> / <sub>4</sub> "	20"
CCEC121615-00_0	CCEC241615-00_0	16"	15 <sup>3</sup> / <sub>8</sub> "

### End Ceiling Enclosure

12" High	24" High	
ECEC120013-00_0	ECEC240013-00_0	12 <sup>3</sup> / <sub>4</sub> "D
ECEC120016-00_0	ECEC240016-00_0	16"D
ECEC120022-00_0	ECEC240022-00_0	22"D


## Pipe Enclosures

### Wall Bench Enclosure



Pipe Enclosures conceal vertical service lines coming down from the ceiling to the worksurface. One side is removable for easy access with no exposed fasteners. Designed to sit on top of curb on worksurface and penetrate drop ceiling. Available in powder coated steel or stainless steel.

### Center Bench Enclosure



**U01C** [ ] [ ] [ ] [ ] -00\_0 Wall Bench Enclosure  
**U02C** [ ] [ ] [ ] [ ] -00\_0 Center Bench Enclosure

Available Heights	Available Depths	Available Lengths	Available Materials
60=60"	05=5"	12=12"	0 = Powder Coated Steel
66=66"	06=6"	18=18"	3 = Stainless Steel
72=72"	07=7"	24=24"	
78=78"	08=8"	30=30"	
84=84"	09=9"	36=36"	
90=90"	10=10"	42=42"	
96=96"		48=48"	

# Fume Hoods & Ventilated Enclosures

## Ducted Ventilated Products

Kewaunee offers a full line of Laboratory Fume Hoods to meet the needs of every lab application.

In addition to our newest line of **Supreme Air Venturi Fume Hoods**, we also offer uniquely designed

hoods for education and specialty hoods for custom applications.

### Supreme Air – General Purpose Fume Hoods



The design is based on rigid frame construction that assures a solid installation and low vibration and sound levels. Access panels are easily removable to access service lines and fittings. Radius corner posts and airfoils, plus, easy operator control of interior baffle settings assure a high level of comfort, safety and efficiency.

### Supreme Air Venturi – Next Generation Fume Hoods

The Supreme Air Venturi utilizes proven containment enhancing techniques coupled with innovative laboratory convenience features, including: adjustable LED lights, electromechanical sash stop, ultra-smooth sash system, and push-button sash management. The Supreme Air Venturi provides the safest user friendly fume hood manufactured today.

### TruView – Teaching Fume Hoods

Designed for today's classroom allowing clear line-of-sight through multiple hoods to allow easier supervision of students by the instructor. May be used as a stand-alone hood or in assemblies, both single-sided for use against the wall, or double-sided for use in the center of the laboratory.



For more information on Kewaunee's **Supreme Air Fume Hoods** contact your local Kewaunee sales representative or visit the Kewaunee website: <http://www.kewaunee.com/LaboratoryHoods/FumeHoods/SupremeAirHoods/>

For more information on Kewaunee's **Supreme Air Venturi Fume Hoods** contact your local Kewaunee sales representative or visit the Kewaunee website: <http://www.kewaunee.com/LaboratoryHoods/FumeHoods/SupremeAirVenturi/>

For more information on Kewaunee's **TruView Teaching Fume Hoods** contact your local Kewaunee sales representative or visit the Kewaunee website: <http://www.kewaunee.com/LaboratoryHoods/FumeHoods/TruViewHood/>


## Fume Hoods & Ventilated Enclosures

### Ductless Ventilated Products

Kewaunee offers a full line of filtered products to meet every protection requirement. **DetectAir** ductless fume hoods provide a safe alternative to conventional chemical fume hoods

#### **DetectAir – Ductless/Filtered Fume Hoods**

With the introduction of DetectAir, the filtered fume hood has come of age. DetectAir combines Kewaunee's traditional fume hood expertise with the filtering technology of AirClean Systems.




For more information on Kewaunee's **DetectAir - Ductless/Filtered Fume Hoods** contact your local Kewaunee sales representative or visit the Kewaunee website: <http://www.kewaunee.com/LaboratoryHoods/DuctlessHoods/>

when total exhaust is not an option. The **Interceptor Biological Safety Cabinets** are offered in both class II A2 and B2 models, (NSF Certified) for both recirculating and direct exhaust

#### **the Interceptor – Biological Safety Cabinet**

The Interceptor Class II, A2 and Class II, B2 Biological Safety Cabinets are Kewaunee's latest design. Developed with comprehensive features created to enhance both safety and ergonomics... the Interceptor leads its class.




For more information on **the Interceptor Biological Safety Cabinet** from Kewaunee, contact your local Kewaunee sales representative for visit the Kewaunee website: <http://www.kewaunee.com/LaboratoryHoods/BioSafetyCabinets/>

requirements. The **Interceptor VLF** provides ISO 14644-1 Standard Class 5 conditions, providing a sterile environment for product protection.

#### **the Interceptor – Vertical Laminar Flow Workstation**

Interceptor Vertical Laminar Flow Workstations feature a front loading HEPA filter, hinged front panel with self-support brackets, 10° angled safety glass, Minihelic gauge and T5 florescent lamp.



For more information on **the Interceptor - Vertical Laminar Flow Workstation** from Kewaunee, contact your local Kewaunee sales representative for visit the Kewaunee website: <http://www.kewaunee.com/LaboratoryHoods/LaminarFlowHoods/>

# Work Tops


Work Top material selections are among the most critical of the laboratory furniture decisions. The intended applications for your work tops dictate material choice. Surfaces could

be exposed to reagents and caustic substances. Heavy equipment may need to be supported or only light duty tasks may be performed.

Kewaunee can recommend and furnish materials fabricated to your individual needs, complete with appropriate sinks, drains and accessories.

## Kemresin® Work Tops

*manufactured exclusively by Kewaunee*



**Appearance:** Ultra-smooth non-porous surface.

**Available Colors:** Black, Grey, Putty, and Slate

**Composition:** 1" thick molded epoxy resin.

**Physical Properties:** Excellent resistance to practically all acids, alkalies and solvents. Highly resistant to heat and normal physical abuse.

**Applications:** Use wherever optimum physical and chemical resisting properties are required.

## Other Work Top Materials

MATERIAL	APPEARANCE	COMPOSITION	CHEMICAL RESISTANCE	APPLICATIONS
STAINLESS STEEL	No.4 Finish	Type 304 furnished as standard, but other types are available when specified. Type 316 supplied for Perchloric Acid fume hoods.	Good chemical resistance to a wide range of reagents.	Recommended for sterile and radiological laboratories.
HARDWOOD, OAK OR MAPLE	Natural varnish, or penetrating oil finish	Solid edge grain Maple or Oak strips, glued together electronically, and finished with baked-on clear synthetic varnish or penetrating oil.	Fair chemical resistance, but not recommended for use in chemical laboratories.	Excellent top material for Physics tables, Electronic benches, and other tables where mechanical or physical work levels are prevalent.
COMPOSITE RESIN	Hard, smooth non-glare surface	Cellulose fiber reinforced phenolic resin core with a polyurethane copolymer surface. Available in multiple colors and patterns.	Good chemical resistance to wide range of acids, alkalies and solvents.	Use wherever a high degree of chemical resistance is required. (Not recommended for fume hood worksurfaces.)
HIGH-PRESSURE LAMINATES	Hard, smooth surface	Plastic laminated surface attached with water resistant adhesive to solid particle board. Available in choice of colors and patterns.	Fair chemical resistance to wide range of acids, alkalies and solvents.	Recommended for light duty purposes.

## Work Tops

### Specifications:


The configurations shown are for Kewaunee base cabinets. Standard top thickness is 1", standard curb height is 4".

To order, select configuration and top material desired. Specify number of linear or square feet per product number.

Figurations are available in all materials. Indented tops available on special order.


Dimensions are nominal. Not all con-

### Wall Bench Tops




**Flat Top no curb** **Style B**

Standard	18"
Sizes:	23" 24" 30" 31"




**Top with Applied Curb** **Style C**

Standard	18"
Sizes:	23" 24" 30" 31"




**Corner Top no curbs** **Style K**

	<b>A</b>	<b>B</b>	<b>C</b>
Standard	33"	24"	9"
Sizes:	39"	30"	9"
	40"	31"	9"



**Flat Top & Ledge no curb** **Style E**

	<b>A</b>	<b>B</b>	<b>C</b>
Standard	30"	23"	7"
Sizes:	31"	23"	8"




**Top & Ledge with Applied curb** **Style F**

	<b>A</b>	<b>B</b>	<b>C</b>
Standard	30"	23"	7"
Sizes:	31"	23"	8"

**Corner Top with Applied Curbs** **Style L**


	<b>A</b>	<b>B</b>	<b>C</b>
Standard	33"	24"	9"
Sizes:	39"	30"	9"
	40"	31"	9"

### Center Bench Tops




**Top without Center Joint** **Style H**

Standard	54"
Sizes:	60" 62"



**Top with Center Joint** **Style I**


	<b>A</b>	<b>B</b>	<b>C</b>
Standard	54"	27"	27"
Sizes:	54"	23"	31"
	60"	30"	30"
	62"	31"	31"



**Top with Center Ledge** **Style J**


	<b>A</b>	<b>B</b>	<b>C</b>	<b>D</b>
Standard	54"	23"	8"	23"
Sizes:	60"	23"	14"	23"
	62"	23"	16"	23"

### Ledges and Shelves



**Ledge Unfinished Edges** **Style A**

Standard	8"
Sizes:	10" 12" 14" 16"



**Shelf Finished Edges** **Style N**

Standard	8"
Sizes:	10" 12" 14" 16"

# Center Bench Shelving Assemblies


## Specifications:

Shelving Support Assemblies are designed to add shelving to center benches. Each assembly consists of four uprights with slots on 1" centers and a top frame. They are 12" deep in lengths as shown.

Each assembly is furnished with mounting brackets for attaching to the counter top and cabinets below but no fasteners. Fasteners must be furnished in field depending on top material and cabinet configuration.

Shelving Assembly Support Struts are especially designed to provide added support to the Shelving Support Assemblies while allowing space for connection of plumbing or electrical connections.

## 12" Deep Shelving Support Assembly




Length	Part Number
24"	<b>F17-4824-00</b>
30"	<b>F17-4830-00</b>
36"	<b>F17-4836-00</b>
42"	<b>F17-4842-00</b>
48"	<b>F17-4848-00</b>
54"	<b>F17-4854-00</b>
60"	<b>F17-4860-00</b>
66"	<b>F17-4866-00</b>
72"	<b>F17-4872-00</b>

## Shelving Assembly Support Strut



Width	Part Number
12"	<b>K17-3612-00</b>
14"	<b>K17-3614-00</b>

## Shelving Uprights



### Powder Coated Steel Uprights

Length	Single Slot Upright	Double Slot Upright	
24"	<b>F13M000224-00</b>	<b>F12M000224-00</b>	each
30"	<b>F13M000230-00</b>	<b>F12M000230-00</b>	each
36"	<b>F13M000236-00</b>	<b>F12M000236-00</b>	each
42"	<b>F13M000242-00</b>	<b>F12M000242-00</b>	each
48"	<b>F13M000248-00</b>	<b>F12M000248-00</b>	each
72"	<b>F13M000272-00</b>	<b>F12M000272-00</b>	each
79"	<b>F13M000279-00</b>	<b>F12M000279-00</b>	each
84"	<b>F13M000284-00</b>	<b>F12M000284-00</b>	each

### Stainless Steel Uprights

Length	Single Slot Upright	Double Slot Upright	
24"	<b>F13S000224-00</b>	<b>F12S000224-00</b>	each
30"	<b>F13S000230-00</b>	<b>F12S000230-00</b>	each
36"	<b>F13S000236-00</b>	<b>F12S000236-00</b>	each
42"	<b>F13S000242-00</b>	<b>F12S000242-00</b>	each
48"	<b>F13S000248-00</b>	<b>F12S000248-00</b>	each
72"	<b>F13S000272-00</b>	<b>F12S000272-00</b>	each
79"	<b>F13S000279-00</b>	<b>F12S000279-00</b>	each
84"	<b>F13S000284-00</b>	<b>F12S000284-00</b>	each

Designed to be used with Adjustable Module Shelving Assemblies to create wall mounted shelving.

# Center Bench Shelving Assemblies

## Specifications:

Pass-Thru Shelf Assemblies for 12" Island Modules are six inches deeper than the Shelf Bracket to allow the Shelf to pass through the module core and form a continuous surface. These shelves may also be moved forward and mounted to form


a 2" cord gap between the shelves. Shelf Brackets are 11 gauge steel, adjustable on 1" centers.

**Painted Steel Shelves** – 1" thick, formed from 16 gauge cold rolled steel.


**Veneer Plywood Shelves** – 1" thick hardwood plywood with matching veneer on all exposed surfaces. Available in Oak, Maple or Beech.

**Phenolic Resin Shelves** – 1" thick phenolic resin.

## Pass-Thru Shelf Assemblies



Pass-Thru Shelves



Pass-Thru Shelves with Cord Gap

	Length	Painted Steel	Veneer Plywood	Phenolic Resin
<b>15" DEEP SHELF</b> 9" Deep Bracket	24"	F21M0924-60	F21W0924-6_	F21T0924-60
	30"	F21M0930-60	F21W0930-6_	F21T0930-60
	36"	F21M0936-60	F21W0936-6_	F21T0936-60
	42"	F21M0942-60	F21W0942-6_	F21T0942-60
	48"	F21M0948-60	F21W0948-6_	F21T0948-60
	54"	F21M0954-60		
	60"	F21M0960-60		
	66"	F21M0966-60		
<b>18" DEEP SHELF</b> 12" Deep Bracket	24"	F21M1224-60	F21W1224-6_	F21T1224-60
	30"	F21M1230-60	F21W1230-6_	F21T1230-60
	36"	F21M1236-60	F21W1236-6_	F21T1236-60
	42"	F21M1242-60	F21W1242-6_	F21T1242-60
	48"	F21M1248-60	F21W1248-6_	F21T1248-60
	54"	F21M1254-60		
	60"	F21M1260-60		
	66"	F21M1266-60		
<b>24" DEEP SHELF</b> 18" Deep Bracket	24"	F21M1824-60	F21W1824-6_	F21T1824-60
	30"	F21M1830-60	F21W1830-6_	F21T1830-60
	36"	F21M1836-60	F21W1836-6_	F21T1836-60
	42"	F21M1842-60	F21W1842-6_	F21T1842-60
	48"	F21M1848-60	F21W1848-6_	F21T1848-60
	60"	F21M1860-60		

To specify Veneer Plywood Species replace blank at end of part number with number that indicates material required.

### Veneer Plywood Material Options

- 1 = Red Oak                      -5 = FSC® Red Oak
- 2 = White Maple                -6 = FSC® White Maple
- 7 = Steamed Beech              -8 = FSC® St'md Beech

## Shelf Retaining Lips



Length	Shelf Lip	
24"	F63M0224-00	each
30"	F63M0230-00	each
36"	F63M0236-00	each
42"	F63M0242-00	each
48"	F63M0248-00	each
54"	F63M0254-00	each
60"	F63M0260-00	each
66"	F63M0266-00	each
72"	F63M0272-00	each

Designed to be added to Pass-thru and Wall Shelf Assemblies. Steel Shelf Lips may be mounted to steel,

veneer plywood, or phenolic resin shelves with supplied screws.



# Wall Shelving Assemblies

## Specifications:

**Vertical Uprights** are 1 1/4" deep and either 79" or 42" high. They are manufactured with a double-slotted face with slots on 1" centers. They are designed to be mounted between two Support Rails to hold accurate spacing.


**Support Rails** are available in 8 foot lengths and may be cut to the required length. Support Rails are notched on 6" centers to accept the Vertical Uprights.

**Support Rail Covers** are required to finish the assembly and must be

ordered in pairs to match the Vertical Upright spacing.

**Optional Center Rails** may be ordered to provide a horizontal member at counter top height. They should be sized to match the Vertical Upright spacing.

## Full Height Wall Suspension System



### Vertical Upright

**F10M0279-00** 79"H pair

### Support Rail

**F11M0196-00** 96"L each

### Support Rail Covers

**F11M0324-00** 24"L pair

**F11M0330-00** 30"L pair

**F11M0336-00** 36"L pair

**F11M0342-00** 42"L pair

**F11M0348-00** 48"L pair

**F11M0354-00** 54"L pair

**F11M0360-00** 60"L pair

**F11M0366-00** 66"L pair

**F11M0372-00** 72"L pair

### Optional Center Rail Assembly

**F11M0224-00** 24"L each

**F11M0230-00** 30"L each

**F11M0236-00** 36"L each

**F11M0242-00** 42"L each

**F11M0248-00** 48"L each


**F11M0254-00** 54"L each

**F11M0260-00** 60"L each

**F11M0266-00** 66"L each

**F11M0272-00** 72"L each

## Above Countertop Wall Suspension System



### Vertical Upright

**F10M0242-00** 42"H pair

### Support Rail

**F11M0196-00** 96"L each

### Support Rail Covers

**F11M0324-00** 24"L pair

**F11M0330-00** 30"L pair

**F11M0336-00** 36"L pair

**F11M0342-00** 42"L pair

**F11M0348-00** 48"L pair

**F11M0354-00** 54"L pair

**F11M0360-00** 60"L pair

**F11M0366-00** 66"L pair

**F11M0372-00** 72"L pair

# Wall Shelving Assemblies

## Specifications:

Adjustable Wall Shelf Assemblies consist of a pair of 11 gauge steel Shelf Brackets and a Shelf. The Shelf Bracket engages the inner slot of the double slotted Wall Vertical Upright and is adjustable on 1" centers. Shelves are available in Painted



Steel, Wood veneer plywood, or Phenolic Resin.

**Painted Steel Shelves** – 1" thick, formed from 16 gauge cold rolled steel.

**Veneer Plywood Shelves** – 1" thick hardwood plywood with matching veneer on all exposed surfaces. Available in Oak, Maple or Beech.

**Phenolic Resin Shelves** – 1" thick phenolic resin.

## Adjustable Wall Shelf Assemblies



### Optional Shelf Locking Tab

Sliding tab mounted on Shelf Bracket locks Shelf into upright to keep it from being accidentally removed.

To specify, add an **K** to the end of the Shelf Assembly part number.

Example: **F21M1236-00,K**

	Length	Painted Steel	Veneer Plywood	Phenolic Resin
<b>6" DEEP SHELF</b> 6" Deep Bracket	24"	<b>F21M0624-00</b>	<b>F21W0624-0_</b>	<b>F21T0624-00</b>
	30"	<b>F21M0630-00</b>	<b>F21W0630-0_</b>	<b>F21T0630-00</b>
	36"	<b>F21M0636-00</b>	<b>F21W0636-0_</b>	<b>F21T0636-00</b>
	42"	<b>F21M0642-00</b>	<b>F21W0642-0_</b>	<b>F21T0642-00</b>
	48"	<b>F21M0648-00</b>	<b>F21W0648-0_</b>	<b>F21T0648-00</b>
	54"	<b>F21M0654-00</b>		
	60"	<b>F21M0660-00</b>		
	66"	<b>F21M0666-00</b>		
<b>9" DEEP SHELF</b> 9" Deep Bracket	24"	<b>F21M0924-00</b>	<b>F21W0924-0_</b>	<b>F21T0924-00</b>
	30"	<b>F21M0930-00</b>	<b>F21W0930-0_</b>	<b>F21T0930-00</b>
	36"	<b>F21M0936-00</b>	<b>F21W0936-0_</b>	<b>F21T0936-00</b>
	42"	<b>F21M0942-00</b>	<b>F21W0942-0_</b>	<b>F21T0942-00</b>
	48"	<b>F21M0948-00</b>	<b>F21W0948-0_</b>	<b>F21T0948-00</b>
	54"	<b>F21M0954-00</b>		
	60"	<b>F21M0960-00</b>		
	66"	<b>F21M0966-00</b>		
<b>12" DEEP SHELF</b> 12" Deep Bracket	24"	<b>F21M1224-00</b>	<b>F21W1224-0_</b>	<b>F21T1224-00</b>
	30"	<b>F21M1230-00</b>	<b>F21W1230-0_</b>	<b>F21T1230-00</b>
	36"	<b>F21M1236-00</b>	<b>F21W1236-0_</b>	<b>F21T1236-00</b>
	42"	<b>F21M1242-00</b>	<b>F21W1242-0_</b>	<b>F21T1242-00</b>
	48"	<b>F21M1248-00</b>	<b>F21W1248-0_</b>	<b>F21T1248-00</b>
	54"	<b>F21M1254-00</b>		
	60"	<b>F21M1260-00</b>		
	66"	<b>F21M1266-00</b>		
<b>18" DEEP SHELF</b> 18" Deep Bracket	24"	<b>F21M1824-00</b>	<b>F21W1824-0_</b>	<b>F21T1824-00</b>
	30"	<b>F21M1830-00</b>	<b>F21W1830-0_</b>	<b>F21T1830-00</b>
	36"	<b>F21M1836-00</b>	<b>F21W1836-0_</b>	<b>F21T1836-00</b>
	42"	<b>F21M1842-00</b>	<b>F21W1842-0_</b>	<b>F21T1842-00</b>
	48"	<b>F21M1848-00</b>	<b>F21W1848-0_</b>	<b>F21T1848-00</b>
	54"	<b>F21M1854-00</b>		
	60"	<b>F21M1860-00</b>		
	<b>24" DEEP SHELF</b> 24" Deep Bracket	24"	<b>F21M2424-00</b>	<b>F21W2424-0_</b>
30"		<b>F21M2430-00</b>	<b>F21W2430-0_</b>	<b>F21T2430-00</b>
36"		<b>F21M2436-00</b>	<b>F21W2436-0_</b>	<b>F21T2436-00</b>
42"		<b>F21M2442-00</b>	<b>F21W2442-0_</b>	<b>F21T2442-00</b>
48"		<b>F21M2448-00</b>	<b>F21W2448-0_</b>	<b>F21T2448-00</b>
60"		<b>F21M2460-00</b>		

To specify Veneer Plywood Species replace blank at end of part number with number that indicates material required.

### Veneer Plywood Material Options

- 1 = Red Oak
- 2 = White Maple
- 7 = Steamed Beech
- 5 = FSC® Red Oak
- 6 = FSC® White Maple
- 8 = FSC® St'md Beech

# Reagent Rack & Shelving Systems

**Reagent Rack Uprights** are available in two heights and two depths. All are 1½" thick and available in either powder coated steel to match or complement the casework color or in type 304 stainless steel. 12" deep uprights are designed for


use on Center Benches with a 14" deep shelf, 10" deep uprights are designed for use on Wall Benches with a 12" deep shelf, but both can be used on either bench.

Installation hardware is included with each upright.


Spacing between upright should not exceed 48".

**Plastic Laminate Shelves** are 1" thick particle board faced on top and all four edges with black plastic laminate. A backer sheet covers the bottom.


## Typical Reagent Rack Assemblies



**Single Tier Center Bench Reagent Rack** using 18" high by 12" deep uprights and 14" deep shelving.




**Double Tier Center Bench Reagent Rack** using 18" high by 12" deep and 12" by 12" deep uprights and 14" deep shelving.



**Single Tier Wall Bench Reagent Rack** using 18" high by 10" deep uprights and 12" deep shelving.

## Reagent Rack Uprights


18" High x 12" Deep



18" H x 12" D x 1½" T

**RRU-1812-0000** P.C. Steel  
**RRU-1812-0030** St. Steel


12" High x 12" Deep



12" H x 12" D x 1½" T

**RRU-1212-0000** P.C. Steel  
**RRU-1212-0030** St. Steel


18" High x 10" Deep



18" H x 10" D x 1½" T

**RRU-1810-0000** P.C. Steel  
**RRU-1810-0030** St. Steel

12" High x 10" Deep




12" H x 10" D x 1½" T

**RRU-1210-0000** P.C. Steel  
**RRU-1210-0030** St. Steel


## Plastic Laminate Shelves

14" Deep Shelf



**PL-4814-00** 14" Deep x 48" Long  
**PL-7214-00** 14" Deep x 72" Long

12" Deep Shelf



**PL-4812-00** 12" Deep x 48" Long  
**PL-7212-00** 12" Deep x 72" Long

# Reagent Rack & Shelving Systems

**Fixed Wall Mounted Shelf Brackets** are designed to mount directly to the wall and furnished without mounting hardware.

**Adjustable Wall Mounted Shelf Brackets** are designed to mount to KF-0170-04 Shelving Uprights and


include mounting bolts and nuts.

**KF-0170-04 Shelving Uprights** are designed to hold a maximum of three shelves.


Brackets and Uprights are finished in Kewaunee's unique corrosion

resistant powder coat finish to match or complement the casework color.


## Typical Reagent Rack Assemblies



Fixed Wall Shelving




Adjustable Wall Shelving



Adjustable Wall Shelving


## Flat Top Shelf Bracket    Lipped Shelf Bracket    Book End Bracket    Shelving Upright



**Fixed Shelf Brackets**

	bracket size	shelf size
<b>M-0031-0A</b>	11" x 8"	12"
<b>M-0030-0A</b>	15" x 10"	16"
<b>M-0032-0A</b>	17" x 12"	18"


part number is for single bracket



**Fixed Shelf Brackets**

	bracket size	shelf size
<b>M-0038-0A</b>	12 <sup>1</sup> / <sub>8</sub> " x 6"	12"
<b>M-0044-0A</b>	16 <sup>1</sup> / <sub>8</sub> " x 6"	16"
<b>M-0039-0A</b>	18 <sup>1</sup> / <sub>8</sub> " x 6"	18"

part number is for single bracket



**Fixed Shelf Brackets**


	bracket size	shelf size
<b>M-0060-0A</b>	12" x 6"	12"
<b>M-0061-0A</b>	16" x 6"	16"
<b>M-0062-0A</b>	18" x 6"	18"

part number is for pair of brackets



**KF-0170-04**                      48" L


High-strength, cold-rolled steel shelf support upright allows adjustable shelf bracket to be positioned along 48" length for maximum use of space. A pair of uprights holds a maximum of three shelf brackets.



**Adjustable Shelf Brackets**

	bracket size	shelf size
<b>M-0031-AD</b>	11" x 8"	12"
<b>M-0030-AD</b>	15" x 10"	16"
<b>M-0032-AD</b>	17" x 12"	18"


part number is for single bracket



**Adjustable Shelf Brackets**

	bracket size	shelf size
<b>M-0038-AD</b>	12 <sup>1</sup> / <sub>8</sub> " x 6"	12"
<b>M-0044-AD</b>	16 <sup>1</sup> / <sub>8</sub> " x 6"	16"
<b>M-0039-AD</b>	18 <sup>1</sup> / <sub>8</sub> " x 6"	18"

part number is for single bracket




**Adjustable Shelf Brackets**

	bracket size	shelf size
<b>M-0060-AD</b>	12" x 6"	12"
<b>M-0061-AD</b>	16" x 6"	16"
<b>M-0062-AD</b>	18" x 6"	18"

part number is for pair of brackets

# Laboratory Pegboards


## Kemresin Laboratory Pegboards



**Pegboard without Drip Trough**  
29" High x 1" Thick x 20" Long

	Pegboard Color	Peg Color
<b>X-020014-BK</b>	Black	Black
<b>X-020014-GR</b>	Grey	White
<b>X-020014-PT</b>	Putty	White
<b>X-020014-SL</b>	Slate	White


**Peg Size:** 1/2" dia. x 6 1/8" long  
**Peg Material:** Polypropylene  
**Number of Pegs:** 32



**Pegboard without Drip Trough**  
29" High x 1" Thick x 32" Long

	Pegboard Color	Peg Color
<b>X-020015-BK</b>	Black	Black
<b>X-020015-GR</b>	Grey	White
<b>X-020015-PT</b>	Putty	White
<b>X-020015-SL</b>	Slate	White


**Peg Size:** 1/2" dia. x 6 1/8" long  
**Peg Material:** Polypropylene  
**Number of Pegs:** 53



**Pegboard with Drip Trough**  
30" High x 1" Thick x 20" Long

	Pegboard Color	Peg Color
<b>X-020016-BK</b>	Black	Black
<b>X-020016-GR</b>	Grey	White
<b>X-020016-PT</b>	Putty	White
<b>X-020016-SL</b>	Slate	White

**Peg Size:** 1/2" dia. x 6 1/8" long  
**Peg Material:** Polypropylene  
**Number of Pegs:** 32  
**Trough Size:** 1" high x 2 1/4" deep



**Pegboard with Drip Trough**  
30" High x 1" Thick x 32" Long

	Pegboard Color	Peg Color
<b>X-020017-BK</b>	Black	Black
<b>X-020017-GR</b>	Grey	White
<b>X-020017-PT</b>	Putty	White
<b>X-020017-SL</b>	Slate	White

**Peg Size:** 1/2" dia. x 6 1/8" long  
**Peg Material:** Polypropylene  
**Number of Pegs:** 53  
**Trough Size:** 1" high x 2 1/4" deep



**Stainless Steel Drip Trough**  
for use on  
**F-2939-00** X-020014 Pegboard


**Size:** 1 1/2"H x 2 1/2"W x 20"L  
**Material:** Type 304 Stainless Steel  
**Includes:** 1/2" OD drain outlet



**Stainless Steel Drip Trough**  
for use on  
**F-2940-00** X-020015 Pegboard

**Size:** 1 1/2"H x 2 1/2"W x 32"L  
**Material:** Type 304 Stainless Steel  
**Includes:** 1/2" OD drain outlet


## Type 304 Stainless Steel Laboratory Pegboards



**Pegboard without Drip Trough**  
30" High x 3/4" Thick x 24" Long

<b>X-020008</b>	with Open Back
<b>X-020010</b>	with Finished Back


**Peg Size:** 1/2" diameter x 6" long  
**Peg Material:** Polypropylene  
**Peg Color:** White  
**Number of Pegs:** 32



**Pegboard without Drip Trough**  
30" High x 3/4" Thick x 30" Long

<b>X-020009</b>	with Open Back
<b>X-020011</b>	with Finished Back


**Peg Size:** 1/2" diameter x 6" long  
**Peg Material:** Polypropylene  
**Peg Color:** White  
**Number of Pegs:** 50



**Pegboard with Drip Trough**  
30" High x 3/4" Thick x 24" Long

<b>X-020004</b>	with Open Back
<b>X-020006</b>	with Finished Back

**Peg Size:** 1/2" diameter x 6" long  
**Peg Material:** Polypropylene  
**Peg Color:** White  
**Number of Pegs:** 32




**Pegboard with Drip Trough**  
30" High x 3/4" Thick x 30" Long

<b>X-020005</b>	with Open Back
<b>X-020007</b>	with Finished Back

**Peg Size:** 1/2" diameter x 6" long  
**Peg Material:** Polypropylene  
**Peg Color:** White  
**Number of Pegs:** 50

# Miscellaneous Accessories


## Writing Boards & Tackboards



**Black Chalkboard**  
**F-5336-00** 18<sup>1</sup>/<sub>4</sub>"H x 24<sup>1</sup>/<sub>4</sub>"L  
**F-5337-00** 24<sup>1</sup>/<sub>4</sub>"H x 36<sup>1</sup>/<sub>4</sub>"L  
**F-5338-00** 36<sup>1</sup>/<sub>4</sub>"H x 48<sup>1</sup>/<sub>4</sub>"L

---


Surface: washable Duracite  
 Frame Material: anodized aluminum  
 Trough Material: anodized aluminum



**White Magnetic Markerboard**  
**F-5342-00** 18<sup>1</sup>/<sub>4</sub>"H x 24<sup>1</sup>/<sub>4</sub>"L  
**F-5343-00** 24<sup>1</sup>/<sub>4</sub>"H x 36<sup>1</sup>/<sub>4</sub>"L  
**F-5344-00** 36<sup>1</sup>/<sub>4</sub>"H x 48<sup>1</sup>/<sub>4</sub>"L

---


Surface: 1/4" white magnetic  
 Frame Material: grey mylar vinyl



**Green Chalkboard**  
**F-5333-00** 18<sup>1</sup>/<sub>4</sub>"H x 24<sup>1</sup>/<sub>4</sub>"L  
**F-5340-00** 24<sup>1</sup>/<sub>4</sub>"H x 36<sup>1</sup>/<sub>4</sub>"L  
**F-5341-00** 36<sup>1</sup>/<sub>4</sub>"H x 48<sup>1</sup>/<sub>4</sub>"L

---

Surface: washable Duracite  
 Frame Material: anodized aluminum  
 Trough Material: anodized aluminum




**Cork Tackboard**  
**F-5418-24** 18"H x 24"L  
**F-5424-36** 24"H x 36"L  
**F-5436-48** 36"H x 48"L

---

Surface: Resilient natural cork  
 Frame Material: anodized aluminum

## Kemresin Balance Table




30" High x 24" Deep x 35" Long  
**K7-302435-BK** Black  
**K7-302435-GR** Grey  
**K7-302435-PT** Putty  
**K7-302435-SL** Slate

---

Work Top: 2" thick Kemresin  
 Work Top Size: 2"T x 24"D x 35"L  
 Pedestal Legs: 2" thick Kemresin  
 Total Weight: 300 pounds

## Marvel Refrigerators & Freezers



**L-110** Refrigerator/Freezer  
 Size: 34<sup>1</sup>/<sub>2</sub>" H x 24" D x 23<sup>7</sup>/<sub>8</sub>" L  
 6.1 cubic feet  
 Door Swing: Right Hand  
 (User changeable with swing reversing kit)  
 Shelves: (3) Adjustable  
 Service: Single-phase,  
 60 hertz, 120 volt AC

---

**L-105** Refrigerator  
 Same as above except  
 without freezer compartment.

---


**L-115** Freezer  
 Same as above except  
 entire unit is a freezer.

## Gillis Stainless Steel Cart



**F-5349-00**  
 Overall Size: 33<sup>3</sup>/<sub>8</sub>"H x 19"D x 31"L  
 Shelf Size: (3) 18"D x 27"L  
 (1" high lip on 3 sides of  
 top & middle shelf)  
 Capacity: 400 pounds  
 Casters: (4) Black rubber  
 4" Diameter – Swivel

## Lund Key Cases




**F-5332-00**  
 Height: 11"  
 Depth: 2<sup>1</sup>/<sub>2</sub>"  
 Length: 17<sup>5</sup>/<sub>16</sub>"  
 Key Capacity: 90 Keys

---

**F-5332-01**  
 Height: 24<sup>3</sup>/<sub>4</sub>"  
 Depth: 2<sup>1</sup>/<sub>2</sub>"  
 Length: 17<sup>5</sup>/<sub>16</sub>"  
 Key Capacity: 200 Keys

# Miscellaneous Accessories

## Pipe Space Kemstruts




Floor mounted support to provide stability to rear overhang of fume hood tops or worksurfaces, and for mounting piping and drainlines. Adjustable to match variances in height.

**K12-**    **-0A**

Available Heights	Available Depths	
<i>adjustable</i>	<i>depth</i>	<i>depth</i>
<b>29</b> = 27" - 31"	<b>05</b> = 5"	<b>13</b> = 13"
<b>32</b> = 31" - 35"	<b>06</b> = 6"	<b>14</b> = 14"
<b>35</b> = 34" - 38"	<b>07</b> = 7"	<b>15</b> = 15"
	<b>08</b> = 8"	<b>16</b> = 16"
	<b>09</b> = 9"	<b>17</b> = 17"
	<b>10</b> = 10"	<b>18</b> = 18"
	<b>11</b> = 11"	<b>19</b> = 19"
	<b>12</b> = 12"	<b>20</b> = 20"

## Corner Kemstruts




**K-0020-1A** 28" - 30" H  
**K-0016-1A** 33" - 36" H

Floor and wall mounted support system for adding stability to corner tops in L-shaped furniture assemblies. Adjustable to match variances in height. Includes two support angles to attach to adjoining base cabinets.


**NOTE:** Kemstrut Corner Top Support Assemblies are not required when a corner base cabinet is used.

## Kemstrut Pipe Hangers



**K-0005-0A** Drainline Support  
**M-0051-0A** Pipe Clamp for  $\frac{3}{8}$ " -  $\frac{1}{2}$ " IPS  
**M-0052-0A** Pipe Clamp for  $\frac{3}{4}$ " - 1" IPS


## Support Angles



**L-2034-06** 6" L  
**L-2034-12** 12" L  
**L-2034-18** 18" L  
**L-2034-24** 24" L  
**L-2034-48** 48" L  
**L-2034-96** 96" L

Formed from 16 gauge cold-rolled steel, with  $1\frac{1}{2}$ " sides.

## LED Task Lighting




Part Number	Length
<b>X-030001-00</b>	300mm - 11 $\frac{7}{8}$ "
<b>X-030002-00</b>	500mm - 19 $\frac{5}{8}$ "
<b>X-030003-00</b>	1000mm - 39 $\frac{3}{8}$ "

**LED Task Lighting** consists of a thin-profile, white, LED strip, 9.5mm high by 33mm deep with a touch dimmer switch. It is UL listed and furnished with 120 VAC to 24 VDC power supply with ten foot power cord. The light includes a magnetic mounting strip for mounting to steel shelving and integral tabs for mounting to wood and plastic laminate.

# Miscellaneous Accessories

## Rod Assembly



**4740-00** Rod Assembly

Consists of two 0553-00 upright rods, one 0557-00 cross rod and two 0550-00 rod clamps.

## Support Rods



- 0549-18** 1/2" Dia. x 18" H
- 0549-24** 1/2" Dia. x 24" H
- 0549-36** 1/2" Dia. x 36" H
- 0553-00** 3/4" Dia. x 36" H

Duraluminum support rods are available in 1/2" and 3/4" diameter and have a tapered end to fit 0548-00 rod sockets.


## Cross Rods



- 0554-00** 1/2" Dia. x 39" L
- 0556-00** 1/2" Dia. x 18" L
- 0557-00** 3/4" Dia. x 39" L

Duraluminum cross rods are available in 1/2" and 3/4" diameter and have rounded ends


## Rod Clamp



**0550-00** Rod Clamp

Made of cast aluminum. Designed to clamp horizontal cross rods on vertical upright rods. Can be used for either 1/2" or 3/4" diameter rods


## Rod Sockets



**0548-00** Rod Socket

Made of chrome-plated brass. Use with 1/2" and 3/4" diameter upright rods. 1" diameter x 1 7/8" long threaded shank. Complete with locknut and washer.


## Label Holder



**F-0260-00** Label Holder

Aluminum with adhesive back for easy installation. Holds 1 1/8" H x 2 1/2" L identification card.

## Base Molding




**F-0862-04** 48" L

4" high black base molding. Made of pliable material to match variances in floor surface. Comes with high-strength adhesive applied to back.

**Note: Specify number of 4' sections required.**


## Corner Guard



**F-1368-00** Corner Guard

4" high black corner guard to finish off corner joints of base molding.

## Wall Cabinet Hanger Rail



- L01-020109-00** for 15" wall cabinet
- L01-020112-00** for 18" wall cabinet
- L01-020115-00** for 21" wall cabinet
- L01-020118-00** for 24" wall cabinet
- L01-020124-00** for 30" wall cabinet
- L01-020130-00** for 36" wall cabinet
- L01-020136-00** for 42" wall cabinet
- L01-020142-00** for 48" wall cabinet

Formed of heavy gauge cold-rolled steel. Designed to capture lip on top rear of wall cabinets for easy installation.



# General Plumbing Information

### Recommended Supply Line Pressures

**Gas** ..... 0.216 psi  
**Air** (Ordinary) ..... 15 psi  
**Air** (Air Powered Stirrer Operation) ..... 40 psi  
**Vacuum** (Ordinary) ..... 18"-22" of Mercury  
**Vacuum** (High) ..... 28.5" of Mercury  
**Hot and Cold Water** ..... Local Pressure  
**Steam** (Low Pressure) ..... 0-15 psi as required  
**Steam** (High Pressure) ..... 15-30 psi as required

### Service & Rough-in Abbreviations

**A** – Air  
**E** – Electric - 120 volt AC  
**G** – Gas  
**N** – Nitrogen  
**O** – Oxygen  
**S** – Steam  
**V** – Vacuum  
**W** – Waste  
**CW** – Cold Water  
**DW** – Distilled Water  
**HV** – Hood Vent  
**HW** – Hot Water  
**WV** – Waste Vent

### Service Fitting Index Color Coding

When ordering service fittings, specify service required, so that the proper index button(s) will be supplied.

**Air** – Orange  
**Cold Water** – Dark Green  
**Distilled Water** – White  
**Gas** – Dark Blue  
**Hot Water** – Red  
**Hydrogen** – Pink  
**Nitrogen** – Brown  
**Oxygen** – Light Green  
**Special Gases** – Light Blue  
**Steam** – Black  
**Vacuum** – Yellow

### Recommended Pipe & Drain Line Sizes

Service	Number of Outlets	Pipe Size	Pipe Material
Steam	1-8	3/8"	In accordance with governing codes
	9-16	1/2"	
	17-24	3/4"	
*Water	1-10	1/2"	In accordance with governing codes
	11-18	3/4"	
	19-24	1"	
Gas	1-12	3/8"	In accordance with governing codes
	13-18	1/2"	
	19-24	3/4"	
	25-30	1"	
Air or Vacuum	1-8	3/8"	In accordance with governing codes
	9-16	1/2"	
	17-24	3/4"	
	<b>Fittings or Filter Pumps</b>	<b>Pipe Size</b>	<b>Pipe Material</b>
Drain Lines	1-6 Fittings or 4 Pumps	1 1/2"	In accordance with governing codes
	7-10 Fittings or 8 Pumps	2"	

\* If table is equipped with a series of filter pumps (aspiration), next larger pipe size should be used.

**NOTE:** For information not supplied, or if additional details are required, contact your local Kewaunee Representative.


### Open Reagent Rack Horizontal Pipe Spacing

Water pipe lines should always be the lowest and a minimum of 9" above the working surface. Electric conduit should always be the uppermost service line. (Insulated pipe will require a greater spacing than indicated below.)

**Minimum Vertical Spacing Between Centerlines of 2 Service Pipe Lines**


**Pipe Diameters**

3/8" and 3/8" ..... 2"  
 3/8" and 1/2" ..... 2"  
 1/2" and 1/2" ..... 2"  
 1/2" and 3/4" ..... 2 1/2"  
 3/4" and 3/4" ..... 2 1/2"  
 3/4" and 1" ..... 3"  
 1" and 1" ..... 3"




### Minimum Spacing Between Fittings

These drawings show the allowable minimum spacing between service and electrical fittings on laboratory bench tops and ledges.




# Laboratory Service Fittings


## Laboratory Ball Valves (for use with air, vacuum or gases – not to be used for liquids, steam or oxygen)




**W-0242-00 Gas, Air or Vacuum Valve** is chrome plated brass. Valve provides on/off control of gas, air or vacuum at pressures up to 75 PSI.  $\frac{3}{8}$ " I.P.S. male inlet. **Specify service desired.**




**W-0286-02 Double Gas, Air or Vacuum Fitting (curb or panel mounted)** is chrome plated brass. Consists of two-way 70° outlet panel flange and two laboratory ball valves. Complete with  $\frac{3}{8}$ " I.P.S. tank nipple, locknut and washer. **Specify service desired.**




**W-0287-00 Gas, Air or Vacuum Fitting (curb or panel mounted)** is chrome plated brass. Consists of a single outlet panel flange and a laboratory ball valve. Complete with a  $\frac{3}{8}$ " I.P.S. tank nipple, locknut and washer. **Specify service desired.**




**W-0260-00 Gas, Air or Vacuum Deck Mounted Fitting** is chrome plated brass. Consists of a single outlet turret with a laboratory ball valve. Complete with a  $\frac{3}{8}$ " I.P.S. tank nipple, locknut and washer. **Specify service desired.**



**W-0263-00 Double Gas, Air or Vacuum Deck Mounted Fitting** Same as No. W-0260-00 except consists of a 2-way 90° outlet turret and two laboratory ball valves. Complete with a  $\frac{3}{8}$ " I.P.S. tank nipple, locknut and washer. **Specify service desired.**




**W-0266-00 Double Gas, Air or Vacuum Deck Mounted Fitting** Same as No. W-0260-00, except consists of a 2-way 180° outlet turret and two laboratory ball valves. Complete with a  $\frac{3}{8}$ " I.P.S. tank nipple, locknut and washer. **Specify service desired.**




**W-0272-00 Quadruple Gas, Air or Vacuum Deck Mounted Fitting** Same as W-0260-00, except consists of a 4-way outlet turret and four valves. Complete with a  $\frac{3}{8}$ " I.P.S. tank nipple, locknut and washer. **Specify service desired.**


## Needle Point Cock Valves (for use with air, vacuum or gases – not to be used for liquids or steam)



**W-0243-00 Gas, Air or Vacuum Valve** is chrome plated brass. Single needle point valve with integral straight serrated hose ends. Brass floating cone and integral brass seat. Valve provides metering control of laboratory gases at pressures up to 60 PSI.  $\frac{3}{8}$ " I.P.S. male inlet. **Specify service desired.**



**W-0400-00 Double Gas, Air or Vacuum Fitting (curb or panel mounted)** with two needle point valves with integral serrated hose ends. Brass floating cone and integral brass seat. Complete with a  $\frac{3}{8}$ " I.P.S. tank nipple, locknut and washer. **Specify service desired.**




**W-0401-00 Gas, Air or Vacuum Fitting (curb or panel mounted)** Same as No. W-0400-00 except single needle point valve. Complete with a  $\frac{3}{8}$ " I.P.S. tank nipple, locknut and washer. **Specify service desired.**

**NOTE:** If valve is to be used for oxygen service, order W-0243-01.


# Laboratory Service Fittings

## Needle Point Cock Valves


(for use with air, vacuum or gases – not to be used for liquids or steam)




**W-0259-00 Single Needle Cock** Deck Mounted Fitting  
Turret Base is chrome plated brass. Single needle point valve with integral straight serrated hose end. Complete with a  $\frac{3}{8}$ " I.P.S. tank nipple, locknut and washer.  
**Specify service desired.**



**W-0402-00 Double Needle Cock** Deck Mounted Fitting  
Turret Base is chrome plated brass. Two needle point valves with integral straight serrated hose ends at 90°. Complete with a  $\frac{3}{8}$ " I.P.S. tank nipple, locknut and washer.  
**Specify service desired.**




**W-0403-00 Double Needle Cock** Deck Mounted Fitting  
Turret Base is chrome plated brass. Same as No. W-0402-00 except integral straight serrated hose ends are at 180°. Complete with a  $\frac{3}{8}$ " I.P.S. tank nipple, locknut and washer.  
**Specify service desired.**




**W-0404-00 Quadruple Needle Cock** Deck Mounted Fitting  
Turret Base is chrome plated brass. Four needle point valves with integral straight serrated hose ends. Complete with a  $\frac{3}{8}$ " I.P.S. tank nipple, locknut and washer.  
**Specify service desired.**

## Water Valve




**W-0306-00 Water Fitting (deck or panel mounted)**  
Made of chrome plated brass. Consists of a single outlet turret and a straight water valve with removable straight serrated hose connector. Complete with a  $\frac{3}{8}$ " I.P.S. tank nipple, locknut and washer.

## Steam Valves




**W-0300-00 Steam Fitting (deck or panel mounted)**  
Made of chrome plated brass. Consists of a single outlet turret and a straight steam valve. Designed for maximum steam pressure of 20 psi. Complete with a  $\frac{3}{8}$ " I.P.S. tank nipple, locknut and washer.




**W-0301-00 Steam Fitting (curb or panel mounted)**  
Made of chrome plated brass. Consists of a single outlet panel flange and a straight steam valve. Designed for maximum steam pressure of 20 psi. Complete with a  $\frac{3}{8}$ " I.P.S. tank nipple, locknut and washer.

## Multiple Service Fitting



**W-0273-02 Deck Mounted Combination Water/Gas Fitting** Made of chrome plated brass. Water faucet has a 6" spread gooseneck with removable aerator outlet and disc type valve. Includes a two-way turret with two laboratory ball valves set at 180°. Complete with two  $\frac{3}{8}$ " I.P.S. tank nipples, locknuts and washers.  
**(Specify service desired for ball valves.)**

**W-0273-VB Deck Mounted Combination Water/Gas Fitting**  
Same as W-0273-02 except furnished with vacuum breaker type gooseneck and removable straight serrated hose connector.  
**(Specify service desired for ball valves.)**




**W-0274-VB Deck Mounted Combination Water/Gas Fitting** Made of chrome plated brass. Water faucet has a 8" spread vacuum breaker type gooseneck with removable straight serrated hose connector and disc type valve. Includes a two-way turret with two laboratory ball valves set at 180°. Complete with two  $\frac{3}{8}$ " I.P.S. tank nipples, locknuts and washers.  
**(Specify service desired for valves.)**


**W-0274-02 Deck Mounted Combination Water/Gas Fitting**  
Same as W-0274-VB except furnished with standard gooseneck and removable aerator outlet.  
**(Specify service desired for valves.)**

# Laboratory Service Fittings


## Single Water Faucet




**W-0337-00 Deck Mounted Gooseneck Water Faucet**  
 Chrome plated brass with an 8" spread rigid gooseneck. Removable straight serrated hose connector. Valve is disc type with a renewable internal assembly. Faucet body is furnished with a 3/8" I.P.S. threaded mounting shank with locknut and washer.




**W-0337-0V Deck Mounted Vacuum Breaker Gooseneck Water Faucet**  
 Same as W-0337-00 except furnished with 8" spread rigid vacuum breaker type gooseneck.




**W-0337-06 Deck Mounted Gooseneck Water Faucet**  
 Same as W-0337-00 except furnished with 6" spread rigid gooseneck




**W-0337-6V Deck Mounted Vacuum Breaker Gooseneck Water Faucet**  
 Same as W-0337-0V except furnished with 6" spread rigid vacuum breaker type gooseneck.




**W-0350-00 Right Hand Curb Mounted Gooseneck Water Faucet**  
 Made of chrome plated brass. Has an 8" spread rigid gooseneck with removable straight serrated hose connector. Valve is disc type with renewable internal assembly. Flanged inlet is complete with 1/2" I.P.S. tank nipple, locknut and washer.




**W-0350-0V Right Hand Curb Mounted Vacuum Breaker Gooseneck Water Faucet**  
 Same as W-0350-00 except supplied with vacuum breaker type gooseneck.



**W-0351-00 Left Hand Curb Mounted Gooseneck Water Faucet**  
 Same as W-0350-00 except supplied with valve handle on left.




**W-0351-0V Left Hand Curb Mounted Vacuum Breaker Gooseneck Water Faucet**  
 Same as W-0350-0V except supplied with valve handle on left.




**W-0405-00 Portable Pump Gooseneck Faucet** for laboratory demonstration table. Made of chrome plated brass with 7<sup>5/8</sup>" spread, rigid, 1/4" I.P.S. gooseneck. Faucet/pump body is complete with an integral threaded mounting shank and flanged locknut. 5/8" O.D. hose connector inlet.

# Laboratory Service Fittings


## Mixing Faucet




**W-0340-00 Deck Mounted Hot and Cold Water Mixing Gooseneck Faucet** Made of chrome plated brass with an 8" spread rigid gooseneck. Removable straight serrated hose connector. Valves are disc type with renewable internal assemblies. Faucet body is complete with an integral threaded mounting shank and a flanged locknut. Furnished with 1/2" O.D. slip joint coupling nuts for 3/8" tubing.




**W-0340-0V Deck Mounted Hot and Cold Water Mixing Vacuum Breaker Gooseneck Faucet** Same as W-0340-00 except supplied with vacuum breaker type gooseneck.




**W-0340-06 Deck Mounted Hot and Cold Water Mixing Gooseneck Faucet** Same as W-0340-00 except supplied with 6" spread rigid gooseneck.




**W-0340-6V Deck Mounted Hot and Cold Water Mixing Vacuum Breaker Gooseneck Faucet** Same as W-0340-06 except supplied with vacuum breaker type gooseneck.




**W-0333-00 Deck Mounted Hot and Cold Water Mixing Swing Spout Faucet** Made of chrome plated brass. Has a 9" spread swinging spout with anti-splash outlet. Valves are disc type with renewable internal assemblies. Faucet body is complete with an integral threaded shank and a flanged locknut. Furnished with 1/2" O.D. slip joint coupling nuts for 3/8" tubing.



**W-0354-00 Curb Mounted Hot and Cold Water Mixing Gooseneck Faucet** Made of chrome plated brass. Has an 8" spread rigid gooseneck with removable straight serrated hose connector. Valves are disc type with renewable internal assemblies. Flanged inlets are complete with 3/8" I.P.S. tank nipples, locknuts and washers on 8" centers.




**W-0354-0V Curb Mounted Hot and Cold Water Mixing Vacuum Breaker Gooseneck Faucet** Same as No. W-0354-00 except supplied with vacuum breaker type gooseneck.



**W-0355-00 Curb Mounted Hot and Cold Water Mixing Swing Spout Faucet** Made of chrome plated brass. Has a 9" spread swinging spout with anti-splash outlet. Valves are disc type with renewable internal assemblies. Flanged inlets are complete with 3/8" I.P.S. tank nipples, locknuts and washers on 8" centers.

# Laboratory Service Fittings

## Mixing Faucet




**W-0361-00 Panel Mounted Hot and Cold Water Mixing Gooseneck Faucet** Chrome plated brass. Designed for mounting on fascia type reagent rack plumbing box uprights. It has an 8" spread rigid gooseneck equipped with a removable straight serrated hose connector. Valves are disc type with renewable internal assemblies. Flanged inlets on 3" centers are complete with  $\frac{3}{8}$ " I.P.S. tank nipples, locknuts and washers.

**W-0361-0V Panel Mounted Hot and Cold Water Mixing Vacuum Breaker Gooseneck Faucet** Same as W-0361-00 except supplied with vacuum breaker type gooseneck.

**W-0825-00 Floor Mounted Foot Control Double Pedal Mixing Valve** is chrome plated brass. Has twin pedals with self-closing valves for mixing hot and cold water.  $\frac{1}{2}$ " I.P.S. female thread inlets and outlet at the rear. Inlets are on 2" centers,  $\frac{7}{8}$ " above floor. Outlet is on center  $1\frac{1}{2}$ " above floor.

**W-0833-00 Deck Mounted Water Gooseneck Fitting** is chrome plated brass. Consisting of a 8" gooseneck with mounting flange and W-0532-00 aerator.  $\frac{3}{8}$ " I.P.S. female thread inlet complete with tank nipple, locknut and washer. **Intended for use as an outlet for W-0825-00.**

## Faucet Accessories




**W-0532-00 Aerator** Anti-splash aerator can be used in lieu of the removable straight serrated hose connectors supplied with the gooseneck water faucet.  $\frac{3}{8}$ " I.P.S. male thread inlet.

**W-0510-00 Filter Pump (Aspirator)**  $5\frac{1}{2}$ " high overall. Has serrated hose connector side and bottom outlets.  $\frac{3}{8}$ " I.P.S. male thread inlet.

**W-0535-00 In-Line Vacuum Breaker** Chrome plated. Mounts between the outlet end of a fitting (gooseneck, flange, etc.) and a removable serrated hose connector.  $\frac{3}{8}$ " I.P.S. male thread inlet and female thread outlet. (Check plumbing codes for required location of vacuum breakers before ordering.)

# Laboratory Service Fittings


## Distilled Water Faucet



**W-0367-00 Right Hand Distilled Water (Curb or Panel Mounted) Faucet** Chrome plated with a hot tin lined interior. It has a self-closing handle and a disc type valve with a renewable internal assembly. The 8" spread rigid gooseneck has a removable serrated hose end. 1/2" I.P.S. flanged female thread inlet complete with tank nipple, locknut and washer.

**W-0368-00 Deck Mounted Distilled Water Gooseneck Faucet** Made of chrome plated brass. Has a hot tin lined interior. It has a self-closing handle and a disc type valve with renewable internal assembly. The 8" spread rigid gooseneck has a removable serrated hose end. Faucet body is complete with an integral threaded mounting shank and a flanged locknut. 1/2" I.P.S. male thread inlet with coupling nut and 1/4" I.P.S. tailpiece.

## Ultra-Pure Deionized Water Faucet




**W-0407-00 Distilled, Deionized and Ultra-Pure Water Faucet (deck mounted)** made of chrome plated brass. All polypropylene internal construction with 6" spread gooseneck. Furnished with 1/2" I.P.S. shank with locknut and washer. Inlet is 3/8" O.D. compression.

**W-0408-00 Distilled, Deionized and Ultra-Pure Water Faucet (panel mounted)** made of chrome plated brass. All polypropylene internal construction with 6" spread gooseneck. Furnished with 1/2" I.P.S. shank with locknut and washer. Inlet is 3/8" O.D. compression.


**W-0409-00 Distilled Deionized and Ultra-Pure Water Faucet (panel mounted)** made of chrome plated brass. All polypropylene internal construction with 8" spread gooseneck. Furnished with 1/2" I.P.S. shank with locknut and washer. Inlet is 3/8" O.D. compression.

# Laboratory Service Fittings


## Vandal-Resistant Fittings




**W-0410-VR Vandal Resistant Deck Mounted Gas, Air or Vacuum Fitting** with single vandal-resistant laboratory ball valve. Complete with a  $\frac{3}{8}$ " I.P.S. tank nipple, locknut and washer. **Specify service desired.**




**W-0411-VR Vandal Resistant Deck Mounted Double Gas, Air or Vacuum Fitting** with two vandal resistant laboratory ball valves with integral straight serrated hose ends at 90°. Complete with a  $\frac{3}{8}$ " I.P.S. tank nipple, locknut and washer. **Specify service desired.**




**W-0412-VR Vandal Resistant Deck Mounted Double Gas, Air or Vacuum Fitting** Same as W-0411-VR except at 180°. Complete with a  $\frac{3}{8}$ " I.P.S. tank nipple, locknut and washer. **Specify service desired.**




**W-0413-VR Vandal Resistant Deck Mounted Hot and Cold Mixing Faucet** with 8" rigid vacuum breaker gooseneck. Valves are disc type with renewable internal assemblies. Faucet body is complete with an integral threaded mounting shank and a flanged locknut.  $\frac{1}{2}$ " I.P.S. male thread inlets with coupling nuts for slip joint connections.




**W-0414-VR Vandal Resistant Panel Mounted Hot and Cold Mixing Faucet** with 8" rigid vacuum breaker gooseneck and mounting shanks.




**W-0415-VR Vandal Resistant Deck Mounted Water Faucet** with 8" rigid vacuum breaker gooseneck. Faucet body is furnished with a  $\frac{3}{8}$ " I.P.S. threaded mounting shank with locknut and washer.



**W-0417-VR Vandal Resistant Deck Mounted Hand-Held Drench Hose** with single spray type outlet heads, squeeze-lever operated valve and flexible stainless steel hose.  $\frac{3}{8}$ " I.P.S. male inlet.



**W-0416-VB Vandal Resistant Deck Mounted Combination Vacuum Breaker Water/Gas Fitting** Made of chrome plated base. Rigid gooseneck with disc type valve and in-line vacuum breaker. Includes two-way turret with two laboratory ball valves. Complete with two  $\frac{3}{8}$ " I.P.S. tank nipples, locknuts and washers. **Specify service desired.**




**W-0416-VR Vandal Resistant Deck Mounted Combination Water/Gas Fitting** Same as W-0416-VB except furnished without in-line vacuum breakers. **Specify service desired.**




# Laboratory Service Fittings


## Emergency Eye Wash/Shower Equipment




**W-0927-00 Panel Mounted Eye Wash Fixture** Consists of stainless steel bowl, two Gentle Spray outlet heads, a "PUSH TO OPERATE" ball valve and mounting bracket. 1/2" I.P.S. female thread inlet and 1 1/4" I.P.S. female thread outlet.




**W-0933-00 Eye/Face Wash (wall mounted)** the same as No. W-0927-00, except with four Gentle Spray outlet heads for complete coverage of eyes and face.




**W-0930-00 BARRIER-FREE (ADA) Eye/Face Wash (wall mounted)** with four Gentle Spray outlet heads, stainless steel skirt and bowl and stay-open ball valve. 1/2" I.P.S. female thread inlet and 1 1/4" I.P.S. female thread outlet.




**W-0931-00 Eye Wash (deck mounted)** with two Gentle Spray outlet heads, stay-open ball valve and 90° swivel feature (right hand mounting).




**W-0932-00 Eye Wash (deck mounted)** the same as No. W-0931-00 but (left hand mounting).




**W-0923-00 Deck Mounted Hose Spray Fitting** Consists of 8 ft. reinforced PVC hose with a squeeze-lever operated valve and spray-type outlet head with protective rubber collar. Flange, shank and locknut for attaching to table top are included. 3/8" I.P.S. male thread swivel inlet.




**W-0934-00 Deck Mounted Dual Purpose Eye Wash/Drench Hose** with two Gentle Spray outlet heads angled at 45° and squeeze valve with locking clip. Flange, shank and locknut for attaching to table top are included. 3/8" I.P.S. male thread swivel inlet.




**W-0935-00 Wall Mounted Dual Purpose Eye Wash/Drench Hose** the same as No. W-0934-00 except wall mounted.



**W-0374-01 Emergency Shower Head Assembly (horizontally mounted).** Includes 10" diameter orange ABS plastic shower head. Operated by a rigid pull-rod and easy to grip handle. Complete with an elbow, connecting pipe and ball valve. 1" I.P.S. female thread inlet.



**W-0928-00 Emergency Shower Head Assembly (vertically mounted).** The same as No. W-0374-01 except vertically mounted.




**W-0926-00 Drench Shower and Eye Wash Combination** A necessity for chemical stockrooms and other areas handling a quantity of corrosive chemicals. Free-standing 10" dia. drench shower made of orange ABS plastic. Has a 1" ball valve operated by a rigid pull-rod and easy to grip handle. Eye wash fixture with stainless steel bowl, two Gentle Spray heads and a "PUSH TO OPERATE" ball valve. Eye wash and shower operate independently. Overall height is 80". 1 1/4" I.P.S. female thread inlet and outlet.

**W-0936-00 BARRIER-FREE (ADA) Drench Shower and Eye Wash Combination** Similar to No. W-0926-00 except this unit is barrier-free and eye wash furnished with four Gentle Spray heads.

# Cupsinks & Drain Accessories

## Cupsinks



**0490-BE** Black Epoxy  
**0490-GE** Grey Epoxy  
**0490-PE** Putty Epoxy  
**0490-SE** Slate Epoxy  
 Complete with removable strainer. 6"x 3" inside dimension. 1 1/2" I.P.S. male straight thread outlet. (Overall height is 7 1/2")



**0491-BP** Black Poly  
**0491-GP** Grey Poly  
 Complete with removable strainer. 6"x 3" inside dimension. 1 1/2" I.P.S. male straight thread outlet. (Overall height is 8 1/2")



**0492-BP** Black Poly  
**0492-GP** Grey Poly  
 Complete with removable strainer. 9"x 3" inside dimension. 1 1/2" I.P.S. male straight thread outlet. (Overall height is 7")



**0493-BE** Black Epoxy  
**0493-GE** Grey Epoxy  
**0493-PE** Putty Epoxy  
**0493-SE** Slate Epoxy  
 Complete with removable type 316 stainless steel wire screen. 13 3/8" x 4 1/8" x 5 1/4" I.D. 1 1/2" I.P.S. male straight thread outlet. (Overall height is 8")

## Stainless Steel Cupsinks




**0950-0F** Stainless Steel  
 Made of type 304 stainless steel and has integral cross bars. 5" I.D. 1 1/2" I.P.S. male straight thread outlet.




**0975-0F** Stainless Steel  
 Made of type 316 stainless steel and has integral cross bars. 6"x 3" I.D. 1 1/2" I.P.S. male straight thread outlet.

## Drain Accessories



**0464-00** Polyethylene Trap and Sink Outlet  
 Adjustable riser (tailpiece) can be cut to the desired length. Sink outlet is complete with plug strainer, gasket and locknut. Has a 1 1/2" I.P.S. male straight thread outlet. Running type trap can be used for either "P" or "S" type drain connections. 1 1/2" I.P.S. male straight thread outlet.



**0469-00** Polyethylene Trap and Sink Outlet  
 Adjustable riser (tailpiece) with 4 1/2" adjustment that does not require cutting. The sink outlet is complete with plug, strainer, gasket and locknut. Has a 1 1/2" I.P.S. male straight thread outlet. Trap is the bottle "P" type, with removable bottom for cleanout. 1 1/2" I.P.S. male straight thread outlet.

## Sink Outlet




**0482-BP** Black Poly  
**0482-GP** Grey Poly  
 Complete with locknut. 1 1/2" I.P.S. male straight thread outlet. For use with epoxy resin tub sinks and troughs.

## Sink Overflow



**0494BP-04** 4" High  
**0494BP-06** 6" High  
**0494BP-08** 8" High  
**0494BP-10** 10" High  
 Sink Overflows are designed to be used with 0482-BP sink outlets. 1 1/2" diameter.

## St. Steel Sink Drain



**0995-00** Stainless Steel  
 Crumb Cup Strainer  
 Complete with gaskets, locknut and basket strainer. 1 1/2" male straight thread outlet.

# Kemresin Tub Sinks & Drain Troughs

Kemresin® Tub Sinks and Drain Troughs are made from a combination of modified epoxy resins, upgraded for maximum properties with carefully selected curing resins and inert fillers to produce a thermosetting material. They are formed in metal molds to provide one-piece units with coved corners and bottoms pitched to the drain outlet. After being removed

from the molds, they are cured at high temperatures for several hours to complete the chemical reaction, and ensure maximum chemical, impact and thermal shock resistance. Drain outlets are sized for 1 1/2" sink outlets unless otherwise specified.

Kemresin Drain Troughs are molded with no longitudinal joints. A slope of 1/8" per foot assures complete drainage, and a 7/8" cove radius

in the bottom corners is provided for ease of cleaning. Five trough sections, or any part of them, are available to be cemented together to provide Open End Troughs up to 21' long, Closed End Troughs with an end drain outlet up to 21'-8" long and Closed End Troughs with a center drain outlet up to 42'-8" long. Section TR-4 is provided with a center drain recess for a standard 1 1/2" sink outlet.

## Epoxy Resin Drain Troughs

## Epoxy Resin Tub Sinks




Closed End Trough

Open End Trough







Double Tub Sink



Center Bench End Sink with Trough Cutout



Single Tub Sink



# Kemresin Tub Sinks & Drain Troughs

## Kemresin Colors

Kemresin® Tub Sinks and Drain Troughs are available in four standard colors: Black, Grey, Putty, and Slate.

To specify color, add color code to the end of the part number.  
**Example: 1005-00-BK**

## Color Codes

**BK** = Black      **PT** = Putty  
**GR** = Grey      **SL** = Slate

## Kemresin Tub Sinks

Catalog Number	Inside Dimensions (at Top of Sink)			Drain Location		Overall Sink Height
	Length	Width	Height	A	B	
<b>1000-00-__</b>	14"	10"	6"	4"	5"	6 <sup>7</sup> / <sub>8</sub> "
<b>1001-00-__</b>	16"	8"	7"	4"	4"	8"
<b>1003-00-__</b>	18"	15"	11 <sup>13</sup> / <sub>16</sub> "	4"	4"	12 <sup>7</sup> / <sub>8</sub> "
(1) <b>1003-DT-__</b>	37"	15"	11"	4 <sup>1</sup> / <sub>2</sub> "	4 <sup>1</sup> / <sub>2</sub> "	12"
<b>1005-00-__</b>	25"	15"	10"	4"	4"	11 <sup>1</sup> / <sub>4</sub> "
<b>1006-00-__</b>	18"	15"	8"	4"	4"	9"
(2) <b>1006-DT-__</b>	36 <sup>3</sup> / <sub>4</sub> "	15"	8 <sup>1</sup> / <sub>4</sub> "	4 <sup>1</sup> / <sub>2</sub> "	4 <sup>1</sup> / <sub>2</sub> "	9"
<b>1007-00-__</b>	24"	16"	8"	4"	8"	9 <sup>1</sup> / <sub>8</sub> "
<b>1008-00-__</b>	16"	12"	8"	4"	6"	9 <sup>1</sup> / <sub>8</sub> "
<b>1011-00-__</b>	18"	6 <sup>1</sup> / <sub>2</sub> "	6"	9"	3 <sup>1</sup> / <sub>4</sub> "	7"
(3) <b>1014-00-__</b>	25"	15"	10"	4"	4"	11 <sup>1</sup> / <sub>4</sub> "
(4) <b>1015-00-__</b>	23 <sup>3</sup> / <sub>8</sub> "	13 <sup>3</sup> / <sub>4</sub> "	13 <sup>1</sup> / <sub>4</sub> "	11 <sup>11</sup> / <sub>16</sub> "	9 <sup>1</sup> / <sub>2</sub> "	17 <sup>3</sup> / <sub>4</sub> "
(5) <b>1016-00-__</b>	23 <sup>3</sup> / <sub>8</sub> "	13 <sup>3</sup> / <sub>4</sub> "	13 <sup>1</sup> / <sub>4</sub> "	11 <sup>11</sup> / <sub>16</sub> "	9 <sup>1</sup> / <sub>2</sub> "	17 <sup>3</sup> / <sub>4</sub> "
<b>1021-00-__</b>	16"	16"	7 <sup>1</sup> / <sub>2</sub> "	4"	4"	8 <sup>3</sup> / <sub>4</sub> "
<b>1022-00-__</b>	21"	18"	10"	4 <sup>1</sup> / <sub>2</sub> "	9"	11"
<b>1023-00-__</b>	21 <sup>1</sup> / <sub>2</sub> "	15 <sup>1</sup> / <sub>2</sub> "	7"	4"	4"	8 <sup>3</sup> / <sub>8</sub> "
<b>1024-00-__</b>	21 <sup>1</sup> / <sub>2</sub> "	15 <sup>1</sup> / <sub>2</sub> "	11"	4"	4"	12 <sup>1</sup> / <sub>8</sub> "
(6) <b>1024-DT-__</b>	44"	15 <sup>1</sup> / <sub>2</sub> "	11"	4"	4"	12 <sup>1</sup> / <sub>2</sub> "
<b>1025-00-__</b>	28"	20"	9 <sup>1</sup> / <sub>8</sub> "	14"	6"	10 <sup>1</sup> / <sub>8</sub> "
<b>1026-00-__</b>	35 <sup>1</sup> / <sub>2</sub> "	19 <sup>1</sup> / <sub>2</sub> "	10"	17 <sup>3</sup> / <sub>4</sub> "	9 <sup>3</sup> / <sub>4</sub> "	11 <sup>1</sup> / <sub>8</sub> "
<b>1027-00-__</b>	42"	20"	10"	21"	10"	11 <sup>1</sup> / <sub>4</sub> "
<b>1030-00-__</b>	16"	12"	8"	8"	6"	9"
<b>1031-00-__</b>	24"	16"	8"	12"	8"	9"
<b>1032-00-__</b>	29 <sup>3</sup> / <sub>4</sub> "	14 <sup>3</sup> / <sub>4</sub> "	6"	14 <sup>7</sup> / <sub>8</sub> "	7 <sup>3</sup> / <sub>8</sub> "	7"
<b>1034-00-__</b>	21 <sup>1</sup> / <sub>2</sub> "	19 <sup>1</sup> / <sub>4</sub> "	8"	10 <sup>3</sup> / <sub>4</sub> "	9 <sup>5</sup> / <sub>8</sub> "	9 <sup>1</sup> / <sub>8</sub> "
<b>1035-00-__</b>	45 <sup>3</sup> / <sub>4</sub> "	17 <sup>3</sup> / <sub>8</sub> "	12 <sup>1</sup> / <sub>2</sub> "	22 <sup>7</sup> / <sub>8</sub> "	8 <sup>3</sup> / <sub>4</sub> "	13 <sup>3</sup> / <sub>8</sub> "
<b>1036-00-__</b>	48 <sup>5</sup> / <sub>8</sub> "	24 <sup>1</sup> / <sub>2</sub> "	5 <sup>3</sup> / <sub>8</sub> "	3 <sup>3</sup> / <sub>4</sub> "	12 <sup>1</sup> / <sub>4</sub> "	6"
(7) <b>1037-00-__</b>	23 <sup>1</sup> / <sub>4</sub> "	14 <sup>1</sup> / <sub>4</sub> "	12 <sup>3</sup> / <sub>4</sub> "	11 <sup>5</sup> / <sub>8</sub> "	9 <sup>1</sup> / <sub>4</sub> "	22 <sup>3</sup> / <sub>8</sub> "
(8) <b>1038-00-__</b>	12" diameter		7 <sup>5</sup> / <sub>8</sub> "	center		8 <sup>3</sup> / <sub>8</sub> "
(9) <b>1039-00-__</b>	18" diameter		8 <sup>1</sup> / <sub>8</sub> "	center		8 <sup>7</sup> / <sub>8</sub> "
<b>SM-01-00-__</b>	12"	6"	6"	6"	3"	7"
<b>SM-02-00-__</b>	19 <sup>1</sup> / <sub>2</sub> "	16 <sup>1</sup> / <sub>2</sub> "	8"	9 <sup>3</sup> / <sub>4</sub> "	8 <sup>1</sup> / <sub>4</sub> "	9"
<b>SM-03-00-__</b>	24"	21"	8"	6"	4 <sup>1</sup> / <sub>2</sub> "	9 <sup>1</sup> / <sub>8</sub> "
<b>SM-04-00-__</b>	25 <sup>1</sup> / <sub>2</sub> "	16 <sup>1</sup> / <sub>2</sub> "	19"	4 <sup>1</sup> / <sub>2</sub> "	4 <sup>1</sup> / <sub>2</sub> "	20 <sup>1</sup> / <sub>8</sub> "
<b>SM-05-00-__</b>	28"	18"	10"	4"	4"	11 <sup>1</sup> / <sub>8</sub> "
<b>SM-06-00-__</b>	28"	18"	10"	24"	4"	11 <sup>1</sup> / <sub>8</sub> "
<b>SM-07-00-__</b>	30"	16"	10"	4"	4"	11 <sup>1</sup> / <sub>8</sub> "
<b>SM-08-00-__</b>	30"	16"	10"	26"	4"	11 <sup>1</sup> / <sub>8</sub> "
<b>SM-09-00-__</b>	30"	18"	18"	4 <sup>1</sup> / <sub>2</sub> "	4 <sup>3</sup> / <sub>4</sub> "	18 <sup>3</sup> / <sub>4</sub> "
<b>SM-10-00-__</b>	36"	21"	4 <sup>3</sup> / <sub>4</sub> "	4"	4"	5 <sup>3</sup> / <sub>4</sub> "
<b>SM-11-00-__</b>	48"	13"	12 <sup>1</sup> / <sub>4</sub> "	24"	6 <sup>1</sup> / <sub>2</sub> "	13 <sup>1</sup> / <sub>4</sub> "
<b>SM-12-00-__</b>	8"	8"	8"	4"	4"	8 <sup>7</sup> / <sub>8</sub> "

- (1) Double tub sink. Distance between centerlines of drain outlets is 19<sup>3</sup>/<sub>8</sub>".
- (2) Double tub sink. Distance between centerlines of drain outlets is 19".
- (3) Sink has 7<sup>1</sup>/<sub>8</sub>" long by 5<sup>5</sup>/<sub>8</sub>" high trough cutout.
- (4) Center table end sink has 7<sup>1</sup>/<sub>8</sub>" long by 5<sup>5</sup>/<sub>8</sub>" high trough cutout and 24<sup>1</sup>/<sub>2</sub>" long by 3<sup>3</sup>/<sub>4</sub>" high back splash. Table side of sink is not tapered.
- (5) Center table end sink has 24<sup>1</sup>/<sub>2</sub>" long by 3<sup>3</sup>/<sub>4</sub>" high back splash without trough cutout. Table side of sink is not tapered.
- (6) Double tub sink. Distance between centerlines of drain outlets is 22<sup>1</sup>/<sub>2</sub>".
- (7) Assembly of tub sink, back splash, ledge and end fillers. Used for T69W4525-00 Scrub Sink Assembly.
- (8) Cylindrical sink –12" inside diameter at top of sink with center drain outlet.
- (9) Cylindrical sink –18" inside diameter at top of sink with center drain outlet.
- (10) Sides have 1<sup>1</sup>/<sub>2</sub>" taper per foot in height.
- (11) Exterior of one 21" long side is not tapered.

Catalog Number	Inside Dimensions (at Top of Sink)			Drain Location		Overall Sink Height
	Length	Width	Height	A	B	
<b>SM-13-00-__</b>	14"	10"	8"	4"	5"	9"
<b>SM-14-00-__</b>	23 <sup>1</sup> / <sub>4</sub> "	14 <sup>1</sup> / <sub>4</sub> "	12 <sup>3</sup> / <sub>4</sub> "	11 <sup>5</sup> / <sub>8</sub> "	9 <sup>1</sup> / <sub>4</sub> "	13 <sup>5</sup> / <sub>8</sub> "
<b>SM-15-00-__</b>	19 <sup>3</sup> / <sub>4</sub> "	18"	24"	9 <sup>7</sup> / <sub>8</sub> "	9"	25"
<b>SM-16-00-__</b>	40"	18"	18 <sup>1</sup> / <sub>4</sub> "	20"	9"	19 <sup>1</sup> / <sub>4</sub> "
<b>SM-17-00-__</b>	29 <sup>3</sup> / <sub>4</sub> "	23 <sup>1</sup> / <sub>2</sub> "	6"	14 <sup>7</sup> / <sub>8</sub> "	11 <sup>3</sup> / <sub>4</sub> "	7"
<b>SM-18-00-__</b>	18"	15"	15 <sup>1</sup> / <sub>4</sub> "	4 <sup>3</sup> / <sub>4</sub> "	4 <sup>3</sup> / <sub>4</sub> "	16 <sup>1</sup> / <sub>2</sub> "
<b>SM-19-00-__</b>	35 <sup>1</sup> / <sub>2</sub> "	19 <sup>1</sup> / <sub>2</sub> "	15 <sup>3</sup> / <sub>8</sub> "	3 <sup>1</sup> / <sub>2</sub> "	4 <sup>1</sup> / <sub>2</sub> "	16 <sup>1</sup> / <sub>2</sub> "
<b>SM-20-00-__</b>	21"	18"	15 <sup>1</sup> / <sub>4</sub> "	4 <sup>1</sup> / <sub>2</sub> "	4 <sup>1</sup> / <sub>2</sub> "	16 <sup>1</sup> / <sub>2</sub> "
<b>SM-21-00-__</b>	60 <sup>1</sup> / <sub>2</sub> "	22"	7 <sup>1</sup> / <sub>4</sub> "	6"	6"	8 <sup>1</sup> / <sub>8</sub> "
<b>SM-22-00-__</b>	24"	16"	15 <sup>3</sup> / <sub>8</sub> "	12"	8"	16 <sup>1</sup> / <sub>8</sub> "
<b>SM-23-00-__</b>	28"	18"	15 <sup>1</sup> / <sub>2</sub> "	5 <sup>1</sup> / <sub>2</sub> "	4"	16 <sup>1</sup> / <sub>2</sub> "
<b>SM-24-00-__</b>	48"	8"	8"	24"	4"	9"
<b>SM-25-00-__</b>	21"	18"	15 <sup>1</sup> / <sub>8</sub> "	10 <sup>1</sup> / <sub>2</sub> "	9"	16 <sup>1</sup> / <sub>8</sub> "
<b>SM-26-00-__</b>	21"	11 <sup>3</sup> / <sub>8</sub> "	15 <sup>1</sup> / <sub>8</sub> "	10 <sup>1</sup> / <sub>2</sub> "	5 <sup>1</sup> / <sub>16</sub> "	16 <sup>1</sup> / <sub>8</sub> "
<b>SM-27-00-__</b>	36"	20"	16 <sup>1</sup> / <sub>4</sub> "	18"	10"	17 <sup>1</sup> / <sub>4</sub> "
<b>SM-28-00-__</b>	30"	18"	4 <sup>1</sup> / <sub>8</sub> "	15"	9"	5 <sup>1</sup> / <sub>4</sub> "
<b>SM-29-00-__</b>	36"	24"	12 <sup>1</sup> / <sub>4</sub> "	18"	12"	13 <sup>1</sup> / <sub>4</sub> "
<b>SM-30-00-__</b>	40"	30"	8 <sup>1</sup> / <sub>8</sub> "	20"	15"	9"
<b>SM-31-00-__</b>	24"	15 <sup>3</sup> / <sub>4</sub> "	12 <sup>1</sup> / <sub>8</sub> "	4 <sup>3</sup> / <sub>4</sub> "	8"	13"
<b>SM-32-00-__</b>	59"	21"	12 <sup>5</sup> / <sub>8</sub> "	29 <sup>1</sup> / <sub>2</sub> "	10 <sup>1</sup> / <sub>2</sub> "	13 <sup>3</sup> / <sub>8</sub> "
<b>SM-33-00-__</b>	24"	18"	11"	4"	9"	11 <sup>7</sup> / <sub>8</sub> "
(10) <b>SM-34-00-__</b>	15"	14"	10 <sup>1</sup> / <sub>2</sub> "	7 <sup>1</sup> / <sub>2</sub> "	7"	11 <sup>1</sup> / <sub>4</sub> "
<b>SM-35-00-__</b>	51"	17"	12 <sup>1</sup> / <sub>2</sub> "	25 <sup>1</sup> / <sub>2</sub> "	8 <sup>1</sup> / <sub>2</sub> "	13 <sup>1</sup> / <sub>4</sub> "
<b>SM-36-00-__</b>	14"	10"	6"	7"	5"	6 <sup>7</sup> / <sub>8</sub> "
(11) <b>SM-37-00-__</b>	24"	20"	10 <sup>5</sup> / <sub>8</sub> "	6 <sup>3</sup> / <sub>4</sub> "	10"	11 <sup>1</sup> / <sub>4</sub> "
<b>SM-39-00-__</b>	35"	6 <sup>5</sup> / <sub>8</sub> "	4 <sup>3</sup> / <sub>4</sub> "	17 <sup>1</sup> / <sub>2</sub> "	3 <sup>5</sup> / <sub>16</sub> "	5 <sup>1</sup> / <sub>2</sub> "
<b>SM-40-00-__</b>	33 <sup>3</sup> / <sub>8</sub> "	29 <sup>1</sup> / <sub>8</sub> "	15 <sup>3</sup> / <sub>4</sub> "	6"	14 <sup>9</sup> / <sub>16</sub> "	16 <sup>1</sup> / <sub>2</sub> "
<b>SM-41-00-__</b>	18"	15"	8 <sup>1</sup> / <sub>4</sub> "	9"	7 <sup>1</sup> / <sub>2</sub> "	9"
<b>SM-42-00-__</b>	25"	15"	10 <sup>1</sup> / <sub>2</sub> "	12 <sup>1</sup> / <sub>2</sub> "	7 <sup>1</sup> / <sub>2</sub> "	11 <sup>1</sup> / <sub>4</sub> "
<b>SM-43-00-__</b>	21"	18"	10 <sup>3</sup> / <sub>8</sub> "	10 <sup>1</sup> / <sub>2</sub> "	9"	11 <sup>1</sup> / <sub>8</sub> "
<b>SM-44-00-__</b>	36"	16"	11"	12"	4"	12"
<b>SM-45-00-__</b>	27 <sup>3</sup> / <sub>4</sub> "	14 <sup>3</sup> / <sub>4</sub> "	8 <sup>1</sup> / <sub>2</sub> "	10 <sup>1</sup> / <sub>2</sub> "	7 <sup>1</sup> / <sub>2</sub> "	9 <sup>3</sup> / <sub>8</sub> "
<b>SM-47-00-__</b>	43"	15 <sup>1</sup> / <sub>2</sub> "	11 <sup>1</sup> / <sub>8</sub> "	21 <sup>1</sup> / <sub>2</sub> "	7 <sup>3</sup> / <sub>4</sub> "	12"
<b>SM-48-00-__</b>	56"	24"	8 <sup>1</sup> / <sub>4</sub> "	28"	12"	9"
<b>SM-49-00-__</b>	36"	16"	8 <sup>3</sup> / <sub>8</sub> "	6"	8"	9"
<b>SM-50-00-__</b>	24"	24"	12 <sup>1</sup> / <sub>2</sub> "	12"	12"	13 <sup>1</sup> / <sub>4</sub> "
<b>SM-51-00-__</b>	30"	15 <sup>7</sup> / <sub>8</sub> "	16 <sup>3</sup> / <sub>8</sub> "	6"	7 <sup>15</sup> / <sub>16</sub> "	17 <sup>1</sup> / <sub>8</sub> "
<b>SM-52-00-__</b>	45 <sup>1</sup> / <sub>8</sub> "	11 <sup>1</sup> / <sub>2</sub> "	8 <sup>1</sup> / <sub>8</sub> "	22 <sup>9</sup> / <sub>16</sub> "	5 <sup>3</sup> / <sub>4</sub> "	9"
<b>SM-53-00-__</b>	31 <sup>1</sup> / <sub>2</sub> "	16"	16 <sup>1</sup> / <sub>2</sub> "	7 <sup>1</sup> / <sub>8</sub> "	8"	17 <sup>1</sup> / <sub>4</sub> "
<b>SM-55-00-__</b>	16"	16"	12 <sup>3</sup> / <sub>8</sub> "	4 <sup>1</sup> / <sub>4</sub> "	4 <sup>1</sup> / <sub>4</sub> "	13"
<b>SM-56-00-__</b>	42 <sup>1</sup> / <sub>4</sub> "	18 <sup>1</sup> / <sub>8</sub> "	12 <sup>5</sup> / <sub>8</sub> "	21 <sup>1</sup> / <sub>8</sub> "	9 <sup>1</sup> / <sub>16</sub> "	13 <sup>3</sup> / <sub>8</sub> "

### General Notes:

- Sinks with center drain outlets, except 1011-00, can be furnished with a drain recess sized for a crumb cup strainer or garbage disposal when specified.
- Any sink may be furnished with trough cutout. Specify location and size. (Trough location must be at least 1<sup>1</sup>/<sub>2</sub>" from inside corner and 2" from bottom.)
- Most sinks can be joined together to form multiple tub sinks provided adjoining sides of sinks are the same size. Sinks longer than 30", wider than 24" or deeper than 12" cannot be factory joined, but they can be factory prepared for field assembly. Assemblies of three or more tub sinks should be no longer than 60" overall. For longer assemblies, one factory joint can be made with a second joint prepared for field assembly.





# Kemresin Drop-in Sinks

Kemresin Drop-in Sinks are formed in metal molds to provide a lipped one-piece tub with covered corners and bottoms pitched to the drain outlet. They are made from a combination of modified epoxy resins, upgraded for maximum properties with carefully selected curing resins and inert fillers to produce a thermosetting


material. After being removed from the molds, they are cured at high temperatures for several hours to complete the chemical reaction, and ensure maximum chemical, impact and thermal shock resistance. Drain outlets are sized for 1/2" sink outlets unless otherwise specified.

They are designed to be mounted in Kemresin tops in a specially rabbetted cutout so that no support from below is required.

## Drop-in Kemresin Sinks



PLAN VIEW  
**1048-DI-**



SECTION VIEW  
**1049-DI-**

Part Number	I.D.		Depth	O.D.		Drain Location	
	A	B	C	D	E	G	H
<b>1000-DI-</b> ___	13 <sup>7</sup> / <sub>8</sub> "	9 <sup>7</sup> / <sub>8</sub> "	6"	16 <sup>1</sup> / <sub>8</sub> "	12 <sup>1</sup> / <sub>8</sub> "	7"	5"
<b>1001-DI-</b> ___	8"	16"	7"	10 <sup>1</sup> / <sub>4</sub> "	18 <sup>1</sup> / <sub>4</sub> "	4"	4"
<b>1003-DI-</b> ___	17 <sup>15</sup> / <sub>16</sub> "	14 <sup>15</sup> / <sub>16</sub> "	11"	20 <sup>3</sup> / <sub>16</sub> "	17 <sup>3</sup> / <sub>16</sub> "	3 <sup>1</sup> / <sub>2</sub> "	3 <sup>1</sup> / <sub>2</sub> "
<b>1005-DI-</b> ___	25"	15"	10"	27 <sup>1</sup> / <sub>4</sub> "	17 <sup>1</sup> / <sub>4</sub> "	4"	4"
<b>1006-DI-</b> ___	17 <sup>7</sup> / <sub>8</sub> "	14 <sup>7</sup> / <sub>8</sub> "	8"	20 <sup>1</sup> / <sub>8</sub> "	17 <sup>1</sup> / <sub>8</sub> "	3 <sup>1</sup> / <sub>2</sub> "	3 <sup>1</sup> / <sub>2</sub> "
<b>1008-DI-</b> ___	16"	12"	8"	18 <sup>1</sup> / <sub>4</sub> "	14 <sup>1</sup> / <sub>4</sub> "	4"	6"
<b>1011-DI-</b> ___	18"	6 <sup>1</sup> / <sub>2</sub> "	6"	20 <sup>1</sup> / <sub>4</sub> "	8 <sup>3</sup> / <sub>4</sub> "	9"	3 <sup>1</sup> / <sub>4</sub> "
<b>1021-DI-</b> ___	16"	16"	7 <sup>1</sup> / <sub>2</sub> "	18 <sup>1</sup> / <sub>4</sub> "	18 <sup>1</sup> / <sub>4</sub> "	4"	4"
<b>1022-DI-</b> ___	21"	18"	10"	23 <sup>1</sup> / <sub>4</sub> "	20 <sup>1</sup> / <sub>4</sub> "	4"	9"
(1) <b>1038-DI-</b> ___	12" diameter		7 <sup>5</sup> / <sub>8</sub> "	14 <sup>1</sup> / <sub>4</sub> " diameter		center center	
(2) <b>1039-DI-</b> ___	18" diameter		8 <sup>1</sup> / <sub>8</sub> "	20 <sup>1</sup> / <sub>4</sub> " diameter		center center	
<b>1040-DI-</b> ___	15 <sup>3</sup> / <sub>8</sub> "	11 <sup>3</sup> / <sub>8</sub> "	8"	17.6"	13.6"	4"	4"
<b>1041-DI-</b> ___	23 <sup>3</sup> / <sub>8</sub> "	15 <sup>3</sup> / <sub>8</sub> "	8.3"	25.6"	17.6"	4"	4"
<b>1042-DI-</b> ___	28"	15"	11.8"	30 <sup>1</sup> / <sub>4</sub> "	17 <sup>1</sup> / <sub>4</sub> "	5"	5"
<b>1043-DI-</b> ___	30"	16"	17.8"	32 <sup>1</sup> / <sub>4</sub> "	18 <sup>1</sup> / <sub>4</sub> "	5"	8"
<b>1044-DI-</b> ___	12"	8"	6"	14 <sup>1</sup> / <sub>4</sub> "	10 <sup>1</sup> / <sub>4</sub> "	6"	4"
<b>1045-DI-</b> ___	21"	17"	10"	23 <sup>1</sup> / <sub>4</sub> "	19 <sup>1</sup> / <sub>4</sub> "	4 <sup>1</sup> / <sub>2</sub> "	8 <sup>1</sup> / <sub>2</sub> "
<b>1046-DI-</b> ___	21 <sup>1</sup> / <sub>2</sub> "	15 <sup>1</sup> / <sub>2</sub> "	11"	23 <sup>3</sup> / <sub>4</sub> "	17 <sup>3</sup> / <sub>4</sub> "	4"	4"
<b>1047-DI-</b> ___	24"	18"	11"	26 <sup>1</sup> / <sub>4</sub> "	20 <sup>1</sup> / <sub>4</sub> "	4"	4"
(3) <b>1048-DI-</b> ___	17 <sup>1</sup> / <sub>2</sub> "	12 <sup>3</sup> / <sub>8</sub> "	6 <sup>1</sup> / <sub>2</sub> "	19 <sup>1</sup> / <sub>2</sub> "	14 <sup>5</sup> / <sub>8</sub> "	8 <sup>3</sup> / <sub>4</sub> "	6 <sup>3</sup> / <sub>16</sub> "
(4) <b>1049-DI-</b> ___	18"	15"	11"	20 <sup>1</sup> / <sub>4</sub> "	17 <sup>1</sup> / <sub>4</sub> "	4 <sup>1</sup> / <sub>2</sub> "	5"
<b>1070-DI-</b> ___	24"	16"	15 <sup>1</sup> / <sub>2</sub> "	26 <sup>1</sup> / <sub>4</sub> "	18 <sup>1</sup> / <sub>4</sub> "	12"	8"
ADA Compliant							
<b>1000-ADA-DI-</b> ___	13 <sup>7</sup> / <sub>8</sub> "	9 <sup>7</sup> / <sub>8</sub> "	5"	16 <sup>1</sup> / <sub>8</sub> "	12 <sup>1</sup> / <sub>8</sub> "	4"	4"
<b>1005-ADA-DI-</b> ___	25"	15"	4.8"	27 <sup>1</sup> / <sub>4</sub> "	17 <sup>1</sup> / <sub>4</sub> "	4"	4"
<b>1006-ADA-DI-</b> ___	17 <sup>7</sup> / <sub>8</sub> "	14 <sup>7</sup> / <sub>8</sub> "	5"	20 <sup>1</sup> / <sub>8</sub> "	17 <sup>1</sup> / <sub>8</sub> "	3 <sup>1</sup> / <sub>2</sub> "	3 <sup>1</sup> / <sub>2</sub> "

(1) Cylindrical sink, 12" inside diameter  
(2) Cylindrical sink, 18" inside diameter

(3) Triangular sink used in X92 and X94 assemblies  
(4) Stepped sink, 4<sup>3</sup>/<sub>4</sub>" deep at front

## Kemresin Colors

Kemresin® Drop-in Sinks are available in four standard colors: Black, Grey, Putty, and Slate.

To specify color, add color code to the end of the part number.  
**Example: 1005-DI-BK**

## Color Codes

**BK** = Black      **PT** = Putty  
**GR** = Grey      **SL** = Slate

# Research Collection

## Stainless Steel Drop-in Sinks


Seamless die-drawn from nickel bearing type 304, 18-8, 18 gauge stainless steel. Allows easy cleaning and stain resistance. Interior and top surfaces polished to a gleaming finish with a highlighted bowl rim. Corners and bottom turns are radiused  $1\frac{3}{4}$ ". Sinks are thoroughly sound deadened with environmentally

safe compounds, to ensure proper sound insulation and to reduce condensation.

All sinks are self-rimming with stainless steel mounting channels for easy, top mount, snap-in clip installation. Sinks are compatible with any type work-surface material from  $\frac{3}{4}$ " to  $2\frac{1}{2}$ " thick.


Each sink is center punched to receive a standard cup strainer. ADA compliant sinks are punched in the center,  $4\frac{1}{2}$ " from the rear of the bowl for extra clearance.

### Stainless Steel Sinks



Part Number	O.D.		I.D.		Sink Depth
	Front to Back	Side to Side	Front to Back	Side to Side	
F-5353-08	15"	15"	12"	12"	7 $\frac{1}{2}$ "
F-5354-08	17"	19"	14"	16"	7 $\frac{1}{2}$ "
F-5355-08	19"	22"	16"	19"	7 $\frac{1}{2}$ "
F-5356-08	19"	25"	16"	22"	7 $\frac{1}{2}$ "
F-5357-08	19"	31"	16"	28"	7 $\frac{1}{2}$ "
F-5353-09	15"	15"	12"	12"	9"
F-5354-11	17"	19"	14"	16"	10 $\frac{1}{2}$ "
F-5355-11	19"	22"	16"	19"	10 $\frac{1}{2}$ "
F-5356-11	19"	25"	16"	22"	10 $\frac{1}{2}$ "
F-5357-11	19"	31"	16"	28"	10 $\frac{1}{2}$ "
F-5354-12	17"	19"	14"	16"	12"
F-5355-12	19"	22"	16"	19"	12"
F-5356-12	19"	25"	16"	22"	12"
F-5357-12	19"	31"	16"	28"	12"


### Stainless Steel Sinks with Fixture Ledge



Part Number	O.D.		I.D.		Sink Depth	Fixture Centers
	Front to Back	Side to Side	Front to Back	Side to Side		
F-5358-08	18"	15"	12"	12"	7 $\frac{1}{2}$ "	4"
F-5359-08	20"	19"	14"	16"	7 $\frac{1}{2}$ "	6"
F-5360-08	22"	22"	16"	19"	7 $\frac{1}{2}$ "	6"
F-5361-08	22"	25"	16"	22"	8"	8"
F-5362-08	22"	31"	16"	28"	7 $\frac{1}{2}$ "	8"
F-5358-09	18"	15"	12"	12"	9"	4"
F-5359-11	20"	19"	14"	16"	10 $\frac{1}{2}$ "	6"
F-5360-11	22"	22"	16"	19"	10 $\frac{1}{2}$ "	6"
F-5361-11	22"	25"	16"	22"	10 $\frac{1}{2}$ "	8"
F-5362-11	22"	31"	16"	28"	10 $\frac{1}{2}$ "	8"
F-5360-12	22"	22"	16"	19"	12"	6"
F-5361-12	22"	25"	16"	22"	12"	8"
F-5362-12	22"	31"	16"	28"	12"	8"


### ADA Compliant

Part Number	O.D.		I.D.		Sink Depth
	Front to Back	Side to Side	Front to Back	Side to Side	
F-5353-06	15"	15"	12"	12"	5 $\frac{1}{2}$ "
F-5354-06	17"	19"	14"	16"	5 $\frac{1}{2}$ "
F-5355-06	19"	22"	16"	19"	5 $\frac{1}{2}$ "
F-5356-06	19"	25"	16"	22"	5 $\frac{1}{2}$ "
F-5357-06	19"	31"	16"	28"	5 $\frac{1}{2}$ "



### ADA Compliant with Fixture Ledge





Part Number	O.D.		I.D.		Sink Depth	Fixture Centers
	Front to Back	Side to Side	Front to Back	Side to Side		
F-5358-06	18"	15"	12"	12"	5 $\frac{1}{2}$ "	4"
F-5359-06	20"	19"	14"	16"	5 $\frac{1}{2}$ "	6"
F-5360-06	22"	22"	16"	19"	5 $\frac{1}{2}$ "	6"
F-5361-06	22"	25"	16"	22"	5 $\frac{1}{2}$ "	8"
F-5362-06	22"	31"	16"	28"	5 $\frac{1}{2}$ "	8"










# Laboratory Electrical Fixtures

## Electrical Receptacles

All flush mount fixtures furnished complete with device, stainless steel face plate and electrical box. Pedestals are furnished with device, wrinkled black epoxy powder coated pedestal box, stainless steel face plate, tank nipple and locknut. All devices and boxes are UL listed. (Ship loose)









Type 2 and Type 3 pedestal boxes furnished with 3/4" IPS tank nipple, Type 7 and Type 8 pedestal boxes furnished with 1" IPS tank nipple.		Flush mount single gang box	Single Wide Pedestal Type 2 – Single Faced Type 3 – Double Faced			Flush mount double gang box	Double Wide Pedestal Type 7 – Single Faced Type 8 – Double Faced	
Electrical Device	Color	Type 1	Type 2	Type 3	Type 6	Type 7	Type 8	
<b>Specification Grade</b> 120 Volt AC, 20 amp, 3-wire polarized duplex receptacle NEMA 5-20 	Black	<b>S581-1K</b>	<b>S581-2K</b>	<b>S581-3K</b>	<b>S581-6K</b>	<b>S581-7K</b>	<b>S581-8K</b>	
	Ivory	<b>S581-1V</b>	<b>S581-2V</b>	<b>S581-3V</b>	<b>S581-6V</b>	<b>S581-7V</b>	<b>S581-8V</b>	
	White	<b>S581-1W</b>	<b>S581-2W</b>	<b>S581-3W</b>	<b>S581-6W</b>	<b>S581-7W</b>	<b>S581-8W</b>	
	Grey	<b>S581-1G</b>	<b>S581-2G</b>	<b>S581-3G</b>	<b>S581-6G</b>	<b>S581-7G</b>	<b>S581-8G</b>	
	Red	<b>S581-1R</b>	<b>S581-2R</b>	<b>S581-3R</b>	<b>S581-6R</b>	<b>S581-7R</b>	<b>S581-8R</b>	
<b>Specification Grade</b> 120 Volt AC, 20 amp, 3-wire polarized duplex receptacle NEMA 5-20 	Black	<b>S681-1K</b>	<b>S681-2K</b>	<b>S681-3K</b>	<b>S681-6K</b>	<b>S681-7K</b>	<b>S681-8K</b>	
	Ivory	<b>S681-1V</b>	<b>S681-2V</b>	<b>S681-3V</b>	<b>S681-6V</b>	<b>S681-7V</b>	<b>S681-8V</b>	
	White	<b>S681-1W</b>	<b>S681-2W</b>	<b>S681-3W</b>	<b>S681-6W</b>	<b>S681-7W</b>	<b>S681-8W</b>	
	Grey	<b>S681-1G</b>	<b>S681-2G</b>	<b>S681-3G</b>	<b>S681-6G</b>	<b>S681-7G</b>	<b>S681-8G</b>	
	Red	<b>S681-1R</b>	<b>S681-2R</b>	<b>S681-3R</b>	<b>S681-6R</b>	<b>S681-7R</b>	<b>S681-8R</b>	
<b>Hospital Grade</b> 120 Volt AC, 20 amp, 3-wire polarized duplex receptacle NEMA 5-20 	Black	<b>H581-1K</b>	<b>H581-2K</b>	<b>H581-3K</b>	<b>H581-6K</b>	<b>H581-7K</b>	<b>H581-8K</b>	
	Ivory	<b>H581-1V</b>	<b>H581-2V</b>	<b>H581-3V</b>	<b>H581-6V</b>	<b>H581-7V</b>	<b>H581-8V</b>	
	White	<b>H581-1W</b>	<b>H581-2W</b>	<b>H581-3W</b>	<b>H581-6W</b>	<b>H581-7W</b>	<b>H581-8W</b>	
	Grey	<b>H581-1G</b>	<b>H581-2G</b>	<b>H581-3G</b>	<b>H581-6G</b>	<b>H581-7G</b>	<b>H581-8G</b>	
	Red	<b>H581-1R</b>	<b>H581-2R</b>	<b>H581-3R</b>	<b>H581-6R</b>	<b>H581-7R</b>	<b>H581-8R</b>	
<b>Specification Grade</b> 240 Volt AC, 20 amp, 3-wire polarized duplex receptacle NEMA 6-20 	Black	<b>S582-1K</b>	<b>S582-2K</b>	<b>S582-3K</b>	<b>S582-6K</b>	<b>S582-7K</b>	<b>S582-8K</b>	
	Ivory	<b>S582-1V</b>	<b>S582-2V</b>	<b>S582-3V</b>	<b>S582-6V</b>	<b>S582-7V</b>	<b>S582-8V</b>	
	White	<b>S582-1W</b>	<b>S582-2W</b>	<b>S582-3W</b>	<b>S582-6W</b>	<b>S582-7W</b>	<b>S582-8W</b>	
	Grey	<b>S582-1G</b>	<b>S582-2G</b>	<b>S582-3G</b>	<b>S582-6G</b>	<b>S582-7G</b>	<b>S582-8G</b>	
	Red	<b>S582-1R</b>	<b>S582-2R</b>	<b>S582-3R</b>	<b>S582-6R</b>	<b>S582-7R</b>	<b>S582-8R</b>	
<b>Ground Fault Protected</b> 120 Volt AC, 20 amp, specification grade, duplex receptacle NEMA 5-20 	Black	<b>S656-1K</b>	<b>S656-2K</b>	<b>S656-3K</b>				
	Ivory	<b>S656-1V</b>	<b>S656-2V</b>	<b>S656-3V</b>				
	White	<b>S656-1W</b>	<b>S656-2W</b>	<b>S656-3W</b>				
	Grey	<b>S656-1G</b>	<b>S656-2G</b>	<b>S656-3G</b>				
	Red	<b>S656-1R</b>	<b>S656-2R</b>	<b>S656-3R</b>				
<b>Ground Fault Protected</b> 120 Volt AC, 20 amp, duplex receptacle and 3-wire polarized duplex receptacle NEMA 5-20 	Black				<b>S657-6K</b>	<b>S657-7K</b>	<b>S657-8K</b>	
	Ivory				<b>S657-6V</b>	<b>S657-7V</b>	<b>S657-8V</b>	
	White				<b>S657-6W</b>	<b>S657-7W</b>	<b>S657-8W</b>	
	Grey				<b>S657-6G</b>	<b>S657-7G</b>	<b>S657-8G</b>	
	Red				<b>S657-6R</b>	<b>S657-7R</b>	<b>S657-8R</b>	
<b>Ground Fault Protected</b> 120 Volt AC, 20 amp, duplex receptacle and 3-wire polarized duplex receptacle NEMA 5-20 	Black				<b>S658-6K</b>	<b>S658-7K</b>	<b>S658-8K</b>	
	Ivory				<b>S658-6V</b>	<b>S658-7V</b>	<b>S658-8V</b>	
	White				<b>S658-6W</b>	<b>S658-7W</b>	<b>S658-8W</b>	
	Grey				<b>S658-6G</b>	<b>S658-7G</b>	<b>S658-8G</b>	
	Red				<b>S658-6R</b>	<b>S658-7R</b>	<b>S658-8R</b>	

# Laboratory Electrical Fixtures

## Electrical Switches and Data Ports

All flush mount fixtures furnished complete with device, stainless steel face plate and electrical box. All devices and boxes are UL listed. (Ship loose)



Electrical Device		Color	Flush mount single gang box		Flush mount double gang box	
			Cat-6	Cat-6a	Cat-6	Cat-6a
<b>Specification Grade</b> Single pole, 120/240 Volt AC, 20 amp toggle switch. 	Black Ivory White Grey Red	<b>S655-1K</b> <b>S655-1V</b> <b>S655-1W</b> <b>S655-1G</b> <b>S655-1R</b>				
<b>RJ-45 Data Port Duplex</b> Cat-6 Data Port 	Black Ivory White Grey Red		<b>0720-1K-6</b> <b>0720-1V-6</b> <b>0720-1W-6</b> <b>0720-1G-6</b> <b>0720-1R-6</b>	<b>0720-1K-6a</b> <b>0720-1V-6a</b> <b>0720-1W-6a</b> <b>0720-1G-6a</b> <b>0720-1R-6a</b>		
<b>RJ-45 Data Port Duplex</b> Cat-6 Data Port with 120 Volt AC, 20 amp, specification grade, duplex receptacle. NEMA 5-20 	Black Ivory White Grey Red				<b>S721-6K-6</b> <b>S721-6V-6</b> <b>S721-6W-6</b> <b>S721-6G-6</b> <b>S721-6R-6</b>	<b>S721-6K-6a</b> <b>S721-6V-6a</b> <b>S721-6W-6a</b> <b>S721-6G-6a</b> <b>S721-6R-6a</b>
<b>RJ-45 Data Port Duplex</b> Cat-6 Data Port with 120 Volt AC, 20 amp, specification grade, GFI duplex receptacle. NEMA 5-20 	Black Ivory White Grey Red				<b>S751-6K-6</b> <b>S751-6V-6</b> <b>S751-6W-6</b> <b>S751-6G-6</b> <b>S751-6R-6</b>	<b>S751-6K-6a</b> <b>S751-6V-6a</b> <b>S751-6W-6a</b> <b>S751-6G-6a</b> <b>S751-6R-6a</b>



# Laboratory Seating

**Specifications:**

The seating shown on this page features chrome-plated metal parts. Upholstery is No. 684 Black Vinyl. Soft

Touch pneumatic mechanism allows easy adjustment to the seating height. Although not shown in the photographs

below, models F-4940-00 thru F-4949-00 come with 2" soft wheel ball bearing casters for hard floors.



- Upholstered seat/backrest
- 5-legged base with 18" fixed footring
- Tubular base

Model No.	Seat Height Adjustment	Footring Height
<b>F-4940-00</b>	16 <sup>3</sup> / <sub>4</sub> "-22"	8 <sup>1</sup> / <sub>2</sub> "
<b>F-4941-00</b>	22 <sup>1</sup> / <sub>4</sub> "-27 <sup>1</sup> / <sub>2</sub> "	8 <sup>1</sup> / <sub>2</sub> "
<b>F-4942-00</b>	26 <sup>1</sup> / <sub>4</sub> "-31 <sup>1</sup> / <sub>2</sub> "	12 <sup>1</sup> / <sub>2</sub> "



- Upholstered seat/backrest
- 5-legged base with 18" fixed foot ring
- Tubular base
- Backrest with lumbar support

Model No.	Seat Height Adjustment	Footring Height
<b>F-4943-00</b>	16 <sup>3</sup> / <sub>4</sub> "-22"	8 <sup>1</sup> / <sub>2</sub> "
<b>F-4944-00</b>	22 <sup>1</sup> / <sub>4</sub> "-27 <sup>1</sup> / <sub>2</sub> "	8 <sup>1</sup> / <sub>2</sub> "
<b>F-4945-00</b>	26 <sup>1</sup> / <sub>4</sub> "-31 <sup>1</sup> / <sub>2</sub> "	12 <sup>1</sup> / <sub>2</sub> "

# Laboratory Seating

**Specifications:**

The seating shown on this page features chrome-plated metal parts. Upholstery is No. 684 Black Vinyl. Soft

Touch pneumatic mechanism allows easy adjustment to the seating height. Although not shown in the photographs

below, models F-4940-00 thru F-4949-00 come with 2" soft wheel ball bearing casters for hard floors.



- 13" round chrome-plated seat
- 5-legged base with 18" fixed foot ring
- Tubular base

Model No.	Seat Height Adjustment	Footring Height
<b>F-4946-00</b>	18 <sup>1</sup> / <sub>2</sub> "-22 <sup>3</sup> / <sub>4</sub> "	8 <sup>1</sup> / <sub>2</sub> "
<b>F-4947-00</b>	22 <sup>1</sup> / <sub>2</sub> "-27 <sup>3</sup> / <sub>4</sub> "	12 <sup>1</sup> / <sub>2</sub> "



- 15" round upholstered seat
- 5-legged base with 18" fixed foot ring
- Tubular base

Model No.	Seat Height Adjustment	Footring Height
<b>F-4948-00</b>	21 <sup>1</sup> / <sub>2</sub> "-24 <sup>3</sup> / <sub>4</sub> "	8 <sup>1</sup> / <sub>2</sub> "
<b>F-4949-00</b>	25 <sup>1</sup> / <sub>2</sub> "-20 <sup>3</sup> / <sub>4</sub> "	12 <sup>1</sup> / <sub>2</sub> "




- 13" round wood seat
- Tubular 4-legged base with 16" fixed foot ring
- Black tubular steel base

Model No.	Seat Height
<b>F-5350-00</b>	18"
<b>F-5351-00</b>	24"
<b>F-5352-00</b>	30"

# Steel Cabinet Technical Information

## Inset Style Wall Bench End Views




**24" Deep Wall Bench**

**30" Deep Wall Bench**

**31" Deep Wall Bench**

## Inset Style Center Bench End Views




**Center Bench End/Section View**


**Center Bench End Sink Cabinet**

# Steel Cabinet Technical Information


## Inset Style Full Height Storage Cabinet End Views



12<sup>3</sup>/<sub>4</sub>" Deep Full Height Cabinet




16" Deep Full Height Cabinet



22" Deep Full Height Cabinet


## Inset Style Center Bench Plan View




Plan View  
Center Bench End Sink Cabinet

# Steel Cabinet Technical Information


## Full Overlay Style Wall Bench End Views



**24" Deep Wall Bench**




**30" Deep Wall Bench**




**31" Deep Wall Bench**

## Full Overlay Style Center Bench End Views




**Center Bench End/Section View**



**Center Bench End Sink Cabinet**

# Steel Cabinet Technical Information

## Full Overlay Style Full Height Storage Cabinet End Views




12<sup>3</sup>/<sub>4</sub>" Deep Full Height Cabinet

16" Deep Full Height Cabinet

22" Deep Full Height Cabinet

## Full Overlay Style Center Bench Plan View



Plan View  
Center Bench End Sink Cabinet

# Lock Schedule Sample and Examples

## Lock Schedule Examples

### Keyed Different & Master-keyed

An example of a simple lock schedule with disk tumbler locks. All locks on the project are keyed different with a single master key that can open all locks. Either disk tumbler or pin tumbler locks could be used for this arrangement.

### LOCK SCHEDULE

Sheet \_\_\_\_ of \_\_\_\_

Project Name: SAMPLE PROJECT				Project/Quote Number		Date:		for Order Dept. use	
Lock Type	<input checked="" type="checkbox"/> Disk Tumbler <input type="checkbox"/> Pin Tumbler	Number of Keys per Lock	2	Number of Master Keys per change	4	Number of Grand-Master Keys per change			
Room Number	Item Number	Cabinet (all cabinets with locks must be listed)	Number of Locks	Keyed Different	Keyed Alike	Key Group	Masterkey Group	Grand-Masterkey Group	
100	1	D30C3518-0101-J1	4	<input checked="" type="checkbox"/>	<input type="checkbox"/>	A			
100	1	D36C3536-0101-J1	6	<input checked="" type="checkbox"/>	<input type="checkbox"/>	A			
100	2	E40C3524-0101-J1	2	<input checked="" type="checkbox"/>	<input type="checkbox"/>	A			
105	6	D30C3524-0101-J1	4	<input checked="" type="checkbox"/>	<input type="checkbox"/>	A			
105	6	E41C3536-0101-J1	3	<input checked="" type="checkbox"/>	<input type="checkbox"/>	A			
105	6	W20C2936-0101-J1	1	<input checked="" type="checkbox"/>	<input type="checkbox"/>	A			

### Keyed Alike in Sets

An example of another simple lock schedule with disk tumbler locks. On this project all the locks in each room are keyed alike, but each room is different from the other. A master key could be used in this type of keying arrangement, but in this example it is not. Either disk tumbler or pin tumbler locks could be used for this arrangement.

### LOCK SCHEDULE

Sheet \_\_\_\_ of \_\_\_\_

Project Name: SAMPLE PROJECT				Project/Quote Number:		Date:		for Order Dept. use	
Lock Type	<input checked="" type="checkbox"/> Disk Tumbler <input type="checkbox"/> Pin Tumbler	Number of Keys per Lock	2	Number of Master Keys per change		Number of Grand-Master Keys per change			
Room Number	Item Number	Cabinet (all cabinets with locks must be listed)	Number of Locks	Keyed Different	Keyed Alike	Key Group	Masterkey Group	Grand-Masterkey Group	
115	5	D30C3518-0101-K1	4	<input type="checkbox"/>	<input checked="" type="checkbox"/>	A			
115	5	E56C3536-0101-K1	4	<input type="checkbox"/>	<input checked="" type="checkbox"/>	A			
115	5	E41C3536-0101-K1	3	<input type="checkbox"/>	<input checked="" type="checkbox"/>	A			
126	16	E40C3524-0101-K1	2	<input type="checkbox"/>	<input checked="" type="checkbox"/>	B			
126	16	E42C3536-0101-K1	5	<input type="checkbox"/>	<input checked="" type="checkbox"/>	B			
126	17	C03C3536-0101-K1	1	<input type="checkbox"/>	<input checked="" type="checkbox"/>	B			
126	17	D30C3524-0101-K1	4	<input type="checkbox"/>	<input checked="" type="checkbox"/>	B			

### Keyed Different – Master-keyed & Grand-Masterkeyed

An example of a more complex lock schedule with pin tumbler locks. All locks are keyed different, with a separate master key for each room, and a grand-master key for the entire project. This keying arrangement could not be used with disk tumbler locks.

### LOCK SCHEDULE

Sheet \_\_\_\_ of \_\_\_\_

Project Name: SAMPLE PROJECT				Project/Quote Number:		Date:		for Order Dept. use	
Lock Type	<input type="checkbox"/> Disk Tumbler <input checked="" type="checkbox"/> Pin Tumbler	Number of Keys per Lock	2	Number of Master Keys per change	5	Number of Grand-Master Keys per change	5		
Room Number	Item Number	Cabinet (all cabinets with locks must be listed)	Number of Locks	Keyed Different	Keyed Alike	Key Group	Masterkey Group	Grand-Masterkey Group	
102	1	E43C3536-0109-J2	6	<input checked="" type="checkbox"/>	<input type="checkbox"/>	A	A		
102	1	E40C3518-0109-J2	2	<input checked="" type="checkbox"/>	<input type="checkbox"/>	A	A		
102	1	D30C3524-0109-J2	4	<input checked="" type="checkbox"/>	<input type="checkbox"/>	A	A		
105	3	D40C3524-0109-J2	2	<input checked="" type="checkbox"/>	<input type="checkbox"/>	B	A		
105	3	D30C3524-0109-J2	4	<input checked="" type="checkbox"/>	<input type="checkbox"/>	B	A		
1113	9	D32C3536-0109-J2	8	<input checked="" type="checkbox"/>	<input type="checkbox"/>	C	A		
1113	9	C03C3536-0109-J2	1	<input checked="" type="checkbox"/>	<input type="checkbox"/>	C	A		
1113	10	D30C3524-0109-J2	4	<input checked="" type="checkbox"/>	<input type="checkbox"/>	C	A		





# Alphabetical Index

<b>A</b>	Acid Storage Cabinet–ADA Height . . . . .	18 & 21		
	Acid Storage Cabinet–Standing Height. . . . .	27 & 32		
	Acid Storage Fume Hood Cabinet. . . . .	21 & 32		
	Acid Storage Cabinet - Full Height . . . . .	57		
	<b>ADA Height Cabinets</b> . . . . .	14-21		
	ADA Compliant Sink . . . . .	94-95		
	Adjustable Slip Scribes . . . . .	65		
	Aerator . . . . .	87		
	Air Valves . . . . .	83-84		
	Applique Panels. . . . .	58-63		
	Apron Rails . . . . .	43		
	Aprons . . . . .	38-43		
	Aspirator . . . . .	87		
<b>B</b>				
	Ball Valves . . . . .	83		
	Blind Corner Apron . . . . .	42		
	Blind Corner Cabinet . . . . .	18 & 29		
	Balance Table . . . . .	79		
	Base Molding. . . . .	81		
<b>C</b>				
	Cabinet Dimensions . . . . .	5 & 100		
	Cabinet Styles . . . . .	3		
	Cart - Stainless Steel . . . . .	79		
	Caster Options. . . . .	6		
	Catalog Number Explanation. . . . .	4		
	Center Bench Shelving . . . . .	72-73		
	Ceiling Enclosure. . . . .	67		
	Chalkboard . . . . .	79		
	Cold Water Faucet. . . . .	85		
	Combination Faucet. . . . .	84		
	Construction Details. . . . .	2		
	Corner Box Apron . . . . .	42		
	Corner Cabinets–ADA Height . . . . .	18		
	Corner Cabinets–Sitting Height. . . . .	13		
	Corner Cabinets–Standing Height. . . . .	29		
	Corner Cabinets–Wall Mounted. . . . .	52		
	Corner Fillers . . . . .	59, 61 & 63		
	Corner Guard. . . . .	81		
	Crumb Cup Strainer. . . . .	91		
	Cupboard/Drawer Cabinets–ADA Height . . . . .	16-17		
	Cupboard/Drawer Cabinets–Sitting Height. . . . .	12-13		
	Cupboard/Drawer Cabinets–Standing Height. . . . .	24-25		
	Cupboard Cabinets–ADA Height . . . . .	14		
	Cupboard Cabinets–Sitting Height . . . . .	11		
	Cupboard Cabinets–Standing Height . . . . .	22		
	Cupsink . . . . .	91		
<b>D</b>				
	Deionized Water Faucet . . . . .	88		
	Diagonal Corner Cabinets–ADA Height . . . . .	18		
	Diagonal Corner Cabinets–Sitting Height . . . . .	13		
	Diagonal Corner Cabinets–Standing Height . . . . .	29		
	Distilled Water Faucet . . . . .	88		
	Double-Sided Wall Cabinet . . . . .	49		
	Drain Outlet . . . . .	91		
	Drain Trough . . . . .	92		
	Drawer Cabinets–ADA Height. . . . .	15		
	Drawer Cabinets–Sitting Height. . . . .	11-12		
	Drawer Cabinets–Standing Height . . . . .	22-23		
	Drench Shower . . . . .	90		
<b>E</b>				
	Electrical Fixtures. . . . .	96-97		
	Emergency Eye Wash/Shower. . . . .	90		
	End Leg Assembly. . . . .	44		
	End View Dimensions . . . . .	100-103		
	Epoxy Resin Cupsink. . . . .	91		
	Epoxy Resin Drain Trough. . . . .	92		
	Epoxy Resin Top . . . . .	70-71		
	Epoxy Resin Tub Sink . . . . .	92-93		
	Eye Wash. . . . .	90		
<b>F</b>				
	Face Wash . . . . .	90		
	Faucet Accessories . . . . .	87		
	Filler Panel . . . . .	58-65		
	Filter Pump. . . . .	87		
	Finished Back . . . . .	64		
	Fixtures–Electrical . . . . .	96-97		
	Foot Valve . . . . .	87		
	Freezer. . . . .	79		
	<b>Full Height Storage Cabinets</b> . . . . .	53-57		
	Full Overlay Applique Panels. . . . .	59, 61 & 63		
	Fume Hood Acid Storage Cabinet–ADA Ht. . . . .	20		
	Fume Hood Acid Storage Cabinet–Standing Ht. . . . .	31		
	Fume Hood Base Cabinet–ADA Height. . . . .	20		
	Fume Hood Base Cabinet–Standing Height . . . . .	30		
	Fume Hoods . . . . .	68-69		
<b>G</b>				
	Gas Valves. . . . .	83-89		
	Gooseneck. . . . .	84-87		
	Gooseneck Coupling. . . . .	87		
<b>H</b>				
	Hoods . . . . .	68-69		
	Hose Spray . . . . .	90		
	Hot & Cold Water Mixing Faucet . . . . .	86-87		
<b>I</b>				
	Intermediate Leg Assembly. . . . .	44		
<b>K</b>				
	Kemresin Work Tops . . . . .	70-71		
	Kemresin Trough . . . . .	92		
	Kemresin Tub Sink. . . . .	92-94		
	Kemstruts. . . . .	80		
	Key Case . . . . .	79		
	Keying Options . . . . .	9-10		
	Knee Space Panel. . . . .	45		
<b>L</b>				
	Label Holder . . . . .	81		
	Laboratory Cart . . . . .	79		
	Laboratory Service Fittings . . . . .	83-90		
	LED Lighting . . . . .	80		
	Leg Assemblies . . . . .	44		
	Lock Keying. . . . .	9-10		
	Lock Panels . . . . .	8		
	Lock Schedule Examples . . . . .	104		
	Lock Schedule Sample . . . . .	105		
	Locks . . . . .	8		

## Alphabetical Index (continued)

<b>M</b>	Markerboard . . . . .	79	Stainless Steel Sink . . . . .	95
	Mobile Cabinets. . . . .	34-37	Stainless Steel Sink Outlet . . . . .	91
	Mixing Faucet . . . . .	86-87	<b>Standing Height Base Cabinets</b> . . . . .	22-33
	Multiple Service Fittings. . . . .	84	Steam Valves . . . . .	84
<b>N</b>			Stools . . . . .	98-99
	Needle Point Cock Valves . . . . .	83-84	Styles . . . . .	3
<b>O</b>			Support Angle . . . . .	80
	Open Cabinets–Sitting Height. . . . .	11	Support Strut . . . . .	72 & 80
	Open Cabinets–ADA Height . . . . .	14	Swinging Glazed Door Full Height Cabinet. . . . .	54
	Open Cabinets–Standing Height. . . . .	22	Swinging Glazed Door Wall Cabinet . . . . .	50
	Open Cabinets–Full Height . . . . .	53	Swinging Panel Door Full Height Cabinet . . . . .	54
	Open Cabinets–Wall Mounted. . . . .	46	Swinging Panel Door Wall Cabinet . . . . .	51
<b>P</b>			Switch . . . . .	97
	Part Number Explanation. . . . .	4	<b>T</b>	
	Pass-Thru Wall Cabinet . . . . .	49	Tackboard . . . . .	79
	Pedestal Legs . . . . .	45	Task Lighting . . . . .	80
	Pegboard . . . . .	78	Tops . . . . .	68-69
	Pipe Enclosure. . . . .	67	Trap . . . . .	91
	Plastic Laminate Work Tops . . . . .	70-71	Trough . . . . .	92
	Plastic Laminate Ledges & Shelves . . . . .	76	Tub Sink. . . . .	92-93
	Plumbing Information. . . . .	82	<b>U</b>	
	Portable Pump Faucet. . . . .	85	Ultra-Pure Water Faucet . . . . .	88
<b>R</b>			Under Counter Refrigerator . . . . .	79
	Reagent Rack Assemblies. . . . .	76	Upright Rod Assembly. . . . .	81
	Reagent Rack Ledges. . . . .	76	<b>V</b>	
	Rear Scribes & Filler . . . . .	59 & 63	Vacuum Breaker . . . . .	87
	Receptacle - Electric . . . . .	96-97	Vacuum Pump Storage Fume Hood Cabinet 20&33	
	Refrigerator . . . . .	79	Vacuum Pump Cart. . . . .	33
	Rods & Accessories . . . . .	81	Vacuum Valves . . . . .	83-84
<b>S</b>			Vandal-Resistant Fittings . . . . .	89
	Safety Shower . . . . .	90	Ventilated Enclosures . . . . .	68-69
	Scribe . . . . .	58-67	<b>W</b>	
	Seating. . . . .	98-99	<b>Wall Cabinets</b> . . . . .	46-52
	Service Fittings . . . . .	82-90	Wall Shelving . . . . .	75-77
	Shelf Brackets . . . . .	77	Wardrobe Storage Cabinet . . . . .	56
	Shelf Lips . . . . .	6 & 73	Waste Storage Cabinet . . . . .	29
	Shelving Assemblies . . . . .	72-75	Water Faucet . . . . .	84-89
	Sink . . . . .	91-95	Water Mixing Faucet . . . . .	86-89
	Sink Cabinets–ADA Height . . . . .	19	Water Valves . . . . .	84
	Sink Cabinets–Sitting Height . . . . .	13	Water/Gas Fitting . . . . .	84
	Sink Cabinets–Standing Height. . . . .	26	Work Tops . . . . .	70-71
	Sink Outlet . . . . .	91	Work Top Support Frame. . . . .	43
	Sink Overflow. . . . .	91	Wrist Blade Handle . . . . .	87
	Sink Trap . . . . .	91		
	<b>Sitting Height Base Cabinets</b> . . . . .	11-13		
	Sloping Tops . . . . .	66		
	Sliding Glazed Door Full Height Cabinet. . . . .	53		
	Sliding Glazed Door Wall Cabinet . . . . .	48		
	Sliding Panel Door Full Height Cabinet . . . . .	53		
	Sliding Panel Door Wall Cabinet . . . . .	48		
	Sliding Plate Glass Door Wall Cabinet. . . . .	47 & 49		
	Solvent Storage Cabinet . . . . .	19-20, 28 & 32		
	Spreader . . . . .	45		
	Spring Clips . . . . .	65		
	Stainless Steel Cart . . . . .	79		
	Stainless Steel Cupsink . . . . .	91		

# Numerical Index

<b>Part Number</b>	<b>Page</b>	<b>Part Number</b>	<b>Page</b>
019623	31	C03C292215L_____	11
0464-00	91	C03C322215L_____	14
0469-00	91	C03C352212L_____	22
0482-BP	91	CCEC121320-00_0	67
0490-BE	91	CSTC081327-00_0	66
0491-BP	91	D00C292215-_____	11
0492-BP	91	D00C322215-_____	15
0493-BE	91	D00C352215-_____	22
0494BP-04	91	D02C322230-_____	15
0548-00	81	D02C352230-_____	22
0720-1K-6	97	D08C352215-_____	22
0950-0F	91	D26C292218-_____	11
0995-00	91	D30C292215-_____	11
1000-00-__	93	D30C322215-_____	15
1000-DI-__	94	D30C352212-_____	23
1000-ADA-DI-__	94	D31C292230-_____	11
4740-00	81	D31C322230-_____	15
A00C022215-00_0	43	D31C352230-_____	23
A00C042215-__0	38	D32C292230-_____	11
A00C042815-__0	39	D32C322230-_____	15
A04C042248-1__	42	D32C352230-_____	23
A06C020215-00_0	43	D36C322230-_____	15
A06C040315-1__0	43	D36C352230-_____	23
A07C292202-00_0	44	D50C322215-_____	15
A09C020230-00_0	45	D50C352215-_____	23
A15C290202-00_0	45	D51C322230-_____	15
A10C250015-00_0	45	D51C352230-_____	23
A11C292202-00_0	45	D52C322230-_____	15
A27C292202-00_0	44	D52C352230-_____	23
A51C042220-1__	40	D81C292218-_____	12
A51C042820-1__	41	D81C352218-_____	23
A56C040236-00_0	43	D83C322215-_____	15
A57C352202-00_0	44	D84C292218-_____	12
A65C350202-00_0	45	D84C322215-_____	15
A77C352202-00_0	44	D84C352218-_____	23
ABSC290012-00_0	65	D85C292218-_____	12
AWSC241312-00_0	65	E05C292230L_____	12
AFSC841312-00_0	65	E05C322230L_____	16
B00C292215-00_0	11	E05C352230L_____	24
B00C322215-00_0	14	E40C292215L_____	12
B00C352215-00_0	22	E40C322215L_____	16
BCAC290003-11_0	59	E40C352212L_____	24
BCFC290001-00_0	59	E41C292230-_____	12
BEAC290003-11_0	58	E41C322230-_____	16
BESC290001-00_0	58	E41C352230-_____	24
BFBC290112-00_0	64	E42C322230L_____	16
BIAC290006-11_0	58	E42C352230L_____	24
BIFC290006-11_0	58	E43C322230L_____	16
BPC290001-00_0	59	E43C352230L_____	24
BRSC290001-00_0	59	E45C292230L_____	13
BPAC290023-16_0	59	E45C322230L_____	16
		E45C352230L_____	24
		E56C322236-_____	16
		E56C352236-_____	24
		E60C352212L_____	25
		E61C352230-_____	25
		E62C322230L_____	17
		E62C352230L_____	25

# Research Collection

## Numerical Index (continued)

Part Number	Page	Part Number	Page
E63C322230L_____ - E63C322248-_____	17	G35C322024L_____ - G35C322248-_____	20
E63C352230L_____ - E63C352248-_____	25	G35C352024L_____ - G35C352248-_____	33
E65C322230L_____ - E65C322248-_____	17	G63C322018L_100 - G63C322248-_____	19
E65C352230L_____ - E65C352248-_____	25	G63C352024L_100 - G63C352248-_____	28
ECEC120013-00_0 - ECEC240022-00_0	67	G68C322024L_100 - G68C322248-_____	20
F-0260-00 - _____	81	G68C352024L_100 - G68C352248-_____	32
F-0862-04 - _____	81	G80C322018L_____ - G80C322248-_____	21
F-1368-00 - _____	81	G80C352018L_____ - G80C352260-_____	31
F-2939-00 - F-2940-00	78	H581-1K - H581-8R	96
F-4940-00 - F-4945-00	98	K-0005-0A - K-0020-1A	80
F-4946-00 - F-5352-00	99	K7-302435-BK - K7-302435-SL	79
F-4995-PK - _____	65	K12-2905-0A - K12-3520-0A	80
F-5332-00 - F-5332-01	79	K17-3612-00 - K17-3614-00	73
F-5333-00 - F-5348-00	79	KF-0170-04 - _____	77
F-5349-00 - _____	79	K-MSS-12 - K-MSS-60	13, 19 & 26
F-5353-06 - F-5362-12	95	L-105 - L-115	79
F-8440-15 - F-8441-48	27	L-2034-06 - L-2034-96	80
F10M0242-00 - F10M0279-00	74	L01-020109-00 - L01-020142-00	47 & 81
F11M0196-00 - F11M0372-00	74	M-0051-0A - M0052-0A	80
F12M000224-00 - F13S000284-00	72	M-0030-0A - M-0044-AD	77
F17-4824-00 - F17-4872-00	72	M-0060-0A - M-0062-AD	77
F21M0624-00 - F21W2448-0_	75	M10C262215L_____ - M10C262224-_____	34
F21M0924-60 - F21W1848-6_	73	M10C332215L_____ - M10C332224-_____	35
F63M0224-00 - F63M0272-00	73	M20C262215L_____ - M20C262224-_____	34
F73M0015-00 - F73S0024-00	37	M20C332215L_____ - M20C332224-_____	35
F94M2215-00 - F94M2348-00	37	M25C262215L_____ - M25C262224-_____	34
F94P2215-00 - F94T2448-00	36	M25C332215L_____ - M25C332224-_____	35
FCAC840003-11_0 - FCAC840012-16_0	63	M30C262215-_____ - M30C262224-_____	34
FCEC120018-00_0 - FCEC240096-00_0	67	M35C262215-_____ - M35C263224-_____	34
FCCFC841301-00_0 - FRSC840024-00_0	63	M40C332215-_____ - M40C332224-_____	35
FEAC840003-11_0 - FEAC840024-16_0	62	M45C262215-_____ - M45C263224-_____	34
FESC841301-11_0 - FESC842224-00_0	62	M46C332215-_____ - M46C332224-_____	35
FFBC840118-00_0 - FFBC840148-00_0	64	M56C332215-_____ - M56C332224-_____	35
FIAC840006-11_0 - FIAC840024-16_0	62	M62C262215-_____ - M62C262224-_____	34
FIFC841301-00_0 - FIFC842224-00_0	62	M66C262215-_____ - M66C262224-_____	34
FPAC840013-16_0 - FPAC840022-16_0	63	M68C332215-_____ - M68C332224-_____	35
FRSC560001-00_0 - FRSC840024-00_0	63	M71C332215-_____ - M71C332224-_____	35
G00C292215L_____ - G00C292248-_____	13	PL-4812-00 - PL-7214-00	76
G00C322212L_____ - G00C322248-_____	19	RRU-1210-0000 - RRU-1812-0030	76
G00C352212L_____ - G00C352260-_____	26	S00C841318-00-0 - S00C842248-00_0	53
G03C322236-_____ - G03C322248-_____	21	S12C841330-_____0 - S12C842248-_____0	53
G03C352230-_____ - G03C352248-_____	26	S25C841318L_____ - S25C842248-_____	54
G04C352252-_____ - G04C352260-_____	26	S35C841330-_____ - S35C842248-_____	56
G07C322236-_____ - G07C322248-_____	21	S40C842218L_____ - S40C842248-_____	56
G08C322012L_____ - G08C322248-_____	20	S41C841330-_____0 - S41C842248-_____0	53
G08C352012-_____ - G08C352260-_____	30	S50C841318L_____ - S50C842248-_____	54
G12C322248-_____ - G13C322248-_____	18	S51C841624L_____ - S51C842248-_____	55
G12C352248-_____ - G13C352248-_____	29	S60C841330-_____ - S60C842248-_____	56
G14C292214-_____ - _____	13	S70C841624L_____ - S70C842248-_____	57
G14C322214-_____ - G15C352214-_____	18		
G14C352214-_____ - G15C352224-_____	29		
G20C352215-_____ - G20C352224-_____	29		
G23C352218-_____ - _____	29		
G21C322018L_____ - G21C322248-_____	18		
G21C352015L_____ - G21C352248-_____	27		
G23C322218-_____ - _____	19		

# Numerical Index (continued)

<b>Part Number</b>	<b>Page</b>	<b>Part Number</b>	<b>Page</b>
S581-1K - S581-8R	96	W00C241315-00_0 - W00C481648-00_0	46
S582-1K - S582-8R	96	W10C241330-00_0 - W10C481648-00_0	47
S655-1K - S655-1R	97	W12C241330-__0 - W12C481648-__0	48
S656-1K - S656-3R	96	W13C241330-00_0 - W13C361648-00-0	49
S657-6K - S657-8R	96	W25C241315-____ - W25C481648-____	51
S658-6K - S658-8R	96	W26C241320-____ - W26C481615-____	52
S681-1K - S681-8R	96	W41C241330-__0 - W41C481648-__0	48
S721-6K-6 - S721-6R-6a	97	W50C241315-____ - W50C481648-____	50
S751-6K-6 - S751-6R-6a	97	W51C241320-____ - W51C481615-____	52
SCAC010103-00_0 - SCAC010196-00_0	65	WCAC240003-11_0 - WCAC480012-16_0	61
SCFC081301-00_0 - SCFC142206-00_0	66	WCFC241301-00_0 - WCFC481612-00_0	61
SLTC081315-00_0 - SLTC142296-00_0	66	WEAC240003-11_0 - WEAC470024-16_0	60
SM-01-00-__ - SM-56-00-__	93	WESC241301-00_0 - WESC481624-00_0	60
STFC081301-00_0 - STFC142212-00_0	66	WFBC240115-00_0 - WFBC480148-00_0	64
TR-01-00-__ - TR-05-00-__	92	WIAC240006-11_0 - WIAC480048-16_0	60
U01C600512-00_0 - U02C961048-00_0	67	WIFC241301-00_0 - WIFC481648-00_0	60
VPCC061414-00_0 - VPCC061644-00_0	33	WPAC240013-16_0 - WPAC480016-16_0	61
W-0242-00 - W-0243-00	83	X-020004 - X-020017	78
W-0259-00	84	X-030001-00 - X-030003-00	80
W-0260-00 - W-0272-00	83		
W-0273-02 - W-0274-VB	84		
W-0286-02 - W-0287-00	83		
W-0300-00 - W-0306-00	84		
W-0333-00	86		
W-0337-00 - W-0337-6V	85		
W-0340-00 - W-0340-0V	86		
W-0350-00 - W-0350-0V	85		
W-0351-00 - W-0351-0V	85		
W-0354-00 - W-0355-00	86		
W-0358-00	87		
W-0361-00 - W-361-0V	87		
W-0367-00 - W-0368-00	88		
W-0374-01	90		
W-0400-00 - W-0401-00	83		
W-0402-00 - W-0404-00	84		
W-0405-00	84		
W-0407-00 - W-0409-00	88		
W-0410-VR - W-0417-VR	89		
W-0510-00	87		
W-0520-00	87		
W-0532-00	87		
W-0535-00	87		
W-0825-00	87		
W-0833-00	87		
W-0923-00	90		
W-0926-00	90		
W-0927-00	90		
W-0928-00	90		
W-0930-00	90		
W-0931-00	90		
W-0932-00	90		
W-0933-00	90		
W-0934-00	90		
W-0935-00	90		
W-0936-00	90		

# Research Collection

# Typical Installations



## Trademarks

Kewaunee, Alpha, Basic Bench, EarthLine, Enterprise, Evolution, Interceptor, Kemresin, Kewaunee Matrix, Research Collection, Signature, Sturdilite, and Supreme Air are registered trademarks of Kewaunee Scientific Corporation.

## Patents:

Products in this catalog are protected by one or more of the following patents:

4,856,420	5,241,788	5,302,030	5,407,389
5,582,472	5,728,626	D388,986	DES. 346,437

Additional patents are pending.

Part numbers, illustrations, descriptions, specifications and other information contained in this catalog are based on the most current design and production information available

at the time of publication. Kewaunee Scientific Corporation reserves the right to make changes at any time without notice.

# Typical Installations



We are pleased to have been selected to supply laboratory furniture for the facilities of some of the best companies, organizations and institutions in the world. We believe the repeat opportunities, over time, to participate in their construction projects is good evidence of our product quality and experience in successfully completing their laboratory projects.

Whether your laboratory projects are large or small, new construction or remodeling, Kewaunee has the qualified people and quality products to transform your visions into successful reality.

From pre-planning to final completion, experienced Kewaunee representatives are available to help you through the intricacies of laboratory design and layout, product selection, ordering, scheduling and installation.



Kewaunee Scientific Corporation is dedicated to manufacturing high quality products for the laboratory marketplace. We offer **Total Laboratory Solutions, Innovative Designs, and Technological Expertise** with a worldwide distribution network to provide laboratories that are truly World Class.

**Laboratory Products from Kewaunee®:**

ALPHA® – Adaptable Casework System

ALPHA Columns – Expandable Workstations

ALPHA Overhead Carrier Systems

BasikBench® – Welded C-Leg Bench

DetectAir – Filtered/Ductless Fume Hoods

EarthLine® – Environmentally Friendly Casework

Enterprise® – Moveable Workstations

Evolution® – Column Based Workstations

Interceptor® – Biological Safety Cabinet

Interceptor® – Laminar Flow Cabinet

Kemresin® – Epoxy Resin Counters & Sinks

Kewaunee Matrix® – Educational Stations

Research Collection® – Steel Casework

Signature® Series – Wood Casework

Spektrum – Laminate Architectural Millwork

Sturdilite® – Flexible Pedestal System

Supreme Air® – Fume Hoods

For our complete product offering, specifications, catalogs, brochures, case studies, and dealer list, go to [www.kewaunee.com](http://www.kewaunee.com), or contact your local sales representative directly.



Kewaunee, Alpha, Basic Bench, EarthLine, Enterprise, Evolution, Interceptor, Kemresin, Kewaunee Matrix, Research Collection, Signature, Sturdilite, and Supreme Air are registered trademarks of Kewaunee Scientific Corporation



P.O. Box 1842 • Statesville, NC 28687-1842  
Phone: (704) 873-7202 • Fax: (704) 873-5160  
E-mail: [kscmarketing@kewaunee.com](mailto:kscmarketing@kewaunee.com)  
Web Site: [www.kewaunee.com](http://www.kewaunee.com)

**Kewaunee Labway India Pvt. Ltd.**

No. CA-9A, 2nd Floor, Jigani Link Road  
Jigani, Bangalore 562 106, India  
Phone: 9180-27826725 • Fax: 9180-27826724  
Web Site: [www.kewaunee.in](http://www.kewaunee.in)

**Kewaunee Labway Asia Pte. Ltd.**

194 Pandan Loop, #06-22  
Pantech Business Hub, Singapore 128 383  
Phone: 65-6773-0288 • Fax: 65-6773-2322

**Kewaunee Scientific (Suzhou) Co Ltd**

333 Xingpu Road, Modern Industrial Square  
Building 4, Unit #2-A • Suzhou Industrial Park • China  
Phone: (0512) 6262 5090 • Fax: (0512) 6262 5089  
Web Site: [china.kewaunee.com](http://china.kewaunee.com)

# Research Collection

Kewaunee Steel Cabinet Construction	2
Kewaunee Steel Cabinet Styles	3
Steel Cabinet Information	4
Steel Cabinet Options	6
Lock Keying Options	9
Steel Sitting Height Base Cabinets	11
Steel ADA Height Base Cabinets	14
Steel Standing Height Base Cabinets	22
Fume Hood Base Cabinets	30
Steel Sitting Height Mobile Cabinets	34
Steel Standing Height Mobile Cabinets	35
Mobile Cabinet Tops & Accessories	36
Steel Apron Assemblies	38
22" Deep Apron	38
22" Deep Apron	38
Steel Box Apron	38
22" Deep Apron	38
Applique Panel	38
Steel Box Apron	38
28" Deep Apron	39
28" Deep Apron	39
Steel Box Apron	39
28" Deep Apron	39
Applique Panel	39
Steel Box Apron	39
Steel Leg Assemblies & Pedestals	44
Steel Apron Rails & Panels	45
Steel Wall Mounted Storage Cabinets	46
Steel Full Height Storage Cabinets	53
Base Cabinet Scribes & Fillers	58
Wall Cabinet Scribes & Fillers	60
Full Height Cabinet Scribes & Fillers	62
Steel Cabinet Finished Backs	64
Steel Adjustable Slip Scribes	65
Steel Sloping Tops	66
Steel Ceiling & Pipe Enclosures	67
Fume Hoods & Ventilated Enclosures	68
Work Tops	70
Center Bench Shelving Assemblies	72
Wall Shelving Assemblies	74
Reagent Rack & Shelving Systems	76
Pegboards	78
Miscellaneous Accessories	79
Miscellaneous Accessories	80
K12-	80
-0A	80
General Plumbing Information	82
Laboratory Service Fittings	83
Cupsinks & Drain Accessories	91
Kemresin Tub Sinks & Drain Troughs	92
Kemresin Drop-in Sinks	94
Stainless Steel Drop-in Sinks	95
Laboratory Electrical Fixtures	96
Laboratory Seating	98
Steel Cabinet Technical Information	100
Lock Schedule Sample and Examples	104
Alphabetical Index	106
Numerical Index	108
Typical Installations	112

# Signature Series

## Wood Laboratory Furniture



KEWAUNEE®

encouraging new discovery... *Worldwide*

# Signature Series

## Wood Laboratory Furniture





## Contents

Construction & Styles .....	2
Catalog Number Description .....	4
Cabinet Options .....	6
Sitting Height Cabinets .....	14
ADA Height Cabinets .....	18
Standing Height Cabinets .....	28
Special Purpose Base Cabinets .....	42
Mobile Cabinets .....	44
Apron & Leg Assemblies.....	48
Wall Mounted Storage Cabinets .....	66
Full Height Storage Cabinets.....	74
Filler, Soffits, & Panels .....	82
Educational Assemblies .....	84
Fume Hoods.....	105
Work Tops.....	112
Shelving & Reagent Racks .....	114
Pegboards.....	120
Miscellaneous Accessories .....	121
Plumbing Information.....	124
Service Fittings .....	125
Sinks & Drain Accessories .....	133
Electrical Fixtures .....	138
Seating .....	140
Technical Information .....	142
Alphabetical Index .....	150
Numerical Index.....	153
Typical Installations.....	148

Signature Series Wood Laboratory Furniture is one of Kewaunee's new generation product lines designed specifically to meet today's quality/value laboratory standards.

The highest standards of wood casework construction, plus multiple choices of premium wood species and beautiful acid-resistant finishes, make the Signature Series furniture the choice of laboratory equipment specifiers and users worldwide.

# Kewaunee Quality Construction

**KEWAUNEE'S SIGNATURE SERIES WOOD FURNITURE** is engineered and fabricated to exacting specifications. From raw material to finished products, our quality assurance program is designed to certify that Kewaunee products exceed established industry standards.

## DOVETAILED DRAWER CONSTRUCTION



**Drawer Sides** - 1/2" nine-ply birch plywood, dovetailed to drawer sub-front and back.

**Drawer Bottom** - 1/4" hardboard set and glued into grooves on all four sides.

**Optional Drawer Bottom:** white melamine on hardboard or birch plywood.



## PULL STYLE OPTIONS



## DOOR HINGE OPTIONS



Full swing opening allows easier access to the cabinet interior. Hinges are rated for 150 lb. load.

## SIGNATURE SERIES WOOD FURNITURE CABINET FEATURES

- **QUALITY CONSTRUCTION** - Joints are multiple doweled, glued, and screwed. Top rails and solid bottom provide complete case rigidity.
- **RIDGED END PANELS** - 3/4" hardwood plywood with front exposed edges faced with solid hardwood.
- **ONE-PIECE CABINET BACK** - 1/4" hardboard, dadoed into end panels and securely fastened.
- **FULL DEPTH SHELVES** - Hardwood plywood, with 1/8" hardwood facing, matching cabinet interior. 3/4" in Base cabinets - 1" in Wall and Tall cabinets.
- **FULL EXTENTION**, ball bearing, drawer slides on all drawers, pullboards, and pullout trays and pans. Load rated at 100 lbs.
- **OPTIONAL ECO-FRIENDLY** - chain of custody FSC® certified plywood, NAUF plywood and particleboards.

# Signature Series Wood Cabinet Styles



## Style-1

### CONVENTIONAL RADIUS LIP

Conventional radius lip doors and drawers. Available in oak, maple, or beech with a wide selection of pull options.



## Style-3

### SILHOUETTE

Radiused edge door and drawer fronts with an integral, recessed pull that is machined into the solid oak hardwood panels. Available in oak only.



## Style-5

### CONTEMPORARY FULL OVERLAY

Square-edged full overlay doors and drawers with vertical matching grain. Available in oak, maple, or beech with a wide selection of pull options.



# Wood Cabinet Information

## Catalog Number Explanation

**E 4 0 W 3 6 2 2 3 6 - 1 0 1 1 AB1C...**

**Cabinet Type**

A = Apron Assemblies and Legs  
 B = Open Cupboard Base Cabinet  
 C = All Cupboard Base Cabinet  
 D = All Drawer Base Cabinet  
 E = Drawer & Cupboard Base Cabinet  
 G = Special Purpose Base Cabinet  
 S = Full Height Storage Cabinet  
 W = Wall Storage Cabinet

**Door & Drawer Configuration**

**Unit Construction and Design**

**Cabinet Height**

**Base Cabinets**  
 31 = 30" High    34 = 32<sup>1</sup>/<sub>2</sub>" High    36 = 35" High

**Aprons**  
 3 = 3" High    5 = 4<sup>17</sup>/<sub>32</sub>" High

**Full Height Cabinets**  
 84 = 84" High

**Wall Cabinet**  
 24 = 24" High    30 = 30" High    36 = 36" High    48 = 48" High

**Cabinet Options**  
(See Cabinet Option list next page)

**Style – Hardware – Material Selections**  
(See Selection list below)

**Special Configurations**  
 L = Reverse Construction/Left Hand Swing  
 A = Adjustable Length (Apron Assemblies)  
 F = Flush with D&D Front (Style 5 only)

**Cabinet Length**  
 Standard Lengths: 15" – 18" – 21" – 24"  
 30" – 36" – 42" – 48"

**Cabinet Depth (Nominal)**

**Base Cabinets**  
 16 = 16"    22 = 22"    28 = 27<sup>1</sup>/<sub>2</sub>"

**Aprons**  
 22 = 22"    28 = 28"

**Wall Cabinets & Full Height Cabinets**  
 12 = 12<sup>3</sup>/<sub>4</sub>"    16 = 16<sup>3</sup>/<sub>4</sub>"    22 = 22"

## Door & Drawer Head Style – Hardware – Material Selections

**E 4 0 W 3 6 2 2 3 6 - 1 0 1 1 AB1C...**

**Door and Drawer Head Style**

1 = Conventional Radius Lip  
 3 = Silhouette  
 5 = Contemporary Full Overlay

**Hinge Type**

0 = Stainless Steel – Five Knuckle  
 2 = Black Powder Coated – Five Knuckle  
 3 = White Powder Coated – Five Knuckle

**Cabinet Exterior Veneer**

1 = Plain Sliced Red Oak  
 2 = Plain Sliced White Maple – available for Style 1 and Style 5 only  
 3 = NAUF Red Oak  
 4 = NAUF White Maple – available for Style 1 and Style 5 only  
 5 = FSC/NAUF Red Oak  
 6 = FSC/NAUF White Maple – available for Style 1 and Style 5 only  
 7 = Steamed Beech – available for Style 1 and Style 5 only  
 8 = FSC/NAUF Steamed Beech – available for Style 1 and Style 5 only

**Pull Type**


\*0 = No Pull  
 1 = Rounded Aluminum  
 2 = Black Wire  
 3 = White Wire  
 4 = Stainless Steel Wire  
 5 = Integral Machined  
 6 = Chrome Plated Wire  
 ‡7 = Beveled Label Holder  
 8 = Square Aluminum

\* pull holes bored on 4" centers, no pull provided  
 ‡ mounts horizontally



# Wood Cabinet Information

## Base Cabinet Dimensional Information



Base cabinets are designed on a 32mm high "unit" system to allow maximum flexibility in the selection of drawer and/or cupboard door sizes.

Standing height cabinets (35" high) have an opening 24 "units" high. Sitting height cabinets (30" high) have an opening 20 "units" high. ADA height cabinets (32 1/2" high) have an opening 22 "units" high. What this means is that any combination of drawer and cupboard door heights equaling 24, 22, or 20 units can be used in a base cabinet.

The adjacent tables indicate the inside clearance dimensions of the base cabinet drawers and cupboards shown in this catalog. Simply match the number code shown on the cabinet to the appropriate code in tables.

The dimensions shown in these tables are minimums. They may be slightly larger on some cabinet types and/or styles.

These tables are intended to show drawer and cupboard clearances, not actual construction dimensions. If actual dimensions are required, contact your Kewaunee representative.

### Wood Drawer Size Codes

Units High	Drawer Opening Height	Inside Dimensions	
		Depth/Height	Length * Front-Back
2	Pullboard		14"
3	2 <sup>13</sup> / <sub>16</sub> "	2 <sup>1</sup> / <sub>8</sub> " †	16 <sup>3</sup> / <sub>4</sub> "
4	4 <sup>1</sup> / <sub>16</sub> "	3 <sup>3</sup> / <sub>8</sub> " †	16 <sup>3</sup> / <sub>4</sub> "
5	5 <sup>5</sup> / <sub>16</sub> "	4 <sup>5</sup> / <sub>8</sub> " †	16 <sup>3</sup> / <sub>4</sub> "
6	6 <sup>9</sup> / <sub>16</sub> "	5 <sup>7</sup> / <sub>8</sub> " †	16 <sup>3</sup> / <sub>4</sub> "
7	7 <sup>13</sup> / <sub>16</sub> "	7 <sup>1</sup> / <sub>8</sub> " †	16 <sup>3</sup> / <sub>4</sub> "
8	9 <sup>1</sup> / <sub>8</sub> "	8 <sup>3</sup> / <sub>8</sub> " †	16 <sup>3</sup> / <sub>4</sub> "
10	11 <sup>5</sup> / <sub>8</sub> "	10 <sup>7</sup> / <sub>8</sub> " †	16 <sup>3</sup> / <sub>4</sub> "
12	14 <sup>1</sup> / <sub>8</sub> "	13 <sup>1</sup> / <sub>8</sub> " †	16 <sup>3</sup> / <sub>4</sub> "

† 1/4" less if top drawer in cabinet  
\* 11 3/4" on 16" deep base cabinets

### Wood Drawer Width I.D. (side-to-side)

Base Cabinet Length	Number of Drawers Across				
	1	2	3	4	6
15"	11 1/2"				
18"	14 1/2"				
21"	17 1/2"				
24"	20 1/2"	8 3/4"			
30"	26 1/2"	11 3/4"			
36"	32 1/2"	14 3/4"	8 3/4"		
42"	38 1/2"	17 3/4"	10 3/4"		
48"	44 1/2"	20 3/4"	12 3/4"	8 7/8"	
60"		26 3/4"	16 3/4"	10 7/8"	
72"		32 3/4"	20 3/4"	14 7/8"	8 7/8"

### Wood Cupboard Size Codes

Units High	Opening Height	Inside Dimensions	
		Height Top-Bottom	Depth * Front-Back
12	14 <sup>3</sup> / <sub>8</sub> "	14 <sup>3</sup> / <sub>8</sub> "	19 <sup>1</sup> / <sub>2</sub> "
15	18 <sup>1</sup> / <sub>8</sub> "	18 <sup>1</sup> / <sub>8</sub> "	19 <sup>1</sup> / <sub>2</sub> "
16	19 <sup>3</sup> / <sub>8</sub> "	19 <sup>3</sup> / <sub>8</sub> "	19 <sup>1</sup> / <sub>2</sub> "
17	20 <sup>5</sup> / <sub>8</sub> "	20 <sup>5</sup> / <sub>8</sub> "	19 <sup>1</sup> / <sub>2</sub> "
18	21 <sup>7</sup> / <sub>8</sub> "	21 <sup>7</sup> / <sub>8</sub> "	19 <sup>1</sup> / <sub>2</sub> "
20	24 <sup>1</sup> / <sub>4</sub> "	24 <sup>1</sup> / <sub>4</sub> "	19 <sup>1</sup> / <sub>2</sub> "
22	26 <sup>3</sup> / <sub>4</sub> "	26 <sup>3</sup> / <sub>4</sub> "	19 <sup>1</sup> / <sub>2</sub> "
24	29 <sup>1</sup> / <sub>4</sub> "	29 <sup>1</sup> / <sub>4</sub> "	19 <sup>1</sup> / <sub>2</sub> "

\* 14 1/2" on 16" deep base cabinets  
26" on 27 1/2" deep base cabinets

### Wood Cupboard Width I.D. (side-to-side)

Base Cabinet Length	Number of Cupboards Across				
	1	2	3	4	6
15"	13 1/2"				
18"	16 1/2"				
21"	19 1/2"				
24"	22 1/2"	10 3/4"			
30"	28 1/2"	13 3/4"			
36"	34 1/2"	16 3/4"	10 3/4"		
42"	40 1/2"	19 3/4"	12 3/4"		
48"	46 1/2"	22 3/4"	14 3/4"	10 7/8"	
60"		28 3/4"	18 3/4"	13 7/8"	
72"		34 3/4"	22 3/4"	16 7/8"	10 7/8"

The dimensions shown in these tables are minimums. They may be slightly larger on some cabinet types and/or styles.

## Cabinet Options – (Add to End of Catalog Number)

- A – Unfinished End Panels
- B11 – 4" Casters w/counter balance & brakes
- B21 – 4" Casters w/counter balance w/o brakes
- B31 – 4" Casters w/counter balance & front brakes
- B41 – 3" Casters w/counter balance & brakes
- B51 – 3" Casters w/counter balance w/o brakes
- B61 – 3" Casters w/counter balance & front .brakes
- D – Full Top Frame
- E – Finished Back
- F – Removable Back (Base cabinet cupboards only)
- G – Toespace Grill
- \*H – Label Holders
- I – 1/4" Laminated Safety Glass
- \*J – Disk Locks w/Lock Panels
- \*K – Disk Locks w/o Lock Panels
- \*L – Padlock Hasps
- M – Levelers for base and full height cabinets

- N1 – Extended End Panel on Left End
- N2 – Extended End Panel on Right End
- N3 – Extended End Panel on Both Ends
- O – Shelf Lips
- P1 – One Pullout Tray - at bottom of cupboard
- P2 – Two Pullout Trays - center & bottom
- P5 – One Pullout Pan - at bottom of cupboard
- P6 – Two Pullout Pans - center & bottom
- R – 1/8" Tempered Safety Glass
- U1 – Toe Space on Left End of Base Cabinet
- U2 – Toe Space on Right End of Base Cabinet
- U3 – Toe Space on Both Ends of Base Cabinet
- V – Epoxy Coated 3/4 Extension Drawer Slide
- W – White Melamine Drawer Bottom
- Y – Intermediate Rails on Style 5 Cabinets
- Z – Birch Plywood Drawer Bottom  
(1/4" thick up to 24" wide – 1/2" over 24" wide)

\* Furnished on all cupboards and drawers on cabinet. — Lock keying and/or numbering must be specified.

# Wood Cabinet Options

## Unfinished End Panels

### Option A

When a cabinet's end panel will not be exposed to view, the cabinet may be furnished with unfinished ends. The A Option provides unfinished ends on both ends of the cabinet.

Unfinished ends may be specified for base cabinets, full height cabinets or wall cabinets.

**Note:** Unfinished ends cannot be specified on Mobile and Suspended Cabinets, on Center Bench End Sink Cabinets or on cabinets with Toe Spaces on the end of the cabinet.

To specify Unfinished Ends, add **A** to the end of the cabinet number.

Example: **E40W362236-1011A**


## Base Cabinet Caster Option

### Option B

The Caster Option adds casters and a finished back to base cabinets when mobility is required. Six caster styles are available and all are furnished with a counter balance to minimize tipping when open.

3" diameter casters add no extra height to the cabinet, 4" diameter casters add 1 1/4" to the overall height.

- B11** – 4" Casters with brakes
- B21** – 4" Casters without brakes
- B31** – 4" Casters with brakes at front only
- B41** – 3" Casters with brakes
- B51**– 3" Casters without brakes
- B61** – 3" Casters with brakes at front only



To specify Casters, add **B** to the end of the cabinet number.


Example: **E40W362236-1011B61**  
(Available on base cabinets only.)

## Full Top Frame

### Option D

Option D adds two side rails to the base cabinet's front and back rails to form a full, four-sided frame at the top of the cabinet.

**Note:** D option not available on base cabinets with pullboards.



To specify a cabinet with Full Top Frame, add **D** to the end of the cabinet number.

Example: **E40W362236-1011D**  
(Available on base cabinets only.)

## Finished Back

### Option E

When a cabinet will be installed away from a wall so that its back is exposed to view, the cabinet can be constructed with a hardwood plywood back flush with the end

panels instead of the standard recessed hardboard back. Finished backs may be specified for base cabinets, full height cabinets or wall cabinets.

**Note:** Finished backs are standard when Option B is specified.

To specify Finished Back, add **E** to the end of the cabinet number.

Example: **E40W362236-1011E**


# Wood Cabinet Options

## Removable Back

## Option F

Removable backs allow easy access to plumbing spaces behind base cabinets. When a removable back is specified, a 16" high access opening is cut in the cupboard back and

covered with a nominal 1/4"-(5.5mm) hardboard back held in place by four screws. A removable back may be added to any base cabinet with a cupboard at least 17 1/2" high.



To specify Removable Back, add **F** to the end of the cabinet number.


Example: **E40W362236-1011F**  
 (Available in base cabinets only.  
 Not available on mobile cabinets.)

## Toespace Grill

## Option G

Toespace Grills are available when a base cabinet or full height cabinet must be installed in front of, or on top of, heating registers or ducts. A

two inch high opening is cut in the cabinet toespace rail and is fitted with a black expanded steel grill.



To specify Toespace Grill, add **G** to the end of the cabinet number.


Example: **E40W362236-1011G**

## Label Holders

## Option H

The H Option adds the F-0260-00 aluminum Label Holder to the upper right hand corner of each drawer head and to the center, 4" down on each door on the cabinet. Overall

size is 3" x 1 5/8" and holds a 2 1/2" x 1 1/8" card. (Label holder is placed on right hand door only of double door cupboards.)



To specify Label Holders, add **H** to the end of the cabinet number.

Example: **E40W362236-1011H**  
 (Added to all cupboards and drawers.)

One-half Actual Size


# Wood Cabinet Options

**Locks with Lock Panels**  
**Locks without Lock Panels**

**Option J**  
**Option K**

All cabinets may be ordered with locks by specifying the K option. In addition, base cabinets may be specified with the J option to add lock panels. Lock panels are

installed between drawers and/or drawers and cupboards to provide additional security when the locks on a cabinet are keyed different.



Lock Panel

**Locks Types**

Cabinet locks are available in two types. Disk Tumbler which offers medium security, and Pin Tumbler which offer a higher level of security. Pin Tumbler also allows more keying arrangements. Both are fabricated with die cast sleeves and cylinders

with a satin nickel face and zinc plated steel cams. Both locks use the same sleeve and cam assembly allowing the lock cylinder to be removed and switched after the lock has been installed. All locks are furnished with two keys unless otherwise specified.

**Disk Tumbler Locks**

**Security:** Medium

**Keying:** Different, Alike, Alike in Sets, Masterkeyed

**Key**

**Changes:** Up to 230

**Master Key**

**Changes:** Up to 5

**Note:** Maximum number of master key changes can vary based on number of key changes.

**Pin Tumbler Locks**

**Security:** Medium-High

**Keying:** Different, Alike, Alike in Sets, Masterkeyed, Grand-masterkeyed, Sub-masterkeyed, Mason Keyed

**Key**

**Changes:** Up to 2000


**Master Key**

**Changes:** Up to 24


**Grand Master Key**

**Changes:** Up to 5

**Note:** Maximum number of master key changes and grand-master key changes can vary based on number of key changes.



**Disk Tumbler Lock**  
(Medium Security)



**Pin Tumbler Lock**  
(Higher Security)

**J** – Disk Locks with Lock Panels  
**K** – Disk Locks w/o Lock Panels

**J2** – Pin Locks with Lock Panels  
**K2** – Pin Locks w/o Lock Panels

**Lock schedule with keying arrangement must be furnished when locks are specified.** See page 12 & 13 for keying arrangement diagrams.

To specify Disk Locks without Lock Panels, add **K** to the end of the cabinet number.

Example: **E40W362236-1011K**  
(Added to all cupboards and drawers.)

To specify Pin Locks with Lock Panels, add **J2** to the end of the cabinet number.

Example: **E40W362236-1011J2**  
(Added to all cupboards and drawers.)

**Padlock Hasps**

**Option L**



Chrome plated padlock hasp mounts in same hose as, and is interchangeable with, tumbler locks. Includes stainless steel wear plate that protects drawer heads and door wear and denting

by the padlock. The L option adds Padlock Hasps and lock panels to all cupboards and drawers on the cabinet.

**(Padlocks not included.)**

To specify Padlock Hasps, add **L** to the end of the cabinet number.

Example: **E40W362236-1011L**  
(Added to all cupboards and drawers.)




# Wood Cabinet Options

## Cabinet Levelers

## Option M

Base cabinets and full height cabinets can be furnished with adjustable levelers when required. Each cabinet is furnished with four (4)

levelers. The levelers are accessible for adjustment through holes in the cabinet bottom that are covered with plug buttons.



To specify Cabinet Leveler, add **M** to the end of the cabinet number.


Example: **E40W362236-1011M**  
(Available on base cabinets and full height cabinets only.)

## Extended End Panel

## Option N

Adds 7" to the width of the base cabinet end panel extending the overall depth from a nominal 22" to 29".

Available for the left, right, or both ends of the cabinet.



- N1** – Extended End Panel on Left End of Cabinet
- N2** – Extended End Panel on Right End of Cabinet
- N3** – Extended End Panel on Both Ends of Cabinet

To specify an Extended End Panel, add **N\_** to the end of the cabinet number.


Example: **E40W362236-1011N1**  
(example adds Extended End to Left End of Cabinet as illustrated)

## Shelf Lips

## Option O

Adds a 1/2" thick solid hardwood lip to the front of all shelves in the cabinet to form a 1/2" high retaining rim. Shelf lips may be added to Base Cabinets, Wall or Counter Mounted

Storage Cabinets or Full Height Storage Cabinets. The overall depth of the shelf remains the same.



To specify Shelf Lips, add **O** to the end of the cabinet number.

Example: **E40W362236-1011O**  
(Added to all shelves in cabinet.)

# Wood Cabinet Options

## Pullout Trays & Pullout Pans


## Option P

- P1** – One Pullout Tray - at bottom of cupboard
- P2** – Two Pullout Trays - center & bottom

Pullout trays are fabricated from 1/2" thick birch plywood with 1" high rims on both sides and back. They are mounted on 100 pound capacity lineal ball bearing full-extension slides. Available in two configurations: one tray, mounted at bottom of

cupboard; or two trays, mounted at bottom and in center of cupboard. When pullout trays are added, the adjustable shelf is omitted.

Pullout trays are available for base cabinets only.



To specify a Pullout Tray, add **P1** or **P2** to the end of the cabinet number.


Example: **E40W362236-1011P2**  
(Available for base cabinets only.)

- P5** – One Pullout Pans - at bottom of cupboard
- P6** – Two Pullout Pans - center & bottom

Pullout Pans are fabricated from 1/2" thick birch plywood with 1" high rims on both sides and back, and 3/4" solid hardwood front. They are mounted on 100 pound capacity lineal ball bearing full-extension slides. Available in two configurations: one tray, mounted

at bottom of cupboard; or two trays, mounted at bottom and in center of cupboard. When pullout trays are added, the adjustable shelf is omitted.

Pullout trays are available for base cabinets only.



To specify a Pullout Pan, add **P5** or **P6** to the end of the cabinet number.


Example: **E40W362236-1011P6**  
(Available for base cabinets only.)

## Laminated Safety Glass Tempered Safety Glass

## Option I Option R

Option I replaces the standard 1/8" float glass in plate glass or glazed door wall, counter mounted, and full height cabinets with 1/4" Laminated Safety Glass

Option R replaces the standard 1/8" float glass in glazed door wall, counter mounted, and full height cabinets with 1/8" Tempered Safety Glass



To specify Laminated Glass, add **I** to the end of the cabinet number.  
Example: **W20W361236-1011I**

To specify Tempered Glass, add **R** to the end of the cabinet number.  
Example: **W20W361236-1011R**

# Wood Cabinet Options


## Toe Space on End of Base Cabinet

## Option U

Base cabinets used at the end of a center bench assembly will usually require a toe space notch at the end

of cabinet as well as the front. Toe Spaces may be specified on the left, the right, or both ends.

- U1** – Toe Space on Left End of Cabinet
- U2** – Toe Space on Right End of Cabinet
- U3** – Toe Space on Both Ends of Cabinet



To specify a Toe Space on End of Base Cabinet, add **U\_** to the end of the cabinet number.


Example: **E40W362236-1011U1**  
(example adds Toe Space to Left End of Cabinet as illustrated)

## Epoxy Coated 3/4 Extention Drawer Slides

## Option V

Epoxy Coated 3/4 Extention drawer slide with positive drawer hold-open feature. Self-closing final 3 inches. Loading capacity of 100 pounds. Replaces standard full extension

slides on all drawers in cabinet except file drawers, when V option is specified.



To specify Epoxy Coated 3/4 Extension Slide, add **V** to the end of the cabinet number.

Example: **E40W362236-1011V**  
(Added to all drawers in cabinet except file drawers.)

## White Drawer Bottoms


## Option W

## Plywood Drawer Bottoms

## Option Z

Standard hardboard drawer bottoms may be replaced with a white melamine coated hardboard bottom or plywood bottoms. White Drawer Bottoms are nominal 1/4"-(5.5mm) thick hardboard

with a white melamine coating. Plywood Drawer Bottoms are 1/4" thick for drawer 24" wide or smaller, 1/2" thick for drawers over 24" wide.



To specify White Drawer Bottoms, add **W** to the end of the cabinet number.

To specify Plywood Drawer Bottoms, add **Z** to the end of the cabinet number.

Example: **E40W362236-1011W**  
(Will be added to all drawers in the cabinet)

Example: **E40W362236-1011Z**  
(Will be added to all drawers in the cabinet)

## Intermediate Rails on Style 5

## Option Y

Style 5 Contemporary Full Overlay base cabinets are furnished standard without intermediate rails between drawers or between cupboards and drawers. They may be added to the cabinet by specifying them with Option Y.

**Note:** Intermediate Rails are automatically added to Style 5 base cabinets when any lock option is selected.

**Note:** Option Y not required for Styles 1 and 3 as they are furnished standard with Intermediate Rails.



To specify Intermediate Rails, add **Y** to the end of the cabinet number.


Example: **D30W362236-1011Y**  
(Available for Style 5 - Contemporary Full Overlay base cabinets only.)

# Lock Keying Options

## Keyed Different

- Disk Tumbler Locks:** Available
- Pin Tumbler Locks:** Available


All locks are keyed different from each other with individual keys for each lock. No key will open more than one lock.



## Keyed Alike

- Disk Tumbler Locks:** Available
- Pin Tumbler Locks:** Available


All locks are keyed alike. One key will open all locks.



## Keyed Alike in Sets

- Disk Tumbler Locks:** Available
- Pin Tumbler Locks:** Available


Locks are grouped into sets, with as many locks as needed in each set. All locks in a set can be opened by a single key, but that key can open no other locks. No key will open all locks.



## Masterkeyed

- Disk Tumbler Locks:** Available
- Pin Tumbler Locks:** Available

All locks are keyed different from each other with individual keys for each lock. Master key provided that can open all locks in the set. Multiple sets are available, but no key can open all locks in all sets.






# Lock Keying Options

## Grand-masterkeyed


**Disk Tumbler Locks:** N/A  
**Pin Tumbler Locks:** Available

All locks are keyed different from each other with individual keys for each lock. Master key provided that can open all locks in the set. Multiple sets are available, with Grand-master key that can open all locks in all sets.



## Sub-masterkeyed & Mason Keyed

**Disk Tumbler Locks:** N/A  
**Pin Tumbler Locks:** Available




### Sub-masterkeyed

All locks are keyed different from each other with individual keys for each lock. Sub-master key provided that can open all locks in sub-set. Multiple sub-sets are available, with Master key that can open all locks in all sub-sets. Multiple sets are available, with Grand-master key that can open all locks in all sets and sub-sets.

### Mason Keyed

Locks are keyed different from each other with individual keys for each lock. Each individual key can also open all common locks in set. Master key provided that can open all locks in set. Multiple sets are available with Grand-master key that can open all locks in all sets.

# Sitting Height Base Cabinets




## Specifications:

These Sitting Height Cabinets are 30" high and a nominal 22" deep. Lengths as shown.


Open and Cupboard cabinets have one adjustable shelf. Open Cubbie cabinets have equal size compartments as shown.

**Note:** Last four digits of part number designates door and drawer head style, hinge, hardwood species, and pull selection. See page 4 for full explanation. Actual product may vary slightly from catalog illustration due to style and hardware selection.

## Open Cabinets



**B00W312215-** 15"L  
**B00W312218-** 18"L  
**B00W312221-** 21"L  
**B00W312224-** 24"L




**B00W312230-** 30"L  
**B00W312236-** 36"L  
**B00W312242-** 42"L  
**B00W312248-** 48"L

## Style 5 Flush Option


Open Cabinets may be ordered 7/8" deeper so that cabinet fronts line up with the doors and drawer fronts on adjacent Style 5 cabinets. To specify the Flush Option, replace the dash in the part number with an **F**.

Example: **B00W312236F**


## Open Cubbie Cabinets



**B01W312212-** 12"L  
**B01W312215-** 15"L  
**B01W312218-** 18"L



**B01W312221-** 21"L  
**B01W312224-** 24"L  
**B01W312230-** 30"L




**B01W312236-** 36"L  
**B01W312242-** 42"L  
**B01W312248-** 48"L

## Style 5 Flush Option


Open Cabinets may be ordered 7/8" deeper so that cabinet fronts line up with the doors and drawer fronts on adjacent Style 5 cabinets. To specify the Flush Option, replace the dash in the part number with an **F**.

Example: **B01W312236F**


## Cupboard Cabinets



**C03W312215L** 15"L  
**C03W312218L** 18"L  
**C03W312221L** 21"L  
**C03W312224L** 24"L




**C03W312215-** 15"L  
**C03W312218-** 18"L  
**C03W312221-** 21"L  
**C03W312224-** 24"L



**C03W312230-** 30"L  
**C03W312236-** 36"L  
**C03W312242-** 42"L  
**C03W312248-** 48"L


## Drawer Cabinets



**D00W312215- 0** 15"L  
**D00W312218- 0** 18"L  
**D00W312221- 0** 21"L  
**D00W312224- 0** 24"L




**D00W312230- 0** 30"L  
**D00W312236- 0** 36"L  
**D00W312242- 0** 42"L  
**D00W312248- 0** 48"L



**D02W312230- 0** 30"L  
**D02W312236- 0** 36"L  
**D02W312242- 0** 42"L  
**D02W312248- 0** 48"L

# Sitting Height Base Cabinets




### Specifications:


These Sitting Height Cabinets are 30" high and a nominal 22" deep. Lengths as shown.

**Note:** Last four digits of part number designates door and drawer head style, hinge, hardwood species, and pull selection. See page 4 for full explanation. Actual product may vary slightly from catalog illustration due to style and hardware selection.


### Drawer Cabinets




- D08W312215- 0 \_\_\_ 15"L
- D08W312218- 0 \_\_\_ 18"L
- D08W312221- 0 \_\_\_ 21"L
- D08W312224- 0 \_\_\_ 24"L



- D10W312215- 0 \_\_\_ 15"L
- D10W312218- 0 \_\_\_ 18"L
- D10W312221- 0 \_\_\_ 21"L
- D10W312224- 0 \_\_\_ 24"L




- D20W312218- 0 \_\_\_ 18"L
- D20W312224- 0 \_\_\_ 24"L




- D26W312218- 0 \_\_\_ 18"L
- D26W312224- 0 \_\_\_ 24"L

Pullboard




- D27W312224- 0 \_\_\_ 24"L

Pullboard




- D28W312218- 0 \_\_\_ 18"L
- D28W312224- 0 \_\_\_ 24"L

Pullboard  
Book Compartment




- D29W312224- 0 \_\_\_ 24"L


Pullboard  
Book Compartment




- D30W312215- 0 \_\_\_ 15"L
- D30W312218- 0 \_\_\_ 18"L
- D30W312221- 0 \_\_\_ 21"L
- D30W312224- 0 \_\_\_ 24"L




- D30W312230- 0 \_\_\_ 30"L
- D30W312236- 0 \_\_\_ 36"L
- D30W312242- 0 \_\_\_ 42"L
- D30W312248- 0 \_\_\_ 48"L



- D31W312230- 0 \_\_\_ 30"L
- D31W312236- 0 \_\_\_ 36"L
- D31W312242- 0 \_\_\_ 42"L
- D31W312248- 0 \_\_\_ 48"L




- D32W312230- 0 \_\_\_ 30"L
- D32W312236- 0 \_\_\_ 36"L
- D32W312242- 0 \_\_\_ 42"L
- D32W312248- 0 \_\_\_ 48"L




- D81W312218- 0 \_\_\_ 18"L
- D81W312221- 0 \_\_\_ 21"L
- D81W312224- 0 \_\_\_ 24"L

2 File Drawers with Hangers



- D85W312218- 0 \_\_\_ 18"L
- D85W312221- 0 \_\_\_ 21"L
- D85W312224- 0 \_\_\_ 24"L


1 File Drawers with Hangers



- D90W312218- 0 \_\_\_ 18"L
- D90W312221- 0 \_\_\_ 21"L
- D90W312224- 0 \_\_\_ 24"L

Pullboard  
1 File Drawer with Hangers

# Sitting Height Base Cabinets




## Specifications:

These Sitting Height Cabinets are 30" high and a nominal 22" deep. Lengths as shown.


E40, E41, and E44 Cabinets have one adjustable shelf.

**Note:** Last four digits of part number designates door and drawer head style, hinge, hardwood species, and pull selection. See page 4 for full explanation. Actual product may vary slightly from catalog illustration due to style and hardware selection.


## Cupboard and Drawer Cabinets




**E24W312224L** \_\_\_\_\_ 24"L  
Pullboard  
Book Compartment




**E24W312224-** \_\_\_\_\_ 24"L  
Pullboard  
Book Compartment




**E25W312218L** \_\_\_\_\_ 18"L  
**E25W312224L** \_\_\_\_\_ 24"L  
Pullboard




**E25W312218-** \_\_\_\_\_ 18"L  
**E25W312224-** \_\_\_\_\_ 24"L  
Pullboard




**E40W312215L** \_\_\_\_\_ 15"L  
**E40W312218L** \_\_\_\_\_ 18"L  
**E40W312221L** \_\_\_\_\_ 21"L  
**E40W312224L** \_\_\_\_\_ 24"L




**E40W312215-** \_\_\_\_\_ 15"L  
**E40W312218-** \_\_\_\_\_ 18"L  
**E40W312221-** \_\_\_\_\_ 21"L  
**E40W312224-** \_\_\_\_\_ 24"L




**E40W312230-** \_\_\_\_\_ 30"L  
**E40W312236-** \_\_\_\_\_ 36"L  
**E40W312242-** \_\_\_\_\_ 42"L  
**E40W312248-** \_\_\_\_\_ 48"L




**E41W312230-** \_\_\_\_\_ 30"L  
**E41W312236-** \_\_\_\_\_ 36"L  
**E41W312242-** \_\_\_\_\_ 42"L  
**E41W312248-** \_\_\_\_\_ 48"L




**E44W312230L** \_\_\_\_\_ 30"L  
**E44W312236L** \_\_\_\_\_ 36"L  
**E44W312242L** \_\_\_\_\_ 42"L  
**E44W312248L** \_\_\_\_\_ 48"L



**E44W312230-** \_\_\_\_\_ 30"L  
**E44W312236-** \_\_\_\_\_ 36"L  
**E44W312242-** \_\_\_\_\_ 42"L  
**E44W312248-** \_\_\_\_\_ 48"L




**E86W312224L** \_\_\_\_\_ 24"L  
Pullboard




**E86W312224-** \_\_\_\_\_ 24"L  
Pullboard

## Diagonal Corner Cabinet




**G14W312215-** \_\_\_\_\_ 32"L



PLAN VIEW OF  
**G14W312215-**

# Sitting Height Base Cabinets




### Specifications:

These Sitting Height Cabinets are 30" high and a nominal 22" deep. Lengths as shown.


Sink Cabinets are provided with lowered back for plumbing access. Sink supports must be ordered separate when required.

**Note:** Last four digits of part number designates door and drawer head style, hinge, hardwood species, and pull selection. See page 4 for full explanation. Actual product may vary slightly from catalog illustration due to style and hardware selection.


### Sink Cabinets



G00W312218L \_\_\_\_\_ 18"L  
G00W312224L \_\_\_\_\_ 24"L



G00W312218- \_\_\_\_\_ 18"L  
G00W312224- \_\_\_\_\_ 24"L




G00W312230- \_\_\_\_\_ 30"L  
G00W312236- \_\_\_\_\_ 36"L  
G00W312242- \_\_\_\_\_ 42"L  
G00W312248- \_\_\_\_\_ 48"L

### Sink Support Kits


K-WSS-18 18"L  
K-WSS-24 24"L  
K-WSS-30 30"L  
K-WSS-36 36"L  
K-WSS-42 42"L  
K-WSS-48 48"L

### Solvent Storage Base Cabinets (available in Style 1 and Style 5 only)




20" Deep

G63W312024L\_0\_0 24"L



20" Deep

G63W312024-\_0\_0 24"L



20" Deep

G63W312030-\_0\_0 30"L  
G63W312036-\_0\_0 36"L  
G63W312048-\_0\_0 48"L

22" Deep

G63W312224L\_0\_0 24"L

22" Deep

G63W312224-\_0\_0 24"L


22" Deep

G63W312230-\_0\_0 30"L  
G63W312236-\_0\_0 36"L  
G63W312248-\_0\_0 48"L


G63 Fume Hood Solvent Storage Cabinets are available as mobile. To add casters to the cabinets above, add **-B7** to the end of the part number.

Example: **G63W312224-5020-B7**  
Mobile – Style 5 – Maple cabinet  
22" deep by 24" long.


### Blind Corner Cabinets



G12W312248- \_\_\_\_\_ 48"L




PLAN VIEW OF  
G12W312248- & G13W312248-



G13W312248- \_\_\_\_\_ 48"L

# ADA Height Base Cabinets




## Specifications:

These ADA Height Cabinets are 32 1/2" high and a nominal 22" deep. Lengths as shown.


Open and Cupboard Cabinets have one adjustable shelf. Open Cubbie Cabinets have equal size compartments as shown.

**Note:** Last four digits of part number designates door and drawer head style, hinge, hardwood species, and pull selection. See page 4 for full explanation. Actual product may vary slightly from catalog illustration due to style and hardware selection.

## Open Cabinets



**B00W342215-**\_\_\_\_\_ 15"L  
**B00W342218-**\_\_\_\_\_ 18"L  
**B00W342221-**\_\_\_\_\_ 21"L  
**B00W342224-**\_\_\_\_\_ 24"L




**B00W342230-**\_\_\_\_\_ 30"L  
**B00W342236-**\_\_\_\_\_ 36"L  
**B00W342242-**\_\_\_\_\_ 42"L  
**B00W342248-**\_\_\_\_\_ 48"L

## Style 5 Flush Option


Open Cabinets may be ordered 7/8" deeper so that cabinet fronts line up with the doors and drawer fronts on adjacent Style 5 cabinets. To specify the Flush Option, replace the dash in the part number with an **F**.

Example: **B00W342236F**\_\_\_\_\_


## Open Cubbie Cabinets



**B01W342212-**\_\_\_\_\_ 12"L  
**B01W342215-**\_\_\_\_\_ 15"L  
**B01W342218-**\_\_\_\_\_ 18"L



**B01W342221-**\_\_\_\_\_ 21"L  
**B01W342224-**\_\_\_\_\_ 24"L  
**B01W342230-**\_\_\_\_\_ 30"L




**B01W342236-**\_\_\_\_\_ 36"L  
**B01W342242-**\_\_\_\_\_ 42"L  
**B01W342248-**\_\_\_\_\_ 48"L

## Style 5 Flush Option


Open Cabinets may be ordered 7/8" deeper so that cabinet fronts line up with the doors and drawer fronts on adjacent Style 5 cabinets. To specify the Flush Option, replace the dash in the part number with an **F**.

Example: **B01W342236F**\_\_\_\_\_


## Cupboard Cabinets



**C03W342215L**\_\_\_\_\_ 15"L  
**C03W342218L**\_\_\_\_\_ 18"L  
**C03W342221L**\_\_\_\_\_ 21"L  
**C03W342224L**\_\_\_\_\_ 24"L




**C03W342215-**\_\_\_\_\_ 15"L  
**C03W342218-**\_\_\_\_\_ 18"L  
**C03W342221-**\_\_\_\_\_ 21"L  
**C03W342224-**\_\_\_\_\_ 24"L




**C03W342230-**\_\_\_\_\_ 30"L  
**C03W342236-**\_\_\_\_\_ 36"L  
**C03W342242-**\_\_\_\_\_ 42"L  
**C03W342248-**\_\_\_\_\_ 48"L

## Drawer Cabinets



**D30W342215-**\_\_\_\_\_ 15"L  
**D30W342218-**\_\_\_\_\_ 18"L  
**D30W342221-**\_\_\_\_\_ 21"L  
**D30W342224-**\_\_\_\_\_ 24"L

# ADA Height Base Cabinets




### Specifications:

These ADA Height Cabinets are 32 1/2" high and a nominal 22" deep. Lengths as shown.


Sink Cabinets are provided with sink supports.

**Note:** Last four digits of part number designates door and drawer head style, hinge, hardwood species, and pull selection. See page 4 for full explanation. Actual product may vary slightly from catalog illustration due to style and hardware selection.


### Cupboard & Drawer Cabinets




- E40W342215L \_\_\_\_\_ 15"L
- E40W342218L \_\_\_\_\_ 18"L
- E40W342221L \_\_\_\_\_ 21"L
- E40W342224L \_\_\_\_\_ 24"L




- E40W342215- \_\_\_\_\_ 15"L
- E40W342218- \_\_\_\_\_ 18"L
- E40W342221- \_\_\_\_\_ 21"L
- E40W342224- \_\_\_\_\_ 24"L




- E40W342230- \_\_\_\_\_ 30"L
- E40W342236- \_\_\_\_\_ 36"L
- E40W342242- \_\_\_\_\_ 42"L
- E40W342248- \_\_\_\_\_ 48"L




- E41W342230- \_\_\_\_\_ 30"L
- E41W342236- \_\_\_\_\_ 36"L
- E41W342242- \_\_\_\_\_ 42"L
- E41W342248- \_\_\_\_\_ 48"L




- E43W342230L \_\_\_\_\_ 30"L
- E43W342236L \_\_\_\_\_ 36"L
- E43W342242L \_\_\_\_\_ 42"L
- E43W342248L \_\_\_\_\_ 48"L



- E43W342230- \_\_\_\_\_ 30"L
- E43W342236- \_\_\_\_\_ 36"L
- E43W342242- \_\_\_\_\_ 42"L
- E43W342248- \_\_\_\_\_ 48"L




- E45W342230L \_\_\_\_\_ 30"L
- E45W342236L \_\_\_\_\_ 36"L
- E45W342242L \_\_\_\_\_ 42"L
- E45W342248L \_\_\_\_\_ 48"L



- E45W342230- \_\_\_\_\_ 30"L
- E45W342236- \_\_\_\_\_ 36"L
- E45W342242- \_\_\_\_\_ 42"L
- E45W342248- \_\_\_\_\_ 48"L

# ADA Height Base Cabinets




## Specifications:

These ADA Height Cabinets are 32 1/2" high and a nominal 22" deep. Lengths as shown.


Sink Cabinets are provided with lowered back for plumbing access. Sink supports must be ordered separate when required.

**Note:** Last four digits of part number designates door and drawer head style, hinge, hardwood species, and pull selection. See page 4 for full explanation. Actual product may vary slightly from catalog illustration due to style and hardware selection.


## Sink Cabinets



- G00W342215L \_\_\_\_\_ 15"L
- G00W342218L \_\_\_\_\_ 18"L
- G00W342221L \_\_\_\_\_ 21"L
- G00W342224L \_\_\_\_\_ 24"L




- G00W342215- \_\_\_\_\_ 15"L
- G00W342218- \_\_\_\_\_ 18"L
- G00W342221- \_\_\_\_\_ 21"L
- G00W342224- \_\_\_\_\_ 24"L




- G00W342230- \_\_\_\_\_ 30"L
- G00W342236- \_\_\_\_\_ 36"L
- G00W342242- \_\_\_\_\_ 42"L
- G00W342248- \_\_\_\_\_ 48"L

## Sink Support Kit for Wood Cabinets


- K-WSS-15 \_\_\_\_\_ 15"L
- K-WSS-18 \_\_\_\_\_ 18"L
- K-WSS-21 \_\_\_\_\_ 21"L
- K-WSS-24 \_\_\_\_\_ 24"L
- K-WSS-30 \_\_\_\_\_ 30"L
- K-WSS-36 \_\_\_\_\_ 36"L
- K-WSS-42 \_\_\_\_\_ 42"L
- K-WSS-48 \_\_\_\_\_ 48"L




## Acid Storage Cabinets



- G21W342215L \_\_\_\_\_ 15"L
- G21W342218L \_\_\_\_\_ 18"L
- G21W342221L \_\_\_\_\_ 21"L
- G21W342224L \_\_\_\_\_ 24"L



- G21W342215- \_\_\_\_\_ 15"L
- G21W342218- \_\_\_\_\_ 18"L
- G21W342221- \_\_\_\_\_ 21"L
- G21W342224- \_\_\_\_\_ 24"L




- G21W342230- \_\_\_\_\_ 30"L
- G21W342236- \_\_\_\_\_ 36"L
- G21W342242- \_\_\_\_\_ 42"L
- G21W342248- \_\_\_\_\_ 48"L

Cabinets are fully lined with 1/8" thick corrosion resistant polyethylene sheets and are provided with one full depth removable shelf.



# ADA Height Base Cabinets




**Specifications:**

These ADA Height Cabinets are 32 1/2" high and a nominal 22" deep. Lengths as shown.


Blind Corner Cabinets have one adjustable shelf.

**Note:** Last four digits of part number designates door and drawer head style, hinge, hardwood species, and pull selection. See page 4 for full explanation. Actual product may vary slightly from catalog illustration due to style and hardware selection.


**Blind Corner Cabinets**



G12W342248- 48"L




PLAN VIEW OF G12W342248- & G13W342248-




G13W342248- 48"L


**Diagonal Corner Cabinets**




G14W342215- \_\_\_\_\_



G15W342215- \_\_\_\_\_  
Provided with Lazy Susan




G16W342215- \_\_\_\_\_  
Provided with Sink Supports



PLAN VIEW OF G14W342215- G15W342215- G16W342215-

# ADA Height Base Cabinets




## Specifications:


These ADA Height Cabinets are 32 1/2" high and a nominal 22" deep. Lengths as shown.

**Note:** Last four digits of part number designates door and drawer head style, hinge, hardwood species, and pull selection. See page 4 for full explanation. Actual product may vary slightly from catalog illustration due to style and hardware selection.

## Wheelchair Accessible Base Cabinet with Doors




The Wheelchair Accessible Base Cabinet with Doors is designed to be used as a sink unit, but may be used anywhere wheelchair accessibility is needed. The cabinet is constructed without a bottom, and with the toe-space rails attached to the doors so it will match adjacent cabinets when the doors are closed. The cabinet incorporates a 3 1/4" top rail to give a minimum of 29 1/2" of clearance from the floor to the bottom of the rail.




- G03W342236-\_\_\_\_\_ 36"L
- G03W342242-\_\_\_\_\_ 42"L
- G03W342248-\_\_\_\_\_ 48"L

## Wheelchair Accessible Base Cabinet without Doors



The Wheelchair Accessible Base Cabinet without Doors is designed to be used as a sink unit, but may be used anywhere wheelchair accessibility is needed. The cabinet is constructed without a bottom, and without doors for easy accessibility. The cabinet incorporates a 5 1/2" top rail to give a minimum of 27" of clearance from the floor to the bottom of the rail and a removable sloping knee space panel to enclose piping and drain lines.

### Style 5 Flush Option




Open Cabinets may be ordered 7/8" deeper so that cabinet fronts line up with the doors and drawer fronts on adjacent Style 5 cabinets. To specify the Flush Option, replace the dash in the part number with an **F**.

- G07W342236-\_\_\_\_\_ 36"L
- G07W342242-\_\_\_\_\_ 42"L
- G07W342248-\_\_\_\_\_ 48"L

Example: **G07W342236F**\_\_\_\_\_

# ADA Height Base Cabinets




### Specifications:

These ADA Height Cabinets are 32 1/2" high and a nominal 22" deep. Lengths as shown.


**Note:** Last four digits of part number designates door and drawer head style, hinge, hardwood species, and pull selection. See page 4 for full explanation. Actual product may vary slightly from catalog illustration due to style and hardware selection.

### Fume Hood Base Cabinets




20" Deep

G08W342015L	15"L
G08W342018L	18"L
G08W342021L	21"L
G08W342024L	24"L



20" Deep

G08W342015-	15"L
G08W342018-	18"L
G08W342021-	21"L
G08W342024-	24"L



20" Deep

G08W342030-	30"L
G08W342036-	36"L
G08W342042-	42"L
G08W342048-	48"L

22" Deep

G08W342215L	15"L
G08W342218L	18"L
G08W342221L	21"L
G08W342224L	24"L


22" Deep

G08W342215-	15"L
G08W342218-	18"L
G08W342221-	21"L
G08W342224-	24"L

22" Deep


G08W342230-	30"L
G08W342236-	36"L
G08W342242-	42"L
G08W342248-	48"L

### Fume Hood Acid Storage Cabinets




20" Deep

G80W342024L	24"L
-------------	------



20" Deep

G80W342024-	24"L
-------------	------



20" Deep

G80W342030-	30"L
G80W342036-	36"L
G80W342048-	48"L

22" Deep

G80W342224L	24"L
-------------	------

22" Deep

G80W342224-	24"L
-------------	------

22" Deep

G80W342230-	30"L
G80W342236-	36"L
G80W342248-	48"L

See Page 36 for Fume Hood Acid Storage Cabinet specification.

# ADA Height Solvent Storage Cabinets

## Specifications:


Wood Solvent Storage Base Cabinets are 32½ or 35" high, and 22" deep. Lengths as shown. Solvent Storage Cabinets are specifically designed for the storage of flammable and combustible liquids. Construction meets the requirements listed by NFPA No. 30-1993 and are UL listed. The cabinet is constructed of 1" exterior grade plywood and is

available in Style 1 and Style 5 only. (Doors on Style 1 are flush-inlay, doors on Style 5 are full-overlay) The right hand door is equipped with a three-point latch. A 2" deep liquid-tight pan covers the entire bottom of the cabinet to contain liquid leaks and spills, and similar pan acts as an adjustable shelf. Two diametrically opposed vents with spark screens


are available when required. The cabinet is labeled: CAUTION – FLAMMABLE – KEEP FIRE AWAY.

**Note:**

Solvent Storage Cabinets may be locked, but can not be keyed or masterkeyed to match other locks in laboratory.




### Vent Locations




**Vent Outlets**  
1½" Male NPT

## Solvent Storage Fume Hood Base Cabinets (available in Style 1 and Style 5 only)




**20" Deep**

**G68W342024L\_0\_0** 24"L



**20" Deep**

**G68W342024-\_0\_0** 24"L



**20" Deep**

**G68W342030-\_0\_0** 30"L

**G68W342036-\_0\_0** 36"L

**G68W342048-\_0\_0** 48"L

G68 Fume Hood Solvent Storage Cabinets are available as mobile. To add casters to the cabinets above, add **-B7** to the end of the part number.

**22" Deep**

**G68W342224L\_0\_0** 24"L

**22" Deep**

**G68W342224-\_0\_0** 24"L

**22" Deep**


**G68W342230-\_0\_0** 30"L

**G68W342236-\_0\_0** 36"L

**G68W342248-\_0\_0** 48"L


Example: **G68W342024-1010-B7**  
Mobile – Style 1 – Oak cabinet  
20" deep by 24" long.

## Solvent Storage Base Cabinets (available in Style 1 and Style 5 only)




**20" Deep**

**G63W342024L\_0\_0** 24"L



**20" Deep**

**G63W342024-\_0\_0** 24"L



**20" Deep**

**G63W342030-\_0\_0** 30"L

**G63W342036-\_0\_0** 36"L

**G63W342048-\_0\_0** 48"L

G63 Fume Hood Solvent Storage Cabinets are available as mobile. To add casters to the cabinets above, add **-B7** to the end of the part number.

**22" Deep**

**G63W342224L\_0\_0** 24"L

**22" Deep**

**G63W342224-\_0\_0** 24"L

**22" Deep**

**G63W342230-\_0\_0** 30"L

**G63W342236-\_0\_0** 36"L

**G63W342248-\_0\_0** 48"L

Example: **G63W342224-5020-B7**  
Mobile – Style 5 – Maple cabinet  
22" deep by 24" long.


# ADA Vacuum Pump Storage Cabinets

## Specifications


Vacuum Pump Storage Fume Hood Cabinets are designed without a bottom to allow vacuum pumps and other equipment to be rolled in or out of the cabinets. The interior is lined with an easy to clean sound absorbant material. Each cabinet is furnished with

a 120 VAC, 20 amp, duplex receptacle mounted on the inside cabinet back and a pilot lighted toggle switch mounted in the top front rail. (Wiring is not included.) Cabinets are also furnished with a 1 1/2" diameter PVC vent pipe for venting or access to the

hood above. (Requires a 2" hole in the hood worksurface.) The toespace rail is attached to the door to allow total access to the cabinet. Cabinet inside clearance at the floor is 14 1/2" or 16 1/2" front-to-back, 25" high, and 3" less than the overall cabinet length.




## Fume Hood Vacuum Pump Storage Cabinets




20" Deep

G35W342024L\_0\_0 24"L



20" Deep

G35W342024-\_0\_0 24"L



20" Deep

G35W342030-\_0\_0 30"L  
G35W342036-\_0\_0 36"L  
G35W342048-\_0\_0 48"L

22" Deep

G35W342224L\_0\_0 24"L

22" Deep


G35W342224-\_0\_0 24"L

22" Deep

G35W342230-\_0\_0 30"L  
G35W342236-\_0\_0 36"L  
G35W342248-\_0\_0 48"L

## Vacuum Pump Carts

Vacuum Pump Carts are designed for use with Kewaunee's Vacuum Pump cabinets shown on this page. The cart consists of a 1" deep pan, fabricated from powder coated or stainless steel and mounted on 4" swivel casters.



### Carts for 20" Deep Cabinets

Part Number	Length
VPCC061414-00_0	14"
VPCC061420-00_0	20"
VPCC061426-00_0	26"
VPCC061432-00_0	32"
VPCC061438-00_0	38"
VPCC061444-00_0	44"

### Carts for 22" Deep Cabinets

Part Number	Length
VPCC061614-00_0	14"
VPCC061620-00_0	20"
VPCC061626-00_0	26"
VPCC061632-00_0	32"
VPCC061638-00_0	38"
VPCC061644-00_0	44"

# ADA Aprons & Legs

## Specifications:

ADA Aprons are 4<sup>17/32</sup>" high with a 2<sup>1/8</sup>" notch in the center of the apron to allow maximum clearance for wheelchairs. They are available in lengths and depths shown.





Dimensions shown are sizes with legs. The actual size of Aprons with

one set of legs is <sup>5</sup>/<sub>8</sub>" less in depth and <sup>5</sup>/<sub>16</sub>" less in length. Aprons with legs on both ends are <sup>5</sup>/<sub>8</sub>" less in depth and length.

Order legs separately in style and height required.





**Note:** Last four digits of part number designates drawer head style, hinge, hardwood species, and pull selection. See page 4 for full explanation. Actual product may vary slightly from catalog illustration due to style and hardware selection.

## ADA Aprons


	No Legs	Legs on One End	Legs on Both Ends
			
	21 <sup>1</sup> / <sub>4</sub> " Deep	21 <sup>1</sup> / <sub>4</sub> " Deep	21 <sup>1</sup> / <sub>4</sub> " Deep
	<b>A05W052236-00_0</b> 36"L	<b>A11W052236-00_0</b> 35 <sup>1</sup> / <sub>2</sub> "L	<b>A13W052235-__0</b> 35"L
	<b>A05W052248-00_0</b> 48"L	<b>A11W052248-00_0</b> 47 <sup>1</sup> / <sub>2</sub> "L	<b>A13W052247-__0</b> 47"L
	<b>A05W052260-00_0</b> 60"L	<b>A11W052260-00_0</b> 59 <sup>1</sup> / <sub>2</sub> "L	<b>A13W052259-__0</b> 59"L
	<b>A05W052272-00_0</b> 72"L	<b>A11W052272-00_0</b> 71 <sup>1</sup> / <sub>2</sub> "L	<b>A13W052271-__0</b> 71"L
	27 <sup>1</sup> / <sub>4</sub> " Deep	27 <sup>1</sup> / <sub>4</sub> " Deep	27 <sup>1</sup> / <sub>4</sub> " Deep
	<b>A05W052836-00_0</b> 36"L	<b>A11W052836-00_0</b> 35 <sup>1</sup> / <sub>2</sub> "L	<b>A13W052835-__0</b> 35"L
	<b>A05W052848-00_0</b> 48"L	<b>A11W052848-00_0</b> 47 <sup>1</sup> / <sub>2</sub> "L	<b>A13W052847-__0</b> 47"L
	<b>A05W052860-00_0</b> 60"L	<b>A11W052860-00_0</b> 59 <sup>1</sup> / <sub>2</sub> "L	<b>A13W052859-__0</b> 59"L
	<b>A05W052872-00_0</b> 72"L	<b>A11W052872-00_0</b> 71 <sup>1</sup> / <sub>2</sub> "L	<b>A13W052871-__0</b> 71"L

## Legs with Leg Shoes

## Legs with Glides

			
<b>A20W3400-00_0</b> 32 <sup>1</sup> / <sub>2</sub> "H	21 <sup>1</sup> / <sub>4</sub> " Deep <b>A21W3422-00_0</b> 32 <sup>1</sup> / <sub>2</sub> "H	<b>A40W3400-00_0</b> 32 <sup>1</sup> / <sub>2</sub> "H	21 <sup>1</sup> / <sub>4</sub> " Deep <b>A41W3422-00_0</b> 32 <sup>1</sup> / <sub>2</sub> "H
	27 <sup>1</sup> / <sub>4</sub> " Deep <b>A21W3428-00_0</b> 32 <sup>1</sup> / <sub>2</sub> "H		27 <sup>1</sup> / <sub>4</sub> " Deep <b>A41W3428-00_0</b> 32 <sup>1</sup> / <sub>2</sub> "H

## Legs with Casters


21 <sup>1</sup> / <sub>4</sub> " Deep <b>A44W3422-00_0</b> 32 <sup>1</sup> / <sub>2</sub> "H
27 <sup>1</sup> / <sub>4</sub> " Deep <b>A44W3428-00_0</b> 32 <sup>1</sup> / <sub>2</sub> "H

# Adjustable Height Tables

## Specifications:

Adjustable Height Tables wrap wood rails around a steel frame for maximum strength. The tables are adjustable in height using a 2½" square solid wood leg with a 1" square stainless steel insert adjustable in 1" increments from 29" to 38" table frame height.

Adjustable Height Table Work Tops are available in three materials: Kemresin, Plastic Laminate and Acid Resistant Plastic Laminate.

Kemresin is available in four colors as shown in the chart to the right. To specify Kemresin color, add the two


letter color code to the end of the Kemresin Work Top part number.

**Example:** X10WT012436-R-BK

Kemresin Color Codes

**BK** = Black      **GR** = Grey  
**SL** = Slate      **PT** = Putty

## Adjustable Height Tables



## Adjustable Height Table Frames


Legs: Adjustable Plastic Swivel Glides

Table Frame	Height	Depth	Length
X10WC302436-00_0	29"-38"	22"	34"
X10WC302442-00_0	29"-38"	22"	40"
X10WC302448-00_0	29"-38"	22"	46"
X10WC302454-00_0	29"-38"	22"	52"
X10WC302460-00_0	29"-38"	22"	58"
X10WC302466-00_0	29"-38"	22"	64"
X10WC302472-00_0	29"-38"	22"	72"
X10WC302496-00_0	29"-38"	22"	96"
X10WC303036-00_0	29"-38"	28"	34"
X10WC303042-00_0	29"-38"	28"	40"
X10WC303048-00_0	29"-38"	28"	46"
X10WC303054-00_0	29"-38"	28"	52"
X10WC303060-00_0	29"-38"	28"	58"
X10WC303066-00_0	29"-38"	28"	64"
X10WC303072-00_0	29"-38"	28"	72"
X10WC303096-00_0	29"-38"	28"	96"

## Adjustable Height Table Work Tops

Kemresin/Epoxy Resin	Black Plastic Laminate	Black Acid Resistant Laminate	Thickness	Depth	Length
X10WT012436-R-___ <i>specify color</i>	X10WT012436-P	X10WT012436-D	1"	24"	36"
X10WT012442-R-___ <i>specify color</i>	X10WT012442-P	X10WT012442-D	1"	24"	42"
X10WT012448-R-___ <i>specify color</i>	X10WT012448-P	X10WT012448-D	1"	24"	48"
X10WT012454-R-___ <i>specify color</i>	X10WT012454-P	X10WT012454-D	1"	24"	54"
X10WT012460-R-___ <i>specify color</i>	X10WT012460-P	X10WT012460-D	1"	24"	60"
X10WT012466-R-___ <i>specify color</i>	X10WT012466-P	X10WT012466-D	1"	24"	66"
X10WT012472-R-___ <i>specify color</i>	X10WT012472-P	X10WT012472-D	1"	24"	72"
X10WT012496-R-___ <i>specify color</i>	X10WT012496-P	X10WT012496-D	1"	24"	96"
X10WT013036-R-___ <i>specify color</i>	X10WT013036-P	X10WT013036-D	1"	30"	36"
X10WT013042-R-___ <i>specify color</i>	X10WT013042-P	X10WT013042-D	1"	30"	42"
X10WT013046-R-___ <i>specify color</i>	X10WT013046-P	X10WT013046-D	1"	30"	48"
X10WT013054-R-___ <i>specify color</i>	X10WT013054-P	X10WT013054-D	1"	30"	54"
X10WT013060-R-___ <i>specify color</i>	X10WT013060-P	X10WT013060-D	1"	30"	60"
X10WT013066-R-___ <i>specify color</i>	X10WT013066-P	X10WT013066-D	1"	30"	66"
X10WT013072-R-___ <i>specify color</i>	X10WT013072-P	X10WT013072-D	1"	30"	72"
X10WT013096-R-___ <i>specify color</i>	X10WT013096-P	X10WT013096-D	1"	30"	96"

# Standing Height Base Cabinets



## Specifications:

These Standing Height Cabinets are 35" high and a nominal 22" deep. Lengths as shown.


Open Cabinets have two adjustable shelves. Open Cubbie Cabinets equal size compartments as shown.

**Note:** Last four digits of part number designates door and drawer head style, hinge, hardwood species, and pull selection. See page 4 for full explanation. Actual product may vary slightly from catalog illustration due to style and hardware selection.

### Open Cabinets



**B00W362215-** 15"L  
**B00W362218-** 18"L  
**B00W362221-** 21"L  
**B00W362224-** 24"L




**B00W362230-** 30"L  
**B00W362236-** 36"L  
**B00W362242-** 42"L  
**B00W362248-** 48"L

### Style 5 Flush Option


Open Cabinets may be ordered 7/8" deeper so that cabinet fronts line up with the doors and drawer fronts on adjacent Style 5 cabinets. To specify the Flush Option, replace the dash in the part number with an "F".

Example: **B00W362236F**\_\_\_\_\_


### Open Cubbie Cabinets



**B01W362212-** 12"L  
**B01W362215-** 15"L  
**B01W362218-** 18"L



**B01W362221-** 21"L  
**B01W362224-** 24"L  
**B01W362230-** 30"L



**B01W362236-** 36"L  
**B01W362242-** 42"L  
**B01W362248-** 48"L


### Style 5 Flush Option

Open Cabinets may be ordered 7/8" deeper so that cabinet fronts line up with the doors and drawer fronts on adjacent Style 5 cabinets. To specify the Flush Option, replace the dash in the part number with an "F".

Example: **B01W362236F**\_\_\_\_\_



# Standing Height Base Cabinets




**Specifications:**

These Standing Height Cabinets are 35" high and a nominal 22" deep. Lengths as shown.


Cupboards Cabinets have one adjustable shelf except for C07 cabinets which have two cupboards.

**Note:** Last four digits of part number designates door and drawer head style, hinge, hardwood species, and pull selection. See page 4 for full explanation. Actual product may vary slightly from catalog illustration due to style and hardware selection.


**Cupboard Cabinets**




- C03W362215L \_\_\_\_\_ 15"L
- C03W362218L \_\_\_\_\_ 18"L
- C03W362221L \_\_\_\_\_ 21"L
- C03W362224L \_\_\_\_\_ 24"L



- C03W362215- \_\_\_\_\_ 15"L
- C03W362218- \_\_\_\_\_ 18"L
- C03W362221- \_\_\_\_\_ 21"L
- C03W362224- \_\_\_\_\_ 24"L




- C03W362230- \_\_\_\_\_ 30"L
- C03W362236- \_\_\_\_\_ 36"L
- C03W362242- \_\_\_\_\_ 42"L
- C03W362248- \_\_\_\_\_ 48"L




- C06W362230- \_\_\_\_\_ 30"L
- C06W362236- \_\_\_\_\_ 36"L
- C06W362242- \_\_\_\_\_ 42"L
- C06W362248- \_\_\_\_\_ 48"L


Divided Cabinet




- C07W362212L \_\_\_\_\_ 12"L
- C07W362215L \_\_\_\_\_ 15"L




- C07W362212- \_\_\_\_\_ 12"L
- C07W362215- \_\_\_\_\_ 15"L



- C07W362218L \_\_\_\_\_ 18"L
- C07W362221L \_\_\_\_\_ 21"L
- C07W362224L \_\_\_\_\_ 24"L




- C07W362218- \_\_\_\_\_ 18"L
- C07W362221- \_\_\_\_\_ 21"L
- C07W362224- \_\_\_\_\_ 24"L



- C15W362230-\_0\_0 \_\_\_\_\_ 30"L
- C15W362236-\_0\_0 \_\_\_\_\_ 36"L
- C15W362242-\_0\_0 \_\_\_\_\_ 42"L
- C15W362248-\_0\_0 \_\_\_\_\_ 48"L

Sliding Panel Doors

# Standing Height Base Cabinets



## Specifications:


These Standing Height Cabinets are 35" high and a nominal 22" deep. Lengths as shown.

**Note:** Last four digits of part number designates door and drawer head style, hinge, hardwood species, and pull selection. See page 4 for full explanation. Actual product may vary slightly from catalog illustration due to style and hardware selection.


## Drawer Cabinets




D00W362212- 0 \_\_ 12"L  
D00W362215- 0 \_\_ 15"L




D00W362218- 0 \_\_ 18"L  
D00W362221- 0 \_\_ 21"L  
D00W362224- 0 \_\_ 24"L




D00W362230- 0 \_\_ 30"L  
D00W362236- 0 \_\_ 36"L  
D00W362242- 0 \_\_ 42"L  
D00W362248- 0 \_\_ 48"L




D02W362230- 0 \_\_ 30"L  
D02W362236- 0 \_\_ 36"L  
D02W362242- 0 \_\_ 42"L  
D02W362248- 0 \_\_ 48"L




D03W362236- 0 \_\_ 36"L  
D03W362248- 0 \_\_ 48"L




D30W362212- 0 \_\_ 12"L  
D30W362215- 0 \_\_ 15"L




D30W362218- 0 \_\_ 18"L  
D30W362221- 0 \_\_ 21"L  
D30W362224- 0 \_\_ 24"L




D30W362230- 0 \_\_ 30"L  
D30W362236- 0 \_\_ 36"L  
D30W362242- 0 \_\_ 42"L  
D30W362248- 0 \_\_ 48"L




D31W362230- 0 \_\_ 30"L  
D31W362236- 0 \_\_ 36"L  
D31W362242- 0 \_\_ 42"L  
D31W362248- 0 \_\_ 48"L




D32W362230- 0 \_\_ 30"L  
D32W362236- 0 \_\_ 36"L  
D32W362242- 0 \_\_ 42"L  
D32W362248- 0 \_\_ 48"L



D33W362236- 0 \_\_ 36"L  
D33W362248- 0 \_\_ 48"L




D35W362230- 0 \_\_ 30"L  
D35W362236- 0 \_\_ 36"L  
D35W362242- 0 \_\_ 42"L  
D35W362248- 0 \_\_ 48"L



D36W362230- 0 \_\_ 30"L  
D36W362236- 0 \_\_ 36"L  
D36W362242- 0 \_\_ 42"L  
D36W362248- 0 \_\_ 48"L

# Standing Height Base Cabinets




### Specifications:

These Standing Height Cabinets are 35" high and a nominal 22" deep. Lengths as shown.


Cupboards have one adjustable shelf.

**Note:** Last four digits of part number designates door and drawer head style, hinge, hardwood species, and pull selection. See page 4 for full explanation. Actual product may vary slightly from catalog illustration due to style and hardware selection.


### Drawer Cabinets




- D37W362230L\_0 \_\_ 30"L
- D37W362236L\_0 \_\_ 36"L
- D37W362242L\_0 \_\_ 42"L
- D37W362248L\_0 \_\_ 48"L




- D37W362230- 0 \_\_ 30"L
- D37W362236- 0 \_\_ 36"L
- D37W362242- 0 \_\_ 42"L
- D37W362248- 0 \_\_ 48"L




- D38W362215- 0 \_\_ 15"L
- D38W362218- 0 \_\_ 18"L
- D38W362221- 0 \_\_ 21"L
- D38W362224- 0 \_\_ 24"L




- D39W362230- 0 \_\_ 30"L
- D39W362236- 0 \_\_ 36"L
- D39W362242- 0 \_\_ 42"L
- D39W362248- 0 \_\_ 48"L




- D40W362230- 0 \_\_ 30"L
- D40W362236- 0 \_\_ 36"L
- D40W362242- 0 \_\_ 42"L
- D40W362248- 0 \_\_ 48"L



- D46W362248- 0 \_\_ 48"L




- D81W362218- 0 \_\_ 18"L
- D81W362221- 0 \_\_ 21"L
- D81W362224- 0 \_\_ 24"L




- D84W362218- 0 \_\_ 18"L
- D84W362221- 0 \_\_ 21"L
- D84W362224- 0 \_\_ 24"L

2 File Drawers with Hangers


Pullboard  
1 File Drawers with Hangers




- D50W362215- 0 \_\_ 15"L
- D50W362218- 0 \_\_ 18"L
- D50W362221- 0 \_\_ 21"L
- D50W362224- 0 \_\_ 24"L




- D51W362230- 0 \_\_ 30"L
- D51W362236- 0 \_\_ 36"L
- D51W362242- 0 \_\_ 42"L
- D51W362248- 0 \_\_ 48"L




- D52W362230- 0 \_\_ 30"L
- D52W362236- 0 \_\_ 36"L
- D52W362242- 0 \_\_ 42"L
- D52W362248- 0 \_\_ 48"L




- D86W362248- 0 \_\_ 48"L



- D87W362218- 0 \_\_ 18"L
- D87W362224- 0 \_\_ 24"L




- D87W362236- 0 \_\_ 36"L
- D87W362248- 0 \_\_ 48"L



- D88W362236- 0 \_\_ 36"L
- D88W362248- 0 \_\_ 48"L

# Standing Height Base Cabinets




### Specifications:

These Standing Height Cabinets are 35" high and a nominal 22" deep. Lengths as shown.


Cupboards have one adjustable shelf.

**Note:** Last four digits of part number designates door and drawer head style, hinge, hardwood species, and pull selection. See page 4 for full explanation. Actual product may vary slightly from catalog illustration due to style and hardware selection.


### Cupboard and Drawer Cabinets




**E00W362212L** \_\_\_\_\_ 12"L  
**E00W362215L** \_\_\_\_\_ 15"L




**E00W362212-** \_\_\_\_\_ 12"L  
**E00W362215-** \_\_\_\_\_ 15"L




**E00W362218L** \_\_\_\_\_ 18"L  
**E00W362221L** \_\_\_\_\_ 21"L  
**E00W362224L** \_\_\_\_\_ 24"L




**E00W362218-** \_\_\_\_\_ 18"L  
**E00W362221-** \_\_\_\_\_ 21"L  
**E00W362224-** \_\_\_\_\_ 24"L




**E03W362230L** \_\_\_\_\_ 30"L  
**E03W362236L** \_\_\_\_\_ 36"L  
**E03W362242L** \_\_\_\_\_ 42"L  
**E03W362248L** \_\_\_\_\_ 48"L




**E03W362230-** \_\_\_\_\_ 30"L  
**E03W362236-** \_\_\_\_\_ 36"L  
**E03W362242-** \_\_\_\_\_ 42"L  
**E03W362248-** \_\_\_\_\_ 48"L




**E04W362236L** \_\_\_\_\_ 36"L  
**E04W362248L** \_\_\_\_\_ 48"L



**E04W362236-** \_\_\_\_\_ 36"L  
**E04W362248-** \_\_\_\_\_ 48"L




**E05W362236L** \_\_\_\_\_ 36"L  
**E05W362248L** \_\_\_\_\_ 48"L



**E05W362236-** \_\_\_\_\_ 36"L  
**E05W362248-** \_\_\_\_\_ 48"L

# Standing Height Base Cabinets




**Specifications:**

These Standing Height Cabinets are 35" high and a nominal 22" deep. Lengths as shown.


Cupboards have one adjustable shelf.

**Note:** Last four digits of part number designates door and drawer head style, hinge, hardwood species, and pull selection. See page 4 for full explanation. Actual product may vary slightly from catalog illustration due to style and hardware selection.


**Cupboard and Drawer Cabinets**




E40W362212L \_\_\_\_\_ 12"L  
E40W362215L \_\_\_\_\_ 15"L




E40W362212- \_\_\_\_\_ 12"L  
E40W362215- \_\_\_\_\_ 15"L




E40W362218L \_\_\_\_\_ 18"L  
E40W362221L \_\_\_\_\_ 21"L  
E40W362224L \_\_\_\_\_ 24"L




E40W362218- \_\_\_\_\_ 18"L  
E40W362221- \_\_\_\_\_ 21"L  
E40W362224- \_\_\_\_\_ 24"L




E40W362230- \_\_\_\_\_ 30"L  
E40W362236- \_\_\_\_\_ 36"L  
E40W362242- \_\_\_\_\_ 42"L  
E40W362248- \_\_\_\_\_ 48"L




E41W362230- \_\_\_\_\_ 30"L  
E41W362236- \_\_\_\_\_ 36"L  
E41W362242- \_\_\_\_\_ 42"L  
E41W362248- \_\_\_\_\_ 48"L




E42W362230L \_\_\_\_\_ 30"L  
E42W362236L \_\_\_\_\_ 36"L  
E42W362242L \_\_\_\_\_ 42"L  
E42W362248L \_\_\_\_\_ 48"L




E42W362230- \_\_\_\_\_ 30"L  
E42W362236- \_\_\_\_\_ 36"L  
E42W362242- \_\_\_\_\_ 42"L  
E42W362248- \_\_\_\_\_ 48"L




E43W362230L \_\_\_\_\_ 30"L  
E43W362236L \_\_\_\_\_ 36"L  
E43W362242L \_\_\_\_\_ 42"L  
E43W362248L \_\_\_\_\_ 48"L




E43W362230- \_\_\_\_\_ 30"L  
E43W362236- \_\_\_\_\_ 36"L  
E43W362242- \_\_\_\_\_ 42"L  
E43W362248- \_\_\_\_\_ 48"L




E44W362230L \_\_\_\_\_ 30"L  
E44W362236L \_\_\_\_\_ 36"L  
E44W362242L \_\_\_\_\_ 42"L  
E44W362248L \_\_\_\_\_ 48"L



E44W362230- \_\_\_\_\_ 30"L  
E44W362236- \_\_\_\_\_ 36"L  
E44W362242- \_\_\_\_\_ 42"L  
E44W362248- \_\_\_\_\_ 48"L




E45W362230L \_\_\_\_\_ 30"L  
E45W362236L \_\_\_\_\_ 36"L  
E45W362242L \_\_\_\_\_ 42"L  
E45W362248L \_\_\_\_\_ 48"L



E45W362230- \_\_\_\_\_ 30"L  
E45W362236- \_\_\_\_\_ 36"L  
E45W362242- \_\_\_\_\_ 42"L  
E45W362248- \_\_\_\_\_ 48"L

# Standing Height Base Cabinets




## Specifications:

These Standing Height Cabinets are 35" high and a nominal 22" deep. Lengths as shown.


Cupboards have one adjustable shelf.

**Note:** Last four digits of part number designates door and drawer head style, hinge, hardwood species, and pull selection. See page 4 for full explanation. Actual product may vary slightly from catalog illustration due to style and hardware selection.


## Cupboard and Drawer Cabinets




**E52W362215L** \_\_\_\_\_ 15"L  
**E52W362218L** \_\_\_\_\_ 18"L  
**E52W362221L** \_\_\_\_\_ 21"L  
**E52W362224L** \_\_\_\_\_ 24"L




**E52W362215-** \_\_\_\_\_ 15"L  
**E52W362218-** \_\_\_\_\_ 18"L  
**E52W362221-** \_\_\_\_\_ 21"L  
**E52W362224-** \_\_\_\_\_ 24"L




**E52W362230-** \_\_\_\_\_ 30"L  
**E52W362236-** \_\_\_\_\_ 36"L  
**E52W362242-** \_\_\_\_\_ 42"L  
**E52W362248-** \_\_\_\_\_ 48"L




**E53W362230-** \_\_\_\_\_ 30"L  
**E53W362236-** \_\_\_\_\_ 36"L  
**E53W362242-** \_\_\_\_\_ 42"L  
**E53W362248-** \_\_\_\_\_ 48"L




**E56W362230-** \_\_\_\_\_ 30"L  
**E56W362236-** \_\_\_\_\_ 36"L  
**E56W362242-** \_\_\_\_\_ 42"L  
**E56W362248-** \_\_\_\_\_ 48"L




**E60W362218L** \_\_\_\_\_ 18"L  
**E60W362224L** \_\_\_\_\_ 24"L




**E60W362218-** \_\_\_\_\_ 18"L  
**E60W362224-** \_\_\_\_\_ 24"L




**E61W362236-** \_\_\_\_\_ 36"L  
**E61W362248-** \_\_\_\_\_ 48"L




**E63W362230L** \_\_\_\_\_ 30"L  
**E63W362236L** \_\_\_\_\_ 36"L  
**E63W362242L** \_\_\_\_\_ 42"L  
**E63W362248L** \_\_\_\_\_ 48"L



**E63W362230-** \_\_\_\_\_ 30"L  
**E63W362236-** \_\_\_\_\_ 36"L  
**E63W362242-** \_\_\_\_\_ 42"L  
**E63W362248-** \_\_\_\_\_ 48"L




**E65W362230L** \_\_\_\_\_ 30"L  
**E65W362236L** \_\_\_\_\_ 36"L  
**E65W362242L** \_\_\_\_\_ 42"L  
**E65W362248L** \_\_\_\_\_ 48"L



**E65W362230-** \_\_\_\_\_ 30"L  
**E65W362236-** \_\_\_\_\_ 36"L  
**E65W362242-** \_\_\_\_\_ 42"L  
**E65W362248-** \_\_\_\_\_ 48"L

# Standing Height Base Cabinets







**Specifications:**

These Standing Height Cabinets are 35" high and a nominal 22" deep. Lengths as shown.

Sink Cabinets are provided with lowered back for plumbing access except as noted. Sink supports must be ordered separate when required.

**Note:** Last four digits of part number designates door and drawer head style, hinge, hardwood species, and pull selection. See page 4 for full explanation. Actual product may vary slightly from catalog illustration due to style and hardware selection.

**Sink Cabinets**





- G00W362215L \_\_\_\_\_ 15"L
- G00W362218L \_\_\_\_\_ 18"L
- G00W362221L \_\_\_\_\_ 21"L
- G00W362224L \_\_\_\_\_ 24"L

- G00W362215- \_\_\_\_\_ 15"L
- G00W362218- \_\_\_\_\_ 18"L
- G00W362221- \_\_\_\_\_ 21"L
- G00W362224- \_\_\_\_\_ 24"L

- G00W362230- \_\_\_\_\_ 30"L
- G00W362236- \_\_\_\_\_ 36"L
- G00W362242- \_\_\_\_\_ 42"L
- G00W362248- \_\_\_\_\_ 48"L

- G01W362260- \_\_\_\_\_ 60"L
- furnished with full height back

**Center Bench End Sink Cabinets**




- G04W362251- \_\_\_\_\_ 50 1/2"L
  - G04W362257- \_\_\_\_\_ 56 1/2"L
  - G04W362259- \_\_\_\_\_ 58 1/2"L
- furnished with full height back


- G04W362252- \_\_\_\_\_ 52"L
  - G04W362258- \_\_\_\_\_ 58"L
  - G04W362260- \_\_\_\_\_ 60"L
- furnished with full height back

**Sink Support Kit for Wood Cabinets**

- K-WSS-15 \_\_\_\_\_ 15"L
- K-WSS-18 \_\_\_\_\_ 18"L
- K-WSS-21 \_\_\_\_\_ 21"L
- K-WSS-24 \_\_\_\_\_ 24"L
- K-WSS-30 \_\_\_\_\_ 30"L
- K-WSS-36 \_\_\_\_\_ 36"L
- K-WSS-42 \_\_\_\_\_ 42"L
- K-WSS-48 \_\_\_\_\_ 48"L
- K-WSS-51 \_\_\_\_\_ 51"L
- K-WSS-52 \_\_\_\_\_ 52"L
- K-WSS-57 \_\_\_\_\_ 57"L
- K-WSS-58 \_\_\_\_\_ 58"L
- K-WSS-59 \_\_\_\_\_ 59"L
- K-WSS-60 \_\_\_\_\_ 60"L



# Standing Height Base Cabinets




**Specifications:**

These Standing Height Cabinets are 35" high and a nominal 22" deep. Lengths as shown.

**Note:** Last four digits of part number designates door and drawer head style, hinge, hardwood species, and pull selection. See page 4 for full explanation. Actual product may vary slightly from catalog illustration due to style and hardware selection.


**Waste Storage**



**G20W362218-00** \_\_\_ 18"L

Furnished with 11" high drawer body with full height drawer front mounted on full extension slides. Includes 6" x 5" (opening) stainless steel waste chute.


**Miscellaneous**




**G23W362218-00** \_\_\_ 18"L

Furnished with 11" high drawer body with full height drawer front mounted on full extension slides.


**Acid Storage Cabinets**



**G21W362215L** \_\_\_\_\_ 15"L  
**G21W362218L** \_\_\_\_\_ 18"L  
**G21W362221L** \_\_\_\_\_ 21"L  
**G21W362224L** \_\_\_\_\_ 24"L




**G21W362215-** \_\_\_\_\_ 15"L  
**G21W362218-** \_\_\_\_\_ 18"L  
**G21W362221-** \_\_\_\_\_ 21"L  
**G21W362224-** \_\_\_\_\_ 24"L



**G21W362230-** \_\_\_\_\_ 30"L  
**G21W362236-** \_\_\_\_\_ 36"L  
**G21W362242-** \_\_\_\_\_ 42"L  
**G21W362248-** \_\_\_\_\_ 48"L

Cabinets are fully lined with 1/8" thick corrosion resistant polyethylene sheets and are provided with one full depth removable shelf.

**Tote Tray Storage**




**G86W362248-** \_\_\_\_\_ 48"L

Tote Tray Storage Cabinet is provided with (24) OP-6551 plastic tote trays. Tote tray size: 19" x 13<sup>3</sup>/<sub>4</sub>" x 3<sup>1</sup>/<sub>2</sub>".



# Standing Height Base Cabinets




**Specifications:**

These Standing Height Cabinets are 35" high and a nominal 22" deep. Lengths as shown.


Blind Corner Cabinets have one adjustable shelf.

**Note:** Last four digits of part number designates door and drawer head style, hinge, hardwood species, and pull selection. See page 4 for full explanation. Actual product may vary slightly from catalog illustration due to style and hardware selection.


**Blind Corner Cabinets**



G12W362248-\_\_\_\_ 48"L




PLAN VIEW OF  
G12W362248- & G13W362248-




G13W362248-\_\_\_\_ 48"L

**Diagonal Corner Cabinets**




G14W362215-\_\_\_\_




G15W362215-\_\_\_\_

Provided with Lazy Susan




G16W362215-\_\_\_\_

Provided with Sink Supports



PLAN VIEW OF  
G14W362215-  
G15W362215-  
G16W362215-

# Fume Hood Base Cabinets




**Specifications:**

Fume Hood Base Cabinets are 35" high and a nominal 22" deep. Lengths as shown.


**Note:** Last four digits of part number designates door and drawer head style, hinge, hardwood species, and pull selection. See page 4 for full explanation. Actual product may vary slightly from catalog illustration due to style and hardware selection.

**Fume Hood Base Cabinets**




**20" Deep**

<b>G08W362012L</b> _____	12"L
<b>G08W362018L</b> _____	18"L
<b>G08W362024L</b> _____	24"L




**20" Deep**

<b>G08W362012-</b> _____	12"L
<b>G08W362018-</b> _____	18"L
<b>G08W362024-</b> _____	24"L



**20" Deep**

<b>G08W362030-</b> _____	30"L
<b>G08W362036-</b> _____	36"L
<b>G08W362048-</b> _____	48"L



**20" Deep**

<b>G08W362060-</b> _____	60"L
--------------------------	------

**22" Deep**

<b>G08W362212L</b> _____	12"L
<b>G08W362218L</b> _____	18"L
<b>G08W362224L</b> _____	24"L

**22" Deep**

<b>G08W362212-</b> _____	12"L
<b>G08W362218-</b> _____	18"L
<b>G08W362224-</b> _____	24"L

**22" Deep**

<b>G08W362230-</b> _____	30"L
<b>G08W362236-</b> _____	36"L
<b>G08W362248-</b> _____	48"L

**22" Deep**

<b>G08W362260-</b> _____	60"L
--------------------------	------

# Fume Hood Base Cabinets


## Specifications:

Acid Storage Fume Hood Base Cabinets are specifically designed for the storage of corrosive chemicals. These cabinets are lined with a molded one piece linear low density polyethylene tub with covered corners and a 1" lip at the bottom front. The cabinet doors are lined

with 1/8" sheet polyethylene and the doors are latched using a nylon roller catch.

Each cabinet is furnished with a PVC vent pipe for venting to the fume hood above. (Requires a 2" hole in the hood worksurface.)


An optional half-depth, removable, phenolic resin shelf is available separately.



## Optional Removable Shelf


- 019624** Shelf for 24" cabinet
- 019625** Shelf for 30" cabinet
- 019626** Shelf for 36" cabinet
- 019627** Shelf for 48" cabinet

## Vent Location




- 24" Cabinet – on Center
- 30" Cabinet – 7 5/8" Right of Center
- 36" Cabinet – 5" Right of Center
- 48" Cabinet – 5" Right of Center
- 60" Cabinet – on Center

## Acid Storage Fume Hood Base Cabinets




20" Deep

G80W362024L \_\_\_\_\_ 24"L




20" Deep

G80W362024- \_\_\_\_\_ 24"L



20" Deep

G80W362030- \_\_\_\_\_ 30"L  
G80W362036- \_\_\_\_\_ 36"L  
G80W362048- \_\_\_\_\_ 48"L



20" Deep

G80W362060- \_\_\_\_\_ 60"L

22" Deep

G80W362224L \_\_\_\_\_ 24"L

22" Deep

G80W362224- \_\_\_\_\_ 24"L

22" Deep

G80W362230- \_\_\_\_\_ 30"L  
G80W362236- \_\_\_\_\_ 36"L  
G80W362248- \_\_\_\_\_ 48"L

22" Deep

G80W362260- \_\_\_\_\_ 60"L

# Solvent Storage Base Cabinets

## Specifications:


Wood Solvent Storage Base Cabinets are 32<sup>1</sup>/<sub>2</sub> or 35" high, and 22" deep. Lengths as shown. Solvent Storage Cabinets are specifically designed for the storage of flammable and combustible liquids. Construction meets the requirements listed by NFPA No. 30-1993 and are UL listed. The cabinet is constructed of 1" exterior grade plywood and is

available in Style 1 and Style 5 only. (Doors on Style 1 are flush-inlay, doors on Style 5 are full-overlay) The right hand door is equipped with a three-point latch. A 2" deep liquid-tight pan covers the entire bottom of the cabinet to contain liquid leaks and spills, and similar pan acts as an adjustable shelf. Two diametrically opposed vents with spark screens


are available when required. The cabinet is labeled: CAUTION – FLAMMABLE – KEEP FIRE AWAY.

**Note:**

Solvent Storage Cabinets may be locked, but can not be keyed or masterkeyed to match other locks in laboratory.




### Vent Locations




**Vent Outlets**  
1 1/2" Male NPT

## Wood Solvent Storage Fume Hood Base Cabinets (available in Style 1 and Style 5 only)




**20" Deep**

G68W362024L\_0\_0 24"L



**20" Deep**

G68W362024-\_0\_0 24"L



**20" Deep**

G68W362030-\_0\_0 30"L  
G68W362036-\_0\_0 36"L  
G68W362048-\_0\_0 48"L

**22" Deep**

G68W362224L\_0\_0 24"L

**22" Deep**

G68W362224-\_0\_0 24"L

**22" Deep**

G68W362230-\_0\_0 30"L  
G68W362236-\_0\_0 36"L  
G68W362248-\_0\_0 48"L

G68 Fume Hood Solvent Storage Cabinets are available as mobile. To add casters to the cabinets above, add **-B7** to the end of the part number.

Example: **G68W362024-1010-B7**  
Mobile – Style 1 – Oak cabinet  
20" deep by 24" long.

## Wood Solvent Storage Base Cabinets (available in Style 1 and Style 5 only)




**20" Deep**

G63W362024L\_0\_0 24"L



**20" Deep**

G63W362024-\_0\_0 24"L



**20" Deep**

G63W362030-\_0\_0 30"L  
G63W362036-\_0\_0 36"L  
G63W362048-\_0\_0 48"L

**22" Deep**

G63W362224L\_0\_0 24"L

**22" Deep**

G63W362224-\_0\_0 24"L

**22" Deep**

G63W362230-\_0\_0 30"L  
G63W362236-\_0\_0 36"L  
G63W362248-\_0\_0 48"L

G63 Fume Hood Solvent Storage Cabinets are available as mobile. To add casters to the cabinets above, add **-B7** to the end of the part number.

Example: **G63W362224-5020-B7**  
Mobile – Style 5 – Maple cabinet  
22" deep by 24" long.


# Vacuum Pump Storage Base Cabinets

## Specifications


Vacuum Pump Storage Fume Hood Cabinets are designed without a bottom to allow vacuum pumps and other equipment to be rolled in or out of the cabinets. The interior is lined with an easy to clean sound absorbant material. Each cabinet is furnished with

a 120 VAC, 20 amp, duplex receptacle mounted on the inside cabinet back and a pilot lighted toggle switch mounted in the top front rail. (Wiring is not included.) Cabinets are also furnished with a 1 1/2" diameter PVC vent pipe for venting or access to the

hood above. (Requires a 2" hole in the hood worksurface.) The toespace rail is attached to the door to allow total access to the cabinet. Cabinet inside clearance at the floor is 14 1/2" or 16 1/2" front-to-back, 25" high, and 3" less than the overall cabinet length.




## Fume Hood Vacuum Pump Storage Cabinets




20" Deep

G35W362024L \_\_\_\_\_ 24"L



20" Deep

G35W362024- \_\_\_\_\_ 24"L



20" Deep

G35W362030- \_\_\_\_\_ 30"L  
 G35W362036- \_\_\_\_\_ 36"L  
 G35W362048- \_\_\_\_\_ 48"L

22" Deep

G35W362224L \_\_\_\_\_ 24"L

22" Deep


G35W362224- \_\_\_\_\_ 24"L

22" Deep

G35W362230- \_\_\_\_\_ 30"L  
 G35W362236- \_\_\_\_\_ 36"L  
 G35W362248- \_\_\_\_\_ 48"L

## Vacuum Pump Carts

Vacuum Pump Carts are designed for use with Kewaunee's Vacuum Pump cabinets shown on this page. The cart consists of a 1" deep pan, fabricated from powder coated or stainless steel and mounted on 4" swivel casters.




### Carts for 20" Deep Cabinets

Part Number	Length
VPCC061414-00_0	14"
VPCC061420-00_0	20"
VPCC061426-00_0	26"
VPCC061432-00_0	32"
VPCC061438-00_0	38"
VPCC061444-00_0	44"

### Carts for 22" Deep Cabinets

Part Number	Length
VPCC061614-00_0	14"
VPCC061620-00_0	20"
VPCC061626-00_0	26"
VPCC061632-00_0	32"
VPCC061638-00_0	38"
VPCC061644-00_0	44"

# Special Purpose Base Cabinets




Also available in ADA height

### Specifications:


Instructor Desk Drawer and Sink Cabinets are designed to be used in Kewaunee's standard Instructor Desk Assemblies. They are 35" high and a nominal 27<sup>1</sup>/<sub>2</sub>" deep. Lengths as shown.  
Sink Cabinets are furnished with sink supports and an electrical box cutout.

**Note:** Last four digits of part number designates door and drawer head style, hinge, hardwood species, and pull selection. See page 4 for full explanation. Actual product may vary slightly from catalog illustration due to style and hardware selection.

### Instructor Desk Drawer Cabinets




- G71W362818-\_\_0\_\_ 18"L
- G71W362824-\_\_0\_\_ 24"L
- G71W342818-\_\_0\_\_ 18"L
- G71W342824-\_\_0\_\_ 24"L




- G73W362836-\_\_0\_\_ 36"L
- G73W362848-\_\_0\_\_ 48"L
- G73W342836-\_\_0\_\_ 36"L
- G73W342848-\_\_0\_\_ 48"L


### Instructor Desk Sink Cabinets




- G79W362818L\_\_\_\_\_ 18"L
- G79W362824L\_\_\_\_\_ 24"L
- G79W342818L\_\_\_\_\_ 18"L
- G79W342824L\_\_\_\_\_ 24"L



- G72W362818-\_\_\_\_\_ 18"L
- G72W362824-\_\_\_\_\_ 24"L
- G72W342818-\_\_\_\_\_ 18"L
- G72W342824-\_\_\_\_\_ 24"L



- G74W362836-\_\_\_\_\_ 36"L
- G74W362848-\_\_\_\_\_ 48"L
- G74W342836-\_\_\_\_\_ 36"L
- G74W342848-\_\_\_\_\_ 48"L



- G78W362836-\_\_\_\_\_ 36"L
- G78W362848-\_\_\_\_\_ 48"L
- G78W342836-\_\_\_\_\_ 36"L
- G78W342848-\_\_\_\_\_ 48"L


Provided with Sink Supports

Provided with Sink Supports

Provided with Sink Supports

Provided with Sink Supports

# Special Purpose Base Cabinets




Also available in ADA height

### Specifications:

These Special Purpose Base Cabinets are 35" high and a nominal 16" deep. Lengths as shown. These cabinets are designed to be used in 18" deep wall benches and 42" deep center benches but can be used anywhere a shallower base cabinet is required.

**Note:** Last four digits of part number designates door and drawer head style, hinge, hardwood species, and pull selection. See page 4 for full explanation. Actual product may vary slightly from catalog illustration due to style and hardware selection.



### Open Storage Cabinet



- G91W361618-00\_0 18"L
- G91W361624-00\_0 24"L
- G91W341618-00\_0 18"L
- G91W341624-00\_0 24"L

Book Comp 2 - Adj. Shelves





### Tote Tray Storage Cabinets



- G94W361618L\_\_\_\_\_ 18"L
- G94W361624L\_\_\_\_\_ 24"L
- G94W361618-\_\_\_\_\_ 18"L
- G94W361624-\_\_\_\_\_ 24"L
- G94W341618L\_\_\_\_\_ 18"L
- G94W341624L\_\_\_\_\_ 24"L
- G94W341618-\_\_\_\_\_ 18"L
- G94W341624-\_\_\_\_\_ 24"L





Book Compartment & 6 - OP-6140 Tote Trays

### Drawer Cabinets






- G93W361618-\_0\_\_ 18"L
- G93W361624-\_0\_\_ 24"L
- G93W341618-\_0\_\_ 18"L
- G93W341624-\_0\_\_ 24"L
- G95W361618-\_0\_\_ 18"L
- G95W361624-\_0\_\_ 24"L
- G95W341618-\_0\_\_ 18"L
- G95W341624-\_0\_\_ 24"L
- G96W361618-\_0\_\_ 18"L
- G96W361624-\_0\_\_ 24"L
- G96W341618-\_0\_\_ 18"L
- G96W341624-\_0\_\_ 24"L
- G98W361636-\_0\_\_ 36"L
- G98W341636-\_0\_\_ 36"L

### Cupboard and Drawer Cabinets




- G90W361618L\_\_\_\_\_ 18"L
- G90W361624L\_\_\_\_\_ 24"L
- G90W341618L\_\_\_\_\_ 18"L
- G90W341624L\_\_\_\_\_ 24"L
- G90W361618-\_\_\_\_\_ 18"L
- G90W361624-\_\_\_\_\_ 24"L
- G90W341618-\_\_\_\_\_ 18"L
- G90W341624-\_\_\_\_\_ 24"L
- G92W361618L\_\_\_\_\_ 18"L
- G92W361624L\_\_\_\_\_ 24"L
- G92W341618L\_\_\_\_\_ 18"L
- G92W341624L\_\_\_\_\_ 24"L
- G92W361618-\_\_\_\_\_ 18"L
- G92W361624-\_\_\_\_\_ 24"L
- G92W341618-\_\_\_\_\_ 18"L
- G92W341624-\_\_\_\_\_ 24"L



- G97W361636-\_\_\_\_\_ 36"L
- G97W341636-\_\_\_\_\_ 36"L
- G99W361618L\_\_\_\_\_ 18"L
- G99W361624L\_\_\_\_\_ 24"L
- G99W341618L\_\_\_\_\_ 18"L
- G99W341624L\_\_\_\_\_ 24"L
- G99W361618-\_\_\_\_\_ 18"L
- G99W361624-\_\_\_\_\_ 24"L
- G99W341618-\_\_\_\_\_ 18"L
- G99W341624-\_\_\_\_\_ 24"L

# Wood Sitting Height Mobile Cabinets



## Specifications:





These Mobile Cabinets are 26<sup>3</sup>/<sub>4</sub>" high and a nominal 22" deep. Lengths as shown. They are designed to fit under 30" high worksurface frames.

Drawer cabinets are furnished with counter weights, multiple-drawer cabinets include an interlock system. All mobile cabinets have finished backs.

**Note:** Last four digits of part number designates door and drawer head style, hinge, hardwood species, and pull selection. See page 42 for full explanation. Actual product may vary slightly from catalog illustration due to style and hardware selection.

### Cupboard Cabinets

### Cupboard Cabinets with Pullboards



M10W272215L \_\_\_\_\_ 15"  
M10W272218L \_\_\_\_\_ 18"  
M10W272221L \_\_\_\_\_ 21"  
M10W272224L \_\_\_\_\_ 24"





M10W272215- \_\_\_\_\_ 15"  
M10W272218- \_\_\_\_\_ 18"  
M10W272221- \_\_\_\_\_ 21"  
M10W272224- \_\_\_\_\_ 24"

M15W272215L \_\_\_\_\_ 15"  
M15W272218L \_\_\_\_\_ 18"  
M15W272221L \_\_\_\_\_ 21"  
M15W272224L \_\_\_\_\_ 24"

M15W272215- \_\_\_\_\_ 15"  
M15W272218- \_\_\_\_\_ 18"  
M15W272221- \_\_\_\_\_ 21"  
M15W272224- \_\_\_\_\_ 24"

### Drawer/Cupboard Cabinets

### Drawer/Cupboard Cabinets with Pullboards



M20W272215L \_\_\_\_\_ 15"  
M20W272218L \_\_\_\_\_ 18"  
M20W272221L \_\_\_\_\_ 21"  
M20W272224L \_\_\_\_\_ 24"




M20W272215- \_\_\_\_\_ 15"  
M20W272218- \_\_\_\_\_ 18"  
M20W272221- \_\_\_\_\_ 21"  
M20W272224- \_\_\_\_\_ 24"

M25W272215L \_\_\_\_\_ 15"  
M25W272218L \_\_\_\_\_ 18"  
M25W272221L \_\_\_\_\_ 21"  
M25W272224L \_\_\_\_\_ 24"

M25W272215- \_\_\_\_\_ 15"  
M25W272218- \_\_\_\_\_ 18"  
M25W272221- \_\_\_\_\_ 21"  
M25W272224- \_\_\_\_\_ 24"

### Drawer Cabinets

### Drawer Cabinets with Pullboards



M30W272215- \_\_\_\_\_ 15"  
M30W272218- \_\_\_\_\_ 18"  
M30W272221- \_\_\_\_\_ 21"  
M30W272224- \_\_\_\_\_ 24"




M40W272215- \_\_\_\_\_ 15"  
M40W272218- \_\_\_\_\_ 18"  
M40W272221- \_\_\_\_\_ 21"  
M40W272224- \_\_\_\_\_ 24"

M35W272215- \_\_\_\_\_ 15"  
M35W272218- \_\_\_\_\_ 18"  
M35W272221- \_\_\_\_\_ 21"  
M35W272224- \_\_\_\_\_ 24"

**Note:** Wood Mobile Cabinets are constructed with finished top panels. Tops may be added but are not required to close the top of the cabinet.

Tops are recommended on wood mobile cabinets with pullboards.

**Tops must be ordered separately.** (see page 43)




M60W272215- \_\_\_\_\_ 15"  
M60W272218- \_\_\_\_\_ 18"  
M60W272221- \_\_\_\_\_ 21"  
M60W272224- \_\_\_\_\_ 24"

M62W272215- \_\_\_\_\_ 15"  
M62W272218- \_\_\_\_\_ 18"  
M62W272221- \_\_\_\_\_ 21"  
M62W272224- \_\_\_\_\_ 24"

M65W272215- \_\_\_\_\_ 15"  
M65W272218- \_\_\_\_\_ 18"  
M65W272221- \_\_\_\_\_ 21"  
M65W272224- \_\_\_\_\_ 24"



# Wood Standing Height Mobile Cabinets



### Specifications:


These Mobile Cabinets are 33" high and a nominal 22" deep. Lengths as shown. They are designed to fit under 36" high worksurface frames.

Drawer cabinets are furnished with counter weights, multiple-drawer cabinets include an interlock system. All mobile cabinets have finished backs.


**Note:** Last four digits of part number designates door and drawer head style, hinge, hardwood species, and pull selection. See page 42 for full explanation. Actual product may vary slightly from catalog illustration due to style and hardware selection.

### Cupboard Cabinets


### Cupboard Cabinets with Pullboards



- M10W332215L \_\_\_\_\_ 15"L
- M10W332218L \_\_\_\_\_ 18"L
- M10W332221L \_\_\_\_\_ 21"L
- M10W332224L \_\_\_\_\_ 24"L



- M10W332215- \_\_\_\_\_ 15"L
- M10W332218- \_\_\_\_\_ 18"L
- M10W332221- \_\_\_\_\_ 21"L
- M10W332224- \_\_\_\_\_ 24"L




- M15W332215L \_\_\_\_\_ 15"L
- M15W332218L \_\_\_\_\_ 18"L
- M15W332221L \_\_\_\_\_ 21"L
- M15W332224L \_\_\_\_\_ 24"L




- M15W332215- \_\_\_\_\_ 15"L
- M15W332218- \_\_\_\_\_ 18"L
- M15W332221- \_\_\_\_\_ 21"L
- M15W332224- \_\_\_\_\_ 24"L

### Drawer/Cupboard Cabinets


### Drawer/Cupboard Cabinets with Pullboards




- M20W332215L \_\_\_\_\_ 15"L
- M20W332218L \_\_\_\_\_ 18"L
- M20W332221L \_\_\_\_\_ 21"L
- M20W332224L \_\_\_\_\_ 24"L



- M20W332215- \_\_\_\_\_ 15"L
- M20W332218- \_\_\_\_\_ 18"L
- M20W332221- \_\_\_\_\_ 21"L
- M20W332224- \_\_\_\_\_ 24"L




- M25W332215L \_\_\_\_\_ 15"L
- M25W332218L \_\_\_\_\_ 18"L
- M25W332221L \_\_\_\_\_ 21"L
- M25W332224L \_\_\_\_\_ 24"L




- M25W332215- \_\_\_\_\_ 15"L
- M25W332218- \_\_\_\_\_ 18"L
- M25W332221- \_\_\_\_\_ 21"L
- M25W332224- \_\_\_\_\_ 24"L

### Drawer Cabinets

### Drawer Cabinets with Pullboards



- M50W332215- \_\_\_\_\_ 15"L
- M50W332218- \_\_\_\_\_ 18"L
- M50W332221- \_\_\_\_\_ 21"L
- M50W332224- \_\_\_\_\_ 24"L




- M45W332215- \_\_\_\_\_ 15"L
- M45W332218- \_\_\_\_\_ 18"L
- M45W332221- \_\_\_\_\_ 21"L
- M45W332224- \_\_\_\_\_ 24"L


**Note:** Wood Mobile Cabinets are constructed with finished top panels. Tops may be added but are not required to close the top of the cabinet.

Tops are recommended on wood mobile cabinets with pullboards.


**Tops must be ordered separately.** (see page 43)




- M62W332215- \_\_\_\_\_ 15"L
- M62W332218- \_\_\_\_\_ 18"L
- M62W332221- \_\_\_\_\_ 21"L
- M62W332224- \_\_\_\_\_ 24"L



- M63W332215- \_\_\_\_\_ 15"L
- M63W332218- \_\_\_\_\_ 18"L
- M63W332221- \_\_\_\_\_ 21"L
- M63W332224- \_\_\_\_\_ 24"L



- M75W332215- \_\_\_\_\_ 15"L
- M75W332218- \_\_\_\_\_ 18"L
- M75W332221- \_\_\_\_\_ 21"L
- M75W332224- \_\_\_\_\_ 24"L



- M67W332215- \_\_\_\_\_ 15"L
- M67W332218- \_\_\_\_\_ 18"L
- M67W332221- \_\_\_\_\_ 21"L
- M67W332224- \_\_\_\_\_ 24"L

# Add-A-Drawer Mobile Cabinets

## Specifications:

**Add-A-Drawer Mobile Cabinets** allow the removal of the top drawer assembly to convert the cabinet from standing height to sitting height. The base section of the add-a-drawer cabinet is 26<sup>3</sup>/<sub>4</sub>" high, the add-on drawer 6<sup>1</sup>/<sub>4</sub>", for an overall height of 33".

All add-a-drawer mobile cabinets have finished backs and counter-weights. The M95 four-drawer cabinet incorporates a drawer interlock system on the bottom three drawers.

The cabinets are a nominal 22" deep with lengths as shown. They are designed to fit under 30" and 36" high worksurface frames.

Cupboards have one adjustable shelf.


Tops of both the base section and the add-on drawer are finished to match the cabinet.

**Note:** Last four digits of part number designates door and drawer head style, hinge, hardwood species, and pull selection. See page 42 for full explanation. Actual product may vary slightly from catalog illustration due to style and hardware selection.



Tops may be added but are not required for a finished cabinet.

**When desired, tops must be ordered separately.** (see next page)

### Wood Add-A-Drawer Mobile Cabinets




### Drawer/Cupboard Cabinets



<b>M90W332215L</b> -----	15"L	<b>M90W332215-</b> -----	15"L
<b>M90W332218L</b> -----	18"L	<b>M90W332218-</b> -----	18"L
<b>M90W332221L</b> -----	21"L	<b>M90W332221-</b> -----	21"L
<b>M90W332224L</b> -----	24"L	<b>M90W332224-</b> -----	24"L

### Drawer Cabinets



<b>M95W332215-</b> -----	15"L
<b>M95W332218-</b> -----	18"L
<b>M95W332221-</b> -----	21"L
<b>M95W332224-</b> -----	24"L


# Mobile Cabinet Tops & Accessories

**Specifications:**

**1" Thick Mobile Cabinet Tops** are available in three depths, 22", 23", or 24" and three materials: Kemresin, Phenolic Resin or Plastic Laminate. They may be used on any Mobile


Cabinet and add 1" to the overall height of the cabinets. All edges of top are finished.

**1" Thick — Mobile Cabinet Top (adds 1" to height of Mobile Cabinet)**




**22" Deep flush to front of case**

Length	Kemresin®	Kemresin Colors	Black Phenolic Resin	Black Plastic Laminate
15"	<b>F94R2215-00-__</b>	<b>BK = Black</b>	<b>F94T2215-00</b>	<b>F94P2215-00</b>
18"	<b>F94R2218-00-__</b>	<b>GR = Grey</b>	<b>F94T2218-00</b>	<b>F94P2218-00</b>
21"	<b>F94R2221-00-__</b>	<b>PT = Putty</b>	<b>F94T2221-00</b>	<b>F94P2221-00</b>
24"	<b>F94R2224-00-__</b>	<b>SL = Slate</b>	<b>F94T2224-00</b>	<b>F94P2224-00</b>
30"	<b>F94R2230-00-__</b>		<b>F94T2230-00</b>	<b>F94P2230-00</b>
36"	<b>F94R2236-00-__</b>		<b>F94T2236-00</b>	<b>F94P2236-00</b>
42"	<b>F94R2242-00-__</b>		<b>F94T2242-00</b>	<b>F94P2242-00</b>
48"	<b>F94R2248-00-__</b>		<b>F94T2248-00</b>	<b>F94P2248-00</b>



**23" Deep flush to overlay door/drawer fronts - 1" overhang of case**


Length	Kemresin®	Kemresin Colors	Black Phenolic Resin	Black Plastic Laminate
15"	<b>F94R2315-00-__</b>	<b>BK = Black</b>	<b>F94T2315-00</b>	<b>F94P2315-00</b>
18"	<b>F94R2318-00-__</b>	<b>GR = Grey</b>	<b>F94T2318-00</b>	<b>F94P2318-00</b>
21"	<b>F94R2321-00-__</b>	<b>PT = Putty</b>	<b>F94T2321-00</b>	<b>F94P2321-00</b>
24"	<b>F94R2324-00-__</b>	<b>SL = Slate</b>	<b>F94T2324-00</b>	<b>F94P2324-00</b>
30"	<b>F94R2330-00-__</b>		<b>F94T2330-00</b>	<b>F94P2330-00</b>
36"	<b>F94R2336-00-__</b>		<b>F94T2336-00</b>	<b>F94P2336-00</b>
42"	<b>F94R2342-00-__</b>		<b>F94T2342-00</b>	<b>F94P2342-00</b>
48"	<b>F94R2348-00-__</b>		<b>F94T2348-00</b>	<b>F94P2348-00</b>




**24" Deep 1" overhang of overlay door/drawer fronts**

Length	Kemresin®	Kemresin Colors	Black Phenolic Resin	Black Plastic Laminate
15"	<b>F94R2415-00-__</b>	<b>BK = Black</b>	<b>F94T2415-00</b>	<b>F94P2415-00</b>
18"	<b>F94R2418-00-__</b>	<b>GR = Grey</b>	<b>F94T2418-00</b>	<b>F94P2418-00</b>
21"	<b>F94R2421-00-__</b>	<b>PT = Putty</b>	<b>F94T2421-00</b>	<b>F94P2421-00</b>
24"	<b>F94R2424-00-__</b>	<b>SL = Slate</b>	<b>F94T2424-00</b>	<b>F94P2424-00</b>
30"	<b>F94R2430-00-__</b>		<b>F94T2430-00</b>	<b>F94P2430-00</b>
36"	<b>F94R2436-00-__</b>		<b>F94T2436-00</b>	<b>F94P2436-00</b>
42"	<b>F94R2442-00-__</b>		<b>F94T2442-00</b>	<b>F94P2442-00</b>
48"	<b>F94R2448-00-__</b>		<b>F94T2448-00</b>	<b>F94P2448-00</b>

**Mobile Cabinet Grabbar**



Length	Painted Steel	Stainless Steel
15"	<b>F73M0015-00</b>	<b>F73S0015-00</b>
18"	<b>F73M0018-00</b>	<b>F73S0018-00</b>
21"	<b>F73M0021-00</b>	<b>F73S0021-00</b>
24"	<b>F73M0024-00</b>	<b>F73S0024-00</b>



# Aprons & Legs

## Specifications:

Aprons are 4<sup>17/32</sup>" high, in lengths and depths shown. They are available with plain rails, drawers or e-book compartments.

Dimensions shown are sizes with legs. The actual size of Standard Aprons with one set of legs is 5/8"


less in depth and 5/16" less in length. Aprons with legs on both ends are 5/8" less in depth and length. Aprons with Flush Panels are the same dimension with or without legs.


Order legs separately in style and height required.


**Note:** Last four digits of part number designates drawer head style, hinge, hardwood species, and pull selection. See page 4 for full explanation. Actual product may vary slightly from catalog illustration due to style and hardware selection.

## Box Aprons for Use Between Two Base Cabinets or Pedestal Legs


**Aprons without Flush Panels** (for use with cabinet styles 1 & 3)



	<b>Plain Rail</b>
<b>20<sup>5/8</sup>" Deep</b>	
<b>A00W052224-00_0</b>	24"L
<b>A00W052230-00_0</b>	30"L
<b>A00W052236-00_0</b>	36"L
<b>A00W052248-00_0</b>	48"L
<b>A00W052260-00_0</b>	60"L
<b>A00W052272-00_0</b>	72"L


	<b>Single Drawer</b>
<b>20<sup>5/8</sup>" Deep</b>	
<b>A01W052224-_0__</b>	24"L
<b>A01W052230-_0__</b>	30"L
<b>A01W052236-_0__</b>	36"L
<b>A01W052248-_0__</b>	48"L


### Single Full Length Drawer


	<b>Single Full Length Drawer</b>
<b>20<sup>5/8</sup>" Deep</b>	
<b>A51W052224-_0__</b>	24"L
<b>A51W052230-_0__</b>	30"L
<b>A51W052236-_0__</b>	36"L
<b>A51W052248-_0__</b>	48"L

## Box Aprons for Use Between Two Base Cabinets or Pedestal Legs

**Aprons with Flush Panels** (for use with cabinet style 5)



	<b>Plain Rail</b>
<b>21<sup>1/2</sup>" Deep</b>	
<b>A00W052224F50_0</b>	24"L
<b>A00W052230F50_0</b>	30"L
<b>A00W052236F50_0</b>	36"L
<b>A00W052248F50_0</b>	48"L
<b>* A00W052260F50_0</b>	60"L
<b>* A00W052272F50_0</b>	72"L

	<b>Single Drawer</b>
<b>21<sup>1/2</sup>" Deep</b>	
<b>A01W052224F50__</b>	24"L
<b>A01W052230F50__</b>	30"L
<b>A01W052236F50__</b>	36"L
<b>A01W052248F50__</b>	48"L

\* two piece panel - split in middle

# Aprons & Legs

## Specifications:

Aprons are 4<sup>17/32</sup>" high, in lengths and depths shown. They are available with plain rails, drawers or e-book compartments.

Dimensions shown are sizes with legs. The actual size of Standard Aprons with one set of legs is <sup>5</sup>/<sub>8</sub>"

less in depth and <sup>5</sup>/<sub>16</sub>" less in length. Aprons with legs on both ends are <sup>5</sup>/<sub>8</sub>" less in depth and length. Aprons with Flush Panels are the same dimension with or without legs.


Order legs separately in style and height required.

**Note:** Last four digits of part number designates drawer head style, hinge, hardwood species, and pull selection. See page 4 for full explanation. Actual product may vary slightly from catalog illustration due to style and hardware selection.

## Box Aprons for Use Between Two Base Cabinets or Pedestal Legs

Aprons without Flush Panels (for use with cabinet styles 1 & 3)


### Double Drawer



20<sup>5</sup>/<sub>8</sub>" Deep

A02W052248-__0__	48"L
A02W052260-__0__	60"L
A02W052272-__0__	72"L


### Single E-Book Compartment



20<sup>5</sup>/<sub>8</sub>" Deep

A03W052224-00_0	24"L
A03W052230-00_0	30"L
A03W052236-00_0	36"L
A03W052248-00_0	48"L

### Double E-Book Compartment




20<sup>5</sup>/<sub>8</sub>" Deep

A04W052248-00_0	48"L
A04W052260-00_0	60"L
A04W052272-00_0	72"L

## Box Aprons for Use Between Two Base Cabinets or Pedestal Legs

Aprons with Flush Panels (for use with cabinet style 5)

### Double Drawer



21<sup>1</sup>/<sub>2</sub>" Deep

A02W052248F50__	48"L
A02W052260F50__	60"L
A02W052272F50__	72"L

# Aprons & Legs

## Specifications:

Aprons are 4<sup>17</sup>/<sub>32</sub>" high, in lengths and depths shown. They are available with plain rails, drawers or e-book compartments.

Dimensions shown are sizes with legs. The actual size of Standard Aprons with one set of legs is 5/8"



less in depth and 5/16" less in length. Aprons with legs on both ends are 5/8" less in depth and length. Aprons with Flush Panels are the same dimension with or without legs.

Order legs separately in style and height required.

**Note:** Last four digits of part number designates drawer head style, hinge, hardwood species, and pull selection. See page 4 for full explanation. Actual product may vary slightly from catalog illustration due to style and hardware selection.

### Apron with Legs on Left End

Aprons without Flush Panels (for use with cabinet styles 1 & 3)




#### 21<sup>1</sup>/<sub>4</sub>" Deep

<b>A10W052224-00_0</b>	23 <sup>1</sup> / <sub>2</sub> "L
<b>A10W052230-00_0</b>	29 <sup>1</sup> / <sub>2</sub> "L
<b>A10W052236-00_0</b>	35 <sup>1</sup> / <sub>2</sub> "L
<b>A10W052248-00_0</b>	47 <sup>1</sup> / <sub>2</sub> "L
<b>A10W052260-00_0</b>	59 <sup>1</sup> / <sub>2</sub> "L
<b>A10W052272-00_0</b>	71 <sup>1</sup> / <sub>2</sub> "L

#### 27<sup>1</sup>/<sub>4</sub>" Deep

<b>A10W052836-00_0</b>	35 <sup>1</sup> / <sub>2</sub> "L
<b>A10W052848-00_0</b>	47 <sup>1</sup> / <sub>2</sub> "L
<b>A10W052860-00_0</b>	59 <sup>1</sup> / <sub>2</sub> "L
<b>A10W052872-00_0</b>	71 <sup>1</sup> / <sub>2</sub> "L

#### Single Drawer



#### 21<sup>1</sup>/<sub>4</sub>" Deep



<b>A15W052224-__0__</b>	23 <sup>1</sup> / <sub>2</sub> "L
<b>A15W052230-__0__</b>	29 <sup>1</sup> / <sub>2</sub> "L
<b>A15W052236-__0__</b>	35 <sup>1</sup> / <sub>2</sub> "L
<b>A15W052248-__0__</b>	47 <sup>1</sup> / <sub>2</sub> "L

#### 27<sup>1</sup>/<sub>4</sub>" Deep

<b>A15W052824-__0__</b>	23 <sup>1</sup> / <sub>2</sub> "L
<b>A15W052830-__0__</b>	29 <sup>1</sup> / <sub>2</sub> "L
<b>A15W052836-__0__</b>	35 <sup>1</sup> / <sub>2</sub> "L
<b>A15W052848-__0__</b>	47 <sup>1</sup> / <sub>2</sub> "L

### Apron with Legs on Left End

Aprons with Flush Panels (for use with cabinet style 5)



#### 22<sup>1</sup>/<sub>4</sub>" Deep

<b>A10W052224F50_0</b>	23 <sup>1</sup> / <sub>2</sub> "L
<b>A10W052230F50_0</b>	29 <sup>1</sup> / <sub>2</sub> "L
<b>A10W052236F50_0</b>	35 <sup>1</sup> / <sub>2</sub> "L
<b>A10W052248F50_0</b>	47 <sup>1</sup> / <sub>2</sub> "L
<b>* A10W052260F50_0</b>	59 <sup>1</sup> / <sub>2</sub> "L
<b>* A10W052272F50_0</b>	71 <sup>1</sup> / <sub>2</sub> "L

\* two piece panel - split in middle

#### 27<sup>3</sup>/<sub>4</sub>" Deep

<b>A10W052836F50_0</b>	35 <sup>1</sup> / <sub>2</sub> "L
<b>A10W052848F50_0</b>	47 <sup>1</sup> / <sub>2</sub> "L
<b>* A10W052860F50_0</b>	59 <sup>1</sup> / <sub>2</sub> "L
<b>* A10W052872F50_0</b>	71 <sup>1</sup> / <sub>2</sub> "L

\* two piece panel - split in middle

#### Single Drawer



#### 22<sup>1</sup>/<sub>4</sub>" Deep

<b>A15W052224F50__</b>	23 <sup>1</sup> / <sub>2</sub> "L
<b>A15W052230F50__</b>	29 <sup>1</sup> / <sub>2</sub> "L
<b>A15W052236F50__</b>	35 <sup>1</sup> / <sub>2</sub> "L
<b>A15W052248F50__</b>	47 <sup>1</sup> / <sub>2</sub> "L

#### 27<sup>3</sup>/<sub>4</sub>" Deep

<b>A15W052824F50__</b>	23 <sup>1</sup> / <sub>2</sub> "L
<b>A15W052830F50__</b>	29 <sup>1</sup> / <sub>2</sub> "L
<b>A15W052836F50__</b>	35 <sup>1</sup> / <sub>2</sub> "L
<b>A15W052848F50__</b>	47 <sup>1</sup> / <sub>2</sub> "L

# Aprons & Legs

## Specifications:

Aprons are 4<sup>17/32</sup>" high, in lengths and depths shown. They are available with plain rails, drawers or e-book compartments.

Dimensions shown are sizes with legs. The actual size of Standard Aprons with one set of legs is 5/8"

less in depth and 5/16" less in length. Aprons with legs on both ends are 5/8" less in depth and length. Aprons with Flush Panels are the same dimension with or without legs.


Order legs separately in style and height required.

**Note:** Last four digits of part number designates drawer head style, hinge, hardwood species, and pull selection. See page 4 for full explanation. Actual product may vary slightly from catalog illustration due to style and hardware selection.

## Apron with Legs on Left End

Aprons without Flush Panels (for use with cabinet styles 1 & 3)

### Double Drawer




#### 21<sup>1/4</sup>" Deep

A16W052248-\_\_0\_\_ 47<sup>1/2</sup>"L  
 A16W052254-\_\_0\_\_ 53<sup>1/2</sup>"L  
 A16W052260-\_\_0\_\_ 59<sup>1/2</sup>"L  
 A16W052272-\_\_0\_\_ 71<sup>1/2</sup>"L

#### 27<sup>1/4</sup>" Deep

A16W052848-\_\_0\_\_ 47<sup>1/2</sup>"L  
 A16W052854-\_\_0\_\_ 53<sup>1/2</sup>"L  
 A16W052860-\_\_0\_\_ 59<sup>1/2</sup>"L  
 A16W052872-\_\_0\_\_ 71<sup>1/2</sup>"L

### Single E-Book Compartment



#### 21<sup>1/4</sup>" Deep

A53W052224-00\_0 23<sup>1/2</sup>"L  
 A53W052230-00\_0 29<sup>1/2</sup>"L  
 A53W052236-00\_0 35<sup>1/2</sup>"L  
 A53W052248-00\_0 47<sup>1/2</sup>"L


#### 27<sup>1/4</sup>" Deep

A53W052824-00\_0 23<sup>1/2</sup>"L  
 A53W052830-00\_0 29<sup>1/2</sup>"L  
 A53W052836-00\_0 35<sup>1/2</sup>"L  
 A53W052848-00\_0 47<sup>1/2</sup>"L

## Apron with Legs on Left End

Aprons with Flush Panels (for use with cabinet style 5)

### Double Drawer



#### 22<sup>1/4</sup>" Deep

A16W052248F50\_\_ 47<sup>1/2</sup>"L  
 A16W052254F50\_\_ 53<sup>1/2</sup>"L  
 A16W052260F50\_\_ 59<sup>1/2</sup>"L  
 A16W052272F50\_\_ 71<sup>1/2</sup>"L

#### 27<sup>3/4</sup>" Deep

A16W052848F50\_\_ 47<sup>1/2</sup>"L  
 A16W052854F50\_\_ 53<sup>1/2</sup>"L  
 A16W052860F50\_\_ 59<sup>1/2</sup>"L  
 A16W052872F50\_\_ 71<sup>1/2</sup>"L

# Aprons & Legs

## Specifications:

Aprons are 4<sup>17/32</sup>" high, in lengths and depths shown. They are available with plain rails, drawers or e-book compartments.

Dimensions shown are sizes with legs. The actual size of Standard Aprons with one set of legs is 5/8" less in depth and 5/16" less in length.




Aprons with legs on both ends are 5/8" less in depth and length. Aprons with Flush Panels are the same dimension with or without legs.

Order legs separately in style and height required.

**Note:** Last four digits of part number designates drawer head style, hinge, hardwood species, and pull selection. See page 4 for full explanation. Actual product may vary slightly from catalog illustration due to style and hardware selection.

### Apron with Legs on Right End

Aprons without Flush Panels (for use with cabinet styles 1 & 3)



#### 21<sup>1/4</sup>" Deep

<b>A10W052224-00_0</b>	23 <sup>1/2</sup> "L
<b>A10W052230-00_0</b>	29 <sup>1/2</sup> "L
<b>A10W052236-00_0</b>	35 <sup>1/2</sup> "L
<b>A10W052248-00_0</b>	47 <sup>1/2</sup> "L
<b>A10W052260-00_0</b>	59 <sup>1/2</sup> "L
<b>A10W052272-00_0</b>	71 <sup>1/2</sup> "L

#### 21<sup>1/4</sup>" Deep

<b>A17W052224-__0__</b>	23 <sup>1/2</sup> "L
<b>A17W052230-__0__</b>	29 <sup>1/2</sup> "L
<b>A17W052236-__0__</b>	35 <sup>1/2</sup> "L
<b>A17W052248-__0__</b>	47 <sup>1/2</sup> "L

#### 27<sup>1/4</sup>" Deep




<b>A10W052836-00_0</b>	35 <sup>1/2</sup> "L
<b>A10W052848-00_0</b>	47 <sup>1/2</sup> "L
<b>A10W052860-00_0</b>	59 <sup>1/2</sup> "L
<b>A10W052872-00_0</b>	71 <sup>1/2</sup> "L

#### 27<sup>1/4</sup>" Deep

<b>A17W052824-__0__</b>	23 <sup>1/2</sup> "L
<b>A17W052830-__0__</b>	29 <sup>1/2</sup> "L
<b>A17W052836-__0__</b>	35 <sup>1/2</sup> "L
<b>A17W052848-__0__</b>	47 <sup>1/2</sup> "L

### Apron with Legs on Right End

Aprons with Flush Panels (for use with cabinet style 5)



#### 22<sup>1/4</sup>" Deep

<b>A10W052224F50_0</b>	23 <sup>1/2</sup> "L
<b>A10W052230F50_0</b>	29 <sup>1/2</sup> "L
<b>A10W052236F50_0</b>	35 <sup>1/2</sup> "L
<b>A10W052248F50_0</b>	47 <sup>1/2</sup> "L
<b>* A10W052260F50_0</b>	59 <sup>1/2</sup> "L
<b>* A10W052272F50_0</b>	71 <sup>1/2</sup> "L

#### 22<sup>1/4</sup>" Deep

<b>A17W052224F50__</b>	23 <sup>1/2</sup> "L
<b>A17W052230F50__</b>	29 <sup>1/2</sup> "L
<b>A17W052236F50__</b>	35 <sup>1/2</sup> "L
<b>A17W052248F50__</b>	47 <sup>1/2</sup> "L

\* two piece panel - split in middle

#### 27<sup>3/4</sup>" Deep

<b>A10W052836F50_0</b>	35 <sup>1/2</sup> "L
<b>A10W052848F50_0</b>	47 <sup>1/2</sup> "L
<b>* A10W052860F50_0</b>	59 <sup>1/2</sup> "L
<b>* A10W052872F50_0</b>	71 <sup>1/2</sup> "L

#### 27<sup>3/4</sup>" Deep

<b>A17W052824F50__</b>	23 <sup>1/2</sup> "L
<b>A17W052830F50__</b>	29 <sup>1/2</sup> "L
<b>A17W052836F50__</b>	35 <sup>1/2</sup> "L
<b>A17W052848F50__</b>	47 <sup>1/2</sup> "L

\* two piece panel - split in middle



# Aprons & Legs

## Specifications:

Aprons are 4<sup>17/32</sup>" high, in lengths and depths shown. They are available with plain rails, drawers or e-book compartments.

Dimensions shown are sizes with legs. The actual size of Standard Aprons with one set of legs is 5/8"

less in depth and 5/16" less in length. Aprons with legs on both ends are 5/8" less in depth and length. Aprons with Flush Panels are the same dimension with or without legs.


Order legs separately in style and height required.

**Note:** Last four digits of part number designates drawer head style, hinge, hardwood species, and pull selection. See page 4 for full explanation. Actual product may vary slightly from catalog illustration due to style and hardware selection.

## Apron with Legs on Right End

Aprons without Flush Panels (for use with cabinet styles 1 & 3)

### Double Drawer




#### 21<sup>1/4</sup>" Deep

A18W052248-__0__	47 <sup>1/2</sup> "L
A18W052254-__0__	53 <sup>1/2</sup> "L
A18W052260-__0__	59 <sup>1/2</sup> "L
A18W052272-__0__	71 <sup>1/2</sup> "L

#### 27<sup>1/4</sup>" Deep

A18W052848-__0__	47 <sup>1/2</sup> "L
A18W052854-__0__	53 <sup>1/2</sup> "L
A18W052860-__0__	59 <sup>1/2</sup> "L
A18W052872-__0__	71 <sup>1/2</sup> "L

### Single E-Book Compartment



#### 21<sup>1/4</sup>" Deep

A54W052224-00_0	23 <sup>1/2</sup> "L
A54W052230-00_0	29 <sup>1/2</sup> "L
A54W052236-00_0	35 <sup>1/2</sup> "L
A54W052248-00_0	47 <sup>1/2</sup> "L


#### 27<sup>1/4</sup>" Deep

A54W052824-00_0	23 <sup>1/2</sup> "L
A54W052830-00_0	29 <sup>1/2</sup> "L
A54W052836-00_0	35 <sup>1/2</sup> "L
A54W052848-00_0	47 <sup>1/2</sup> "L

## Apron with Legs on Right End

Aprons with Flush Panels (for use with cabinet style 5)

### Double Drawer



#### 22<sup>1/4</sup>" Deep

A18W052248F50__	47 <sup>1/2</sup> "L
A18W052254F50__	53 <sup>1/2</sup> "L
A18W052260F50__	59 <sup>1/2</sup> "L
A18W052272F50__	71 <sup>1/2</sup> "L

#### 27<sup>3/4</sup>" Deep

A18W052848F50__	47 <sup>1/2</sup> "L
A18W052854F50__	53 <sup>1/2</sup> "L
A18W052860F50__	59 <sup>1/2</sup> "L
A18W052872F50__	71 <sup>1/2</sup> "L

# Aprons & Legs

## Specifications:

Aprons are 4<sup>17/32</sup>" high, in lengths and depths shown. They are available with plain rails, drawers or e-book compartments.

Dimensions shown are sizes with legs. The actual size of Standard Aprons with one set of legs is 5/8" less in depth and 5/16" less in length.


Aprons with legs on both ends are 5/8" less in depth and length. Aprons with Flush Panels are the same dimension with or without legs.


Order legs separately in style and height required.


**Note:** Last four digits of part number designates drawer head style, hinge, hardwood species, and pull selection. See page 4 for full explanation. Actual product may vary slightly from catalog illustration due to style and hardware selection.

## Free-Standing Table Aprons

**Aprons without Flush Panels** (for use with cabinet styles 1 & 3)





<b>Plain Rail</b>	
	
<b>21<sup>1/4</sup>" Deep</b>	
<b>A30W052235-00_0</b>	35"L
<b>A30W052247-00_0</b>	47"L
<b>A30W052259-00_0</b>	59"L
<b>A30W052271-00_0</b>	71"L
<b>27<sup>1/4</sup>" Deep</b>	
<b>A30W052835-00_0</b>	35"L
<b>A30W052847-00_0</b>	47"L
<b>A30W052859-00_0</b>	59"L
<b>A30W052871-00_0</b>	71"L


<b>Single Drawer</b>	
	
<b>21<sup>1/4</sup>" Deep</b>	
<b>A31W052235-__0__</b>	35"L
<b>A31W052247-__0__</b>	47"L
<b>27<sup>1/4</sup>" Deep</b>	
<b>A31W052835-__0__</b>	35 <sup>1/2</sup> "L
<b>A31W052847-__0__</b>	47 <sup>1/2</sup> "L

## Free-Standing Table Aprons

**Aprons with Flush Panels** (for use with cabinet style 5)



<b>Plain Rail</b>	
	
<b>21<sup>1/4</sup>" Deep</b>	
<b>A30W052235F50_0</b>	35"L
<b>A30W052247F50_0</b>	47"L
<b>* A30W052259F50_0</b>	59"L
<b>* A30W052271F50_0</b>	71"L
* two piece panel - split in middle	
<b>27<sup>3/4</sup>" Deep</b>	
<b>A30W052835F50_0</b>	35"L
<b>A30W052847F50_0</b>	47"L
<b>* A30W052859F50_0</b>	59"L
<b>* A30W052871F50_0</b>	71"L
* two piece panel - split in middle	

<b>Single Drawer</b>	
	
<b>21<sup>1/4</sup>" Deep</b>	
<b>A31W052235F50__</b>	35"L
<b>A31W052247F50__</b>	47"L
<b>27<sup>3/4</sup>" Deep</b>	
<b>A31W052835F50__</b>	35 <sup>1/2</sup> "L
<b>A31W052847F50__</b>	47 <sup>1/2</sup> "L

# Aprons & Legs

## Specifications:

Aprons are 4<sup>17/32</sup>" high, in lengths and depths shown. They are available with plain rails, drawers or e-book compartments.

Dimensions shown are sizes with legs. The actual size of Standard Aprons with one set of legs is 5/8"

less in depth and 5/16" less in length. Aprons with legs on both ends are 5/8" less in depth and length. Aprons with Flush Panels are the same dimension with or without legs.


Order legs separately in style and height required.

**Note:** Last four digits of part number designates drawer head style, hinge, hardwood species, and pull selection. See page 4 for full explanation. Actual product may vary slightly from catalog illustration due to style and hardware selection.

## Free-Standing Table Aprons

**Aprons without Flush Panels** (for use with cabinet styles 1 & 3)

### Double Drawer




#### 21<sup>1/4</sup>" Deep

A32W052247-__0__	47"L
A32W052253-__0__	53"L
A32W052259-__0__	59"L
A32W052271-__0__	71"L

#### 27<sup>1/4</sup>" Deep

A32W052847-__0__	47"L
A32W052853-__0__	53"L
A32W052859-__0__	59"L
A32W052871-__0__	71"L

### Single E-Book Compartment




#### 21<sup>1/4</sup>" Deep

A33W052235-00_0	35"L
A33W052247-00_0	47"L

#### 27<sup>1/4</sup>" Deep

A33W052835-00_0	35 <sup>1/2</sup> "L
A33W052847-00_0	47 <sup>1/2</sup> "L

### Double E-Book Compartment



#### 21<sup>1/4</sup>" Deep

A34W052247-00_0	47"L
A34W052253-00_0	53"L
A34W052259-00_0	59"L
A34W052271-00_0	71"L


#### 27<sup>1/4</sup>" Deep

A34W052847-00_0	47"L
A34W052853-00_0	53"L
A34W052859-00_0	59"L
A34W052871-00_0	71"L

## Free-Standing Table Aprons

**Aprons with Flush Panels** (for use with cabinet style 5)

### Double Drawer



#### 21<sup>1/4</sup>" Deep

A32W052247F50__	47"L
A32W052253F50__	53"L
A32W052259F50__	59"L
A32W052271F50__	71"L

#### 27<sup>3/4</sup>" Deep

A32W052847F50__	47"L
A32W052853F50__	53"L
A32W052859F50__	59"L
A32W052871F50__	71"L

# Aprons & Legs

## Specifications:

Aprons are  $4\frac{17}{32}$ " high, in lengths and depths shown. They are available with plain rails, drawers or e-book compartments.

Dimensions shown are sizes with legs. The actual size of Standard Aprons with one set of legs is  $\frac{5}{8}$ "


less in depth and  $\frac{5}{16}$ " less in length. Aprons with legs on both ends are  $\frac{5}{8}$ " less in depth and length. Aprons with Flush Panels are the same dimension with or without legs.

Order legs separately in style and height required.


**Note:** Last four digits of part number designates drawer head style, hinge, hardwood species, and pull selection. See page 4 for full explanation. Actual product may vary slightly from catalog illustration due to style and hardware selection.

## Double-sided Free-Standing Table Aprons

**Aprons without Flush Panels** (for use with cabinet styles 1 & 3)



### Plain Rail



### 32 $\frac{1}{2}$ " Deep


<b>A35W053447-00_0</b>	47"L
<b>A35W053453-00_0</b>	53"L
<b>A35W053459-00_0</b>	59"L
<b>A35W053471-00_0</b>	71"L

### 38 $\frac{1}{2}$ " Deep


<b>A35W054047-00_0</b>	47"L
<b>A35W054053-00_0</b>	53"L
<b>A35W054059-00_0</b>	59"L
<b>A35W054071-00_0</b>	71"L

## Double-sided Free-Standing Table Aprons

**Aprons with Flush Panels** (for use with cabinet style 5)



### Plain Rail



### 32 $\frac{1}{2}$ " Deep

<b>A35W053447F50_0</b>	47"L
<b>* A35W053453F50_0</b>	53"L
<b>* A35W053459F50_0</b>	59"L
<b>* A35W053471F50_0</b>	71"L

\* two piece panel - split in middle

### 38 $\frac{1}{2}$ " Deep

<b>A35W054047F50_0</b>	47"L
<b>* A35W054053F50_0</b>	53"L
<b>* A35W054059F50_0</b>	59"L
<b>* A35W054071F50_0</b>	71"L

\* two piece panel - split in middle

# Aprons & Legs

## Specifications:

Aprons are 4<sup>17/32</sup>" high, in lengths and depths shown. They are available with plain rails, drawers or e-book compartments.

Dimensions shown are sizes with legs. The actual size of Standard Aprons with one set of legs is <sup>5</sup>/<sub>8</sub>"

less in depth and <sup>5</sup>/<sub>16</sub>" less in length. Aprons with legs on both ends are <sup>5</sup>/<sub>8</sub>" less in depth and length. Aprons with Flush Panels are the same dimension with or without legs.


Order legs separately in style and height required.

**Note:** Last four digits of part number designates drawer head style, hinge, hardwood species, and pull selection. See page 4 for full explanation. Actual product may vary slightly from catalog illustration due to style and hardware selection.

## Double-sided Free-Standing Table Aprons

Aprons without Flush Panels *(for use with cabinet styles 1 & 3)*

### Single E-Book Compartment




#### 32<sup>1</sup>/<sub>2</sub>" Deep

<b>A38W053447-00_0</b>	47"L
<b>A38W053453-00_0</b>	53"L
<b>A38W053459-00_0</b>	59"L
<b>A38W053471-00_0</b>	71"L

#### 38<sup>1</sup>/<sub>2</sub>" Deep

<b>A38W054047-00_0</b>	47"L
<b>A38W054053-00_0</b>	53"L
<b>A38W054059-00_0</b>	59"L
<b>A38W054071-00_0</b>	71"L

### Double E-Book Compartment



#### 32<sup>1</sup>/<sub>2</sub>" Deep

<b>A39W053447-00_0</b>	47"L
<b>A39W053453-00_0</b>	53"L
<b>A39W053459-00_0</b>	59"L
<b>A39W053471-00_0</b>	71"L

#### 38<sup>1</sup>/<sub>2</sub>" Deep

<b>A39W054047-00_0</b>	47"L
<b>A39W054053-00_0</b>	53"L
<b>A39W054059-00_0</b>	59"L
<b>A39W054071-00_0</b>	71"L

# Aprons & Legs

## Specifications:

Aprons are 4<sup>17</sup>/<sub>32</sub>" high, in lengths and depths shown. They are available with plain rails, drawers or e-book compartments.

Dimensions shown are sizes with legs. The actual size of Standard Aprons with one set of legs is 5/8"


less in depth and 5/16" less in length. Aprons with legs on both ends are 5/8" less in depth and length. Aprons with Flush Panels are the same dimension with or without legs.

Order legs separately in style and height required.

**Note:** Last four digits of part number designates drawer head style, hinge, hardwood species, and pull selection. See page 4 for full explanation. Actual product may vary slightly from catalog illustration due to style and hardware selection.

## Apron with Single Leg on One End

**Aprons without Flush Panels** (for use with cabinet styles 1 & 3)



### Plain Rail



### 20<sup>5</sup>/<sub>8</sub>" Deep


<b>A60W052212-00_0</b>	12"L
<b>A60W052218-00_0</b>	18"L
<b>A60W052224-00_0</b>	24"L
<b>A60W052230-00_0</b>	30"L
<b>A60W052236-00_0</b>	36"L

### 26<sup>5</sup>/<sub>8</sub>" Deep

<b>A60W052812-00_0</b>	12"L
<b>A60W052818-00_0</b>	18"L
<b>A60W052824-00_0</b>	24"L
<b>A60W052830-00_0</b>	30"L
<b>A60W052836-00_0</b>	36"L

## Apron with Single Leg on One End

**Aprons with Flush Panels** (for use with cabinet style 5)



### Plain Rail



### 22" Deep

<b>A60W052212F50_0</b>	12"L
<b>A60W052218F50_0</b>	18"L
<b>A60W052224F50_0</b>	24"L
<b>A60W052230F50_0</b>	30"L
<b>A60W052236F50_0</b>	36"L

### 28" Deep

<b>A60W052812F50_0</b>	12"L
<b>A60W052818F50_0</b>	18"L
<b>A60W052824F50_0</b>	24"L
<b>A60W052830F50_0</b>	30"L
<b>A60W052836F50_0</b>	36"L

# Aprons & Legs

## Specifications:

Aprons are 4<sup>17</sup>/<sub>32</sub>" high, in lengths and depths shown. They are available with plain rails, drawers or e-book compartments.

Dimensions shown are sizes with legs. The actual size of Standard Aprons with one set of legs is 5/8" less in depth and 5/16" less in length.


Aprons with legs on both ends are 5/8" less in depth and length. Aprons with Flush Panels are the same dimension with or without legs.

Order legs separately in style and height required.


**Note:** Last four digits of part number designates drawer head style, hinge, hardwood species, and pull selection. See page 4 for full explanation. Actual product may vary slightly from catalog illustration due to style and hardware selection.

### Apron with Single Leg on Left End

Aprons without Flush Panels (for use with cabinet styles 1 & 3)



#### Single E-Book Compartment



##### 21 1/4" Deep


A63W052224-00_0	24"L
A63W052230-00_0	30"L
A63W052236-00_0	36"L

##### 27 1/4" Deep


A63W052824-00_0	24"L
A63W052830-00_0	30"L
A63W052836-00_0	36"L

### Apron with Single Leg on Right End

Aprons without Flush Panels (for use with cabinet styles 1 & 3)



#### Single E-Book Compartment



##### 21 1/4" Deep

A64W052224-00_0	24"L
A64W052230-00_0	30"L
A64W052236-00_0	36"L

##### 27 1/4" Deep

A64W052824-00_0	24"L
A64W052830-00_0	30"L
A64W052836-00_0	36"L

# Aprons & Legs

## Specifications:

Aprons are 4<sup>17</sup>/<sub>32</sub>" high, in lengths and depths shown. They are available with plain rails, drawers or e-Book compartments.

Dimensions shown are sizes with legs. The actual size of Standard Aprons with one set of legs is 5/8" less in depth and 5/16" less in length.


Aprons with legs on both ends are 5/8" less in depth and length. Aprons with Flush Panels are the same dimension with or without legs.

Order legs separately in style and height required.


**Note:** Last four digits of part number designates drawer head style, hinge, hardwood species, and pull selection. See page 4 for full explanation. Actual product may vary slightly from catalog illustration due to style and hardware selection.

## Double-sided Apron with Single Leg on One End

Aprons without Flush Panels (for use with cabinet styles 1 & 3)



### Plain Rail




#### 32<sup>1</sup>/<sub>2</sub>" Deep

<b>A65W053418-00_0</b>	18"L
<b>A65W053424-00_0</b>	24"L
<b>A65W053430-00_0</b>	30"L
<b>A65W053436-00_0</b>	36"L

#### 38<sup>1</sup>/<sub>2</sub>" Deep

<b>A65W054018-00_0</b>	18"L
<b>A65W054024-00_0</b>	24"L
<b>A65W054030-00_0</b>	30"L
<b>A65W054036-00_0</b>	36"L

### Single E-Book Compartment



#### 32<sup>1</sup>/<sub>2</sub>" Deep


<b>A68W053418-00_0</b>	18"L
<b>A68W053424-00_0</b>	24"L
<b>A68W053430-00_0</b>	30"L
<b>A68W053436-00_0</b>	36"L

#### 38<sup>1</sup>/<sub>2</sub>" Deep

<b>A68W054018-00_0</b>	18"L
<b>A68W054024-00_0</b>	24"L
<b>A68W054030-00_0</b>	30"L
<b>A68W054036-00_0</b>	36"L

## Double-sided Apron with Single Leg on One End

Aprons with Flush Panels (for use with cabinet style 5)



### Plain Rail



#### 34" Deep

<b>A65W053418F50_0</b>	18"L
<b>A65W053424F50_0</b>	24"L
<b>A65W053430F50_0</b>	30"L
<b>A65W053436F50_0</b>	36"L

#### 40" Deep

<b>A65W054018F50_0</b>	18"L
<b>A65W054024F50_0</b>	24"L
<b>A65W054030F50_0</b>	30"L
<b>A65W054036F50_0</b>	36"L



# Aprons & Legs


## Specifications:

Apron Rails and Adjustable Apron Assemblies can be cut to fit in the field. All Apron Rails and Adjustable Apron Assemblies are furnished with hardware for attachment to adjacent units.

**Note:** Last four digits of part number designates drawer head style, hinge, hardwood species, and pull selection. See page 4 for full explanation. Actual product may vary slightly from catalog illustration due to style and hardware selection.


## Apron Rails

### Apron Rails without Flush Panels *(for use with cabinet styles 1 & 3)*




<b>A00W0524-00_0</b>	24"L	<b>A00W0324-00_0</b>	24"L
<b>A00W0530-00_0</b>	30"L	<b>A00W0330-00_0</b>	30"L
<b>A00W0536-00_0</b>	36"L	<b>A00W0336-00_0</b>	36"L
<b>A00W0548-00_0</b>	48"L	<b>A00W0348-00_0</b>	48"L
<b>A00W0560-00_0</b>	60"L	<b>A00W0360-00_0</b>	60"L
<b>A00W0572-00_0</b>	72"L	<b>A00W0372-00_0</b>	72"L

### Apron Rails with Flush Panels *(for use with cabinet style 5)*




<b>A00W0524F50_0</b>	24"L	<b>A00W0324F50_0</b>	24"L
<b>A00W0530F50_0</b>	30"L	<b>A00W0330F50_0</b>	30"L
<b>A00W0536F50_0</b>	36"L	<b>A00W0336F50_0</b>	36"L
<b>A00W0548F50_0</b>	48"L	<b>A00W0348F50_0</b>	48"L
<b>A00W0560F50_0</b>	60"L	<b>A00W0360F50_0</b>	60"L
<b>A00W0572F50_0</b>	72"L	<b>A00W0372F50_0</b>	72"L

## Adjustable Apron Assemblies




### 4 1/2" High Plain Rail

<b>A00W052224A00_0</b>	24"L
<b>A00W052248A00_0</b>	48"L
<b>A00W052272A00_0</b>	72"L




### 3" High Plain Rail

<b>A00W032224A00_0</b>	24"L
<b>A00W032248A00_0</b>	48"L
<b>A00W032272A00_0</b>	72"L



### Single Drawer

<b>A01W052224A_0</b>	24"L
<b>A01W052248A_0</b>	48"L
<b>A01W052272A_0</b>	72"L



### Single E-Book Compartment

<b>A08W052224A00_0</b>	24"L
<b>A08W052248A00_0</b>	48"L
<b>A08W052272A00_0</b>	72"L

# Aprons & Legs


## Specifications:

Legs are 2 1/2" x 2 1/2" square, in lengths shown. They are available with leg shoes and holddown angles for permanent mounting to the floor, or with nylon glides or casters for mobility.

Leg Assemblies are complete with hanger bolt, nut and washer for attaching to apron assembly.


**Note:** Last four digits of part number designates drawer head style, hinge, hardwood species, and pull selection. See page 4 for full explanation. Actual product may vary slightly from catalog illustration due to style and hardware selection.

### Single Leg with Leg Shoe



A24W3100-00\_0 30"H  
A24W3400-00\_0 32 1/2"H  
A24W3600-00\_0 35"H


### Leg Pair with Leg Shoes



**For use with Standard Aprons without Flush Panels (cabinet style 1 & 3)**

A20W3100-00\_0 30"H  
A20W3400-00\_0 32 1/2"H  
A20W3600-00\_0 35"H

### Leg Pair with Glides




A40W3100-00\_0 30"H  
A40W3400-00\_0 32 1/2"H  
A40W3600-00\_0 35"H

**For use with Aprons with Flush Panels (cabinet style 5)**

A20W3100F00\_0 30"H  
A20W3400F00\_0 32 1/2"H  
A20W3600F00\_0 35"H

A40W3100F00\_0 30"H  
A40W3400F00\_0 32 1/2"H  
A40W3600F00\_0 35"H


## Spreaders



A09W0335-00\_0 35"L  
A09W0347-00\_0 47"L  
A09W0353-00\_0 53"L  
A09W0359-00\_0 59"L  
A09W0371-00\_0 71"L

For use with A30 Free-standing Table Aprons and either A21, A41, or A44 Leg Assemblies. Furnished with all needed hardware. Select Spreader to match apron length.

## Vibration Isolator



### KTS-070

Can be installed on any open frame table to minimize floor vibrations for balance procedures. Consists of a rubber pad formed around a threaded stud and a metal bracket with a T-nut. Pads can be adjusted to any desired height. Specify quantity needed.

# Aprons & Legs

## Specifications:

Legs are 2 1/2" x 2 1/2" square, in lengths shown. They are available with leg shoes and holddown angles for permanent mounting to the floor, or with nylon glides or casters for mobility.




Leg Assemblies are complete with

hanger bolt, nut and washer for attaching to apron assembly.

Pedestal Legs are 1 1/2" thick. Depths and heights shown. They may be used instead of base cabinets to support aprons, tops or frames.

**Note:** Last four digits of part number designates drawer head style, hinge, hardwood species, and pull selection. See page 4 for full explanation. Actual product may vary slightly from catalog illustration due to style and hardware selection.

Leg Assemblies w/Shoes	Leg Assemblies w/Glides	Leg Assemblies w/Casters
------------------------	-------------------------	--------------------------



**For use with Standard Aprons without Flush Panels (cabinet style 1 & 3)**

Depth	Part Number	Height	Part Number	Height	Part Number	Height
21 1/4" Deep	A21W3122-00_0	30"H	A41W3122-00_0	30"H	A44W3122-00_0	30"H
	A21W3422-00_0	32 1/2"H	A41W3422-00_0	32 1/2"H	A44W3422-00_0	32 1/2"H
	A21W3622-00_0	35"H	A41W3622-00_0	35"H	A44W3622-00_0	35"H
27 3/4" Deep	A21W3128-00_0	30"H	A41W3128-00_0	30"H	A44W3128-00_0	30"H
	A21W3428-00_0	32 1/2"H	A41W3428-00_0	32 1/2"H	A44W3428-00_0	32 1/2"H
	A21W3628-00_0	35"H	A41W3628-00_0	35"H	A44W3628-00_0	35"H
32 1/2" Deep	A21W3128-00_0	30"H	A41W3128-00_0	30"H	A44W3128-00_0	30"H
	A21W3428-00_0	32 1/2"H	A41W3428-00_0	32 1/2"H	A44W3428-00_0	32 1/2"H
	A21W3628-00_0	35"H	A41W3628-00_0	35"H	A44W3628-00_0	35"H
38 1/2" Deep	A21W3128-00_0	30"H	A41W3128-00_0	30"H	A44W3128-00_0	30"H
	A21W3428-00_0	32 1/2"H	A41W3428-00_0	32 1/2"H	A44W3428-00_0	32 1/2"H
	A21W3628-00_0	35"H	A41W3628-00_0	35"H	A44W3628-00_0	35"H

**For use with Aprons with Flush Panels (cabinet style 5)**

Depth	Part Number	Height	Part Number	Height	Part Number	Height
22 1/4" Deep	A21W3122F00_0	30"H	A41W3122F00_0	30"H	A44W3122F00_0	30"H
	A21W3422F00_0	32 1/2"H	A41W3422F00_0	32 1/2"H	A44W3422F00_0	32 1/2"H
	A21W3622F00_0	35"H	A41W3622F00_0	35"H	A44W3622F00_0	35"H
27 3/4" Deep	A21W3128F00_0	30"H	A41W3128F00_0	30"H	A44W3128F00_0	30"H
	A21W3428F00_0	32 1/2"H	A41W3428F00_0	32 1/2"H	A44W3428F00_0	32 1/2"H
	A21W3628F00_0	35"H	A41W3628F00_0	35"H	A44W3628F00_0	35"H
34" Deep	A21W3134F00_0	30"H	A41W3134F00_0	30"H	A44W3134F00_0	30"H
	A21W3434F00_0	32 1/2"H	A41W3434F00_0	32 1/2"H	A44W3434F00_0	32 1/2"H
	A21W3634F00_0	35"H	A41W3634F00_0	35"H	A44W3634F00_0	35"H
40" Deep	A21W3140F00_0	30"H	A41W3140F00_0	30"H	A44W3140F00_0	30"H
	A21W3440F00_0	32 1/2"H	A41W3440F00_0	32 1/2"H	A44W3440F00_0	32 1/2"H
	A21W3640F00_0	35"H	A41W3640F00_0	35"H	A44W3640F00_0	35"H

# Aprons & Legs



## Specifications:

Legs are 2 1/2" x 2 1/2" square, in lengths shown. They are available with leg shoes and holddown angles for permanent mounting to the floor, or with nylon glides or casters for mobility.

Leg Assemblies are complete with hanger bolt, nut and washer for attaching to apron assembly.

**Note:** Last four digits of part number designates drawer head style, hinge, hardwood species, and pull selection. See page 4 for full explanation. Actual product may vary slightly from catalog illustration due to style and hardware selection.

## Pedestal Legs



### For use with Standard Aprons without Flush Panels (cabinet style 1 & 3)

<b>21 1/4" Deep</b>		<b>50 1/2" Deep</b>	
A25W3122-00_0	30"H	A26W3151-00_0	30"H
A25W3422-00_0	33"H	A26W3451-00_0	33"H
A25W3622-00_0	35"H	A26W3651-00_0	35"H
<b>22 1/4" Deep</b>		<b>56 1/2" Deep</b>	
A25W3123-00_0	30"H	A26W3157-00_0	30"H
A25W3423-00_0	33"H	A26W3457-00_0	33"H
A25W3623-00_0	35"H	A26W3657-00_0	35"H
<b>29 1/4" Deep</b>		<b>58 1/2" Deep</b>	
A25W3130-00_0	30"H	A26W3159-00_0	30"H
A25W3430-00_0	33"H	A26W3459-00_0	33"H
A25W3630-00_0	35"H	A26W3659-00_0	35"H

Back edge of single-depth Pedestal Legs are unfinished. To specify the back edge to be banded and finished, add an **E** to the end of the part number.

Example: **A25W3622-0010E**


### For use with Aprons with Flush Panels (cabinet style 5)

<b>22 1/8" Deep</b>		<b>52 1/4" Deep</b>	
A25W3122F00_0	30"H	A26W3151F00_0	30"H
A25W3422F00_0	33"H	A26W3451F00_0	33"H
A25W3622F00_0	35"H	A26W3651F00_0	35"H
<b>23 1/8" Deep</b>		<b>58 1/4" Deep</b>	
A25W3123F00_0	30"H	A26W3157F00_0	30"H
A25W3423F00_0	33"H	A26W3457F00_0	33"H
A25W3623F00_0	35"H	A26W3657F00_0	35"H
<b>30 1/8" Deep</b>		<b>60 1/4" Deep</b>	
A25W3130F00_0	30"H	A26W3159F00_0	30"H
A25W3430F00_0	33"H	A26W3459F00_0	33"H
A25W3630F00_0	35"H	A26W3659F00_0	35"H

Back edge of single-depth Pedestal Legs are unfinished. To specify the back edge to be banded and finished, add an **E** to the end of the part number.

Example: **A25W3622F5010E**

# Work Top Support and Keyboard Tray




**Specifications:**

Work Top Support Frames may be used anywhere maximum knee space is wanted but top support is required. Fabricated of heavy gauge steel channels with a solid wood facing, it is only 7/8" high and is designed to mount between base units or pedestal legs.


The keyboard tray may be added to any of the frames or may be used independently.

**Work Top Support Frame**




- A80W002230-00\_0 30"L
- A80W002236-00\_0 36"L
- A80W002242-00\_0 42"L
- A80W002248-00\_0 48"L

**Keyboard Tray**



- F-6853-00 22"L
- 22 1/4" Long x 10 1/2" Deep x 3 1/2" High

**Knee Space Panels**



**36" High**


- |               |       |               |        |
|---------------|-------|---------------|--------|
| X-WK2436-00_0 | 24" L | X-WK4836-00_0 | 48" L  |
| X-WK3036-00_0 | 30" L | X-WK5436-00_0 | 54" L* |
| X-WK3636-00_0 | 36" L | X-WK6036-00_0 | 60" L* |
| X-WK4236-00_0 | 42" L | X-WK7236-00_0 | 72" L* |

3/4" thick plywood panels finished on face side only. Designed to be field-cut to exact height required. Furnished with two cleats.

**Note:** To specify panel to be finished on both sides, add an **-S** to the part number. Example: **X-WK4236-0010-S**

\* Knee Space Panels over 48" long are furnished in two sections.

**Knee Space Panels with Base Rail**



**36" High**


- |                |       |                |        |
|----------------|-------|----------------|--------|
| X-WKR2436-00_0 | 24" L | X-WKR4836-00_0 | 48" L  |
| X-WKR3036-00_0 | 30" L | X-WKR5436-00_0 | 54" L* |
| X-WKR3636-00_0 | 36" L | X-WKR6036-00_0 | 60" L* |
| X-WKR4236-00_0 | 42" L | X-WKR7236-00_0 | 72" L* |

3/4" thick plywood panels with separate 4" high base rail, finished on face side only. Designed to be field-cut to exact height required. Furnished with two cleats.

**Note:** To specify panel to be finished on both sides, add an **-S** to the part number. Example: **X-WKR4236-0010-S**

\* Knee Space Panels over 48" long are furnished in two sections.

# Wall Mounted Storage Cabinets











## Specifications:


Open Wall Mounted Storage Cabinets are available in four heights, 24", 30", 36" and 48" and two depths; 12" and 16". Lengths are as shown. All are furnished with one, two or three adjustable shelves as shown.

**Note:** Last four digits of part number designate door and drawer head style, hinge, hardwood species, and pull selection. See page 4 for full explanation. Actual product may vary slightly from catalog illustration due to style and hardware selection.

## Open Storage Cabinets

<b>24" High</b>  	<b>12" Deep</b> <hr/> W00W241215-00_0 15"L W00W241218-00_0 18"L W00W241221-00_0 21"L W00W241224-00_0 24"L	<b>24" High</b>  	<b>12" Deep</b> <hr/> W00W241230-00_0 30"L W00W241236-00_0 36"L W00W241242-00_0 42"L W00W241248-00_0 48"L
	<b>16" Deep</b> <hr/> W00W241615-00_0 15"L W00W241618-00_0 18"L W00W241621-00_0 21"L W00W241624-00_0 24"L		<b>16" Deep</b> <hr/> W00W241630-00_0 30"L W00W241636-00_0 36"L W00W241642-00_0 42"L W00W241648-00_0 48"L
<b>30" High</b>  	<b>12" Deep</b> <hr/> W00W301215-00_0 15"L W00W301218-00_0 18"L W00W301221-00_0 21"L W00W301224-00_0 24"L	<b>30" High</b>  	<b>12" Deep</b> <hr/> W00W301230-00_0 30"L W00W301236-00_0 36"L W00W301242-00_0 42"L W00W301248-00_0 48"L
	<b>16" Deep</b> <hr/> W00W301615-00_0 15"L W00W301618-00_0 18"L W00W301621-00_0 21"L W00W301624-00_0 24"L		<b>16" Deep</b> <hr/> W00W301630-00_0 30"L W00W301636-00_0 36"L W00W301642-00_0 42"L W00W301648-00_0 48"L
<b>36" High</b>  	<b>12" Deep</b> <hr/> W00W361215-00_0 15"L W00W361218-00_0 18"L W00W361221-00_0 21"L W00W361224-00_0 24"L	<b>36" High</b>  	<b>12" Deep</b> <hr/> W00W361230-00_0 30"L W00W361236-00_0 36"L W00W361242-00_0 42"L W00W361248-00_0 48"L
	<b>16" Deep</b> <hr/> W00W361615-00_0 15"L W00W361618-00_0 18"L W00W361621-00_0 21"L W00W361624-00_0 24"L		<b>16" Deep</b> <hr/> W00W361630-00_0 30"L W00W361636-00_0 36"L W00W361642-00_0 42"L W00W361648-00_0 48"L
<b>48" High</b>  	<b>12" Deep</b> <hr/> W00W481215-00_0 15"L W00W481218-00_0 18"L W00W481221-00_0 21"L W00W481224-00_0 24"L	<b>48" High</b>  	<b>12" Deep</b> <hr/> W00W481230-00_0 30"L W00W481236-00_0 36"L W00W481242-00_0 42"L W00W481248-00_0 48"L
	<b>16" Deep</b> <hr/> W00W481615-00_0 15"L W00W481618-00_0 18"L W00W481621-00_0 21"L W00W481624-00_0 24"L		<b>16" Deep</b> <hr/> W00W481630-00_0 30"L W00W481636-00_0 36"L W00W481642-00_0 42"L W00W481648-00_0 48"L

# Wall Mounted Storage Cabinets




### Specifications:

Open Wall Mounted Cubbie Storage Cabinets are available in four heights, 24", 30", 36" and 48" and two depths; 12" and 16". Lengths are as shown. Shelves and vertical dividers are fixed to form multiple storage compartments as shown.

**Note:** Last four digits of part number designates door and drawer head style, hinge, hardwood species, and pull selection. See page 4 for full explanation. Actual product may vary slightly from catalog illustration due to style and hardware selection.

## Open Cubbie Storage Cabinets

### 24" High




#### 12" Deep

W01W241212-00_0	12"L
W01W241215-00_0	15"L
W01W241218-00_0	18"L

#### 16" Deep

W01W241612-00_0	12"L
W01W241615-00_0	15"L
W01W241618-00_0	18"L

### 24" High




#### 12" Deep

W01W241221-00_0	21"L
W01W241224-00_0	24"L
W01W241230-00_0	30"L

#### 16" Deep

W01W241621-00_0	21"L
W01W241624-00_0	24"L
W01W241630-00_0	30"L

### 30" High




#### 12" Deep

W01W301212-00_0	12"L
W01W301215-00_0	15"L
W01W301218-00_0	18"L

#### 16" Deep

W01W301612-00_0	12"L
W01W301615-00_0	15"L
W01W301618-00_0	18"L

### 30" High




#### 12" Deep

W01W301221-00_0	21"L
W01W301224-00_0	24"L
W01W301230-00_0	30"L

#### 16" Deep

W01W301621-00_0	21"L
W01W301624-00_0	24"L
W01W301630-00_0	30"L

### 36" High




#### 12" Deep

W01W361212-00_0	12"L
W01W361215-00_0	15"L
W01W361218-00_0	18"L

#### 16" Deep

W01W361612-00_0	12"L
W01W361615-00_0	15"L
W01W361618-00_0	18"L

### 36" High




#### 12" Deep

W01W361221-00_0	21"L
W01W361224-00_0	24"L
W01W361230-00_0	30"L

#### 16" Deep

W01W361621-00_0	21"L
W01W361624-00_0	24"L
W01W361630-00_0	30"L

### 48" High




#### 12" Deep

W01W481212-00_0	12"L
W01W481215-00_0	15"L
W01W481218-00_0	18"L

#### 16" Deep

W01W481612-00_0	12"L
W01W481615-00_0	15"L
W01W481618-00_0	18"L

### 48" High




#### 12" Deep

W01W481221-00_0	21"L
W01W481224-00_0	24"L
W01W481230-00_0	30"L

#### 16" Deep

W01W481621-00_0	21"L
W01W481624-00_0	24"L
W01W481630-00_0	30"L

# Wall Mounted Storage Cabinets




## Specifications:

Open Wall Mounted Cubbie Storage Cabinets are available in four heights, 24", 30", 36" and 48" and two depths; 12" and 16". Lengths are as shown. Shelves and vertical dividers are fixed to form multiple storage compartments as shown.

**Note:** Last four digits of part number designates door and drawer head style, hinge, hardwood species, and pull selection. See page 4 for full explanation. Actual product may vary slightly from catalog illustration due to style and hardware selection.

## Open Cubbie Storage Cabinets

### 24" High




#### 12" Deep

W01W241236-00_0	36"L
W01W241242-00_0	42"L
W01W241248-00_0	48"L

#### 16" Deep

W01W241636-00_0	36"L
W01W241642-00_0	42"L
W01W241648-00_0	48"L

### 30" High




#### 12" Deep

W01W301236-00_0	36"L
W01W301242-00_0	42"L
W01W301248-00_0	48"L

#### 16" Deep

W01W301636-00_0	36"L
W01W301642-00_0	42"L
W01W301648-00_0	48"L

### 36" High




#### 12" Deep

W01W361236-00_0	36"L
W01W361242-00_0	42"L
W01W361248-00_0	48"L

#### 16" Deep

W01W361636-00_0	36"L
W01W361642-00_0	42"L
W01W361648-00_0	48"L

### 48" High



#### 12" Deep


W01W481236-00_0	36"L
W01W481242-00_0	42"L
W01W481248-00_0	48"L

#### 16" Deep

W01W481636-00_0	36"L
W01W481642-00_0	42"L
W01W481648-00_0	48"L



# Wall Mounted Storage Cabinets




**Specifications:**

Plate Glass Door Wall Mounted Storage Cabinets are available in four heights, 24", 30", 36" and 48" and two depths; 12" and 16". Lengths are as shown. All are furnished with one, two or three adjustable shelves as shown. Doors are 1/4" float glass and are mounted on friction-free roller glides.

**Note:** Last four digits of part number designates door and drawer head style, hinge, hardwood species, and pull selection. See page 4 for full explanation. Actual product may vary slightly from catalog illustration due to style and hardware selection.

**Sliding Plate Glass Door Cabinets**                      **Double Sided — Pass-Thru Cabinets**

**24" High**




**12" Deep**

W10W241230-00_0	30"L
W10W241236-00_0	36"L
W10W241242-00_0	42"L
W10W241248-00_0	48"L

**16" Deep**

W10W241630-00_0	30"L
W10W241636-00_0	36"L
W10W241642-00_0	42"L
W10W241648-00_0	48"L

**30" High**




**12" Deep**

W10W301230-00_0	30"L
W10W301236-00_0	36"L
W10W301242-00_0	42"L
W10W301248-00_0	48"L

**16" Deep**

W10W301630-00_0	30"L
W10W301636-00_0	36"L
W10W301642-00_0	42"L
W10W301648-00_0	48"L

**30" High**




**12" Deep**

W13W301230-00_0	30"L
W13W301236-00_0	36"L
W13W301242-00_0	42"L
W13W301248-00_0	48"L

**16" Deep**

W13W301630-00_0	30"L
W13W301636-00_0	36"L
W13W301642-00_0	42"L
W13W301648-00_0	48"L

**36" High**




**12" Deep**

W10W361230-00_0	30"L
W10W361236-00_0	36"L
W10W361242-00_0	42"L
W10W361248-00_0	48"L

**16" Deep**

W10W361630-00_0	30"L
W10W361636-00_0	36"L
W10W361642-00_0	42"L
W10W361648-00_0	48"L

**36" High**




**12" Deep**

W13W361230-00_0	30"L
W13W361236-00_0	36"L
W13W361242-00_0	42"L
W13W361248-00_0	48"L

**16" Deep**

W13W361630-00_0	30"L
W13W361636-00_0	36"L
W13W361642-00_0	42"L
W13W361648-00_0	48"L

**48" High**




**12" Deep**

W10W481230-00_0	30"L
W10W481236-00_0	36"L
W10W481242-00_0	42"L
W10W481248-00_0	48"L

**16" Deep**

W10W481630-00_0	30"L
W10W481636-00_0	36"L
W10W481642-00_0	42"L
W10W481648-00_0	48"L

# Wall Mounted Storage Cabinets



## Specifications:


Sliding Door Wall Mounted Storage Cabinets are available in four heights, 24", 30", 36" and 48" and two depths; 12" and 16". Lengths are as shown. All are furnished with one, two or three adjustable shelves as shown. Sliding doors are suspended by rollers from an upper track. Glazed doors are fitted with

$\frac{1}{8}$ " float glass unless otherwise specified. **Note:** Last four digits of part number designates door and drawer head style, hinge, hardwood species, and pull selection. See page 4 for full explanation. Actual product may vary slightly from catalog illustration due to style and hardware selection.

### Sliding Glazed Door Cabinets

### Sliding Panel Door Cabinets

#### 24" High




#### 12" Deep

W11W241230-00_0	30"L
W11W241236-00_0	36"L
W11W241242-00_0	42"L
W11W241248-00_0	48"L

#### 16" Deep

W11W241630-00_0	30"L
W11W241636-00_0	36"L
W11W241642-00_0	42"L
W11W241648-00_0	48"L

#### 24" High




#### 12" Deep

W12W241230-00_0	30"L
W12W241236-00_0	36"L
W12W241242-00_0	42"L
W12W241248-00_0	48"L

#### 16" Deep

W12W241630-00_0	30"L
W12W241636-00_0	36"L
W12W241642-00_0	42"L
W12W241648-00_0	48"L

#### 30" High




#### 12" Deep

W11W301230-00_0	30"L
W11W301236-00_0	36"L
W11W301242-00_0	42"L
W11W301248-00_0	48"L

#### 16" Deep

W11W301630-00_0	30"L
W11W301636-00_0	36"L
W11W301642-00_0	42"L
W11W301648-00_0	48"L

#### 30" High




#### 12" Deep

W12W301230-00_0	30"L
W12W301236-00_0	36"L
W12W301242-00_0	42"L
W12W301248-00_0	48"L

#### 16" Deep

W12W301630-00_0	30"L
W12W301636-00_0	36"L
W12W301642-00_0	42"L
W12W301648-00_0	48"L

#### 36" High




#### 12" Deep

W11W361230-00_0	30"L
W11W361236-00_0	36"L
W11W361242-00_0	42"L
W11W361248-00_0	48"L

#### 16" Deep

W11W361630-00_0	30"L
W11W361636-00_0	36"L
W11W361642-00_0	42"L
W11W361648-00_0	48"L

#### 36" High




#### 12" Deep

W12W361230-00_0	30"L
W12W361236-00_0	36"L
W12W361242-00_0	42"L
W12W361248-00_0	48"L

#### 16" Deep

W12W361630-00_0	30"L
W12W361636-00_0	36"L
W12W361642-00_0	42"L
W12W361648-00_0	48"L

#### 48" High




#### 12" Deep

W11W481230-00_0	30"L
W11W481236-00_0	36"L
W11W481242-00_0	42"L
W11W481248-00_0	48"L

#### 16" Deep

W11W481630-00_0	30"L
W11W481636-00_0	36"L
W11W481642-00_0	42"L
W11W481648-00_0	48"L

#### 48" High




#### 12" Deep

W12W481230-00_0	30"L
W12W481236-00_0	36"L
W12W481242-00_0	42"L
W12W481248-00_0	48"L

#### 16" Deep

W12W481630-00_0	30"L
W12W481636-00_0	36"L
W12W481642-00_0	42"L
W12W481648-00_0	48"L

# Wall Mounted Storage Cabinets




**Specifications:**

Swinging Glazed Door Wall Mounted Storage Cabinets are available in four heights, 24", 30", 36" and 48" and two nominal depths; 12<sup>3</sup>/<sub>4</sub>" and 16<sup>3</sup>/<sub>4</sub>". Lengths are as shown. Cabinets are glazed with 1/8" float glass unless otherwise specified and are furnished with one, two or three adjustable shelves as shown.

**Note:** Last four digits of part number designates door style, hinge, hardwood species, and pull selection. See page 4 for full explanation. Actual product may vary slightly from catalog illustration due to style and hardware selection.


**Swinging Glazed Door Cabinets**

**24" High**




Left Hand	Right Hand	12 <sup>3</sup> / <sub>4</sub> " D
W20W241215L	W20W241215-	15"L
W20W241218L	W20W241218-	18"L
W20W241221L	W20W241221-	21"L
W20W241224L	W20W241224-	24"L
		<b>16<sup>3</sup>/<sub>4</sub>" D</b>
W20W241615L	W20W241615-	15"L
W20W241618L	W20W241618-	18"L
W20W241621L	W20W241621-	21"L
W20W241624L	W20W241624-	24"L

**24" High**




12 <sup>3</sup> / <sub>4</sub> " Deep
W20W241230- 30"L
W20W241236- 36"L
W20W241242- 42"L
W20W241248- 48"L
<b>16<sup>3</sup>/<sub>4</sub>" Deep</b>
W20W241630- 30"L
W20W241636- 36"L
W20W241642- 42"L
W20W241648- 48"L

**30" High**




Left Hand	Right Hand	12 <sup>3</sup> / <sub>4</sub> " D
W20W301215L	W20W301215-	15"L
W20W301218L	W20W301218-	18"L
W20W301221L	W20W301221-	21"L
W20W301224L	W20W301224-	24"L
		<b>16<sup>3</sup>/<sub>4</sub>" D</b>
W20W301615L	W20W301615-	15"L
W20W301618L	W20W301618-	18"L
W20W301621L	W20W301621-	21"L
W20W301624L	W20W301624-	24"L

**30" High**




12 <sup>3</sup> / <sub>4</sub> " Deep
W20W301230- 30"L
W20W301236- 36"L
W20W301242- 42"L
W20W301248- 48"L
<b>16<sup>3</sup>/<sub>4</sub>" Deep</b>
W20W301630- 30"L
W20W301636- 36"L
W20W301642- 42"L
W20W301648- 48"L

**36" High**




Left Hand	Right Hand	12 <sup>3</sup> / <sub>4</sub> " D
W20W361215L	W20W361215-	15"L
W20W361218L	W20W361218-	18"L
W20W361221L	W20W361221-	21"L
W20W361224L	W20W361224-	24"L
		<b>16<sup>3</sup>/<sub>4</sub>" D</b>
W20W361615L	W20W361615-	15"L
W20W361618L	W20W361618-	18"L
W20W361621L	W20W361621-	21"L
W20W361624L	W20W361624-	24"L

**36" High**




12 <sup>3</sup> / <sub>4</sub> " Deep
W20W361230- 30"L
W20W361236- 36"L
W20W361242- 42"L
W20W361248- 48"L
<b>16<sup>3</sup>/<sub>4</sub>" Deep</b>
W20W361630- 30"L
W20W361636- 36"L
W20W361642- 42"L
W20W361648- 48"L

**48" High**




Left Hand	Right Hand	12 <sup>3</sup> / <sub>4</sub> " D
W20W481215L	W20W481215-	15"L
W20W481218L	W20W481218-	18"L
W20W481221L	W20W481221-	21"L
W20W481224L	W20W481224-	24"L
		<b>16<sup>3</sup>/<sub>4</sub>" D</b>
W20W481615L	W20W481615-	15"L
W20W481618L	W20W481618-	18"L
W20W481621L	W20W481621-	21"L
W20W481624L	W20W481624-	24"L

**48" High**



12 <sup>3</sup> / <sub>4</sub> " Deep
W20W481230- 30"L
W20W481236- 36"L
W20W481242- 42"L
W20W481248- 48"L
<b>16<sup>3</sup>/<sub>4</sub>" Deep</b>
W20W481630- 30"L
W20W481636- 36"L
W20W481642- 42"L
W20W481648- 48"L

# Wall Mounted Storage Cabinets




## Specifications:

Swinging Panel Door Wall Mounted Storage Cabinets are available in four heights, 24", 30", 36" and 48" and two nominal depths; 12<sup>3</sup>/<sub>4</sub>" and 16<sup>3</sup>/<sub>4</sub>". Lengths are as shown. Cabinets are furnished with one, two or three adjustable shelves as shown.

**Note:** Last four digits of part number designates door style, hinge, hardwood species, and pull selection. See page 4 for full explanation. Actual product may vary slightly from catalog illustration due to style and hardware selection.


## Swinging Panel Door Cabinets

### 24" High




Left Hand	Right Hand	12 <sup>3</sup> / <sub>4</sub> " D
W25W241215L	W25W241215-	15"L
W25W241218L	W25W241218-	18"L
W25W241221L	W25W241221-	21"L
W25W241224L	W25W241224-	24"L
		16 <sup>3</sup> / <sub>4</sub> " D
W25W241615L	W25W241615-	15"L
W25W241618L	W25W241618-	18"L
W25W241621L	W25W241621-	21"L
W25W241624L	W25W241624-	24"L

### 24" High




12 <sup>3</sup> / <sub>4</sub> " Deep
W25W241230- 30"L
W25W241236- 36"L
W25W241242- 42"L
W25W241248- 48"L
16 <sup>3</sup> / <sub>4</sub> " Deep
W25W241630- 30"L
W25W241636- 36"L
W25W241642- 42"L
W25W241648- 48"L

### 30" High




Left Hand	Right Hand	12 <sup>3</sup> / <sub>4</sub> " D
W25W301215L	W25W301215-	15"L
W25W301218L	W25W301218-	18"L
W25W301221L	W25W301221-	21"L
W25W301224L	W25W301224-	24"L
		16 <sup>3</sup> / <sub>4</sub> " D
W25W301615L	W25W301615-	15"L
W25W301618L	W25W301618-	18"L
W25W301621L	W25W301621-	21"L
W25W301624L	W25W301624-	24"L

### 30" High




12 <sup>3</sup> / <sub>4</sub> " Deep
W25W301230- 30"L
W25W301236- 36"L
W25W301242- 42"L
W25W301248- 48"L
16 <sup>3</sup> / <sub>4</sub> " Deep
W25W301630- 30"L
W25W301636- 36"L
W25W301642- 42"L
W25W301648- 48"L

### 36" High




Left Hand	Right Hand	12 <sup>3</sup> / <sub>4</sub> " D
W25W361215L	W25W361215-	15"L
W25W361218L	W25W361218-	18"L
W25W361221L	W25W361221-	21"L
W25W361224L	W25W361224-	24"L
		16 <sup>3</sup> / <sub>4</sub> " D
W25W361615L	W25W361615-	15"L
W25W361618L	W25W361618-	18"L
W25W361621L	W25W361621-	21"L
W25W361624L	W25W361624-	24"L

### 36" High




12 <sup>3</sup> / <sub>4</sub> " Deep
W25W361230- 30"L
W25W361236- 36"L
W25W361242- 42"L
W25W361248- 48"L
16 <sup>3</sup> / <sub>4</sub> " Deep
W25W361630- 30"L
W25W361636- 36"L
W25W361642- 42"L
W25W361648- 48"L

### 48" High




Left Hand	Right Hand	12 <sup>3</sup> / <sub>4</sub> " D
W25W481215L	W25W481215-	15"L
W25W481218L	W25W481218-	18"L
W25W481221L	W25W481221-	21"L
W25W481224L	W25W481224-	24"L
		16 <sup>3</sup> / <sub>4</sub> " D
W25W481615L	W25W481615-	15"L
W25W481618L	W25W481618-	18"L
W25W481621L	W25W481621-	21"L
W25W481624L	W25W481624-	24"L

### 48" High



12 <sup>3</sup> / <sub>4</sub> " Deep
W25W481230- 30"L
W25W481236- 36"L
W25W481242- 42"L
W25W481248- 48"L
16 <sup>3</sup> / <sub>4</sub> " Deep
W25W481630- 30"L
W25W481636- 36"L
W25W481642- 42"L
W25W481648- 48"L

# Wall Mounted Storage Cabinets




**Specifications:**

Corner Wall Mounted Storage Cabinets are available in four heights, 24", 30", 36" and 48", two depths; 12" and 16", and with either glazed or panel swinging doors. Cabinets are furnished with one, two or three adjustable shelves as shown.


**Note:** Last four digits of part number designates door style, hinge, hardwood species, and pull selection. See page 4 for full explanation. Actual product may vary slightly from catalog illustration due to style and hardware selection.

**Swinging Glazed Door Corner Cabinets      Swinging Panel Door Corner Cabinets**


**24" High**




W21W241217L  
W21W241623L




W21W241217-  
W21W241623-



W26W241217L  
W26W241623L




W26W241217-  
W26W241623-




12<sup>3</sup>/<sub>4</sub>" D  
16<sup>3</sup>/<sub>4</sub>" D

PLAN VIEW OF  
12<sup>3</sup>/<sub>4</sub>" Deep  
Corner Cabinets


**30" High**




W21W301217L  
W21W301623L



W21W301217-  
W21W301623-




W26W301217L  
W26W301623L




W26W301217-  
W26W301623-

12<sup>3</sup>/<sub>4</sub>" D  
16<sup>3</sup>/<sub>4</sub>" D


**36" High**




W21W361217L  
W21W361623L




W21W361217-  
W21W361623-



W26W361217L  
W26W361623L




W26W361217-  
W26W361623-




12<sup>3</sup>/<sub>4</sub>" D  
16<sup>3</sup>/<sub>4</sub>" D

PLAN VIEW OF  
16<sup>3</sup>/<sub>4</sub>" Deep  
Corner Cabinets


**48" High**




W21W481217L  
W21W481623L



W21W481217-  
W21W481623-




W26W481217L  
W26W481623L



W26W481217-  
W26W481623-

12<sup>3</sup>/<sub>4</sub>" D  
16<sup>3</sup>/<sub>4</sub>" D

# Full Height Storage Cabinets



**Specifications:**


Full Height Storage Cabinets are 84" high. Depths and lengths as shown.  
Open Cubbie Cabinets have five fixed shelves creating a cabinet with six, twelve, or eighteen equal sized compartments as shown.

**Caution:**  
Full Height Storage Cabinets may present a potential tipping over problem. Proper installation requires

these cabinets to be anchored to a wall or other cabinets to improve stability and ensure smooth door movement.

**Note:** Last four digits of part number designates door style, hinge, hardwood species, and pull selection. See page 4 for full explanation. Actual product may vary slightly from catalog illustration due to style and hardware selection.

## Open Cubbie Storage



**12" Deep**


S01W841212-00_0	12"L
S01W841215-00_0	15"L
S01W841218-00_0	18"L

**16" Deep**

S01W841612-00_0	12"L
S01W841615-00_0	15"L
S01W841618-00_0	18"L

**21 1/4" Deep**

S01W842212-00_0	12"L
S01W842215-00_0	15"L
S01W842218-00_0	18"L



**12" Deep**


S01W841221-00_0	21"L
S01W841224-00_0	24"L
S01W841230-00_0	30"L

**16" Deep**

S01W841621-00_0	21"L
S01W841624-00_0	24"L
S01W841630-00_0	30"L

**21 1/4" Deep**

S01W842221-00_0	21"L
S01W842224-00_0	24"L
S01W842230-00_0	30"L



**12" Deep**

S01W841236-00_0	36"L
S01W841242-00_0	42"L
S01W841248-00_0	48"L


**16" Deep**

S01W841636-00_0	36"L
S01W841642-00_0	42"L
S01W841648-00_0	48"L

**21 1/4" Deep**

S01W842236-00_0	36"L
S01W842242-00_0	42"L
S01W842248-00_0	48"L

# Full Height Storage Cabinets



**Specifications:**

Full Height Storage Cabinets are 84" high. Depths and lengths as shown.

These cabinets have four adjustable shelves and one fixed shelf.

Sliding Glazed and Panel Doors are suspended by rollers from an aluminum track at the top for easy movement.


**Caution:**

Full Height Storage Cabinets may present a potential tipping over


problem. Proper installation requires these cabinets to be anchored to a wall or other cabinets to improve stability and ensure smooth door movement.

**Note:** Last four digits of part number designates door style, hinge, hardwood species, and pull selection. See page 4 for full explanation. Actual product may vary slightly from catalog illustration due to style and hardware selection.

**Open Storage**




<b>12" Deep</b>	
S00W841218-00_0	18"L
S00W841224-00_0	24"L
<b>16" Deep</b>	
S00W841618-00_0	18"L
S00W841624-00_0	24"L
<b>21 1/4" Deep</b>	
S00W842218-00_0	18"L
S00W842224-00_0	24"L




<b>12" Deep</b>	
S00W841230-00_0	30"L
S00W841236-00_0	36"L
S00W841242-00_0	42"L
S00W841248-00_0	48"L
<b>16" Deep</b>	
S00W841630-00_0	30"L
S00W841636-00_0	36"L
S00W841642-00_0	42"L
S00W841648-00_0	48"L
<b>21 1/4" Deep</b>	
S00W842230-00_0	30"L
S00W842236-00_0	36"L
S00W842242-00_0	42"L
S00W842248-00_0	48"L

**Sliding Glazed Door      Sliding Panel Door**




<b>12" Deep</b>	
S11W841230-00_0	30"L
S11W841236-00_0	36"L
S11W841242-00_0	42"L
S11W841248-00_0	48"L
<b>16" Deep</b>	
S11W841630-00_0	30"L
S11W841636-00_0	36"L
S11W841642-00_0	42"L
S11W841648-00_0	48"L
<b>21 1/4" Deep</b>	
S11W842230-00_0	30"L
S11W842236-00_0	36"L
S11W842242-00_0	42"L
S11W842248-00_0	48"L



<b>12" Deep</b>	
S12W841230-00_0	30"L
S12W841236-00_0	36"L
S12W841242-00_0	42"L
S12W841248-00_0	48"L
<b>16" Deep</b>	
S12W841630-00_0	30"L
S12W841636-00_0	36"L
S12W841642-00_0	42"L
S12W841648-00_0	48"L
<b>21 1/4" Deep</b>	
S12W842230-00_0	30"L
S12W842236-00_0	36"L
S12W842242-00_0	42"L
S12W842248-00_0	48"L

# Full Height Storage Cabinets



**Specifications:**

Full Height Storage Cabinets are 84" high. Depths and lengths as shown.


These cabinets have four adjustable shelves and one fixed shelf.

**Caution:**

Full Height Storage Cabinets may present a potential tipping over problem. Proper installation requires these cabinets to be anchored to a wall or other cabinets to improve stability and ensure smooth door movement.

**Note:** Last four digits of part number designates door style, hinge, hardwood species, and pull selection. See page 4 for full explanation. Actual product may vary slightly from catalog illustration due to style and hardware selection.

**Swinging Glazed Doors**



**12<sup>3</sup>/<sub>4</sub>" Deep**


S20W841218L	_____	18"L
S20W841224L	_____	24"L

**16<sup>3</sup>/<sub>4</sub>" Deep**

S20W841618L	_____	18"L
S20W841624L	_____	24"L

**22" Deep**

S20W842218L	_____	18"L
S20W842224L	_____	24"L



**12<sup>3</sup>/<sub>4</sub>" Deep**


S20W841218-	_____	18"L
S20W841224-	_____	24"L

**16<sup>3</sup>/<sub>4</sub>" Deep**

S20W841618-	_____	18"L
S20W841624-	_____	24"L

**22" Deep**

S20W842218-	_____	18"L
S20W842224-	_____	24"L



**12<sup>3</sup>/<sub>4</sub>" Deep**

S20W841230-	_____	30"L
S20W841236-	_____	36"L
S20W841242-	_____	42"L
S20W841248-	_____	48"L

**16<sup>3</sup>/<sub>4</sub>" Deep**


S20W841630-	_____	30"L
S20W841636-	_____	36"L
S20W841642-	_____	42"L
S20W841648-	_____	48"L

**22" Deep**

S20W842230-	_____	30"L
S20W842236-	_____	36"L
S20W842242-	_____	42"L
S20W842248-	_____	48"L



# Full Height Storage Cabinets



**Specifications:**

Full Height Storage Cabinets are 84" high. Depths and lengths as shown.


These cabinets have four adjustable shelves and one fixed shelf.

**Caution:**

Full Height Storage Cabinets may present a potential tipping over problem. Proper installation requires these cabinets to be anchored to a wall or other cabinets to improve stability and ensure smooth door movement.

**Note:** Last four digits of part number designates door style, hinge, hardwood species, and pull selection. See page 4 for full explanation. Actual product may vary slightly from catalog illustration due to style and hardware selection.

**Swinging Panel Doors**



**12<sup>3</sup>/<sub>4</sub>" Deep**

S25W841218L	18"L
S25W841224L	24"L

**16<sup>3</sup>/<sub>4</sub>" Deep**

S25W841618L	18"L
S25W841624L	24"L

**22" Deep**

S25W842218L	18"L
S25W842224L	24"L



**12<sup>3</sup>/<sub>4</sub>" Deep**

S25W841218-	18"L
S25W841224-	24"L


**16<sup>3</sup>/<sub>4</sub>" Deep**

S25W841618-	18"L
S25W841624-	24"L

**22" Deep**

S25W842218-	18"L
S25W842224-	24"L

**Swinging Panel Doors      Display Cabinet**



**12<sup>3</sup>/<sub>4</sub>" Deep**


S25W841230-	30"L
S25W841236-	36"L
S25W841242-	42"L
S25W841248-	48"L

**16<sup>3</sup>/<sub>4</sub>" Deep**

S25W841630-	30"L
S25W841636-	36"L
S25W841642-	42"L
S25W841648-	48"L

**22" Deep**

S25W842230-	30"L
S25W842236-	36"L
S25W842242-	42"L
S25W842248-	48"L



**16" Deep**


S72W841648-_0_0	48"L
-----------------	------

**21<sup>1</sup>/<sub>4</sub>" Deep**

S72W842248-_0_0	48"L
-----------------	------

Two adjustable plate glass shelves are provided in the upper portion behind plate glass sliding doors. One adjustable shelf is provided in the lower portion behind sliding panel doors. A fluorescent light fixture (bulbs not included) with reflector and switch wired to a plug is furnished at the top of the cabinet.

# Full Height Storage Cabinets



**Specifications:**

Full Height Storage Cabinets are 84" high. Depths and lengths as shown.

Sliding Door Tote Tray Storage Cabinets have sliding panel doors suspended by rollers from an aluminum track at the top for easy movement.


**Caution:**

Full Height Storage Cabinets may present a potential tipping over problem. Proper installation requires these cabinets to be anchored to

a wall or other cabinets to improve stability and ensure smooth door movement.

**Note:** Last four digits of part number designates door style, hinge, hardwood species, and pull selection. See page 4 for full explanation. Actual product may vary slightly from catalog illustration due to style and hardware selection.

**Split Swinging Glazed Doors**



**12<sup>3</sup>/<sub>4</sub>" Deep**

<b>S30W841230-</b> _____	30"L
<b>S30W841236-</b> _____	36"L
<b>S30W841242-</b> _____	42"L
<b>S30W841248-</b> _____	48"L


**16<sup>3</sup>/<sub>4</sub>" Deep**

<b>S30W841630-</b> _____	30"L
<b>S30W841636-</b> _____	36"L
<b>S30W841642-</b> _____	42"L
<b>S30W841648-</b> _____	48"L

**22" Deep**

<b>S30W842230-</b> _____	30"L
<b>S30W842236-</b> _____	36"L
<b>S30W842242-</b> _____	42"L
<b>S30W842248-</b> _____	48"L

**Split Swinging Panel Doors**



**12<sup>3</sup>/<sub>4</sub>" Deep**

<b>S35W841230-</b> _____	30"L
<b>S35W841236-</b> _____	36"L
<b>S35W841242-</b> _____	42"L
<b>S35W841248-</b> _____	48"L


**16<sup>3</sup>/<sub>4</sub>" Deep**

<b>S35W841630-</b> _____	30"L
<b>S35W841636-</b> _____	36"L
<b>S35W841642-</b> _____	42"L
<b>S35W841648-</b> _____	48"L

**22" Deep**

<b>S35W842230-</b> _____	30"L
<b>S35W842236-</b> _____	36"L
<b>S35W842242-</b> _____	42"L
<b>S35W842248-</b> _____	48"L

**Sliding Door Tote Tray Storage**




**22" Deep**

<b>S62W842236-</b> _____	36"L
<b>S62W842248-</b> _____	48"L

Provided with one fixed and 10 adjustable shelves. Consecutive number plates furnished for shelves and case bottom. The 36" cabinet holds (36) OP-6141 tote trays. The 48" cabinet holds (48). Order tote trays separately.

**Swinging Door Tote Tray Storage**




**22" Deep**

<b>S64W842236-</b> _____	36"L
<b>S64W842248-</b> _____	48"L

Provided with one fixed and 10 adjustable shelves. Consecutive number plates furnished for shelves and case bottom. The 36" cabinet holds (36) OP-6141 tote trays. The 48" cabinet holds (48). Order tote trays separately.

# Full Height Storage Cabinets



**Specifications:**

Full Height Storage Cabinets are 84" high. Depths and lengths as shown.

Microscope Storage Cabinets have five adjustable shelves and 10, 15, or 25 vertical partitions to create the number of storage compartments indicated.


**Caution:**

Full Height Storage Cabinets may present a potential tipping over problem. Proper installation requires these cabinets to be anchored to

a wall or other cabinets to improve stability and ensure smooth door movement.

**Note:** Last four digits of part number designates door style, hinge, hardwood species, and pull selection. See page 4 for full explanation. Actual product may vary slightly from catalog illustration due to style and hardware selection.

**Microscope Storage Cabinets — Swinging Glazed Doors with Lock**




16<sup>3</sup>/<sub>4</sub>" Deep

S21W841624L \_\_\_\_\_ 24"L

22" Deep

S21W842224L \_\_\_\_\_ 24"L

Provides 15 storage compartments




16<sup>3</sup>/<sub>4</sub>" Deep

S21W841624- \_\_\_\_\_ 24"L

22" Deep

S21W842224- \_\_\_\_\_ 24"L

Provides 15 storage compartments




16<sup>3</sup>/<sub>4</sub>" Deep

S21W841636- \_\_\_\_\_ 36"L

22" Deep

S21W842236- \_\_\_\_\_ 36"L

Provides 20 storage compartments.



16<sup>3</sup>/<sub>4</sub>" Deep

S21W841648- \_\_\_\_\_ 48"L


22" Deep

S21W842248- \_\_\_\_\_ 48"L

Provides 30 storage compartments.

All Microscope Storage Cabinet shelves and vertical partitions are adjustable.

**Instructor Storage Cabinet**




22" Deep

S41W842236L \_\_\_\_\_ 36"L

S41W842248L \_\_\_\_\_ 48"L

Provided with a fixed full length upper storage shelf, three fixed half length storage shelves, an extension slide garment hanger, two hanging file drawers, and a tote tray. 36" cabinet, OP-6140 tray 48" cabinet, OP-661 tray




22" Deep

S41W842236- \_\_\_\_\_ 36"L

S41W842248- \_\_\_\_\_ 48"L

Provided with a fixed full length upper storage shelf, three fixed half length storage shelves, an extension slide garment hanger, two hanging file drawers, and a tote tray. 36" cabinet, OP-6140 tray 48" cabinet, OP-661 tray

# Full Height Storage Cabinets



**Specifications:**

Full Height Storage Cabinets are 84" high. Depths and lengths as shown.

Acid Storage Cabinets have a polyethylene corrosion resistant lining and five full-width corrosion resistant shelves. Ventilation provided by louvers at door bottom and a 3 1/2" O.D. vent outlet at the top of the case.

**Caution:**





Full Height Storage Cabinets may present a potential tipping over problem. Proper installation requires

these cabinets to be anchored to a wall or other cabinets to improve stability and ensure smooth door movement.

**Note:** Last four digits of part number designates door style, hinge, hardwood species, and pull selection. See page 4 for full explanation. Actual product may vary slightly from catalog illustration due to style and hardware selection.

**Garment Locker Cabinets – 16 3/4" Deep**




**Garment Locker Cabinets – 22" Deep**



**S40W841624L** \_\_\_\_\_ 24"L    **S40W841624-** \_\_\_\_\_ 24"L  
Provided with a hat shelf and an extension slide garment hanger.

**S40W842236L** \_\_\_\_\_ 36"L    **S40W842236-** \_\_\_\_\_ 36"L  
**S40W842248L** \_\_\_\_\_ 48"L    **S40W842248-** \_\_\_\_\_ 48"L  
Provided with a fixed full length upper storage shelf, three fixed half length storage shelves, and an extension slide garment hanger.

**Acid Storage Cabinets**



**S70W841624L** \_\_\_\_\_ 16 3/4"D x 24"L

**S70W841624-** \_\_\_\_\_ 16 3/4"D x 24"L

**S70W841636-** \_\_\_\_\_ 16 3/4"D x 36"L

**S70W841648-** \_\_\_\_\_ 16 3/4"D x 48"L


**S70W842224L** \_\_\_\_\_ 22"D x 24"L

**S70W842224-** \_\_\_\_\_ 22"D x 24"L

**S70W842236-** \_\_\_\_\_ 22"D x 36"L

**S70W842248-** \_\_\_\_\_ 22"D x 48"L

# Full Height Storage Cabinets



**Specifications:**

Full Height Storage Cabinets are 84" high. Depths and lengths as shown.





**Caution:**

Full Height Storage Cabinets may present a potential tipping over problem. Proper installation requires these cabinets to be anchored to a wall or other cabinets to improve stability and ensure smooth door movement.

**Note:** Last four digits of part number designates door style, hinge, hardwood species, and pull selection. See page 4 for full explanation. Actual product may vary slightly from catalog illustration due to style and hardware selection.

**Skeleton Storage Cabinet – 24" Long**

**Torso Storage Cabinet – 24" Long**



S61W841624L \_\_\_\_\_ 16<sup>3</sup>/<sub>4</sub>"D  
S61W842224L \_\_\_\_\_ 22"D

S61W841624- \_\_\_\_\_ 16<sup>3</sup>/<sub>4</sub>"D  
S61W842224- \_\_\_\_\_ 22"D

S63W841624L \_\_\_\_\_ 16<sup>3</sup>/<sub>4</sub>"D  
S63W842224L \_\_\_\_\_ 22"D





S63W841624- \_\_\_\_\_ 16<sup>3</sup>/<sub>4</sub>"D  
S63W842224- \_\_\_\_\_ 22"D

Provided with an extension slide carrier.

Provided with one fixed full extension pull out shelf, and two full depth adjustable shelves.

**Glass Tubing Storage Cabinet – 24" Long**

**14 Drawer Storage Cabinet – 22" Deep**



S60W841624L \_\_\_\_\_ 16<sup>3</sup>/<sub>4</sub>"D

S60W841624- \_\_\_\_\_ 16<sup>3</sup>/<sub>4</sub>"D

S85W842236- \_\_\_\_\_ 36"L

S85W842248- \_\_\_\_\_ 48"L


Provided with fixed slanting compartments for storage of long lengths of glass tubing.

Provided with 14 box drawers: 14"L x 18<sup>1</sup>/<sub>8</sub>"D x 7"H I.D.

Provided with 14 box drawers: 20"L x 18<sup>1</sup>/<sub>8</sub>"D x 7"H I.D.

# Fillers, Soffits & Panels

## Face Fillers



### 36" High (H)

<b>X-WP0236-00_0</b>	2" W
<b>X-WP0436-00_0</b>	4" W
<b>X-WP0636-00_0</b>	6" W
<b>X-WP1236-00_0</b>	12" W


### 48" High (H)

<b>X-WP0248-00_0</b>	2" W
<b>X-WP0448-00_0</b>	4" W
<b>X-WP0648-00_0</b>	6" W
<b>X-WP1248-00_0</b>	12" W

### 84" High (H)


<b>X-WP0284-00_0</b>	2" W
<b>X-WP0484-00_0</b>	4" W
<b>X-WP0684-00_0</b>	6" W
<b>X-WP1284-00_0</b>	12" W

## Toe Space Fillers



<b>X-WTP412-0000</b>	12" W
<b>X-WTP424-0000</b>	24" W
<b>X-WTP436-0000</b>	36" W
<b>X-WTP448-0000</b>	48" W


## Rear Fillers



### 36" High (H)

<b>X-WP0136-00_0</b>	1" W
<b>X-WP0736-00_0</b>	7" W
<b>X-WP0836-00_0</b>	8" W
<b>X-WP1436-00_0</b>	14" W
<b>X-WP1636-00_0</b>	16" W

## Corner Fillers




### 36" High (H)

<b>X-WC0236-00_0</b>	2" W
<b>X-WC0436-00_0</b>	4" W
<b>X-WC0636-00_0</b>	6" W

All fillers designed to be field-cut to required length and fit. All fabricated from 3/4" thick hardwood plywood with edges cut square for flush installation. Fillers finished face side only, rear coated with clear sealer. Toe Space Filler finished in clear sealer only.

**Note:** To specify panel to be finished on both sides, add an **S** to the part number.  
Example: **X-WP0236-0010S**

## Finished Panels Standard Grain



### 36"W x 48"H Panel

<b>X-WS3648-00_0</b>	1/4" Thick
<b>X-WP3648-00_0</b>	3/4" Thick


### 24"W x 96"H Panel

<b>X-WS2496-00_0</b>	1/4" Thick
<b>X-WP2496-00_0</b>	3/4" Thick

### 48"W x 96"H Panel

<b>X-WS4896-00_0</b>	1/4" Thick
<b>X-WP4896-00_0</b>	3/4" Thick


## Finished Panels Cross Grain



### 48"W x 36"H Panel

<b>X-WS4836-00_0</b>	1/4" Thick
<b>X-WP4836-00_0</b>	3/4" Thick

## Corner Soffits



### for 12" Deep Cabinets

<b>X-WSF212-00_0</b>	2" W
<b>X-WSF412-00_0</b>	4" W
<b>X-WSF612-00_0</b>	6" W

### for 16" Deep Cabinets

<b>X-WSF216-00_0</b>	2" W
<b>X-WSF416-00_0</b>	4" W
<b>X-WSF616-00_0</b>	6" W

Panels can be field-cut to size for fabricating finished backs, soffits, miscellaneous fillers, etc. Finished front face only, rear coated with clear sealer.


**Note:** To specify panel to be finished on both sides, add an **S** to the part number. Example: **X-WP0236-0010S**

1/4" plywood panels used to close opening between two adjacent wall cabinets on opposite legs of L-shaped furniture assemblies. Order two to close off top and bottom

# Fillers, Soffits & Panels


## Base Cabinet Fillers with Flush Panels ( for use with Style 5 Cabinets)

### Left Hand Fillers




Sitting Height	30" High (H)
X-WBP0431L50_0	4" W
X-WBP0631L50_0	6" W
X-WBP1231L50_0	12" W
ADA Height	32 1/2" High (H)
X-WBP0434L50_0	4" W
X-WBP0634L50_0	6" W
X-WBP1234L50_0	12" W
Standing Height	35" High (H)
X-WBP0436L50_0	4" W
X-WBP0636L50_0	6" W
X-WBP1236L50_0	12" W

### Right Hand Fillers




Sitting Height	30" High (H)
X-WBP0431-50_0	4" W
X-WBP0631-50_0	6" W
X-WBP1231-50_0	12" W
ADA Height	32 1/2" High (H)
X-WBP0434-50_0	4" W
X-WBP0634-50_0	6" W
X-WBP1234-50_0	12" W
Standing Height	35" High (H)
X-WBP0436-50_0	4" W
X-WBP0636-50_0	6" W
X-WBP1236-50_0	12" W

### Intermediate Fillers



Sitting Height	30" High (H)
X-WBP0431I50_0	4" W
X-WBP0631I50_0	6" W
X-WBP1231I50_0	12" W
ADA Height	32 1/2" High (H)
X-WBP0434I50_0	4" W
X-WBP0634I50_0	6" W
X-WBP1234I50_0	12" W
Standing Height	35" High (H)
X-WBP0436I50_0	4" W
X-WBP0636I50_0	6" W
X-WBP1236I50_0	12" W


### Corner Fillers



Sitting Height	30" High (H)
X-WBC0431-50_0	4" W
X-WBC0631-50_0	6" W
ADA Height	32 1/2" High (H)
X-WBC0434-50_0	4" W
X-WBC0634-50_0	6" W
Standing Height	35" High (H)
X-WBC0436-50_0	4" W
X-WBC0636-50_0	6" W


## Wall & Full Height Cabinet Fillers with Flush Panels ( for use with Style 5 Cabinets)

### Left Hand Fillers




24" High (H)	
X-WWP0424L50_0	4" W
X-WWP0624L50_0	6" W
X-WWP1224L50_0	12" W
30" High (H)	
X-WWP0430L50_0	4" W
X-WWP0630L50_0	6" W
X-WWP1230L50_0	12" W
36" High (H)	
X-WWP0436L50_0	4" W
X-WWP0636L50_0	6" W
X-WWP1236L50_0	12" W
48" High (H)	
X-WWP0448L50_0	4" W
X-WWP0648L50_0	6" W
X-WWP1248L50_0	12" W

### Right Hand Fillers




24" High (H)	
X-WWP0424-50_0	4" W
X-WWP0624-50_0	6" W
X-WWP1224-50_0	12" W
30" High (H)	
X-WWP0430-50_0	4" W
X-WWP0630-50_0	6" W
X-WWP1230-50_0	12" W
36" High (H)	
X-WWP0436-50_0	4" W
X-WWP0636-50_0	6" W
X-WWP1236-50_0	12" W
48" High (H)	
X-WWP0448-50_0	4" W
X-WWP0648-50_0	6" W
X-WWP1248-50_0	12" W

### Intermediate Fillers



24" High (H)	
X-WWP0424I50_0	4" W
X-WWP0624I50_0	6" W
X-WWP1224I50_0	12" W
30" High (H)	
X-WWP0430I50_0	4" W
X-WWP0630I50_0	6" W
X-WWP1230I50_0	12" W
36" High (H)	
X-WWP0436I50_0	4" W
X-WWP0636I50_0	6" W
X-WWP1236I50_0	12" W
48" High (H)	
X-WWP0448I50_0	4" W
X-WWP0648I50_0	6" W
X-WWP1248I50_0	12" W

### Full Height Fillers



Left Hand	84" High (H)
X-WFP0484L50_0	4" W
X-WFP0684L50_0	6" W
X-WFP1284L50_0	12" W
Right Hand	84" High (H)
X-WFP0484-50_0	4" W
X-WFP0684-50_0	6" W
X-WFP1284-50_0	12" W
Intermediate	84" High (H)
X-WFP0484I50_0	4" W
X-WFP0684I50_0	6" W
X-WFP1284I50_0	12" W

# Educational Assemblies

To further speed your progress, we offer a Wood Assemblies Program. It's a special collection of furniture, tops and accessories that have proven popular with clients whose needs may be similar to your own.

Physics labs. Biology labs. We have many different configurations ready for you to use.

Each configuration is based on thorough research. We've looked at literally hundreds of actual labs in public schools, colleges and universities. And talked with numerous designers and architects throughout the country. The result is our assemblies give you highly-detailed, real-world information.

For instance, if you're building a school science lab, you can choose from an assortment of tables from instructor's desk to multiple work tables. Each with different configurations of sinks, cabinets, outlets and gas fittings.

Before you plan your next lab, ask about our Wood Assemblies Program.





# Educational Assemblies



# Kewaunee Matrix Instructor Station

## Specifications


**Kewaunee Matrix Instructor Station** provides the space required to conduct effective table top demonstrations. Stations are available in either left or right hand configurations. The Cabinet Assembly includes a cupboard, drawer, and multi-media storage, finished back, and filler panels.

The Work Top Assembly includes a 1" thick Kemresin® worksurface with a drop-in sink and sink outlet.

The Fixture Package includes: fittings, fixtures, rod sockets, upright rod assembly, keyboard tray, monitor arm, and CPU holder.

**Note:** Last four digits of part number designates door style, hinge, hardwood species, and pull selection. See page 4 of Signature Series Wood catalog for full explanation. Actual product may vary slightly from catalog illustration due to style and hardware selection.

## Kewaunee Matrix Instructor Station



**X23 (Back View)**

**X23 Instructors Station** 36" High x 30" Deep x 140" Long (Sink on Right)

**X33 Instructors Station** 36" High x 30" Deep x 140" Long (Sink on Left)

**Cabinet Assembly** (Sink on Right)

**Cabinet Assembly** (Sink on Left)

**Cabinets:** **G40W312838**  
**G77W362848**  
**G83W362848**

**X23WC363000-**\_\_\_\_\_ *specify style*

**X33WC363000-**\_\_\_\_\_ *specify style*

Desk Assembly  
Sink/Multi-media Storage Cabinet  
Drawer Cupboard Cabinet

**Work Top Assembly** (Sink on Right)

**Work Top Assembly**(Sink on Left)

**Work Top:** 1" thick Kemresin - 30" deep x 140" long  
**Sink:** (1) **1003-DI-**\_\_\_\_ Drop-in Kemresin® - 11" D x 15" W x 18" L  
(1) **0482-**\_\_\_\_ 1 1/2" I.P.S. Sink Outlet

**X23WT013000R-**\_\_\_ *specify color*

**X33WT013000R-**\_\_\_ *specify color*

**Fixture Package** (Sink on Right)


**Fixture Package** (Sink on Left)

**Electrical Fixtures:** (3) **S656-1**\_\_\_\_ 120 VAC GFI Duplex Receptacle  
(1) **0720-1**\_\_\_\_ Duplex RJ45 Data Jack  
(1) **S721-6**\_\_\_\_ 120 VAC Duplex/Duplex RJ45 Data  
**Fittings:** (1) **W-0340-0V** Water Mixing Faucet w/VB Gooseneck  
(1) **W-0263-0G** 90° Double Outlet Gas Fitting  
(2) **0548-00** Rod Sockets  
(1) **4740-00** Upright Rod Assembly  
(1) **F72M1220-00** Keyboard Tray  
(1) **F72M0808-0T** Mouse Platform  
(1) **F-5435-00** Monitor Arm  
(1) **F-5436-00** CPU Holder


**X23FP-**\_\_\_,\_\_\_ *specify receptacle color, service*

**X33FP-**\_\_\_,\_\_\_ *specify receptacle color, service*


# Kewaunee Matrix Instructor Station



**X33** (Back View)



**X23** (Front View)



**X23** (End View)

**For Example Classroom Layouts see pages 147-149**

# Kewaunee Matrix Student Stations

## Specifications

**Kewaunee Matrix Power Adjustable Height Student Stations** provide a movable, adaptable teaching space on a steel tubular frame with hardwood veneer plywood modesty panels. The stations are adjustable in height from 31" to 38" by means of an electronically controlled lift


mechanism with up and down push button controls and three programmable memory location buttons.

Stations are available on casters or glides and are furnished with a six foot, 3-wire, 20 amp pigtail.

Student Station Work Tops are available in Kemresin, Black Phenolic Resin, and Black Plastic Laminate.

An optional six outlet power strip is available to provide power to the student work space.

## Kewaunee Matrix Power Adjustable Height Student Station



### X-020018 Optional Power Strip

120 VAC, 20 amp, six outlet power strip with six foot 12/3 cord. May be mounted to either side of the table frame to provide power for student work areas.

### X91 Power Adjustable Height Student Station 36" High x 48" Deep x 101" Long

Student Station Tubular Frame with Glides

**X91WC364800-00\_0** *specify wood species*

Student Station Tubular Frame with Casters

**X91WC364800-00\_0-B9** *specify wood species*

**Support Frame:** Steel Tubular Frame - Hardwood Veneer Plywood Modesty Panels  
**Work Top Height:** Adjustable from 31" to 38"  
**Height Adjustment Type:** Electric Motor with Three (3) position Memory Switch  
**Dynamic Load Capacity:** 1200 pounds (evenly distributed)

Student Station Work Top - Kemresin

**X91WT014800R-\_\_\_** *specify color*


Student Station Work Top - Black Phenolic Resin

**X91WT014800T**

Student Station Work Top - Black Plastic Laminate

**X91WT014800P**

**Work Top Size:** 1" thick x 48" deep x 101" long



For Example Classroom Layouts see pages 141-143

# Kewaunee Matrix Student Stations

## Specifications

**Kewaunee Matrix Manual Adjustable Height Student Stations** provide a movable, adaptable teaching space on a steel tubular frame with hardwood veneer plywood modesty panels. The stations are adjustable in height from 31" to 38" in 1" increments.


Stations are available on casters or glides.

Student Station Work Tops are available in Kemresin, Black Phenolic Resin, and Black Plastic Laminate.

An optional six outlet power strip is available to provide power to the student work space.

## Kewaunee Matrix Manual Adjustable Height Student Station

**X-020018 Optional Power Strip**  
120 VAC, 20 amp, six outlet power strip with six foot 12/3 cord. May be mounted to either side of the table frame to provide power for student work areas.



### X90 Manual Adjustable Height Student Station 36" High x 48" Deep x 101" Long

Student Station Tubular Frame with Glides  
Student Station Tubular Frame with Casters

**X90WC364800-00\_0** *specify wood species*

**X90WC364800-00\_0-B9** *specify wood species*

**Support Frame:** Steel Tubular Frame - Hardwood Veneer Plywood Modesty Panels  
**Work Top Height:** Adjustable from 31" to 38"  
**Height Adjustment Type:** Bolt-thru on 1" centers  
**Dynamic Load Capacity:** 1200 pounds (evenly distributed)

Student Station Work Top - Kemresin

**X90WT014800R-\_\_\_** *specify color*


Student Station Work Top - Black Phenolic Resin

**X90WT014800T**

Student Station Work Top - Black Plastic Laminate

**X90WT014800P**

**Work Top Size:** 1" thick x 48" deep x 101" long



For Example Classroom Layouts see pages 141-143

# Kewaunee Matrix Pedestals

## Specifications


**Kewaunee Matrix Sink Pedestals** are designed to be used in conjunction with the Kewaunee Matrix Student Stations. The three 27" long edges match the 27" long edges of the student station allowing multiple configuration options.

Sink Pedestals are available in three standard heights, with or without, three 120 VAC GFI Duplex Receptacles.

Work Tops are available in Kemresin® with a drop-in sink and sink outlet.

The standard fixture package includes a cold water faucet with vacuum breaker gooseneck.

## Kewaunee Matrix Sink Pedestal



### X94 Sink Pedestal with Electric Outlets – 22" Deep x 27" Long

#### Pedestal Base

Standing Height	35" High	<b>X94WC362727-00_0</b> <i>specify wood species</i>
ADA Height	32 <sup>1</sup> / <sub>2</sub> " High	<b>X94WC342727-00_0</b> <i>specify wood species</i>
Sitting Height	30" High	<b>X94WC312727-00_0</b> <i>specify wood species</i>

#### Pedestal Work Top – Kemresin

**X94WT362727R-**\_\_\_ *specify color*

**Work Top:** 1" thick x 22" deep x 27" long

**Sink:** (1) **1048-DI-BK** Drop-in Kemresin - 6<sup>1</sup>/<sub>2</sub>"H 12<sup>3</sup>/<sub>8</sub>"W x 15<sup>3</sup>/<sub>4</sub>"L  
(1) **0482-BP** 1<sup>1</sup>/<sub>2</sub>" I.P.S. Sink Outlet

#### Pedestal Fixture Package

**X94FP-**\_\_\_ *specify receptacle color*

**Electrical Fixtures:** (1) **W-0337-0V** Cold Water Faucet with VB Gooseneck  
**Fittings:** (3) **S656-1\_** 120 VAC GFI Duplex Receptacle

### X92 Sink Pedestal without Electric Outlets – 22" Deep x 27" Long

#### Pedestal Base

Standing Height	35" High	<b>X92WC362727-00_0</b> <i>specify wood species</i>
ADA Height	32 <sup>1</sup> / <sub>2</sub> " High	<b>X92WC342727-00_0</b> <i>specify wood species</i>
Sitting Height	30" High	<b>X92WC312727-00_0</b> <i>specify wood species</i>

#### Pedestal Work Top – Kemresin

**X92WT362727R-**\_\_\_ *specify color*

**Work Top:** 1" thick x 22" deep x 27" long

**Sink:** (1) **1048-DI-BK** Drop-in Kemresin - 6<sup>1</sup>/<sub>2</sub>"H 12<sup>3</sup>/<sub>8</sub>"W x 15<sup>3</sup>/<sub>4</sub>"L  
(1) **0482-BP** 1<sup>1</sup>/<sub>2</sub>" I.P.S. Sink Outlet


#### Pedestal Fixture Package

**X92FP**


**Electrical Fixtures:** (1) **W-0337-0V** Cold Water Faucet with VB Gooseneck

## Example Classroom Layouts


Layout 1



Layout 2



Layout 3



Drawings based on 28' x 34' classroom

For more Example Classroom Layouts see pages 141-143

# Kewaunee Matrix Pedestals

## Specifications

**Kewaunee Matrix Fixed Pedestals** are designed to be used in conjunction with the Kewaunee Matrix Student Stations. The three 27" long edges match the 27" long edges of the student station allowing multiple configuration options.

Fixed Pedestals are available in three standard heights, with or without, three 120 VAC GFI Duplex Receptacles.

Work Tops are available in Kemresin, Black Phenolic Resin, and Black Plastic Laminate.

**Kewaunee Matrix Adjustable Height Pedestals** use the same size work top as the Fixed Pedestal, but is mounted on an adjustable height, three-legged steel base. The pedestal is adjustable in height from 31" to 38".

### Kewaunee Matrix Fixed Pedestal

#### X95 Fixed Pedestal with Electric Outlets – 22" Deep x 27" Long

##### Pedestal Base

Standing Height	35" High	<b>X95WC362727-00_0</b> <i>specify wood species</i>
ADA Height	32 1/2" High	<b>X95WC342727-00_0</b> <i>specify wood species</i>
Sitting Height	30" High	<b>X95WC312727-00_0</b> <i>specify wood species</i>

##### Pedestal Work Top


Kemresin	<b>X95WT012727R-___</b> <i>specify color</i>
Black Phenolic Resin	<b>X95WT012727T</b>
Black Plastic Laminate	<b>X95WT012727P</b>

**Work Top Size:** 1" thick x 22" deep x 27" long

##### Pedestal Fixture Package

**X95FP-\_\_** *specify receptacle color*

**Fittings:** (3) **S656-1\_** 120 VAC GFI Duplex Receptacle



#### X93 Fixed Pedestal without Electric Outlets – 22" Deep x 27" Long

##### Pedestal Base


Standing Height	35" High	<b>X93WC362727-00_0</b> <i>specify wood species</i>
ADA Height	32 1/2" High	<b>X93WC342727-00_0</b> <i>specify wood species</i>
Sitting Height	30" High	<b>X93WC312727-00_0</b> <i>specify wood species</i>

##### Pedestal Work Top

Kemresin	<b>X93WT012727R-___</b> <i>specify color</i>
Black Phenolic Resin	<b>X93WT012727T</b>
Black Plastic Laminate	<b>X93WT012727P</b>

**Work Top Size:** 1" thick x 22" deep x 27" long

### Kewaunee Matrix Adjustable Height Pedestal



#### X96 Adjustable Height Pedestal – 22" Deep x 27" Long

**Pedestal Base** 31-38" High **X96WC362727**

##### Pedestal Work Top

Kemresin	<b>X96WT362727R-___</b> <i>specify color</i>
Black Phenolic Resin	<b>X96WT362727T</b>
Black Plastic Laminate	<b>X96WT362727P</b>

**Work To Size:** 1" thick x 22" deep x 27" long


# Instructor and Demonstration Benches

**Instructor's Demonstration Benches** provide the space required to conduct effective table top demonstrations. Benches are available in two sizes and either left or right hand configurations. Both have cupboard and drawer storage areas and include a finished back, knee space filler panel, base molding

and molding corner guards. The Work Top includes a 1" thick black Kemresin® work top, sink, and sink outlet. The standard Fixture Package includes: service fittings, electrical fixtures, rod sockets, and an upright rod assembly.

**Note:** Last four digits of part number designates door style, hinge, hardwood species, and pull selection. See page 4 for full explanation. Actual product may vary slightly from catalog illustration due to style and hardware selection.

## 60" Instructor's Demonstration Bench



**T20 Instructors Bench** 30" Deep x 60" Long (Sink on Right)

**T30 Instructors Bench** 30" Deep x 60" Long (Sink on Left)

### Cabinet Assembly

Standing Height	(Sink on Right)	<b>T20WC363060-</b> _____ <i>specify style</i>
ADA Height	(Sink on Right)	<b>T20WC343060-</b> _____ <i>specify style</i>
Standing Height	(Sink on Left)	<b>T30WC363060-</b> _____ <i>specify style</i>
ADA Height	(Sink on Left)	<b>T30WC343060-</b> _____ <i>specify style</i>

Cabinets: (1) **\*G72W362818** Sink Cabinet (*right hand*)  
 or **\*G77W362818** Sink Cabinet (*left hand*)  
 (1) **\*G71W362818** Four Drawer Cabinet

\* ADA height bench uses ADA height cabinets

<b>Work Top Assembly</b>	(Sink on Right)	<b>T20WT013060R-</b> ___ <i>specify color</i>
<b>Work Top Assembly</b>	(Sink on Left)	<b>T30WT013060R-</b> ___ <i>specify color</i>

Work Top: **Kemresin®** 1" thick x 30" deep x 60" long  
 Sink: (1) **1000-00-**\_\_\_ Kemresin® – 6" deep x 10" wide x 14" long  
 (1) **0482-**\_\_\_ 1½" I.P.S. Sink Outlet

<b>Fixture Package</b>	(Sink on Right)	<b>T20FP-</b> ___ <i>specify receptacle color, service</i>
<b>Fixture Package</b>	(Sink on Left)	<b>T30FP-</b> ___ <i>specify receptacle color, service</i>

Electric Fixtures: (1) **S656-1**\_\_\_ 120 VAC GFCI Duplex Receptacle  
 Fittings: (1) **W-0337-6V** Cold Water w/VB Gooseneck  
 (1) **W-0260-0G** Single Outlet Service Fitting  
 (2) **0548-00** Rod Socket  
 (1) **4740-00** Upright Rod Assembly

**No service piping or electric conduit included.**



# Instructor and Demonstration Benches

## Instructor's Demonstration

**Benches** provide the space required to conduct effective table top demonstrations. Benches are available in two sizes and either left or right hand configurations. Both have cupboard and drawer storage areas and include a finished back, knee space filler panel, base molding


and molding corner guards.

The Work Top includes a 1" thick black Kemresin® work top, sink, and sink outlet.

The standard Fixture Package includes: service fittings, electrical fixtures, rod sockets, and an upright rod assembly.

**Note:** Last four digits of part number designates door style, hinge, hardwood species, and pull selection. See page 4 for full explanation. Actual product may vary slightly from catalog illustration due to style and hardware selection.

## 96" Instructor's Demonstration Bench



**T21 Instructors Bench** 30" Deep x 96" Long (Sink on Right)

**T31 Instructors Bench** 30" Deep x 96" Long (Sink on Left)

### Cabinet Assembly

Standing Height	(Sink on Right)	<b>T21WC363060-</b> _____ <i>specify style</i>
ADA Height	(Sink on Right)	<b>T21WC343060-</b> _____ <i>specify style</i>
Standing Height	(Sink on Left)	<b>T31WC363060-</b> _____ <i>specify style</i>
ADA Height	(Sink on Left)	<b>T31WC343060-</b> _____ <i>specify style</i>

<b>Cabinets:</b>	(1) <b>*G74W362836</b>	Sink Cabinet ( <i>right hand</i> )
	or	<b>*G78W362836</b>
	(1) <b>*G73W362836</b>	Sink Cabinet ( <i>left hand</i> )
		Five Drawer Cabinet

\* ADA height bench uses ADA height cabinets

<b>Work Top Assembly</b>	(Sink on Right)	<b>T21WT013060R-</b> ___ <i>specify color</i>
<b>Work Top Assembly</b>	(Sink on Left)	<b>T31WT013060R-</b> ___ <i>specify color</i>

<b>Work Top:</b>	<b>Kemresin®</b>	1" thick x 30" deep x 96" long
<b>Sink:</b>	(1) <b>1003-00-</b> ___	Kemresin® – 11" deep x 15" wide x 18" long I.D.
	(1) <b>0482-</b> ___	1½" I.P.S. Sink Outlet

<b>Fixture Package</b>	(Sink on Right)	<b>T21FP-</b> ___ <i>specify receptacle color, service</i>
<b>Fixture Package</b>	(Sink on Left)	<b>T31FP-</b> ___ <i>specify receptacle color, service</i>

<b>Electric Fixtures:</b>	(1) <b>S656-1</b> ___	120 VAC GFCI Duplex Receptacle
<b>Fittings:</b>	(2) <b>W-0337-0V</b>	Water w/VB Gooseneck
	(1) <b>W-0263-0G</b>	Double Outlet Service Fitting
	(2) <b>0548-00</b>	Rod Socket
	(1) <b>4740-00</b>	Upright Rod Assembly


**No service piping or electric conduit included.**

# Instructor and Demonstration Benches

**Instructor's Bench Desk** is designed to be added to the right or left side of the Instructor's Demonstration Benches creating additional working and storage space for the instructor.

**Note:** Last four digits of part number designates door style, hinge, hardwood species, and pull selection. See page 4 for full explanation. Actual product may vary slightly from catalog illustration due to style and hardware selection.

## Instructor's Bench Desk with Apron Drawer



**T24 Instructors Bench Desk** 31" High x 30" Deep x 48" Long **(Right Side)**  
**T34 Instructors Bench Desk** 31" High x 30" Deep x 60" Long **(Left Side)**

**Cabinet-Apron Assembly** (Right Side) **T24WC313048-**\_\_\_\_\_ *specify style*  
**Cabinet-Apron Assembly** (Left Side) **T34WC313048-**\_\_\_\_\_ *specify style*

**Includes: Cabinets** (1) **D90W301118-**\_\_\_\_\_ 3-Drawer Cabinet w/Pull Board  
 (1) **A01W052229-0** Apron with Drawer – 29" Long  
 (1) **G70W240647-00\_0** Book Shelf 6" D x 26" H x 47" L

### Right Side Work Tops

Kemresin **T24WT013048R-**\_\_\_ *specify color*  
 Black Plastic Laminate **T24WT013048P**  
 Black Acid Resistant Laminate **T24WT013048D**


**Work Top Size:** 1" thick x 30" deep x 48" long

### Left Side Work Tops

Kemresin **T34WT013048R-**\_\_\_ *specify color*  
 Black Plastic Laminate **T34WT013048P**  
 Black Acid Resistant Laminate **T34WT013048D**

**Work Top Size:** 1" thick x 30" deep x 48" long

## Instructor's Bench Desk with Keyboard Tray



**T25 Instructors Bench Desk** 31" High x 30" Deep x 48" Long **(Right Side)**  
**T35 Instructors Bench Desk** 31" High x 30" Deep x 60" Long **(Left Side)**

**Cabinet-Tray Assembly** (Right Side) **T25WC313048-**\_\_\_\_\_ *specify style*  
**Cabinet-Tray Assembly** (Left Side) **T35WC313048-**\_\_\_\_\_ *specify style*

**Includes: Cabinets** (1) **D90W301118-**\_\_\_\_\_ 3-Drawer Cabinet w/Pull Board  
 (1) **A80W002229-00\_0** Frame with Keyboard Tray – 29" Long  
 (1) **G70W240647-00\_0** Book Shelf 6" D x 26" H x 47" L

### Right Side Work Tops

Kemresin **T25WT013048R-**\_\_\_ *specify color*  
 Black Plastic Laminate **T25WT013048P**  
 Black Acid Resistant Laminate **T25WT013048D**

**Work Top Size:** 1" thick x 30" deep x 48" long

### Left Side Work Tops

Kemresin **T35WT013048R-**\_\_\_ *specify color*  
 Black Plastic Laminate **T35WT013048P**  
 Black Acid Resistant Laminate **T35WT013048D**


**Work Top Size:** 1" thick x 30" deep x 48" long

# Laboratory Carts

**Wood Laboratory Carts** are designed for efficient interaction where sturdy, dependable work surface height mobility is needed. Both the open and enclosed carts are mounted on 3" diameter rubber tired casters. Casters have brakes for stationary use.

**Note:** Last four digits of part number designates door style, hinge, hardwood species, and pull selection. See page 4 for full explanation. Actual product may vary slightly from catalog illustration due to style and hardware selection.

## Wood Laboratory Carts



### T00 Laboratory Cart 36" High x 25" Deep x 37" Long


**Table Assembly** **T00WC362537-00\_0** *specify species*

**Apron/Leg Assembly Size:** 35" High x 22" Deep x 35" Long

#### Work Top

Kemresin **T00WT012537R-\_\_** *specify color*  
 Black Plastic Laminate Top **T00WT012537P**  
 Black Acid Resistant Laminate **T00WT012537D**

**Work Top Size:** 1" thick x 25" deep x 37" long



### T01 Laboratory Cart 36" High x 30" Deep x 37" Long

**Cabinet Assembly** **T01WC363037-00\_0** *specify species*

**Cabinet Size:** 35" High x 26<sup>3</sup>/<sub>4</sub>" Deep x 35" Long

#### Work Top

Kemresin **T01WT013037R-\_\_** *specify color*  
 Black Plastic Laminate Top **T01WT013037P**  
 Black Acid Resistant Laminate **T01WT013037D**

**Work Top Size:** 1" thick x 30" deep x 37" long

# Mobile Demonstration Bench

**Mobile Demonstration Bench** is a totally self-contained movable laboratory bench with its own services as well as storage space for supplies. It is furnished complete with water supply and waste containers.

**Note:** Last four digits of part number designates door style, hinge, hardwood species, and pull selection. See page 4 for full explanation. Actual product may vary slightly from catalog illustration due to style and hardware selection.

## Mobile Demonstration Bench

**T56 Mobile Demo Bench** 36" High x 30" Deep x 66" Long

**Cabinet Assembly** **T56WC363066-**\_\_\_\_\_ *specify style*

**Cabinet:** 35" H x 27 1/2" D x 45" L

### Work Top

Black Plastic Laminate **T56WT013066-P**  
Black Acid Resistant Laminate **T56WT013066-D**

**Main Work Top:** 1" T x 30" D x 48" L

**Drop Leaf:** 1" T x 30" D x 18" L

### Fixture Package


**T56FP-**\_\_\_\_\_ *specify receptacle color*

**Sink:** (1) 3 1/2" dia. Stainless Steel Crumb Cup


**Electrical Fixtures:** (1) **S581-1**\_\_\_\_\_ 120 VAC Duplex Receptacle with 10 foot long extension cord

**Note: assembly is not pre-wired**

**Fittings:** (1) **W-0405-00** Cold Water Pump Gooseneck  
(2) One-Gallon Plastic Carboys  
(2) **0548-00** Rod Sockets  
(2) **0549-18** 1/2" Dia Upright Rods – 18" Long  
(1) **4740-00** Upright Rod Assembly



(front view)



(rear view)

# ADA Height Mobile Workstation

**ADA Height Mobile Workstation** designed to meet the height requirements of the Americans with Disabilities Act. Totally self-contained with its own services as well as storage space for supplies, it is furnished complete with water supply and waste containers. Electrical

pedestal is furnished with an eight-foot drop cord. Gas turret is furnished with hose bib on the supply side to connect to laboratory gas supply.

**Note:** Last four digits of part number designates door style, hinge, hardwood species, and pull selection. See page 4 for full explanation. Actual product may vary slightly from catalog illustration due to style and hardware selection.

## ADA Height Mobile Workstation

### T57 Mobile Workstation 31<sup>3</sup>/<sub>4</sub>" High x 30" Deep x 60" Long

**Cabinet Assembly** **T57WC323060-**\_\_\_\_\_ *specify style*

Cabinet: 30<sup>3</sup>/<sub>4</sub>" H x 27<sup>1</sup>/<sub>2</sub>" D x 58" L

**Work Top - Kemresin** **T57WT013060R-**\_\_\_ *specify color*


Work Top Size: 1" T x 30" D x 60" L

**Fixture Package** **T57FP-**\_\_\_ *specify rctp color, service*


**Sink:** (1) 3<sup>1</sup>/<sub>2</sub>" dia. Stainless Steel Crumb Cup

**Electrical Fixtures:** (1) **S581-2** 120 VAC Duplex Receptacle with 8 foot long extension cord  
**Note: assembly is not pre-wired**

**Fittings:** (1) **W-0405-00** Cold Water Pump Gooseneck  
(2) One-Gallon Plastic Carboys  
(1) **W-0260-0** Deck Mounted Service Outlet  
(2) **0548-00** Rod Sockets  
(2) **0549-18** 1/2" Dia Upright Rods – 18" Long  
(1) **4740-00** Upright Rod Assembly



(front view)



(rear view)

# Perimeter Student Table Assemblies

Perimeter Student Tables are ideally suited for the teaching laboratory. Select the exact combination of storage space needed for the students.


The X75 Assembly is furnished with two pair of 16" deep x 24" long cabinets. one pair at the wall and

the other under the peninsula. The X80 Assembly is the same, except 16" deep x 36" long cabinets are furnished under the peninsula.

To specify the 16" deep base cabinets you want, insert the cabinet's two digit code number into the assembly catalog number.

**Note:** Last four digits of part number designates drawer head style, hinge, hardwood species, and pull selection. See page 4 for full explanation. Actual product may vary slightly from catalog illustration due to style and hardware selection.

## Perimeter Student Tables



### X75W369296-3015

Shown with G92W361624-3015 at wall and G96W361624-3015 under peninsula.

### X75 Perimeter Student Table 84" Deep x 90" Long

#### Cabinet Assembly

**ADA Height** X75WC34 -\_\_\_\_\_ specify cabinets & style

**Standing Height** X75WC36 -\_\_\_\_\_ specify cabinets & style

Cabinets: (2) pair as specified in part number

#### Work Top Assembly

**Kemresin** X75WT018490-R-\_\_\_\_\_ specify color

**Work Top:** 1" thick - 18" x 90" wall - 42" x 66" peninsula

**Sink:** (1) 1011-00-\_\_ Kemresin -6" H x 18" L x 6 1/2" W

(1) 0482-\_\_ 1 1/2" I.P.S. Sink Outlet


#### Fixture Package X75FP-\_\_ specify receptacle color

**Fixtures:** (2) S656-1\_ 120 VAC GFCI Duplex Rctp.

**Fittings:** (2) W-0274-VB Combination Water & Gas

(4) 0548-00 Rod Socket

**No service piping or electric conduit included.**



### X80W369398-3015\_

Shown with G93W361624-3015 at wall and G98W361636-3015 under peninsula.

### X80 Perimeter Student Table 84" Deep x 90" Long

#### Cabinet Assembly

**ADA Height** X80WC34 -\_\_\_\_\_ specify cabinets & style

**Standing Height** X80WC36 -\_\_\_\_\_ specify cabinets & style

Cabinets: (2) pair as specified in part number

#### Work Top Assembly

**Kemresin** X80WT018490-R-\_\_\_\_\_ specify color

**Work Top:** 1" thick - 18" x 90" wall - 42" x 66" peninsula

**Sink:** (1) 1011-00-\_\_ Kemresin -6" H x 18" L x 6 1/2" W

(1) 0482-\_\_ 1 1/2" I.P.S. Sink Outlet

#### Fixture Package X80FP-\_\_ specify receptacle color

**Fixtures:** (2) S656-1\_ 120 VAC GFCI Duplex Rctp.

**Fittings:** (2) W-0274-VB Combination Water & Gas

(4) 0548-00 Rod Socket

**No service piping or electric conduit included.**

# Perimeter Student Table Assemblies


Each Perimeter Student Table Cabinet Assembly is furnished with filler panels, scribe strips, sink supports, base molding and molding corner guards. Work Top Assemblies are furnished with 1" Kemresin top, sink, and sink outlet. Illustrated at the bottom of this page

are examples of only some of the perimeter tables that can be created. Any combination of the cabinets shown below may be used.


**Note:** All cabinets are furnished in pairs. 24" long cupboard cabinets are furnished with one left hand door swing and one right hand.

**Note:** Last four digits of part number designates drawer head style, hinge, hardwood species, and pull selection. See page 4 for full explanation. Actual product may vary slightly from catalog illustration due to style and hardware selection.


## Perimeter Student Table - Cabinet Options (Nominal 16" Deep)




G92W341624- 24" L  
G92W361624- 24" L  
Two Digit Code **92**




G93W341624- 24" L  
G93W361624- 24" L  
Two Digit Code **93**




G94W341624- 24" L  
G94W361624- 24" L  
Two Digit Code **94**  
(6) OP-6140 Trays Included




G95W341624- 24" L  
G95W361624- 24" L  
Two Digit Code **95**




G96W361624- 24" L  
Two Digit Code **96**



G99W361624- 24" L  
Two Digit Code **99**




G97W361636- 36" L  
Two Digit Code **97**




G98W361636- 36" L  
Two Digit Code **98**


## Cabinet Combination Examples




X75W369293-




X75W369292-




X75W369394-




X75W369494-




X75W369294-




X75W369296-




X75W369393-




X75W369495-




X75W369599-



X75W369699-



X80W369397-



X80W369698-


# All-Science Student Table Assemblies

**Octagon All-Science Tables** have been designed to serve as an "all purpose" science table. It is completely free-standing and offers up to four students equal work space with the option of four suspended base cabinets. It consists of a 24" by 24" pedestal with a service access panel and an apron to support the 1"

thick epoxy resin work top. Eight optional Suspended Base Cabinets are available for the Octagon All-Science Table. All are 15" deep by 21" long and either 18" or 20" high.

**Note:** Last four digits of part number designates drawer head style, hinge, hardwood species, and pull selection. See page 4 for full explanation. Actual product may vary slightly from catalog illustration due to style and hardware selection.

## Octagon All-Science Table



### X02 Octagon Science Table 55 1/2" x 55 1/2"

**Cabinet Assembly**

ADA Height **X02WC345656-00\_0** *specify material*  
Standing Height **X02WC365656-00\_0** *specify material*

**Pedestal Size:** 23 1/2" x 23 1/2"  
**Cabinets:** Optional – Order separately

**Work Top Assembly**

Kemresin **X02WT015656-R-\_\_\_** *specify color*

**Work Top:** 1" Thick Kemresin  
55 1/2" X 55 1/2" — 66 1/2" X 66 1/2" diagonally


**Fixture Package** **X02FP-\_\_** *specify receptacle color*

**Electrical Fixtures:** (2) **S656-1** 120 VAC GFCI Duplex Rctpt.  
**Fittings:** (4) **0548-00** Rod Socket  
(4) **0553-00** Upright Rod

**No service piping or electric conduit included.**

## Perimeter Annex

**Perimeter Annex** assembly is available to connect the Octagon All-Science Table to a perimeter wall assembly. It is available with a Kemresin work top and filler panels to hide electrical and service piping.



### X05 Perimeter Annex 13" Wide x 41 1/2" Long

**Cabinet Assembly**

ADA Height **X05WC341342-00\_0** *specify material*  
Standing Height **X05WC361342-00\_0** *specify material*

**Pedestal Size:** 11" Wide x 41 1/2" Long

**Work Top Assembly**

Kemresin **X05WT011342-R-\_\_\_** *specify color*

**Work Top:** 1" Thick Kemresin - 13" Wide x 18" Long



# All-Science Student Table Assemblies

## Octagon All-Science Tables


have been designed to serve as an "all purpose" science table. It is completely free-standing and offers up to four students equal work space with the option of four suspended base cabinets. It consists of a 24" by 24" pedestal with a service access panel and an apron to support the 1"

thick epoxy resin work top.

Eight optional Suspended Base Cabinets are available for the Octagon All-Science Table. All are 15" deep by 21" long and either 18" or 20" high.

**Note:** Last four digits of part number designates drawer head style, hinge, hardwood species, and pull selection. See page 4 for full explanation. Actual product may vary slightly from catalog illustration due to style and hardware selection.

## Octagon All-Science Table



### X01 Octagon Science Table

55 $\frac{1}{2}$ " x 55 $\frac{1}{2}$ "

#### Cabinet Assembly

ADA Height **X01WC345656-00\_0** *specify material*  
 Standing Height **X01WC365656-00\_0** *specify material*

**Pedestal Size:** 23 $\frac{1}{2}$ " x 23 $\frac{1}{2}$ "

**Cabinets:** Optional – Order separately

#### Work Top Assembly

Kemresin **X01WT015656-R-\_\_** *specify color*

**Work Top:** 1" Thick Kemresin

55 $\frac{1}{2}$ " X 55 $\frac{1}{2}$ " — 66 $\frac{1}{2}$ " X 66 $\frac{1}{2}$ " diagonally

**Sink:** (1) **SM-39-00-\_\_** 4 $\frac{3}{4}$ " H x 6 $\frac{3}{8}$ " W x 35" L I.D.  
 (1) **0482-BP** 1 $\frac{1}{2}$ " I.P.S. Sink Outlet

#### Fixture Package

**X01FP-\_\_** *specify receptacle color*


**Electrical Fixtures:** (2) **S656-1\_** 120 VAC GFCI Duplex Receptacle

**Fittings:** (2) **W-0416-VB** Combination Water & Gas

(4) **0548-00** Rod Socket

(4) **0553-00** Upright Rod

**No service piping or electric conduit included.**



### X03 Octagon Science Table

55 $\frac{1}{2}$ " x 55 $\frac{1}{2}$ "

#### Cabinet Assembly

ADA Height **X03WC345656-00\_0** *specify material*  
 Standing Height **X03WC365656-00\_0** *specify material*

**Pedestal Size:** 23 $\frac{1}{2}$ " x 23 $\frac{1}{2}$ "

**Cabinets:** Optional – Order separately

#### Work Top Assembly

Kemresin **X01WT015656-R-\_\_** *specify color*

**Work Top:** 1" Thick Kemresin

55 $\frac{1}{2}$ " X 55 $\frac{1}{2}$ " — 66 $\frac{1}{2}$ " X 66 $\frac{1}{2}$ " diagonally

**Sink:** (1) **1039-00-\_\_** 8" High x 18" Dia. I.D.  
 (1) **0482-\_P** 1 $\frac{1}{2}$ " I.P.S. Sink Outlet

#### Fixture Package

**X01FP-\_\_** *specify receptacle color*

**Electrical Fixtures:** (2) **S656-1\_** 120 VAC GFCI Duplex Receptacle

**Fittings:** (2) **W-0416-VB** Combination Water & Gas

(4) **0548-00** Rod Socket

(4) **0553-00** Upright Rod

**No service piping or electric conduit included.**

# All-Science Student Table Assemblies

**Optional Suspended Base Cabinets** for the Octagon All-Science Table are available either with or without an open book compartment at the top of the cabinet.

Cabinets without the book compartment are 18" high, 15" deep, and 21" long. Cabinets with book


compartments are 20" high, 15" deep, and 21" long.

Both are available in four different door and drawer configurations and the three standard Kewaunee drawer head styles.

**Note:** Last four digits of part number designates drawer head style, hinge, hardwood species, and pull selection. See page 4 for full explanation. Actual product may vary slightly from catalog illustration due to style and hardware selection.


## Optional Suspended Base Cabinets

### Cupboard Cabinet




**G41W181521-**\_\_\_\_\_ 21" L  
(1) 12u Cupboard

### 3 Drawer Cabinet




**G42W181521-**\_\_\_\_\_ 21" L  
(3) 4u Drawers

### 2 Drawer Cabinet



**G43W181521-**\_\_\_\_\_ 21" L  
(1) 4u Drawer  
(1) 8u Drawer


### 6 Drawer Cabinet



**G44W181521-**\_\_\_\_\_ 21" L  
(6) 4u Drawers


## Optional Suspended Base Cabinets with Book Compartment

### Cupboard Cabinet




**G51W201521-**\_\_\_\_\_ 21" L  
(1) 12u Cupboard

### 3 Drawer Cabinet




**G52W201521-**\_\_\_\_\_ 21" L  
(3) 4u Drawers

### 2 Drawer Cabinet




**G53W201521-**\_\_\_\_\_ 21" L  
(1) 4u Drawer  
(1) 8u Drawer

### 6 Drawer Cabinet




**G54W201521-**\_\_\_\_\_ 21" L  
(6) 4u Drawers

Order Suspended Base Cabinets separately. Each Octagon Table can hold up to four cabinets.



**X02W365656-0010** with  
(4) **G41W181521-5018** Cabinets



**X02W365656-0010**  
with (4) **G52W181521-5018** Cabinets

# Multi-Science Tables

## Sitting Height Multi-Student Tables:


Free-standing multi-student tables that can be used for combination of sciences. Designed for laboratories where permanent individual student storage space at the table is not required. Each student is provided with an open book/tote tray storage area at his seat. All tables include

cupboard space to store equipment. Cupboard includes one F-5320 Tote Tray (12" x 14" x 4"). Available in single-sided and double-sided designs.

**Note:** Last four digits of part number designates drawer head style, hinge, hardwood species, and pull selection. See page 4 for full explanation. Actual product may vary slightly from catalog illustration due to style and hardware selection.

## 24" Deep Single-sided Table

### T44 Two Student Table 31" High x 24" Deep x 60" Long



#### Cabinet Assembly **T44WC312460-**\_\_\_\_\_ *specify style*

Cabinet: (1) **G75W242114-**\_\_\_\_\_ Cupboard Cabinet with Tote Tray

#### Work Top - Kemresin **T44WT012460R-**\_\_\_ *specify color*

Work Top: 1" T x 24" D x 60" L


#### Fixture Package **T44FP-**\_\_\_,\_\_\_ *specify color & service*

Electric Fixtures: (1) **S656-1**\_ 120 VAC GFI Duplex Receptacle  
(specify color)  
Fittings: (1) **W-0263-0**\_ Double Outlet Service Fitting  
(specify service)  
(2) **0548-00** Rod Sockets

**No service piping or electrical conduit included.**

## 38" Deep Double-sided Table

### T84 Four Student Table 31" High x 24" Deep x 60" Long



#### Cabinet Assembly **T84WC313860-**\_\_\_\_\_ *specify style*

Cabinet: (1) **G76W243714-**\_\_\_\_\_ Cupboard Cabinet with Tote Tray

#### Work Top - Kemresin **T84WT013860R-**\_\_\_ *specify color*

Work Top: 1" T x 24" D x 60" L

#### Fixture Package **T84FP-**\_\_\_,\_\_\_ *specify color & service*

Electric Fixtures: (2) **S656-1**\_ 120 VAC GFI Duplex Receptacle  
(specify color)  
Fittings: (1) **W-0272-0**\_ Four-way Outlet Service Fitting  
(specify service)  
(2) **0548-00** Rod Sockets

**No service piping or electrical conduit included.**

# Safety Cabinets

**Safety Cabinet:** a valuable addition for any laboratory to concentrate many safety equipment items in one location. Two half-depth fixed shelves and the lower cupboard provide storage space for fire blankets, fire extinguisher, first aid kit and supplies. The cabinet is equipped with an Emergency Shower Head Assembly and a Deck


Mounted Eyewash Assembly.

The **ADA Safety Cabinet** is the same except the Emergency Shower is fitted with a longer handle, and a Dual Purpose Eye Wash/Drench Hose is substituted for the Deck Mounted Eye Wash.

Both are furnished with base molding and corner clips.

**Note:** Last four digits of part number designates door and drawer head style, hinge, hardwood species, and pull selection. See page 4 for full explanation. Actual product may vary slightly from catalog illustration due to style and hardware selection.

## Safety Cabinet



### G81 Safety Cabinet 90" High x 21 1/4" Deep x 36" Long

#### Cabinet Assembly **G81WC363066-**\_\_\_\_\_ *specify style*

**Cabinet:** 90" H x 21 1/4" D x 36" L with two fixed shelves and double cupboard


#### Fixture Package **G81FP**

**Eyewash:** Deck Mounted with two Gentle Spray Outlet Heads, 1 1/2" Dia. ABS Plastic Bowl and Stay-Open Ball Valve with Epoxy Coated Flag Handle

**Shower:** 10" Dia. ABS Plastic Shower Head with 1" I.P.S. Stay-Open Ball Valve with 29" long Pull Rod.

**No service piping or electrical conduit included.**

## ADA Safety Cabinet



### G85 ADA Safety Cabinet 90" High x 31 5/8" Deep x 36" Long

#### Cabinet Assembly **G85WC902236-00\_0** *specify species*

**Cabinet:** 90" H x 21 1/4" D x 36" L with two fixed shelves and apron

#### Work Top - Kemresin **G85WT012236R-**\_\_\_ *specify color*

**Cabinet:** 1" T x 28 5/8" D x 34 1/2" L with two fixed shelves and apron

#### Fixture Package **G85FP**

**Eyewash:** Deck Mounted Stainless Steel Bowl with four Gentle Spray Outlet Heads and plunger operated ball valve.

**Shower:** 10" Dia. ABS Plastic Shower Head with 1" I.P.S. Stay-Open Ball Valve with 59 1/2" long Pull Rod.

**No service piping or electrical conduit included.**

# Demonstration & Teaching Hoods

## Duo-Vision Demonstration Hood

Double-sided teaching and demonstration hood



### T95 Duo-vision Demonstration Hood

#### Hood Assembly

Glass End Panels	<b>T95G483058-GE</b>
Glass End Panels/One Sash Fixed	<b>T95G483058-FS</b>
Solid End Panels	<b>T95G483058-SE</b>

**Exterior Material:** Powder Coated Cold Rolled Steel

**Exterior Size:** 48" H x 30" D x 58" L

**Interior Liner:** Kemglass - Fiberglass Reinforced Polyester

**Interior Size:** 40" H x 25" D x 53" L

#### Cabinet Assembly

ADA Height **T95WC343059-\_\_\_** *specify style*

Standing Height **T95WC363059-00\_0** *specify style*

**Cabinet Size:** 34<sup>3</sup>/<sub>4</sub>"/32<sup>1</sup>/<sub>2</sub>" H x 30" D x 59" L

**Cabinets:** Optional – Order separately

#### Work Top Assembly

Kemresin **T95WT0132606R-\_\_\_** *specify color*

**Work Top:** 1<sup>1</sup>/<sub>4</sub>" Thick Dished Kemresin - 30" D x 60" L

**Sink:** (1) **1038-00-\_\_\_** 7<sup>5</sup>/<sub>8</sub>" High x 12" Dia. I.D.

(1) **0482-BP** 1<sup>1</sup>/<sub>2</sub>" I.P.S. Sink Outlet

#### Fixture Package **T95FP-\_\_\_** *specify rpt. color & service*

**Electrical Fixtures:** (2) **S656-1\_** 120 VAC GFCI Duplex Receptacle  
 (2) 3-way Light Switch  
 (1) Blower Switch with Warning Light  
 (1) Blower Warning Light  
 (2) T8 Light Fixtures (tubes not included)

**Fittings:** (1) **W-0337-0V** Cold Water VB Gooseneck

(1) **W-0260-0\_** Single Outlet Service Fitting

**No service piping or electric conduit included.**

## Features:

- Dual working sashes with glass end panels for excellent visibility into hood.
- Available with one sash in fixed closed position for higher access control.
- Furnished with 12" diameter sink with water faucet and service fitting.
- Furnished with 3-way light switch, blower switch with warning light and gfi duplex receptacles.
- Available with or without base cabinet, work top, fittings and fixtures.

## Specifications:

Recommended Face Velocity: 80 feet/minute

Sash Opening: 53" x 30" \*

Duct Collar: 12" Dia

Exhaust Volume: 940 cfm

Static Pressure: 0.30" @ 940 c.f.m.

\* based on only one sash open at a time.

# TruView Teaching Hoods



## Features:




- Available in 4-foot, 5-foot, and 6-foot lengths, either single-sided for use against walls, or double-sided for use in center of laboratory.
- Glass end-walls and center baffles provide unobstructed view through hood.
- Uni-baffle design allows both halves of the double-sided hood to function independently as a complete fume hood.
- Sash-stops provided at 18" for economical hood operation.
- Designed for use as a stand-alone hood or in assemblies of multiple hoods sharing common endwalls.
- Frameless sash with full-length handle for neat, clean appearance and streamline air flow.
- Flush-sill design allows easy access and excellent airflow.
- Designed and tested to meet ASHRAE 110-2016 standards.

## Typical Standing Height Configurations

4'-0"

5'-0"

6'-0"



## Typical ADA High Configurations






4'-0"

5'-0"

4'-0"

5'-0"

6'-0"




# TruView Teaching Hoods

## Ordering Instructions:

1. Select Teaching Hood superstructure. When planning a multiple-hood assembly, select both a Right End and Left End hood, and as many Middle hoods as needed to finish the assembly.
2. Select a Hood Work Top for each hood. Work Tops are available in either Black or Grey dished Kemresin.
3. Select the required number of cupsinks for each hood.
  - 1 for 48" Single-sided
  - 2 for 48" Double-sided
  - 2 for 60" and 72" Single-sided
  - 4 for 60" and 72" Double-sided
4. Select Service Fittings and Electrical Fixtures. Each hood can accept up to four service fittings at each cupsink location. Hoods ship with removable plug buttons in service fitting holes.
5. Select base cabinets to fill the area below hood. Sink cabinet is needed on right end of 48" hoods and both ends of 60" and 72" hoods to house cupsink and service fittings.
6. Select Filler Panels as required.


## Single-sided Hood Superstructure



Stand Alone	<b>T90M3048-SA</b>	<b>T90M3060-SA</b>	<b>T90M3072-SA</b>
Left End	<b>T90M3048-LE</b>	<b>T90M3060-LE</b>	<b>T90M3072-LE</b>
Right End	<b>T90M3048-RE</b>	<b>T90M3060-RE</b>	<b>T90M3072-RE</b>
Middle	<b>T90M3048-AD</b>	<b>T90M3060-AD</b>	<b>T90M3072-AD</b>

Overall Hood Length	Duct Collar Dimensions	Sash Opening Sq. Ft. *	Total CFM and Static Pressure			
			80 FPM	Static Pressure	100 FPM	Static Pressure
4'-0" / 48"	5" x 20"	5.5	440	0.10	550	0.15
5'-0" / 60"	5" x 25"	7.0	560	0.10	700	0.15
6'-0" / 72"	5" x 30"	8.5	680	0.10	850	0.15

## Double-sided Hood Superstructure



Stand Alone	<b>T90M6048-SA</b>	<b>T90M6060-SA</b>	<b>T90M6072-SA</b>
Left End	<b>T90M6048-LE</b>	<b>T90M6060-LE</b>	<b>T90M6072-LE</b>
Right End	<b>T90M6048-RE</b>	<b>T90M6060-RE</b>	<b>T90M6072-RE</b>
Middle	<b>T90M6048-AD</b>	<b>T90M6060-AD</b>	<b>T90M6072-AD</b>

Overall Hood Length	Duct Collar Dimensions	Sash Opening Sq. Ft. *	Total CFM and Static Pressure			
			80 FPM	Static Pressure	100 FPM	Static Pressure
4'-0" / 48"	9" x 20"	11.0	880	0.10	1100	0.20
5'-0" / 60"	9" x 25"	14.0	1120	0.15	1400	0.20
6'-0" / 72"	9" x 30"	17.0	1360	0.15	1700	0.20

\* Based on 18" high sash opening

# TruView Teaching Hood Accessories


**Cabinet Options:**

Shown below are the most common cabinets used below TruView fume hoods. Most other base cabinets shown in this catalog may be used.

**Note:** Last four digits of part number designates door and drawer head style, hinge, hardwood species, and pull selection. See page 4 for full explanation. Actual product may


vary slightly from catalog illustration due to style and hardware selection.

**Wood Standing Height Base Cabinets 35" high**




- |                               |                               |                                 |                               |
|-------------------------------|-------------------------------|---------------------------------|-------------------------------|
| <b>G10W362212L</b> _____ 12"L | <b>G10W362212-</b> _____ 12"L | <b>D30W362218- 0</b> _____ 18"L | <b>C03W362230-</b> _____ 30"L |
| <b>G10W362218L</b> _____ 18"L | <b>G10W362218-</b> _____ 18"L | <b>D30W362224- 0</b> _____ 24"L | <b>C03W362236-</b> _____ 36"L |
| without fixture cutouts       | without fixture cutouts       |                                 |                               |
| <b>G00W362212L</b> _____ 12"L | <b>G00W362212-</b> _____ 12"L |                                 |                               |
| <b>G00W362218L</b> _____ 18"L | <b>G00W362218-</b> _____ 18"L |                                 |                               |

**Wood ADA Height Base Cabinets 32½" high Pedestal Legs**




- |                               |                               |   |   |
|-------------------------------|-------------------------------|---|---|
| <b>G10W342212L</b> _____ 12"L | <b>G10W342212-</b> _____ 12"L | <b>A25W3430-00_0</b> 33" H  | <b>A26W3459-00_0</b> 33" H  |
| <b>G10W342218L</b> _____ 18"L | <b>G10W342218-</b> _____ 18"L | <b>A25W3630-00_0</b> 35" H  | <b>A26W3659-00_0</b> 35" H  |
| without fixture cutouts       | without fixture cutouts       | Pedestal Legs are 1½" thick and 29¾" deep. They may be used instead of base cabinets to support the hood. They should be used with a Knee Space Panel and Apron Rail for stability. | Pedestal Legs are 1½" thick and 58½" deep. They may be used instead of base cabinets to support the hood. They should be used with a Knee Space Panel and Apron Rail for stability. |
| <b>G00W342212L</b> _____ 12"L | <b>G00W342212-</b> _____ 12"L |   |   |
| <b>G00W342218L</b> _____ 18"L | <b>G00W342218-</b> _____ 18"L |   |   |


**Apron Rails Knee Space Panels Filler Panels**



**A00W0330-00\_0** 30"L  
**A00W0336-00\_0** 36"L  
Apron Rail ¾" thick, 3" high, can be cut to length in field. Furnished with angles for attachment to adjacent units or pedestal legs.



**X-WK3036-00\_0** 30" L  
**X-WK3636-00\_0** 36" L  
¾" thick plywood, 36" W. Designed to be field-cut to exact width and length required. Furnished with two cleats




**X-WP0836-00\_0** 8"W x 36"H  
**X-WP1436-00\_0** 14"W x 36"H  
Fabricated from ¾" thick hardwood plywood with edges cut square for flush installation.



# TruView Teaching Hood Accessories

## Single-sided Hood Work Tops




Black Kemresin

**T90R3048-BK** 48"L  
**T90R3060-BK** 60"L  
**T90R3072-BK** 72"L

Six foot and five foot tops are furnished with two cupsink cutouts, one each on the left and the right. Four foot tops are furnished with only one, on the right.

## Double-sided Hood Work Tops



Black Kemresin

**T90R6048-BK** 48"L  
**T90R6060-BK** 60"L  
**T90R6072-BK** 72"L

Six foot and five foot tops are furnished with four cupsink cutouts, two on each side, one on the left and one on the right. Four foot tops are furnished with only two, on the right of each side.

## Cupsinks



**0491-BP** Black Poly

Complete with gasket and removable strainer. 6"x 3" inside dimension. 1½" IPS male straight thread outlet. (Overall height is 8")


## Service Fittings



**W-0744-A** Air  
**W-0744-G** Gas  
**W-0744-V** Vacuum

Front-load remote control service valves mount in top rail of sink base cabinet. Air, Gas, and Vacuum valves are furnished with color-coded, angled, nylon hose connector.


## Water Fittings




**W-0744-W** Water

Front load remote control water valves mount in top rail of sink base cabinet. Water valves are furnished with color-coded epoxy-coated brass hose tip.

## Electrical Fixtures




**S656-1V** Receptacle  
 120 volt, GFI, specification grade, 20 amp, ground fault protected, single duplex receptacle with electric box and stainless steel face plate.



**S655-1V** Light Switch  
 single pole, 120/240 volt AC, 20 amp. with electric box and stainless steel face plate.

## Air Alert 300 Face Velocity Airflow Alarm



**X-018528-S Air Alert 300/AirGard 200**

consists of a therm-istor sensor mounted through the end wall of the hood, and a control monitor that gives both a visual and audible alarm. The alarm monitors the fume hood face velocity and sounds an alarm when the airflow falls below safe levels. A glowing green light signals when conditions are again safe. The control monitor also contains a test/reset button that allows the hood user to verify alarm readiness.

**Double-sided TruView hoods require two alarms, one for each sash opening.**

# Fume Hoods & Ventilated Enclosures

## Ducted Ventilated Products

Kewaunee offers a full line of Laboratory Fume Hoods to meet the needs of every lab application.

In addition to our newest line of **Supreme Air Venturi Fume Hoods**, we also offer uniquely designed

hoods for education and specialty hoods for custom applications.

### **Supreme Air – General Purpose Fume Hoods**



The design is based on rigid frame construction that assures a solid installation and low vibration and sound levels. Access panels are easily removable to access service lines and fittings. Radius corner posts and airfoils, plus, easy operator control of interior baffle settings assure a high level of comfort, safety and efficiency.

### **Supreme Air Venturi – Next Generation Fume Hoods**

The Supreme Air Venturi utilizes proven containment enhancing techniques coupled with innovative laboratory convenience features, including: adjustable LED lights, electromechanical sash stop, ultra-smooth sash system, and push-button sash management. The Supreme Air Venturi provides the safest user friendly fume hood manufactured today.

### **TruView – Teaching Fume Hoods**

Designed for today's classroom allowing clear line-of-sight through multiple hoods to allow easier supervision of students by the instructor. May be used as a stand-alone hood or in assemblies, both single-sided for use against the wall, or double-sided for use in the center of the laboratory.



For more information on Kewaunee's **Supreme Air Fume Hoods** contact your local Kewaunee sales representative or visit the Kewaunee website: <http://www.kewaunee.com/LaboratoryHoods/FumeHoods/SupremeAirHoods/>

For more information on Kewaunee's **Supreme Air Venturi Fume Hoods** contact your local Kewaunee sales representative or visit the Kewaunee website: <http://www.kewaunee.com/LaboratoryHoods/FumeHoods/SupremeAirVenturi/>

For more information on Kewaunee's **TruView Teaching Fume Hoods** contact your local Kewaunee sales representative or visit the Kewaunee website: <http://www.kewaunee.com/LaboratoryHoods/FumeHoods/TruViewHood/>


# Fume Hoods & Ventilated Enclosures

## Ductless Ventilated Products

Kewaunee offers a full line of filtered products to meet every protection requirement. **DetectAir** ductless fume hoods provide a safe alternative to conventional chemical fume hoods

### **DetectAir – Ductless/Filtered Fume Hoods**

With the introduction of DetectAir, the filtered fume hood has come of age. DetectAir combines Kewaunee's traditional fume hood expertise with the filtering technology of AirClean Systems.



For more information on Kewaunee's **DetectAir - Ductless/Filtered Fume Hoods** contact your local Kewaunee sales representative or visit the Kewaunee website: <http://www.kewaunee.com/LaboratoryHoods/DuctlessHoods/>

when total exhaust is not an option. The **Interceptor Biological Safety Cabinets** are offered in both class II A2 and B2 models, (NSF Certified) for both recirculating and direct exhaust

### **the Interceptor – Biological Safety Cabinet**

The Interceptor Class II, A2 and Class II, B2 Biological Safety Cabinets are Kewaunee's latest design. Developed with comprehensive features created to enhance both safety and ergonomics... the Interceptor leads its class.




For more information on **the Interceptor Biological Safety Cabinet** from Kewaunee, contact your local Kewaunee sales representative for visit the Kewaunee website: <http://www.kewaunee.com/LaboratoryHoods/BioSafetyCabinets/>

requirements. The **Interceptor VLF** provides ISO 14644-1 Standard Class 5 conditions, providing a sterile environment for product protection.

### **the Interceptor – Vertical Laminar Flow Workstation**

Interceptor Vertical Laminar Flow Workstations feature a front loading HEPA filter, hinged front panel with self-support brackets, 10° angled safety glass, Minihelic gauge and T5 fluorescent lamp.



For more information on **the Interceptor - Vertical Laminar Flow Workstation** from Kewaunee, contact your local Kewaunee sales representative for visit the Kewaunee website: <http://www.kewaunee.com/LaboratoryHoods/LaminarFlowHoods/>


# Work Tops

Work Top material selections are among the most critical of the laboratory furniture decisions. The intended applications for your work tops dictate material choice.

Surfaces could be exposed to reagents and caustic substances. Heavy equipment may need to be supported or only light duty tasks may be performed.

Kewaunee can recommend and furnish materials fabricated to your individual needs, complete with appropriate sinks, drains and accessories.

**Kemresin® Work Tops manufactured exclusively by Kewaunee**



**Appearance:** Ultra-smooth non-porous surface.

**Available Colors:** Grey, Slate, Putty, and Black

**Composition:** 1" and thick molded epoxy resin.

**Physical Properties:** Excellent resistance to practically all acids, alkalis and solvents. Highly resistant to heat and normal physical abuse.

**Applications:** Use wherever optimum physical and chemical resisting properties are required.

**Other Work Top Materials**

MATERIAL	APPEARANCE	COMPOSITION	CHEMICAL RESISTANCE	APPLICATIONS
STAINLESS STEEL	No.4 Finish	Type 304 furnished as standard, but other types are available when specified. Type 316 supplied for Perchloric Acid fume hoods.	Good chemical resistance to a wide range of reagents.	Recommended for sterile and radiological laboratories.
HARDWOOD, OAK OR MAPLE	Natural varnish, or penetrating oil finish	Solid edge grain Maple or Oak strips, glued together electronically, and finished with baked-on clear synthetic varnish or penetrating oil.	Fair chemical resistance, but not recommended for use in chemical laboratories.	Excellent top material for Physics tables, Electronic benches, and other tables where mechanical or physical work levels are prevalent.
COMPOSITE RESIN	Hard, smooth non-glare surface	Cellulose fiber reinforced phenolic resin core with a polyurethane copolymer surface. Available in multiple colors and patterns.	Good chemical resistance to wide range of acids, alkalis and solvents.	Use wherever a high degree of chemical resistance is required. (Not recommended for fume hood worksurfaces.)
HIGH-PRESSURE LAMINATES	Hard, smooth surface	Plastic laminated surface attached with water resistant adhesive to solid particle board. Available in choice of colors and patterns.	Fair chemical resistance to wide range of acids, alkalis and solvents.	Recommended for light duty purposes.

# Work Tops

**Specifications:**

The configurations shown are for Kewaunee base cabinets. Standard top thickness is 1", standard curb height is 4".

To order, select configuration and top material desired. Specify number of linear or square feet per product number.

Dimensions are nominal. Not all configurations are available in all materials. Indented tops available on special order.


**Wall Bench Tops**

**Flat Top no curb Style B**

Standard	18"
Sizes:	23" 24" 30" 31"

**Top with Applied Curb Style C**

Standard	18"
Sizes:	23" 24" 30" 31"



**Corner Top no curbs Style K**

	A	B	C
Standard	33"	24"	9"
Sizes:	39"	30"	9"
	40"	31"	9"

**Flat Top & Ledge no curb Style E**

	A	B	C
Standard	30"	23"	7"
Sizes:	31"	23"	8"

**Top & Ledge with Applied curb Style F**

	A	B	C
Standard	30"	23"	7"
Sizes:	31"	23"	8"

**Corner Top with Applied Curbs Style L**

	A	B	C
Standard	33"	24"	9"
Sizes:	39"	30"	9"
	40"	31"	9"

**Center Bench Tops**

**Top without Center Joint Style H**

Standard	54"
Sizes:	60" 62"

**Top with Center Joint Style I**

	A	B	C
Standard	54"	27"	27"
Sizes:	54"	23"	31"
	60"	30"	30"
	62"	31"	31"

**Top with Center Ledge Style J**

	A	B	C	D
Standard	54"	23"	8"	23"
Sizes:	60"	23"	14"	23"
	62"	23"	16"	23"

**Ledges and Shelves**

**Ledge Unfinished Edges Style A**

Standard	8"
Sizes:	10" 12" 14" 16"

**Shelf Finished Edges Style N**

Standard	8"
Sizes:	10" 12" 14" 16"

# Center Bench Shelving Assemblies


## Specifications:

Shelving Support Assemblies are designed to add shelving to center benches. Each assembly consists of four uprights with slots on 1" centers and a top frame. They are 12" deep in lengths as shown.

Each assembly is furnished with mounting brackets for attaching to the counter top and cabinets below but no fasteners. Fasteners must be furnished in field depending on top material and cabinet configuration.


Shelving Assembly Support Struts are especially designed to provide added support to the Shelving Support Assemblies while allowing space for connection of plumbing or electrical connections.

### 12" Deep Shelving Support Assembly




Length	Part Number
24"	<b>F17-4824-00</b>
30"	<b>F17-4830-00</b>
36"	<b>F17-4836-00</b>
42"	<b>F17-4842-00</b>
48"	<b>F17-4848-00</b>
54"	<b>F17-4854-00</b>
60"	<b>F17-4860-00</b>
66"	<b>F17-4866-00</b>
72"	<b>F17-4872-00</b>

### Shelving Assembly Support Strut



Width	Part Number
12"	<b>K17-3612-00</b>
14"	<b>K17-3614-00</b>

### Shelving Uprights



#### Powder Coated Steel Uprights

Length	Single Slot Upright	Double Slot Upright	
24"	<b>F13M000224-00</b>	<b>F12M000224-00</b>	each
30"	<b>F13M000230-00</b>	<b>F12M000230-00</b>	each
36"	<b>F13M000236-00</b>	<b>F12M000236-00</b>	each
42"	<b>F13M000242-00</b>	<b>F12M000242-00</b>	each
48"	<b>F13M000248-00</b>	<b>F12M000248-00</b>	each
72"	<b>F13M000272-00</b>	<b>F12M000272-00</b>	each
79"	<b>F13M000279-00</b>	<b>F12M000279-00</b>	each
84"	<b>F13M000284-00</b>	<b>F12M000284-00</b>	each

#### Stainless Steel Uprights

Length	Single Slot Upright	Double Slot Upright	
24"	<b>F13S000224-00</b>	<b>F12S000224-00</b>	each
30"	<b>F13S000230-00</b>	<b>F12S000230-00</b>	each
36"	<b>F13S000236-00</b>	<b>F12S000236-00</b>	each
42"	<b>F13S000242-00</b>	<b>F12S000242-00</b>	each
48"	<b>F13S000248-00</b>	<b>F12S000248-00</b>	each
72"	<b>F13S000272-00</b>	<b>F12S000272-00</b>	each
79"	<b>F13S000279-00</b>	<b>F12S000279-00</b>	each
84"	<b>F13S000284-00</b>	<b>F12S000284-00</b>	each

Designed to be used with Adjustable Module Shelving Assemblies to create wall mounted shelving.

# Center Bench Shelving Assemblies

## Specifications:

Pass-Thru Shelf Assemblies for 12" Island Modules are six inches deeper than the Shelf Bracket to allow the Shelf to pass through the module core and form a continuous surface. These shelves may also be moved forward and mounted to form


a 2" cord gap between the shelves. Shelf Brackets are 11 gauge steel, adjustable on 1" centers.

**Painted Steel Shelves** – 1" thick, formed from 16 gauge cold rolled steel.


**Veneer Plywood Shelves** – 1" thick hardwood plywood with matching veneer on all exposed surfaces. Available in Oak, Maple or Beech.

**Phenolic Resin Shelves** – 1" thick phenolic resin.

## Pass-Thru Shelf Assemblies



Pass-Thru Shelves



Pass-Thru Shelves with Cord Gap


	Length	Painted Steel	Veneer Plywood	Phenolic Resin
<b>15" DEEP SHELF</b> 9" Deep Bracket	24"	<b>F21M0924-60</b>	<b>F21W0924-6_</b>	<b>F21T0924-60</b>
	30"	<b>F21M0930-60</b>	<b>F21W0930-6_</b>	<b>F21T0930-60</b>
	36"	<b>F21M0936-60</b>	<b>F21W0936-6_</b>	<b>F21T0936-60</b>
	42"	<b>F21M0942-60</b>	<b>F21W0942-6_</b>	<b>F21T0942-60</b>
	48"	<b>F21M0948-60</b>	<b>F21W0948-6_</b>	<b>F21T0948-60</b>
	54"	<b>F21M0954-60</b>		
	60"	<b>F21M0960-60</b>		
	66"	<b>F21M0966-60</b>		
<b>18" DEEP SHELF</b> 12" Deep Bracket	24"	<b>F21M1224-60</b>	<b>F21W1224-6_</b>	<b>F21T1224-60</b>
	30"	<b>F21M1230-60</b>	<b>F21W1230-6_</b>	<b>F21T1230-60</b>
	36"	<b>F21M1236-60</b>	<b>F21W1236-6_</b>	<b>F21T1236-60</b>
	42"	<b>F21M1242-60</b>	<b>F21W1242-6_</b>	<b>F21T1242-60</b>
	48"	<b>F21M1248-60</b>	<b>F21W1248-6_</b>	<b>F21T1248-60</b>
	54"	<b>F21M1254-60</b>		
	60"	<b>F21M1260-60</b>		
	66"	<b>F21M1266-60</b>		
<b>24" DEEP SHELF</b> 18" Deep Bracket	24"	<b>F21M1824-60</b>	<b>F21W1824-6_</b>	<b>F21T1824-60</b>
	30"	<b>F21M1830-60</b>	<b>F21W1830-6_</b>	<b>F21T1830-60</b>
	36"	<b>F21M1836-60</b>	<b>F21W1836-6_</b>	<b>F21T1836-60</b>
	42"	<b>F21M1842-60</b>	<b>F21W1842-6_</b>	<b>F21T1842-60</b>
	48"	<b>F21M1848-60</b>	<b>F21W1848-6_</b>	<b>F21T1848-60</b>
	60"	<b>F21M1860-60</b>		

To specify Veneer Plywood Species replace blank at end of part number with number that indicates material required.

### Veneer Plywood Material Options

- 1 = Red Oak                      5 = FSC® Red Oak
- 2 = White Maple                6 = FSC® White Maple
- 7 = Steamed Beech              8 = FSC® St'md Beech

## Shelf Retaining Lips



Length	Shelf Lip	
24"	<b>F63M0224-00</b>	each
30"	<b>F63M0230-00</b>	each
36"	<b>F63M0236-00</b>	each
42"	<b>F63M0242-00</b>	each
48"	<b>F63M0248-00</b>	each
54"	<b>F63M0254-00</b>	each
60"	<b>F63M0260-00</b>	each
66"	<b>F63M0266-00</b>	each
72"	<b>F63M0272-00</b>	each

Designed to be added to Pass-thru and Wall Shelf Assemblies. Steel Shelf Lips may be mounted to steel,

veneer plywood, or phenolic resin shelves with supplied screws.

# Wall Shelving Assemblies

## Specifications:

**Vertical Uprights** are 1 1/4" deep and either 79" or 42" high. They are manufactured with a double-slotted face with slots on 1" centers. They are designed to be mounted between two Support Rails to hold accurate spacing.


**Support Rails** are available in 8 foot lengths and may be cut to the required length. Support Rails are notched on 6" centers to accept the Vertical Uprights.

**Support Rail Covers** are required to finish the assembly and must be

ordered in pairs to match the Vertical Upright spacing.

**Optional Center Rails** may be ordered to provide a horizontal member at counter top height. They should be sized to match the Vertical Upright spacing.

## Full Height Wall Suspension System



### Vertical Upright

<b>F10M0279-00</b>	79"H	pair
--------------------	------	------

### Support Rail

<b>F11M0196-00</b>	96"L	each
--------------------	------	------


### Support Rail Covers

<b>F11M0324-00</b>	24"L	pair
<b>F11M0330-00</b>	30"L	pair
<b>F11M0336-00</b>	36"L	pair
<b>F11M0342-00</b>	42"L	pair
<b>F11M0348-00</b>	48"L	pair
<b>F11M0354-00</b>	54"L	pair
<b>F11M0360-00</b>	60"L	pair
<b>F11M0366-00</b>	66"L	pair
<b>F11M0372-00</b>	72"L	pair

### Optional Center Rail Assembly

<b>F11M0224-00</b>	24"L	each
<b>F11M0230-00</b>	30"L	each
<b>F11M0236-00</b>	36"L	each
<b>F11M0242-00</b>	42"L	each
<b>F11M0248-00</b>	48"L	each
<b>F11M0254-00</b>	54"L	each
<b>F11M0260-00</b>	60"L	each
<b>F11M0266-00</b>	66"L	each
<b>F11M0272-00</b>	72"L	each

## Above Countertop Wall Suspension System



### Vertical Upright

<b>F10M0242-00</b>	42"H	pair
--------------------	------	------

### Support Rail

<b>F11M0196-00</b>	96"L	each
--------------------	------	------

### Support Rail Covers

<b>F11M0324-00</b>	24"L	pair
<b>F11M0330-00</b>	30"L	pair
<b>F11M0336-00</b>	36"L	pair
<b>F11M0342-00</b>	42"L	pair
<b>F11M0348-00</b>	48"L	pair
<b>F11M0354-00</b>	54"L	pair
<b>F11M0360-00</b>	60"L	pair
<b>F11M0366-00</b>	66"L	pair
<b>F11M0372-00</b>	72"L	pair



# Wall Shelving Assemblies

## Specifications:

Adjustable Wall Shelf Assemblies consist of a pair of 11 gauge steel Shelf Brackets and a Shelf. The Shelf Bracket engages the inner slot of the double slotted Wall Vertical Upright and is adjustable on 1" centers. Shelves are available in Painted



Steel, Wood veneer plywood, or Phenolic Resin.

**Painted Steel Shelves** – 1" thick, formed from 16 gauge cold rolled steel.

**Veneer Plywood Shelves** – 1" thick hardwood plywood with matching veneer on all exposed surfaces. Available in Oak, Maple or Beech.

**Phenolic Resin Shelves** – 1" thick phenolic resin.

## Adjustable Wall Shelf Assemblies



### Optional Shelf Locking Tab

Sliding tab mounted on Shelf Bracket locks Shelf into upright to keep it from being accidentally removed.

To specify, add an **K** to the end of the Shelf Assembly part number.

Example: **F21M1236-00K**

	Length	Painted Steel	Veneer Plywood	Phenolic Resin
<b>6" DEEP SHELF</b> 6" Deep Bracket	24"	<b>F21M0624-00</b>	<b>F21W0624-0_</b>	<b>F21T0624-00</b>
	30"	<b>F21M0630-00</b>	<b>F21W0630-0_</b>	<b>F21T0630-00</b>
	36"	<b>F21M0636-00</b>	<b>F21W0636-0_</b>	<b>F21T0636-00</b>
	42"	<b>F21M0642-00</b>	<b>F21W0642-0_</b>	<b>F21T0642-00</b>
	48"	<b>F21M0648-00</b>	<b>F21W0648-0_</b>	<b>F21T0648-00</b>
	54"	<b>F21M0654-00</b>		
	60"	<b>F21M0660-00</b>		
	66"	<b>F21M0666-00</b>		
<b>9" DEEP SHELF</b> 9" Deep Bracket	24"	<b>F21M0924-00</b>	<b>F21W0924-0_</b>	<b>F21T0924-00</b>
	30"	<b>F21M0930-00</b>	<b>F21W0930-0_</b>	<b>F21T0930-00</b>
	36"	<b>F21M0936-00</b>	<b>F21W0936-0_</b>	<b>F21T0936-00</b>
	42"	<b>F21M0942-00</b>	<b>F21W0942-0_</b>	<b>F21T0942-00</b>
	48"	<b>F21M0948-00</b>	<b>F21W0948-0_</b>	<b>F21T0948-00</b>
	54"	<b>F21M0954-00</b>		
	60"	<b>F21M0960-00</b>		
	66"	<b>F21M0966-00</b>		
<b>12" DEEP SHELF</b> 12" Deep Bracket	24"	<b>F21M1224-00</b>	<b>F21W1224-0_</b>	<b>F21T1224-00</b>
	30"	<b>F21M1230-00</b>	<b>F21W1230-0_</b>	<b>F21T1230-00</b>
	36"	<b>F21M1236-00</b>	<b>F21W1236-0_</b>	<b>F21T1236-00</b>
	42"	<b>F21M1242-00</b>	<b>F21W1242-0_</b>	<b>F21T1242-00</b>
	48"	<b>F21M1248-00</b>	<b>F21W1248-0_</b>	<b>F21T1248-00</b>
	54"	<b>F21M1254-00</b>		
	60"	<b>F21M1260-00</b>		
	66"	<b>F21M1266-00</b>		
<b>18" DEEP SHELF</b> 18" Deep Bracket	24"	<b>F21M1824-00</b>	<b>F21W1824-0_</b>	<b>F21T1824-00</b>
	30"	<b>F21M1830-00</b>	<b>F21W1830-0_</b>	<b>F21T1830-00</b>
	36"	<b>F21M1836-00</b>	<b>F21W1836-0_</b>	<b>F21T1836-00</b>
	42"	<b>F21M1842-00</b>	<b>F21W1842-0_</b>	<b>F21T1842-00</b>
	48"	<b>F21M1848-00</b>	<b>F21W1848-0_</b>	<b>F21T1848-00</b>
	54"	<b>F21M1854-00</b>		
	60"	<b>F21M1860-00</b>		
	<b>24" DEEP SHELF</b> 24" Deep Bracket	24"	<b>F21M2424-00</b>	<b>F21W2424-0_</b>
30"		<b>F21M2430-00</b>	<b>F21W2430-0_</b>	<b>F21T2430-00</b>
36"		<b>F21M2436-00</b>	<b>F21W2436-0_</b>	<b>F21T2436-00</b>
42"		<b>F21M2442-00</b>	<b>F21W2442-0_</b>	<b>F21T2442-00</b>
48"		<b>F21M2448-00</b>	<b>F21W2448-0_</b>	<b>F21T2448-00</b>
54"		<b>F21M2454-00</b>		
60"		<b>F21M2460-00</b>		

To specify Veneer Plywood Species replace blank at end of part number with number that indicates material required.

### Veneer Plywood Material Options

- 1 = Red Oak                      5 = FSC® Red Oak
- 2 = White Maple                6 = FSC® White Maple
- 7 = Steamed Beech              8 = FSC® St'md Beech

# Reagent Rack & Shelving Systems

**Reagent Rack Uprights** are available in two heights and two depths. All are 1½" thick and available in either powder coated steel to match or complement the casework color or in type 304 stainless steel. 12" deep uprights are designed for


use on Center Benches with a 14" deep shelf, 10" deep uprights are designed for use on Wall Benches with a 12" deep shelf, but both can be used on either bench.

Installation hardware is included with each upright.


Spacing between upright should not exceed 48".

**Plastic Laminate Shelves** are 1" thick particle board faced on top and all four edges with black plastic laminate. A backer sheet covers the bottom.


## Typical Reagent Rack Assemblies



**Single Tier Center Bench Reagent Rack** using 18" high by 12" deep uprights and 14" deep shelving.




**Double Tier Center Bench Reagent Rack** using 18" high by 12" deep and 12" by 12" deep uprights and 14" deep shelving.



**Single Tier Wall Bench Reagent Rack** using 18" high by 10" deep uprights and 12" deep shelving.

## Reagent Rack Uprights


**18" High x 12" Deep**



18" H x 12" D x 1½" T

**RRU-1812-0000** P.C. Steel  
**RRU-1812-0030** St. Steel


**12" High x 12" Deep**



12" H x 12" D x 1½" T

**RRU-1212-0000** P.C. Steel  
**RRU-1212-0030** St. Steel


**18" High x 10" Deep**



18" H x 10" D x 1½" T

**RRU-1810-0000** P.C. Steel  
**RRU-1810-0030** St. Steel

**12" High x 10" Deep**




12" H x 10" D x 1½" T

**RRU-1210-0000** P.C. Steel  
**RRU-1210-0030** St. Steel


## Plastic Laminate Shelves

**14" Deep Shelf**



**PL-4814-00** 14" Deep x 48" Long  
**PL-7214-00** 14" Deep x 72" Long

**12" Deep Shelf**



**PL-4812-00** 12" Deep x 48" Long  
**PL-7212-00** 12" Deep x 72" Long

# Reagent Rack & Shelving Systems

**Fixed Wall Mounted Shelf Brackets** are designed to mount directly to the wall and furnished without mounting hardware.

**Adjustable Wall Mounted Shelf Brackets** are designed to mount to KF-0170-04 Shelving Uprights and


include mounting bolts and nuts.

**KF-0170-04 Shelving Uprights** are designed to hold a maximum of three shelves.


Brackets and Uprights are finished in Kewaunee's unique corrosion

resistant powder coat finish to match or complement the casework color.


## Typical Reagent Rack Assemblies



Fixed Wall Shelving




Adjustable Wall Shelving



Adjustable Wall Shelving


## Flat Top Shelf Bracket    Lipped Shelf Bracket    Book End Bracket    Shelving Upright



**Fixed Shelf Brackets**

	bracket size	shelf size
<b>M-0031-0A</b>	11" x 8"	12"
<b>M-0030-0A</b>	15" x 10"	16"
<b>M-0032-0A</b>	17" x 12"	18"


part number is for single bracket



**Fixed Shelf Brackets**

	bracket size	shelf size
<b>M-0038-0A</b>	12 <sup>1</sup> / <sub>8</sub> " x 6"	12"
<b>M-0044-0A</b>	16 <sup>1</sup> / <sub>8</sub> " x 6"	16"
<b>M-0039-0A</b>	18 <sup>1</sup> / <sub>8</sub> " x 6"	18"

part number is for single bracket



**Fixed Shelf Brackets**


	bracket size	shelf size
<b>M-0060-0A</b>	12" x 6"	12"
<b>M-0061-0A</b>	16" x 6"	16"
<b>M-0062-0A</b>	18" x 6"	18"

part number is for pair of brackets



**KF-0170-04**                      48" L


High-strength, cold-rolled steel shelf support upright allows adjustable shelf bracket to be positioned along 48" length for maximum use of space. A pair of uprights holds a maximum of three shelf brackets.



**Adjustable Shelf Brackets**

	bracket size	shelf size
<b>M-0031-AD</b>	11" x 8"	12"
<b>M-0030-AD</b>	15" x 10"	16"
<b>M-0032-AD</b>	17" x 12"	18"


part number is for single bracket



**Adjustable Shelf Brackets**

	bracket size	shelf size
<b>M-0038-AD</b>	12 <sup>1</sup> / <sub>8</sub> " x 6"	12"
<b>M-0044-AD</b>	16 <sup>1</sup> / <sub>8</sub> " x 6"	16"
<b>M-0039-AD</b>	18 <sup>1</sup> / <sub>8</sub> " x 6"	18"

part number is for single bracket




**Adjustable Shelf Brackets**

	bracket size	shelf size
<b>M-0060-AD</b>	12" x 6"	12"
<b>M-0061-AD</b>	16" x 6"	16"
<b>M-0062-AD</b>	18" x 6"	18"

part number is for pair of brackets

# Pegboards


## Kemresin Laboratory Pegboards



**Pegboard without Drip Trough**  
29" High x 1" Thick x 20" Long

	Pegboard Color	Peg Color
<b>X-020014-BK</b>	Black	Black
<b>X-020014-GR</b>	Grey	White
<b>X-020014-PT</b>	Putty	White
<b>X-020014-SL</b>	Slate	White


**Peg Size:** 1/2" diameter x 5" long  
**Peg Material:** Polypropylene  
**Number of Pegs:** 32



**Pegboard without Drip Trough**  
29" High x 1" Thick x 32" Long

	Pegboard Color	Peg Color
<b>X-020015-BK</b>	Black	Black
<b>X-020015-GR</b>	Grey	White
<b>X-020015-PT</b>	Putty	White
<b>X-020015-SL</b>	Slate	White


**Peg Size:** 1/2" diameter x 5" long  
**Peg Material:** Polypropylene  
**Number of Pegs:** 53



**Pegboard with Drip Trough**  
30" High x 1" Thick x 20" Long

	Pegboard Color	Peg Color
<b>X-020016-BK</b>	Black	Black
<b>X-020016-GR</b>	Grey	White
<b>X-020016-PT</b>	Putty	White
<b>X-020016-SL</b>	Slate	White

**Peg Size:** 1/2" diameter x 5" long  
**Peg Material:** Polypropylene  
**Number of Pegs:** 32  
**Trough Size:** 1" high x 2 1/4" deep



**Pegboard with Drip Trough**  
30" High x 1" Thick x 32" Long

	Pegboard Color	Peg Color
<b>X-020017-BK</b>	Black	Black
<b>X-020017-GR</b>	Grey	White
<b>X-020017-PT</b>	Putty	White
<b>X-020017-SL</b>	Slate	White

**Peg Size:** 1/2" diameter x 5" long  
**Peg Material:** Polypropylene  
**Number of Pegs:** 53  
**Trough Size:** 1" high x 2 1/4" deep



**Stainless Steel Drip Trough**  
for use on  
**F-2939-00** X-020014 Pegboard


**Size:** 1 1/2"H x 2 1/2"W x 20"L  
**Material:** Type 304 Stainless Steel  
**Includes:** 1/2" OD drain outlet



**Stainless Steel Drip Trough**  
for use on  
**F-2940-00** X-020015 Pegboard

**Size:** 1 1/2"H x 2 1/2"W x 32"L  
**Material:** Type 304 Stainless Steel  
**Includes:** 1/2" OD drain outlet


## Type 304 Stainless Steel Laboratory Pegboards



**Pegboard without Drip Trough**  
30" High x 3/4" Thick x 24" Long

<b>X-020008</b>	with Open Back	
<b>X-020010</b>	with Finished Back	


**Peg Size:** 1/2" diameter x 6" long  
**Peg Material:** Polypropylene  
**Peg Color:** White  
**Number of Pegs:** 32



**Pegboard without Drip Trough**  
30" High x 3/4" Thick x 30" Long

<b>X-020009</b>	with Open Back	
<b>X-020011</b>	with Finished Back	


**Peg Size:** 1/2" diameter x 6" long  
**Peg Material:** Polypropylene  
**Peg Color:** White  
**Number of Pegs:** 50



**Pegboard with Drip Trough**  
30" High x 3/4" Thick x 24" Long

<b>X-020004</b>	with Open Back	
<b>X-020006</b>	with Finished Back	

**Peg Size:** 1/2" diameter x 6" long  
**Peg Material:** Polypropylene  
**Peg Color:** White  
**Number of Pegs:** 32




**Pegboard with Drip Trough**  
30" High x 3/4" Thick x 30" Long

<b>X-020009</b>	with Open Back	
<b>X-020011</b>	with Finished Back	

**Peg Size:** 1/2" diameter x 6" long  
**Peg Material:** Polypropylene  
**Peg Color:** White  
**Number of Pegs:** 50

# Miscellaneous Accessories


## Writing Boards & Tackboards



**Black Chalkboard**  
**F-5336-00** 18<sup>1</sup>/<sub>4</sub>"H x 24<sup>1</sup>/<sub>4</sub>"L  
**F-5337-00** 24<sup>1</sup>/<sub>4</sub>"H x 36<sup>1</sup>/<sub>4</sub>"L  
**F-5338-00** 36<sup>1</sup>/<sub>4</sub>"H x 48<sup>1</sup>/<sub>4</sub>"L

---


**Surface:** washable Duracite  
**Frame Material:** anodized aluminum  
**Trough Material:** anodized aluminum



**White Magnetic Markerboard**  
**F-5342-00** 18<sup>1</sup>/<sub>4</sub>"H x 24<sup>1</sup>/<sub>4</sub>"L  
**F-5343-00** 24<sup>1</sup>/<sub>4</sub>"H x 36<sup>1</sup>/<sub>4</sub>"L  
**F-5344-00** 36<sup>1</sup>/<sub>4</sub>"H x 48<sup>1</sup>/<sub>4</sub>"L

---


**Surface:** 1/4" white magnetic  
**Frame Material:** grey mylar vinyl



**Green Chalkboard**  
**F-5333-00** 18<sup>1</sup>/<sub>4</sub>"H x 24<sup>1</sup>/<sub>4</sub>"L  
**F-5340-00** 24<sup>1</sup>/<sub>4</sub>"H x 36<sup>1</sup>/<sub>4</sub>"L  
**F-5341-00** 36<sup>1</sup>/<sub>4</sub>"H x 48<sup>1</sup>/<sub>4</sub>"L

---

**Surface:** washable Duracite  
**Frame Material:** anodized aluminum  
**Trough Material:** anodized aluminum




**Cork Tackboard**  
**F-5418-24** 18"H x 24"L  
**F-5424-36** 24"H x 36"L  
**F-5436-48** 36"H x 48"L

---

**Surface:** Resilient natural cork  
**Frame Material:** anodized aluminum

## Kemresin Balance Table




30" High x 24" Deep x 35" Long  
**K7-302435-BK** Black  
**K7-302435-GR** Grey  
**K7-302435-PT** Putty  
**K7-302435-SL** Slate

---

**Work Top:** 2" thick Kemresin  
**Work Top Size:** 2"T x 24"D x 35"L  
**Pedestal Legs:** 2" thick Kemresin  
**Total Weight:** 300 pounds

## Marvel Refrigerators & Freezers



**L-110** Refrigerator/Freezer  
**Size:** 34<sup>1</sup>/<sub>2</sub>" H x 24" D x 23<sup>7</sup>/<sub>8</sub>" L  
 6.1 cubic feet  
**Door Swing:** Right Hand  
 (User changeable with swing reversing kit)  
**Shelves:** (3) Adjustable  
**Service:** Single-phase,  
 60 hertz, 120 volt AC

---

**L-105** Refrigerator  
 Same as above except  
 without freezer compartment.


---

**L-115** Freezer  
 Same as above except  
 entire unit is a freezer.

## Gillis Stainless Steel Cart Lund Key Cases



**F-5349-00**  
**Overall Size:** 33<sup>3</sup>/<sub>8</sub>"H x 19"D x 31"L  
**Shelf Size:** (3) 18"D x 27"L  
 (1" high lip on 3 sides of  
 top & middle shelf)  
**Capacity:** 400 pounds  
**Casters:** (4) Black rubber  
 4" Diameter – Swivel




**F-5332-00**  
**Height:** 11"  
**Depth:** 2<sup>1</sup>/<sub>2</sub>"  
**Length:** 17<sup>5</sup>/<sub>16</sub>"  
**Key Capacity:** 90 Keys

---

**F-5332-01**  
**Height:** 24<sup>3</sup>/<sub>4</sub>"  
**Depth:** 2<sup>1</sup>/<sub>2</sub>"  
**Length:** 17<sup>5</sup>/<sub>16</sub>"  
**Key Capacity:** 200 Keys

# Miscellaneous Accessories

## Pipe Space Kemstruts




Floor mounted support to provide stability to rear overhang of fume hood tops or worksurfaces, and for mounting piping and drainlines. Adjustable to match variances in height.

**K12-**    **-0A**

Available Heights		Available Depths	
<i>adjustable</i>		<i>depth</i>	<i>depth</i>
<b>29</b> = 27" - 31"		<b>05</b> = 5"	<b>13</b> = 13"
<b>32</b> = 31" - 35"		<b>06</b> = 6"	<b>14</b> = 14"
<b>35</b> = 34" - 38"		<b>07</b> = 7"	<b>15</b> = 15"
		<b>08</b> = 8"	<b>16</b> = 16"
		<b>09</b> = 9"	<b>17</b> = 17"
		<b>10</b> = 10"	<b>18</b> = 18"
		<b>11</b> = 11"	<b>19</b> = 19"
		<b>12</b> = 12"	<b>20</b> = 20"

## Corner Kemstruts




**K-0020-1A** 28" - 30" H  
**K-0016-1A** 33" - 36" H

Floor and wall mounted support system for adding stability to corner tops in L-shaped furniture assemblies. Adjustable to match variances in height. Includes two support angles to attach to adjoining base cabinets.


**NOTE:** Kemstrut Corner Top Support Assemblies are not required when a corner base cabinet is used.

## Kemstrut Pipe Hangers



**K-0005-0A** Drainline Support  
**M-0051-0A** Pipe Clamp for  $\frac{3}{8}$ " -  $\frac{1}{2}$ " IPS  
**M-0052-0A** Pipe Clamp for  $\frac{3}{4}$ " - 1" IPS

## Support Angles




**L-2034-06** 6" L  
**L-2034-12** 12" L  
**L-2034-18** 18" L  
**L-2034-24** 24" L  
**L-2034-48** 48" L  
**L-2034-96** 96" L

Formed from 16 gauge cold-rolled steel, with  $1\frac{1}{2}$ " sides.

## Pipe Enclosures

### Wall Bench Enclosure




Pipe Enclosures conceal vertical service lines coming down from the ceiling to the worksurface. One side is removable for easy access with no exposed fasteners. Designed to sit on top of curb on worksurface and penetrate drop ceiling.

**U01W**     **-00\_0** Wall Bench Enclosure  
**U02W**     **-00\_0** Center Bench Enclosure


Available Heights	Available Depths	Available Lengths	Available Materials
<b>60</b> = 60"	<b>05</b> = 5"	<b>12</b> = 12"	<b>1</b> = Plain Sliced Red Oak
<b>66</b> = 66"	<b>06</b> = 6"	<b>18</b> = 18"	<b>2</b> = Plain Sliced White Maple
<b>72</b> = 72"	<b>07</b> = 7"	<b>24</b> = 24"	<b>3</b> = NAUF Red Oak
<b>78</b> = 78"	<b>08</b> = 8"	<b>30</b> = 30"	<b>4</b> = NAUF White Maple
<b>84</b> = 84"	<b>09</b> = 9"	<b>36</b> = 36"	<b>5</b> = FSC/NAUF Red Oak
<b>90</b> = 90"	<b>10</b> = 10"	<b>42</b> = 42"	<b>6</b> = FSC/NAUF White Maple
<b>96</b> = 96"		<b>48</b> = 48"	<b>7</b> = Steamed Beech
			<b>8</b> = FSC/NAUF Steamed Beech

### Center Bench Enclosure



# Miscellaneous Accessories


## Rod Assembly



**4740-00** Rod Assembly

Consists of two 0553-00 upright rods, one 0557-00 cross rod and two 0550-00 rod clamps.


## Support Rods



- 0549-18** 1/2" Dia. x 18" H
- 0549-24** 1/2" Dia. x 24" H
- 0549-36** 1/2" Dia. x 36" H
- 0553-00** 3/4" Dia. x 36" H

Duralumin support rods are available in 1/2" and 3/4" diameter and have a tapered end to fit 0548-00 rod sockets.


## Cross Rods



- 0554-00** 1/2" Dia. x 39" L
- 0556-00** 1/2" Dia. x 18" L
- 0557-00** 3/4" Dia. x 39" L

Duralumin cross rods are available in 1/2" and 3/4" diameter and have rounded ends


## Rod Clamp



**0550-00** Rod Clamp

Made of cast aluminum. Designed to clamp horizontal cross rods on vertical upright rods. Can be used for either 1/2" or 3/4" diameter rods


## Rod Sockets



**0548-00** Rod Socket

Made of chrome-plated brass. Use with 1/2" and 3/4" diameter upright rods. 1" diameter x 17/8" long threaded shank. Complete with locknut and washer.

## Label Holder



**F-0260-00** Label Holder

Aluminum with adhesive back for easy installation. Holds 1 1/8" H x 2 1/2" L identification card.

## Base Molding



**F-0862-04** 48" L

4" high black base molding. Made of pliable material to match variances in floor surface. Comes with high-strength adhesive applied to back.

**Note: Specify number of 4' sections required.**


## Corner Guard



**F-1368-00** Corner Guard

4" high black corner guard to finish off corner joints of base molding.


## LED Task Lighting



- | Part Number        | Length           |
|--------------------|------------------|
| <b>X-030001-00</b> | 300mm – 11 7/8"  |
| <b>X-030002-00</b> | 500mm – 19 5/8"  |
| <b>X-030003-00</b> | 1000mm – 39 3/8" |

**LED Task Lighting** consists of a thin-profile, white, LED strip, 9.5mm high by 33mm deep with a touch dimmer switch. It is UL listed and furnished with 120 VAC to 24 VDC power supply with ten foot power cord. The light includes a magnetic mounting strip for mounting to steel shelving and integral tabs for mounting to wood and plastic laminate.

## Tote Trays



- OP-6140** 14" x 13 1/4" x 3 1/8" for use in G94 base cabinets
- OP-6141** 16 3/8" x 11" x 2 1/2" for use in S62 & S64 storage cabinets
- OP-6551** 19" x 13 3/4" x 3 1/2" for use in G86 base cabinets

Made of tan high strength plastic.

# General Plumbing Information

## Recommended Supply Line Pressures

Gas ..... 0.216 psi  
 Air (Ordinary) ..... 15 psi  
 Air (Air Powered Stirrer Operation) ..... 40 psi  
 Vacuum (Ordinary) ..... 18"-22" of Mercury  
 Vacuum (High) ..... 28.5" of Mercury  
 Hot and Cold Water ..... Local Pressure  
 Steam (Low Pressure) ..... 0-15 psi as required  
 Steam (High Pressure) ..... 15-30 psi as required

## Service & Rough-in Abbreviations

**A** – Air  
**E** – Electric - 120 volt AC  
**G** – Gas  
**N** – Nitrogen  
**O** – Oxygen  
**S** – Steam  
**V** – Vacuum  
**W** – Waste  
**CW** – Cold Water  
**DW** – Distilled Water  
**HV** – Hood Vent  
**HW** – Hot Water  
**WV** – Waste Vent

## Service Fitting Index Color Coding

When ordering service fittings, specify service required, so that the proper index button(s) will be supplied.

**Air** – Orange  
**Cold Water** – Dark Green  
**Distilled Water** – White  
**Gas** – Dark Blue  
**Hot Water** – Red  
**Hydrogen** – Pink  
**Nitrogen** – Brown  
**Oxygen** – Light Green  
**Special Gases** – Light Blue  
**Steam** – Black  
**Vacuum** – Yellow

## Recommended Pipe & Drain Line Sizes

Service	Number of Outlets	Pipe Size	Pipe Material
Steam	1-8	3/8"	In accordance with governing codes
	9-16	1/2"	
	17-24	3/4"	
*Water	1-10	1/2"	In accordance with governing codes
	11-18	3/4"	
	19-24	1"	
Gas	1-12	3/8"	In accordance with governing codes
	13-18	1/2"	
	19-24	3/4"	
	25-30	1"	
Air or Vacuum	1-8	3/8"	In accordance with governing codes
	9-16	1/2"	
	17-24	3/4"	
	<b>Fittings or Filter Pumps</b>	<b>Pipe Size</b>	<b>Pipe Material</b>
Drain Lines	1-6 Fittings or 4 Pumps	1 1/2"	In accordance with governing codes
	7-10 Fittings or 8 Pumps	2"	


\* If table is equipped with a series of filter pumps (aspiration), next larger pipe size should be used.

**NOTE:** For information not supplied, or if additional details are required, contact your local Kewaunee Representative.

## Open Reagent Rack Horizontal Pipe Spacing


Water pipe lines should always be the lowest and a minimum of 9" above the working surface. Electric conduit should always be the uppermost service line. (Insulated pipe will require a greater spacing than indicated below.)

Pipe Diameters	Minimum Vertical Spacing Between Centerlines of 2 Service Pipe Lines
3/8" and 3/8"	2"
3/8" and 1/2"	2"
1/2" and 1/2"	2"
1/2" and 3/4"	2 1/2"
3/4" and 3/4"	2 1/2"
3/4" and 1"	3"
1" and 1"	3"



## Minimum Spacing Between Fittings


These drawings show the allowable minimum spacing between service and electrical fittings on laboratory bench tops and ledges.






# Laboratory Service Fittings


## Laboratory Ball Valves (for use with air, vacuum or gases – not to be used for liquids, steam or oxygen)




**W-0242 Gas, Air or Vacuum Valve** is chrome plated brass. Valve provides on/off control of gas, air or vacuum at pressures up to 75 PSI.  $\frac{3}{8}$ " I.P.S. male inlet. **Specify service desired.**




**W-0286 Double Gas, Air or Vacuum Fitting (curb or panel mounted)** is chrome plated brass. Consists of two-way 70° outlet panel flange and two laboratory ball valves. Complete with  $\frac{3}{8}$ " I.P.S. tank nipple, locknut and washer. **Specify service desired.**




**W-0287 Gas, Air or Vacuum Fitting (curb or panel mounted)** is chrome plated brass. Consists of a single outlet panel flange and a laboratory ball valve. Complete with a  $\frac{3}{8}$ " I.P.S. tank nipple, locknut and washer. **Specify service desired.**




**W-0260 Gas, Air or Vacuum Deck Mounted Fitting** is chrome plated brass. Consists of a single outlet turret with a laboratory ball valve. Complete with a  $\frac{3}{8}$ " I.P.S. tank nipple, locknut and washer. **Specify service desired.**



**W-0263 Double Gas, Air or Vacuum Deck Mounted Fitting** Same as No. W-0260-00 except consists of a 2-way 90° outlet turret and two laboratory ball valves. Complete with a  $\frac{3}{8}$ " I.P.S. tank nipple, locknut and washer. **Specify service desired.**




**W-0266 Double Gas, Air or Vacuum Deck Mounted Fitting** Same as No. W-0260-00, except consists of a 2-way 180° outlet turret and two laboratory ball valves. Complete with a  $\frac{3}{8}$ " I.P.S. tank nipple, locknut and washer. **Specify service desired.**




**W-0272 Quadruple Gas, Air or Vacuum Deck Mounted Fitting** Same as W-0260-00, except consists of a 4-way outlet turret and four valves. Complete with a  $\frac{3}{8}$ " I.P.S. tank nipple, locknut and washer. **Specify service desired.**


## Needle Point Cock Valves (for use with air, vacuum or gases – not to be used for liquids or steam)



**W-0243 Gas, Air or Vacuum Valve** is chrome plated brass. Single needle point valve with integral straight serrated hose ends. Stainless steel floating cone and removable brass seat. Valve provides metering control of laboratory gases at pressures up to 75 PSI.  $\frac{3}{8}$ " I.P.S. male inlet. **Specify service desired.**




**W-0400 Double Gas, Air or Vacuum Fitting (curb or panel mounted)** with two needle point valves with integral serrated hose ends. Stainless steel floating cone and removable brass seat. Complete with a  $\frac{3}{8}$ " I.P.S. tank nipple, locknut and washer. **Specify service desired.**




**W-0401 Gas, Air or Vacuum Fitting (curb or panel mounted)** Same as No. W-0400-0\_ except single needle point valve. Complete with a  $\frac{3}{8}$ " I.P.S. tank nipple, locknut and washer. **Specify service desired.**

# Laboratory Service Fittings


## Needle Point Cock Valves (for use with air, vacuum or gases – not to be used for liquids or steam)




**W-0259 Single Needle Cock Deck Mounted Fitting** Turret Base is chrome plated brass. Single needle point valve with integral straight serrated hose end. Complete with a 3/8" I.P.S. tank nipple, locknut and washer.  
**Specify service desired.**



**W-0402 Double Needle Cock Deck Mounted Fitting** Turret Base is chrome plated brass. Two needle point valves with integral straight serrated hose ends at 90°. Complete with a 3/8" I.P.S. tank nipple, locknut and washer.  
**Specify service desired.**




**W-0403 Double Needle Cock Deck Mounted Fitting** Turret Base is chrome plated brass. Same as No. W-0402-00 except integral straight serrated hose ends are at 180°. Complete with a 3/8" I.P.S. tank nipple, locknut and washer.  
**Specify service desired.**




**W-0404 Quadruple Needle Cock Deck Mounted Fitting** Turret Base is chrome plated brass. Four needle point valves with integral straight serrated hose ends. Complete with a 3/8" I.P.S. tank nipple, locknut and washer.  
**Specify service desired.**

## Water Valve




**W-0306 Water Fitting (deck or panel mounted)** Made of chrome plated brass. Consists of a single outlet turret and a straight water valve with removable straight serrated hose connector. Complete with a 3/8" I.P.S. tank nipple, locknut and washer.

## Steam Valves




**W-0300 Steam Fitting (deck or panel mounted)** Made of chrome plated brass. Consists of a single outlet turret and a straight steam valve. Designed for maximum steam pressure of 20 psi. Complete with a 3/8" I.P.S. tank nipple, locknut and washer.




**W-0301 Steam Fitting (curb or panel mounted)** Made of chrome plated brass. Consists of a single outlet panel flange and a straight steam valve. Designed for maximum steam pressure of 20 psi. Complete with a 3/8" I.P.S. tank nipple, locknut and washer.

## Multiple Service Fitting



**W-0273 Deck Mounted Combination Water/Gas Fitting** Made of chrome plated brass. Water faucet has a 6" spread gooseneck with removable aerator outlet and disc type valve. Includes a two-way turret with two laboratory ball valves set at 180°. Complete with two 3/8" I.P.S. tank nipples, locknuts and washers.  
**(Specify service desired for ball valves.)**

**W-0273-VB Deck Mounted Combination Water/Gas Fitting** Same as W-0273 except furnished with vacuum breaker type gooseneck and removable straight serrated hose connector.  
**(Specify service desired for ball valves.)**




**W-0274-VB Deck Mounted Combination Water/Gas Fitting** Made of chrome plated brass. Water faucet has a 8" spread vacuum breaker type gooseneck with removable straight serrated hose connector and disc type valve. Includes a two-way turret with two laboratory ball valves set at 180°. Complete with two 3/8" I.P.S. tank nipples, locknuts and washers.  
**(Specify service desired for valves.)**


**W-0274 Deck Mounted Combination Water/Gas Fitting** Same as W-0274-8VB except furnished with standard gooseneck and removable aerator outlet.  
**(Specify service desired for valves.)**

# Laboratory Service Fittings


## Single Water Faucet




**W-0337-00 Deck Mounted Gooseneck Water Faucet**  
 Chrome plated brass with an 8" spread rigid gooseneck. Removable straight serrated hose connector. Valve is disc type with a renewable internal assembly. Faucet body is complete with 3/8" I.P.S. threaded mounting shank with locknut and washer.




**W-0337-0V Deck Mounted Vacuum Breaker Gooseneck Water Faucet**  
 Same as W-0337-00 except furnished with 8" spread rigid vacuum breaker type gooseneck.




**W-0337-06 Deck Mounted Gooseneck Water Faucet**  
 Same as W-0337-00 except furnished with 6" spread rigid gooseneck




**W-0337-6V Deck Mounted Vacuum Breaker Gooseneck Water Faucet**  
 Same as W-0337-0V except furnished with 6" spread rigid vacuum breaker type gooseneck.




**W-0350-00 Right Hand Curb Mounted Gooseneck Water Faucet**  
 Made of chrome plated brass. Has an 8" spread rigid gooseneck with removable straight serrated hose connector. Valve is disc type with renewable internal assembly. Flanged inlet is complete with 1/2" I.P.S. tank nipple, locknut and washer.




**W-0350-0V Right Hand Curb Mounted Vacuum Breaker Gooseneck Water Faucet**  
 Same as W-0350-00 except supplied with vacuum breaker type gooseneck.



**W-0351-00 Left Hand Curb Mounted Gooseneck Water Faucet**  
 Same as W-0350-00 except supplied with valve handle on left.




**W-0351-0V Left Hand Curb Mounted Vacuum Breaker Gooseneck Water Faucet**  
 Same as W-0350-0V except supplied with valve handle on left.




**W-0405-00 Portable Pump Gooseneck Faucet** for laboratory demonstration table. Made of chrome plated brass with 7<sup>5</sup>/<sub>8</sub>" spread, rigid, 1/4" I.P.S. gooseneck. Faucet/pump body is complete with an integral threaded mounting shank and flanged locknut. 5/8" O.D. hose connector inlet.

# Laboratory Service Fittings


## Mixing Faucet




**W-0340-00 Deck Mounted Hot and Cold Water Mixing Gooseneck Faucet** Made of chrome plated brass with an 8" spread rigid gooseneck. Removable straight serrated hose connector. Valves are disc type with renewable internal assemblies. Faucet body is complete with an integral threaded mounting shank and a flanged locknut. Furnished with 1/2" O.D. slip joint coupling nuts for 3/8" tubing.




**W-0340-0V Deck Mounted Hot and Cold Water Mixing Vacuum Breaker Gooseneck Faucet** Same as W-0340-00 except supplied with vacuum breaker type gooseneck.




**W-0340-06 Deck Mounted Hot and Cold Water Mixing Vacuum Breaker Gooseneck Faucet** Same as W-0340-00 except furnished with 6" spread rigid gooseneck




**W-0340-6V Deck Mounted Hot and Cold Water Mixing Vacuum Breaker Gooseneck Faucet** Same as W-0340-0V except furnished with 6" spread rigid vacuum breaker type gooseneck.




**W-0333-00 Deck Mounted Hot and Cold Water Mixing Swing Spout Faucet** Made of chrome plated brass. Has a 9" spread swinging spout with anti-splash outlet. Valves are disc type with renewable internal assemblies. Faucet body is complete with an integral threaded shank and a flanged locknut. Furnished with 1/2" O.D. slip joint coupling nuts for 3/8" tubing.



**W-0354-00 Curb Mounted Hot and Cold Water Mixing Gooseneck Faucet** Made of chrome plated brass. Has an 8" spread rigid gooseneck with removable straight serrated hose connector. Valves are disc type with renewable internal assemblies. Flanged inlets are complete with 3/8" I.P.S. tank nipples, locknuts and washers on 8" centers.




**W-0354-0V Curb Mounted Hot and Cold Water Mixing Vacuum Breaker Gooseneck Faucet** Same as No. W-0354-00 except supplied with vacuum breaker type gooseneck.



**W-0355-00 Curb Mounted Hot and Cold Water Mixing Swing Spout Faucet** Made of chrome plated brass. Has a 9" spread swinging spout with anti-splash outlet. Valves are disc type with renewable internal assemblies. Flanged inlets are complete with 3/8" I.P.S. tank nipples, locknuts and washers on 8" centers.

# Laboratory Service Fittings

## Mixing Faucet




**W-0361-00 Panel Mounted Hot and Cold Water Mixing Gooseneck Faucet** Chrome plated brass. Designed for mounting on facia type reagent rack plumbing box uprights. It has an 8" spread rigid gooseneck equipped with a removable straight serrated hose connector. Valves are disc type with renewable internal assemblies. Flanged inlets on 3" centers are complete with 3/8" I.P.S. tank nipples, locknuts and washers.

**W-0361-0V Panel Mounted Hot and Cold Water Mixing Vacuum Breaker Gooseneck Faucet** Same as W-0361-00 except supplied with vacuum breaker type gooseneck.

**W-0825-00 Floor Mounted Foot Control Double Pedal Mixing Valve** is chrome plated brass. Has twin pedals with self-closing valves for mixing hot and cold water. 1/2" I.P.S. female thread inlets and outlet at the rear. Inlets are on 2" centers, 7/8" above floor. Outlet is on center 1 1/2" above floor.

**W-0833-00 Deck Mounted Water Gooseneck Fitting** is chrome plated brass. Consisting of a 8" gooseneck with mounting flange and W-0532-00 aerator. 3/8" I.P.S. female thread inlet complete with tank nipple, locknut and washer. **Intended for use as an outlet for W-0825-00.**

## Faucet Accessories/Options




**W-0532-00 Aerator**  
Anti-splash aerator can be used in lieu of the removable straight serrated hose connectors supplied with the gooseneck water faucet. 3/8" I.P.S. male thread inlet.

**W-0510-00 Filter Pump (Aspirator)** 5 1/2" high overall. Has serrated hose connector side and bottom outlets. 3/8" I.P.S. male thread inlet.

**W-0535-00 In-Line Vacuum Breaker**  
Chrome plated. Mounts between the outlet end of a fitting (gooseneck, flange, etc.) and a removable serrated hose connector. 3/8" I.P.S. male thread inlet and female thread outlet. (Check plumbing codes for required location of vacuum breakers before ordering.)


# Laboratory Service Fittings

## Distilled Water Faucet



**W-0367-00 Right Hand Distilled Water (Curb or Panel Mounted) Faucet**


Chrome plated with a hot tin lined interior. It has a self-closing handle and a disc type valve with a renewable internal assembly. The 8" spread rigid gooseneck has a removable serrated hose end. 1/2" I.P.S. flanged female thread inlet complete with tank nipple, locknut and washer.




**W-0368-00 Deck Mounted Distilled Water Gooseneck Faucet**

Made of chrome plated brass. Has a hot tin lined interior. It has a self-closing handle and a disc type valve with renewable internal assembly. The 8" spread rigid gooseneck has a removable serrated hose end. Faucet body is complete with an integral threaded mounting shank and a flanged locknut. 1/2" I.P.S. male thread inlet with coupling nut and 1/4" I.P.S. tailpiece.


## Ultra-Pure Deionized Water Faucet



**W-0407-00 Distilled, Deionized and Ultra-Pure Water Faucet (deck mounted)** made of chrome plated brass. All polypropylene internal construction with 6" spread gooseneck. Furnished with 1/2" I.P.S. shank with locknut and washer. Inlet is 3/8" O.D. compression.




**W-0408-00 Distilled, Deionized and Ultra-Pure Water Faucet (panel mounted)** made of chrome plated brass. All polypropylene internal construction with 6" spread gooseneck. Furnished with 1/2" I.P.S. shank with locknut and washer. Inlet is 3/8" O.D. compression.




**W-0409-00 Distilled Deionized and Ultra-Pure Water Faucet (panel mounted)** made of chrome plated brass. All polypropylene internal construction with 8" spread gooseneck. Furnished with 1/2" I.P.S. shank with locknut and washer. Inlet is 3/8" O.D. compression.

# Laboratory Service Fittings


## Vandal-Resistant Fittings




**W-0410-VR Vandal Resistant Deck Mounted Gas, Air or Vacuum Fitting** with single vandal-resistant laboratory ball valve. Complete with a 3/8" I.P.S. tank nipple, locknut and washer. **Specify service desired.**




**W-0411-VR Vandal Resistant Deck Mounted Double Gas, Air or Vacuum Fitting** with two vandal resistant laboratory ball valves with integral straight serrated hose ends at 90°. Complete with a 3/8" I.P.S. tank nipple, locknut and washer. **Specify service desired.**




**W-0412-VR Vandal Resistant Deck Mounted Double Gas, Air or Vacuum Fitting** Same as W-0411-VR except at 180°. Complete with a 3/8" I.P.S. tank nipple, locknut and washer. **Specify service desired.**




**W-0413-VR Vandal Resistant Deck Mounted Hot and Cold Mixing Faucet** with 8" rigid vacuum breaker gooseneck. Valves are disc type with renewable internal assemblies. Faucet body is complete with an integral threaded mounting shank and a flanged locknut. 1/2" I.P.S. male thread inlets with coupling nuts for slip joint connections.




**W-0414-VR Vandal Resistant Panel Mounted Hot and Cold Mixing Faucet** with 8" rigid vacuum breaker gooseneck and mounting shanks.




**W-0415-VR Vandal Resistant Deck Mounted Water Faucet** with 8" rigid vacuum breaker gooseneck. Furnished with 3/8" IPS mounting shank with locknut and washer.



**W-0417-VR Vandal Resistant Deck Mounted Hand-Held Drench Hose** with single spray type outlet heads, squeeze-lever operated valve and flexible stainless steel hose. 3/8" I.P.S. male inlet.




**W-0416-VB Vandal Resistant Deck Mounted Combination Vacuum Breaker Water/Gas Fitting** Made of chrome plated base. Rigid gooseneck with disc type valve and in-line vacuum breaker. Includes two-way turret with two laboratory ball valves. Complete with two 3/8" I.P.S. tank nipples, locknuts and washers. **Specify service desired.**




**W-0416-VR Vandal Resistant Deck Mounted Combination Water/Gas Fitting** Same as W-416-VB except furnished without in-line vacuum breakers. **Specify service desired.**

# Laboratory Service Fittings


## Emergency Eye Wash/Shower Equipment




**W-0927-00 Panel Mounted Eye Wash Fixture** Consists of stainless steel bowl, two Gentle Spray outlet heads, a "PUSH TO OPERATE" ball valve and mounting bracket. 1/2" I.P.S. female thread inlet and 1 1/4" I.P.S. female thread outlet.




**W-0933-00 Eye/Face Wash (wall mounted)** the same as No. W-0927-00, except with four Gentle Spray outlet heads for complete coverage of eyes and face.




**W-0930-00 BARRIER-FREE (ADA) Eye/Face Wash (wall mounted)** with four Gentle Spray outlet heads, stainless steel skirt and bowl and stay-open ball valve. 1/2" I.P.S. female thread inlet and 1 1/4" I.P.S. female thread outlet.




**W-0931-00 Eye Wash (deck mounted)** with two Gentle Spray outlet heads, stay-open ball valve and 90° swivel feature (right hand mounting).




**W-0932-00 Eye Wash (deck mounted)** the same as No. W-0931-00 but (left hand mounting).




**W-0923-00 Deck Mounted Hose Spray Fitting** Consists of 8 ft. reinforced PVC hose with a squeeze-lever operated valve and spray-type outlet head with protective rubber collar. Flange, shank and locknut for attaching to table top are included. 3/8" I.P.S. male thread swivel inlet.




**W-0934-00 Deck Mounted Dual Purpose Eye Wash/Drench Hose** with two Gentle Spray outlet heads angled at 45° and squeeze valve with locking clip. Flange, shank and locknut for attaching to table top are included. 3/8" I.P.S. male thread swivel inlet.




**W-0935-00 Wall Mounted Dual Purpose Eye Wash/Drench Hose** the same as No. W-0934-00 except wall mounted.



**W-0374-01 Emergency Shower Head Assembly (horizontally mounted).** Includes 10" diameter orange ABS plastic shower head. Operated by a rigid pull-rod and easy to grip handle. Complete with an elbow, connecting pipe and ball valve. 1" I.P.S. female thread inlet.



**W-0928-00 Emergency Shower Head Assembly (vertically mounted).** The same as No. W-0374-01 except vertically mounted.



**W-0926-00 Drench Shower and Eye Wash Combination** A necessity for chemical stockrooms and other areas handling a quantity of corrosive chemicals. Free-standing 10" dia. drench shower made of orange ABS plastic. Has a 1" ball valve operated by a rigid pull-rod and easy to grip handle. Eye wash fixture with stainless steel bowl, two Gentle Spray heads and a "PUSH TO OPERATE" ball valve. Eye wash and shower operate independently. Overall height is 7'6". 1 1/4" I.P.S. female thread inlet and outlet.

**W-0936-00**

**W-0936-00 BARRIER-FREE (ADA) Drench Shower and Eye Wash Combination** Similar to No. W-0926-00 except this unit is barrier-free and eye wash furnished with four Gentle Spray heads.



# Cupsinks & Drain Accessories

## Cupsinks



**0490-BE** Black Epoxy  
**0490-GE** Grey Epoxy  
**0490-PE** Putty Epoxy  
**0490-SE** Slate Epoxy  
 Complete with removable strainer. 6"x 3" inside dimension. 1 1/2" I.P.S. male straight thread outlet. (Overall height is 7 1/2")



**0491-BP** Black Poly  
**0491-GP** Grey Poly  
 Complete with removable strainer. 6"x 3" inside dimension. 1 1/2" I.P.S. male straight thread outlet. (Overall height is 8 1/2")



**0492-BP** Black Poly  
**0492-GP** Grey Poly  
 Complete with removable strainer. 9"x 3" inside dimension. 1 1/2" I.P.S. male straight thread outlet. (Overall height is 7")



**0493-BE** Black Epoxy  
**0493-GE** Grey Epoxy  
**0493-PE** Putty Epoxy  
**0493-SE** Slate Epoxy  
 Complete with removable type 316 stainless steel wire screen. 13 3/8" x 4 1/8" x 5 1/4" I.D. 1 1/2" I.P.S. male straight thread outlet. (Overall height is 8")

## Stainless Steel Cupsinks




**0950-0F** Stainless Steel  
 Made of type 304 stainless steel and has integral cross bars. 5" I.D. 1 1/2" I.P.S. male straight thread outlet.




**0975-0F** Stainless Steel  
 Made of type 316 stainless steel and has integral cross bars. 6"x 3" I.D. 1 1/2" I.P.S. male straight thread outlet.

## Drain Accessories



**0464-00** Polyethylene Trap and Sink Outlet  
 Adjustable riser (tailpiece) can be cut to the desired length. Sink outlet is complete with plug strainer, gasket and locknut. Has a 1 1/2" I.P.S. male straight thread outlet. Running type trap can be used for either "P" or "S" type drain connections. 1 1/2" I.P.S. male straight thread outlet.



**0469-00** Polyethylene Trap and Sink Outlet  
 Adjustable riser (tailpiece) with 4 1/2" adjustment that does not require cutting. The sink outlet is complete with plug, strainer, gasket and locknut. Has a 1 1/2" I.P.S. male straight thread outlet. Trap is the bottle "P" type, with removable bottom for cleanout. 1 1/2" I.P.S. male straight thread outlet.

## Sink Outlet




**0482-BP** Black Poly  
**0482-GP** Grey Poly  
**0482-PP** Putty Poly  
**0482-SP** Slate Poly  
 Complete with locknut. 1 1/2" I.P.S. male straight thread outlet. For use with epoxy resin tub sinks and troughs.

## Sink Overflow



**0494BP-04** 4" High  
**0494BP-06** 6" High  
**0494BP-08** 8" High  
**0494BP-10** 10" High  
 Sink Overflows are designed to be used with 0482-BP sink outlets. 1 1/2" diameter.

## St. Steel Sink Drain



**0995-00** Stainless Steel  
 Crumb Cup Strainer  
 Complete with gaskets, locknut and basket strainer. 1 1/2" male straight thread outlet.

# Kemresin Tub Sinks & Drain Troughs

Kemresin® Tub Sinks and Drain Troughs are made from a combination of modified epoxy resins, upgraded for maximum properties with carefully selected curing resins and inert fillers to produce a thermosetting material. They are formed in metal molds to provide one-piece units with coved corners and bottoms pitched to the drain outlet. After being removed

from the molds, they are cured at high temperatures for several hours to complete the chemical reaction, and ensure maximum chemical, impact and thermal shock resistance. Drain outlets are sized for 1½" sink outlets unless otherwise specified.

Kemresin Drain Troughs are molded with no longitudinal joints. A slope of 1/8" per foot assures complete drainage, and a 7/8" cove radius

in the bottom corners is provided for ease of cleaning. Five trough sections, or any part of them, are available to be cemented together to provide Open End Troughs up to 21' long, Closed End Troughs with an end drain outlet up to 21'-8" long and Closed End Troughs with a center drain outlet up to 42'-8" long. Section TR-4 is provided with a center drain recess for a standard 1½" sink outlet.

## Epoxy Resin Drain Troughs

## Epoxy Resin Tub Sinks


**Closed End Trough**



**Open End Trough**











**Double Tub Sink**



**Center Bench End Sink with Trough Cutout**



**Single Tub Sink**



# Kemresin Tub Sinks & Drain Troughs

### Kemresin Colors

Kemresin® Tub Sinks and Drain Troughs are available in three standard colors: Black, Grey, and Slate.

To specify color, add color code to the end of the part number.  
**Example: 1005-00-BK**

### Color Codes

**BK** = Black      **PT** = Putty  
**GR** = Grey      **SL** = Slate

### Kemresin Tub Sinks

Catalog Number	Inside Dimensions (at Top of Sink)			Drain Location		Overall Sink Height	Catalog Number	Inside Dimensions (at Top of Sink)			Drain Location		Overall Sink Height
	Length	Width	Height	A	B			Length	Width	Height	A	B	
1000-00-__	14"	10"	6"	4"	5"	6 <sup>7</sup> / <sub>8</sub> "	SM-13-00-__	14"	10"	8"	4"	5"	9"
1001-00-__	16"	8"	7"	4"	4"	8"	SM-14-00-__	23 <sup>1</sup> / <sub>4</sub> "	14 <sup>1</sup> / <sub>4</sub> "	12 <sup>3</sup> / <sub>4</sub> "	11 <sup>5</sup> / <sub>8</sub> "	9 <sup>1</sup> / <sub>4</sub> "	13 <sup>5</sup> / <sub>8</sub> "
1003-00-__	18"	15"	11 <sup>13</sup> / <sub>16</sub> "	4"	4"	12 <sup>7</sup> / <sub>8</sub> "	SM-15-00-__	19 <sup>3</sup> / <sub>4</sub> "	18"	24"	9 <sup>7</sup> / <sub>8</sub> "	9"	25"
(1) 1003-DT-__	37"	15"	11"	4 <sup>1</sup> / <sub>2</sub> "	4 <sup>1</sup> / <sub>2</sub> "	12"	SM-16-00-__	40"	18"	18 <sup>1</sup> / <sub>4</sub> "	20"	9"	19 <sup>1</sup> / <sub>4</sub> "
1005-00-__	25"	15"	10"	4"	4"	11 <sup>1</sup> / <sub>4</sub> "	SM-17-00-__	29 <sup>3</sup> / <sub>4</sub> "	23 <sup>1</sup> / <sub>2</sub> "	6"	14 <sup>7</sup> / <sub>8</sub> "	11 <sup>3</sup> / <sub>4</sub> "	7"
1006-00-__	18"	15"	8"	4"	4"	9"	SM-18-00-__	18"	15"	15 <sup>1</sup> / <sub>4</sub> "	4 <sup>3</sup> / <sub>4</sub> "	4 <sup>3</sup> / <sub>4</sub> "	16 <sup>1</sup> / <sub>2</sub> "
(2) 1006-DT-__	36 <sup>3</sup> / <sub>4</sub> "	15"	8 <sup>1</sup> / <sub>4</sub> "	4 <sup>1</sup> / <sub>2</sub> "	4 <sup>1</sup> / <sub>2</sub> "	9"	SM-19-00-__	35 <sup>1</sup> / <sub>2</sub> "	19 <sup>1</sup> / <sub>2</sub> "	15 <sup>3</sup> / <sub>8</sub> "	3 <sup>1</sup> / <sub>2</sub> "	4 <sup>1</sup> / <sub>2</sub> "	16 <sup>1</sup> / <sub>2</sub> "
1007-00-__	24"	16"	8"	4"	8"	9 <sup>1</sup> / <sub>8</sub> "	SM-20-00-__	21"	18"	15 <sup>1</sup> / <sub>4</sub> "	4 <sup>1</sup> / <sub>2</sub> "	4 <sup>1</sup> / <sub>2</sub> "	16 <sup>1</sup> / <sub>2</sub> "
1008-00-__	16"	12"	8"	4"	6"	9 <sup>1</sup> / <sub>8</sub> "	SM-21-00-__	60 <sup>1</sup> / <sub>2</sub> "	22"	7 <sup>1</sup> / <sub>4</sub> "	6"	6"	8 <sup>1</sup> / <sub>8</sub> "
1011-00-__	18"	6 <sup>1</sup> / <sub>2</sub> "	6"	9"	3 <sup>1</sup> / <sub>4</sub> "	7"	SM-22-00-__	24"	16"	15 <sup>3</sup> / <sub>8</sub> "	12"	8"	16 <sup>1</sup> / <sub>8</sub> "
(3) 1014-00-__	25"	15"	10"	4"	4"	11 <sup>1</sup> / <sub>4</sub> "	SM-23-00-__	28"	18"	15 <sup>1</sup> / <sub>2</sub> "	5 <sup>1</sup> / <sub>2</sub> "	4"	16 <sup>1</sup> / <sub>2</sub> "
(4) 1015-00-__	23 <sup>3</sup> / <sub>8</sub> "	13 <sup>3</sup> / <sub>4</sub> "	13 <sup>1</sup> / <sub>4</sub> "	11 <sup>1</sup> / <sub>16</sub> "	9 <sup>1</sup> / <sub>2</sub> "	17 <sup>3</sup> / <sub>4</sub> "	SM-24-00-__	48"	8"	8"	24"	4"	9"
(5) 1016-00-__	23 <sup>3</sup> / <sub>8</sub> "	13 <sup>3</sup> / <sub>4</sub> "	13 <sup>1</sup> / <sub>4</sub> "	11 <sup>1</sup> / <sub>16</sub> "	9 <sup>1</sup> / <sub>2</sub> "	17 <sup>3</sup> / <sub>4</sub> "	SM-25-00-__	21"	18"	15 <sup>1</sup> / <sub>8</sub> "	10 <sup>1</sup> / <sub>2</sub> "	9"	16 <sup>1</sup> / <sub>8</sub> "
1021-00-__	16"	16"	7 <sup>1</sup> / <sub>2</sub> "	4"	4"	8 <sup>3</sup> / <sub>4</sub> "	SM-26-00-__	21"	11 <sup>3</sup> / <sub>8</sub> "	15 <sup>1</sup> / <sub>8</sub> "	10 <sup>1</sup> / <sub>2</sub> "	5 <sup>1</sup> / <sub>16</sub> "	16 <sup>1</sup> / <sub>8</sub> "
1022-00-__	21"	18"	10"	4 <sup>1</sup> / <sub>2</sub> "	9"	11"	SM-27-00-__	36"	20"	16 <sup>1</sup> / <sub>4</sub> "	18"	10"	17 <sup>1</sup> / <sub>4</sub> "
1023-00-__	21 <sup>1</sup> / <sub>2</sub> "	15 <sup>1</sup> / <sub>2</sub> "	7"	4"	4"	8 <sup>3</sup> / <sub>8</sub> "	SM-28-00-__	30"	18"	4 <sup>1</sup> / <sub>8</sub> "	15"	9"	5 <sup>1</sup> / <sub>4</sub> "
1024-00-__	21 <sup>1</sup> / <sub>2</sub> "	15 <sup>1</sup> / <sub>2</sub> "	11"	4"	4"	12 <sup>1</sup> / <sub>8</sub> "	SM-29-00-__	36"	24"	12 <sup>1</sup> / <sub>4</sub> "	18"	12"	13 <sup>1</sup> / <sub>4</sub> "
(6) 1024-DT-__	44"	15 <sup>1</sup> / <sub>2</sub> "	11"	4"	4"	12 <sup>1</sup> / <sub>2</sub> "	SM-30-00-__	40"	30"	8 <sup>1</sup> / <sub>8</sub> "	20"	15"	9"
1025-00-__	28"	20"	9 <sup>1</sup> / <sub>8</sub> "	14"	6"	10 <sup>1</sup> / <sub>8</sub> "	SM-31-00-__	24"	15 <sup>3</sup> / <sub>4</sub> "	12 <sup>1</sup> / <sub>8</sub> "	4 <sup>3</sup> / <sub>4</sub> "	8"	13"
1026-00-__	35 <sup>1</sup> / <sub>2</sub> "	19 <sup>1</sup> / <sub>2</sub> "	10"	17 <sup>3</sup> / <sub>4</sub> "	9 <sup>3</sup> / <sub>4</sub> "	11 <sup>1</sup> / <sub>8</sub> "	SM-32-00-__	59"	21"	12 <sup>5</sup> / <sub>8</sub> "	29 <sup>1</sup> / <sub>2</sub> "	10 <sup>1</sup> / <sub>2</sub> "	13 <sup>3</sup> / <sub>8</sub> "
1027-00-__	42"	20"	10"	21"	10"	11 <sup>1</sup> / <sub>4</sub> "	SM-33-00-__	24"	18"	11"	4"	9"	11 <sup>7</sup> / <sub>8</sub> "
1030-00-__	16"	12"	8"	8"	6"	9"	(10) SM-34-00-__	15"	14"	10 <sup>1</sup> / <sub>2</sub> "	7 <sup>1</sup> / <sub>2</sub> "	7"	11 <sup>1</sup> / <sub>4</sub> "
1031-00-__	24"	16"	8"	12"	8"	9"	SM-35-00-__	51"	17"	12 <sup>1</sup> / <sub>2</sub> "	25 <sup>1</sup> / <sub>2</sub> "	8 <sup>1</sup> / <sub>2</sub> "	13 <sup>1</sup> / <sub>4</sub> "
1032-00-__	29 <sup>3</sup> / <sub>4</sub> "	14 <sup>3</sup> / <sub>4</sub> "	6"	14 <sup>7</sup> / <sub>8</sub> "	7 <sup>3</sup> / <sub>8</sub> "	7"	SM-36-00-__	14"	10"	6"	7"	5"	6 <sup>7</sup> / <sub>8</sub> "
1034-00-__	21 <sup>1</sup> / <sub>2</sub> "	19 <sup>1</sup> / <sub>4</sub> "	8"	10 <sup>3</sup> / <sub>4</sub> "	9 <sup>5</sup> / <sub>8</sub> "	9 <sup>1</sup> / <sub>8</sub> "	(11) SM-37-00-__	24"	20"	10 <sup>5</sup> / <sub>8</sub> "	6 <sup>3</sup> / <sub>4</sub> "	10"	11 <sup>1</sup> / <sub>4</sub> "
1035-00-__	45 <sup>3</sup> / <sub>4</sub> "	17 <sup>3</sup> / <sub>8</sub> "	12 <sup>1</sup> / <sub>2</sub> "	22 <sup>7</sup> / <sub>8</sub> "	8 <sup>3</sup> / <sub>4</sub> "	13 <sup>3</sup> / <sub>8</sub> "	SM-39-00-__	35"	6 <sup>5</sup> / <sub>8</sub> "	4 <sup>3</sup> / <sub>4</sub> "	17 <sup>1</sup> / <sub>2</sub> "	3 <sup>5</sup> / <sub>16</sub> "	5 <sup>1</sup> / <sub>2</sub> "
1036-00-__	48 <sup>5</sup> / <sub>8</sub> "	24 <sup>1</sup> / <sub>2</sub> "	5 <sup>3</sup> / <sub>8</sub> "	3 <sup>3</sup> / <sub>4</sub> "	12 <sup>1</sup> / <sub>4</sub> "	6"	SM-40-00-__	33 <sup>3</sup> / <sub>8</sub> "	29 <sup>1</sup> / <sub>8</sub> "	15 <sup>3</sup> / <sub>4</sub> "	6"	14 <sup>9</sup> / <sub>16</sub> "	16 <sup>1</sup> / <sub>2</sub> "
(7) 1037-00-__	23 <sup>1</sup> / <sub>4</sub> "	14 <sup>1</sup> / <sub>4</sub> "	12 <sup>3</sup> / <sub>4</sub> "	11 <sup>5</sup> / <sub>8</sub> "	9 <sup>1</sup> / <sub>4</sub> "	22 <sup>3</sup> / <sub>8</sub> "	SM-41-00-__	18"	15"	8 <sup>1</sup> / <sub>4</sub> "	9"	7 <sup>1</sup> / <sub>2</sub> "	9"
(8) 1038-00-__	12" diameter		7 <sup>5</sup> / <sub>8</sub> "	center		8 <sup>3</sup> / <sub>8</sub> "	SM-42-00-__	25"	15"	10 <sup>1</sup> / <sub>2</sub> "	12 <sup>1</sup> / <sub>2</sub> "	7 <sup>1</sup> / <sub>2</sub> "	11 <sup>1</sup> / <sub>4</sub> "
(9) 1039-00-__	18" diameter		8 <sup>1</sup> / <sub>8</sub> "	center		8 <sup>7</sup> / <sub>8</sub> "	SM-43-00-__	21"	18"	10 <sup>3</sup> / <sub>8</sub> "	10 <sup>1</sup> / <sub>2</sub> "	9"	11 <sup>1</sup> / <sub>8</sub> "
SM-01-00-__	12"	6"	6"	6"	3"	7"	SM-44-00-__	36"	16"	11"	12"	4"	12"
SM-02-00-__	19 <sup>1</sup> / <sub>2</sub> "	16 <sup>1</sup> / <sub>2</sub> "	8"	9 <sup>3</sup> / <sub>4</sub> "	8 <sup>1</sup> / <sub>4</sub> "	9"	SM-45-00-__	27 <sup>3</sup> / <sub>4</sub> "	14 <sup>3</sup> / <sub>4</sub> "	8 <sup>1</sup> / <sub>2</sub> "	10 <sup>1</sup> / <sub>2</sub> "	7 <sup>1</sup> / <sub>2</sub> "	9 <sup>3</sup> / <sub>8</sub> "
SM-03-00-__	24"	21"	8"	6"	4 <sup>1</sup> / <sub>2</sub> "	9 <sup>1</sup> / <sub>8</sub> "	SM-47-00-__	43"	15 <sup>1</sup> / <sub>2</sub> "	11 <sup>1</sup> / <sub>8</sub> "	21 <sup>1</sup> / <sub>2</sub> "	7 <sup>3</sup> / <sub>4</sub> "	12"
SM-04-00-__	25 <sup>1</sup> / <sub>2</sub> "	16 <sup>1</sup> / <sub>2</sub> "	19"	4 <sup>1</sup> / <sub>2</sub> "	4 <sup>1</sup> / <sub>2</sub> "	20 <sup>1</sup> / <sub>8</sub> "	SM-48-00-__	56"	24"	8 <sup>1</sup> / <sub>4</sub> "	28"	12"	9"
SM-05-00-__	28"	18"	10"	4"	4"	11 <sup>1</sup> / <sub>8</sub> "	SM-49-00-__	36"	16"	8 <sup>3</sup> / <sub>8</sub> "	6"	8"	9"
SM-06-00-__	28"	18"	10"	24"	4"	11 <sup>1</sup> / <sub>8</sub> "	SM-50-00-__	24"	24"	12 <sup>1</sup> / <sub>2</sub> "	12"	12"	13 <sup>1</sup> / <sub>4</sub> "
SM-07-00-__	30"	16"	10"	4"	4"	11 <sup>1</sup> / <sub>8</sub> "	SM-51-00-__	30"	15 <sup>7</sup> / <sub>8</sub> "	16 <sup>3</sup> / <sub>8</sub> "	6"	7 <sup>15</sup> / <sub>16</sub> "	17 <sup>1</sup> / <sub>8</sub> "
SM-08-00-__	30"	16"	10"	26"	4"	11 <sup>1</sup> / <sub>8</sub> "	SM-52-00-__	45 <sup>1</sup> / <sub>8</sub> "	11 <sup>1</sup> / <sub>2</sub> "	8 <sup>1</sup> / <sub>8</sub> "	22 <sup>9</sup> / <sub>16</sub> "	5 <sup>3</sup> / <sub>4</sub> "	9"
SM-09-00-__	30"	18"	18"	4 <sup>1</sup> / <sub>2</sub> "	4 <sup>3</sup> / <sub>4</sub> "	18 <sup>3</sup> / <sub>4</sub> "	SM-53-00-__	31 <sup>1</sup> / <sub>2</sub> "	16"	16 <sup>1</sup> / <sub>2</sub> "	7 <sup>1</sup> / <sub>8</sub> "	8"	17 <sup>1</sup> / <sub>4</sub> "
SM-10-00-__	36"	21"	4 <sup>3</sup> / <sub>4</sub> "	4"	4"	5 <sup>3</sup> / <sub>4</sub> "	SM-55-00-__	16"	16"	12 <sup>3</sup> / <sub>8</sub> "	4 <sup>1</sup> / <sub>4</sub> "	4 <sup>1</sup> / <sub>4</sub> "	13"
SM-11-00-__	48"	13"	12 <sup>1</sup> / <sub>4</sub> "	24"	6 <sup>1</sup> / <sub>2</sub> "	13 <sup>1</sup> / <sub>4</sub> "	SM-56-00-__	42 <sup>1</sup> / <sub>4</sub> "	18 <sup>1</sup> / <sub>8</sub> "	12 <sup>5</sup> / <sub>8</sub> "	21 <sup>1</sup> / <sub>8</sub> "	9 <sup>1</sup> / <sub>16</sub> "	13 <sup>3</sup> / <sub>8</sub> "
SM-12-00-__	8"	8"	8"	4"	4"	8 <sup>7</sup> / <sub>8</sub> "							

- (1) Double tub sink. Distance between centerlines of drain outlets is 19<sup>3</sup>/<sub>8</sub>".
- (2) Double tub sink. Distance between centerlines of drain outlets is 19".
- (3) Sink has 7<sup>1</sup>/<sub>8</sub>" long by 5<sup>5</sup>/<sub>8</sub>" high trough cutout.
- (4) Center table end sink has 7<sup>1</sup>/<sub>8</sub>" long by 5<sup>5</sup>/<sub>8</sub>" high trough cutout and 24<sup>1</sup>/<sub>2</sub>" long by 3<sup>3</sup>/<sub>4</sub>" high back splash. Table side of sink is not tapered.
- (5) Center table end sink has 24<sup>1</sup>/<sub>2</sub>" long by 3<sup>3</sup>/<sub>4</sub>" high back splash without trough cutout. Table side of sink is not tapered.
- (6) Double tub sink. Distance between centerlines of drain outlets is 22<sup>1</sup>/<sub>2</sub>".
- (7) Assembly of tub sink, back splash, ledge and end fillers. Used for T69W4525-00 Scrub Sink Assembly.
- (8) Cylindrical sink – 12" inside diameter at top of sink with center drain outlet.
- (9) Cylindrical sink – 18" inside diameter at top of sink with center drain outlet.
- (10) Sides have 1<sup>1</sup>/<sub>2</sub>" taper per foot in height.
- (11) Exterior of one 21" long side is not tapered.

**General Notes:**

- Sinks with center drain outlets, except 1011-00, can be furnished with a drain recess sized for a crumb cup strainer or garbage disposal when specified.
- Any sink may be furnished with trough cutout. Specify location and size. (Trough location must be at least 1<sup>1</sup>/<sub>2</sub>" from inside corner and 2" from bottom.)
- Most sinks can be joined together to form multiple tub sinks provided adjoining sides of sinks are the same size. Sinks longer than 30", wider than 24" or deeper than 12" cannot be factory joined, but they can be factory prepared for field assembly. Assemblies of three or more tub sinks should be no longer than 60" overall. For longer assemblies, one factory joint can be made with a second joint prepared for field assembly.




# Kemresin Drop-in Sinks

Kemresin Drop-in Sinks are formed in metal molds to provide a lipped one-piece tub with covered corners and bottoms pitched to the drain outlet. They are made from a combination of modified epoxy resins, upgraded for maximum properties with carefully selected curing resins and inert fillers to produce a thermosetting


material. After being removed from the molds, they are cured at high temperatures for several hours to complete the chemical reaction, and ensure maximum chemical, impact and thermal shock resistance. Drain outlets are sized for 1½" sink outlets unless otherwise specified.

They are designed to be mounted in Kemresin tops in a specially rabbetted cutout so that no support from below is required.


## Drop-in Kemresin Sinks



Part Number	I.D.		Depth	O.D.		Drain Location	
	A	B	C	D	E	G	H
<b>1000-DI</b> ---	13 <sup>7</sup> / <sub>8</sub> "	9 <sup>7</sup> / <sub>8</sub> "	6"	16 <sup>1</sup> / <sub>8</sub> "	12 <sup>1</sup> / <sub>8</sub> "	7"	5"
<b>1001-DI</b> ---	8"	16"	7"	10 <sup>1</sup> / <sub>4</sub> "	18 <sup>1</sup> / <sub>4</sub> "	4"	4"
<b>1003-DI</b> ---	17 <sup>15</sup> / <sub>16</sub> "	14 <sup>15</sup> / <sub>16</sub> "	11"	20 <sup>3</sup> / <sub>16</sub> "	17 <sup>3</sup> / <sub>16</sub> "	3 <sup>1</sup> / <sub>2</sub> "	3 <sup>1</sup> / <sub>2</sub> "
<b>1005-DI</b> ---	25"	15"	10"	27 <sup>1</sup> / <sub>4</sub> "	17 <sup>1</sup> / <sub>4</sub> "	4"	4"
<b>1006-DI</b> ---	17 <sup>7</sup> / <sub>8</sub> "	14 <sup>7</sup> / <sub>8</sub> "	8"	20 <sup>1</sup> / <sub>8</sub> "	17 <sup>1</sup> / <sub>8</sub> "	3 <sup>1</sup> / <sub>2</sub> "	3 <sup>1</sup> / <sub>2</sub> "
<b>1008-DI</b> ---	16"	12"	8"	18 <sup>1</sup> / <sub>4</sub> "	14 <sup>1</sup> / <sub>4</sub> "	4"	6"
<b>1011-DI</b> ---	18"	6 <sup>1</sup> / <sub>2</sub> "	6"	20 <sup>1</sup> / <sub>4</sub> "	8 <sup>3</sup> / <sub>4</sub> "	9"	3 <sup>1</sup> / <sub>4</sub> "
<b>1021-DI</b> ---	16"	16"	7 <sup>1</sup> / <sub>2</sub> "	18 <sup>1</sup> / <sub>4</sub> "	18 <sup>1</sup> / <sub>4</sub> "	4"	4"
<b>1022-DI</b> ---	21"	18"	10"	23 <sup>1</sup> / <sub>4</sub> "	20 <sup>1</sup> / <sub>4</sub> "	4"	9"
(1) <b>1038-DI</b> ---	12" diameter		7 <sup>5</sup> / <sub>8</sub> "	14 <sup>1</sup> / <sub>4</sub> " diameter		center	
(2) <b>1039-DI</b> ---	18" diameter		8 <sup>1</sup> / <sub>8</sub> "	20 <sup>1</sup> / <sub>4</sub> " diameter		center	
<b>1040-DI</b> ---	15 <sup>3</sup> / <sub>8</sub> "	11 <sup>3</sup> / <sub>8</sub> "	8"	17.6"	13.6"	4"	4"
<b>1041-DI</b> ---	23 <sup>3</sup> / <sub>8</sub> "	15 <sup>3</sup> / <sub>8</sub> "	8.3"	25.6"	17.6"	4"	4"
<b>1042-DI</b> ---	28"	15"	11.8"	30 <sup>1</sup> / <sub>4</sub> "	17 <sup>1</sup> / <sub>4</sub> "	5"	5"
<b>1043-DI</b> ---	30"	16"	17.8"	32 <sup>1</sup> / <sub>4</sub> "	18 <sup>1</sup> / <sub>4</sub> "	5"	8"
<b>1044-DI</b> ---	12"	8"	6"	14 <sup>1</sup> / <sub>4</sub> "	10 <sup>1</sup> / <sub>4</sub> "	6"	4"
<b>1045-DI</b> ---	21"	17"	10"	23 <sup>1</sup> / <sub>4</sub> "	19 <sup>1</sup> / <sub>4</sub> "	4 <sup>1</sup> / <sub>2</sub> "	8 <sup>1</sup> / <sub>2</sub> "
<b>1046-DI</b> ---	21 <sup>1</sup> / <sub>2</sub> "	15 <sup>1</sup> / <sub>2</sub> "	11"	23 <sup>3</sup> / <sub>4</sub> "	17 <sup>3</sup> / <sub>4</sub> "	4"	4"
<b>1047-DI</b> ---	24"	18"	11"	26 <sup>1</sup> / <sub>4</sub> "	20 <sup>1</sup> / <sub>4</sub> "	4"	4"
(3) <b>1048-DI</b> ---	17 <sup>1</sup> / <sub>2</sub> "	12 <sup>3</sup> / <sub>8</sub> "	6 <sup>1</sup> / <sub>2</sub> "	19 <sup>1</sup> / <sub>2</sub> "	14 <sup>5</sup> / <sub>8</sub> "	8 <sup>3</sup> / <sub>4</sub> "	6 <sup>3</sup> / <sub>16</sub> "
(4) <b>1049-DI</b> ---	18"	15"	11"	20 <sup>1</sup> / <sub>4</sub> "	17 <sup>1</sup> / <sub>4</sub> "	4 <sup>1</sup> / <sub>2</sub> "	5"
<b>1070-DI</b> ---	24"	16"	15 <sup>1</sup> / <sub>2</sub> "	26 <sup>1</sup> / <sub>4</sub> "	18 <sup>1</sup> / <sub>4</sub> "	12"	8"



PLAN VIEW  
**1048-DI**



SECTION VIEW  
**1049-DI**

### ADA Compliant

<b>1000-ADA-DI</b> ---	13 <sup>7</sup> / <sub>8</sub> "	9 <sup>7</sup> / <sub>8</sub> "	5"	16 <sup>1</sup> / <sub>8</sub> "	12 <sup>1</sup> / <sub>8</sub> "	4"	4"
<b>1005-ADA-DI</b> ---	25"	15"	4.8"	27 <sup>1</sup> / <sub>4</sub> "	17 <sup>1</sup> / <sub>4</sub> "	4"	4"
<b>1006-ADA-DI</b> ---	17 <sup>7</sup> / <sub>8</sub> "	14 <sup>7</sup> / <sub>8</sub> "	5"	20 <sup>1</sup> / <sub>8</sub> "	17 <sup>1</sup> / <sub>8</sub> "	3 <sup>1</sup> / <sub>2</sub> "	3 <sup>1</sup> / <sub>2</sub> "

(1) Cylindrical sink, 12" inside diameter  
(2) Cylindrical sink, 18" inside diameter

(3) Triangular sink used in X92 and X94 assemblies  
(4) Stepped sink, 4<sup>3</sup>/<sub>4</sub>" deep at front

### Kemresin Colors

Kemresin® Drop-in Sinks are available in four standard colors: Black, Grey, Putty, and Slate.

To specify color, add color code to the end of the part number.

**Example: 1005-DI-BK**

### Color Codes

**BK** = Black      **PT** = Putty  
**GR** = Grey      **SL** = Slate

# Stainless Steel Drop-in Sinks


Seamless die-drawn from nickel bearing type 304, 18-8, 18 gauge stainless steel. Allows easy cleaning and stain resistance. Interior and top surfaces polished to a gleaming finish with a highlighted bowl rim. Corners and bottom turns are radiused  $1\frac{3}{4}$ ". Sinks are thoroughly sound deadened with environmentally

safe compounds, to ensure proper sound insulation and to reduce condensation.

All sinks are self-rimming with stainless steel mounting channels for easy, top mount, snap-in clip installation. Sinks are compatible with any type work-surface material from  $\frac{3}{4}$ " to  $2\frac{1}{2}$ " thick.


Each sink is center punched to receive a standard cup strainer. ADA compliant sinks are punched in the center,  $4\frac{1}{2}$ " from the rear of the bowl for extra clearance.

## Stainless Steel Sinks



Part Number	O.D.		I.D.		Sink Depth
	Front to Back	Side to Side	Front to Back	Side to Side	
<b>F-5353-08</b>	15"	15"	12"	12"	7 $\frac{1}{2}$ "
<b>F-5354-08</b>	17"	19"	14"	16"	7 $\frac{1}{2}$ "
<b>F-5355-08</b>	19"	22"	16"	19"	7 $\frac{1}{2}$ "
<b>F-5356-08</b>	19"	25"	16"	22"	7 $\frac{1}{2}$ "
<b>F-5357-08</b>	19"	31"	16"	28"	7 $\frac{1}{2}$ "
<b>F-5353-09</b>	15"	15"	12"	12"	9"
<b>F-5354-11</b>	17"	19"	14"	16"	10 $\frac{1}{2}$ "
<b>F-5355-11</b>	19"	22"	16"	19"	10 $\frac{1}{2}$ "
<b>F-5356-11</b>	19"	25"	16"	22"	10 $\frac{1}{2}$ "
<b>F-5357-11</b>	19"	31"	16"	28"	10 $\frac{1}{2}$ "
<b>F-5354-12</b>	17"	19"	14"	16"	12"
<b>F-5355-12</b>	19"	22"	16"	19"	12"
<b>F-5356-12</b>	19"	25"	16"	22"	12"
<b>F-5357-12</b>	19"	31"	16"	28"	12"


## Stainless Steel Sinks with Fixture Ledge



Part Number	O.D.		I.D.		Sink Depth	Fixture Centers
	Front to Back	Side to Side	Front to Back	Side to Side		
<b>F-5358-08</b>	18"	15"	12"	12"	7 $\frac{1}{2}$ "	4"
<b>F-5359-08</b>	20"	19"	14"	16"	7 $\frac{1}{2}$ "	6"
<b>F-5360-08</b>	22"	22"	16"	19"	7 $\frac{1}{2}$ "	6"
<b>F-5361-08</b>	22"	25"	16"	22"	8"	8"
<b>F-5362-08</b>	22"	31"	16"	28"	7 $\frac{1}{2}$ "	8"
<b>F-5358-09</b>	18"	15"	12"	12"	9"	4"
<b>F-5359-11</b>	20"	19"	14"	16"	10 $\frac{1}{2}$ "	6"
<b>F-5360-11</b>	22"	22"	16"	19"	10 $\frac{1}{2}$ "	6"
<b>F-5361-11</b>	22"	25"	16"	22"	10 $\frac{1}{2}$ "	8"
<b>F-5362-11</b>	22"	31"	16"	28"	10 $\frac{1}{2}$ "	8"
<b>F-5360-12</b>	22"	22"	16"	19"	12"	6"
<b>F-5361-12</b>	22"	25"	16"	22"	12"	8"
<b>F-5362-12</b>	22"	31"	16"	28"	12"	8"


## ADA Compliant

Part Number	O.D.		I.D.		Sink Depth
	Front to Back	Side to Side	Front to Back	Side to Side	
<b>F-5353-06</b>	15"	15"	12"	12"	5 $\frac{1}{2}$ "
<b>F-5354-06</b>	17"	19"	14"	16"	5 $\frac{1}{2}$ "
<b>F-5355-06</b>	19"	22"	16"	19"	5 $\frac{1}{2}$ "
<b>F-5356-06</b>	19"	25"	16"	22"	5 $\frac{1}{2}$ "
<b>F-5357-06</b>	19"	31"	16"	28"	5 $\frac{1}{2}$ "



## ADA Compliant with Fixture Ledge





Part Number	O.D.		I.D.		Sink Depth	Fixture Centers
	Front to Back	Side to Side	Front to Back	Side to Side		
<b>F-5358-06</b>	18"	15"	12"	12"	5 $\frac{1}{2}$ "	4"
<b>F-5359-06</b>	20"	19"	14"	16"	5 $\frac{1}{2}$ "	6"
<b>F-5360-06</b>	22"	22"	16"	19"	5 $\frac{1}{2}$ "	6"
<b>F-5361-06</b>	22"	25"	16"	22"	5 $\frac{1}{2}$ "	8"
<b>F-5362-06</b>	22"	31"	16"	28"	5 $\frac{1}{2}$ "	8"















# Laboratory Electrical Fixtures

## Electrical Receptacles

All flush mount fixtures furnished complete with device, stainless steel face plate and electrical box. Pedestals are furnished with device, wrinkled black epoxy powder coated pedestal box, stainless steel face plate, tank nipple and locknut. All devices and boxes are UL listed. (Ship loose)









Type 2 and Type 3 pedestal boxes furnished with 3/4" IPS tank nipple, Type 7 and Type 8 pedestal boxes furnished with 1" IPS tank nipple.		Flush mount single gang box	Single Wide Pedestal Type 2 – Single Faced Type 3 – Double Faced			Flush mount double gang box	Double Wide Pedestal Type 7 – Single Faced Type 8 – Double Faced	
Electrical Device	Color	Type 1	Type 2	Type 3	Type 6	Type 7	Type 8	
<b>Specification Grade</b> 120 Volt AC, 20 amp, 3-wire polarized duplex receptacle NEMA 5-20  	Black Ivory White Grey Red	<b>S581-1K</b> <b>S581-1V</b> <b>S581-1W</b> <b>S581-1G</b> <b>S581-1R</b>	<b>S581-2K</b> <b>S581-2V</b> <b>S581-2W</b> <b>S581-2G</b> <b>S581-2R</b>	<b>S581-3K</b> <b>S581-3V</b> <b>S581-3W</b> <b>S581-3G</b> <b>S581-3R</b>	<b>S581-6K</b> <b>S581-6V</b> <b>S581-6W</b> <b>S581-6G</b> <b>S581-6R</b>	<b>S581-7K</b> <b>S581-7V</b> <b>S581-7W</b> <b>S581-7G</b> <b>S581-7R</b>	<b>S581-8K</b> <b>S581-8V</b> <b>S581-8W</b> <b>S581-8G</b> <b>S581-8R</b>	
<b>Specification Grade</b> 120 Volt AC, 20 amp, 3-wire polarized duplex receptacle NEMA 5-20 	Black Ivory White Grey Red	<b>S681-1K</b> <b>S681-1V</b> <b>S681-1W</b> <b>S681-1G</b> <b>S681-1R</b>	<b>S681-2K</b> <b>S681-2V</b> <b>S681-2W</b> <b>S681-2G</b> <b>S681-2R</b>	<b>S681-3K</b> <b>S681-3V</b> <b>S681-3W</b> <b>S681-3G</b> <b>S681-3R</b>	<b>S681-6K</b> <b>S681-6V</b> <b>S681-6W</b> <b>S681-6G</b> <b>S681-6R</b>	<b>S681-7K</b> <b>S681-7V</b> <b>S681-7W</b> <b>S681-7G</b> <b>S681-7R</b>	<b>S681-8K</b> <b>S681-8V</b> <b>S681-8W</b> <b>S681-8G</b> <b>S681-8R</b>	
<b>Hospital Grade</b> 120 Volt AC, 20 amp, 3-wire polarized duplex receptacle NEMA 5-20  	Black Ivory White Grey Red	<b>H581-1K</b> <b>H581-1V</b> <b>H581-1W</b> <b>H581-1G</b> <b>H581-1R</b>	<b>H581-2K</b> <b>H581-2V</b> <b>H581-2W</b> <b>H581-2G</b> <b>H581-2R</b>	<b>H581-3K</b> <b>H581-3V</b> <b>H581-3W</b> <b>H581-3G</b> <b>H581-3R</b>	<b>H581-6K</b> <b>H581-6V</b> <b>H581-6W</b> <b>H581-6G</b> <b>H581-6R</b>	<b>H581-7K</b> <b>H581-7V</b> <b>H581-7W</b> <b>H581-7G</b> <b>H581-7R</b>	<b>H581-8K</b> <b>H581-8V</b> <b>H581-8W</b> <b>H581-8G</b> <b>H581-8R</b>	
<b>Specification Grade</b> 240 Volt AC, 20 amp, 3-wire polarized duplex receptacle NEMA 6-20  	Black Ivory White Grey Red	<b>S582-1K</b> <b>S582-1V</b> <b>S582-1W</b> <b>S582-1G</b> <b>S582-1R</b>	<b>S582-2K</b> <b>S582-2V</b> <b>S582-2W</b> <b>S582-2G</b> <b>S582-2R</b>	<b>S582-3K</b> <b>S582-3V</b> <b>S582-3W</b> <b>S582-3G</b> <b>S582-3R</b>	<b>S582-6K</b> <b>S582-6V</b> <b>S582-6W</b> <b>S582-6G</b> <b>S582-6R</b>	<b>S582-7K</b> <b>S582-7V</b> <b>S582-7W</b> <b>S582-7G</b> <b>S582-7R</b>	<b>S582-8K</b> <b>S582-8V</b> <b>S582-8W</b> <b>S582-8G</b> <b>S582-8R</b>	
<b>Ground Fault Protected</b> 120 Volt AC, 20 amp, specification grade, duplex receptacle NEMA 5-20 	Black Ivory White Grey Red	<b>S656-1K</b> <b>S656-1V</b> <b>S656-1W</b> <b>S656-1G</b> <b>S656-1R</b>	<b>S656-2K</b> <b>S656-2V</b> <b>S656-2W</b> <b>S656-2G</b> <b>S656-2R</b>	<b>S656-3K</b> <b>S656-3V</b> <b>S656-3W</b> <b>S656-3G</b> <b>S656-3R</b>				
<b>Ground Fault Protected</b> 120 Volt AC, 20 amp, duplex receptacle and 3-wire polarized duplex receptacle NEMA 5-20  	Black Ivory White Grey Red				<b>S657-6K</b> <b>S657-6V</b> <b>S657-6W</b> <b>S657-6G</b> <b>S657-6R</b>	<b>S657-7K</b> <b>S657-7V</b> <b>S657-7W</b> <b>S657-7G</b> <b>S657-7R</b>	<b>S657-8K</b> <b>S657-8V</b> <b>S657-8W</b> <b>S657-8G</b> <b>S657-8R</b>	
<b>Ground Fault Protected</b> 120 Volt AC, 20 amp, duplex receptacle and 3-wire polarized duplex receptacle NEMA 5-20  	Black Ivory White Grey Red				<b>S658-6K</b> <b>S658-6V</b> <b>S658-6W</b> <b>S658-6G</b> <b>S658-6R</b>	<b>S658-7K</b> <b>S658-7V</b> <b>S658-7W</b> <b>S658-7G</b> <b>S658-7R</b>	<b>S658-8K</b> <b>S658-8V</b> <b>S658-8W</b> <b>S658-8G</b> <b>S658-8R</b>	

# Laboratory Electrical Fixtures

## Electrical Switches and Data Ports

All flush mount fixtures furnished complete with device, stainless steel face plate and electrical box. All devices and boxes are UL listed. (Ship loose)



Electrical Device	Color	Flush mount single gang box		Flush mount double gang box		
		Cat-6	Cat-6a	Cat-6	Cat-6a	
<b>Specification Grade</b> Single pole, 120/240 Volt AC, 20 amp toggle switch. 	Black Ivory White Grey Red	<b>S655-1K</b> <b>S655-1V</b> <b>S655-1W</b> <b>S655-1G</b> <b>S655-1R</b>				
<b>RJ-45 Data Port Duplex</b> Cat-6 Data Port 	Black Ivory White Grey Red		<b>0720-1K-6</b> <b>0720-1V-6</b> <b>0720-1W-6</b> <b>0720-1G-6</b> <b>0720-1R-6</b>	<b>0720-1K-6a</b> <b>0720-1V-6a</b> <b>0720-1W-6a</b> <b>0720-1G-6a</b> <b>0720-1R-6a</b>		
<b>RJ-45 Data Port Duplex</b> Cat-6 Data Port with 120 Volt AC, 20 amp, specification grade, duplex receptacle. NEMA 5-20 	Black Ivory White Grey Red				<b>S721-6K-6</b> <b>S721-6V-6</b> <b>S721-6W-6</b> <b>S721-6G-6</b> <b>S721-6R-6</b>	<b>S721-6K-6a</b> <b>S721-6V-6a</b> <b>S721-6W-6a</b> <b>S721-6G-6a</b> <b>S721-6R-6a</b>
<b>RJ-45 Data Port Duplex</b> Cat-6 Data Port with 120 Volt AC, 20 amp, specification grade, GFI duplex receptacle. NEMA 5-20 	Black Ivory White Grey Red				<b>S751-6K-6</b> <b>S751-6V-6</b> <b>S751-6W-6</b> <b>S751-6G-6</b> <b>S751-6R-6</b>	<b>S751-6K-6a</b> <b>S751-6V-6a</b> <b>S751-6W-6a</b> <b>S751-6G-6a</b> <b>S751-6R-6a</b>

# Laboratory Seating

**Specifications:**

The seating shown on this page features chrome-plated metal parts. Upholstery is No. 684 Black Vinyl. Soft

Touch pneumatic mechanism allows easy adjustment to the seating height. Although not shown in the photographs

below, models F-4940-00 thru F-4949-00 come with 2" soft wheel ball bearing casters for hard floors.



- Upholstered seat/backrest
- 5-legged base with 18" fixed footring
- Tubular base

Model No.	Seat Height Adjustment	Footring Height
<b>F-4940-00</b>	16 <sup>3</sup> / <sub>4</sub> "-22"	8 <sup>1</sup> / <sub>2</sub> "
<b>F-4941-00</b>	22 <sup>1</sup> / <sub>4</sub> "-27 <sup>1</sup> / <sub>2</sub> "	8 <sup>1</sup> / <sub>2</sub> "
<b>F-4942-00</b>	26 <sup>1</sup> / <sub>4</sub> "-31 <sup>1</sup> / <sub>2</sub> "	12 <sup>1</sup> / <sub>2</sub> "



- Upholstered seat/backrest
- 5-legged base with 18" fixed foot ring
- Tubular base
- Backrest with lumbar support

Model No.	Seat Height Adjustment	Footring Height
<b>F-4943-00</b>	16 <sup>3</sup> / <sub>4</sub> "-22"	8 <sup>1</sup> / <sub>2</sub> "
<b>F-4944-00</b>	22 <sup>1</sup> / <sub>4</sub> "-27 <sup>1</sup> / <sub>2</sub> "	8 <sup>1</sup> / <sub>2</sub> "
<b>F-4945-00</b>	26 <sup>1</sup> / <sub>4</sub> "-31 <sup>1</sup> / <sub>2</sub> "	12 <sup>1</sup> / <sub>2</sub> "



# Laboratory Seating

**Specifications:**

The seating shown on this page features chrome-plated metal parts. Upholstery is No. 684 Black Vinyl. Soft

Touch pneumatic mechanism allows easy adjustment to the seating height. Although not shown in the photographs

below, models F-4940-00 thru F-4949-00 come with 2" soft wheel ball bearing casters for hard floors.



- 13" round chrome-plated seat
- 5-legged base with 18" fixed foot ring
- Tubular base

Model No.	Seat Height Adjustment	Footring Height
<b>F-4946-00</b>	18 <sup>1</sup> / <sub>2</sub> "-22 <sup>3</sup> / <sub>4</sub> "	8 <sup>1</sup> / <sub>2</sub> "
<b>F-4947-00</b>	22 <sup>1</sup> / <sub>2</sub> "-27 <sup>3</sup> / <sub>4</sub> "	12 <sup>1</sup> / <sub>2</sub> "



- 15" round upholstered seat
- 5-legged base with 18" fixed foot ring
- Tubular base

Model No.	Seat Height Adjustment	Footring Height
<b>F-4948-00</b>	21 <sup>1</sup> / <sub>2</sub> "-24 <sup>3</sup> / <sub>4</sub> "	8 <sup>1</sup> / <sub>2</sub> "
<b>F-4949-00</b>	25 <sup>1</sup> / <sub>2</sub> "-20 <sup>3</sup> / <sub>4</sub> "	12 <sup>1</sup> / <sub>2</sub> "




- 13" round wood seat
- Tubular 4-legged base with 16" fixed foot ring
- Black tubular steel base


Model No.	Seat Height
<b>F-5350-00</b>	18"
<b>F-5351-00</b>	24"
<b>F-5352-00</b>	30"

# Technical Information


## Wall Bench End Views



**24" Deep Wall Bench**




**30" Deep Wall Bench**




**31" Deep Wall Bench**

## Center Bench End Views




**Center Bench End/Section View**



**Center Bench End Sink Cabinet**

# Technical Information

## Full Height Storage Cabinet End Views




12<sup>3</sup>/<sub>4</sub>" Deep Full Height Cabinet

16<sup>3</sup>/<sub>4</sub>" Deep Full Height Cabinet

22" Deep Full Height Cabinet


## Center Bench Plan View




**Plan View**  
Center Bench End Sink Cabinet

# Technical Information

## X01 - X02 - X03 Octagon All-science Table Dimensions




## X05 Perimeter Annex Dimensions




# Technical Information

## X75 Perimeter Student Table Dimensions





## X80 Perimeter Student Table Dimensions




# Technical Information

## Student Station Work Top Dimensions

## Pedestal Work Top Dimensions




## Instructor Station Work Top Dimensions




# Technical Information

## Example Classroom Layouts

Layout 1




Layout 2




Drawings based on 28' x 34' classroom

Layout 3



Layout 4




Drawings based on 28' x 34' classroom


# Technical Information

## Example Classroom Layouts

Layout 5




Layout 6




Drawings based on 28' x 34' classroom

Layout 7



Layout 8




Drawings based on 28' x 34' classroom




# Technical Information

## Example Classroom Layouts - Multiple Configurations of Fixed Sink Locations

Layout 9a




Layout 9b




Drawings based on 32' x 38' classroom

Layout 10a



Layout 10b



Drawings based on 32' x 38' classroom

# Alphabetical Index

Description	Page	Description	Page
<b>A</b>		<b>C</b>	
<b>Acid Storage Cabinet</b>		<b>Cubbie Cabinets</b>	
ADA Height . . . . .	23	ADA Height . . . . .	18
Full Height . . . . .	80	Full Height . . . . .	74
Standing Height . . . . .	39	Sitting Height . . . . .	14
ADA Compliant Sink . . . . .	136-137	Standing Height . . . . .	28
ADA Height Cabinet . . . . .	18-27	Wall Mounted . . . . .	68
ADA Mobile Workstation . . . . .	97	<b>Cupboard Cabinets</b>	
ADA Safety Cabinet . . . . .	104	ADA Height . . . . .	18
ADA Aprons . . . . .	26	Mobile . . . . .	44-45
Adjustable Height Tables . . . . .	27	Sitting Height . . . . .	14
Aerator . . . . .	129	Standing Height . . . . .	29
Air Valves . . . . .	125-126	<b>Cupboard/Drawer Cabinets</b>	
All-Science Student Table . . . . .	100-102	ADA Height . . . . .	19
<b>Aprons</b>		Mobile . . . . .	44-45
ADA Aprons . . . . .	26	Sitting Height . . . . .	16
Adjustable Aprons . . . . .	61	Standing Height . . . . .	32-34
Apron Rails . . . . .	61	Cupsinks . . . . .	133
Box Aprons . . . . .	48-49	<b>D</b>	
Four-legged Table Aprons . . . . .	54-57	Deionized Water Faucet . . . . .	130
Legs on Left End . . . . .	50-51	Demonstration Hood . . . . .	105
Legs on Right End . . . . .	52-53	<b>Diagonal Corner Cabinets</b>	
One-legged Aprons . . . . .	58-60	ADA Height . . . . .	21
Aspirator . . . . .	129	Sitting Height . . . . .	17
Assemblies . . . . .	84	Standing Height . . . . .	37
<b>B</b>		Display Cabinet . . . . .	77
Balance Table . . . . .	121	Distilled Water Faucet . . . . .	130
Ball Valve . . . . .	125	Drain Outlet . . . . .	133
Base Molding . . . . .	123	Drain Trough . . . . .	134
<b>Blind Corner Cabinets</b>		Drawer Bottom Option . . . . .	11
ADA Height . . . . .	21	<b>Drawer Cabinets</b>	
Sitting Height . . . . .	17	ADA Height . . . . .	18
Standing Height . . . . .	37	Sitting Height . . . . .	14-15
<b>C</b>		Standing Height . . . . .	30-31
Cabinet & Assembly Dimensions . . . . .	142-147	Drawer Storage Cabinet - Full Height . . . . .	81
Cabinet Levelers . . . . .	9	Drawer Slide Option . . . . .	11
Cabinet Options . . . . .	4-11	Drench Shower . . . . .	132
Cabinet Styles . . . . .	3	Duo-Vision Demonstration Hood . . . . .	106-109
<b>Carts</b>		Duralumin Rods . . . . .	121
Stainless Steel . . . . .	121	<b>E</b>	
Vacuum Pump . . . . .	24-40	Educational Assemblies . . . . .	84-109
Wood . . . . .	95	Eight-Sided Tables . . . . .	100-102
Casters . . . . .	6	Electrical Fixtures . . . . .	138-139
Catalog Number Explanation . . . . .	4	Emergency Eye Wash/Shower . . . . .	132
Chalkboard . . . . .	121	End View Dimensions . . . . .	142-143
Cold Water Faucet . . . . .	127	Epoxy Resin Cupsink . . . . .	133
Combination Faucet . . . . .	126	Epoxy Resin Drain Trough . . . . .	134
Construction Details . . . . .	2	Epoxy Resin Top . . . . .	112
Contemporary Full Overlay . . . . .	3	Epoxy Resin Drop-in Sink . . . . .	136
Conventional Radius Lip . . . . .	3	Epoxy Resin Tub Sink . . . . .	135
<b>Corner Cabinets</b>		Extended End Panels . . . . .	9
ADA Height . . . . .	21	Eye Wash . . . . .	132
Sitting Height . . . . .	17	<b>F</b>	
Standing Height . . . . .	37	Face Wash . . . . .	132
Counter Mounted Cabinets . . . . .	68-73	Faucet Accessories . . . . .	129
Crumb Cup Strainer . . . . .	133	Filler Panels . . . . .	82-83
		Filter Pump . . . . .	129
		Finished Back . . . . .	6
		Finished Panel . . . . .	82
		Fixtures - Electrical . . . . .	138-139

# Alphabetical Index (continued)

Description	Page	Description	Page
<b>F</b> Foot Valve . . . . .	129	<b>L</b> Lock Keying . . . . .	12-13
Freezer . . . . .	121	Lock Panels . . . . .	8
		Locks . . . . .	8
<b>Full Height Storage Cabinets</b>		<b>M</b> <b>Matrix</b>	
Acid Storage . . . . .	80	Instructor Station . . . . .	86
Cubbie Storage . . . . .	74	Student Stations . . . . .	88-89
Display Cabinet . . . . .	77	Pedestals . . . . .	90-91
Drawer Storage . . . . .	81	Markerboard . . . . .	121
Garment Locker . . . . .	80	Matched Vertical Grain . . . . .	3
Glass Tubing Storage . . . . .	81	Microscope Storage Cabinet . . . . .	79
Instructor Storage . . . . .	79	Mixing Faucet . . . . .	128-129
Microscope Storage . . . . .	79	<b>Mobile</b>	
Open Storage . . . . .	75	Cabinets . . . . .	44-47
Skeleton Storage . . . . .	81	Demonstration Bench . . . . .	96
Sliding Glazed Doors . . . . .	75	Workstations . . . . .	97
Sliding Panel Doors . . . . .	75	Multi-Science Tables . . . . .	102
Split Swinging Doors . . . . .	78	Multiple Service Fittings . . . . .	126
Swinging Glazed Doors . . . . .	76	<b>N</b> Needle Point Cock Valves . . . . .	125
Swinging Panel Doors . . . . .	77	<b>O</b> Octagon Table . . . . .	100-102
Torso Storage . . . . .	81	<b>Open Cabinets</b>	
Tote Tray Storage . . . . .	78	ADA Height . . . . .	18
Full Radius Edge . . . . .	3	Sitting Height . . . . .	14
Full Top Frame Option . . . . .	6	Standing Height . . . . .	28
<b>Fume Hood Base Cabinets</b>		<b>P</b> Padlock Hasps . . . . .	8
ADA Height . . . . .	23-25	Panels . . . . .	82
Standing Height . . . . .	38-41	Pedestal Legs . . . . .	64
Fume Hoods . . . . .	105-111	Pegboard . . . . .	120
<b>G</b> Garment Locker Cabinet . . . . .	80	Perimeter Annex . . . . .	100
Gas Valves . . . . .	125-126	Perimeter Bench . . . . .	98-99
Glass Tubing Storage . . . . .	81	Pipe Enclosure . . . . .	122
Gooseneck . . . . .	129	Plastic Laminate Top . . . . .	112
<b>H</b> Hardwood Top . . . . .	112	Plumbing Information . . . . .	124
Hinge Options . . . . .	2	Plywood Drawer Bottom . . . . .	11
Hoods . . . . .	105-111	Portable Pump Faucet . . . . .	127
Hose Spray . . . . .	132	Pullout Tray . . . . .	10
Hot & Cold Water Mixing Faucet . . . . .	128-129	Pull Styles . . . . .	2
<b>I</b> Instructor Desk . . . . .	86 & 92	<b>R</b> Reagent Rack Assemblies . . . . .	118
Instructor Desk Cabinets . . . . .	42	Rear Filler . . . . .	82
Instructor Storage Cabinet . . . . .	79	Receptacle - Electric . . . . .	138-139
Intermediate Rail Option . . . . .	11	Refrigerator . . . . .	121
<b>K</b> <b>Kemresin</b>		Removable Back . . . . .	7
Cupsink . . . . .	133	Rod Clamps . . . . .	123
Drop-in Sink . . . . .	136	Rod Socket . . . . .	123
Top . . . . .	112	Rods . . . . .	123
Trough . . . . .	134	<b>S</b> Safety Cabinet . . . . .	104
Tub Sink . . . . .	135	Safety Shower . . . . .	132
Kemstrut . . . . .	122	Seating . . . . .	140-141
Key Case . . . . .	121	Service Fittings . . . . .	125-131
Keyboard Tray . . . . .	65	Shelf Brackets . . . . .	119
Keying Options . . . . .	12-13	Shelf Lip . . . . .	9
Knee Space Panel . . . . .	65	Shelving . . . . .	114-119
<b>L</b> Label Holder . . . . .	123	Silhouette . . . . .	3
Label Holder Option . . . . .	7	Single Water Faucet . . . . .	127
Laboratory Cart . . . . .	95 & 121	Sink . . . . .	134-137
Laboratory Service Fittings . . . . .	125-132		
Leg Assemblies . . . . .	62-64		
Levelers . . . . .	9		

# Alphabetical Index

<b>Description</b>	<b>Page</b>	<b>Description</b>	<b>Page</b>
<b>S</b>			
<b>Sink Cabinets</b>		Tote Tray Storage - Full Height . . . . .	78
ADA Height . . . . .	20 & 22	Trap . . . . .	133
Sitting Height . . . . .	17	Trough . . . . .	134
Standing Height . . . . .	35	TruView Fume Hood & Accessories . . . . .	106-109
Sink Supports . . . . .	17 / 20 / 35	Tub Sink . . . . .	134
Sink Outlet . . . . .	133	<b>U</b>	
Sink Overflow . . . . .	133	Ultra-Pure Water Faucet . . . . .	130
Sink Trap . . . . .	133	Under Counter Refrigerator . . . . .	121
Sitting Height Base Cabinets . . . . .	14-17	Upright Rod Assembly . . . . .	123
Skeleton Storage Cabinet . . . . .	81	<b>V</b>	
Soffit . . . . .	82	Vacuum Breaker . . . . .	129
<b>Solvent Storage Cabinet</b>		<b>Vacuum Pump Cabinets</b>	
ADA Height . . . . .	24	ADA Height . . . . .	25
Standing Height . . . . .	40	Carts . . . . .	25 & 41
Special Purpose Base Cabinets . . . . .	42-43	Standing Height . . . . .	41
Spreaders . . . . .	62	Vacuum Valves . . . . .	125-126
<b>Stainless Steel</b>		Vandal-Resistant Fittings . . . . .	131
Carts . . . . .	121	Vibration Pads . . . . .	62
Cupsinks . . . . .	133	<b>W</b>	
Sinks . . . . .	137	<b>Wall Cabinets</b>	
Tops . . . . .	112	Corner . . . . .	73
Standing Height Base Cabinets . . . . .	28-43	Cubbie . . . . .	67-68
Steam Valves . . . . .	126	Open . . . . .	66
Stools . . . . .	140-141	Sliding Glazed Door . . . . .	70
Support Angle . . . . .	122	Sliding Panel Door . . . . .	70
Support Strut . . . . .	122	Sliding Plate Glass Door . . . . .	69
Switch . . . . .	139	Swinging Glazed Door . . . . .	71
<b>T</b>		Swinging Panel Door . . . . .	72
Tables - Adjustable Height . . . . .	27	Waste Storage Cabinet . . . . .	36
Tackboard . . . . .	121	Water Faucet . . . . .	126-131
Teaching Hood . . . . .	105-109	Water Mixing Faucet . . . . .	128-129
Tempered Glass . . . . .	10	Water Valves . . . . .	126
Toe space Filler . . . . .	82	Water/Gas Fitting . . . . .	126
Toe space Grill . . . . .	7	White Drawer Bottom . . . . .	11
Torso Storage Cabinet . . . . .	81	Work Top Support Frame . . . . .	65
Tote Tray . . . . .	123	Work Tops . . . . .	112
Tote Tray Storage - Base Cabinet . . . . .	36	Wheelchair Accessible Base Cabinets . . . . .	22

# Numerical Index

Part Number	Page	Part Number	Page	Part Number	Page	
0464-00	- 0469-00 . . . . .	133	A64W052224	- A64W052836 . . . . .	59	
0482-BP	- 0482-GP . . . . .	133	A65W053418	- A65W054036F . . . . .	60	
0490-BE	- 0490-SE . . . . .	133	A80W002230	- A80W002248 . . . . .	65	
0491-BP	. . . . .	109				
0491-BP	- 0491-GP . . . . .	133	B00W312215	- B00W312248 . . . . .	14	
0492-BP	- 0492-GP . . . . .	133	B00W342215	- B00W342248 . . . . .	18	
0493-BE	- 0493-SE . . . . .	133	B00W362215	- B00W362248 . . . . .	28	
0494BP-04	- 0494BP-10 . . . . .	133	B01W312212	- B01W312248 . . . . .	14	
0548-00	. . . . .	123	B01W342212	- B00W342248 . . . . .	18	
0549-18	- 0557-00 . . . . .	123	B01W362212	- B01W362248 . . . . .	28	
0720-1G-6	- 0720-1W-6a . . . . .	139				
0950-0F	- 0975-0F . . . . .	133	C03W312215	- C03W312248 . . . . .	14	
0995-00	. . . . .	133	C03W342215	- C03W342248 . . . . .	18	
1000-00	- 1039-00 . . . . .	135	C03W362215	- C03W362248 . . . . .	29	
1000-ADA-DI	- 1006-ADA-DI . . . . .	136	C03W362230	- C03W362236 . . . . .	108	
1000-DI	- 1070-DI . . . . .	136	C06W362230	- C06W362248 . . . . .	29	
4740-00	. . . . .	123	C07W362212	- C07W362224 . . . . .	29	
A00W032224A	- A00W052272A . . . . .	61	C15W362230	- C15W362248 . . . . .	29	
A00W0324	- A00W0572F . . . . .	61				
A00W0330	- A00W0336 . . . . .	108	D00W312215	- D00W312248 . . . . .	14	
A00W052224	- A00W052272F . . . . .	48	D00W362212	- D00W362248 . . . . .	30	
A01W052224	- A01W052248F . . . . .	48	D02W312230	- D02W312248 . . . . .	14	
A01W052224A	- A01W052272A . . . . .	61	D02W362230	- D02W362248 . . . . .	30	
A02W052248	- A02W052272F . . . . .	49	D03W362236	- D03W362248 . . . . .	30	
A03W052224	- A03W052248 . . . . .	49	D08W312215	- D08W312224 . . . . .	15	
A04W052248	- A04W052272 . . . . .	49	D10W312215	- D10W312224 . . . . .	15	
A05W052236	- A05W052872 . . . . .	26	D20W312218	- D20W312224 . . . . .	15	
A08W052224A	- A08W052272A . . . . .	61	D26W312218	- D26W312224 . . . . .	15	
A09W0335	- A09W0371 . . . . .	62	D27W312224	. . . . .	15	
A10W052224	- A10W052872F . . . . .	50	D28W312218	- D28W312224 . . . . .	15	
A11W052236	- A11W052872 . . . . .	26	D29W312224	. . . . .	15	
A13W052236	- A13W052872 . . . . .	26	D30W312215	- D30W312248 . . . . .	15	
A15W052224	- A15W052848F . . . . .	50	D30W342215	- D30W342224 . . . . .	18	
A16W052248	- A16W052872F . . . . .	51	D30W362212	- D30W362248 . . . . .	30	
A17W052224-	A17W052848F . . . . .	52	D30W362218	- D30W362224 . . . . .	108	
A18W052248	- A18W052872F . . . . .	53	D31W312230	- D31W312248 . . . . .	15	
A20W3100	- A20W3600F . . . . .	62	D31W362230	- D31W362248 . . . . .	30	
A20W3400	. . . . .	26	D32W312230	- D32W312248 . . . . .	15	
A21W3122	- A21W3640F . . . . .	63	D32W362230	- D32W362248 . . . . .	30	
A21W3422	- A21W3428 . . . . .	26	D33W362236	- D33W362248 . . . . .	30	
A24W3100	- A24W3600 . . . . .	62	D35W362230	- D35W362248 . . . . .	30	
A25W3122	- A25W3630F . . . . .	64	D36W362230	- D36W362248 . . . . .	30	
A25W3430	- A25W3636 . . . . .	108	D37W362230	- D37W362248 . . . . .	31	
A26W3151	- A26W3659F . . . . .	64	D38W362215	- D38W362224 . . . . .	31	
A26W3459	- A26W3659 . . . . .	108	D39W362230	- D39W362248 . . . . .	31	
A30W052235	- A30W052871F . . . . .	54	D40W362230	- D40W362248 . . . . .	31	
A31W052235	- A31W052847F . . . . .	54	D46W362248	. . . . .	31	
A32W052247	- A32W052871F . . . . .	55	D50W362215	- D50W362224 . . . . .	31	
A33W052235-	A33W052847 . . . . .	55	D51W362230	- D51W362248 . . . . .	31	
A34W052247	- A34W052871 . . . . .	55	D52W362230	- D52W362248 . . . . .	31	
A35W053447	- A35W054071F . . . . .	56	D81W312218	- D81W312224 . . . . .	15	
A38W053447	- A38W054071 . . . . .	57	D81W362218	- D81W362224 . . . . .	31	
A39W053447	- A39W054071 . . . . .	57	D84W362218	- D84W362224 . . . . .	31	
A40W3100	- A40W3600F . . . . .	62	D85W312218	- D85W312224 . . . . .	15	
A40W3400	. . . . .	26	D86W362248	. . . . .	31	
A41W3122	- A41W3640F . . . . .	63	D87W362218	- D87W362248 . . . . .	31	
A41W3422	- A41W3428 . . . . .	26	D88W362236	- D88W362248 . . . . .	31	
A44W3122	- A44W3640F . . . . .	63	D90W312218	- D90W312224 . . . . .	15	
A44W3422	- A44W3428 . . . . .	26				
A53W052224	- A53W052848 . . . . .	51	E00W362212	- E00W362224 . . . . .	32	
A54W052224	- A54W052848 . . . . .	53	E03W362230	- E03W362248 . . . . .	32	
A60W052212	- A60W052836F . . . . .	58	E05W362236	- E05W362248 . . . . .	32	
A63W052224	- A63W052836 . . . . .	59	E24W312224	. . . . .	16	
				E25W312218	- E25W312224 . . . . .	16
				E40W312215	- E40W312248 . . . . .	16
				E40W342215	- E40W342248 . . . . .	19
				E40W362212	- E40W362248 . . . . .	33
				E40W362236	- E40W362248 . . . . .	32
				E41W312230	- E41W312248 . . . . .	16
				E41W342230	- E41W342248 . . . . .	19
				E41W362230	- E41W362248 . . . . .	33
				E42W362230	- E42W362248 . . . . .	33
				E43W342230	- E43W342248 . . . . .	19
				E43W362230	- E43W362248 . . . . .	33
				E44W312230	- E44W312248 . . . . .	16
				E44W362230	- E44W362248 . . . . .	33
				E45W342230	- E45W342248 . . . . .	19
				E45W362230	- E45W362248 . . . . .	33
				E52W362215	- E52W362248 . . . . .	34
				E53W362230	- E53W362248 . . . . .	34
				E56W362230	- E56W362248 . . . . .	34
				E60W362218	- E60W362224 . . . . .	34
				E61W362236	- E61W362248 . . . . .	34
				E63W362230	- E63W362248 . . . . .	34
				E65W362230	- E65W362248 . . . . .	34
				E86W312224	. . . . .	16
				F-0260-00	. . . . .	123
				F0862-04	. . . . .	123
				F10M0242-00	- F10M0279-00 . . . . .	116
				F11M0196-00	. . . . .	116
				F11M0224-00	- F11M0272-00 . . . . .	116
				F11M0324-00	- F11M0372-00 . . . . .	116
				F12M000224	- F12S000284 . . . . .	114
				F-1368-00	. . . . .	123
				F13M000224	- F13S000284 . . . . .	114
				F17-3612-00	- F17-3612 . . . . .	114
				F17-4824-00	- F17-4872 . . . . .	114
				F21M0624-00	- F21W2448-08 . . . . .	117
				F21M0924-60	- F21W1848-68 . . . . .	115
				F-2939-00	- F-2940-00 . . . . .	120
				F-4940-00	- F-4945-00 . . . . .	140
				F-4946-00	- F-4949-00 . . . . .	141
				F-5332-00	- F-5332-01 . . . . .	121
				F-5336-00	- F-5344-00 . . . . .	121
				F-5349-00	. . . . .	121
				F-5350-00	- F-5352-00 . . . . .	141
				F-5353-08	- F-5362-12 . . . . .	137
				F-5418-24	- F-5436-48 . . . . .	121
				F63M0224-00	- F63M0272-00 . . . . .	115
				F-6853-00	. . . . .	65
				F73M0015	- F73S0024 . . . . .	47
				F94P2215	- F94T2448 . . . . .	47
				G00W312218	- G00W312248 . . . . .	17
				G00W342212	- G00W362218 . . . . .	108
				G00W342215	- G00W342248 . . . . .	20
				G00W362215	- G00W362248 . . . . .	35

# Numerical Index (continued)

<b>Part Number</b>	<b>Page</b>	<b>Part Number</b>	<b>Page</b>	<b>Part Number</b>	<b>Page</b>
G01W362260	35	K7-302435-BK - K7302435-SL	121	S655-1G - S655-1W	139
G03W342236	22			S655-1V	109
G04W362251	35	KF-0170-04	119	S656-1G - S656-3W	138
G07W342236	22	KTS-070	62	S656-1V	109
G08W342015	23			S657-6G - S657-8W	138
G08W362012	38	K-WSS-15 - K-WSS-48	17	S658-6G - S658-8W	138
		K-WSS-15 - K-WSS-48	20	S681-1G - S681-8W	138
G10W342212	108	K-WSS-15 - K-WSS-60	35		
				S70W841624 - S70W842248	80
G12W312248	17	L-105 - L-115	121	S721-6G-6 - S721-6W-6a	139
G12W342248	21	L-2034-06 - L-2034-96	122		
G12W362248	37			S72W841648 - S72W842248	77
G13W312248	17	M-0031-0A - M-0062-AD	119	S751-6G-6 - S751-6W-6a	139
G13W342248	21	M-0051-0A - M-0052-0A	122		
G13W362248	37			S85W842236 - S85W842248	81
G14W312215	17	M10W272215 - M10W272224	44	SM-01-00 - SM-56-00	135
G14W342215	21	M10W332215 - M10W332224	45		
G14W362215	37	M15W272215 - M15W272224	44	T00WC - T00WT	95
G15W342215	21	M15W332215 - M15W332224	45	T01WC - T01WT	95
G15W362215	37	M20W272215 - M20W272224	44	T20WC - T20WT	92
G16W342215	21	M20W332215 - M20W332224	45	T21WC - T21WT	93
G16W362215	37	M25W272215 - M25W272224	44	T24WC - T24WT	94
G20W362218	36	M25W332215 - M25W332224	45	T25WC - T25WT	94
		M30W272215 - M30W272224	44	T30WC - T30WT	92
G21W342215	20	M35W272215 - M35W272224	44	T31WC - T31WT	93
G21W362215	36	M40M272215 - M40W272224	44	T34WC - T34WT	94
G23W362218	36	M45W332216 - M45W332224	45	T35WC - T34WT	94
G35W342024	25	M50W332215 - M50W332224	45	T44WC - T44WT	103
G35W362024	41	M60W272215 - M60W272224	44	T56WC - T56WT	96
G41W181521	102	M62W272215 - M62W272224	44	T57WC - T57WT	97
G51W201521	102	M62W332215 - M62W332224	45	T84WC - T84WT	103
G63W342024	24	M63W332215 - M63W332224	45	T90M3048-AD - T90M6072-SA	107
G63W362024	40	M65W272215 - M65W272224	44	T90R3048 - T90R6072	109
G68W342024	24	M67W332215 - M67W332224	45	T95G483058-FS - T94G483058-SE	105
G68W362024	40	M75W332215 - M75W332224	45	T95WC - T95WT	105
G71W342818	42	M90W332215 - M90W332224	46		
G72W342818	42	M95W332215 - M95W332224	46	TR-01-00 - TR-05-00	134
G73W342836	42			U01W600512 - U02W961048	122
G74W342836	42	PL-4812-00 - PL-7214-00	118		
G78W342836	42	RRU-1210-0000 - RRU-1812-0030	118	VPCC061414 - VPCC061644	25
G79W342818	42			VPCC061414 - VPCC061644	41
G80W342024	23			W00W241215 - W00W481648	66
G80W362024	39	S00W841218 - S00W842248	75	W01W241212 - W01W481630	67
		S01W841212 - S01W842230	74	W01W241236 - W01W481648	68
G81WC - G81FP	104	S11W841230 - S11W842248	75		
G85WC - G85WT	104	S12W841230 - S12W842248	75	W-0242-00	125
G86W362248	36	S20W841218 - S20W842248	76	W-0243-00	125
		S21W841624 - S21W842248	79	W-0259-00	126
G90W341618	43	S25W841218 - S25W842248	77	W-0260-00	125
G91W341618	43	S30W841230 - S30W842248	78	W-0263-00	125
G92W341618	43	S35W841230 - S35W842248	78	W-0266-00	125
G93W341618	43	S40W841624 - S40W842248	80	W-0272-00	125
G94W341618	43	S41W842236 - S41W842248	79	W-0273-02 - W-0274-VB	126
G95W341618	43			W-0286-02	125
G96W341618	43	S581-1G - S581-8W	138	W-0287-00	125
G97W341636	43	S582-1G - S582-8W	138	W-0300-00	126
G98W341636	43			W-0301-00	126
G99W341618	43	S60W841624	81	W-0306-00	126
		S61W841624 - S61W842224	81		
H581-1G - H581-8W	138	S62W842236 - S62W842248	78		
		S63W841624 - S63W842224	81		
K-0005-0A	122	S64W842236 - S64W842248	78		
K12-2905-0A - K12-3520-0A	122				

# Numerical Index (continued)

Part Number	Page	Part Number	Page	Part Number	Page
W-0333-00	128	W-0928-00	132	X94WC	- X94WT . . . . .90
W-0337-00	- W-0337-6V . . . . .127	W-0930-00	132	X95WC	- X95WT . . . . .91
W-0340-00	- W-0340-6V . . . . .128	W-0931-00	132	X96WC	- X96WT . . . . .91
W-0350-00	- W-351-0V . . . . .127	W-0932-00	132		
W-0354-00	- W-0354-0V . . . . .128	W-0933-00	132	X-WBC0431	- X-WBC0636 . . . . .83
W-0355-00	128	W-0934-00	132	X-WBP0431	- X-WBP1236L . . . . .83
W-0361-00	- W-0361-0V . . . . .129	W-0935-00	132	X-WC0236	- X-WC0636 . . . . .82
W-0367-00	130	W-0936-00	132	X-WFP0484	- X-WFP1284L . . . . .83
W-0368-00	130			X-WK2436	- X-WKR7236 . . . . .65
W-0374-01	132	W10W241230	- W10W481648 . . . . .69	X-WK3036	- X-WK3636 . . . . .108
W-0400-00	125	W11W241230	- W11W481648 . . . . .70	X-WP0136	- X-WP1636 . . . . .82
W-0401-00	125	W12W241230	- W12W481648 . . . . .70	X-WP0236	- X-WP1284 . . . . .82
W-0402-00	126	W13W301230	- W13W361648 . . . . .69	X-WP0836	- X-WP1436 . . . . .108
W-0403-00	126	W20W241215	- W20W481648 . . . . .71	X-WP3648	- X-WP4896 . . . . .82
W-0404-00	126	W21W241217	- W21W481623 . . . . .73	X-WS3648	- X-WS4896 . . . . .82
W-0405-00	127	W25W241215	- W25W491648 . . . . .72	X-WTP412	- X-WTP448 . . . . .82
W-0407-00	130	W26W241217	- W26W481623 . . . . .73	X-WWP0424	- X-WWP1248L . . . . .83
W-0408-00	130				
W-0409-00	130	X-018528-S	109		
W-0410-VR	131	X01WC	- X01WT . . . . .101		
W-0411-VR	131	X-020008	- X-020017 . . . . .120		
W-0412-VR	131	X02WC	- X02WT . . . . .100		
W-0413-VR	131	X-030001-00	- X-030003-00 . . . . .123		
W-0414-VR	131	X03WC	- X03WT . . . . .101		
W-0415-VR	131	X05WC	- X05WT . . . . .100		
W-0416-VB	- W-0416-VR . . . . .131				
W-0417-VR	131	X10WC302436	- X10WT013096 . . . . .27		
W-0510-00	129				
W-0532-00	129	X23WC	- X23WT . . . . .86		
W-0535-00	129	X33WC	- X33WT . . . . .86		
W-0744-A	- W-0744-W . . . . .109	X75WC	- X75WT . . . . .98		
W-0825-00	129	X80WC	- X80WT . . . . .98		
W-0833-00	129	X90WC	- X90WT . . . . .89		
W-0923-00	132	X91WC	- X91WT . . . . .88		
W-0926-00	132	X92WC	- X92WT . . . . .90		
W-0927-00	132	X93WC	- X93WT . . . . .91		

## Trademarks

Kemresin, Kemrock, Signature, Silhouette and Supreme Air are registered trademarks of Kewaunee Scientific Corporation.

Part numbers, illustrations, descriptions, specifications and other information contained in this catalog are based on the most current design and production information available

at the time of publication. Kewaunee Scientific Corporation reserves the right to make changes at any time without notice.

# Typical Installations





# Typical Installations



Kewaunee Scientific Corporation is dedicated to manufacturing high quality products for the laboratory marketplace. We offer **Total Laboratory Solutions, Innovative Designs, and Technological Expertise** with a worldwide distribution network to provide laboratories that are truly World Class.

**Laboratory Products from Kewaunee®:**

ALPHA® – Adaptable Casework System

ALPHA Columns – Expandable Workstations

ALPHA Overhead Carrier Systems

BasikBench® – Welded C-Leg Bench

DetectAir – Filtered/Ductless Fume Hoods

EarthLine® – Environmentally Friendly Casework

Enterprise® – Moveable Workstations

Evolution® – Column Based Workstations

Interceptor® – Biological Safety Cabinet

Interceptor® – Laminar Flow Cabinet

Kemresin® – Epoxy Resin Counters & Sinks

Kewaunee Matrix® – Educational Stations

Research Collection® – Steel Casework

Signature® Series – Wood Casework

Spektrum – Laminate Architectural Millwork

Sturdilite® – Flexible Pedestal System

Supreme Air® Venturi – Fume Hoods

For our complete product offering, specifications, catalogs, brochures, case studies, and dealer list, go to [www.kewaunee.com](http://www.kewaunee.com), or contact your local sales representative directly.



The mark of responsible forestry

Kewaunee, Alpha, BasikBench, EarthLine, Enterprise, Evolution, Interceptor, Kemresin, Kewaunee Matrix, Research Collection, Signature, Sturdilite, and Supreme Air are registered trademarks of Kewaunee Scientific Corporation



P.O. Box 1842 • Statesville, NC 28687-1842  
Phone: (704) 873-7202 • Fax: (704) 873-5160  
E-mail: [kscmarketing@kewaunee.com](mailto:kscmarketing@kewaunee.com)  
Web Site: [www.kewaunee.com](http://www.kewaunee.com)

**Kewaunee Labway India Pvt. Ltd.**

Plot No. 186, Jigani Link Road  
Jigani, Anekal Taluk, Bangalore - 562106, India  
Phone: 9180-27826725 • Fax: 9180-27826724  
Web Site: [www.kewaunee.in](http://www.kewaunee.in)

**Kewaunee Labway Asia Pte. Ltd.**

194 Pandan Loop, #06-22  
Pantech Business Hub, Singapore 128 383  
Phone: 65-6773-0288 • Fax: 65-6773-2322